

BenchTOF2[™] mass spectrometers

See what you've been missing





An introduction to SepSolve Analytical

Experts in analytical chemistry





Why choose BenchTOF2?

Enhanced sensitivity

Lower detection limits (IDL <20 fg) for a wider range of compounds in a single analysis

Superior spectral quality

Confident identification of analytes, through improved match factors and improved isotope abundances

Powerful mass accuracy

Mass accuracy of **<50ppm** and an easy-to-use mass-to-formula calculator for improved identification of unknowns



Extended dynamic range

Accurate quantification (across 5 orders of magnitude) of high-concentration compounds, while still reaching the required low detection limits. Removing the need for dilutions/repeat analyses.

Hydrogen carrier

Fast chromatographic separations with reduced running costs and faster ROI

Tandem Ionisation®

Patented technology for simultaneous hard and soft EI. Providing uniquely confident identification, through all new streamlined workflows



Time-of-flight (TOF) mass spectrometry

Proprietary BenchTOF2 technology





Improved stability and simplified maintenance

Redesigned ion source for 2025





- A redesigned ion source boosts robustness and prolongs cleanliness
- Enhanced performance for high boilers

Designed for minimal downtime

Tool-free ion source installation and removal





Get MORE from your mass spectrometer

BenchTOF2 provides unparalleled non-target screening...

- ✓ MORE sensitivity
- ✓ MORE confidence
- ✓ MORE productivity
- ✓ MORE efficiency





More confidence

Be sure you're identifying the correct compound





More confidence

Be sure you're identifying the correct compound







More confidence

Be sure you're identifying the correct compound



 User-friendly mass-to-formula calculator in ChromSpace[®] software simplifies identification

ldv	Form	ula		RDB		Delta PPM		_
1	CaH			3		9 4027		
2	C ₃ H	1004		-1		147.9911		
3	C _e H	14		2		-321,147		
4	C ₆ H	5 O 2		4		339.9525		
Element settin	igs:	Basic					~ E0	dit
otope patter	,							
Somula:		C ₂ H ₂₀ O	1					
Dattern fit		0.070						
Pattern fit:		0.978	3					
Pattern fit: Mass		0.978	Abundance		Unknown		f(x)	
Pattern fit: Mass 110.073		0.978	Abundance 100		Unknown 100		f(x)	
Pattern fit: Mass 110.073 111.077		0.978	Abundance 100 7.72412		Unknown 100 8.33857		f(x)	
Pattern fit: Mass 110.073 111.077 112.079		0.978	Abundance 100 7.72412 0.462852		Unknown 100 8.33857 0.263725		f(x)	
Pattern fit: Mass 110.073 111.077 112.079 113.081		0.978	Abundance 100 7.72412 0.462852 0.020607		Unknown 100 8.33857 0.263725 0		f(x)	
Pattern fit: Mass 110.073 111.077 112.079 113.081 114.084		0.978	Abundance 100 7.72412 0.462852 0.020607 0.000577663		Unknown 100 8.33857 0.263725 0 0		f(x)	
Pattern fit: Mass 110.073 111.077 112.079 113.081 114.084 Parametric Exp	pressio	0.978 0.978	Abundance 100 7.72412 0.462852 0.020607 0.000577663		Unknown 100 8.33857 0.263725 0 0		f(x)	
Pattern fit: Mass 110.073 111.077 112.079 113.081 114.084 Parametric Exp Tolerance %	pressi	on 20 🖨	Abundance 100 7.72412 0.462852 0.020607 0.000577663		Unknown 100 8.33857 0.263725 0 0		f(x) ♥ □ □ □	



More productivity

Extended Dynamic Range (EDR) across 5 orders of magnitude





- Confidently detect and quantify high- & lowconcentration compounds in a single calibration series
- Increased sample throughput by reducing the number of dilutions and repeat analyses required



More efficiency

BenchTOF2 is fully certified for use with H₂ carrier gas



- **Economic benefits**: lower running costs than He, for savings of > 50% on your GC gas supply & faster ROI
- Time-savings increase lab productivity: faster throughput and improved chromatographic resolution
- Environmentally-friendly: eliminate He cylinder delivery to reduce your lab's carbon footprint



Patented technology

Tandem Ionisation®



Acquire hard and soft EI simultaneously

 Complementary spectra to confirm compound identity

Confirm isomer identity



Patented technology

Tandem Ionisation®



Benefit from the excellent spectral quality of BenchTOF2 at 70 eV while also acquiring complementary soft EI data Reduced fragmentation and enhanced structurally-significant ions for confirmation of compound identity



Gaining greater insight into sample composition

BenchTOF2 is the ideal partner to GC×GC separations





Why are fast acquisition speeds important?



- 10-20 data points are needed to accurately define a chromatographic peak
- Fast acquisitions speeds are needed to tackle the narrow peak widths in GC×GC and hyper-fast GC





Enhanced separation by GC×GC–TOF MS





Data courtesy of Nadin Boegelsack, Mount Royal University – see webinar "Food, fire and fuel: Recent advances in the application of gas chromatography" for full details

Evaluating cross-contamination

 Glove box study can evaluate cross-contamination by passive transmission (e.g. a leaking or incorrectly sealed container)





Data courtesy of Nadin Boegelsack, Mount Royal University – see webinar "Food, fire and fuel: Recent advances in the application of gas chromatography" for full details

Summary

The BenchTOF2 provides:

- Excellent spectral quality and improved mass accuracy (<50 ppm) to eliminate guesswork in identifications.
- Trace-level sensitivity for targets and unknowns, revealing previously undetected components.
- **Tandem Ionisation** to enhance confidence in identifying challenging compounds, such as isomers.
- Fast acquisition speeds, making it an ideal detector for GC×GC and hyper-fast GC.
- Futureproof technology certified for use with hydrogen carrier gas.
- Simplified workflows and meaningful results through easy-to-use acquisition and processing software







Contact SepSolve

- Email: hello@sepsolve.com
- Tel.: +44 (0)1733 669222 (UK) +49 (0)69 668 108 920 (DE) +1 519 206 0055 (Canada/US)
- Web: www.sepsolve.com
- Twitter: @SepSolve
- LinkedIn: www.linkedin.com/company/sepsolve-analytical



