

# SIMPLIFY AND ACCELERATE YOUR COMPLEX SCREENING METHODS

The Measure of Confidence



## Agilent GC/MSD Forensic/Toxicology Analyzer

### Reduce the strain on your resources by implementing new technologies for analyzing difficult matrices

Based on Agilent's 5977A Series GC/MSD and 7890B GC System, our user-friendly GC/MSD Forensic/Toxicology Analyzer quickly screens and quantitates large numbers of target compounds in complex matrices – all within a single analysis. Its full-scan EI methods give you many advantages for broad-range screening, such as unlimited targets, full-spectrum identity confirmation, and library searching for non-target identification. Additionally, the Analyzer is pre-tested for forensic/toxicology analysis – with inlet, column, capillary flow technology, and software tools factory installed and configured – to save you method development time.

### Screen *more* target compounds... in *less* time

Agilent's GC/MSD Forensic/Toxicology Analyzer makes use of productivity-boosting GC/MS technologies that allow you to:

- Increase the number of targets screened
- Differentiate target compounds from matrix interference
- Reduce the analysis time required per sample
- Perform a complete data review including screening and quantitation in less than 10 minutes
- Produce consistent, high-quality results immediately after installation

### The following components are included – saving you time and money:

- Retention Time Locked application-specific column, ensuring reliable database matching
- Video training tutorials facilitate learning of more advanced Analyzer features
- Quick-start guide and Application Note that demonstrate how to run the screening method provided with the Analyzer
- CD-ROM with analysis methods, data files, and reports



**Agilent Technologies**


These built-in features make it *faster* and *easier* to screen large numbers of target compounds in complex matrices



**Nitrogen-phosphorus detector (NPD) with advanced Blos bead technology** provides sensitive nitrogen detection, ensures stable operation, and more than doubles bead lifetime.



**Retention Time Locking (RTL)** for consistent retention times after column maintenance and easy matching with the 725-compound forensic/toxicology library.




**Deconvolution Reporting Software (DRS)** for fast data review in 10 minutes per sample, with screening and quantitation in one run.



**Capillary Flow Technology (CFT) and backflush** shorten cycle times, reduce chemical background, extend column life and reduce frequency of source cleaning.

See how our GC/MSD Forensic/Toxicology Analyzer can get your lab on the *fast track* to better broad-range screening



**Fast acquisition of multiple signals:** In this toxicology screen of blood extract, Scan, SIM and NPD signals were collected simultaneously in one 9.75-minute run (injection to injection) instead of three 30-minute runs.

**MSD Deconvolution Report**  
 Sample Name: CA5995  
 Data File: C:\msdchem\1\Appnote\FT5\_4x10m\_SamplesSimScan\CA5995\_nss.D  
 Date/Time: 11:39 AM Wednesday, Apr 2 2008

The NIST library was not searched for the components that were found in the AMDIS target library.

R.T.	Cas #	Compound Name	Agilent	AMDIS	
			ChemStation Amount (~ng)	Match	R.T. Diff sec.
1.539	54115	Nicotine	0.03	59	-0.5
1.6446	98920	Nicotinamide	0.27	93	-0.9
2.1631	999401024	Carisoprodol artifact	64.87	93	-0.5
2.6367	486566	Cotinine	1	96	-0.4
2.928	57534	Meprobamate	4.11	99	0.0
3.033	59062	Caffeine	0.04	82	-0.5
3.1832	78444	Carisoprodol	127.4	96	1.0
3.8653	76993	Methadone	0.39	74	-0.1
4.2279	7199293	Cyheptamide	22.5	98	0.1
4.8014	76426	Oxycodone	2.37	82	0.0
5.850	57885	Cholesterol	922.73	97	3.4

**A DRS report (with locked retention times)** of the blood extract toxicology screen. DRS works with Agilent's forensic/toxicology library to quickly and accurately identify target compounds in high-matrix samples.

### Put your lab on the productivity fast track.

Contact your local Agilent Representative or Agilent Authorized Distributor

Call **800-227-9770** (in the U.S. or Canada) or visit [www.agilent.com/chem/appkits](http://www.agilent.com/chem/appkits)

This information is subject to change without notice.  
 © Agilent Technologies, Inc. 2013  
 Printed in U.S.A., April 10, 2013  
 5990-6300EN

### Ordering information:

Order an Agilent **5977A Series GC/MSD** along with an Agilent **7890B GC system** using the following Part Number:

- **G3445B#471: Forensic/Toxicology DRS Screening GC/MSD Analyzer**



**Agilent Technologies**