



Washing Adherent Cells in Agilent Seahorse XF96 Cell Culture Microplates

Basic Procedure

Before performing an Agilent Seahorse XF Assay, growth medium must be replaced with a suitable assay medium (generally this means medium without bicarbonate buffer or serum and with low/no phenol red content).

This procedure describes replacing the growth medium with assay medium for adherent cells grown in XF96 Cell Culture Microplates prior to being assayed using an Agilent Seahorse XF^e96/XF96 Analyzer.

For choosing and preparing the appropriate assay medium, please see <http://www.agilent.com/cs/library/selectionguide/public/5991-7878EN.pdf> and http://www.agilent.com/cs/library/usermanuals/public/XFe96_DAY_OF_MEDIA_PREP.pdf or http://www.agilent.com/cs/library/usermanuals/public/XF96_DAY_OF_MEDIA_PREP.pdf

1. Warm the assay medium to 37°C.
2. Retrieve the XF Cell Culture Microplate plate from the CO₂ incubator.
3. Look at the cells under the microscope to:
 - a. Confirm cell health, morphology, seeding uniformity and purity (no contamination).
 - b. Ensure cells are adhered, showing a consistent monolayer.
 - c. Make sure no cells were plated in the background correction wells.



4. Wash cells with the appropriate assay medium
 - a. Using an Agilent Seahorse XF Prep Station
 - i. Attach bottle of assay medium to XF Prep Station. Open the XF Prep Station software. On the "Media Change" tab, select "Do Prime", set final volume to 180 μ L of assay medium, and unselect "Do Rinse".
 - ii. Place the cell plate vertically onto the tray and remove the lid.
 - iii. Press "Start".
 - b. Without using an Agilent Seahorse XF Prep Station
 - i. Remove all but 20 μ L of the culture medium from each well.
 - ii. Rinse cells two times with 200 μ L of assay medium, leaving 20 μ L behind after each wash.
 - iii. Add 160 μ L of assay medium to each well for a final volume of 180 μ L/well.
5. Look at cells under the microscope to ensure that cells were not disturbed or washed away.
6. Place the plate in a 37°C incubator without CO₂ for 45-60 minutes prior to the assay.

NOTE: Incubating the cell plates without CO₂ allows outgassing from the plate and is required for accurate ECAR measurements.

Learn more
www.agilent.com/en-us/promotions/seahorse-xf-technology

For Research Use Only. Not for use in diagnostic procedures.

Product specifications and descriptions in this document
are subject to change without notice.

© Agilent Technologies, Inc., 2014, 2017
Published in the USA, June 30, 2017



Agilent Technologies