



CHEM-SERVICE^{INC}

Pesticide and Metabolite Standards Catalog



Table of Contents

GENERAL INFORMATION	2
Ordering Information	
Fax Order Form	
Making Dilutions	
Miscibility Table	
How To Read a Chem Service Label	
PESTICIDE AND METABOLITE STANDARDS - ALPHABETICAL LISTING	7
PESTICIDE STANDARDS LISTING BY CAS NUMBER WITH STRUCTURES	28
METABOLITE STANDARDS LISTING BY PARENT COMPOUND WITH STRUCTURES ..	161
GOVERNMENT/STATE/INTERNATIONAL AGENCY STANDARDS	178
EPA 500	180
EPA 600	190
EPA 1600	197
EPA 8000	199
EPA 8100	200
EPA 8200	203
EPA 8300	204
CLP	205
USP	207
ISO	208
Canadian Regulatory Standards	208
DIN (German Regulatory Standards)	209
European Regulatory Standards	209
Cannabis Standards	211
PESTICIDE AND METABOLITE STANDARDS LISTING BY PART NUMBER	213
PESTICIDE AND METABOLITE STANDARDS LISTING BY CAS NUMBER	227
OTHER INFORMATION	229
Product Documentation	
Certificate of Analysis	
Replacement/Return Policy	
Shelf Life Information	
Grade Information	
Terms and Conditions	

Chem Service is Committed to Quality and to Our Customers



QUALITY RECOGNITION

Chem Service is recognized worldwide for the quality of the materials we provide. We have maintained ISO Guide 34, ISO 17025 accreditation and ISO 9001 certification. This stringent level of compliance applies to all of our internal design, development, production, distribution, and servicing of organic neat and synthetic reference materials which, we supply to our customers.



QUALITY STANDARDS

Chem Service provides Certificates of Analysis for our Standards Grade pesticide and metabolite standards in accordance with our internal policy requirements, and ISO Guide 34, ISO 17025, and ISO 9001 Quality System Requirements. Our neat chemical certification policy states that purity and identification must be established independently by performing three or more (when feasible) separate analyses before we certify any of our neat products. More than 95% of our neat Standards Grade materials have a certified purity of 98.0% or greater. Our pesticide standards, whether sold as neat products or as solutions, are certified for purity, identity, stability and are traceable by lot.



QUALITY ASSURANCE

Our internal quality program utilizes our company's years of experience in providing laboratories with highly purified materials for reference standards. Our experience has enabled us to identify "problem" materials that need special attention. For example, any potential peroxide former is tested frequently for peroxide formation to assure both purity and safety. Moreover, with our expertise in analyte compatibility, solubility, solution stability and analysis, Chem Service is uniquely qualified to provide your customized solutions and mixture standards.



QUALITY TESTING

Our on-going stability testing program for Standards Grade materials has additional benefits for our customers. We will supply neat standards only when they have a minimum of the shelf life remaining on their lot certification. If our stability testing shows the material has a purity or stability problem which affects its expiration date, we notify our customers and offer a replacement. Chem Service is committed to supplying quality products and services to our customers.

How to Place Orders

You have several options on how to place an order with Chem Service.

STOCK AND CUSTOM ORDERS

- **Online:** <http://www.chemservice.com>
- **Phone Orders:** 1-800-452-9994 or 1-610-692-3026
- **Mail Orders:** P.O. Box 599, West Chester, PA 19381-0599
- **Fax Orders:** 1-610-692-8729
- **Email orders:** info@chemservice.com

We have included forms for general orders and for custom orders on the following pages. You can copy these forms and fill in the indicated information.

Ordering Information – Phone Orders

ORDERING INFORMATION

Chem Service's business hours are Monday - Thursday 8:00 AM to 6:00 PM ET and Friday 8:00AM - 5:00PM ET for ordering and support of all products.

To facilitate the rapid placement of an order, please have your Chem Service customer account number available. We recommend placing this number on the CAN# line at the top of this page.

To avoid errors, please state the part number and associated catalog name for each item being purchased. At the end of each order, your Chem Service representative will state the order reference number. Please document this reference number for future order status. Generally a mail or fax confirming order is not necessary for phone orders. If a confirming order is sent, please indicate that it is confirming and indicate our order reference number.

Prices are F.O.B. West Chester, PA, net 30 days. Sizes and Prices are subject to change without notice.

RUSH ORDERS (SAME DAY SHIPMENT)

Chem Service will accept rush orders for any in-stock Chem Service item until 4 PM ET, for next day delivery by Federal Express, UPS, or DHL. This service is provided free of charge to all Chem Service customers. Please inform your sales representative whether you are placing a rush order and your preferred carrier. Chem Service, Inc. assumes no liability except for due diligence in handling and shipping of rush orders.

NEW CUSTOMERS

Until credit has been established, we request that first orders over \$100 be shipped C.O.D., be held for check in advance, or be paid for by credit card. A credit application must be completed and received by Chem Service, Inc. prior to the shipment of all first time orders.

INTERNATIONAL CUSTOMERS

Company policy on payment is check in advance, wire transfer or credit card, covering all chemicals, shipping and documentation fees. Check to be drawn on a US Bank in US Funds. A documentation fee may be charged for any additional paperwork required for shipment. Please advise us of any special documentation or handling requirements in advance of order. All applicable duties and taxes are the responsibility of the recipient.

REPLACEMENT POLICY

Prior to the expiration date shown on the label and exclusive of any customer contamination, Chem Service will replace any product showing purity or concentration degradation. You must advise our Customer Service Department immediately to receive approval and instructions on either disposal or return of the product. Call 1-800-452-9994.

PRODUCT DOCUMENTATION

CERTIFICATE OF ANALYSIS (COA) - With the purchase of any Chem Service "Standards Grade" environmental, pesticide, organic chemicals, solutions or mixtures, a COA will be provided free of charge, where available. Please note that "Standards Grade" products assigned a purity of tech, will not receive a COA. COAs are not available for "Tech Grade" and "No Grade" products, but, where available, these products will come with an IR.

SAFETY DATA SHEETS (SDS) - English language SDS's are shipped in the box with the chemicals, exceeding OSHA's requirements. Additional copies, replacement copies, updated versions, SDS in advance of an order, or additional languages, can be downloaded from our website www.chemservice.com.

WARNING - READ YOUR SAFETY DATA SHEET BEFORE HANDLING OR USING ANY CHEMICALS.

WE ACCEPT...



**P.O. BOX 599
WEST CHESTER, PA 19381, USA
PHONE: 610-692-3026
FAX: 610-692-8729**

DATE _____ PURCHASE ORDER NO.# _____

CONTACT NAME _____ COMPANY NAME _____

CHEM SERVICE ACCOUNT NO.# _____ DEPT. _____

PHONE

ATTENTION: _____

ADDRESS

CITY/STATE/ZIP _____

PHONE _____ FAX _____ EMAIL _____

ADDRESS _____

CITY/STATE/ZIP _____

PHONE _____ FAX _____ EMAIL _____

DATE NEEDED _____ SHIP VIA: _____

QUOTE NO.# UPS or FEDEX ACCT. NO.#:

[illegible]

CREDIT CARD: ACCT.#: EXP. DATE: / /

CARD HOLDER BILLING ZIP CODE#:

NOTICE: All orders are provided subject to the Terms and Conditions of Sale which appear in the Chem Service Catalog, on the back of each Packaging List, and on our web site, www.chemservice.com. THERE ARE LIMITATIONS OF WARRANTIES IN THESE TERMS INCLUDING BUT NOT LIMITED TO THE CONDITION THAT THERE IS NO EXPRESS NOR IMPLIED WARRANTY OF MERCHANTABILITY NOR FITNESS FOR A PRACTICAL PURPOSE. There is also no warranty as to the analysis of the product nor the result that will be obtained from its use except as may be set forth in the Terms and Conditions of Sale.

It is often necessary to make a dilution of your chemicals in order to analyze them. Dilutions can be made in one of three manners: weight/weight, weight/volume and volume/volume. The most common of these is weight/volume. The units of measurement used are commonly expressed in parts per million (ppm), or parts per billion (ppb). For weight/weight, it is expressed as 1 ppm = 1 ug/g, for weight/volume, 1 ppm = 1 ug/mL, and for volume/volume, 1 ppm = 1 ul/L. Chem Service supplies most of their solutions at a concentration of 0.1 mg/mL (100 ug/mL or 100 ppm). In order to make a solution at a concentration of 100 ppm, starting with a pure chemical, you need to know your final volume. Then if the material is a solid, the calculations are:

$$0.1 \text{ mg/mL (final concentration)} \times 100 \text{ mL (final volume)} = \text{weight of chemical needed (10 mg)}$$

$$\text{FC} \times \text{FV} = \text{W}$$

If the material is a liquid, it is usually easier to convert to a liquid measurement (i.e. ul), for handling purposes. Therefore, the density is needed and the calculations are:

$$\frac{0.1 \text{ mg/mL (final conc.)} \times 100 \text{ mL (final volume)}}{\text{density (mg/ul)}} = \text{volume (number of ul) of chemical needed}$$

$$\frac{\text{FC} \times \text{FV}}{d} = V$$

If you need to make a further dilution from a prepared solution, the calculations would be basically the same format. For example, to dilute 100 ppm to 1 ppm, take a 1 mL aliquot of 0.1 mg/mL (100 ppm) solution, and dilute to your final volume (100 mL).

$$\frac{0.1 \text{ mg/mL (initial concentration)} \times 1 \text{ mL (volume injected)} \times 1000 \text{ ug (conversion factor)}}{100 \text{ mL (final volume)} \times 1 \text{ mg}} = 1 \text{ ug/mL (final concentration)}$$

$$\frac{\text{IC} \times \text{VI} \times \text{CF}}{\text{FV}} = \text{FC}$$

The above formula gives the final concentration for the solution. The formula to obtain the volume to inject in order to make the dilution for this final volume is:

$$\frac{\text{FC} \times \text{FV} \times \text{CF}}{\text{IC}} = \text{VI}$$

As one can see, these calculations can be manipulated to allow you to achieve any final volume or concentration you need. They can also be used to convert to other types of concentrations (weight/weight or volume/volume).

LEGEND

FC = final concentration

FV = final volume

W = weight needed

d = density

IC = initial concentration

VI = volume injected

CF = conversion factor

V = volume needed

Acetic acid

Acetone

Acetonitrile

Benzene

2-Butanol

Butyl alcohol

tert-Butylmethyl ether (MTBE)

Carbon tetrachloride

Chloroform

Cyclohexane

Cyclopentane

Dichloroethane

N,N-Dimethylformamide (DMF)

1,4-Dioxane

Dipropyl ether

Ethyl acetate

Ethyl alcohol

Ethyl ether

n-Heptane

n-Hexane

Isooctane

Isopropyl alcohol

Methanol

Methylene chloride

Methyl sulfoxide (DMSO)

n-Pentane

1,1,2,2-Tetrachloroethane

Tetrahydrofuran (THF)

Toluene

Trichloroethane

Water

Xylene

Legend



Miscible



Immiscible

Read down column and over
for solvent miscibility

Acetonitrile



Cyclohexane

are immiscible

Neat Standard Label

a
↓
N-11971-100MG
b
↓
Exp 11/30/2018

CHEM SERVICE
 PO BOX 599, WEST CHESTER, PA 19381 800-452-9994

c → Pipronil
 d → Lot: 5848400 Purity: 99.5% ← e
 f → TOXIC-HAZARDOUS TO THE ENVIRONMENT

WARNING! READ your SDS before HANDLING chemical. Full Label information is provided on the outside package.

a	Part number	d	Lot Number	Chem Service Part Number	
b	Expiration Date	e	Purity	N-11971-100MG	
c	Product Description	f	Hazard Warning	Product Number	Size of item.
				N = Neat NG = No Grade MET = Metabolite BZ = PCB	MG = Milligrams G = Grams

Solution and Mixture Standard Label

a
↓
S-13189A1-1ML
b
↓
Exp 1/31/2019

CHEM SERVICE
 PO BOX 599, WEST CHESTER, PA 19381 800-452-9994

c → S-Hydroprene Solution
 100ug/ml in Acetonitrile ← e
 d → Lot: 5827200
 f → FLAMMABLE-POISON-LACHRYMATOR

WARNING! READ your SDS before HANDLING chemical. Full Label information is provided on the outside package.

a	Part number	d	Lot Number	Chem Service Part Number	
b	Expiration Date	e	Solvent and Concentration	S-13189A1-1ML	
c	Product Description	f	Hazard Warning	Product Number	Size of item.
				S = Solution M = Mixture BZ = PCB	ML = Milliliters

Solvent and Concentration Key

Code	Solvent	Code	Solvent	Code	Concentration
A	Acetonitrile	V	Transformer Oil	0	50 ug/ml
B	Acetone	W	Methanol:Acetonitrile (80:20)	1	100 ug/ml
C	Benzene	X	Methylene Chloride	2	200 ug/ml
D	Carbon Disulfide	Y	Isooctane:Acetone (90:10)	3	500 ug/ml
E	Cyclohexane	Z	Acetonitrile:Acetone (90:10)	4	1000 ug/ml
F	Water	AA	Acetonitrile:Water (75:25)	5	2000 ug/ml
G	Ethanol	AB	Benzene:MeCl2:Acetonitrile (40:40:20)	6	2500 ug/ml
H	Ethyl Acetate	AC	Toluene:Hexane (50:50)	7	5000 ug/ml
J	Hexane	AD	CH2Cl2:Benzene (50:50)	8	10000 ug/ml
K	Isooctane	AE	Hexane:Methylene chloride (90:10)	9	20000 ug/ml
L	Isopropanol	AF	Hexane:Acetone (90:10)	10	10 ug/ml
M	Methanol	AG	Hexane:Acetone (80:20)	11	250 ug/ml
N	Methanol:Water	AH	Methanol:Acetonitrile (50:50)	12	4000 ug/ml
P	n-Nonane	AJ	Methylene chloride:Benzene (50:50)	13	1970 ug/ml
Q	n-Pentane	AK	Methanol:Acetone (50:50)	14	3000 ug/ml
R	P&T Methanol	AL	Acetone:Toluene (50:50)	15	400 ug/ml
S	Propanol	AM	2-Methoxyethanol	16	1300 ug/ml
T	t-Butylmethyl ether	AN	Dimethyl sufoxide	17	1770 ug/ml
U	Toluene	AP	Methanol:Water (90:10)	18	25 ug/ml
				19	20 ug/ml
				20	30 ug/ml



Pesticide and Metabolite Standards Alphabetical Listing

a	b	c	d
			
Aastar(TM)	70124-77-5	N-11984-100MG	100MG

a Chem Service Product Description

c Part Number

b Chemical Abstracts Registry Number

d Size

Description	CAS#	Part #	Size	Description	CAS#	Part #	Size
Aastar(TM)	70124-77-5	N-11984-100MG	100MG	Allyl trimethylammonium bromide	3004-51-1	N-11059-250MG	250MG
Aastar(TM) Solution	70124-77-5	S-11984A1-1ML	1ML	Allyl trimethylammonium bromide solution	3004-51-1	S-11059F1-1ML	1ML
Aastar(TM) Solution	70124-77-5	S-11984U1-1ML	1ML	Allyl trimethylammonium bromide solution	3004-51-1	S-11059U1-1ML	1ML
Abamectin	71751-41-2	N-10995-100MG	100MG	Alon	34123-59-6	N-12279-100MG	100MG
Abamectin solution	71751-41-2	S-10995A1-1ML	1ML	Alon Solution	34123-59-6	S-12279A1-1ML	1ML
Abate	3383-96-8	N-10996-100MG	100MG	Alsyntine(TM)	64628-44-0	N-13687-100MG	100MG
Abate Solution	3383-96-8	S-10996A1-1ML	1ML	Alsyntine(TM) Solution	64628-44-0	S-13687A1-1ML	1ML
Abate Solution	3383-96-8	S-10996U1-1ML	1ML	Alsyntine(TM) Solution	64628-44-0	S-13687U1-1ML	1ML
AC 263, 222(TM)	104098-48-8	N-12202-250MG	250MG	Altima(TM)	79622-59-6	N-11980-100MG	100MG
AC 263, 222(TM) Solution	104098-48-8	S-12202A1-1ML	1ML	Altima(TM) Solution	79622-59-6	S-11980A1-1ML	1ML
Acaban(TM)	134098-61-6	N-11962-100MG	100MG	Altima(TM) Solution	79622-59-6	S-11980U1-1ML	1ML
Acaban(TM) Solution	134098-61-6	S-11962A1-1ML	1ML	Amdro	67485-29-4	N-11063-100MG	100MG
Acaban(TM) Solution	134098-61-6	S-11962U1-1ML	1ML	Amdro solution	67485-29-4	S-11063A1-1ML	1ML
Acarol (TM)	18181-80-1	N-10997-250MG	250MG	Amdro solution	67485-29-4	S-11063J1-1ML	1ML
Acarol (TM) Solution	18181-80-1	S-10997A1-1ML	1ML	Amex(TM)	33629-47-9	N-11355-500MG	500MG
Accent (TM)	111991-09-4	N-10998-100MG	100MG	Amex(TM) Solution	33629-47-9	S-11355A1-1ML	1ML
Acclaim(TM)	71283-80-2	N-11958-100MG	100MG	Amex(TM) Solution	33629-47-9	S-11355K1-1ML	1ML
Acclaim(TM) Solution	71283-80-2	S-11958A1-1ML	1ML	Amicarbazone	129909-90-6	N-13075-100MG	100MG
Acclaim(TM) Solution	71283-80-2	S-11958U1-1ML	1ML	Amicarbazone solution	129909-90-6	S-13075A1-1ML	1ML
Acephate	30560-19-1	N-11002-250MG	250MG	Amidosulfuron	120923-37-7	N-13819-100MG	100MG
Acephate solution	30560-19-1	S-11002A1-1ML	1ML	2-Aminobenzimidazole	934-32-7	MET-11138A-250MG	250MG
Acephate solution	30560-19-1	S-11002B4-1ML	1ML	Aminocarb	2032-59-9	N-11065-250MG	250MG
Acephate solution	30560-19-1	S-11002B4-5ML	5ML	Aminocarb solution	2032-59-9	S-11065A4-1ML	1ML
Acequinocyl	57960-19-7	N-11003-250MG	250MG	Aminocarb solution	2032-59-9	S-11065A4-5ML	5ML
Acequinocyl solution	57960-19-7	S-11003A1-1ML	1ML	Aminocarb solution	2032-59-9	S-11065U1-1ML	1ML
Acequinocyl solution	57960-19-7	S-11003U1-1ML	1ML	Aminocarb-phenol	14143-25-0	MET-11065-100MG	100MG
Acetaldehyde polymerized	9002-91-9	N-12381-1G	1G	2-Amino-4-chloro-6-ethylamino-s-triazine	1007-28-9	MET-13213A-50MG	50MG
Acetamidrid	135410-20-7	N-11009-100MG	100MG	4-Amino-2-chloro-6-isopropylamino-s-triazine	6190-65-4	MET-11106B-250MG	250MG
Acetamidrid solution	135410-20-7	S-11009A1-1ML	1ML	2-Amino-3-chloro-1,4-naphthoquinone	2797-51-5	N-13825-100MG	100MG
Acetamidrid solution	135410-20-7	S-11009U1-1ML	1ML	2-Amino-5-chlorotoluene	95-69-2	MET-11426A-250MG	250MG
Acetamidrid-N-desmethyl	190604-92-3	MET-111009A-10MG	10MG	2-Amino-5-chlorotoluene	95-69-2	MET-11426A1-1ML	1ML
Acetochlor	34256-82-1	N-11013-100MG	100MG	Aminocyclopyrachlor	858956-08-8	N-13826-250MG	250MG
Acetochlor solution	34256-82-1	S-11013M1-1ML	1ML	2-Amino-1,3-diethylbenzene	579-66-8	MET-11043C-500MG	500MG
Acetochlor ESA sodium salt	187022-11-3	MET-11013A-25MG	25MG	4-Amino-2-hydroxy-6-isopropylamino-s-triazine	19988-24-0	MET-11106D-50MG	50MG
Acetochlor ESA sodium salt Solution	187022-11-3	MET-11013AM1-1ML	1ML	4-Amino-2-hydroxy-6-isopropylamino-s-triazine	19988-24-0	MET-11106DAM1-1ML	1ML
Acetochlor OA	194992-44-4	MET-11013B-10MG	10MG	2-Amino-N-isopropyl benzamide	30391-89-0	MET-11142C-50MG	50MG
Acetochlor OA Solution	194992-44-4	MET-11013BM1-1ML	1ML	Aminomethyl phosphonic acid	1066-51-9	MET-12133A-100MG	100MG
sec-Acetochlor OA solution	152019-74-4	S-13946M0-1ML	1ML	Aminopyralid	150114-71-9	N-11066-100MG	100MG
3-(a-Acetyl benzyl)-4-hydroxy-coumarin	7779-90-0	N-13758-250MG	250MG	4-Aminopyridine	504-24-5	N-10806-1G	1G
Acibenzolar-S-methyl	135158-54-2	N-11026-250MG	250MG	4-Aminopyridine solution	504-24-5	S-10806M1-1ML	1ML
Acibenzolar-S-methyl solution	135158-54-2	S-11026A1-1ML	1ML	Aminotriazole	61-82-5	N-11067-250MG	250MG
Acibenzolar-S-methyl solution	135158-54-2	S-11026U1-1ML	1ML	Aminotriazole solution	61-82-5	S-11067A1-1ML	1ML
Acifluorfen	50594-66-6	N-11027-250MG	250MG	Aminotriazole solution	61-82-5	S-11067B1-1ML	1ML
Acifluorfen solution	50594-66-6	S-11027A1-1ML	1ML	Amisulbrom solution	348635-87-0	S-13960M1-1ML	1ML
Acifluorfen solution	50594-66-6	S-11027B1-1ML	1ML	Amitraz	33089-61-1	N-11068-250MG	250MG
Acifluorfen solution	50594-66-6	S-11027B1-5ML	5ML	Amitraz solution	33089-61-1	S-11068A1-1ML	1ML
Acifluorfen,amino	50594-66-6	MET-11027A-10MG	10MG	Amitraz Metabolite BTS 27271	33089-74-6	MET-11068B-10MG	10MG
Acifluorfen methyl ester	50594-67-7	N-11028-100MG	100MG	Amitrol	504-24-5	N-10806-1G	1G
Acifluorfen methyl ester solution	50594-67-7	S-11028Y1-1ML	1ML	Amitrol Solution	504-24-5	S-10806M1-1ML	1ML
Acifluorfen methyl ester solution	50594-67-7	S-11028Y1-5ML	5ML	Amitrole	61-82-5	N-11067-250MG	250MG
Acionifen	74070-46-5	N-13991-250MG	250MG	Amitrole Solution	61-82-5	S-11067A1-1ML	1ML
Acionifen solution	74070-46-5	S-13991A1-1ML	1ML	Amitrole Solution	61-82-5	S-11067B1-1ML	1ML
Acrinathrin	101007-06-1	N-13928-100MG	100MG	Ammate (TM)	7773-06-0	N-11070-500MG	500MG
Acrinathrin solution	101007-06-1	S-13928A1-1ML	1ML	Ammonium ethylcarbamoylphosphonate Solution	25954-13-6	S-12018F1-1ML	1ML
Acrobat(TM)	110488-70-5	N-11759-100MG	100MG	Ammonium fluosilicate(Technical)	16919-19-0	N-11069-500MG	500MG
Acrobat(TM) Solution	110488-70-5	S-11759A1-1ML	1ML	Ammonium sulfamate(Technical)	7773-06-0	N-11070-500MG	500MG
Actara(TM)	153719-23-4	N-13562-100MG	100MG	Amobam	3566-10-700	N-11071-100MG	100MG
Actara(TM) Solution	153719-23-4	S-13562M1-1ML	1ML	Amobam solution	3566-10-700	S-11071M1-1ML	1ML
Admiral(TM)	95737-68-1	N-13160-100MG	100MG	Anchor (TM)	77732-09-3	N-12732-100MG	100MG
Admiral(TM) Solution	95737-68-1	S-13160M1-1ML	1ML	Anchor (TM) Solution	77732-09-3	S-12732A1-1ML	1ML
Advance(TM)	69377-81-7	N-12000-100MG	100MG	Anchor (TM) Solution	77732-09-3	S-12732U1-1ML	1ML
Advance(TM) Solution	69377-81-7	S-12000M1-1ML	1ML	Ancymidol	12771-68-5	N-11073-100MG	100MG
Afos(TM)	2595-54-2	N-12359-50MG	50MG	Ancymidol solution	12771-68-5	S-11073M1-1ML	1ML
Afos(TM) Solution	2595-54-2	S-12359A1-1ML	1ML	Anilazine	101-05-3	N-11075-250MG	250MG
Afos(TM) Solution	2595-54-2	S-12359U1-1ML	1ML	Anilazine solution	101-05-3	S-11075A1-1ML	1ML
Afugan(TM)	13457-18-6	N-13146-100MG	100MG	Anilazine solution	101-05-3	S-11075U1-1ML	1ML
Afugan(TM) Solution	13457-18-6	S-13146A1-1ML	1ML	Anilazine solution	101-05-3	S-11075U1-5ML	5ML
Afugan(TM) Solution	13457-18-6	S-13146U1-1ML	1ML	Anilofos	64249-01-0	N-13131-100MG	100MG
Agribon(TM)	122-11-2	N-13242-250MG	250MG	Anilofos solution	64249-01-0	S-13131K1-1ML	1ML
Agribon(TM) Solution	122-11-2	S-13242A1-1ML	1ML	Aniloguard	64249-01-0	N-13131-100MG	100MG
Agribrom	16079-88-2	N-11041-500MG	500MG	Aniloguard Solution	64249-01-0	S-13131K1-1ML	1ML
Agrothrin(TM)	65731-84-2	N-11192-250MG	250MG	Anthio™	2540-82-1	N-12017-100MG	100MG
Agrothrin(TM) Solution	65731-84-2	S-11192A1-1ML	1ML	Anthio™ Solution	2540-82-1	S-12017A1-1ML	1ML
Agrothrin(TM) Solution	65731-84-2	S-11192T1-1ML	1ML	Anthio™ Solution	2540-82-1	S-12017J1-1ML	1ML
Aim(TM)	128639-02-1	N-11411-10MG	10MG	9,10-Anthraquinone	84-65-1	N-10970-250MG	250MG
Aim(TM) Solution	128639-02-1	S-11411A1-1ML	1ML	9,10-Anthraquinone solution	84-65-1	S-10970M1-1ML	1ML
Aim(TM) Solution	128639-02-1	S-11411U1-1ML	1ML	ANTU	86-88-4	N-10983-250MG	250MG
Akton (TM)	1757-18-2	N-11042-100MG	100MG	ANTU Solution	86-88-4	S-10983A1-1ML	1ML
Akton (TM) solution	1757-18-2	S-11042A1-1ML	1ML	ANTU Solution	86-88-4	S-10983U1-1ML	1ML
Akton (TM) solution	1757-18-2	S-11042U1-1ML	1ML	Anvil(TM)	79983-71-4	N-12165-100MG	100MG
Alachlor	15972-60-8	N-11043-250MG	250MG	Anvil(TM) Solution	79983-71-4	S-12165A1-1ML	1ML
Alachlor (ring 13C6) Solution	15972-60-8	S-FC2072S-1.2ML	1.2ML	Anvil(TM) Solution	79983-71-4	S-12165U1-1ML	1ML
Alachlor solution	15972-60-8	S-11043M1-1ML	1ML	Apache(TM)	95465-99-9	N-11392-100MG	100MG
Alachlor solution	15972-60-8	S-11043M1-5ML	5ML	Apache(TM) Solution	95465-99-9	S-11392A1-1ML	1ML
Alachlor acetate	166407-15-4	MET-11043B-100MG	100MG	Apache(TM) Solution	95465-99-9	S-11392U1-1ML	1ML
Alachlor ESA sodium salt	140939-15-7	MET-11043D-25MG	25MG	Apollo	74115-24-5	N-11490-100MG	100MG
Alachlor ESA sodium salt Solution	140939-15-7	MET-11043DM1-1ML	1ML	Apollo Solution	74115-24-5	S-11490M1-1ML	1ML
Alachlor OA	171262-17-2	MET-11043E-10MG	10MG	Apron	57837-19-1	N-12380-100MG	100MG
Alachlor OA Solution	171262-17-2	MET-11043EM1-1ML	1ML	Apron Solution	57837-19-1	S-12380A1-1ML	1ML
Albendazole	54965-21-8	N-12959-10MG	10MG	Apron Solution	57837-19-1	S-12380U1-1ML	1ML
Aldicarb	116-06-3	N-11044-100MG	100MG	Aquaclyne(TM)	2058-46-0	N-12743-250MG	250MG
Aldicarb solution	116-06-3	S-11044M1-1ML	1ML	Aqua-Kleen (TM)	1929-73-3	N-10610-250MG	250MG
Aldicarb oxime	1646-75-9	MET-11044A-10MG	10MG	Aqua-Kleen (TM) Solution	1929-73-3	S-10610A1-1ML	1ML
Aldicarb sulfone	1646-88-4	N-11045-10MG	10MG	Aqua-Kleen (TM) Solution	1929-73-3	S-10610A4-1ML	1ML
Aldicarb sulfone solution	1646-88-4	S-11045M1-1ML	1ML	Aqua-Kleen (TM) Solution	1929-73-3	S-10610A4-5ML	5ML
Aldicarb sulfoxide	1646-87-3	N-11046-10MG	10MG	Aracid	80-38-6	N-13812-250MG	250MG
Aldicarb sulfoxide solution	1646-87-3	S-11046A1-1ML	1ML	Aramite	140-57-8	N-11085-10MG	10MG
Aldrin	309-00-2	N-11049-250MG	250MG	Aramite solution	140-57-8	S-11085J1-1ML	1ML
Aldrin (TM) (13C12, 99%) Solution	309-00-2	S-FC89S-1.2ML	1.2ML	Aresin(TM)	1746-81-2	N-12494-250MG	250MG
Aldrin solution	309-00-2	S-11049M1-1ML	1ML	Aresin(TM) Solution	1746-81-2	S-12494A1-1ML	1ML
Aldrin solution	309-00-2	S-11049M1-5ML	5ML	Aresin(TM) Solution	1746-81-2	S-12494U1-1ML	1ML
Alfa-4(TM)	15879-93-3	N-10978-250MG	250MG	Armour(TM)	76674-21-0	N-12005-100MG	100MG
Allethrin	584-79-2	N-11051-100MG	100MG	Armour(TM) Solution	76674-21-0	S-12005A1-1ML	1ML
Allethrin solution	584-79-2	S-11051X1-1ML	1ML	Armour(TM) Solution	76674-21-0	S-12005U1-1ML	1ML
Allethrin solution	584-79-2	S-11051X1-1ML	1ML	Arpege(TM)	112281-77-3	N-13538-100MG	100MG
Allethrin solution	584-79-2	S-11051X1-5ML	5ML	Arpege(TM) Solution	112281-77-3	S-13538A1-1ML	1ML
d-trans-Allethrin	28057-48-9	N-11205-250MG	250MG	Arpege(TM) Solution	112281-77-3	S-13538U1-1ML	1ML
d-trans-Allethrin Solution	28057-48-9	S-11205A1-1ML	1ML	Arsenal(TM)	81334-34-1	N-12203-100MG	100MG
d-trans-Allethrin Solution	28057-48-9	S-11205T1-1ML	1ML	Arsenal(TM) Solution	81334-34-1	S-12203A1-1ML	1ML
Allidochlor	93-71-0	N-12617-250MG	250MG	Asana(TM)	66230-04-4	N-11102-100MG	100MG
Allidochlor Solution	93-71-0	S-12617M1-1ML	1ML	Asana(TM) Solution	66230-04-4	S-11102A1-1ML	1ML
Alloxydim-Sodium	66003-55-2	N-11052-50MG	50MG	Asana(TM) Solution	66230-04-4	S-11102T1-1ML	1ML
4-Allyl-2-methoxyphenol	97-53-0	N-10803-250MG	250MG	Ascurit(TM)	67747-09-5	N-13093-250MG	250MG
4-Allyl-2-methoxyphenol solution	97-53-0	S-10803M1-1ML	1ML	Ascurit(TM) Solution	67747-09-5	S-13093A1-1ML	1ML

Description	CAS#	Part #	Size	Description	CAS#	Part #	Size
Ascurit(TM) Solution	67747-09-5	S-13093U1-1ML	1ML	Bendiocarb solution	22781-23-3	S-11135U1-1ML	1ML
Asomate	3586-60-5	N-13904-50MG	50MG	Benfluralin	1861-40-1	N-11136-250MG	250MG
Aspon (TM)	3244-90-4	N-11103-10MG	10MG	Benfluralin solution	1861-40-1	S-11136A1-1ML	1ML
Aspon (TM) solution	3244-90-4	S-11103A1-1ML	1ML	Benfluralin solution	1861-40-1	S-11136J4-1ML	1ML
Aspon (TM) solution	3244-90-4	S-11103J4-1ML	1ML	Benfluralin solution	1861-40-1	S-11136J4-5ML	5ML
Aspon (TM) solution	3244-90-4	S-11103J4-5ML	5ML	Benfuracarb	82560-54-1	N-11137-100MG	100MG
Assert(TM)	81405-85-8	N-12200-100MG	100MG	Benodanil	15310-01-7	N-12991-10MG	10MG
Assert(TM) Solution	81405-85-8	S-12200A1-1ML	1ML	Benomyl	17804-35-2	N-11138-100MG	100MG
Assert(TM) Solution	81405-85-8	S-12200U1-1ML	1ML	Benoxacor	98730-04-2	N-12957-10MG	10MG
Assure	76578-14-8	N-13174-100MG	100MG	Bensulfuron-methyl	83055-99-6	N-11139-100MG	100MG
Assure Solution	76578-14-8	S-13174A1-1ML	1ML	Bensulfuron-methyl solution	83055-99-6	S-11139A1-1ML	1ML
Assure Solution	76578-14-8	S-13174U1-1ML	1ML	Bensulide	741-58-2	N-11140-250MG	250MG
Assure II(TM)	100646-51-3	N-13175-100MG	100MG	Bensulide solution	741-58-2	S-11140M4-1ML	1ML
Assure II(TM) Solution	100646-51-3	S-13175A1-1ML	1ML	Bensulide solution	741-58-2	S-11140M4-5ML	5ML
Assure II(TM) Solution	100646-51-3	S-13175U1-1ML	1ML	Bensultap	17606-31-4	N-11141-250MG	250MG
Astix(TM)	16484-77-8	N-12363-250MG	250MG	Bensultap solution	17606-31-4	S-11141M1-1ML	1ML
Astix(TM) Solution	16484-77-8	S-12363A1-1ML	1ML	Bentazon	25057-89-0	N-11142-250MG	250MG
Astix(TM) Solution	16484-77-8	S-12363T1-1ML	1ML	Bentazon solution	25057-89-0	S-11142A4-1ML	1ML
Asulam	3337-71-1	N-11104-250MG	250MG	Bentazon solution	25057-89-0	S-11142A4-5ML	5ML
Asulam solution	3337-71-1	S-11104M1-1ML	1ML	Bentazon solution	25057-89-0	S-11142B1-1ML	1ML
Atabron(TM)	71422-67-8	N-11430-50MG	50MG	Bentazon methyl derivative	61592-45-8	N-11143-50MG	50MG
Atabron(TM) Solution	71422-67-8	S-11430A1-1ML	1ML	Bentazon methyl derivative solution	61592-45-8	S-11143Y1-1ML	1ML
Atraton	1610-17-9	N-11105-250MG	250MG	Bentazon methyl derivative solution	61592-45-8	S-11143Y1-5ML	5ML
Atraton solution	1610-17-9	S-11105A1-1ML	1ML	Bentazone metabolite	30391-89-0	MET-11142C-50MG	50MG
Atraton solution	1610-17-9	S-11105T1-1ML	1ML	Benthiavalicarb-isopropyl	177406-68-7	N-14076-50MG	50MG
Atraton solution	1610-17-9	S-11105T1-5ML	5ML	Benthiocarb	28249-77-6	N-11144-250MG	250MG
Atrazine	1912-24-9	N-11106-250MG	250MG	Benthiocarb solution	28249-77-6	S-11144M1-1ML	1ML
Atrazine (13C3) Solution		S-FC2208S-1.2ML	1.2ML	2-(1,3-dimethylbutyl)-Benzenamine	203448-76-4	MET-14035A-10MG	10MG
Atrazine (ethylamine-d5) Solution		S-FD2208S-1.2ML	1.2ML	Benzoic acid-benzyl ester	120-51-4	N-11182-100MG	100MG
Atrazine desethyl	6190-65-4	MET-11106B-250MG	250MG	Benzoic acid-benzyl ester Solution	120-51-4	S-11182A1-1ML	1ML
Atrazine desethyl desisopropyl	3397-62-4	MET-13213C-100MG	100MG	Benzoic acid-benzyl ester Solution	120-51-4	S-11182J5-1ML	1ML
Atrazine desethyl desisopropyl-2-hydroxy	645-92-1	MET-11106E-50MG	50MG	Benzoic acid-benzyl ester Solution	120-51-4	S-11182J7-1ML	1ML
Atrazine-desethyl-2-hydroxy Solution	19988-24-0	MET-11106DAM1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	S-11182J7-5ML	5ML
Atrazine desethyl-2-hydroxy	19988-24-0	MET-11106D-50MG	50MG	Benzoic acid-benzyl ester Solution	120-51-4	MET-13750B-250MG	250MG
Atrazine desisopropyl	1007-28-9	MET-13213A-50MG	50MG	Benzoic acid-benzyl ester Solution	120-51-4	S-14010X1-1ML	1ML
Atrazine (ethylamine-d5)	163165-75-1	N-FD2208-5-5MG	5MG	Benzoic acid-benzyl ester Solution	120-51-4	S-11177A4-1ML	1ML
Atrazine-2-hydroxy	2163-68-0	MET-11106F-100MG	100MG	Benzoic acid-benzyl ester Solution	120-51-4	S-11177A4-5ML	5ML
Authority*	122836-35-5	N-13253-250MG	250MG	Benzoic acid-benzyl ester Solution	120-51-4	S-11177U1-1ML	1ML
Authority* Solution	122836-35-5	S-13253M1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	N-10956-250MG	250MG
Avermectin B1	71751-41-2	N-10995-100MG	100MG	Benzoic acid-benzyl ester Solution	120-51-4	S-10956A1-1ML	1ML
Avermectin B1 Solution	71751-41-2	S-10995A1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	N-11182-100MG	100MG
Axial(TM)	243973-20-8	N-13057-100MG	100MG	Benzoic acid-benzyl ester Solution	120-51-4	S-11182A1-1ML	1ML
Axiom	1757-18-2	N-11042-100MG	100MG	Benzoic acid-benzyl ester Solution	120-51-4	S-11182J5-1ML	1ML
Axiom Solution	1757-18-2	S-11042A1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	S-11182J7-1ML	1ML
Axiom Solution	1757-18-2	S-11042U1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	S-11182J7-5ML	5ML
Azaconazole solution	60207-31-0	S-13976A1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	N-12676-100MG	100MG
Azadirachtin(Technical)	11141-17-6	N-11107-10MG	10MG	Benzoic acid-benzyl ester Solution	120-51-4	S-12676M1-1ML	1ML
Azadirachtin	11141-17-6	N-11107-10MG	10MG	Benzoic acid-benzyl ester Solution	120-51-4	S-11182J7-1ML	1ML
Azak	1918-11-200	N-13513-100MG	100MG	Benzoic acid-benzyl ester Solution	120-51-4	S-11182J7-5ML	5ML
Azimsulfuron solution	120162-55-2	S-13975A1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	N-12676-100MG	100MG
Azinphos-ethyl	2642-71-9	N-12140-100MG	100MG	Benzoic acid-benzyl ester Solution	120-51-4	S-12676M1-1ML	1ML
Azinphos-methyl	86-50-0	N-12139-250MG	250MG	Benzoic acid-benzyl ester Solution	120-51-4	MET-13176A-100MG	100MG
Azinphos-methyl solution	86-50-0	S-12139A1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	N-10001-1G	1G
Azinphos-methyl solution	86-50-0	S-12139U1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	N-10437-250MG	250MG
Azinphos-methyl solution	86-50-0	S-12139U1-5ML	5ML	Benzoic acid-benzyl ester Solution	120-51-4	S-10437M1-1ML	1ML
Azinphos-methyl oxon Solution	961-22-8	MET-12139AU1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	S-10437T1-1ML	1ML
Azobenzene	103-33-3	N-11110-500MG	500MG	Benzoic acid-benzyl ester Solution	120-51-4	S-FC102S-1.2ML	1.2ML
Azobenzene solution	103-33-3	S-11110M1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	S-FC103S-2X1.2ML	2X1.2ML
Azobenzene solution	103-33-3	S-11110M1-5ML	5ML	Benzoic acid-benzyl ester Solution	120-51-4	S-FC105S-1.2ML	1.2ML
Azocyclopin (Technical)	41083-11-8	N-11111-250MG	250MG	Benzoic acid-benzyl ester Solution	120-51-4	N-11194-100MG	100MG
Azocyclopin (Technical) solution	41083-11-8	S-11111A1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	S-11194U1-1ML	1ML
Azomat	3586-60-5	N-13904-50MG	50MG	Benzoic acid-benzyl ester Solution	120-51-4	N-11195-100MG	100MG
Azoxystrobin	131860-33-8	N-11113-100MG	100MG	Benzoic acid-benzyl ester Solution	120-51-4	S-11195U1-1ML	1ML
Azoxystrobin solution	131860-33-8	S-11113A1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	S-11195U1-5ML	5ML
Azoxystrobin solution	131860-33-8	S-11113U1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	N-11196-50MG	50MG
Balance	141112-29-0	N-12287-100MG	100MG	Benzoic acid-benzyl ester Solution	120-51-4	S-11196U1-1ML	1ML
Balance Solution	141112-29-0	S-12287U1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	S-11196U1-5ML	5ML
Banol	671-04-5	N-11122-100MG	100MG	Benzoic acid-benzyl ester Solution	120-51-4	N-11198-250MG	250MG
Banol Solution	671-04-5	S-11122A1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	S-11198A1-1ML	1ML
Banol Solution	671-04-5	S-11122U1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	S-11198U1-1ML	1ML
Banvel (TM)	1918-00-9	N-11656-250MG	250MG	Benzoic acid-benzyl ester Solution	120-51-4	N-11201-50MG	50MG
Banvel (TM) Solution	1918-00-9	S-11656A1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	S-11201A1-1ML	1ML
Banvel (TM) Solution	1918-00-9	S-11656B1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	S-11201U1-1ML	1ML
Banvel (TM) Solution	1918-00-9	S-11656B1-5ML	5ML	Benzoic acid-benzyl ester Solution	120-51-4	N-11202-50MG	250MG
BAP	1214-39-7	N-10956-250MG	250MG	Benzoic acid-benzyl ester Solution	120-51-4	S-11202M1-1ML	1ML
BAP Solution	1214-39-7	S-10956A1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	N-11203-100MG	100MG
Barban	101-27-9	N-11123-100MG	100MG	Benzoic acid-benzyl ester Solution	120-51-4	S-11203A1-1ML	1ML
Barban solution	101-27-9	S-11123M1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	S-11203T1-1ML	1ML
Barban solution	101-27-9	S-11123M1-5ML	5ML	Benzoic acid-benzyl ester Solution	120-51-4	N-11204-100MG	100MG
Barricade(tm)	29091-21-2	N-13096-100MG	100MG	Benzoic acid-benzyl ester Solution	120-51-4	S-11204A1-1ML	1ML
Basitac(TM)	55814-41-0	N-12372-100MG	100MG	Benzoic acid-benzyl ester Solution	120-51-4	S-11204T1-1ML	1ML
Basitac(TM) Solution	55814-41-0	S-12372M1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	N-11205-250MG	250MG
Baycarb (TM)	3766-81-2	N-11126-250MG	250MG	Benzoic acid-benzyl ester Solution	120-51-4	N-13187-250MG	250MG
Baycarb (TM) solution	3766-81-2	S-11126M1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	S-1205A1-1ML	1ML
Baycor	55179-31-2	N-11127-100MG	100MG	Benzoic acid-benzyl ester Solution	120-51-4	S-11205T1-1ML	1ML
Baycor Solution	55179-31-2	S-11127A1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	S-13187A1-1ML	1ML
Baycor Solution	55179-31-2	S-11127U1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	S-13187T1-1ML	1ML
Baygon	114-26-1	N-11128-250MG	250MG	Benzoic acid-benzyl ester Solution	120-51-4	N-10965-250MG	250MG
Baygon Solution	114-26-1	S-11128A1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	S-10965M1-1ML	1ML
Bayothrin	118712-89-3	N-13626-100MG	100MG	Benzoic acid-benzyl ester Solution	120-51-4	N-11206-250MG	250MG
Bayrusil	13593-03-8	N-13168-250MG	250MG	Benzoic acid-benzyl ester Solution	120-51-4	S-11206A1-1ML	1ML
Bayrusil Solution	13593-03-8	S-13168A1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	S-11206T1-1ML	1ML
Bayrusil Solution	13593-03-8	S-13168U1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	N-11207-500MG	500MG
Baytan(TM)	55219-65-3	N-11129-250MG	250MG	Benzoic acid-benzyl ester Solution	120-51-4	S-11207J1-1ML	1ML
Baytan(TM) Solution	55219-65-3	S-11129A1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	S-11207J1-5ML	5ML
Baytan(TM) Solution	55219-65-3	S-11129U1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	S-11207M4-1ML	1ML
Beta-Baythroid(TM)	68359-37-5	N-11191-250MG	250MG	Benzoic acid-benzyl ester Solution	120-51-4	S-11207M4-5ML	5ML
BDMC	672-99-1	N-11131-100MG	100MG	Benzoic acid-benzyl ester Solution	120-51-4	N-10280-250MG	250MG
BDMC solution	672-99-1	S-11131M1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	S-10280A1-1ML	1ML
BDMC solution	672-99-1	S-11131M1-5ML	5ML	Benzoic acid-benzyl ester Solution	120-51-4	S-10280M5-1ML	1ML
Beacon(TM)	86209-51-0	N-13090-100MG	100MG	Benzoic acid-benzyl ester Solution	120-51-4	N-11322-250MG	250MG
Beacon(TM) Solution	86209-51-0	S-13090A1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	MET-12708A-100MG	100MG
Beam™	41814-78-2	N-11132-50MG	50MG	Benzoic acid-benzyl ester Solution	120-51-4	MET-12497A-500MG	500MG
Beam™ Solution	41814-78-2	S-11132A1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	MET-12135A-500MG	500MG
Beam™ Solution	41814-78-2	S-11132T1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	MET-13800B-25MG	25MG
Beam™ Solution	41814-78-2	S-11132T1-5ML	5ML	Benzoic acid-benzyl ester Solution	120-51-4	N-12960-50MG	50MG
Benalaxyl	71626-11-4	N-11133-250MG	250MG	Benzoic acid-benzyl ester Solution	120-51-4	N-11503-500MG	500MG
Benalaxyl solution	71626-11-4	S-11133M1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	N-11323-100MG	100MG
Benazolin	3813-05-600	N-12901-100MG	100MG	Benzoic acid-benzyl ester Solution	120-51-4	N-11127-100MG	100MG
Benazolin-ethyl	25059-80-7	N-11134-250MG	250MG	Benzoic acid-benzyl ester Solution	120-51-4	S-11127A1-1ML	1ML
Benazolin-ethyl solution	25059-80-7	S-11134A1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	S-11127U1-1ML	1ML
Benazolin-ethyl solution	25059-80-7	S-11134T1-1ML	1ML	Benzoic acid-benzyl ester Solution	120-51-4	N-11324-100MG	100MG
Bendiocarb	22781-23-3	N-11135-250MG	250MG	Benzoic acid-benzyl ester Solution	120-51-4	S-11324U1-1ML	1ML

Description	CAS#	Part #	Size	Description	CAS#	Part #	Size
Bladex solution	21725-46-2	S-11327T4-5ML	5ML	Cabrio EG(TM)	175013-18-0	N-13143-100MG	100MG
Blazer methyl ester	50594-67-7	N-11028-100MG	100MG	Cabrio EG(TM) Solution	175013-18-0	S-13143A1-1ML	1ML
Blazer methyl ester Solution	50594-67-7	S-11028Y1-1ML	1ML	Cacodylic acid	75-60-5	N-11779-500MG	500MG
Blazer methyl ester Solution	50594-67-7	S-11028Y1-5ML	5ML	Cacodylic acid Solution	75-60-5	S-11779F4-1ML	1ML
Boltage	89784-60-1	N-13142-10MG	10MG	Cadmium succinate (Technical)	141-00-4	N-11391-250MG	250MG
Boltage Solution	89784-60-1	S-13142A1-1ML	1ML	Cadusafos	95465-99-9	N-11392-100MG	100MG
Boscalid	188425-85-6	N-11328-100MG	100MG	Cadusafos solution	95465-99-9	S-11392A1-1ML	1ML
Boscalid solution	188425-85-6	S-11328A1-1ML	1ML	Cadusafos solution	95465-99-9	S-11392U1-1ML	1ML
BPMC	3766-81-2	N-11126-250MG	250MG	Calixin(TM)	81412-43-3	N-13665-100MG	100MG
BPMC Solution	3766-81-2	S-11126M1-1ML	1ML	Calixin(TM) Solution	81412-43-3	S-13665A1-1ML	1ML
Brace	42509-80-8	N-12222-100MG	100MG	Calixin(TM) Solution	81412-43-3	S-13665U1-1ML	1ML
Brace Solution	42509-80-8	S-12222U1-1ML	1ML	caloxydim Solution	149979-41-9	S-13970A1-1ML	1ML
Brodifacoum	56073-10-0	N-11329-50MG	50MG	Camphechlor	8001-35-2	N-13586-250MG	250MG
Brodifacoum solution	56073-10-0	S-11329A1-1ML	1ML	Camphechlor solution	8001-35-2	S-13586B6-1ML	1ML
Brofenprox	111872-58-3	N-13079-10MG	10MG	Camphechlor solution	8001-35-2	S-13586J4-1ML	1ML
Brofenprox Solution	111872-58-3	S-13079M1-1ML	1ML	Camphechlor solution	8001-35-2	S-13586J4-5ML	5ML
Bromacil	314-40-9	N-11330-250MG	250MG	Camphechlor solution	8001-35-2	S-13586M1-1ML	1ML
Bromacil solution	314-40-9	S-11330A1-1ML	1ML	Camphechlor solution	8001-35-2	S-13586M1-5ML	5ML
Bromacil solution	314-40-9	S-11330T1-1ML	1ML	Capsaicin, Natural	404-86-4	N-11398-100MG	100MG
Bromadiolol	28772-56-7	N-13259-100MG	100MG	Capsaicin, Natural solution	404-86-4	S-11398A1-1ML	1ML
Bromadiolol Solution	28772-56-7	S-13259A1-1ML	1ML	Capsaicin, Natural solution	404-86-4	S-11398U1-1ML	1ML
Bromadiolol Solution	28772-56-7	S-13259U1-1ML	1ML	Capstar	150824-47-8	N-12657-100MG	100MG
Bromoacetyl bromide	598-21-0	N-11333-500MG	500MG	Captafol	2425-06-100	N-11399-250MG	250MG
Bromoacetyl bromide solution	598-21-0	S-11333A1-1ML	1ML	Captan	133-06-2	N-11400-250MG	250MG
Bromoacetyl bromide solution	598-21-0	S-11333U1-1ML	1ML	Captan solution	133-06-2	S-11400U1-1ML	1ML
Bromobutide solution	74712-19-9	S-13977M1-1ML	1ML	Captan solution	133-06-2	S-11400U1-5ML	5ML
4-Bromo-2-chlorophenol	3964-56-5	MET-13097A-500MG	500MG	Capture	82657-04-3	N-11203-100MG	100MG
Bromoconazole	116255-48-2	N-11347-100MG	100MG	Capture Solution	82657-04-3	S-11203A1-1ML	1ML
(2-Bromoethyl)trimethyl ammonium bromide	2758-06-700	N-10001-1G	1G	Capture Solution	82657-04-3	S-11203T1-1ML	1ML
Bromofenoxim	13181-17-4	N-12961-10MG	10MG	Carbamult (TM)	2631-37-0	N-13101-100MG	100MG
De-Bromo-halifenprox		MET-13079A-10MG	10MG	Carbamult (TM) Solution	2631-37-0	S-13101M4-1ML	1ML
p-Bromophenoxy acetic acid	1878-91-7	N-12757-250MG	250MG	Carbanolate	671-04-5	N-11122-100MG	100MG
p-Bromophenoxy acetic acid solution	1878-91-7	S-12757A1-1ML	1ML	Carbanolate solution	671-04-5	S-11122A1-1ML	1ML
p-Bromophenoxy acetic acid solution	1878-91-7	S-12757T1-1ML	1ML	Carbanolate solution	671-04-5	S-11122U1-1ML	1ML
Bromophos	2104-96-3	N-11343-50MG	50MG	Carbaryl	63-25-2	N-11402-250MG	250MG
Bromophos Solution	2104-96-3	S-11343A1-1ML	1ML	Carbaryl (ring-13C6) Solution		S-FC2005S-1.2ML	1.2ML
Bromophos Solution	2104-96-3	S-11343J1-1ML	1ML	Carbaryl solution	63-25-2	S-11402A4-1ML	1ML
Bromophos ethyl	4824-78-6	N-11342-100MG	100MG	Carbaryl solution	63-25-2	S-11402A4-5ML	5ML
Bromophos ethyl solution	4824-78-6	S-11342J1-1ML	1ML	Carbendazim	10605-21-7	N-11404-100MG	100MG
Bromophos methyl	2104-96-3	N-11343-50MG	50MG	Carbendazim solution	10605-21-7	S-11404M1-1ML	1ML
Bromophos methyl solution	2104-96-3	S-11343A1-1ML	1ML	Carbendazim solution	10605-21-7	S-11404M1-5ML	5ML
Bromophos methyl solution	2104-96-3	S-11343J1-1ML	1ML	Carbetamide	16118-49-3	N-12971-10MG	10MG
Bromopropylate	18181-80-1	N-10997-250MG	250MG	Carbetoxy malaoxon	1634-78-2	MET-12346C-100MG	100MG
Bromopropylate solution	18181-80-1	S-10997A1-1ML	1ML	Carbofuran	1563-66-2	N-11405-250MG	250MG
Bromoxynil	1689-84-5	N-11345-250MG	250MG	Carbofuran (ring 13C6) Solution		S-FC2006S-1.2ML	1.2ML
Bromoxynil (ring 13C6) Solution		S-FC2093S-1.2ML	1.2ML	Carbofuran solution	1563-66-2	S-11405M1-1ML	1ML
Bromoxynil solution	1689-84-5	S-11345M1-1ML	1ML	Carbofuran-3-keto	16709-30-1	MET-11405E-100MG	100MG
Bromoxynil-heptanoate	56634-95-8	MET-11345A-50MG	50MG	Carbofuranphenol-3-keto	17781-16-7	MET-11405F-50MG	50MG
Bromoxynil-methyl ether	3336-39-8	MET-11345B-100MG	100MG	Carbofuranphenol-3-keto Solution	17781-16-7	MET-11405FM1-1ML	1ML
Bromoxynil octanoic acid ester	1689-99-2	N-11346-250MG	250MG	Carbophenothion	786-19-6	N-11408-250MG	250MG
Bromoxynil octanoic acid ester solution	1689-99-2	S-11346A1-1ML	1ML	Carbophenothion solution	786-19-6	S-11408M1-1ML	1ML
Bromoxynil octanoic acid ester solution	1689-99-2	S-11346K4-1ML	1ML	Carbophenothion methyl Solution	953-17-3	MET-11408AM1-1ML	1ML
Bromoxynil octanoic acid ester solution	1689-99-2	S-11346K4-5ML	5ML	Carbophenothion methyl-O-analog	7332-32-3	MET-11408B-50MG	50MG
Bromuconazole	116255-48-2	N-11347-100MG	100MG	Carbophenothion methyl-O-analog Solution	7173-84-4	MET-11408BH10-1ML	1ML
Bronopol	52-51-7	N-11348-250MG	250MG	Carbophenothion sulfone	16662-85-4	MET-11408D-25MG	25MG
Bronopol solution	52-51-7	S-11348M1-1ML	1ML	Carbophenothion sulfoxide	17297-40-4	MET-11408C-100MG	100MG
Bucid(TM)	10004-44-1	N-12196-25MG	25MG	Carbosulfan	55285-14-8	N-11409-250MG	250MG
Bucid(TM) Solution	10004-44-1	S-12196M1-1ML	1ML	Carbosulfan solution	55285-14-8	S-11409A1-1ML	1ML
Bupirimate	41483-43-6	N-11349-100MG	100MG	Carboxin	5234-68-4	N-11410-250MG	250MG
Buprofezin	69327-76-0	N-11350-100MG	100MG	Carboxin solution	5234-68-4	S-11410A1-1ML	1ML
Buprofezin solution	69327-76-0	S-11350A1-1ML	1ML	Cardomec(TM)	70288-86-7	N-12289-100MG	100MG
Buprofezin solution	69327-76-0	S-11350U1-1ML	1ML	Cardomec(TM) Solution	70288-86-7	S-12289A1-1ML	1ML
Busan 30A	21564-17-0	N-13500-10MG	10MG	Cardomec(TM) Solution	70288-86-7	S-12289U1-1ML	1ML
Busan 30A Solution	21564-17-0	S-13500A1-1ML	1ML	Carfentrazone-ethyl	128639-02-1	N-11411-10MG	10MG
Busan 30A Solution	21564-17-0	S-13500X4-1ML	1ML	Carfentrazone-ethyl solution	128639-02-1	S-11411A1-1ML	1ML
Busan 30A Solution	21564-17-0	S-13500X4-5ML	5ML	Carfentrazone-ethyl solution	128639-02-1	S-11411U1-1ML	1ML
Busan 40	51026-28-9	N-11351-100MG	100MG	Carproamid Solution Solution	104030-54-8	S-13974A1-1ML	1ML
Busan 40 solution	51026-28-9	S-11351A1-1ML	1ML	Cartap hydrochloride	15263-52-2	N-11412-250MG	250MG
Busan 40 solution	51026-28-9	S-11351F1-1ML	1ML	Cascade(TM)	101463-69-8	N-11987-50MG	50MG
Busan 85	128-03-0	N-11352-100MG	100MG	Cascade(TM) Solution	101463-69-8	S-11987A1-1ML	1ML
Busan 85 solution	128-03-0	S-11352A1-1ML	1ML	Celest(TM)	131341-86-1	N-11985-100MG	100MG
Busan 85 solution	128-03-0	S-11352F1-1ML	1ML	Celest(TM) Solution	131341-86-1	S-11985A1-1ML	1ML
Busan 85 solution	128-03-0	S-11352M1-1ML	1ML	Celest(TM) Solution	131341-86-1	S-11985U1-1ML	1ML
Butachlor	23184-66-9	N-11353-100MG	100MG	Chinomethionate	2439-01-200	N-13172-100MG	100MG
Butachlor solution	23184-66-9	S-11353A1-1ML	1ML	a-Chloralose(Technical)	15879-93-3	N-10978-250MG	250MG
Butachlor solution	23184-66-9	S-11353T1-1ML	1ML	b-Chloralose(Technical)	16376-36-6	N-11117-100MG	100MG
Butachlor solution	23184-66-9	S-11353T1-5ML	5ML	Chloramben	133-90-4	N-11418-100MG	100MG
Butachlor ESA sodium salt	1173022-75-7	MET-11353A-25MG	25MG	Chloramben solution	133-90-4	S-11418A1-1ML	1ML
Butachlor ESA sodium salt Solution	1173022-75-7	MET-11353AM1-1ML	1ML	Chloramben solution	133-90-4	S-11418B1-1ML	1ML
Butam(TM)	35256-85-0	N-13504-250MG	250MG	Chloramben solution	133-90-4	S-11418B1-5ML	5ML
Butam(TM) Solution	35256-85-0	S-13504A1-1ML	1ML	Chloramben ethyl ester	60541-86-8	N-11916-100MG	100MG
Butam(TM) Solution	35256-85-0	S-13504U1-1ML	1ML	Chloramben ethyl ester Solution	60541-86-8	S-11916A1-1ML	1ML
Butisan S(TM)	67129-08-2	N-12385-100MG	100MG	Chloramben ethyl ester Solution	60541-86-8	S-11916T1-1ML	1ML
Butocarbixim	34681-10-2	N-12963-100MG	100MG	Chloramben methyl ester	7286-84-2	N-11419-100MG	100MG
Butocarbixim solution	34681-10-2	S-12963M1-1ML	1ML	Chloramben methyl ester solution	7286-84-2	S-11419Y1-1ML	1ML
Butoxycarbixim	34681-23-7	N-13921-100MG	100MG	Chloramben methyl ester solution	7286-84-2	S-11419Y1-5ML	5ML
Butoxycarbixim solution	34681-23-7	S-13921M1-1ML	1ML	Chloranil	118-75-2	N-11420-250MG	250MG
2-[2-Butoxyethoxy]ethyl thiocyanate	112-56-1	N-10260-250MG	250MG	Chloranil solution	118-75-2	S-11420M1-1ML	1ML
2-[2-Butoxyethoxy]ethyl thiocyanate solution	112-56-1	S-10260A1-1ML	1ML	Chloranil solution	118-75-2	S-11420U1-1ML	1ML
2-[2-Butoxyethoxy]ethyl thiocyanate solution	112-56-1	S-10260J4-1ML	1ML	Chloranocryl		N-11421-250MG	250MG
2-[2-Butoxyethoxy]ethyl thiocyanate solution	112-56-1	S-10260J4-5ML	5ML	Chlorantraniliprole	500008-45-7	N-11422-100MG	100MG
Butoxyethyl-2,4,5-trichlorophenoxyacetate Solution	2545-59-7	S-10648A4-1ML	1ML	Chlorbenside	103-17-3	N-11423-100MG	100MG
Butoxyethyl-2,4,5-trichlorophenoxyacetate Solution	2545-59-7	S-10648A4-5ML	5ML	Chlorbenside solution	103-17-3	S-11423A1-1ML	1ML
Butoxyethyl-2,4,5-trichlorophenoxyacetate Solution	2545-59-7	S-10648U1-1ML	1ML	Chlorbenside solution	103-17-3	S-11423U1-1ML	1ML
Butoxypolypropylene glycol(Technical)	9003-13-8	N-11354-250MG	250MG	Chlorbenside sulfone	7082-99-7	MET-11423A-25MG	25MG
Butralin	33629-47-9	N-11355-500MG	500MG	Chlorbenside sulfone Solution	7082-99-7	MET-11423AJ1-1ML	1ML
Butralin solution	33629-47-9	S-11355A1-1ML	1ML	Chlorbenside sulfoxide	7047-28-1	MET-11423B-100MG	100MG
Butralin solution	33629-47-9	S-11355K1-1ML	1ML	Chlorbenside sulfoxide Solution	7047-28-1	MET-11423BX1-1ML	1ML
Buturon	3766-60-7	N-12944-10MG	10MG	Chlorbromuron	13360-45-7	N-11424-250MG	250MG
Butylate	2008-41-5	N-11383-250MG	250MG	Chlorbromuron solution	13360-45-7	S-11424A1-1ML	1ML
Butylate solution	2008-41-5	S-11383A1-1ML	1ML	Chlorbromuron solution	13360-45-7	S-11424U1-1ML	1ML
Butylate solution	2008-41-5	S-11383T1-1ML	1ML	Chlorbufam	1967-16-4	N-12970-10MG	10MG
Butylate solution	2008-41-5	S-11383T1-5ML	5ML	a-Chlordan	5103-71-9	N-11480-10MG	10MG
4-tert-Butyl-2-chlorophenol	98-28-2	MET-13186A-10MG	10MG	a-Chlordan Solution	5103-71-9	S-11480M1-1ML	1ML
1-n-Butyl-3-[3,4-dichlorophenyl]-1-methylurea Sol	555-37-3	S-12647A4-1ML	1ML	a-Chlordan Solution	5103-71-9	S-11480M1-5ML	5ML
1-n-Butyl-3-[3,4-dichlorophenyl]-1-methylurea Sol	555-37-3	S-12647A4-5ML	5ML	Chlordane	57-74-9	N-11425-250MG	250MG
2-(sec-Butyl)-4,6-dinitrophenol, ethanolamine salt		N-10888-1G	1G	trans-Chlordane (13C10) Solution		S-FC908S-1.2ML	1.2ML
3-Butyl-2,4-dioxo(1,2-a)-s-triazinobenzimidazole	132827-25-9	MET-11138B-500MG	500MG	Chlordane (random-13C1) Solution		S-FC91S-1.2ML	1.2ML
2-Butyl-2-ethyl-1,3-propanediol solution	115-84-4	N-10306-1G	1G	Chlordane solution	57-74-9	S-11425J4-1ML	1ML
2-Butyl-2-ethyl-1,3-propanediol solution	115-84-4	S-10306M1-1ML	1ML	Chlordane solution	57-74-9	S-11425J4-5ML	5ML
Butyl fluorenl	2314-09-200	N-13814-250MG	250MG	Chlordane solution	57-74-9	S-11425M1-1ML	1ML
Butylmesityl oxide oxalate	532-34-3	N-11385-250MG	250MG	beta-Chlordane	5103-74-2	N-13615-10MG	10MG
2-(4-tert-Butylphenoxy)cyclohexanol	1942-71-8	MET-12727A-100MG	100MG	beta-Chlordane Solution	5103-74-2	S-13615M1-1ML	1ML
1-(4-tert-Butylphenoxy)propan-2-ol	2416-30-0	MET-11085A-100MG	100MG	beta-Chlordane Solution	5103-74-2	S-13615M1-5ML	5ML

Description	CAS#	Part #	Size	Description	CAS#	Part #	Size
cis-Chlordane	5103-71-9	N-11480-10MG	10MG	Chlorotoluron solution	15545-48-9	S-11455U1-1ML	1ML
cis-Chlordane solution	5103-71-9	S-11480M1-1ML	1ML	Chlorotoluron, N-desmethyl	22175-22-0	MET-11455A-100MG	100MG
cis-Chlordane solution	5103-71-9	S-11480M1-5ML	5ML	4-Chloro-o-tolylloxycetic acid	94-74-6	N-10818-250MG	250MG
trans-Chlordane	5103-74-2	N-13615-10MG	10MG	4-Chloro-o-tolylloxycetic acid solution	94-74-6	S-10818B1-1ML	1ML
trans-Chlordane solution	5103-74-2	S-13615M1-1ML	1ML	4-Chloro-o-tolylloxycetic acid solution	94-74-6	S-10818B1-5ML	5ML
trans-Chlordane solution	5103-74-2	S-13615M1-5ML	5ML	4-Chloro-o-tolylloxycetic acid solution	94-74-6	S-10818M1-1ML	1ML
Chlordecone	143-50-0	N-12291-100MG	100MG	4-Chloro-o-tolylloxycetic acid methyl ester	2436-73-9	N-10819-100MG	100MG
Chlordecone solution	143-50-0	S-12291A1-1ML	1ML	4-Chloro-o-tolylloxycetic acid methyl ester solution	2436-73-9	S-10819M1-1ML	1ML
Chlordecone solution	143-50-0	S-12291J1-1ML	1ML	4-Chloro-o-tolylloxycetic acid methyl ester solution	2436-73-9	S-10819M1-5ML	5ML
Chlordecone solution	143-50-0	S-12291J1-5ML	5ML	2-Chloro-6-(trichloromethyl)pyridine Solution	1929-82-4	S-12659A1-1ML	1ML
Chlordecone solution	143-50-0	S-12291T4-1ML	1ML	Chloroxuron	1982-47-4	N-11456-250MG	250MG
Chlordecone solution	143-50-0	S-12291T4-5ML	5ML	Chloroxuron solution	1982-47-4	S-11456A1-1ML	1ML
Chlordecone	3734-48-3	MET-11425D-250MG	250MG	Chloroxynil	1891-95-8	N-11457-500MG	500MG
Chlordimeform	6164-98-3	N-11426-250MG	250MG	Chloroxynil solution	1891-95-8	S-11457M1-1ML	1ML
Chlordimeform solution	6164-98-3	S-11426A1-1ML	1ML	Chloropham	101-21-3	N-11458-250MG	250MG
Chlordimeform solution	6164-98-3	S-11426U1-1ML	1ML	Chloropham solution	101-21-3	S-11458A1-1ML	1ML
Chlorethoxyfos	54593-83-8	N-11427-250MG	250MG	Chloropham solution	101-21-3	S-11458T1-1ML	1ML
Chlorethoxyfos solution	54593-83-8	S-11427A1-1ML	1ML	Chloropham solution	101-21-3	S-11458T1-5ML	5ML
Chlorethoxyfos solution	54593-83-8	S-11427U1-1ML	1ML	Chlorpyrifos	2921-88-2	N-11459-250MG	250MG
Chlorfenac	85-34-7	N-11945-250MG	250MG	Chlorpyrifos (diethyl-d10) Solution	285138-81-0	S-FD207S-1.2ML	1.2ML
Chlorfenac Solution	85-34-7	S-11945A1-1ML	1ML	Chlorpyrifos solution	2921-88-2	S-11459A1-1ML	1ML
Chlorfenac Solution	85-34-7	S-11945T1-1ML	1ML	Chlorpyrifos solution	2921-88-2	S-11459K1-1ML	1ML
Chlorfenapyr	122453-73-0	N-11428-100MG	100MG	Chlorpyrifos solution	2921-88-2	S-11459K1-5ML	5ML
Chlorfenapyr solution	122453-73-0	S-11428A1-1ML	1ML	Chlorpyrifos Methyl	5598-13-0	N-11460-250MG	250MG
Chlorfenapyr solution	122453-73-0	S-11428U1-1ML	1ML	Chlorpyrifos Methyl solution	5598-13-0	S-11460A1-1ML	1ML
Chlorfenprop-methyl	14437-17-3	N-12898-50MG	50MG	Chlorpyrifos Methyl solution	5598-13-0	S-11460J4-1ML	1ML
Chlorfenvinphos	470-90-6	N-11429-250MG	250MG	Chlorpyrifos Methyl-O-analog	5598-52-7	MET-11460B-50MG	50MG
Chlorfenvinphos solution	470-90-6	S-11429K1-1ML	1ML	Chlorpyrifos Oxon	5598-15-2	MET-11459B-100MG	100MG
Chlorfluzazuron	71422-67-8	N-11430-50MG	50MG	Chlorpyrifos-o-analog	5598-15-2	MET-11459B-100MG	100MG
Chlorfluzazuron solution	71422-67-8	S-11430A1-1ML	1ML	Chlorosulfuron	64902-72-3	N-11461-100MG	100MG
Chlorflurecol-methyl ester	2536-31-4	N-11431-100MG	100MG	Chlorthal-dimethyl	1861-32-1	N-11462-250MG	250MG
Chlorflurecol-methyl ester solution	2536-31-4	S-11431M1-1ML	1ML	Chlorthal-dimethyl solution	1861-32-1	S-11462B1-1ML	1ML
Chlorimuron ethyl	90982-32-4	N-11432-100MG	100MG	Chlorthal-dimethyl solution	1861-32-1	S-11462B1-5ML	5ML
Chlorimuron ethyl solution	90982-32-4	S-11432A1-1ML	1ML	Chlorthal-dimethyl solution	1861-32-1	S-11462M1-1ML	1ML
Chlormephos	24934-91-6	N-11434-250MG	250MG	Chlorthiamid	1918-13-4	N-12992-10MG	10MG
Chlormephos solution	24934-91-6	S-11434A1-1ML	1ML	Chlorthiophos	60238-56-4	N-12969-100MG	100MG
Chlormephos solution	24934-91-6	S-11434U1-1ML	1ML	Cholecalciferol	67-97-0	N-11463-100MG	100MG
Chlormequat chloride	999-81-5	N-11435-250MG	250MG	Cholecalciferol solution	67-97-0	S-11463A1-1ML	1ML
Chlormequat chloride solution	999-81-5	S-11435F1-1ML	1ML	Cholecalciferol solution	67-97-0	S-11463U1-1ML	1ML
Chloroacetic acid solution	79-11-8	S-11438T1-1ML	1ML	Chromafenozide solution	143807-66-3	S-13973A1-1ML	1ML
Chloroacetic acid	79-11-8	N-11438-250MG	250MG	Cierto(TM)	98886-44-3	N-12020-30MG	30MG
Chlorobenzilate	510-15-6	N-11442-100MG	100MG	Cierto(TM) Solution	98886-44-3	S-12020A1-1ML	1ML
Chlorobenzilate solution	510-15-6	S-11442A1-1ML	1ML	Cierto(TM) Solution	98886-44-3	S-12020U1-1ML	1ML
Chlorobenzilate solution	510-15-6	S-11442J1-1ML	1ML	Cinnamethylin solution	87818-31-3	S-14050M1-1ML	1ML
Chlorobenzilate solution	510-15-6	S-11442J1-5ML	5ML	Cinnacure(TM)	104-55-2	N-11468-500MG	500MG
2-Chlorobenzonitrile	873-32-5	MET-11490A-250MG	250MG	Cinnacure(TM) Solution	104-55-2	S-11468A1-1ML	1ML
6-Chlorobenzoxazol-2(3H)-one	19932-84-4	MET-11957A-100MG	100MG	Cinnacure(TM) Solution	104-55-2	S-11468U1-1ML	1ML
1-Chloro-2,2-bis(p-chlorophenyl)ethene	1022-22-6	MET-12810B-10MG	10MG	Cinnamaldehyde	104-55-2	N-11468-500MG	500MG
4-Chloro-o-cresol	1570-64-5	MET-10818A-1G	1G	Cinnamaldehyde solution	104-55-2	S-11468A1-1ML	1ML
2-Chloro-4,6-diamino-1,3,5-triazine	3397-62-4	MET-13213C-100MG	100MG	Cinnamaldehyde solution	104-55-2	S-11468U1-1ML	1ML
2-Chloro-2',6'-diethylacetanilide	6967-29-9	MET-11043A-100MG	100MG	Cinosulfuron	94593-91-6	N-11470-100MG	100MG
6-Chloro-2,3-dihydrobenzoxazol-2-one	19932-84-4	MET-11957A-100MG	100MG	Clairst(TM)	594-42-3	N-12845-250MG	250MG
(2-Chloroethyl)phosphonic acid	16672-87-0	N-10002-250MG	250MG	Clethodim	99129-21-2	N-11488-100MG	100MG
Chloroethylmercury	107-27-7	N-11940-100MG	100MG	Clipper	219714-96-2	N-12823-50MG	50MG
2-Chloro-N-(2-ethyl-6-methylphenyl) acetamide	32428-71-0	MET-11013-50MG	50MG	Clodinafop-propargyl	105512-06-9	N-11489-100MG	100MG
Chlorofenson (TM)	80-33-1	N-12730-250MG	250MG	Clodinafop-propargyl solution	105512-06-9	S-11489A1-1ML	1ML
Chlorofenson (TM) Solution	80-33-1	S-12730A1-1ML	1ML	Clodinafop-propargyl solution	105512-06-9	S-11489U1-1ML	1ML
Chlorofenson (TM) Solution	80-33-1	S-12730J1-1ML	1ML	Clofentazine	74115-24-5	N-11490-100MG	100MG
p-Chloromandelic acid	492-86-4	N-12761-250MG	250MG	Clofentazine solution	74115-24-5	S-11490M1-1ML	1ML
p-Chloromandelic acid solution	492-86-4	S-12761A1-1ML	1ML	Clomazone	81777-89-1	N-11495-100MG	100MG
p-Chloromandelic acid solution	492-86-4	S-12761T1-1ML	1ML	Clomazone Solution	81777-89-1	S-11495M1-1ML	1ML
4-Chloro-2-methylaniline	95-69-2	MET-11426A-250MG	250MG	Clopropr	101-10-0	N-10254-100MG	100MG
4-Chloro-2-methylaniline Solution	95-69-2	MET-11426AM1-1ML	1ML	Clopropr Solution	101-10-0	S-10254T1-1ML	1ML
4-Chloro-2-methylphenol	1570-64-5	MET-10818A-1G	1G	Clopropr Solution	101-10-0	S-10254T1-1ML	1ML
4-Chloro-2-methylphenoxyacetic acid sodium salt	3653-48-3	MET-12348-250MG	250MG	Clopyralid	1702-17-6	N-12324-250MG	250MG
1-(3-Chloro-4-methylphenyl)-3-methylurea	22175-22-0	MET-11455A-100MG	100MG	Clopyralid solution	1702-17-6	S-12324A1-1ML	1ML
2-Chloro-4-(methylsulfonyl)benzoic acid	53250-83-2	MET-12890A-100MG	100MG	Clopyralid solution	1702-17-6	S-12324T1-1ML	1ML
Chloroneb	2675-77-6	N-11449-250MG	250MG	Cloquintocet	88349-88-6	N-13820-50MG	50MG
Chloroneb	2675-77-6	N-11449-250MG	250MG	Cloquintocet-mexyl	99607-70-2	N-11491-100MG	100MG
Chloroneb Solution	2675-77-6	S-11449A1-1ML	1ML	Cloransulam-methyl	147150-35-4	N-11492-100MG	100MG
Chloroneb Solution	2675-77-6	S-11449A1-1ML	1ML	Cloransulam-methyl solution	147150-35-4	S-11492A1-1ML	1ML
Chloroneb Solution	2675-77-6	S-11449T1-1ML	1ML	Cloransulam-methyl solution	147150-35-4	S-11492U1-1ML	1ML
Chloroneb solution	2675-77-6	S-11449T1-1ML	1ML	Clothianidin	210880-92-5	N-11493-100MG	100MG
Chloroneb Solution	2675-77-6	S-11449T1-5ML	5ML	Clothianidin solution	210880-92-5	S-11493A1-1ML	1ML
Chloroneb solution	2675-77-6	S-11449T1-5ML	5ML	Command	81777-89-1	N-11495-100MG	100MG
6-Chloronicotinic acid	5326-23-8	MET-12066B-100MG	100MG	Command solution	81777-89-1	S-11495M1-1ML	1ML
6-Chloronicotinic acid Solution	5326-23-8	MET-12066BM1-1ML	1ML	Compete(TM)	77501-90-7	N-11999-100MG	100MG
2-Chloro-4-nitrophenol	619-08-9	MET-11658A-250MG	250MG	Comply(TM)	72490-01-8	N-11959-100MG	100MG
Chlorophacinone	3691-35-8	N-11451-250MG	250MG	Comply(TM) Solution	72490-01-8	S-11959A1-1ML	1ML
Chlorophacinone solution	3691-35-8	S-11451A1-1ML	1ML	Confirm(TM)	112410-23-8	N-13501-100MG	100MG
o-Chlorophenoxy acetic acid	614-61-9	N-12679-250MG	250MG	Confirm(TM) Solution	112410-23-8	S-13501A1-1ML	1ML
o-Chlorophenoxy acetic acid solution	614-61-9	S-12679M1-1ML	1ML	Consult(TM)	86479-06-3	N-12167-250MG	250MG
o-Chlorophenoxy acetic acid solution	614-61-9	S-12679T1-1ML	1ML	Consult(TM) Solution	86479-06-3	S-12167A1-1ML	1ML
p-Chlorophenoxyacetic acid	122-88-3	N-12763-250MG	250MG	Consult(TM) Solution	86479-06-3	S-12167U1-1ML	1ML
p-Chlorophenoxyacetic acid solution	122-88-3	S-12763M1-1ML	1ML	Copper Naphthenates(Technical)	1338-02-9	N-11501-250MG	250MG
p-Chlorophenoxyacetic acid solution	122-88-3	S-12763T1-1ML	1ML	Copper oxychloride(Technical)	1332-65-6	N-11502-250MG	250MG
2-[3-Chlorophenoxy]propionamide	5825-87-6	N-10261-1G	1G	Copper-8-Quinolinate	10380-28-6	N-11503-500MG	500MG
2-[3-Chlorophenoxy]propionamide solution	5825-87-6	S-10261M1-1ML	1ML	Corbit	84-65-1	N-10970-250MG	250MG
2-[3-Chlorophenoxy]propionic acid	101-10-0	N-10254-100MG	100MG	Corbit Solution	84-65-1	S-10970M1-1ML	1ML
2-[3-Chlorophenoxy]propionic acid solution	101-10-0	S-10254A1-1ML	1ML	Coumaphos	56-72-4	N-11507-100MG	100MG
2-[3-Chlorophenoxy]propionic acid solution	101-10-0	S-10254T1-1ML	1ML	Coumaphos solution	56-72-4	S-11507A1-1ML	1ML
1-(4-Chlorophenyl)-3-(2,6-difluorobenzoyl)urea	35367-38-5	N-11722-250MG	250MG	Coumaphos solution	56-72-4	S-11507U1-1ML	1ML
1-(4-Chlorophenyl)-3-(2,6-difluorobenzoyl)urea Sol	35367-38-5	S-11722A1-1ML	1ML	Coumaphos solution	56-72-4	S-11507U1-5ML	5ML
2-(p-Chlorophenyl)-isovaleric acid	2012-74-0	MET-13201A-100MG	100MG	Coumaphos-O-analog	321-54-0	MET-11507D-50MG	50MG
Chlorophenyl mercury	100-56-1	N-13013-250MG	250MG	Coumatetralyl	5836-29-3	N-11509-250MG	250MG
2-(p-Chlorophenyl)-3-methyl butyric acid	2012-74-0	MET-13201A-100MG	100MG	Coumatetralyl solution	5836-29-3	S-11509M1-1ML	1ML
3-(4-Chlorophenyl)methyl urea	5352-88-5	MET-11827B-100MG	100MG	4-CPA	122-88-3	N-12763-250MG	250MG
p-Chlorophenyl-2,4,5-trichlorophenylsulfone Sol.	116-29-0	S-13539A1-1ML	1ML	4-CPA Solution	122-88-3	S-12763M1-1ML	1ML
p-Chlorophenyl-2,4,5-trichlorophenylsulfone Sol.	116-29-0	S-13539B3-1ML	1ML	4-CPA Solution	122-88-3	S-12763T1-1ML	1ML
p-Chlorophenyl-2,4,5-trichlorophenylsulfone Sol.	116-29-0	S-13539B3-5ML	5ML	Creatine	57-00-1	MET-11840A-250MG	250MG
1-(p-Chlorophenyl)urea	140-38-5	MET-11722A-100MG	100MG	Croneton sulfone	53380-22-6	MET-11872A-10MG	10MG
1-(4-Chlorophenyl)-1H-pyrazol-4-ol	77458-30-1	MET-13142A-100MG	100MG	Croneton sulfone	53380-23-7	MET-11872B-10MG	10MG
Chlorophenylurea	140-38-5	MET-11722A-100MG	100MG	Croneton sulfone	53380-23-7	MET-11872BA1-1ML	1ML
6-Chloro-2-picolinic acid methyl ester Solution	66356-55-1	MET-12659B1-1ML	1ML	Croneton*	29973-13-5	N-11872-100MG	100MG
Chloropicrin	76-06-2	N-11452-100MG	100MG	Crossbow(TM)	55335-06-3	N-13659-500MG	500MG
Chloropicrin Solution	76-06-2	S-11452B1-1ML	1ML	Crossbow(TM) Solution	55335-06-3	S-13659A1-1ML	1ML
Chloropicrin Solution	76-06-2	S-11452B1-5ML	5ML	Crossbow(TM) Solution	55335-06-3	S-13659T1-1ML	1ML
Chloropropylate	5836-10-200	N-11453-50MG	50MG	Crotoxphos	7700-17-6	N-11513-50MG	50MG
(1E)-N-[(6-Chloro-3-pyridinyl)methyl]-N'	190604-92-3	MET-111009A-10MG	10MG	Crotoxphos solution	7700-17-6	S-11513K1-1ML	1ML
Chlorothalonil	1897-45-6	N-11454-250MG	250MG	Crotoxphos solution	7700-17-6	S-11513K1-5ML	5ML
Chlorothalonil solution	1897-45-6	S-11454A1-1ML	1ML	Cruformate	299-86-5	N-13186-50MG	50MG
Chlorothalonil solution	1897-45-6	S-11454T1-1ML	1ML	Cruformate Solution	299-86-5	S-13186A1-1ML	1ML
Chlorotoluron	15545-48-9	N-11455-250MG	250MG	Cruformate Solution	299-86-5	S-13186U1-1ML	1ML
Chlorotoluron solution	15545-48-9	S-11455A1-1ML	1ML	Cutless*	56425-91-3	N-12002-100MG	100MG

Description	CAS#	Part #	Size	Description	CAS#	Part #	Size
Cutless* Solution	56425-91-3	S-12002M1-1ML	1ML	2,4-D butoxyethyl ester	1929-73-3	N-10610-250MG	250MG
Cyhalofop-amide		MET-11543A-50MG	50MG	2,4-D butoxyethyl ester solution	1929-73-3	S-10610A1-1ML	1ML
Cyhalofop-diacid	252564-94-6	MET-14038A-100MG	100MG	2,4-D butoxyethyl ester solution	1929-73-3	S-10610A4-1ML	1ML
Cyanazine	21725-46-2	N-11327-250MG	250MG	2,4-D butoxyethyl ester solution	1929-73-3	S-10610A4-5ML	5ML
Cyanazine Solution	21725-46-2	S-11327A1-1ML	1ML	2,4-D butyl ester	94-80-4	N-10611-250MG	250MG
Cyanazine Solution	21725-46-2	S-11327T4-1ML	1ML	2,4-D sec-butyl ester	94-79-1	N-10621-250MG	250MG
Cyanazine Solution	21725-46-2	S-11327T4-5ML	5ML	2,4-D butyl ester solution	94-80-4	S-10611A1-1ML	1ML
Cyanofenphos	13067-93-1	N-12881-50MG	50MG	2,4-D butyl ester solution	94-80-4	S-10611U1-1ML	1ML
Cyanofenphos solution	13067-93-1	S-12881M1-1ML	1ML	2,4-D sec-butyl ester solution	94-75-7	S-10621A4-1ML	1ML
RS-a-Cyanophenoxybenzyl-1RS-cis	52315-07-8	N-11545-100MG	100MG	2,4-D sec-butyl ester solution	94-75-7	S-10621A4-5ML	5ML
RS-a-Cyanophenoxybenzyl-1RS-cis Solution	52315-07-8	S-11545A1-1ML	1ML	2,4-D Choline salt	94-75-7	S-10621U1-1ML	1ML
RS-a-Cyanophenoxybenzyl-1RS-cis Solution	52315-07-8	S-11545U1-1ML	1ML	2,4-D Choline salt solution	1048373-72-3	N-13985-100MG	100MG
RS-a-Cyanophenoxybenzyl-1RS-cis Solution	52315-07-8	S-13754A1-1ML	1ML	2,4-D Choline salt solution	1048373-72-3	S-13985F1-1ML	1ML
RS-a-Cyanophenoxybenzyl-1RS-cis Solution	52315-07-8	S-13754U1-1ML	1ML	2,4-D dimethylamine salt	2008-39-1	N-10612-250MG	250MG
O-4-Cyanophenyl	2636-26-2	N-11517-100MG	100MG	2,4-D dimethylamine salt solution	2008-39-1	S-10612M1-1ML	1ML
O-4-Cyanophenyl Solution	2636-26-2	S-11517A1-1ML	1ML	2,4-D ethanolamine salt	3599-58-4	N-10613-250MG	250MG
O-4-Cyanophenyl Solution	2636-26-2	S-11517U1-1ML	1ML	2,4-D ethanolamine salt solution	3599-58-4	S-10613M1-1ML	1ML
Cyanophos	2636-26-2	N-11517-100MG	100MG	2,4-D ethyl ester	533-23-3	N-10614-250MG	250MG
Cyanophos solution	2636-26-2	S-11517A1-1ML	1ML	2,4-D ethyl ester solution	533-23-3	S-10614A1-1ML	1ML
Cyanophos solution	2636-26-2	S-11517U1-1ML	1ML	2,4-D ethyl ester solution	533-23-3	S-10614J1-1ML	1ML
Cyantraniliprole	736994-63-1	N-12886-10MG	10MG	2,4-D 2-ethylhexyl ester solution	1928-43-4	S-10531A4-1ML	1ML
Cyazofamid	120116-88-3	N-11519-100MG	100MG	2,4-D 2-ethylhexyl ester solution	1928-43-4	S-10531A4-5ML	5ML
Cybutryne	28159-98-0	N-11520-100MG	100MG	2,4-D isobutyl ester	1713-15-1	N-10532-250MG	250MG
Cybutryne solution	28159-98-0	S-11520A1-1ML	1ML	2,4-D isobutyl ester solution	1713-15-1	S-10532A1-1ML	1ML
Cybutryne solution	28159-98-0	S-11520U1-1ML	1ML	2,4-D isobutyl ester solution	1713-15-1	S-10532U1-1ML	1ML
Cyclanilide	113136-77-9	N-11521-100MG	100MG	2,4-D isooctyl ester	25168-26-7	N-10615-250MG	250MG
Cyclanilide solution	113136-77-9	S-11521A1-1ML	1ML	2,4-D isooctyl ester solution	25168-26-7	S-10615A1-1ML	1ML
Cyclanilide solution	113136-77-9	S-11521U1-1ML	1ML	2,4-D isooctyl ester solution	25168-26-7	S-10615T1-1ML	1ML
Cycloate	1134-23-2	N-11522-250MG	250MG	2,4-D isopropanol amine salt	6365-72-6	N-10616-250MG	250MG
Cycloate solution	1134-23-2	S-11522A1-1ML	1ML	2,4-D isopropanol amine salt solution	6365-72-6	S-10616M1-1ML	1ML
Cycloate solution	1134-23-2	S-11522T1-1ML	1ML	2,4-D isopropylamine salt	5742-17-6	N-10618-250MG	250MG
Cycloate solution	1134-23-2	S-11522T1-5ML	5ML	2,4-D isopropylamine salt solution	5742-17-6	S-10618M1-1ML	1ML
1,2-Cyclohexanediol (cis & trans mixture)	931-17-9	MET-12727B-250MG	250MG	2,4-D isopropyl ester	94-11-1	N-10617-250MG	250MG
cis-4-Cyclohexene-1,2-dicarboximide	1469-48-3	MET-11399A-250MG	250MG	2,4-D isopropyl ester solution	94-11-1	S-10617A1-1ML	1ML
Cycloheximide	66-81-9	N-11534-100MG	100MG	2,4-D isopropyl ester solution	94-11-1	S-10617U1-1ML	1ML
Cycloheximide solution	66-81-9	S-11534A1-1ML	1ML	2,4-D methyl ester	1928-38-7	N-10619-250MG	250MG
Cycloheximide solution	66-81-9	S-11534U1-1ML	1ML	2,3-D methyl ester	103-26-4	N-12861-100MG	100MG
Cyclopentadiene dimer	77-73-6	N-11686-500MG	500MG	2,4-D methyl ester solution	1928-38-7	S-10619M1-1ML	1ML
Cyclopentadiene dimer Solution	77-73-6	N-11686K1-1ML	1ML	2,4-D methyl ester solution	1928-38-7	S-10619M1-5ML	5ML
Cycloprate	54460-46-7	N-11541-500MG	500MG	2,4-D propylene glycol ester	103564-78-9	N-10620-250MG	250MG
Cycloprate solution	54460-46-7	S-11541A1-1ML	1ML	2,4-D propylene glycol ester solution	103564-78-9	S-10620A1-1ML	1ML
Cycloprate solution	54460-46-7	S-11541T4-1ML	1ML	2,4-D propylene glycol ester solution	103564-78-9	S-10620J1-1ML	1ML
Cycloprate solution	54460-46-7	S-11541T4-5ML	5ML	Daimuron	42609-52-9	N-13955-100MG	100MG
Cyclosulfamuron	136849-15-5	N-11542-250MG	250MG	Daimuron solution	42609-52-9	S-13955A1-1ML	1ML
Cyclosulfamuron solution	136849-15-5	S-11542A1-1ML	1ML	Dalapon	75-99-0	N-11562-250MG	250MG
Cyclosulfamuron solution	136849-15-5	S-11542U1-1ML	1ML	Dalapon solution	75-99-0	S-11562B1-1ML	1ML
Cycloxydimide	101205-02-1	N-13916-10MG	10MG	Dalapon solution	75-99-0	S-11562M1-1ML	1ML
Cyflufenamid	180409-60-3	N-13132-25MG	25MG	Dalapon methyl ester	17640-02-7	N-11563-100MG	100MG
Cyflumetofen	400882-07-7	N-13972-50MG	50MG	Dalapon methyl ester solution	17640-02-7	S-11563M1-1ML	1ML
Cyfluthrin	68359-37-5	N-11130-250MG	250MG	Dalapon Sodium	127-20-8	N-11564-250MG	250MG
beta-Cyfluthrin	68359-37-5	N-11191-250MG	250MG	Damfin(TM)	30864-28-9	N-12388-100MG	100MG
Cyfluthrin solution	68359-37-5	S-11130A4-1ML	1ML	Damfin(TM) Solution	30864-28-9	S-12388A1-1ML	1ML
Cyfluthrin solution	68359-37-5	S-11130A4-5ML	5ML	Damfin(TM) Solution	30864-28-9	S-12388U1-1ML	1ML
Cyfluthrin solution	68359-37-5	S-11130U1-1ML	1ML	Daminozide	1596-84-5	N-11565-250MG	250MG
beta-Cyfluthrin solution	68359-37-5	S-11191A1-1ML	1ML	Daminozide solution	1596-84-5	S-11565A1-1ML	1ML
beta-Cyfluthrin solution	68359-37-5	S-11191T1-1ML	1ML	Daminozide solution	1596-84-5	S-11565T1-1ML	1ML
Cyfoxylate	68359-37-5	N-11130-250MG	250MG	Danienon(TM)	148477-71-8	N-13223-100MG	100MG
Cyfoxylate Solution	68359-37-5	S-11130A4-1ML	1ML	Danisaraba	400882-07-7	N-13972-50MG	50MG
Cyfoxylate Solution	68359-37-5	S-11130A4-5ML	5ML	Danitol	39515-41-8	N-11960-250MG	250MG
Cyfoxylate Solution	68359-37-5	S-11130U1-1ML	1ML	Danitol Solution	39515-41-8	S-11960A1-1ML	1ML
Cyfoxylate Solution	68359-37-5	S-11191A1-1ML	1ML	Danitol Solution	39515-41-8	S-11960T1-1ML	1ML
Cyfoxylate Solution	68359-37-5	S-11191T1-1ML	1ML	Dazomet	533-74-4	N-11566-250MG	250MG
Cyhalofop	122008-78-0	N-14038-25MG	25MG	Dazomet solution	533-74-4	S-11566M1-1ML	1ML
Cyhalofop solution	122008-78-0	S-14038A1-1ML	1ML	2,4-DB	94-82-6	N-10622-250MG	250MG
Cyhalofop-butyl	122008-85-9	N-11543-50MG	50MG	2,4-DB solution	94-82-6	S-10622A1-1ML	1ML
Cyhalofop-butyl solution	122008-85-9	S-11543A1-1ML	1ML	2,4-DB solution	94-82-6	S-10622B1-1ML	1ML
Cyhalofop-butyl solution	122008-85-9	S-11543U1-1ML	1ML	2,4-DB solution	94-82-6	S-10622B1-5ML	5ML
lambda-Cyhalothrin	91465-08-6	N-12307-100MG	100MG	2,4-DB, dimethyl ammonium Solution	2758-42-1	S-10623M1-1ML	1ML
lambda-Cyhalothrin solution	91465-08-6	S-12307A1-1ML	1ML	2,4-DB dimethylamine salt	2758-42-1	N-10623-250MG	250MG
lambda-Cyhalothrin solution	91465-08-6	S-12307J1-1ML	1ML	2,4-DB dimethylamine salt solution	2758-42-1	S-10623M1-1ML	1ML
Cyhexatin	13121-70-5	N-13662-250MG	250MG	2,4-DB methyl ester	18625-12-2	N-10533-100MG	100MG
Cyhexatin Solution	13121-70-5	S-13662M1-1ML	1ML	2,4-DB methyl ester solution	18625-12-2	S-10533M1-1ML	1ML
Cymetrin	1014-70-6	N-13802-50MG	50MG	2,4-DB methyl ester solution	18625-12-2	S-10533M1-5ML	5ML
Cymetrin Solution	1014-70-6	S-13802A1-1ML	1ML	DBCP	96-12-8	N-10149-250MG	250MG
Cymetrin Solution	1014-70-6	S-13802T1-1ML	1ML	DBCP Solution	96-12-8	S-10149M1-1ML	1ML
Cymetrin Solution	1014-70-6	S-13802T1-5ML	5ML	DBCP Solution	96-12-8	S-10149M1-5ML	5ML
Cymiazole solution	61676-87-7	S-14051A1-1ML	1ML	DBP Solution	84-74-2	S-11589A1-1ML	1ML
Cymoxanil	57966-95-7	N-11544-100MG	100MG	DBP Solution	84-74-2	S-11589J1-1ML	1ML
Cymoxanil solution	57966-95-7	S-11544A1-1ML	1ML	DBP Solution	84-74-2	S-11589J1-5ML	5ML
Cymoxanil solution	57966-95-7	S-11544U1-1ML	1ML	DCNA	99-30-9	N-11678-250MG	250MG
Cypermethrin	52315-07-8	N-11545-100MG	100MG	DCNA Solution	99-30-9	S-11678A1-1ML	1ML
Cypermethrin solution	52315-07-8	S-11545A1-1ML	1ML	DCNA Solution	99-30-9	S-11678U4-1ML	1ML
Cypermethrin solution	52315-07-8	S-11545U1-1ML	1ML	DCPA	1861-32-1	N-11462-250MG	250MG
alpha-Cypermethrin	67375-30-8	N-11061-100MG	100MG	DCPA Solution	1861-32-1	S-11462B1-1ML	1ML
alpha-Cypermethrin solution	67375-30-8	S-11061A1-1ML	1ML	DCPA Solution	1861-32-1	S-11462B1-5ML	5ML
alpha-Cypermethrin solution	67375-30-8	S-11061J1-1ML	1ML	DCPA Solution	1861-32-1	S-11462M1-1ML	1ML
beta-Cypermethrin	65731-84-2	N-11192-250MG	250MG	p,p'-DDA	83-05-6	MET-12708A-100MG	100MG
beta-Cypermethrin solution	65731-84-2	S-11192A1-1ML	1ML	DDCA	53179-78-5	MET-11579A-10MG	10MG
beta-Cypermethrin solution	65731-84-2	S-11192T1-1ML	1ML	DDCA	53179-78-5	MET-11579BM1-1ML	1ML
zeta-Cypermethrin	52315-07-8	N-13754-250MG	250MG	4,4'-DDD (ring-d8) Solution	72-54-8	S-FD94S-1.2ML	1.2ML
zeta-Cypermethrin solution	52315-07-8	S-13754A1-1ML	1ML	4,4'-DDD	72-54-8	N-10874-250MG	250MG
d-trans-Cyphenothrin	39515-40-7	N-11561-100MG	100MG	4,4'-DDD solution	72-54-8	S-10874A1-1ML	1ML
d-trans-Cyphenothrin solution	39515-40-7	S-11561A1-1ML	1ML	4,4'-DDD solution	72-54-8	S-10874K1-1ML	1ML
d-trans-Cyphenothrin solution	39515-40-7	S-11561T1-1ML	1ML	4,4'-DDD solution	72-54-8	S-10874K1-5ML	5ML
Cyprazine	22936-86-3	N-11546-50MG	50MG	o,p'-DDD	53-19-0	N-12706-250MG	250MG
Cyprazine solution	22936-86-3	S-11546A1-1ML	1ML	o,p'-DDD solution	53-19-0	S-12706M1-1ML	1ML
Cyprazine solution	22936-86-3	S-11546B1-1ML	1ML	o,p'-DDD olefin	3424-82-6	MET-12706A-10MG	10MG
Cyproconazole	94361-06-5	N-11547-50MG	50MG	o,p'-DDE	3424-82-6	N-12707-50MG	50MG
Cyproconazole solution	94361-06-5	S-11547A1-1ML	1ML	4,4'-DDE (13C12) Solution		S-FC93S-1.2ML	1.2ML
Cyproconazole solution	94361-06-5	S-11547E1-1ML	1ML	4,4'-DDE (13C12)		N-FC93-5-5MG	5MG
Cyprodinil	121552-61-2	N-11548-250MG	250MG	4,4'-DDE solution	72-55-9	S-10875M1-1ML	1ML
Cyprodinil solution	121552-61-2	S-11548A1-1ML	1ML	4,4'-DDE solution	72-55-9	S-10875M1-5ML	5ML
Cyprodinil solution	121552-61-2	S-11548U1-1ML	1ML	o,p'-DDE (13C12) Solution		S-FC2306S-1.2ML	1.2ML
Cyrosulfamide	221667-31-8	N-11549-100MG	100MG	4,4'-DDE	72-55-9	N-10875-100MG	100MG
Cyromazine	66215-27-8	N-11550-250MG	250MG	o,p'-DDE solution	3424-82-6	S-12707M1-5ML	5ML
Cyromazine solution	66215-27-8	S-11550A1-1ML	1ML	o,p'-DDE solution	3424-82-6	S-12707M1-1ML	1ML
Cyromazine solution	66215-27-8	S-11550B1-1ML	1ML	p,p'-DDE	72-55-9	N-10875-100MG	100MG
Cythioate	115-93-5	N-11551-50MG	50MG	p,p'-DDE Solution	72-55-9	S-10875M1-1ML	1ML
Cythioate solution	115-93-5	S-11551A1-1ML	1ML	p,p'-DDE Solution	72-55-9	S-10875M1-5ML	5ML
Cythioate solution	115-93-5	S-11551U1-1ML	1ML	p,p'-DDMU	1022-22-6	MET-12810B-10MG	10MG
2,4-D	94-75-7	N-10609-250MG	250MG	o,p'-DDT (13C12) Solution		S-FC2307S-1.2ML	1.2ML
2,6-D Acid	575-90-6	N-12996-10MG	10MG	4,4'-DDT (13C12) Solution		S-FC92S-1.2ML	1.2ML
2,4-D solution	94-75-7	S-10609B1-1ML	1ML	DDT (mixture p,p' & o,p')	8017-34-3	N-11567-250MG	250MG
2,4-D solution	94-75-7	S-10609M1-1ML	1ML	DDT (mixture p,p' & o,p') solution	8017-34-3	S-11567M1-1ML	1ML

Description	CAS#	Part #	Size	Description	CAS#	Part #	Size
4,4'-DDT	50-29-3	N-10876-100MG	100MG	2,4-Dichlorobenzhydrol	43171-49-9	MET-11948A-100MG	100MG
4,4'-DDT (13C12)	50-29-3	N-FC92-5-5MG	5MG	2,6-Dichlorobenzoic acid	50-30-6	MET-10688A-250MG	250MG
4,4'-DDT solution	789-02-6	S-10876J1-1ML	1ML	3,5-Dichlorobenzoic acid	51-36-5	MET-10770-250MG	250MG
o,p'-DDT	789-02-6	N-12708-50MG	50MG	3,5-Dichlorobenzoic acid solution	51-36-5	S-10770A1-1ML	1ML
o,p'-DDT solution	789-02-6	S-12708M1-1ML	1ML	3,5-Dichlorobenzoic acid solution	51-36-5	S-10770T5-1ML	1ML
o,p'-DDT solution	789-02-6	S-12708M1-5ML	5ML	3,5-Dichlorobenzoic acid solution	51-36-5	S-10770T5-5ML	5ML
cis-Decamethrinic acid methyl ester	113830-50-5	MET-11579C-10MG	10MG	2,6-Dichlorobenzoic acid methyl ester	14920-87-7	MET-10688C-250MG	250MG
DEET	134-62-3	N-12618-250MG	250MG	3,5-Dichlorobenzoic acid methyl ester	2905-67-1	N-10791-100MG	100MG
DEET Solution	134-62-3	S-12618A1-1ML	1ML	3,5-Dichlorobenzoic acid methyl ester solution	2905-67-1	S-10791A1-1ML	1ML
DEET Solution	134-62-3	S-12618B5-1ML	1ML	3,5-Dichlorobenzoic acid methyl ester solution	2905-67-1	S-10791T5-1ML	1ML
DEET Solution	134-62-3	S-12618B5-5ML	5ML	3,5-Dichlorobenzoic acid methyl ester solution	2905-67-1	S-10791T5-5ML	5ML
DEF	78-48-8	N-13194-250MG	250MG	2,6-Dichlorobenzonitrile	1194-65-6	N-10688-250MG	250MG
DEF Solution	78-48-8	S-13194A1-1ML	1ML	2,6-Dichlorobenzonitrile solution	1194-65-6	S-10688M1-1ML	1ML
DEF Solution	78-48-8	S-13194T1-1ML	1ML	4,4'-Dichlorobenzophenone	90-98-2	MET-11442B-100MG	100MG
Deflor(TM)	19937-59-8	N-12480-250MG	250MG	2,4'-Dichlorobenzophenone	85-29-0	MET-11948-100MG	100MG
Deflor(TM) Solution	19937-59-8	S-12480A1-1ML	1ML	3,4-Dichlorocarbonyl acid methyl ester Solution	1918-18-9	S-13260A4-1ML	1ML
Deflor(TM) Solution	19937-59-8	S-12480U1-1ML	1ML	3,4-Dichlorocarbonyl acid methyl ester Solution	1918-18-9	S-13260A4-5ML	5ML
Deltamethrin	52918-63-5	N-11579-250MG	250MG	3,4-Dichlorocarbonyl acid methyl ester Solution	1918-18-9	S-13260U1-1ML	1ML
Deltamethrin solution	52918-63-5	S-11579A1-1ML	1ML	Dichlorodiphenyltrichloroethane	50-29-3	N-10876-100MG	100MG
Deltamethrin solution	52918-63-5	S-11579T1-1ML	1ML	5,7-Dichloro-4-(4-fluorophenoxy)quinoline	124495-18-7	N-13173-100MG	100MG
Demeton-S-methyl	919-86-8	N-12384-50MG	50MG	5,7-Dichloro-4-(4-fluorophenoxy)quinoline Solution	124495-18-7	S-13173A1-1ML	1ML
Demeton-S-methyl Solution	919-86-8	S-12384A1-1ML	1ML	5,7-Dichloro-4-(4-fluorophenoxy)quinoline Solution	124495-18-7	S-13173U1-1ML	1ML
Demeton-S-methyl Solution	919-86-8	S-12384U1-1ML	1ML	3,6-Dichloro-2-hydroxy benzoic acid	3401-80-7	MET-11656B-50MG	50MG
Demeton O	298-03-3	N-11580-100MG	100MG	2,3-Dichloro-1,4-naphthoquinone	117-80-6	N-11661-250MG	250MG
Demeton O solution	298-03-3	S-11580U1-1ML	1ML	2,3-Dichloro-1,4-naphthoquinone Solution	117-80-6	S-11661M1-1ML	1ML
Demeton O&S	8065-48-3	N-11581-100MG	100MG	2,3-Dichloro-1,4-naphthoquinone Solution	117-80-6	S-11661U1-1ML	1ML
Demeton O&S solution	8065-48-3	S-11581U1-1ML	1ML	2,3-Dichloro-1,4-naphthoquinone Solution	117-80-6	S-11661U1-5ML	5ML
Demeton S	126-75-0	N-11582-100MG	100MG	2,3-Dichloronitrobenzene	3209-22-1	N-10584-250MG	250MG
Demeton S solution	126-75-0	S-11582A1-1ML	1ML	2,4-Dichloronitrobenzene	611-06-3	N-10625-250MG	250MG
Demeton S solution	126-75-0	S-11582U1-1ML	1ML	1,4-Dichloro-2-nitrobenzene	89-61-2	N-10670-250MG	250MG
Demeton-O-methyl	867-27-6	N-11583-50MG	50MG	2,5-Dichloronitrobenzene	89-61-2	N-10670-250MG	250MG
Denatonium benzoate	3734-33-6	N-11324-100MG	100MG	1,2-Dichloro-4-nitrobenzene	99-54-7	N-10785-250MG	250MG
Denatonium benzoate Solution	3734-33-6	S-11324F1-1ML	1ML	3,4-Dichloronitrobenzene	99-54-7	N-10785-250MG	250MG
Denatonium benzoate Solution	3734-33-6	S-11324U1-1ML	1ML	3,5-Dichloronitrobenzene	618-62-2	N-10792-250MG	250MG
Desmedipham	13684-56-5	N-11584-250MG	250MG	2,3-Dichloronitrobenzene solution	3209-22-1	S-10584A1-1ML	1ML
Desmedipham solution	13684-56-5	S-11584M1-1ML	1ML	2,3-Dichloronitrobenzene solution	3209-22-1	S-10584U1-1ML	1ML
N-Desmethylthiamethoxam	171103-04-1	MET-13562A-50MG	50MG	2,4-Dichloronitrobenzene solution	611-06-3	S-10625A1-1ML	1ML
N-Desmethylthiamethoxam Solution	171103-04-1	MET-13562AA1-1ML	1ML	2,4-Dichloronitrobenzene solution	611-06-3	S-10625U1-1ML	1ML
Desmetryn	1014-69-3	N-12973-10MG	10MG	2,5-Dichloronitrobenzene solution	89-61-2	S-10670A1-1ML	1ML
Dexon	140-56-7	N-11586-250MG	250MG	2,5-Dichloronitrobenzene solution	89-61-2	S-10670U1-1ML	1ML
Dexon Solution	140-56-7	S-11586M1-1ML	1ML	3,4-Dichloronitrobenzene solution	99-54-7	S-10785A1-1ML	1ML
Di-n-butyl phthalate solution	84-74-2	S-11589A1-1ML	1ML	3,4-Dichloronitrobenzene solution	99-54-7	S-10785U1-1ML	1ML
Di-n-butyl phthalate solution	84-74-2	S-11589J1-1ML	1ML	3,5-Dichloronitrobenzene solution	618-62-2	S-10792A1-1ML	1ML
Di-n-butyl phthalate solution	84-74-2	S-11589J1-5ML	5ML	3,5-Dichloronitrobenzene solution	618-62-2	S-10792U1-1ML	1ML
Diafenthion	80060-09-9	N-11612-250MG	250MG	1,4-Dichloro-2-nitrobenzene Solution	89-61-2	S-10670A1-1ML	1ML
Diafenthion solution	80060-09-9	S-11612A1-1ML	1ML	1,4-Dichloro-2-nitrobenzene Solution	89-61-2	S-10670U1-1ML	1ML
Diafenthion solution	80060-09-9	S-11612U1-1ML	1ML	1,2-Dichloro-4-nitrobenzene Solution	99-54-7	S-10785A1-1ML	1ML
Di-allyl	2303-16-4	N-11587-100MG	100MG	1,2-Dichloro-4-nitrobenzene Solution	99-54-7	S-10785U1-1ML	1ML
Di-allyl solution	2303-16-4	S-11587M1-1ML	1ML	1,1-Dichloro-1-nitroethane	594-72-9	N-10122-1G	1G
Di-allyl solution	2303-16-4	S-11587M1-5ML	5ML	1,1-Dichloro-1-nitroethane solution	594-72-9	S-10122M1-1ML	1ML
N,N-Diallyl-2-chloroacetamide	93-71-0	N-12617-250MG	250MG	1,1-Dichloro-1-nitroethane solution	594-72-9	S-10122X1-1ML	1ML
N,N-Diallyl-2-chloroacetamide solution	93-71-0	S-12617M1-1ML	1ML	Dichlorophen	97-23-4	N-11667-250MG	250MG
N,N-Diallyl-2,2-dichloroacetamide	37764-25-3	N-11662-100MG	100MG	Dichlorophen solution	97-23-4	S-11667M4-1ML	1ML
N,N-Diallyl-2,2-dichloroacetamide Solution	37764-25-3	S-11662A1-1ML	1ML	Dichlorophen solution	97-23-4	S-11667M4-5ML	5ML
N,N-Diallyl-2,2-dichloroacetamide Solution	37764-25-3	S-11662U1-1ML	1ML	2,3-Dichlorophenoxyacetic acid	2976-74-1	N-12953-10MG	10MG
Diammonium ethylenebis(Dithiocarbamate)	3566-10-700	N-11071-100MG	100MG	(2,4-Dichlorophenoxy)acetic acid - ethyl ester	533-23-3	N-10614-250MG	250MG
Diammonium ethylenebis(Dithiocarbamate) Sol.	3566-10-700	S-11071M1-1ML	1ML	(2,4-Dichlorophenoxy)acetic acid - ethyl ester Sol	533-23-3	S-10614A1-1ML	1ML
Diazil	57-68-1	N-13245-250MG	250MG	(2,4-Dichlorophenoxy)acetic acid - ethyl ester Sol	533-23-3	S-10614U1-1ML	1ML
Diazinon	333-41-5	N-11621-250MG	250MG	2,4-Dichlorophenoxy acetic acid-methyl ester Sol	1928-38-7	S-10619M1-1ML	1ML
Diazinon (diethyl-d10) Solution	333-41-5	S-FD2060S-1.2ML	1.2ML	2,4-Dichlorophenoxy acetic acid-methyl ester Sol	1928-38-7	S-10619M1-5ML	5ML
Diazinon solution	333-41-5	S-11621U1-1ML	1ML	N-(2,4-Dichlorophenoxy)acetyl-DL-methionine	62625-13-2	N-12505-10MG	10MG
Diazinon-O-analog	962-58-3	MET-11621A-100MG	100MG	N-(2,4-Dichlorophenoxy)acetyl-DL-methionine	62625-13-2	S-12505A1-1ML	1ML
Diazoxon	962-58-3	MET-11621A-100MG	100MG	N-(2,4-Dichlorophenoxy)acetyl-DL-methionine	62625-13-2	S-12505T1-1ML	1ML
Dibrom	128-04-1	N-12955-10MG	10MG	4-(2,4-Dichlorophenoxy)butyric acid methyl ester	18625-12-2	N-10533-100MG	100MG
Dibrom Solution	300-76-5	N-12840-250MG	250MG	4-(2,4-Dichlorophenoxy)butyric acid methyl ester	18625-12-2	S-10533M1-1ML	1ML
1,2-Dibromo-3-chloropropane	96-12-8	S-12640U1-1ML	1ML	4-(2,4-Dichlorophenoxy)butyric acid methyl ester	18625-12-2	S-10533M1-5ML	5ML
1,2-Dibromo-3-chloropropane solution	96-12-8	N-10149-250MG	250MG	2-(4-(2,4-Dichlorophenoxy)phenoxy)propanoic acid	40843-25-2	MET-11677B-100MG	100MG
1,2-Dibromo-3-chloropropane solution	96-12-8	S-10149M1-1ML	1ML	2,4-Dichlorophenyl benzenesulfonate	97-16-5	N-10628-250MG	250MG
1,2-Dibromo-3-chloropropane solution	96-12-8	S-10149M1-5ML	5ML	2,4-Dichlorophenyl benzenesulfonate solution	97-16-5	S-10628A1-1ML	1ML
2,6-Dibromo-4-cyanophenyl methyl ether	3336-39-8	MET-11345B-100MG	100MG	2,4-Dichlorophenyl benzenesulfonate solution	97-16-5	S-10628U1-1ML	1ML
2,6-Dibromo-4-cyanophenylheptanoate	56634-95-8	MET-11345A-50MG	50MG	3-(3,4-Dichlorophenyl)-1,1-dimethylurea Solution	330-54-1	S-11827A4-1ML	1ML
1,2-Dibromoethane	106-93-4	N-10150-1G	1G	3-(3,4-Dichlorophenyl)-1,1-dimethylurea Solution	330-54-1	S-11827U1-1ML	1ML
1,2-Dibromoethane solution	106-93-4	S-10150M1-1ML	1ML	3-(3,4-Dichlorophenyl)-1-methoxy-1-methylurea Sol	330-55-2	S-13232A6-1ML	1ML
Dibutyl succinate	141-03-7	N-11644-500MG	500MG	3-(3,4-Dichlorophenyl)-1-methoxy-1-methylurea Sol	330-55-2	S-13232A6-5ML	5ML
Dibutyl succinate solution	141-03-7	S-11644A1-1ML	1ML	1-(3,4-Dichlorophenyl)-3-methylurea	3567-62-2	MET-11827C-250MG	250MG
Dibutyl succinate solution	141-03-7	S-11644U1-1ML	1ML	N-(3,4-Dichlorophenyl)-N-methylurea	3567-62-2	MET-11827C-250MG	250MG
Dicamba	1918-00-9	N-11656-250MG	250MG	1-(2,4-Dichlorophenyl)-2-(1H-1,2,4-triazole-1-yl)-N,N-Bis(3,4-dichlorophenyl)urea	58905-16-1	MET-13576B-100MG	100MG
Dicamba solution	1918-00-9	S-11656A1-1ML	1ML	N,N-Bis(3,4-dichlorophenyl)urea	4300-43-0	MET-11827A-100MG	100MG
Dicamba solution	1918-00-9	S-11656B1-1ML	1ML	1,2-Dichloropropane	78-87-5	N-10156-1G	1G
Dicamba solution	1918-00-9	S-11656B1-5ML	5ML	1,2-Dichloropropane solution	78-87-5	S-10156M1-1ML	1ML
Dicamba methyl ester	6597-78-0	N-11657-100MG	100MG	1,2-Dichloropropane solution	78-87-5	S-10156M1-5ML	5ML
Dicamba methyl ester solution	6597-78-0	S-11657M1-1ML	1ML	1,3-Dichloropropane	542-75-6	N-10193-500MG	500MG
Dicamba methyl ester solution	6597-78-0	S-11657M1-5ML	5ML	1,3-Dichloropropane solution	542-75-6	S-10193M1-1ML	1ML
Dicapthion	2463-84-5	N-11658-100MG	100MG	3,7-Dichloroquinoline-8-carboxylic acid Solution	84087-01-4	S-13169A1-1ML	1ML
Dicapthion solution	2463-84-5	S-11658A1-1ML	1ML	3,7-Dichloroquinoline-8-carboxylic acid Solution	84087-01-4	S-13169T1-1ML	1ML
Dicapthion solution	2463-84-5	S-11658U1-1ML	1ML	3,6-Dichlorosalicylic acid	3401-80-7	MET-11656B-50MG	50MG
S-(1,2-Dicarboxyethyl)-O,S-dimethyl thiophosphate	3344-12-500	MET-12346AJ1-1ML	1ML	Dichloroprop	120-36-5	N-11671-250MG	250MG
Dichlobenil	1194-65-6	N-10688-250MG	250MG	Dichloroprop (ring-13C6) Solution	79270-78-3	S-FC962S-1.2ML	1.2ML
Dichlobenil Solution	1194-65-6	S-10688M1-1ML	1ML	Dichloroprop 2-ethylhexyl ester	120-36-5	N-13821-250MG	250MG
Dichlofenthion	97-17-6	N-11659-250MG	250MG	Dichloroprop solution	120-36-5	S-11671B1-1ML	1ML
Dichlofenthion solution	97-17-6	S-11659A1-1ML	1ML	Dichloroprop solution	120-36-5	S-11671B1-5ML	5ML
Dichlofenthion solution	97-17-6	S-11659J4-1ML	1ML	Dichloroprop solution	120-36-5	S-11671M1-1ML	1ML
Dichlofenthion solution	97-17-6	S-11659J4-5ML	5ML	Dichloroprop-butoxyethyl ester	53404-31-2	N-11673-50MG	50MG
Dichlofluandil	1085-98-9	N-11660-250MG	250MG	Dichloroprop methyl ester	57153-17-0	N-11672-100MG	100MG
Dichlofluandil solution	1085-98-9	S-11660U1-1ML	1ML	Dichloroprop methyl ester solution	57153-17-0	S-11672M1-1ML	1ML
Dichlone	117-80-6	N-11661-250MG	250MG	Dichloroprop methyl ester solution	57153-17-0	S-11672M1-5ML	5ML
Dichlone solution	117-80-6	S-11661M1-1ML	1ML	Dichloroprop-P	15165-67-0	N-11674-250MG	250MG
Dichlone solution	117-80-6	S-11661U1-1ML	1ML	Dichloroprop-P solution	15165-67-0	S-11674A1-1ML	1ML
Dichlone solution	117-80-6	S-11661U1-5ML	5ML	Dichloroprop-P solution	15165-67-0	S-11674U1-1ML	1ML
Dichlorimid	37764-25-3	N-11662-100MG	100MG	Dichlorosol	62-73-7	N-11675-250MG	250MG
Dichlorimid solution	37764-25-3	S-11662A1-1ML	1ML	Dichlorosol solution	62-73-7	S-11675J4-1ML	1ML
Dichlorimid solution	37764-25-3	S-11662U1-1ML	1ML	Dichlorosol solution	62-73-7	S-11675U1-1ML	1ML
2,4-Dichloroaniline	554-00-7	MET-11521A-1G	1G	Dichlorosol (dimethyl-d6)	203645-53-8	N-FD2061-E-0.01G	0.01G
3,5-Dichloroaniline	626-43-7	MET-12220B-250MG	250MG	Diclofop acid	40843-25-2	MET-11677B-100MG	100MG
2,6-Dichlorobenzamide	2008-58-4	MET-10688B-1G	1G	Diclofop methyl	51338-27-3	N-11677-250MG	250MG
1,2-Dichlorobenzene	95-50-1	N-10152-1G	1G	Diclofop methyl solution	51338-27-3	S-11677A1-1ML	1ML
1,4-Dichlorobenzene	106-46-7	N-10216-1G	1G	Diclofop methyl solution	51338-27-3	S-11677T1-1ML	1ML
1,2-Dichlorobenzene solution	95-50-1	S-10152M1-1ML	1ML	Dicloran	99-30-9	N-11678-250MG	250MG
1,2-Dichlorobenzene solution	95-50-1	S-10152M1-5ML	5ML	Dicloran solution	99-30-9	S-11678A1-1ML	1ML
1,4-Dichlorobenzene solution	106-46-7	S-10216M1-1ML	1ML	Dicloran solution	99-30-9	S-11678U4-1ML	1ML
1,4-Dichlorobenzene solution	106-46-7	S-10216M1-5ML	5ML	Diclosulam	145701-21-9	N-11679-250MG	250MG
2,4-Dichlorobenzenamine	554-00-7	MET-11521A-1G	1G	Diclosulam solution	14570		

Description	CAS#	Part #	Size	Description	CAS#	Part #	Size
Dicrotaphos	141-66-2	N-11680-100MG	100MG	2,2-Dimethyl-1,3-benzodioxole-4-ol		MET-11135A-500MG	500MG
Dicrotaphos solution	141-66-2	S-11680U1-1ML	1ML	Dimethyl chlorothiophosphate	2524-03-0	N-12704-250MG	250MG
Dicrotaphos solution	141-66-2	S-11680U1-5ML	5ML	Dimethyl chlorothiophosphate Solution	2524-03-0	S-12704A1-1ML	1ML
Dicryl	2164-09-200	N-11421-250MG	250MG	Dimethyl chlorothiophosphate Solution	2524-03-0	S-12704U1-1ML	1ML
Dicuran*	15545-48-9	N-11455-250MG	250MG	O,O-dimethyl-O-[4-nitro-m-tolyl]phosphate	2255-17-6	MET-11955AJ1-1ML	1ML
Dicuran* Solution	15545-48-9	S-11455A1-1ML	1ML	Dimethyl-p-nitrophenylphosphate	950-35-6	N-11775-100MG	100MG
Dicuran* Solution	15545-48-9	S-11455U1-1ML	1ML	Dimethyl-p-nitrophenylphosphate solution	950-35-6	S-11775A1-1ML	1ML
Dicyclopentadiene	77-73-6	N-11686-500MG	500MG	Dimethyl-p-nitrophenylphosphate solution	950-35-6	S-11775T1-1ML	1ML
Dicyclopentadiene solution	77-73-6	S-11686K1-1ML	1ML	Dimethyl-p-nitrophenylphosphate solution	950-35-6	S-11775T1-5ML	5ML
Dieldrin	60-57-1	N-11688-100MG	100MG	N-(2,4-Dimethylphenyl)formamide	60397-77-5	MET-11068A-50MG	50MG
Dieldrin (13C12, 99%) Solution		S-FC90AS-1.2ML	1.2ML	2-(2,4-Dimethylphenylimino)-3-methylthiazoline	61676-87-7	S-14051A1-1ML	1ML
Dieldrin solution	60-57-1	S-11688M1-1ML	1ML	1,1-Dimethyl-3-phenylurea trichloroacetate	4482-55-7	N-11968-50MG	50MG
Dieldrin solution	60-57-1	S-11688M1-5ML	5ML	O,O-Dimethyl phosphorochloridothioate	2524-03-0	N-12704-250MG	250MG
Dienochlor	2227-17-0	N-12824-250MG	250MG	O,O-Dimethyl phosphorochloridothioate solution	2524-03-0	S-12704A1-1ML	1ML
Dienochlor	2227-17-0	N-12824-250MG	250MG	O,O-Dimethyl phosphorochloridothioate solution	2524-03-0	S-12704U1-1ML	1ML
Dienochlor Solution	2227-17-0	S-12824M1-1ML	1ML	Dimethyl phthalate	131-11-3	N-11770-500MG	500MG
Dienochlor solution	2227-17-0	S-12824M1-1ML	1ML	Dimethyl phthalate solution	131-11-3	S-11770H1-1ML	1ML
Diethyl ethyl	38727-55-8	N-11690-100MG	100MG	Dimethyl phthalate solution	131-11-3	S-11770H1-5ML	5ML
Diethyl ethyl solution	38727-55-8	S-11690M1-1ML	1ML	Dimethyl phthalate solution	131-11-3	S-11770M1-1ML	1ML
Diethofencarb	87130-20-9	N-11691-250MG	250MG	Dimethyl phthalate solution	131-11-3	S-11770M1-5ML	5ML
Diethofencarb solution	87130-20-9	S-11691M1-1ML	1ML	N-(4,6-Dimethylpyrimidin-2-yl)aniline Solution	53112-28-0	S-13159A1-1ML	1ML
2-Diethylamino-6-methyl-4-hydroxy pyrimidine	42487-72-9	MET-13064A-100MG	100MG	N-(4,6-Dimethylpyrimidin-2-yl)aniline Solution	53112-28-0	S-13159U1-1ML	1ML
2-Diethylamino-6-methyl pyrimidin-4-ol	42487-72-9	MET-13064A-100MG	100MG	N,N-Dimethylsuccinamic acid	2564-95-6	N-12630-250MG	250MG
2,6-Diethylaniline	579-66-8	MET-11043C-500MG	500MG	N,N-Dimethylsuccinamic acid solution	2564-95-6	S-12630A1-1ML	1ML
Diethyl dithiobis(thionoformate)	502-55-6	N-11696-250MG	250MG	N,N-Dimethylsuccinamic acid solution	2564-95-6	S-12630T1-1ML	1ML
Diethyl dithiobis(thionoformate) solution	502-55-6	S-11696A1-1ML	1ML	Dimethylvinphos solution	2274-67-1	S-14052U1-1ML	1ML
Diethyl dithiobis(thionoformate) solution	502-55-6	S-11696K1-1ML	1ML	Dimoxystrobin solution	149961-52-4	S-13926A1-1ML	1ML
Diethyldithiocarbamic acid 2-Chloroallyl ester	95-06-7	N-13243-250MG	250MG	Diniconazole	83657-24-3	N-11782-100MG	100MG
Diethyldithiocarbamic acid 2-Chloroallyl ester	95-06-7	S-13243M1-1ML	1ML	Diniconazole solution	83657-24-3	S-11782M1-1ML	1ML
Diethyldithiocarbamic acid 2-Chloroallyl ester	95-06-7	S-13243M1-5ML	5ML	Dinitramine	29091-05-2	N-12976-10MG	10MG
O,O-Diethyl-S-(ethylthio)methyl]phosphorothioate	2600-69-3	MET-13035AU1-1ML	1ML	4,6-Dinitro-2-sec-butylphenol ethanolamine salt	87674-68-8	N-10888-1G	1G
N,N-Diethyl-4-methoxy-6-methyl-2-Pyrimidinamine	111697-06-4	MET-13064E-100MG	100MG	4,6-Dinitro-2-sec-butylphenol ethanolamine salt	87674-68-8	S-10888F1-1ML	1ML
O,O-Diethyl-O-P-nitrophenyl phosphorothioate	2164-09-200	N-FD98E-E-0.01G	0.01G	4,6-Dinitro-2-sec-butylphenol ethanolamine salt	87674-68-8	S-10888U1-1ML	1ML
O,O-Diethyl-O-2-parazinyl phosphate	7359-55-9	MET-13760A-50MG	50MG	4,6-Dinitro-o-cresol (contains ~10% water)	534-52-1	N-10889-250MG	250MG
Diethyl phosphate	598-02-2	MET-11621C-100MG	100MG	4,6-Dinitro-o-cresol solution	534-52-1	S-10889M4-1ML	1ML
O,O-Diethyl phosphorochloridothioate	2524-04-100	N-12703-250MG	250MG	2-sec-Butyl-4,6-dinitrophenol triethanolamine salt	6420-47-9	N-10517-500MG	500MG
O,O-Diethyl phosphorochloridothioate solution	2524-04-100	S-12703A1-1ML	1ML	2-sec-Butyl-4,6-dinitrophenol triethanolamine salt	6420-47-9	S-10517F1-1ML	1ML

Description	CAS#	Part #	Size	Description	CAS#	Part #	Size
DNOC	534-52-1	N-10889-250MG	250MG	Endrin solution	72-20-8	S-11854M1-5ML	5ML
DNOC Solution	534-52-1	S-10889M4-1ML	1ML	Endrin aldehyde	7421-93-4	N-11855-10MG	10MG
n-Dodecylguanidin acetate	2439-10-300	N-11840-250MG	250MG	Endrin aldehyde solution	7421-93-4	S-11855M1-1ML	1ML
Dodermorph	1593-77-7	N-11839-250MG	250MG	Endrin aldehyde solution	7421-93-4	S-11855M1-5ML	5ML
Dodermorph solution	1593-77-7	S-11839M1-1ML	1ML	Endrin ketone	53494-70-5	N-11856-10MG	10MG
Dodermorph acetate	31717-87-0	N-12948-100MG	100MG	Endrin ketone solution	53494-70-5	S-11856A1-1ML	1ML
Dodine (Technical)	2439-10-300	N-11840-250MG	250MG	Endrin ketone solution	53494-70-5	S-11856U1-1ML	1ML
Dotan(TM)	24934-91-6	N-11434-250MG	250MG	Envoke(TM)	199119-58-9	N-13685-100MG	100MG
Dotan(TM) Solution	24934-91-6	S-11434A1-1ML	1ML	EPN	2104-64-5	N-11858-250MG	250MG
Dotan(TM) Solution	24934-91-6	S-11434U1-1ML	1ML	EPN solution	2104-64-5	S-11858A1-1ML	1ML
Dowcide 1	90-43-7	N-10280-250MG	250MG	EPN solution	2104-64-5	S-11858J1-1ML	1ML
Dowcide 1 Solution	90-43-7	S-10280A1-1ML	1ML	EPN solution	2104-64-5	S-11858J1-5ML	5ML
Dowcide 1 Solution	90-43-7	S-10280M5-1ML	1ML	EPN solution	2104-64-5	S-11858J4-1ML	1ML
Dowco 217	5598-52-7	MET-11460B-50MG	50MG	EPN solution	2104-64-5	S-11858J4-5ML	5ML
Dowcide 2S(TM)	88-06-2	N-10657-250MG	250MG	Epoxidan	133220-30-1	N-12209-100MG	100MG
Dowcide 2S(TM) Solution	88-06-2	S-10657M1-1ML	1ML	Epoxiconazole	133855-98-8	N-11859-100MG	100MG
Dowcide 2S(TM) Solution	88-06-2	S-10657M1-5ML	5ML	EPTC	759-94-4	N-13188-500MG	500MG
Dowcide 6(TM)	58-90-2	N-10599-50MG	50MG	EPTC Solution	759-94-4	S-13188A1-1ML	1ML
Dowcide 6(TM) Solution	58-90-2	S-10599M1-1ML	1ML	EPTC Solution	759-94-4	S-13188T1-1ML	1ML
Dowpon S(TM)	127-20-8	N-11564-250MG	250MG	EPTC Solution	759-94-4	S-13188T1-5ML	5ML
2,4-DP methyl ester	57153-17-0	N-11672-100MG	100MG	Esbiothrin	84030-86-4	N-11864-100MG	100MG
2,4-DP methyl ester Solution	57153-17-0	S-11672M1-1ML	1ML	Escuri(TM)	83055-99-6	N-11139-100MG	100MG
2,4-DP methyl ester Solution	57153-17-0	S-11672M1-5ML	5ML	Escuri(TM) Solution	83055-99-6	S-11139A1-1ML	1ML
Drop(m)	51707-55-2	N-13564-250MG	250MG	Esfenvalerate	66230-04-4	N-11102-100MG	100MG
Drop(m) Solution	51707-55-2	S-13564A1-1ML	1ML	Esfenvalerate solution	66230-04-4	S-11102A1-1ML	1ML
Drop(m) Solution	51707-55-2	S-13564U1-1ML	1ML	Esfenvalerate solution	66230-04-4	S-11102T1-1ML	1ML
DSMA Hexahydrate	144-21-8	N-11817-500MG	500MG	Esteron 44 Solution	94-11-1	S-10617A1-1ML	1ML
DSMA Hexahydrate Solution	144-21-8	S-11817F1-1ML	1ML	Esteron 44 Solution	94-11-1	S-10617U1-1ML	1ML
Duplosan(TM)	15165-67-0	N-11674-250MG	250MG	Esteron 45	94-11-1	N-10617-250MG	250MG
Duplosan(TM) Solution	15165-67-0	S-11674A1-1ML	1ML	Etaconazole	60207-93-4	N-12985-100MG	100MG
Duplosan(TM) Solution	15165-67-0	S-11674U1-1ML	1ML	Ethalfuralin	55283-68-6	N-11865-250MG	250MG
cis-Duraphos(TM)	26718-65-0	N-11481-100MG	100MG	Ethalfuralin solution	55283-68-6	S-11865A1-1ML	1ML
cis-Duraphos(TM) Solution	26718-65-0	S-11481A1-1ML	1ML	Ethalfuralin solution	55283-68-6	S-11865J4-1ML	1ML
cis-Duraphos(TM) Solution	26718-65-0	S-11481T1-1ML	1ML	Ethalfuralin solution	55283-68-6	S-11865J4-5ML	5ML
Duraset (TM)	85-72-3	N-12557-100MG	100MG	Ethametsulfuron-methyl	97780-06-8	N-11866-100MG	100MG
Duraset (TM) Solution	85-72-3	S-12557A1-1ML	1ML	Ethametsulfuron-methyl solution	97780-06-8	S-11866A1-1ML	1ML
Duraset (TM) Solution	85-72-3	S-12557T1-1ML	1ML	Ethanedial dioxime	557-30-2	N-11867-100MG	100MG
Dyfonate	944-22-9	N-11842-100MG	100MG	Ethanedial dioxime	557-30-2	N-11867-250MG	250MG
Dyfonate Solution	944-22-9	S-11842A1-1ML	1ML	Ethanedial dioxime solution	557-30-2	S-11867M1-1ML	1ML
Dyfonate Solution	944-22-9	S-11842J4-1ML	1ML	Ethephon	16672-87-0	N-10002-250MG	250MG
Dyfonate Solution	944-22-9	S-11842J4-5ML	5ML	Ethide (TM)	594-72-9	N-10122-1G	1G
Dylox (TM)	52-68-6	N-11843-250MG	250MG	Ethide (TM) Solution	594-72-9	S-10122M1-1ML	1ML
Dymron Solution	42609-52-9	S-13955A1-1ML	1ML	Ethide (TM) Solution	594-72-9	S-10122X1-1ML	1ML
Ecothal	2303-17-5	N-13628-250MG	250MG	Ethidimuron	30043-49-3	N-11871-100MG	100MG
Ecothal Solution	2303-17-5	S-13628A1-1ML	1ML	Ethidimuron solution	30043-49-3	S-11871A1-1ML	1ML
Ecothal Solution	2303-17-5	S-13628K1-1ML	1ML	Ethidimuron solution	30043-49-3	S-11871U1-1ML	1ML
Edifenphos	17109-49-8	N-11844-100MG	100MG	Ethidimuron solution	30043-49-3	N-13913-10MG	10MG
Edifenphos solution	17109-49-8	S-11844A1-1ML	1ML	1-Ethynyl-2-methyl-2-penten-1-yl-chrysanthemate	54406-48-3	N-11872-100MG	100MG
Edifenphos solution	17109-49-8	S-11844U1-1ML	1ML	Ethiofencarb	29973-13-5	S-11872-100MG	100MG
EGT	2514-53-6	N-11845-100MG	100MG	Ethiofencarb solution	29973-13-5	S-11872A1-1ML	1ML
Elevate(TM)	126833-17-8	N-11954-250MG	250MG	Ethiofencarb solution	29973-13-5	S-11872U1-1ML	1ML
Elevate(TM) Solution	126833-17-8	S-11954A1-1ML	1ML	Ethiofencarb sulfone	53380-23-7	MET-11872B-10MG	10MG
Elgetol	6420-47-9	N-10517-500MG	500MG	Ethiofencarb sulfone Solution	53380-23-7	MET-11872BA1-1ML	1ML
Elgetol Solution	6420-47-9	S-10517F1-1ML	1ML	Ethiofencarb sulfoxide	53380-22-6	MET-11872A-10MG	10MG
Elgetol Solution	6420-47-9	S-10517U1-1ML	1ML	Ethion	563-12-2	N-11873-250MG	250MG
Emamectin	119791-41-2	N-11847-100MG	100MG	Ethion solution	563-12-2	S-11873M1-1ML	1ML
Emamectin solution	119791-41-2	S-11847A1-1ML	1ML	Ethion solution	563-12-2	S-11873M1-5ML	5ML
Emamectin benzoate	155569-91-8	N-12920-100MG	100MG	Ethiozin	64529-56-2	N-13731-100MG	100MG
Embutox	94-82-6	N-10622-250MG	250MG	Ethiozin solution	64529-56-2	S-13731M1-1ML	1ML
Embutox Solution	94-82-6	S-10622A1-1ML	1ML	Ethiprole	181587-01-9	N-12870-100MG	100MG
Embutox Solution	94-82-6	S-10622B1-1ML	1ML	Ethirimol	23947-60-6	N-12986-10MG	10MG
Embutox Solution	94-82-6	S-10622B1-5ML	5ML	Ethofumesate	26225-79-6	N-11874-250MG	250MG
Emmi (Technical)	2597-93-5	N-10224-100MG	100MG	Ethofumesate solution	26225-79-6	S-11874A1-1ML	1ML
Empenthrin	54406-48-3	N-13913-10MG	10MG	Ethofumesate solution	26225-79-6	S-11874U1-1ML	1ML
a-Endosulfan (13C9) Solution		S-FC202S-1.2ML	1.2ML	Ethohexadiol (TM)	94-96-2	N-10346-500MG	500MG
b-Endosulfan (13C9) Solution		S-FC203S-1.2ML	1.2ML	Ethohexadiol (TM) Solution	94-96-2	S-10346B6-1ML	1ML
a-Endosulfan (d4) Solution		S-FD202S-1.2ML	1.2ML	Ethohexadiol (TM) Solution	94-96-2	S-10346B6-5ML	5ML
b-Endosulfan (d4) Solution		S-FD203S-1.2ML	1.2ML	Ethoprophos (TM)	13194-48-4	N-13110-250MG	250MG
a-Endosulfan	959-98-8	N-10979-100MG	100MG	Ethoprophos (TM) Solution	13194-48-4	S-13110A1-1ML	1ML
a-Endosulfan (d4)		N-FD202-E-0.01G	0.01G	Ethoprophos (TM) Solution	13194-48-4	S-13110U1-1ML	1ML
a-Endosulfan solution	959-98-8	S-10979A1-1ML	1ML	Ethoprophos (TM) Solution	13194-48-4	S-13110U1-5ML	5ML
a-Endosulfan solution	959-98-8	S-10979J1-1ML	1ML	Ethoxyquin	91-53-2	N-11877-100MG	100MG
a-Endosulfan solution	959-98-8	S-10979J1-5ML	5ML	Ethoxysulfuron	126801-58-9	N-11878-100MG	100MG
a-Endosulfan solution	959-98-8	S-10979T1-1ML	1ML	2-Ethylamino-6-methyl-4(3H)-Pyrimidinone	5734-69-0	MET-13064C-100MG	100MG
a-Endosulfan solution	959-98-8	S-10979T1-5ML	5ML	Ethyl-3-amino-2,5-dichlorobenzoate	60541-86-8	N-11916-100MG	100MG
a,b-Endosulfan	115-29-7	N-10992-500MG	500MG	Ethyl-3-amino-2,5-dichlorobenzoate solution	60541-86-8	S-11916A1-1ML	1ML
a,b-Endosulfan solution	115-29-7	S-10992A1-1ML	1ML	Ethyl-3-amino-2,5-dichlorobenzoate solution	60541-86-8	S-11916T1-1ML	1ML
a,b-Endosulfan solution	115-29-7	S-10992J1-1ML	1ML	Ethylan	72-56-0	N-12849-250MG	250MG
a,b-Endosulfan solution	115-29-7	S-10992J1-5ML	5ML	Ethylan solution	72-56-0	S-12849A1-1ML	1ML
Endosulfan alcohol	2157-19-9	MET-11848A-50MG	50MG	Ethylan solution	72-56-0	S-12849K4-1ML	1ML
b-Endosulfan	33213-65-9	N-11118-100MG	100MG	Ethylan solution	72-56-0	S-12849K4-5ML	5ML
b-Endosulfan (d4)		N-FD203-5-5MG	5MG	2-Ethyl-2-butyl-1,3-propanediol	115-84-4	N-10306-1G	1G
b-Endosulfan solution	33213-65-9	S-11118A1-1ML	1ML	2-Ethyl-2-butyl-1,3-propanediol Solution	115-84-4	S-10306M1-1ML	1ML
b-Endosulfan solution	33213-65-9	S-11118J1-1ML	1ML	O-Ethyl-S,S-diphenyl phosphorodithioate	17109-49-8	N-11844-100MG	100MG
b-Endosulfan solution	33213-65-9	S-11118J1-5ML	5ML	O-Ethyl-S,S-diphenyl phosphorodithioate Solution	17109-49-8	S-11844A1-1ML	1ML
beta-Endosulfan Solution	33213-65-9	S-11118A1-1ML	1ML	O-Ethyl-S,S-diphenyl phosphorodithioate Solution	17109-49-8	S-11844U1-1ML	1ML
beta-Endosulfan Solution	33213-65-9	S-11118J1-1ML	1ML	S-Ethyl dipropylthiocarbamate	759-94-4	N-13188-500MG	500MG
beta-Endosulfan Solution	33213-65-9	S-11118J1-5ML	5ML	S-Ethyl dipropylthiocarbamate solution	759-94-4	S-13188A1-1ML	1ML
Endosulfan ether	3369-52-6	MET-11848B-50MG	50MG	S-Ethyl dipropylthiocarbamate solution	759-94-4	S-13188T1-1ML	1ML
Endosulfan lactone	3868-61-9	MET-11848C-50MG	50MG	S-Ethyl dipropylthiocarbamate solution	759-94-4	S-13188T1-5ML	5ML
Endosulfan sulfate	1031-07-8	N-11851-50MG	50MG	Ethylene	3424-82-6	MET-12706A-10MG	10MG
Endosulfan sulfate solution	1031-07-8	S-11851A1-1ML	1ML	Ethylene bis(dithiocarbamic acid)-manganese salt	12427-38-2	N-12355-250MG	250MG
Endosulfan sulfate solution	1031-07-8	S-11851J1-1ML	1ML	Ethylene bis(trichloroacetate)	2514-53-6	N-11845-100MG	100MG
Endosulfan sulfate solution	1031-07-8	S-11851J1-5ML	5ML	Ethylene dibromide	106-93-4	N-10150-1G	1G
Endothal (unstated stereochemistry)	62059-43-2	N-11852-100MG	100MG	Ethylene dibromide Solution	106-93-4	S-10150M1-1ML	1ML
Endothal (unstated stereochemistry) Solution	62059-43-2	S-11852A1-1ML	1ML	2-Ethyl-1,3-hexanediol	94-96-2	N-10346-500MG	500MG
Endothal (unstated stereochemistry) Solution	62059-43-2	S-11852A1-5ML	5ML	2-Ethyl-1,3-hexanediol solution	94-96-2	S-10346B6-1ML	1ML
Endothal (unstated stereochemistry) Solution	62059-43-2	S-11852F0-1ML	1ML	2-Ethyl-1,3-hexanediol solution	94-96-2	S-10346B6-5ML	5ML
Endothal (unstated stereochemistry) Solution	62059-43-2	S-11852F10-1ML	1ML	2-Ethylhexyl-2,4-dichlorophenoxy acetate	1928-43-4	N-10356-1G	1G
Endothal (unstated stereochemistry) Solution	62059-43-2	S-11852U1-1ML	1ML	2-Ethylhexyl-2,4-dichlorophenoxy acetate solution	1928-43-4	S-10356A1-1ML	1ML
Endothal monohydrate	62059-43-2	N-11852-100MG	100MG	2-Ethylhexyl-2,4-dichlorophenoxy acetate solution	1928-43-4	S-10356J1-1ML	1ML
Endothal monohydrate solution	62059-43-2	S-11852A1-1ML	1ML	Ethylmercuric chloride(Technical)	107-27-7	N-11940-100MG	100MG
Endothal monohydrate solution	62059-43-2	S-11852A1-5ML	5ML	2-[(2-Ethyl-6-methylphenyl)amino]-1-propanol	61520-53-4	MET-12478A-50MG	50MG
Endothal monohydrate solution	62059-43-2	S-11852F0-1ML	1ML	Ethyl-1-naphthalene acetate	2122-70-5	N-11914-250MG	250MG
Endothal monohydrate solution	62059-43-2	S-11852F10-1ML	1ML	Ethyl-1-naphthalene acetate solution	2122-70-5	S-11914A1-1ML	1ML
Endothal monohydrate solution	62059-43-2	S-11852F10-5ML	5ML	Ethyl-1-naphthalene acetate solution	2122-70-5	S-11914T1-1ML	1ML
Endothal monohydrate solution	62059-43-2	S-11852U1-1ML	1ML	4-(Ethylthio)-6-methyl-2-(1-methylethyl)-pyrimidine	77738-92-2	MET-11621D-10MG	10MG
Endothal monohydrate solution	62059-43-2	S-11852U1-5ML	5ML	N-Ethylcyclohexylamine	5459-93-8	MET-11522A-250MG	250MG
Endothal-Pentafluorophenyl hydrazine Solution		S-11853T1-1ML	1ML	Etofenprox	80844-07-1	N-11941-100MG	100MG
Endothal-Pentafluorophenyl hydrazine Solution		S-11853T1-5ML	5ML	Etofenprox solution	80844-07-1	S-11941A1-1ML	1ML
Endothal-PFPH solution		S-11853T1-1ML	1ML	Etofenprox solution	80844-07-1	S-11941U1-1ML	1ML
Endothal-PFPH solution		S-11853T1-5ML	5ML	Etofenprox-desethyl Solution	80854-21-3	MET-11941AM1-1ML	1ML
Endrin	72-20-8	N-11854-100MG	100MG	Etoazole	153233-91-1	N-11942-50MG	50MG
Endrin solution	72-20-8	S-11854M1-1ML	1ML	Etoazole solution	153233-91-1	S-11942A1-1ML	1ML
				Etoazole solution	153233-91-1	S-11942U1-1ML	1ML

Description	CAS#	Part #	Size	Description	CAS#	Part #	Size
Etridiazole	2593-15-9	N-13517-250MG	250MG	Fenuron solution	101-42-8	S-11967A4-5ML	5ML
Etridiazole	2593-15-9	N-13517-250MG	250MG	Fenuron solution	101-42-8	S-11967U1-1ML	1ML
Etridiazole Solution	2593-15-9	S-13517A1-1ML	1ML	Fenuron TCA	4482-55-7	N-11968-50MG	50MG
Etridiazole solution	2593-15-9	S-13517A1-1ML	1ML	Fenvalerate	51630-58-1	N-13201-250MG	250MG
Etridiazole Solution	2593-15-9	S-13517T1-1ML	1ML	Fenvalerate solution	51630-58-1	S-13201A4-1ML	1ML
Etridiazole solution	2593-15-9	S-13517T1-1ML	1ML	Fenvalerate solution	51630-58-1	S-13201A4-5ML	5ML
Etrifos	38260-54-7	N-12954-10MG	10MG	Fenvalerate solution	51630-58-1	S-13201U1-1ML	1ML
Eugenol	97-53-0	N-10803-250MG	250MG	Ferbam	14484-64-1	N-11970-250MG	250MG
Eugenol Solution	97-53-0	S-10803M1-1ML	1ML	Ferbam solution	14484-64-1	S-11970M1-1ML	1ML
Euparen M(TM)	731-27-1	N-13584-250MG	250MG	Ferbam solution	14484-64-1	S-11970U1-1ML	1ML
Evergol™	494793-67-8	N-14035-100MG	100MG	Ferimzone	89269-64-7	N-14030-100MG	100MG
Evergol™ Solution	494793-67-8	S-14035A1-1ML	1ML	Ferimzone solution	89269-64-7	S-14030A1-1ML	1ML
EXD	502-55-6	N-11696-250MG	250MG	(E)-Ferimzone solution	77359-18-3	S-14031A1-1ML	1ML
EXD Solution	502-55-6	S-11696A1-1ML	1ML	Fermate	14484-64-1	N-11970-250MG	250MG
EXD Solution	502-55-6	S-11696K1-1ML	1ML	Fermate Solution	14484-64-1	S-11970M1-1ML	1ML
Famoxadone	131807-57-3	N-11943-100MG	100MG	Fermate Solution	14484-64-1	S-11970U1-1ML	1ML
Famoxadone solution	131807-57-3	S-11943A1-1ML	1ML	Fipronil	120068-37-3	N-11971-100MG	100MG
Famoxadone solution	131807-57-3	S-11943U1-1ML	1ML	Fipronil solution	120068-37-3	S-11971A1-1ML	1ML
Famphur	52-85-7	N-11944-250MG	250MG	Fipronil solution	120068-37-3	S-11971U1-1ML	1ML
Famphur solution	52-85-7	S-11944K1-1ML	1ML	Fipronil desulfinyl solution	205650-65-3	S-13920B1-1ML	1ML
Famphur-O-analog	960-25-8	MET-11944A-50MG	50MG	Fipronil sulfide	120067-83-6	MET-11971B-25MG	25MG
Fenac	85-34-7	N-11945-250MG	250MG	Fipronil sulfide Solution	120067-83-6	MET-11971BA1-1ML	1ML
Fenac solution	85-34-7	S-11945A1-1ML	1ML	Fipronil sulfone	120068-36-2	MET-11971A-50MG	50MG
Fenac solution	85-34-7	S-11945T1-1ML	1ML	Fipronil sulfone Solution	120068-36-2	MET-11971AA1-1ML	1ML
Fenamidone	161326-34-7	N-11946-100MG	100MG	FirstRate	147150-35-4	N-11492-100MG	100MG
Fenamidone solution	161326-34-7	S-11946A1-1ML	1ML	FirstRate Solution	147150-35-4	S-11492A1-1ML	1ML
Fenamidone solution	161326-34-7	S-11946U1-1ML	1ML	FirstRate Solution	147150-35-4	S-11492U1-1ML	1ML
Fenamidone Metabolite	332855-88-6	MET-11946-100MG	100MG	Flamprop isopropyl	52756-22-6	N-11972-100MG	100MG
Fenaminosulf	140-56-7	N-11586-250MG	250MG	Flamprop-M-isopropyl	57973-67-8	N-13939-250MG	250MG
Fenaminosulf solution	140-56-7	S-11586M1-1ML	1ML	Flamprop isopropyl solution	52756-22-6	S-11972A1-1ML	1ML
Fenamiphos	22224-92-6	N-12854-250MG	250MG	Flamprop isopropyl solution	52756-22-6	S-11972T1-1ML	1ML
Fenamiphos solution	22224-92-6	S-12854A1-1ML	1ML	Flamprop-M-isopropyl solution	57973-67-8	S-13939M1-1ML	1ML
Fenamiphos solution	22224-92-6	S-12854T1-1ML	1ML	Flamprop-methyl	52756-25-9	N-11973-100MG	100MG
Fenamiphos solution	22224-92-6	S-12854T1-5ML	5ML	Flamprop-methyl solution	52756-25-9	S-11973A1-1ML	1ML
Fenamiphos sulfone	31972-44-8	MET-11947B-100MG	100MG	Flamprop-methyl solution	52756-25-9	S-11973T1-1ML	1ML
Fenamiphos sulfoxide	31972-43-7	MET-11947A-50MG	50MG	Flazasulfuron	104040-78-0	N-11974-50MG	50MG
Fenarimol	60168-88-9	N-11948-100MG	100MG	Flint(TM)	141517-21-7	N-13684-100MG	100MG
Fenarimol solution	60168-88-9	S-11948A1-1ML	1ML	Flint(TM) Solution	141517-21-7	S-13684A1-1ML	1ML
Fenarimol solution	60168-88-9	S-11948T1-1ML	1ML	Flocumafen	90035-08-8	N-11975-10MG	10MG
Fenarimol solution	60168-88-9	S-11948T1-5ML	5ML	Flocumafen solution	90035-08-8	S-11975A1-1ML	1ML
Fenazaquin	120928-09-8	N-11949-100MG	100MG	Flonicamid	158062-67-0	N-11976-50MG	50MG
Fenazaquin solution	120928-09-8	S-11949A1-1ML	1ML	Flonicamid solution	158062-67-0	S-11976A1-1ML	1ML
Fenbuconazole	114369-43-6	N-11950-100MG	100MG	Flonicamid solution	158062-67-0	S-11976U1-1ML	1ML
Fenbuconazole solution	114369-43-6	S-11950M1-1ML	1ML	Floramite(TM)	149877-41-8	N-1201-50MG	50MG
Fenbutatin-oxide(Technical)	13356-08-6	N-11951-250MG	250MG	Floramite(TM) Solution	149877-41-8	S-11201A1-1ML	1ML
Fenchlorim	3740-92-9	N-13265-250MG	250MG	Floramite(TM) Solution	149877-41-8	S-1201U1-1ML	1ML
Fenchlorphos	299-84-3	N-11952-100MG	100MG	Florasulam	145701-23-1	N-13922-10MG	10MG
Fenchlorphos solution	299-84-3	S-11952U1-1ML	1ML	Fluacrypyrim solution	229977-93-9	S-13980M1-1ML	1ML
Fenchlorphos-oxon	299-84-3	S-11952U1-5ML	5ML	Fluazaindolizine solution	1254304-22-7	S-14042A1-1ML	1ML
Fenethanil	3983-45-7	MET-11952A-100MG	100MG	Fluazifop	69335-91-7	N-11977-10MG	10MG
Fenethanil Solution	114369-43-6	N-11950-100MG	100MG	Fluazifop solution	69335-91-7	S-11977A1-1ML	1ML
Fenfuram	24691-80-3	S-11950M1-1ML	1ML	Fluazifop solution	69335-91-7	S-11977A1-5ML	5ML
Fenfuram solution	24691-80-3	N-14013-250MG	250MG	Fluazifop-butyl	69806-50-4	N-11978-250MG	250MG
Fenhexamid	126833-17-8	S-14013M1-1ML	1ML	Fluazifop-butyl solution	69806-50-4	S-11978A1-1ML	1ML
Fenhexamid solution	126833-17-8	N-11954-250MG	250MG	Fluazifop-butyl solution	69806-50-4	S-11978U1-1ML	1ML
Fenitrothion	122-14-5	S-11954A1-1ML	1ML	Fluazifop-methyl	69335-90-6	N-13906-50MG	50MG
Fenitrothion solution	122-14-5	N-11955-250MG	250MG	Fluazifop-P	83066-88-0	N-13907-50MG	50MG
Fenitrothion solution	122-14-5	S-11955A1-1ML	1ML	Fluazifop-P-butyl	79241-46-6	N-11979-100MG	100MG
Fenitrothion-O-analog Solution	122-14-5	S-11955J4-1ML	1ML	Fluazifop-P-butyl solution	79241-46-6	S-11979A1-1ML	1ML
Fenoprop (TM)	93-72-1	MET-11955A1-1ML	1ML	Fluazifop-P-butyl solution	79241-46-6	S-11979U1-1ML	1ML
Fenoprop (TM) Solution	93-72-1	N-13212-250MG	250MG	Fluazifop-P-methyl	69335-90-6	N-13908-50MG	50MG
Fenoprop (TM) Solution	93-72-1	S-13212A1-1ML	1ML	Fluazinam	79622-59-6	N-11980-100MG	100MG
Fenoprop (TM) Solution	93-72-1	S-13212B1-1ML	1ML	Fluazinam solution	79622-59-6	S-11980A1-1ML	1ML
Fenoprop (TM) Solution	93-72-1	S-13212B1-5ML	5ML	Fluazinam solution	79622-59-6	S-11980U1-1ML	1ML
Fenothiocarb solution	62850-32-2	S-13962A1-1ML	1ML	Fluazinam solution	79622-59-6	N-11981-100MG	100MG
Fenoxanil	115852-48-7	N-11956-100MG	100MG	Fluazinam solution	79622-59-6	N-11982-250MG	250MG
Fenoxaprop	95617-09-7	N-13830-50MG	50MG	Fluazinam solution	79622-59-6	S-11982F1-1ML	1ML
Fenoxaprop ethyl	66441-23-4	N-11957-100MG	100MG	Fluazinam solution	79622-59-6	N-11983-50MG	50MG
Fenoxaprop ethyl (racemate CAS # 82110-72-3)	66441-23-4	N-11957-100MG	100MG	Fluazinam solution	79622-59-6	S-11983A1-1ML	1ML
Fenoxaprop ethyl (racemate CAS # 82110-72-3)	66441-23-4	S-11957M1-1ML	1ML	Fluazinam solution	79622-59-6	S-11983J4-1ML	1ML
Fenoxaprop ethyl solution	66441-23-4	S-11957M1-1ML	1ML	Fluazinam solution	79622-59-6	S-11983J4-5ML	5ML
Fenoxaprop-P	113158-40-0	MET-13802-100MG	100MG	Fluazinam solution	79622-59-6	N-11984-100MG	100MG
Fenoxaprop-P-ethyl	71283-80-2	N-11958-100MG	100MG	Fluazinam solution	79622-59-6	S-11984A1-1ML	1ML
Fenoxaprop-P-ethyl solution	71283-80-2	S-11958A1-1ML	1ML	Fluazinam solution	79622-59-6	S-11984U1-1ML	1ML
Fenoxaprop-P-ethyl solution	71283-80-2	S-11958U1-1ML	1ML	Fluazinam solution	79622-59-6	N-11985-100MG	100MG
Fenoxycarb	72490-01-8	N-11959-100MG	100MG	Fluazinam solution	79622-59-6	S-11985A1-1ML	1ML
Fenoxycarb solution	72490-01-8	S-11959A1-1ML	1ML	Fluazinam solution	79622-59-6	S-11985U1-1ML	1ML
Fenpropathrin	39515-41-8	N-11960-250MG	250MG	Fluazinam solution	79622-59-6	N-14084-100MG	100MG
Fenpropathrin solution	39515-41-8	S-11960A1-1ML	1ML	Fluazinam solution	79622-59-6	N-11986-100MG	100MG
Fenpropathrin solution	39515-41-8	S-11960T1-1ML	1ML	Fluazinam solution	79622-59-6	S-11986M1-1ML	1ML
Fenproximate	134098-61-6	N-11962-100MG	100MG	Fluazinam solution	79622-59-6	MET-11986A-25MG	25MG
Fenproximate solution	134098-61-6	S-11962A1-1ML	1ML	Fluazinam solution	79622-59-6	MET-11986AM1-1ML	1ML
Fenproximate solution	134098-61-6	S-11962U1-1ML	1ML	Fluazinam solution	79622-59-6	MET-11986B-10MG	10MG
Fenson	80-38-6	N-13812-250MG	250MG	Fluazinam solution	79622-59-6	MET-11986BM1-1ML	1ML
Fensulfothion	115-90-2	N-11963-100MG	100MG	Fluazinam solution	79622-59-6	MET-11986B-10MG	10MG
Fensulfothion solution	115-90-2	S-11963M1-1ML	1ML	Fluazinam solution	79622-59-6	MET-11986BM1-1ML	1ML
Fensulfothion solution	115-90-2	S-11963M1-5ML	5ML	Fluazinam solution	79622-59-6	N-11987-50MG	50MG
Fensulfothion Oxon Solution	6552-21-2	MET-11963DA1-1ML	1ML	Fluazinam solution	79622-59-6	S-11987A1-1ML	1ML
Fensulfothion oxon sulfone Solution	6132-17-8	MET-11963CA1-1ML	1ML	Fluazinam solution	79622-59-6	N-13139-100MG	100MG
Fensulfothion-sulfide	3070-15-3	MET-11963B-100MG	100MG	Fluazinam solution	79622-59-6	S-13139A1-1ML	1ML
Fensulfothion-sulfone	14255-72-2	MET-11963A-50MG	50MG	Fluazinam solution	79622-59-6	N-11988-250MG	250MG
Fensulfoxide	3761-41-9	MET-11964BX1-1ML	1ML	Fluazinam solution	79622-59-6	S-11988A1-1ML	1ML
Fensulfoxide	3761-41-9	MET-11964D-100MG	100MG	Fluazinam solution	79622-59-6	S-11988B1-1ML	1ML
Fenthion	55-38-9	N-11964-250MG	250MG	Fluazinam solution	79622-59-6	S-11988B1-5ML	5ML
Fenthion solution	55-38-9	S-11964S1-1ML	1ML	Fluazinam solution	79622-59-6	S-11988K1-1ML	1ML
Fenthion solution	55-38-9	S-11964S1-5ML	5ML	Fluazinam solution	79622-59-6	N-11989-100MG	100MG
Fenthion-O-analog Solution	6552-12-100	MET-11964CK1-1ML	1ML	Fluazinam solution	79622-59-6	S-11989A1-1ML	1ML
Fenthion-ethyl	1716-09-2	N-11965-500MG	500MG	Fluazinam solution	79622-59-6	S-11989U1-1ML	1ML
Fenthion-ethyl solution	1716-09-2	S-11965A1-1ML	1ML	Fluazinam solution	79622-59-6	N-11990-500MG	500MG
Fenthion-ethyl solution	1716-09-2	S-11965U1-1ML	1ML	Fluazinam solution	79622-59-6	S-11990A1-1ML	1ML
Fenthion-oxon-sulfone	14086-35-2	MET-11964E-100MG	100MG	Fluazinam solution	79622-59-6	S-11990T1-1ML	1ML
Fenthion oxon sulfide Solution	6552-13-2	MET-11964FA1-1ML	1ML	Fluazinam solution	79622-59-6	N-11991-100MG	100MG
Fenthion-sulfone	3761-42-0	MET-11964A-100MG	100MG	Fluazinam solution	79622-59-6	S-11991A1-1ML	1ML
Fenthion-sulfide	3761-41-9	MET-11964D-100MG	100MG	Fluazinam solution	79622-59-6	S-11991U1-1ML	1ML
Fenthion sulfide Solution	3761-41-9	MET-11964BX1-1ML	1ML	Fluazinam solution	79622-59-6	N-13280-25MG	25MG
Fentin acetate	900-95-8	N-13007-250MG	250MG	Flumorph	211867-47-9	S-13280A1-1ML	1ML
Fentin acetate Solution	900-95-8	S-13007A1-1ML	1ML	Flumorph solution	211867-47-9	N-11991-100MG	100MG
Fentin acetate Solution	900-95-8	S-13007U1-1ML	1ML	Flumyzin(TM)	103361-09-7	S-11991A1-1ML	1ML
Fentin hydroxide	76-87-9	N-11966-250MG	250MG	Flumyzin(TM) Solution	103361-09-7	S-11991U1-1ML	1ML
Fentin hydroxide solution	76-87-9	S-11966A1-1ML	1ML	Flumeturon	2164-17-2	N-11992-250MG	250MG
Fentrazamide	158237-07-1	N-14059-100MG	100MG	Flumeturon solution	2164-17-2	S-11992A4-1ML	1ML
Fentrazamide solution	158237-07-1	N-14059A1-1ML	1ML	Fluopicolide	239110-15-7	N-11993-100MG	100MG
Fenuron	101-42-8	N-11967-250MG	250MG	Fluopyram	658066-35-4	N-13988-50MG	50MG
Fenuron solution	101-42-8	S-11967A4-1ML	1ML	Fluopyram solution	658066-35-4	S-13988A1-1ML	1ML

Description	CAS#	Part #	Size	Description	CAS#	Part #	Size
2-Fluoroacetamide	640-19-7	N-10358-100MG	100MG	Genite (TM) Solution	97-16-5	S-10628U1-1ML	1ML
2-Fluoroacetamide solution	640-19-7	S-10358M1-1ML	1ML	Gestamin	1610-17-9	N-11105-250MG	250MG
Fluoroacetic acid amide	640-19-7	N-10358-100MG	100MG	Gestamin Solution	1610-17-9	S-11105A1-1ML	1ML
Fluoroacetic acid amide Solution	640-19-7	S-10358M1-1ML	1ML	Gestamin Solution	1610-17-9	S-11105T1-1ML	1ML
Fluoroacetic acid-sodium salt	34014-18-1	N-13505-250MG	250MG	Gestamin Solution	1610-17-9	S-11105T1-5ML	5ML
Fluorodifen	15457-05-3	N-11998-250MG	250MG	Gibberellic acid	77-06-5	N-12108-100MG	100MG
Fluorodifen solution	15457-05-3	S-11998M1-1ML	1ML	Gibberellic acid solution	77-06-5	S-12108A1-1ML	1ML
Fluoroglycofen-ethyl	77501-90-7	N-11999-100MG	100MG	Gibberellic acid solution	77-06-5	S-12108T1-1ML	1ML
Fluoxastrobil	361377-29-9	N-13816-100MG	100MG	Gibberellin X	77-06-5	N-12108-100MG	100MG
Flupyradifurone	951659-40-8	N-13941-100MG	100MG	Gibberellin X Solution	77-06-5	S-12108A1-1ML	1ML
Flupyradifurone solution	951659-40-8	S-13941A1-1ML	1ML	Gibberellin X Solution	77-06-5	S-12108T1-1ML	1ML
Flurenol	467-69-6	N-13815-250MG	250MG	Glean	64902-72-3	N-11461-100MG	100MG
Flurenol-butyl	2314-09-200	N-13814-250MG	250MG	Glufosinate-ammonium	77182-82-2	N-12111-100MG	100MG
Flurenol-methyl	1216-44-0	N-13813-250MG	250MG	Glycel(TM)	38641-94-0	N-12134-250MG	250MG
Flurenol-methyl	1216-44-0	N-13813-250MG	250MG	Glycel(TM) Solution	38641-94-0	S-12134F1-1ML	1ML
Fluridone	59756-60-4	N-13217-100MG	100MG	Glycel(TM) Solution	38641-94-0	S-12134U1-1ML	1ML
Fluridone solution	59756-60-4	S-13217T1-1ML	1ML	Glyodin(Technical)	556-22-9	N-12131-250MG	250MG
Fluridone solution	59756-60-4	S-13217T1-5ML	5ML	Glyoxime	557-30-2	N-11867-250MG	250MG
Flurochloridone	61213-25-0	N-12987-100MG	100MG	Glyoxime Solution	557-30-2	S-11867M1-1ML	1ML
Flurochloridone solution	61213-25-0	S-12987M1-1ML	1ML	Glyphosate	1071-83-6	N-12133-250MG	250MG
Fluroxypyr	69377-81-7	N-12000-100MG	100MG	Glyphosate (2-13CP, 15N) Solution		S-FCN104S-1.2ML	1.2ML
Fluroxypyr solution	69377-81-7	S-12000M1-1ML	1ML	Glyphosate-isopropylammonium	38641-94-0	N-12134-250MG	250MG
Fluroxypyr-meptyl	81406-37-3	N-12001-100MG	100MG	Glyphosate-isopropylammonium solution	38641-94-0	S-12134F1-1ML	1ML
Fluroxypyr-meptyl solution	81406-37-3	S-12001M1-1ML	1ML	Glyphosate-isopropylammonium solution	38641-94-0	S-12134U1-1ML	1ML
Flurprimidol	56425-91-3	N-12002-100MG	100MG	Glyphosine	2439-99-8	N-12135-100MG	100MG
Flurprimidol solution	56425-91-3	S-12002M1-1ML	1ML	Glyphosine solution	2439-99-8	S-12135F1-1ML	1ML
Flurtamone solution	96525-23-4	S-13981A1-1ML	1ML	Glyphosine solution	2439-99-8	S-12135U1-1ML	1ML
Flusilazole	85509-19-9	N-12003-100MG	100MG	Gokilaht(TM)	39515-40-7	N-11561-100MG	100MG
Flusilazole solution	85509-19-9	S-12003A1-1ML	1ML	Gokilaht(TM) Solution	39515-40-7	S-11561A1-1ML	1ML
Flusilazole solution	85509-19-9	S-12003U1-1ML	1ML	Gokilaht(TM) Solution	39515-40-7	S-11561T1-1ML	1ML
Flusulfamide	106917-52-6	N-14029-100MG	100MG	Grasp(TM)	87820-88-0	N-13587-500MG	500MG
Flusulfamide solution	106917-52-6	S-14029A1-1ML	1ML	Grasp(TM) Solution	87820-88-0	S-13587A1-1ML	1ML
Fluthiacet-methyl	117337-19-6	N-13129-25MG	25MG	Grasp(TM) Solution	87820-88-0	S-13587U1-1ML	1ML
Fluthiacet-methyl solution	117337-19-6	S-13129A1-1ML	1ML	Gratil	120923-37-7	N-13819-100MG	100MG
Fluthiamide	142459-58-3	N-11986-100MG	100MG	N-Guanyl-N-methyl glycine	57-00-1	MET-11840A-250MG	250MG
Fluthiamide Solution	142459-58-3	S-11986M1-1ML	1ML	Guthion	86-50-0	N-12139-250MG	250MG
Flutolanil	66332-96-5	N-12004-250MG	250MG	Guthion Solution	86-50-0	S-12139A1-1ML	1ML
Flutolanil solution	66332-96-5	S-12004M1-1ML	1ML	Guthion Solution	86-50-0	S-12139U1-1ML	1ML
Flutriafol	76674-21-0	N-12005-100MG	100MG	Guthion Solution	86-50-0	S-12139U1-5ML	5ML
Flutriafol solution	76674-21-0	S-12005A1-1ML	1ML	Guthion ethyl	2642-71-9	N-12140-100MG	100MG
Flutriafol solution	76674-21-0	S-12005U1-1ML	1ML	Guthion ethyl Solution	2642-71-9	S-12140A1-1ML	1ML
tau-Fluvalinate	102851-06-9	N-13263-100MG	100MG	Guthion Ethyl solution	2642-71-9	S-12140A1-1ML	1ML
tau-Fluvalinate solution	102851-06-9	S-13263A1-1ML	1ML	Gypsine	7784-40-9	N-12310-100MG	100MG
tau-Fluvalinate solution	102851-06-9	S-13263K1-1ML	1ML	Halfenprox	111872-58-3	N-13079-10MG	10MG
Fluxaproxad	907204-31-3	N-12921-100MG	100MG	Halfenprox solution	111872-58-3	S-13079M1-1ML	1ML
Fluxofenim	88485-37-4	N-13940-250MG	250MG	Halofenozide	112226-61-6	N-12141-500MG	500MG
Fluxofenim solution	88485-37-4	S-13940M1-1ML	1ML	Halofenozide solution	112226-61-6	S-12141M1-1ML	1ML
Folex (TM)	150-50-5	N-13648-250MG	250MG	Halosulfuron-methyl	100784-20-1	N-12142-100MG	100MG
Folicur(TM)	107534-96-3	N-12006-250MG	250MG	Halosulfuron-methyl solution	100784-20-1	S-12142A1-1ML	1ML
Folicur(TM) Solution	107534-96-3	S-12006A1-1ML	1ML	Haloxypof (free acid)	69806-34-4	N-12143-10MG	10MG
Folicur(TM) Solution	107534-96-3	S-12006K1-1ML	1ML	Haloxypof-P (free acid)	95977-29-0	N-12145-50MG	50MG
Folpet	133-07-3	N-12007-250MG	250MG	Haloxypof-P methyl	72619-32-0	N-12146-50MG	50MG
Folpet solution	133-07-3	S-12007M1-1ML	1ML	Haloxypof-R (free acid)	95977-29-0	N-12145-50MG	50MG
Folpet solution	133-07-3	S-12007U1-1ML	1ML	Haloxypof-R methyl	72619-32-0	N-12146-50MG	50MG
Fomesafen	72178-02-0	N-12008-500MG	500MG	Harmony(TM)	79277-27-3	N-13056-100MG	100MG
Fomesafen solution	72178-02-0	S-12008M1-1ML	1ML	Harmony(TM) Solution	79277-27-3	S-13056A1-1ML	1ML
Fongarid	57646-30-7	N-12978-10MG	10MG	Harness(TM)	34256-82-1	N-11013-100MG	100MG
Fonofos	944-22-9	N-11842-100MG	100MG	Harness(TM) Solution	34256-82-1	S-11013M1-1ML	1ML
Fonofos solution	944-22-9	S-11842A1-1ML	1ML	Heptachlor	76-44-8	N-12147-100MG	100MG
Fonofos solution	944-22-9	S-11842J4-1ML	1ML	Heptachlor (13C10) Solution		S-FC100S-1.2ML	1.2ML
Fonofos solution	944-22-9	S-11842J4-5ML	5ML	Heptachlor solution	76-44-8	S-12147M1-1ML	1ML
Foramsulfuron	173159-57-4	N-12009-100MG	100MG	Heptachlor solution	76-44-8	S-12147M1-5ML	5ML
Foramsulfuron solution	173159-57-4	S-12009A1-1ML	1ML	trans-Heptachlor epoxide	28044-83-9	N-13618-50MG	50MG
Forbid(TM)	283594-90-1	N-13224-100MG	100MG	Heptachlor epoxide (Isomer B)	1024-57-3	N-12148-50MG	50MG
Forchlorfenuron	68157-60-8	N-12010-100MG	100MG	Heptachlor epoxide (Isomer B) (13C1) Solution		S-FC101S-1.2ML	1.2ML
Forchlorfenuron solution	68157-60-8	S-12010A1-1ML	1ML	Heptachlor epoxide (Isomer B) solution	1024-57-3	S-12148M1-1ML	1ML
Formaldehyde (in water)	50-00-0	N-12012-1G	1G	Heptachlor epoxide (Isomer B) solution	1024-57-3	S-12148M1-5ML	5ML
Formalin	50-00-0	N-12012-1G	1G	trans-Heptachlor epoxide solution	28044-83-9	S-13618M1-1ML	1ML
Formetanate hydrochloride	23422-53-9	N-12015-250MG	250MG	cis-Heptachlor epoxide	1024-57-3	N-12148-50MG	50MG
Formetanate hydrochloride solution	23422-53-9	S-12015A1-1ML	1ML	cis-Heptachlor epoxide Solution	1024-57-3	S-12148M1-1ML	1ML
Formothion	2540-82-1	N-12017-100MG	100MG	cis-Heptachlor epoxide Solution	1024-57-3	S-12148M1-5ML	5ML
Formothion solution	2540-82-1	S-12017A1-1ML	1ML	trans-Heptachlor epoxide (Isomer A)	28044-83-9	N-13618-50MG	50MG
Formothion solution	2540-82-1	S-12017J1-1ML	1ML	trans-Heptachlor epoxide (Isomer A) Solution	28044-83-9	S-13618M1-1ML	1ML
Fortress(TM)	54593-83-8	N-11427-250MG	250MG	Heptachlor-tetrahydro-methanoindene	76-44-8	N-12147-100MG	100MG
Fortress(TM) Solution	54593-83-8	S-11427A1-1ML	1ML	Heptachlor-tetrahydro-methanoindene Solution	76-44-8	S-12147M1-1ML	1ML
Fortress(TM) Solution	54593-83-8	S-11427U1-1ML	1ML	Heptachlor-tetrahydro-methanoindene Solution	76-44-8	S-12147M1-5ML	5ML
Fosamine Ammonium solution	25954-13-6	N-12018F1-1ML	1ML	Heptenofos	23560-59-0	N-12153-100MG	100MG
Fosetyl-Al	39148-24-8	N-12019-100MG	100MG	Heptenofos Solution	23560-59-0	S-12153A1-1ML	1ML
Fosetyl-Aluminum(Technical)	39148-24-8	N-12019-100MG	100MG	Heptenofos solution	23560-59-0	S-12153U1-1ML	1ML
Fosthiazate	98886-44-3	N-12020-30MG	30MG	Heptenophos	23560-59-0	N-12153-100MG	100MG
Fosthiazate solution	98886-44-3	S-12020A1-1ML	1ML	Heptenophos solution	23560-59-0	S-12153A1-1ML	1ML
Fosthiazate solution	98886-44-3	S-12020U1-1ML	1ML	Heptenophos solution	23560-59-0	S-12153U1-1ML	1ML
Frontier(TM)	87674-68-8	N-11755-100MG	100MG	Herbazolin	3813-05-600	N-12901-100MG	100MG
Frontier(TM) Solution	87674-68-8	S-11755A1-1ML	1ML	Heritage(TM)	131860-33-8	N-11113-100MG	100MG
Frontier(TM) Solution	87674-68-8	S-11755U1-1ML	1ML	Heritage(TM) Solution	131860-33-8	S-11113A1-1ML	1ML
Frunax-DS(TM)	56073-07-5	N-11719-100MG	100MG	Heritage(TM) Solution	131860-33-8	S-11113U1-1ML	1ML
Frunax-DS(TM) Solution	56073-07-5	S-11719M1-1ML	1ML	Heteroauxin	87-51-4	N-10727-250MG	250MG
Fthalide	27355-22-2	N-13085-10MG	10MG	Heteroauxin Solution	87-51-4	S-10727A1-1ML	1ML
Fuberidazole	3878-19-1	N-12021-100MG	100MG	Heteroauxin Solution	87-51-4	S-10727T1-1ML	1ML
Fuberidazole solution	3878-19-1	S-12021M1-1ML	1ML	Hexachloroacetone	116-16-5	N-12158-500MG	500MG
Furacn(TM)	82560-54-1	N-11137-100MG	100MG	Hexachloroacetone solution	116-16-5	S-12158A1-1ML	1ML
Furalaxyl	57646-30-7	N-12978-10MG	10MG	Hexachloroacetone solution	116-16-5	S-12158J1-1ML	1ML
2-Furanacrylic acid	539-47-9	N-10363-100MG	100MG	Hexachlorobenzene	118-74-1	N-12159-250MG	250MG
2-Furanacrylic acid solution	539-47-9	S-10363A1-1ML	1ML	Hexachlorobenzene solution	118-74-1	S-12159K4-1ML	1ML
2-Furanacrylic acid solution	539-47-9	S-10363T1-1ML	1ML	Hexachlorobenzene solution	118-74-1	S-12159M1-1ML	1ML
6-Furfurylamino-9H-purine	525-79-1	N-12294-250MG	250MG	Hexachlorobenzene solution	118-74-1	S-12159M1-5ML	5ML
Furilazole	121776-33-8	N-13817-10MG	10MG	Hexachlorodimethyl sulfone	3064-70-8	N-12161-250MG	250MG
Furmecyclox	60568-05-0	N-12979-10MG	10MG	Hexachlorodimethyl sulfone	3064-70-8	N-12161-250MG	250MG
Furmecyclox solution	60568-05-0	S-12979M1-1ML	1ML	Hexachlorodimethyl sulfone solution	3064-70-8	S-12161A1-1ML	1ML
Furmetamid	60568-05-0	N-12979-10MG	10MG	Hexachlorodimethyl sulfone solution	3064-70-8	S-12161U1-1ML	1ML
Furmetamid Solution	60568-05-0	S-12979M1-1ML	1ML	Hexachlorophene	70-30-4	N-12163-250MG	250MG
Fury(TM)	52315-07-8	N-13754-250MG	250MG	Hexachlorophene solution	70-30-4	S-12163M1-1ML	1ML
2-(2-Furyl)benzimidazole Solution	3878-19-1	S-12021M1-1ML	1ML	Hexachlorophene solution	70-30-4	S-12163M1-5ML	5ML
Fusilade(TM)	69806-50-4	N-11978-250MG	250MG	Hexachloro-2-propanone Solution	116-16-5	S-12158A1-1ML	1ML
Fusilade(TM) Solution	69806-50-4	S-11978A1-1ML	1ML	Hexaconazole	79983-71-4	N-12165-100MG	100MG
Fusilade(TM) Solution	69806-50-4	S-11978U1-1ML	1ML	Hexaconazole solution	79983-71-4	S-12165A1-1ML	1ML
Galben(TM)	71626-11-4	N-11133-250MG	250MG	Hexaconazole solution	79983-71-4	S-12165U1-1ML	1ML
Galben(TM) Solution	71626-11-4	S-11133M1-1ML	1ML	Hexaflumuron	86479-06-3	N-12167-250MG	250MG
Galtak(TM)	25059-80-7	N-11134-250MG	250MG	Hexaflumuron solution	86479-06-3	S-12167A1-1ML	1ML
Galtak(TM) Solution	25059-80-7	S-11134A1-1ML	1ML	Hexaflumuron solution	86479-06-3	S-12167U1-1ML	1ML
Galtak(TM) Solution	25059-80-7	S-11134T1-1ML	1ML	Hexazinone	51235-04-2	N-13743-250MG	250MG
gamma-Cyhalothrin	76703-62-3	N-12106-100MG	100MG	Hexazinone Solution	51235-04-2	S-13743A1-1ML	1ML
Genite (TM)	97-16-5	N-10628-250MG	250MG	Hexazinone Solution	51235-04-2	S-13743T1-1ML	1ML
Genite(TM) Solution	97-16-5	S-10628A1-1ML	1ML				

Description	CAS#	Part #	Size	Description	CAS#	Part #	Size
Hexazinone Solution	51235-04-2	S-13743T1-5ML	5ML	Irgarol(TM) Solution	28159-98-0	S-11520U1-1ML	1ML
Hexide	70-30-4	N-12163-250MG	250MG	Isazophos	42509-80-8	N-12222-100MG	100MG
Hexide Solution	70-30-4	S-12163M1-1ML	1ML	Isazophos solution	42509-80-8	S-12222U1-1ML	1ML
Hexide Solution	70-30-4	S-12163M1-5ML	5ML	Isobenzan	297-78-9	N-12231-10MG	10MG
Hexythiazox	78587-05-0	N-12184-100MG	100MG	Isobenzan solution	297-78-9	S-12231M1-1ML	1ML
Hexythiazox solution	78587-05-0	S-12184A1-1ML	1ML	Isobornyl thiocyanacetate	115-31-1	N-12232-250MG	250MG
Hinocloa	73250-68-7	N-12364-250MG	250MG	Isobornyl thiocyanacetate solution	115-31-1	S-12232A1-1ML	1ML
Hinocloa Solution	73250-68-7	S-12364M1-1ML	1ML	Isobornyl thiocyanacetate solution	115-31-1	S-12232U1-1ML	1ML
Hostathion(TM)	24017-47-8	S-13639-100MG	100MG	Isobutyl 2,4-dichlorophenoxyacetate	1713-15-1	N-10532-250MG	250MG
Hostathion(TM) Solution	24017-47-8	S-13639A1-1ML	1ML	Isobutyl 2,4-dichlorophenoxyacetate Solution	1713-15-1	S-10532A1-1ML	1ML
Hostathion(TM) Solution	24017-47-8	S-13639U1-1ML	1ML	Isobutyl 2,4-dichlorophenoxyacetate Solution	1713-15-1	S-10532U1-1ML	1ML
Hostathion(TM) Solution	24017-47-8	S-13639U1-5ML	5ML	Isocarbophos solution	24353-61-5	S-13933A1-1ML	1ML
Hydramethylinon	67485-29-4	N-11063-100MG	100MG	Isodrin	465-73-6	N-12248-100MG	100MG
Hydramethylinon Solution	67485-29-4	S-11063A1-1ML	1ML	Isodrin solution	465-73-6	S-12248M1-1ML	1ML
Hydramethylinon Solution	67485-29-4	S-11063J1-1ML	1ML	Isodrin solution	465-73-6	S-12248M1-5ML	5ML
2-Hydrazino-4,6-dimethylpyrimidine	23906-13-0	MET-14030A-100MG	100MG	Isophenphos	25311-71-1	N-12249-100MG	100MG
2-Hydrazino-4,6-dimethylpyrimidine Solution	23906-13-0	MET-14030AM1-1ML	1ML	Isophenphos solution	25311-71-1	S-12249A1-1ML	1ML
(7S)-Hydroprene	65733-18-8	N-13189-50MG	50MG	Isophenphos solution	25311-71-1	S-12249U1-1ML	1ML
S-Hydroprene	65733-18-8	N-13189-50MG	50MG	Isophenphos-des-N-isopropyl	25205-08-7	MET-12249C-100MG	100MG
(7S)-Hydroprene Solution	65733-18-8	S-13189A1-1ML	1ML	Isophenphos-des-N-isopropyl Solution	25205-08-7	MET-12249AJ1-1ML	1ML
(7S)-Hydroprene solution	65733-18-8	S-13189A1-1ML	1ML	Isophenphos-des-N-isopropyl-O-analogue Solution	31120-83-9	MET-12249BJ1-1ML	1ML
(7S)-Hydroprene Solution	65733-18-8	S-13189J1-1ML	1ML	Isophenphos-Methyl	99675-03-3	N-12250-50MG	50MG
S-Hydroprene solution	65733-18-8	S-13189J1-1ML	1ML	Isotetamid	875915-78-9	N-14079-100MG	100MG
7-Hydroxycoumarin	93-35-6	MET-13750A-250MG	250MG	Isomalathion Solution	3344-12-500	MET-12346AJ1-1ML	1ML
4-Hydroxycoumarin	1076-38-6	MET-13750B-250MG	250MG	Isoprocab	2631-40-5	N-12264-250MG	250MG
p-Hydroxy-N,N-dimethylbenzene sulfonamide	15020-57-2	MET-11944B-10MG	10MG	Isoprocab solution	2631-40-5	S-12264A1-1ML	1ML
2-Hydroxyethyl-n-octyl sulphide	3547-33-9	N-10378-100MG	100MG	Isoprocab solution	2631-40-5	S-12264U1-1ML	1ML
2-Hydroxyethyl-n-octyl sulphide solution	3547-33-9	S-10378M1-1ML	1ML	Isopropalin	33820-53-0	N-12265-250MG	250MG
9-Hydroxyfluorene-9-carboxylic acid	467-69-6	N-13815-250MG	250MG	Isopropalin solution	33820-53-0	S-12265A1-1ML	1ML
3-(Hydroxymethyl)diphenyl ether	13826-35-2	MET-12848B-100MG	100MG	Isopropalin solution	33820-53-0	S-12265J4-1ML	1ML
Hydroxyphenyl mercury	100-57-2	N-13014-250MG	250MG	Isopropalin solution	33820-53-0	S-12265J4-5ML	5ML
3-Hydroxy-5-phenylisoxazole	939-05-9	MET-13932A-10MG	10MG	2-Isopropoxyphenol	4812-20-8	MET-11128A-250MG	250MG
2-Hydroxypyrazine	6270-63-9	MET-13760BM1-1ML	1ML	Isopropyl 3,4-diethoxycarbanilate	87130-20-9	N-11691-250MG	250MG
2-Hydroxypyrazine Solution	6270-63-9	MET-13760BM1-1ML	1ML	Isopropyl 3,4-diethoxycarbanilate Solution	87130-20-9	S-11691M1-1ML	1ML
8-Hydroxyquinoline	148-24-3	N-10965-250MG	250MG	2-Isopropyl-4-methyl-6-hydroxypyrimidine	2814-20-2	MET-11621B-250MG	250MG
8-Hydroxyquinoline solution	148-24-3	S-10965M1-1ML	1ML	2-Isopropyl-6-methyl-4-pyrimidinol	2814-20-2	MET-11621B-250MG	250MG
2-Hydroxy-1,3,5-triazine-4,6-diamine	645-92-1	MET-11106E-50MG	50MG	Isoprothiolane	50512-35-1	N-12278-250MG	250MG
Hymexazole	10004-44-1	N-12196-25MG	25MG	Isoprothiolane solution	50512-35-1	S-12278A1-1ML	1ML
Hymexazole solution	10004-44-1	S-12196M1-1ML	1ML	Isoprothiolane solution	50512-35-1	S-12278T1-1ML	1ML
Icaridin	119515-38-7	N-12197-100MG	100MG	Isoproturon	34123-59-6	N-12279-100MG	100MG
Icaridin solution	119515-38-7	S-12197A1-1ML	1ML	Isoproturon solution	34123-59-6	S-12279A1-1ML	1ML
Imazalil	35554-44-0	N-12198-100MG	100MG	Isoproturon monodesmethyl	34123-57-4	MET-12279A-250MG	250MG
Imazalil solution	35554-44-0	S-12198A1-1ML	1ML	Isoqramaz solution	881685-58-1	S-13996A1-1ML	1ML
Imazalil sulfate	35554-44-0	S-12198J1-1ML	1ML	Isoqramaz solution	55861-78-4	S-13971A1-1ML	1ML
Imazalil sulfate solution	35554-44-0	N-12199-100MG	100MG	2-Isovaleryl-1,3-indanedione	83-28-3	N-10384-250MG	250MG
Imazamethabenz (free acid) Solution	58594-72-2	S-12199F1-1ML	1ML	2-Isovaleryl-1,3-indanedione solution	83-28-3	S-10384A1-1ML	1ML
Imazamethabenz-methyl	100728-84-5	MET-12200AA0-1ML	1ML	2-Isovaleryl-1,3-indanedione solution	83-28-3	S-10384U1-1ML	1ML
Imazamethabenz-methyl solution	81405-85-8	N-12200-100MG	100MG	Isoxaben	82558-50-7	N-12285-250MG	250MG
Imazamethabenz-methyl solution	81405-85-8	S-12200A1-1ML	1ML	Isoxaben solution	82558-50-7	S-12285M1-1ML	1ML
Imazamox	81405-85-8	S-12200U1-1ML	1ML	Isoxadifen-ethyl	163520-33-0	N-12286-100MG	100MG
Imazamox solution	114311-32-9	N-12201-250MG	250MG	Isoxafutole	141112-29-0	N-12287-100MG	100MG
Imazamox solution	114311-32-9	S-12201A1-1ML	1ML	Isoxafutole solution	141112-29-0	S-12287U1-1ML	1ML
Imazapic	104098-48-8	S-12201T1-1ML	1ML	Isoxathion	18854-01-8	N-13932-25MG	25MG
Imazapic solution	104098-48-8	N-12202-250MG	250MG	Itraconazole	84625-61-6	N-13958-100MG	100MG
Imazapyr	104098-48-8	S-12202A1-1ML	1ML	Ivermectin	70288-86-7	N-12289-100MG	100MG
Imazapyr solution	81334-34-1	N-12203-100MG	100MG	Ivermectin solution	70288-86-7	S-12289A1-1ML	1ML
Imazaquin	81334-34-1	S-12203A1-1ML	1ML	Ivermectin solution	70288-86-7	S-12289U1-1ML	1ML
Imazaquin solution	81335-37-7	N-12204-100MG	100MG	Kadethrine	58769-20-3	N-12947-10MG	10MG
Imazaquin solution	81335-37-7	S-12204A1-1ML	1ML	Kanemite(TM)	57960-19-7	N-11003-250MG	250MG
Imazaquin solution	81335-37-7	S-12204T1-1ML	1ML	Kanemite(TM) Solution	57960-19-7	S-11003A1-1ML	1ML
Imazethapyr	81335-77-5	N-12205-100MG	100MG	Kanemite(TM) Solution	57960-19-7	S-11003U1-1ML	1ML
Imazethapyr solution	81335-77-5	S-12205M1-1ML	1ML	Karate (TM)	91465-08-6	N-12307-100MG	100MG
Imazosulfuron	122548-33-8	N-13818-100MG	100MG	Karate (TM) Solution	91465-08-6	S-12307A1-1ML	1ML
Imazosulfuron Solution	122548-33-8	S-13818A1-1ML	1ML	Karate (TM) Solution	91465-08-6	S-12307J1-1ML	1ML
Imibenconazole solution	86598-92-7	S-13968A1-1ML	1ML	Karbutilate	4849-32-5	N-12845-250MG	250MG
Imidacloprid	138261-41-3	N-12206-100MG	100MG	Kelthane (T.M.)	115-32-2	N-12290-100MG	100MG
Imidacloprid solution	138261-41-3	S-12206M1-1ML	1ML	Kepone (TM) (13C10) Solution		S-FC2441S-1.2ML	1.2ML
Imidacloprid olefin	115086-54-9	MET-12206-10MG	10MG	Kerb (TM)	23950-58-5	N-13136-250MG	250MG
Imidacloprid urea	120868-66-8	MET-12206A-50MG	50MG	Kerb (TM) Solution	23950-58-5	S-13136M1-1ML	1ML
Imidan (T.M.)	732-11-6	N-12207-100MG	100MG	Kerb (TM) Solution	23950-58-5	S-13136M1-5ML	5ML
Imidan (T.M.) Solution	732-11-6	S-12207A1-1ML	1ML	Kilval(TM)	2275-23-2	N-13741-100MG	100MG
Imidan (T.M.) Solution	732-11-6	S-12207J1-1ML	1ML	Kilval(TM) Solution	2275-23-2	S-13741A1-1ML	1ML
Imidan (T.M.) Solution	732-11-6	S-12207J1-5ML	5ML	Kinetin	525-79-1	N-12294-250MG	250MG
2-Imidazolidinethione	96-45-7	N-10381-250MG	250MG	Kinoprene	42588-37-4	N-12295-100MG	100MG
2-Imidazolidinethione solution	96-45-7	S-10381A7-1ML	1ML	Kinoprene solution	42588-37-4	S-12295M1-1ML	1ML
2-Imidazolidinethione solution	96-45-7	S-10381A7-5ML	5ML	Kinoprene solution	42588-37-4	S-12295T4-1ML	1ML
2-Imidazolidinethione solution	96-45-7	S-10381H1-1ML	1ML	Kinoprene solution	42588-37-4	S-12295T4-5ML	5ML
2-Imidazolidinethione solution	96-45-7	S-10381H1-5ML	5ML	Kresoxim-methyl	143390-89-0	N-12296-250MG	250MG
Iminobis (methyl)phosphonic acid)	6711-48-4	MET-12135A-500MG	500MG	Kresoxim-methyl solution	143390-89-0	S-12296A1-1ML	1ML
Imiprothrin	72963-72-5	N-14006-50MG	50MG	Kresoxim-methyl solution	143390-89-0	S-12296U1-1ML	1ML
Imiprothrin solution	72963-72-5	S-14006A1-1ML	1ML	Lactofen	77501-63-4	N-12304-100MG	100MG
Impact(TM)	210631-68-8	N-13585-100MG	100MG	Lactofen solution	77501-63-4	S-12304A1-1ML	1ML
IMPRELIS™	858956-08-8	N-13826-250MG	250MG	Lactofen solution	77501-63-4	S-12304U1-1ML	1ML
Indanofan	133220-30-1	N-12209-100MG	100MG	[Lactoxyloxy]phenylmercury	122-66-4	N-13015-250MG	250MG
Indaziflam	950782-86-2	N-13074-100MG	100MG	Lead arsenate(Technical)	7784-40-9	N-12310-100MG	100MG
Indaziflam solution	950782-86-2	S-13074A1-1ML	1ML	Legacy(TM)	83164-33-4	N-11723-100MG	100MG
b-Indeppropionic acid	830-96-6	N-10729-250MG	250MG	Lenacil	2164-08-100	N-12314-100MG	100MG
b-Indeppropionic acid Solution	830-96-6	S-10729A1-1ML	1ML	Lenacil solution	2164-08-100	S-12314M1-1ML	1ML
b-Indeppropionic acid Solution	830-96-6	S-10729T1-1ML	1ML	Leptophos	21609-90-5	N-12315-100MG	100MG
3-Indoleacetic acid	87-51-4	N-10727-250MG	250MG	Leptophos solution	21609-90-5	S-12315A1-1ML	1ML
3-Indoleacetic acid solution	87-51-4	S-10727A1-1ML	1ML	Leptophos solution	21609-90-5	S-12315U1-1ML	1ML
3-Indoleacetic acid solution	87-51-4	S-10727T1-1ML	1ML	Leptophos oxon	25006-32-0	MET-12315A-50MG	50MG
3-Indolebutyric acid	133-32-4	N-10728-100MG	100MG	Lethane (TM)	112-56-1	N-10260-250MG	250MG
3-Indolebutyric acid solution	133-32-4	S-10728A1-1ML	1ML	Lethane 384 (TM) Solution	112-56-1	S-10260A1-1ML	1ML
3-Indolebutyric acid solution	133-32-4	S-10728T1-1ML	1ML	Lethane 384 (TM) Solution	112-56-1	S-10260J4-1ML	1ML
3-Indolepropionic acid	830-96-6	N-10729-250MG	250MG	Lethane 384 (TM) Solution	112-56-1	S-10260J4-5ML	5ML
3-Indolepropionic acid solution	830-96-6	S-10729A1-1ML	1ML	Lindane (BHC gamma isomer)	58-89-9	N-12319-250MG	250MG
3-Indolepropionic acid solution	830-96-6	S-10729T1-1ML	1ML	Lindane (BHC gamma isomer) (13C6) Solution		S-FC104S-1.2ML	1.2ML
S-Indoxacarb	173584-44-6	N-12213-100MG	100MG	Lindane (BHC gamma isomer) (13C6, D6) Solution		S-FCD104S-1.2ML	1.2ML
S-Indoxacarb solution	173584-44-6	S-12213A1-1ML	1ML	Lindane (BHC gamma isomer) solution		S-12319A1-1ML	1ML
S-Indoxacarb solution	173584-44-6	S-12213U1-1ML	1ML	Linuron	330-55-2	N-12322-250MG	250MG
Iodosulfuron-methyl-sodium	144550-36-7	N-13923-10MG	10MG	Linuron solution	330-55-2	S-12322A6-1ML	1ML
Ioxynil octanoate	3861-47-0	N-12218-100MG	100MG	Linuron solution	330-55-2	S-12322A6-5ML	5ML
Ioxynil octanoate solution	3861-47-0	S-12218A1-1ML	1ML	Lirinox	94-80-4	N-10611-250MG	250MG
Ioxynil octanoate solution	3861-47-0	S-12218U1-1ML	1ML	Lirinox Solution	94-80-4	S-10611A1-1ML	1ML
Ioxynil-methyl	3336-40-1	MET-10794A-100MG	100MG	Lirinox Solution	94-80-4	S-10611U1-1ML	1ML
Iprobenfos	26087-47-8	N-12219-250MG	250MG	Lontrel™	1702-17-6	N-12324-250MG	250MG
Iprobenfos solution	26087-47-8	S-12219A1-1ML	1ML	Lontrel™ Solution	1702-17-6	S-12324A1-1ML	1ML
Iprobenfos solution	26087-47-8	S-12219U1-1ML	1ML	Lontrel™ Solution	1702-17-6	S-12324T1-1ML	1ML
Iprodione	36734-19-7	N-12220-100MG	100MG	Lufenuron	103055-07-8	N-12325-100MG	100MG
Iprodione solution	36734-19-7	S-12220A1-1ML	1ML	Lufenuron solution	103055-07-8	S-12325A1-1ML	1ML
Iprovalicarb	140923-17-7	N-13924-10MG	10MG	Mach 2(TM)	112228-61-6	N-12141-500MG	500MG
Irgarol(TM)	28159-98-0	N-11520-100MG	100MG	Mach 2(TM) Solution	112228-61-6	S-12141M1-1ML	1ML
Irgarol(TM) Solution	28159-98-0	S-11520A1-1ML	1ML	Malathion	121-75-5	N-12346-100MG	100MG

Description	CAS#	Part #	Size	Description	CAS#	Part #	Size
Malathion (d10) Solution		S-FD2118S-1.2ML	1.2ML	Methaphenamiphos Solution	22224-92-6	S-12854T1-1ML	1ML
Malathion solution	121-75-5	S-12346A1-1ML	1ML	Methaphenamiphos Solution	22224-92-6	S-12854T1-5ML	5ML
Malathion solution	121-75-5	S-12346U1-1ML	1ML	Methasulfocarb	66952-49-6	N-13283-25MG	25MG
Malathion solution	121-75-5	S-12346U1-5ML	5ML	Methasulfocarb solution	66952-49-6	S-13283A1-1ML	1ML
Malathion monocarboxylic acid	35884-76-5	MET-12346D-100MG	100MG	Methidathion	950-37-8	N-12397-250MG	250MG
Malathion monocarboxylic acid Solution	35884-76-5	MET-12346BT1-1ML	1ML	Methidathion solution	950-37-8	S-12397A1-1ML	1ML
Malathion-O-analog	1634-78-2	MET-12346C-100MG	100MG	Methidathion solution	950-37-8	S-12397U1-1ML	1ML
Maleic hydrazide	123-33-1	N-12350-250MG	250MG	Methiocarb	2032-65-7	N-12398-100MG	100MG
Maleic hydrazide solution	123-33-1	S-12350M1-1ML	1ML	Methiocarb solution	2032-65-7	S-12398M1-1ML	1ML
Mandoczeb (Technical)	8018-01-700	N-12353-250MG	250MG	Methiocarb sulfone	2179-25-1	MET-12398A-50MG	50MG
Mandipropamid	374726-62-2	N-12354-100MG	100MG	Methiocarb sulfoxide Solution	2635-10-100	MET-12398BU1-1ML	1ML
Maneb(Technical)	12427-38-2	N-12355-250MG	250MG	Methomyl	16752-77-5	N-12399-100MG	100MG
Maverick(TM) Solution	141776-32-1	N-12355-100MG	100MG	Methomyl solution	16752-77-5	S-12399A1-1ML	1ML
Maverick(TM) Solution	141776-32-1	S-12355A1-1ML	1ML	Methoprene	40596-69-8	N-12400-100MG	100MG
Maverick(TM) Solution	141776-32-1	S-12355U1-1ML	1ML	(S)-Methoprene	65733-16-6	N-13832-100MG	100MG
Mavrik (TM)	102851-06-9	N-13263-100MG	100MG	S-Methoprene	65733-16-6	N-13832-100MG	100MG
Mavrik (TM) Solution	102851-06-9	S-12363A1-1ML	1ML	Methoprene solution	40596-69-8	S-12400T4-1ML	1ML
Mavrik (TM) Solution	102851-06-9	S-12363K1-1ML	1ML	Methoprene solution	40596-69-8	S-12400T4-5ML	5ML
MCPA (TM)	94-74-6	N-10818-250MG	250MG	Methoprotyrine	841-06-5	N-12401-100MG	100MG
MCPA (TM) Solution	94-74-6	S-10818B1-1ML	1ML	Methoprotyrine solution	841-06-5	S-12401A1-1ML	1ML
MCPA (TM) Solution	94-74-6	S-10818B1-5ML	5ML	Methoprotyrine solution	841-06-5	S-12401T1-1ML	1ML
MCPA (TM) Solution	94-74-6	S-10818M1-1ML	1ML	Methothrin	34388-29-9	N-14037-100MG	100MG
MCPA methyl ester	2436-73-9	N-10819-100MG	100MG	o,p'-Methoxychlor	30667-99-3	MET-12403A-100MG	100MG
MCPA methyl ester Solution	2436-73-9	S-10819M1-1ML	1ML	Methoxychlor	72-43-5	N-12403-100MG	100MG
MCPA methyl ester Solution	2436-73-9	S-10819M1-5ML	5ML	o,p'-Methoxychlor Solution	72-43-5	MET-12403AM1-1ML	1ML
MCPA sodium salt	3653-48-3	MET-12348-250MG	250MG	Methoxychlor solution	72-43-5	S-12403M1-1ML	1ML
MCPB	94-81-5	N-12357-100MG	100MG	Methoxychlor solution	72-43-5	S-12403M1-5ML	5ML
MCPB solution	94-81-5	S-12357A1-1ML	1ML	p,p'-Methoxychlor-olefin	2132-70-9	MET-12403B-50MG	50MG
MCPB solution	94-81-5	S-12357T1-1ML	1ML	p,p'-Methoxychlor-olefin Solution	2132-70-9	MET-12403BM1-1ML	1ML
MCPB methyl ester	57153-18-1	N-12358-250MG	250MG	2-Methoxyethyl 2-[4-(tert-butyl)phenyl]cyanoacetate	400882-57-7	MET-13987A-50MG	50MG
MCPB methyl ester solution	57153-18-1	S-12358A1-1ML	1ML	Methoxyfenozide	161050-58-4	N-12406-100MG	100MG
MCPB methyl ester solution	57153-18-1	S-12358T1-1ML	1ML	Methoxyfenozide solution	161050-58-4	S-12406A1-1ML	1ML
MCPB	7085-19-0	N-12360-100MG	100MG	Methyl anthranilate	134-20-3	N-12414-100MG	100MG
MCPB Solution	7085-19-0	S-12360A1-1ML	1ML	Methyl anthranilate solution	134-20-3	S-12414M1-1ML	1ML
MCPB Solution	7085-19-0	S-12360B1-1ML	1ML	Methylarsonic acid	124-58-3	N-12922-100MG	100MG
MCPB methyl ester	23844-56-6	N-12361-100MG	100MG	1-Methyl-3,5-diphenyl-1H-pyrazolium		MET-11721C-100MG	100MG
MCPB methyl ester Solution	23844-56-6	S-12361M1-1ML	1ML	S,S-Methylene bis[O,O-diethyl]phosphorodithioate	563-12-2	N-11873-250MG	250MG
MCPB methyl ester Solution	23844-56-6	S-12361M1-5ML	5ML	S,S-Methylene bis[O,O-diethyl]phosphorodithioate	563-12-2	S-11873M1-1ML	1ML
MCPB methyl ester Solution	23844-56-6	S-12362B1-1ML	1ML	S,S-Methylene bis[O,O-diethyl]phosphorodithioate	563-12-2	S-11873M1-5ML	5ML
MCPB methyl ester Solution	23844-56-6	S-12362B1-5ML	5ML	N-(1-Methylethyl)-N-(phenyl)oxalamic acid	70628-36-3	MET-13104B-10MG	10MG
Mecarbam	2595-54-2	N-12359-50MG	50MG	N-(1-Methylethyl)-N-(phenyl)oxalamic acid	70628-36-3	MET-13104BM1-1ML	1ML
Mecarbam solution	2595-54-2	S-12359A1-1ML	1ML	2-Methyl-3-furic acid	6947-94-0	MET-14013A-100MG	100MG
Mecarbam solution	2595-54-2	S-12359U1-1ML	1ML	Methyl isothiocyanate	556-61-6	MET-12392A-250MG	250MG
Mecoprop	7085-19-0	N-12360-100MG	100MG	N-Methyl-methamidophos	28167-49-9	N-12558-50MG	50MG
Mecoprop solution	7085-19-0	S-12360A1-1ML	1ML	N-Methyl-methamidophos solution	28167-49-9	S-12558A1-1ML	1ML
Mecoprop solution	7085-19-0	S-12360B1-1ML	1ML	N-Methyl-methamidophos solution	28167-49-9	S-12558U1-1ML	1ML
Mecoprop methyl ester	23844-56-6	N-12361-100MG	100MG	2-Methyl-2-(methylthio)-propionaldehyde oxime	1646-75-9	MET-11044A-10MG	10MG
Mecoprop methyl ester solution	23844-56-6	S-12361M1-1ML	1ML	Methyl mustard oil	556-61-6	MET-12392A-250MG	250MG
Mecoprop methyl ester solution	23844-56-6	S-12361M1-5ML	5ML	Methyl-1-naphthalene acetate	2876-78-0	N-12464-500MG	500MG
Mecoprop-2,4,4-trimethylpentyl ester solution	23844-56-6	S-12362B1-1ML	1ML	Methyl-1-naphthalene acetate solution	2876-78-0	S-12464M1-1ML	1ML
Mecoprop-2,4,4-trimethylpentyl ester solution	23844-56-6	S-12362B1-5ML	5ML	Methyl-1-naphthalene acetate solution	2876-78-0	S-12464U1-1ML	1ML
Mecoprop-P	16484-77-8	N-12363-250MG	250MG	N-Methyl-N-1-naphthyl acetamide	5903-13-9	N-12559-100MG	100MG
Mecoprop-P solution	16484-77-8	S-12363A1-1ML	1ML	N-Methyl-N-1-naphthyl acetamide solution	5903-13-9	S-12559M1-1ML	1ML
Mecoprop-P solution	16484-77-8	S-12363T1-1ML	1ML	3-Methyl-4-nitrophenol	2581-34-2	MET-11955B-100MG	100MG
Mefenacet	73250-68-7	N-12364-250MG	250MG	Methyl nonyl ketone	112-12-9	N-10520-1G	1G
Mefenacet solution	73250-68-7	S-12364M1-1ML	1ML	Methyl nonyl ketone Solution	112-12-9	S-10520N1-1ML	1ML
Mefenoxam	70630-17-0	N-12365-100MG	100MG	Methyl paraoxon	950-35-6	N-11775-100MG	100MG
Mefenoxam solution	70630-17-0	S-12365A1-1ML	1ML	Methyl paraoxon Solution	950-35-6	S-11775A1-1ML	1ML
Mefenoxam solution	70630-17-0	S-12365T1-1ML	1ML	Methyl paraoxon Solution	950-35-6	S-11775T1-1ML	1ML
Mefenpyr-diethyl	135590-91-9	N-12366-100MG	100MG	Methyl paraoxon Solution	950-35-6	S-11775T1-5ML	5ML
Mefenpyr-diethyl solution	135590-91-9	S-12366A1-1ML	1ML	Methyl parathion	298-00-0	N-12452-100MG	100MG
Mefenpyr-diethyl solution	135590-91-9	S-12366U1-1ML	1ML	Methyl parathion solution	298-00-0	S-12452A1-1ML	1ML
Mefluidide	53780-34-0	N-12367-100MG	100MG	Methyl parathion solution	298-00-0	S-12452U1-1ML	1ML
Mefluidide solution	53780-34-0	S-12367M1-1ML	1ML	1-(1-Methyl-1-phenylethyl)-3-p-tolylurea	42609-52-9	N-13955-100MG	100MG
Mepanipyrin	110235-47-7	N-12369-100MG	100MG	N-Methylpiperidine	626-67-5	MET-12371A-250MG	250MG
Mephosfolan	950-10-7	N-12370-100MG	100MG	N-(4-methyl-6-prop-1-ynylpyrimidin-2-yl)aniline	110235-47-7	N-12369-100MG	100MG
Mephosfolan solution	950-10-7	S-12370A1-1ML	1ML	2-Methylsulfonyl-4-trifluoromethyl benzoic acid	142994-06-7	MET-12287B-100MG	100MG
Mephosfolan solution	950-10-7	S-12370U1-1ML	1ML	4-Methylsulfonylphenol	14763-60-1	MET-11964G-50MG	50MG
Mepiquat chloride	24307-26-4	N-12371-100MG	100MG	N-Methyl-Tamaron(TM)	28167-49-9	N-12558-50MG	50MG
Mepiquat chloride solution	24307-26-4	S-12371F1-1ML	1ML	N-Methyl-Tamaron(TM) Solution	28167-49-9	S-12558A1-1ML	1ML
Mepronil	55814-41-0	N-12372-100MG	100MG	N-Methyl-Tamaron(TM) Solution	28167-49-9	S-12558U1-1ML	1ML
Mepronil solution	55814-41-0	S-12372M1-1ML	1ML	Methyl o-trifluoromethylbenzoate	344-96-7	MET-12004A-100MG	100MG
Meptyldinocap	131-72-6	N-12388-25MG	25MG	Metiram(Technical)	9006-42-2	N-13086-250MG	250MG
Meptyldinocap solution	131-72-6	S-12388M1-1ML	1ML	Metobromuron	3060-89-7	N-12477-250MG	250MG
2-Mercaptobenzothiazole monoethanolamine salt	5902-85-2	N-13742-250MG	250MG	Metobromuron solution	3060-89-7	S-12477A1-1ML	1ML
2-Mercaptobenzothiazole monoethanolamine salt	5902-85-2	S-13742A1-1ML	1ML	Metobromuron solution	3060-89-7	S-12477U1-1ML	1ML
2-Mercaptobenzothiazole monoethanolamine salt	5902-85-2	S-13742U1-1ML	1ML	Metolachlor	51218-45-2	N-12478-100MG	100MG
Mercusol	28086-13-7	N-13033-250MG	250MG	S-Metolachlor	87392-12-9	N-13191-100MG	100MG
Mercusol Solution	28086-13-7	S-13033A1-1ML	1ML	Metolachlor (ring 13C6) Solution		S-FC2241S-1.2ML	1.2ML
Mercusol Solution	28086-13-7	S-13033U1-1ML	1ML	Metolachlor solution	51218-45-2	S-12478A1-1ML	1ML
Merphenyl nitrate	55-68-5	N-13032-100MG	100MG	Metolachlor solution	51218-45-2	S-12478T1-1ML	1ML
Mesoprazine	1824-09-5	S-13296-50MG	50MG	Metolachlor solution	51218-45-2	S-12478T1-5ML	5ML
Mesoprazine solution	1824-09-5	S-13296A1-1ML	1ML	S-Metolachlor solution	87392-12-9	S-13191A1-1ML	1ML
Mesosulfuron-methyl	208465-21-8	N-12376-100MG	100MG	S-Metolachlor solution	87392-12-9	S-13191T1-1ML	1ML
Mesosulfuron-methyl solution	208465-21-8	S-12376A1-1ML	1ML	Metolachlor ESA sodium salt	947601-85-6	MET-12478C-25MG	25MG
Mesotrione	104206-82-8	N-12377-100MG	100MG	Metolachlor ESA sodium salt Solution	947601-85-6	MET-12478CM1-1ML	1ML
Metallumzone	139968-49-3	N-12379-100MG	100MG	Metolachlor OA	152019-73-3	MET-12478D-10MG	10MG
Metalaxyl	57837-19-1	N-12380-100MG	100MG	Metolachlor OA Solution	152019-73-3	MET-12478DM1-1ML	1ML
Metalaxyl solution	57837-19-1	S-12380A1-1ML	1ML	Metolachlor oxanilic acid	152019-73-3	MET-12478D-10MG	10MG
Metalaxyl solution	57837-19-1	S-12380U1-1ML	1ML	Metolachlor oxanilic acid	152019-73-3	MET-12478DM1-1ML	1ML
Metalaxyl-M	70630-17-0	N-12365-100MG	100MG	Metolcarb	1129-41-5	N-12479-100MG	100MG
Metalaxyl-M Solution	70630-17-0	S-12365A1-1ML	1ML	Metolcarb solution	1129-41-5	S-12479A1-1ML	1ML
Metalaxyl-M Solution	70630-17-0	S-12365T1-1ML	1ML	Metolcarb solution	1129-41-5	S-12479U1-1ML	1ML
Metaldehyde(Technical)	9002-91-9	N-12381-1G	1G	Metometuron	101200-48-0	N-13640-250MG	250MG
Metam sodium	6734-80-1	N-12382-250MG	250MG	Metometuron Solution	101200-48-0	S-13640A1-1ML	1ML
Metamifop	256412-89-2	N-13084-25MG	25MG	Metoxuron	19937-59-8	N-12480-250MG	250MG
Metamitron	41394-05-2	N-12383-250MG	250MG	Metoxuron solution	19937-59-8	S-12480A1-1ML	1ML
Metamitron solution	41394-05-2	S-12383M1-1ML	1ML	Metoxuron solution	19937-59-8	S-12480U1-1ML	1ML
Metasystox I (TM)	919-86-8	N-12384-50MG	50MG	Metrafenone	220899-03-6	N-12897-100MG	100MG
Metasystox I (TM) solution	919-86-8	S-12384A1-1ML	1ML	Metribuzin	21087-64-9	N-12481-250MG	250MG
Metasystox I (TM) solution	919-86-8	S-12384U1-1ML	1ML	Metribuzin solution	21087-64-9	S-12481A1-1ML	1ML
Metazachlor	67129-08-2	N-12385-100MG	100MG	Metribuzin solution	21087-64-9	S-12481T1-1ML	1ML
Metconazole	125116-23-6	N-12386-100MG	100MG	Metribuzin solution	21087-64-9	S-12481T1-5ML	5ML
Methabenzthiazuron	18691-97-9	N-12387-100MG	100MG	Metsulfuron methyl	74223-64-6	N-12482-100MG	100MG
Methabenzthiazuron solution	18691-97-9	S-12387A1-1ML	1ML	Metsulfuron methyl solution	74223-64-6	S-12482A1-1ML	1ML
Methabenzthiazuron solution	18691-97-9	S-12387U1-1ML	1ML	E-Mevinphos	26718-65-0	N-11481-100MG	100MG
Methacrifos	30864-28-9	N-12388-100MG	100MG	Z-Mevinphos	26718-65-0	N-12906-50MG	50MG
Methacrifos solution	30864-28-9	S-12388A1-1ML	1ML	Mevinphos	7786-34-7	N-13037-250MG	250MG
Methacrifos solution	30864-28-9	S-12388U1-1ML	1ML	E-Mevinphos solution	26718-65-0	S-11481A1-1ML	1ML
Methamidophos	10265-92-6	N-12393-100MG	100MG	E-Mevinphos solution	26718-65-0	S-11481T1-1ML	1ML
Methamidophos solution	10265-92-6	S-12393A1-1ML	1ML	Mevinphos solution	7786-34-7	S-13037A1-1ML	1ML
Methamidophos solution	10265-92-6	S-12393U1-1ML	1ML	Mevinphos solution	7786-34-7	S-13037U1-1ML	1ML
Methaphenamiphos	22224-92-6	N-12854-250MG	250MG	Mevinphos solution	7786-34-7	S-13037U1-5ML	5ML
Methaphenamiphos Solution	22224-92-6	S-12854A1-1ML	1ML	MGK 264 (TM)	113-48-4	N-12483-250MG	250MG

Description	CAS#	Part #	Size	Description	CAS#	Part #	Size
MGK 264 (TM) solution	113-48-4	S-12483A1-1ML	1ML	Nicrazine	330-95-0	N-12651-100MG	100MG
MGK 264 (TM) solution	113-48-4	S-12483T1-1ML	1ML	Nimbin	5945-86-8	N-12912-5MG	5MG
MGK 264 (TM) solution	113-48-4	S-12483T1-5ML	5ML	Nimrod(TM)	41483-43-6	N-11349-100MG	100MG
MGK Repellant 326	136-45-8	N-12484-250MG	250MG	Nirothale	10552-74-6	N-12665-100MG	100MG
Mipcin*	2631-40-5	N-12264-250MG	250MG	Nitenpyram	150824-47-8	N-12657-100MG	100MG
Mipcin* Solution	2631-40-5	S-12264A1-1ML	1ML	Nitenpyram solution	150824-47-8	S-12657M1-1ML	1ML
Mipcin* Solution	2631-40-5	S-12264U1-1ML	1ML	Nitrain	4726-14-1	N-12658-250MG	250MG
Mirex	2385-85-5	N-12486-100MG	100MG	Nitrain solution	4726-14-1	S-12658A1-1ML	1ML
Mirex (13C8) Solution		S-FC2121S-1ML	1ML	Nitrain solution	4726-14-1	S-12658J1-1ML	1ML
Mirex solution	2385-85-5	S-12486M1-1ML	1ML	Nitrapyrin	1929-82-4	N-12659-100MG	100MG
mitotane	53-19-0	N-12706-250MG	250MG	Nitrapyrin solution	1929-82-4	S-12659A1-1ML	1ML
mitotane Solution	53-19-0	S-12706M1-1ML	1ML	asym-Nitro-m-dichlorobenzene	611-06-3	N-10625-250MG	250MG
Molinate	2212-67-1	N-12487-250MG	250MG	asym-Nitro-m-dichlorobenzene Solution	611-06-3	S-10625A1-1ML	1ML
Molinate solution	2212-67-1	S-12487A1-1ML	1ML	asym-Nitro-m-dichlorobenzene Solution	611-06-3	S-10625U1-1ML	1ML
Molinate solution	2212-67-1	S-12487T1-1ML	1ML	Nitrofen	1836-75-5	N-12663-100MG	100MG
Molinate solution	2212-67-1	S-12487T1-5ML	5ML	Nitrofen solution	1836-75-5	S-12663M1-1ML	1ML
Monalide	7287-36-7	N-13917-10MG	10MG	N-Nitrosopendimethalin	68897-50-7	MET-13140A-100MG	100MG
Monceren	66063-05-6	N-12488-100MG	100MG	Nitrothal isopropyl	10552-74-6	N-12665-100MG	100MG
Monceren Solution	66063-05-6	S-12488A1-1ML	1ML	Nominee(TM)	125401-92-5	N-11323-100MG	100MG
Monceren Solution	66063-05-6	S-12488U1-1ML	1ML	Nomolt(TM)	83121-18-0	N-13506-250MG	250MG
Moncut (TM)	66332-96-5	N-12004-250MG	250MG	Nomolt(TM) Solution	83121-18-0	S-13506A1-1ML	1ML
Moncut (TM) Solution	66332-96-5	S-12004M1-1ML	1ML	Norflurazon	27314-13-2	N-12668-100MG	100MG
Monocrotophos	6923-22-4	N-12493-100MG	100MG	Norflurazon solution	27314-13-2	S-12668A1-1ML	1ML
Monocrotophos solution	6923-22-4	S-12493A1-1ML	1ML	Norflurazon solution	27314-13-2	S-12668T1-1ML	1ML
Monocrotophos solution	6923-22-4	S-12493U1-1ML	1ML	Norflurazon solution	27314-13-2	S-12668T1-5ML	5ML
Monocrotophos solution	6923-22-4	S-12493U1-5ML	5ML	Noruron	18530-56-8	N-12669-250MG	250MG
Monodesmethyl Isoproturon	34123-57-4	MET-12279A-250MG	250MG	Noruron solution	18530-56-8	S-12669A1-1ML	1ML
Monolinuron	1746-81-2	N-12494-250MG	250MG	Noruron solution	18530-56-8	S-12669U1-1ML	1ML
Monolinuron solution	1746-81-2	S-12494A1-1ML	1ML	Novaluron	116714-46-6	N-12670-100MG	100MG
Monolinuron solution	1746-81-2	S-12494U1-1ML	1ML	Novaluron solution	116714-46-6	S-12670A1-1ML	1ML
Monomethyl tetrachloroterephthalate	887-54-7	MET-11462A-100MG	100MG	Novaluron solution	116714-46-6	S-12670U1-1ML	1ML
Monosodium acid methane arsonate sesquihydrate	2163-80-6	N-12495-100MG	100MG	Nuairimol	63284-71-9	N-13942-100MG	100MG
Monosodium acid methane arsonate sesquihydrate	2163-80-6	S-12495F1-1ML	1ML	Nuairimol solution	63284-71-9	S-13942M1-1ML	1ML
Monothiono TEPP	645-78-3	N-12496-50MG	50MG	Nustar(tm)	85509-19-9	N-12003-100MG	100MG
Monothiono TEPP solution	645-78-3	S-12496U1-1ML	1ML	Nustar(tm) Solution	85509-19-9	S-12003A1-1ML	1ML
Monuron	150-68-5	N-12497-250MG	250MG	Nustar(tm) Solution	85509-19-9	S-12003U1-1ML	1ML
Monuron solution	150-68-5	S-12497A4-1ML	1ML	Octachloro-hexahydro-methano-1H-indene	57-74-9	N-11425-250MG	250MG
Monuron solution	150-68-5	S-12497U1-1ML	1ML	Octachloro-hexahydro-methano-1H-indene Solution	57-74-9	S-11425J4-1ML	1ML
Monuron TCA	140-41-0	N-12498-100MG	100MG	Octachloro-hexahydro-methano-1H-indene Solution	57-74-9	S-11425J4-5ML	5ML
Mospilan(TM)	135410-20-7	N-11009-100MG	100MG	Octachloro-hexahydro-methano-1H-indene Solution	57-74-9	S-11425M1-1ML	1ML
Mospilan(TM) Solution	135410-20-7	S-11009A1-1ML	1ML	Octamethylpyrophosphoramide	152-16-9	N-12715-50MG	50MG
Mospilan(TM) Solution	135410-20-7	S-11009U1-1ML	1ML	Ochililine	26530-20-1	N-14012-250MG	250MG
MPMC	2425-10-700	N-13274-10MG	10MG	Ochililine solution	26530-20-1	S-14012X1-1ML	1ML
MMSA sesquihydrate	2163-80-6	N-12495-100MG	100MG	2-Octyl-4-isothiazolin-3-one Solution	26530-20-1	S-14012X1-1ML	1ML
MMSA sesquihydrate Solution	2163-80-6	S-12495F1-1ML	1ML	2[Octylthio]ethanol Solution	3547-33-9	S-10378M1-1ML	1ML
MTMC	1129-41-5	N-12479-100MG	100MG	Ofir(TM)	66246-88-6	N-12822-100MG	100MG
MTMC Solution	1129-41-5	S-12479A1-1ML	1ML	Ofir(TM) Solution	66246-88-6	S-12822A1-1ML	1ML
MTMC Solution	1129-41-5	S-12479U1-1ML	1ML	Ofir(TM) Solution	66246-88-6	S-12822U1-1ML	1ML
Muster(TM)	97780-06-8	N-11866-100MG	100MG	Ofunack(TM)	119-12-0	N-13154-250MG	250MG
Muster(TM) Solution	97780-06-8	S-11866A1-1ML	1ML	Ofunack(TM) Solution	119-12-0	S-13154A1-1ML	1ML
Myclobutanil	88671-89-0	N-13261-100MG	100MG	Ofunack(TM) Solution	119-12-0	S-13154U1-1ML	1ML
Myclobutanil solution	88671-89-0	S-13261A1-1ML	1ML	Ofurace solution	58810-48-3	S-12988K10-1ML	1ML
Myclobutanil solution	88671-89-0	S-13261U1-1ML	1ML	Ohric	24096-53-5	N-11753-500MG	500MG
Na2AsH4	62-74-8	N-13216-250MG	250MG	Ohric Solution	24096-53-5	S-11753M1-1ML	1ML
NAA(TM)	86-87-3	N-10075-1G	1G	Omethoate	1113-02-6	N-12726-100MG	100MG
NAA(TM) Solution	86-87-3	S-10075M1-1ML	1ML	Omethoate solution	1113-02-6	S-12726U1-1ML	1ML
NAA(TM) Solution	86-87-3	S-10075T1-1ML	1ML	Omite(TM)	2312-35-8	N-12727-100MG	100MG
Nabam	142-59-6	N-12639-250MG	250MG	Omite(TM) Solution	2312-35-8	S-12727A1-1ML	1ML
Nabam solution	142-59-6	S-12639F1-1ML	1ML	Omite(TM) Solution	2312-35-8	S-12727U1-1ML	1ML
Nabam solution	142-59-6	S-12639U1-1ML	1ML	Orbencarb	34622-58-7	N-12951-10MG	10MG
Naled	300-76-5	N-12640-250MG	250MG	Ormetoprim	6981-18-6	N-12728-250MG	250MG
Naled solution	300-76-5	S-12640U1-1ML	1ML	Ormetoprim solution	6981-18-6	S-12728A1-1ML	1ML
1-Naphthaleneacetamide	86-86-2	N-10074-1G	1G	Ormetoprim solution	6981-18-6	S-12728U1-1ML	1ML
1-Naphthaleneacetamide solution	86-86-2	S-10074M1-1ML	1ML	Orthobencarb	34622-58-7	N-12951-10MG	10MG
1-Naphthaleneacetamide solution	86-86-2	S-10074U1-1ML	1ML	Orthosulfamuron solution	213464-77-8	S-13944B1-1ML	1ML
1-Naphthaleneacetic acid	86-87-3	N-10075-1G	1G	Oryzalin	19044-88-3	N-12729-100MG	100MG
1-Naphthaleneacetic acid solution	86-87-3	S-10075M1-1ML	1ML	Oryzalin solution	19044-88-3	S-12729A4-1ML	1ML
1-Naphthaleneacetic acid solution	86-87-3	S-10075T1-1ML	1ML	Oryzemat(TM)	27605-76-1	N-13092-500MG	500MG
1-Naphthaleneacetic acid, ethyl ester	2122-70-5	N-11914-250MG	250MG	Oryzemat(TM) Solution	27605-76-1	S-13092M1-1ML	1ML
1-Naphthaleneacetic acid, ethyl ester Solution	2122-70-5	S-11914A1-1ML	1ML	Outlook(TM)	163515-14-8	N-11756-100MG	100MG
1-Naphthaleneacetic acid, ethyl ester Solution	2122-70-5	S-11914T1-1ML	1ML	Outlook(TM) Solution	163515-14-8	S-11756A1-1ML	1ML
1-8-Naphthalic anhydride	81-84-5	N-10012-500MG	500MG	Outlook(TM) Solution	163515-14-8	S-11756U1-1ML	1ML
1-8-Naphthalic anhydride solution	81-84-5	S-10012A1-1ML	1ML	Overture(TM)	179101-81-6	N-13153-50MG	50MG
1-8-Naphthalic anhydride solution	81-84-5	S-10012U1-1ML	1ML	Ovex	80-33-1	N-12730-250MG	250MG
Naphthenic acid. copper(2+) salt	1338-02-9	N-11501-250MG	250MG	Ovex solution	80-33-1	S-12730A1-1ML	1ML
2-Naphthoxyacetic acid	120-23-0	N-10437-250MG	250MG	Ovex solution	80-33-1	S-12730J1-1ML	1ML
2-Naphthoxyacetic acid solution	120-23-0	S-10437M1-1ML	1ML	Oxadialryl	39807-15-3	N-14011-100MG	100MG
2-Naphthoxyacetic acid solution	120-23-0	S-10437T1-1ML	1ML	Oxadialryl solution	39807-15-3	S-14011A1-1ML	1ML
2-[a-Naphthoxy]-N,N-diethylpropionamide	15299-99-7	N-11585-250MG	250MG	Oxadiazon	19666-30-9	N-12731-250MG	250MG
2-[a-Naphthoxy]-N,N-diethylpropionamide Solution	15299-99-7	S-11585A1-1ML	1ML	Oxadiazon solution	19666-30-9	S-12731A1-1ML	1ML
2-[a-Naphthoxy]-N,N-diethylpropionamide Solution	15299-99-7	S-11585T1-1ML	1ML	Oxadiazon solution	19666-30-9	S-12731U1-1ML	1ML
2-[a-Naphthoxy]-N,N-diethylpropionamide Solution	15299-99-7	S-11585T1-5ML	5ML	Oxadiazon-hydroxy	39807-19-7	MET-12731A-50MG	50MG
N-1-Naphthylacetamide solution	575-36-0	S-12506M1-1ML	1ML	Oxadiazon-hydroxy Solution	39807-19-7	MET-12731AU10-1ML	1ML
N-1-Naphthylacetamide	575-36-0	N-12506-500MG	500MG	Oxadixyl	77732-09-3	N-12732-100MG	100MG
N-1-Naphthylphthalamic acid	132-66-1	N-12507-250MG	250MG	Oxadixyl solution	77732-09-3	S-12732A1-1ML	1ML
N-1-Naphthylphthalamic acid solution	132-66-1	S-12507A1-1ML	1ML	Oxadixyl solution	77732-09-3	S-12732U1-1ML	1ML
N-1-Naphthylphthalamic acid solution	132-66-1	S-12507T1-1ML	1ML	Oxamyl solution	23135-22-0	S-12735A1-1ML	1ML
N-(1-Naphthyl)phthalimide	5333-99-3	MET-12507A-250MG	250MG	Oxamyl oxime	30558-43-1	MET-12735-10MG	10MG
N-alpha-Naphthyl phthalimide	5333-99-3	MET-12507A-250MG	250MG	Oxolinic acid	14698-29-4	N-12738-250MG	250MG
a-Naphthylthiourea	86-88-4	N-10983-250MG	250MG	Oxolinic acid solution	14698-29-4	S-12738A1-1ML	1ML
a-Naphthylthiourea solution	86-88-4	S-10983A1-1ML	1ML	Oxymyl	23135-22-0	N-12740-100MG	100MG
a-Naphthylthiourea solution	86-88-4	S-10983U1-1ML	1ML	Oxycarboxin	5259-88-1	S-12740A1-1ML	1ML
Napropamide	15299-99-7	N-11585-250MG	250MG	Oxycarboxin solution	5259-88-1	S-12740U1-1ML	1ML
Napropamide solution	15299-99-7	S-11585A1-1ML	1ML	Oxychlorodane Solution	27304-13-8	MET-11425AM1-1ML	1ML
Napropamide solution	15299-99-7	S-11585T1-1ML	1ML	Oxy-chlordane solution	33386-84-4	S-14054J1-1ML	1ML
Napropamide solution	15299-99-7	S-11585T1-5ML	5ML	Oxydemeton-methyl	301-12-2	N-12741-50MG	50MG
Naptalam	132-66-1	N-12507-250MG	250MG	Oxyfluorfen	42874-03-3	N-12742-250MG	250MG
Naptalam Solution	132-66-1	S-12507A1-1ML	1ML	Oxyfluorfen solution	42874-03-3	S-12742M1-1ML	1ML
Naptalam Solution	132-66-1	S-12507T1-1ML	1ML	Oxytetracycline hydrochloride	2058-46-0	N-12743-250MG	250MG
Neburon	555-37-3	N-12647-250MG	250MG	Paclobutrazol	76738-62-0	N-12811-500MG	500MG
Neburon solution	555-37-3	S-12647A4-1ML	1ML	Paclobutrazol solution	76738-62-0	S-12811A1-1ML	1ML
Neburon solution	555-37-3	S-12647A4-5ML	5ML	Paclobutrazol solution	76738-62-0	S-12811U1-1ML	1ML
Neomycin trisulfate hydrate(Technical)	1405-10-3	N-12649-250MG	250MG	Paraaxon	311-45-5	N-12816-100MG	100MG
Neomycin trisulfate salt hydrate	1405-10-3	N-12649-250MG	250MG	Paraaxon solution	311-45-5	S-12816A1-1ML	1ML
Nevirol(TM)	4727-29-1	N-13043-250MG	250MG	Paraaxon solution	311-45-5	S-12816U1-1ML	1ML
Nevirol(TM) Solution	4727-29-1	S-13043A1-1ML	1ML	Paraquat dichloride solution	1910-42-5	S-12818F1-1ML	1ML
Nevirol(TM) Solution	4727-29-1	S-13043T1-1ML	1ML	Paraquat dichloride hydrate	1910-42-5	N-12818-500MG	500MG
Nicarbazine	330-95-0	N-12651-100MG	100MG	Paraquat dichloride hydrate solution	1910-42-5	S-12818F17-1ML	1ML
Niclosamide	50-65-7	N-12653-250MG	250MG	Paraquat dichloride hydrate solution	1910-42-5	S-12818F17-5ML	5ML
Nicobifen	188425-85-6	N-11328-100MG	100MG	Parathion	56-38-2	N-12819-100MG	100MG
Nicobifen Solution	188425-85-6	S-11328A1-1ML	1ML	Parathion (TM) (diethyl-d10) Solution	350820-04-1	S-FD98S-1.2ML	1.2ML
Nicosulfuron	111991-09-4	N-10998-100MG	100MG	Parathion solution	56-38-2	S-12819A1-1ML	1ML
(S)-Nicotine	54-11-5	N-12655-1G	1G	Parathion solution	56-38-2	S-12819U1-1ML	1ML
Nicotine solution	54-11-5	S-12655M1-1ML	1ML				

Description	CAS#	Part #	Size	Description	CAS#	Part #	Size
Parathion O-analog	311-45-5	N-12816-100MG	100MG	Phenylmercuric oleate	104-60-9	N-13030-100MG	100MG
Parathion O-analog Solution	311-45-5	S-12816A1-1ML	1ML	Phenylmercuric oleate Solution	104-60-9	S-13030A1-1ML	1ML
Parathion O-analog Solution	311-45-5	S-12816U1-1ML	1ML	Phenylmercuric oleate solution	104-60-9	S-13030A1-1ML	1ML
Parathion (TM) (diethyl-d10)		N-FD998-E-0.01G	0.01G	Phenylmercuric oleate Solution	104-60-9	S-13030U1-1ML	1ML
PCNB	82-68-8	N-12830-250MG	250MG	Phenylmercuric oleate solution	104-60-9	S-13030U1-1ML	1ML
PCNB Solution	82-68-8	S-12830H5-1ML	1ML	Phenylmercuric oleate(Technical)	104-60-9	N-13030-100MG	100MG
PCNB Solution	82-68-8	S-12830H5-5ML	5ML	Phenyl mercuric phthalate	84-70-8	N-13016-100MG	100MG
PCNB Solution	82-68-8	S-12830K4-1ML	1ML	Phenyl mercuric phthalate	84-70-8	N-13016-100MG	100MG
PCNB Solution	82-68-8	S-12830K4-5ML	5ML	Phenyl mercuric phthalate Solution	84-70-8	S-13016A1-1ML	1ML
PCNB Solution	82-68-8	S-12830K7-1ML	1ML	Phenyl mercuric phthalate solution	84-70-8	S-13016A1-1ML	1ML
PCNB Solution	82-68-8	S-12830K7-5ML	5ML	Phenyl mercuric phthalate Solution	84-70-8	S-13016U1-1ML	1ML
PCNB Solution	82-68-8	S-12830M1-1ML	1ML	Phenyl mercuric phthalate solution	84-70-8	S-13016U1-1ML	1ML
PCNB Solution	82-68-8	S-12830M1-5ML	5ML	Phenylmercury borate(Technical)	102-98-7	N-13031-100MG	100MG
Pebulate Solution	1114-71-2	S-13192A1-1ML	1ML	Phenylmercury nitrate(Technical)	55-68-5	N-13032-100MG	100MG
Pebulate Solution	1114-71-2	S-13192T1-1ML	1ML	Phenylmercury salicylate solution	28086-13-7	S-13033A1-1ML	1ML
Pegasus(TM)	80060-09-9	N-11612-250MG	250MG	Phenylmercury salicylate solution	28086-13-7	S-13033U1-1ML	1ML
Pegasus(TM) Solution	80060-09-9	S-11612A1-1ML	1ML	Phenylmercury salicylate(Technical)	28086-13-7	N-13033-250MG	250MG
Pegasus(TM) Solution	80060-09-9	S-11612U1-1ML	1ML	5-(Phenylmethyl)-3-furanmethanol	20416-09-5	MET-13176A-100MG	100MG
Penconazole	66246-88-6	N-12822-100MG	100MG	Phenylthiocarbamide	103-85-5	N-13020-250MG	250MG
Penconazole solution	66246-88-6	S-12822A1-1ML	1ML	Phenylthiocarbamide Solution	103-85-5	S-13020A1-1ML	1ML
Penconazole solution	66246-88-6	S-12822U1-1ML	1ML	Phenylthiocarbamide Solution	103-85-5	S-13020U1-1ML	1ML
Pencycuron	66063-05-6	N-12488-100MG	100MG	Phenyl-2-thiourea	103-85-5	N-13020-250MG	250MG
Pencycuron solution	66063-05-6	S-12488A1-1ML	1ML	Phenyl-2-thiourea Solution	103-85-5	S-13020A1-1ML	1ML
Pencycuron solution	66063-05-6	S-12488U1-1ML	1ML	Phenyl-2-thiourea Solution	103-85-5	S-13020U1-1ML	1ML
Pendimethalin	40487-42-1	N-13140-100MG	100MG	Phorate	298-02-2	N-13035-100MG	100MG
Pendimethalin solution	40487-42-1	S-13140A1-1ML	1ML	Phorate solution	298-02-2	S-13035A1-1ML	1ML
Pendimethalin solution	40487-42-1	S-13140K4-1ML	1ML	Phorate solution	298-02-2	S-13035J1-1ML	1ML
Penflufen	494793-67-8	N-14035-100MG	100MG	Phorate sulfone	2588-04-700	MET-13035B-100MG	100MG
Penflufen solution	494793-67-8	S-14035A1-1ML	1ML	Phorate sulfone Solution	2588-04-700	MET-13035B1-1ML	1ML
Penoxsulam	219714-96-2	N-12823-50MG	50MG	Phorate sulfone Solution	2588-04-700	MET-13035C-100MG	100MG
Penta(TM)	87-86-5	N-12831-500MG	500MG	Phorate sulfone Solution	2588-04-700	MET-13035AU1-1ML	1ML
Penta(TM) Solution	87-86-5	S-12831M1-1ML	1ML	Phosalone	2600-69-3	N-13036-100MG	100MG
Penta(TM) Solution	87-86-5	S-12831M1-5ML	5ML	Phosalone solution	2310-17-0	N-13036M1-1ML	1ML
Pentachloroaniline	527-20-8	MET-12383A-100MG	100MG	Phosdrin	7786-34-7	N-13037-250MG	250MG
Pentachloroanisole	1825-21-4	N-12826-100MG	100MG	Phosdrin Solution	7786-34-7	S-13037A1-1ML	1ML
Pentachloroanisole solution	1825-21-4	S-12826M1-1ML	1ML	Phosdrin Solution	7786-34-7	S-13037U1-1ML	1ML
2,3,4,5,6-Pentachlorobenzamine	527-20-8	MET-12383A-100MG	100MG	Phosdrin Solution	7786-34-7	S-13037U1-5ML	5ML
Pentachlorobenzene	133-49-3	MET-12159B-100MG	100MG	trans-Phosdrin	338-45-4	N-12906-50MG	50MG
Pentachloronitrobenzene	82-68-8	N-12830-250MG	250MG	Phosfolan	947-02-4	N-13038-100MG	100MG
Pentachloronitrobenzene solution	82-68-8	S-12830H5-1ML	1ML	Phosfolan solution	947-02-4	S-13038A1-1ML	1ML
Pentachloronitrobenzene solution	82-68-8	S-12830H5-5ML	5ML	Phosfolan solution	947-02-4	S-13038U1-1ML	1ML
Pentachloronitrobenzene solution	82-68-8	S-12830K4-1ML	1ML	Phosmat (dimethyl-d6) (TM) Solution		S-FD2129S-1.2ML	1.2ML
Pentachloronitrobenzene solution	82-68-8	S-12830K4-5ML	5ML	Phosmet	732-11-6	N-12207-100MG	100MG
Pentachloronitrobenzene solution	82-68-8	S-12830K7-1ML	1ML	Phosmet solution	732-11-6	S-12207A1-1ML	1ML
Pentachloronitrobenzene solution	82-68-8	S-12830K7-5ML	5ML	Phosmet solution	732-11-6	S-12207J1-1ML	1ML
Pentachloronitrobenzene solution	82-68-8	S-12830M1-1ML	1ML	Phosmet solution	732-11-6	S-12207J1-5ML	5ML
Pentachloronitrobenzene solution	82-68-8	S-12830M1-5ML	5ML	Phosmet-oxon Solution	3735-33-9	MET-12207AU1-1ML	1ML
Pentachlorophenol	87-86-5	N-12831-500MG	500MG	Phosphamidon	13171-21-6	N-13039-250MG	250MG
Pentachlorophenol solution	87-86-5	S-12831M1-1ML	1ML	Phosphamidon solution	13171-21-6	S-13039A1-1ML	1ML
Pentachlorophenol solution	87-86-5	S-12831M1-5ML	5ML	Phosphamidon solution	13171-21-6	S-13039U1-1ML	1ML
Pentachlorophenyl acetate	1441-02-7	N-13824-100MG	100MG	Phosphamidon solution	13171-21-6	S-13039U1-5ML	5ML
Pentachlorothiophenol	133-49-3	MET-12159B-100MG	100MG	Phosphan	115-78-6	N-13040-50MG	50MG
Pentachlor	2307-68-8	N-12844-100MG	100MG	Phosphan solution	115-78-6	S-13040M1-1ML	1ML
Pentachlor solution	2307-68-8	S-12844A1-1ML	1ML	D,L-Phosphinothricin	77182-82-2	N-12111-100MG	100MG
Pentachlor solution	2307-68-8	S-12844U1-1ML	1ML	N-(Phosphonomethyl)glycine	1071-83-6	N-12133-250MG	250MG
Penthiopyrad	183675-82-3	N-17917-50MG	50MG	Phosphoric acid diethyl ester	598-02-7	MET-11621C-100MG	100MG
Perchlorobenzene	118-74-1	N-12159-250MG	250MG	Phosvel	21609-90-5	N-12315-100MG	100MG
Perchlorobenzene Solution	118-74-1	S-12159K4-1ML	1ML	Phosvel Solution	21609-90-5	S-12315A1-1ML	1ML
Perchlorobenzene Solution	118-74-1	S-12159M1-1ML	1ML	Phosvel Solution	21609-90-5	S-12315U1-1ML	1ML
Perchlorobenzene Solution	118-74-1	S-12159M1-5ML	5ML	Phoxim	14816-18-3	N-13041-100MG	100MG
Perchloromethyl mercaptan	594-42-3	N-12845-250MG	250MG	Phoxim solution	14816-18-3	S-13041A1-1ML	1ML
Perchloromethyl mercaptan solution	594-42-3	S-12845A1-1ML	1ML	Phoxim solution	14816-18-3	S-13041U1-1ML	1ML
Perchloromethyl mercaptan solution	594-42-3	S-12845U1-1ML	1ML	Phthalamic acid	88-97-1	MET-12007A-100MG	100MG
Perfluorooctanesulfonic acid solution	1763-23-1	S-14002M1-1ML	1ML	Phthalanilic acid	4727-29-1	N-13043-250MG	250MG
Permethrin	52645-53-1	N-12848-250MG	250MG	Phthalanilic acid solution	4727-29-1	S-13043A1-1ML	1ML
cis-Permethrin (phenoxy-13C6) Solution		S-FC2216AS-1.2ML	1.2ML	Phthalanilic acid solution	4727-29-1	S-13043T1-1ML	1ML
trans-Permethrin (phenoxy-13C6) Solution		S-FC2216BS-1.2ML	1.2ML	Phthalic acid - dimethyl ester	131-11-3	N-11770-500MG	500MG
Permethrin solution	52645-53-1	S-12848A1-1ML	1ML	Phthalic acid - dimethyl ester Solution	131-11-3	S-11770H1-1ML	1ML
Permethrin solution	52645-53-1	S-12848T1-1ML	1ML	Phthalic acid - dimethyl ester Solution	131-11-3	S-11770H1-5ML	5ML
cis-Permethrin	61949-76-6	N-11483-50MG	50MG	Phthalic acid - dimethyl ester Solution	131-11-3	S-11770M1-1ML	1ML
cis-Permethrin solution	61949-76-6	S-11483A1-1ML	1ML	Phthalic acid - dimethyl ester Solution	131-11-3	S-11770M1-5ML	5ML
cis-Permethrin solution	61949-76-6	S-11483T1-1ML	1ML	Phthalide	27355-22-2	N-13085-10MG	10MG
trans-Permethrin	61949-77-7	N-13620-50MG	50MG	Picardin	119515-38-7	N-12197-100MG	100MG
trans-Permethrin solution	61949-77-7	S-13620A1-1ML	1ML	Picardin Solution	119515-38-7	S-12197A1-1ML	1ML
trans-Permethrin solution	61949-77-7	S-13620T1-1ML	1ML	Picloram	1918-02-100	N-13050-250MG	250MG
Permethrinic acid	55701-05-8	MET-12848A-10MG	10MG	Picloram solution	1918-02-100	S-13050A1-1ML	1ML
Permethrinic acid	55701-05-8	MET-12848AM1-1ML	1ML	Picloram solution	1918-02-100	S-13050A4-1ML	1ML
PFO Solution	1763-23-1	S-14002M1-1ML	1ML	Picloram solution	1918-02-100	S-13050A4-5ML	5ML
Phenmedipham	13684-63-4	N-12859-250MG	250MG	Picloram solution	1918-02-100	S-13050B1-1ML	1ML
Phenmedipham solution	13684-63-4	S-12859A1-1ML	1ML	Picloram solution	1918-02-100	S-13050B1-5ML	5ML
Phenmedipham solution	13684-63-4	S-12859U1-1ML	1ML	Picloram methyl ester	14143-55-6	N-13051-100MG	100MG
d-(cis-trans)-Phenothrin	26002-80-2	N-11554-50MG	50MG	Picloram methyl ester solution	14143-55-6	S-13051M1-1ML	1ML
d-(cis-trans)-Phenothrin solution	26002-80-2	S-11554A4-1ML	1ML	Piclorafen solution	137641-05-5	S-13982M1-1ML	1ML
d-(cis-trans)-Phenothrin solution	26002-80-2	S-11554U1-1ML	1ML	Picoxystrobin	117428-22-5	N-13927-100MG	100MG
Phenoxyacetic acid	122-59-8	N-13004-250MG	250MG	Picoxystrobin solution	117428-22-5	S-13927A1-1ML	1ML
Phenoxyacetic acid solution	122-59-8	S-13004M1-1ML	1ML	Pindone	83-26-1	N-10512-250MG	250MG
Phenoxyacetic acid solution	122-59-8	S-13004T1-1ML	1ML	Pindone Solution	83-26-1	S-10512A1-1ML	1ML
3-Phenoxybenzoic acid	3739-38-6	MET-12848C-100MG	100MG	Pindone Solution	83-26-1	S-10512U1-1ML	1ML
m-Phenoxybenzoic acid	3739-38-6	MET-12848C-100MG	100MG	Pinoxaden	243973-20-8	N-13057-100MG	100MG
m-Phenoxybenzyl alcohol	13826-35-2	MET-12848B-100MG	100MG	Piperalin	3478-94-2	N-13058-100MG	100MG
2-Phenoxypropionic acid	940-31-8	N-10507-500MG	500MG	Piperalin solution	3478-94-2	S-13058A1-1ML	1ML
2-Phenoxypropionic acid solution	940-31-8	S-10507M1-1ML	1ML	Piperalin solution	3478-94-2	S-13058U1-1ML	1ML
2-Phenoxypropionic acid solution	940-31-8	S-10507T1-1ML	1ML	Piperonyl butoxide	51-03-6	N-13061-100MG	100MG
Phenthoate	2597-03-700	N-13006-100MG	100MG	Piperonyl butoxide solution	51-03-6	S-13061A1-1ML	1ML
Phenthoate solution	2597-03-700	S-13006A1-1ML	1ML	Piperonyl butoxide solution	51-03-6	S-13061K1-1ML	1ML
Phenthoate solution	2597-03-700	S-13006U1-1ML	1ML	Piperophos	24151-93-7	N-12989-250MG	250MG
Phentinoacetate	900-95-8	N-13007-250MG	250MG	Pipron	3478-94-2	N-13058-100MG	100MG
Phentinoacetate solution	900-95-8	S-13007A1-1ML	1ML	Pipron Solution	3478-94-2	S-13058A1-1ML	1ML
Phentinoacetate solution	900-95-8	S-13007U1-1ML	1ML	Pipron Solution	3478-94-2	S-13058U1-1ML	1ML
Phenyl benzene	92-52-4	N-11207-500MG	500MG	Pirimicarb	23103-98-2	N-13062-250MG	250MG
Phenyl benzene Solution	92-52-4	S-11207J1-1ML	1ML	Pirimicarb solution	23103-98-2	S-13062A1-1ML	1ML
Phenyl benzene Solution	92-52-4	S-11207J1-5ML	5ML	Pirimicarb solution	23103-98-2	S-13062U1-1ML	1ML
Phenyl benzene Solution	92-52-4	S-11207M4-1ML	1ML	Pirimicarb-desmethyl	30614-22-3	MET-13801A-50MG	50MG
Phenyl benzene Solution	92-52-4	S-11207M4-5ML	5ML	Pirimiphos-ethyl	23505-41-1	N-13063-250MG	250MG
4-(2-phenylethoxy)-quinazoline	124427-60-7	N-12908-10MG	10MG	Pirimiphos-ethyl solution	23505-41-1	S-13063U1-1ML	1ML
Phenylm (TM)	122-59-8	N-13004-250MG	250MG	Pirimiphos-ethyl-oxon	36378-61-7	MET-13063A-250MG	250MG
Phenylm (TM) Solution	122-59-8	S-13004M1-1ML	1ML	Pirimiphos-methyl	29232-93-7	N-13064-250MG	250MG
Phenylm (TM) Solution	122-59-8	S-13004T1-1ML	1ML	Pirimiphos-methyl solution	29232-93-7	S-13064A1-1ML	1ML
Phenyl mercuric acetate(Technical)	62-38-4	N-13011-250MG	250MG	Pirimiphos-methyl solution	29232-93-7	S-13064U1-1ML	1ML
Phenyl mercuric benzoate	94-43-9	N-13012-100MG	100MG	Pirimiphos-methyl-N-desethyl	67018-59-1	MET-13064D-100MG	100MG
Phenyl mercuric benzoate	94-43-9	N-13012-100MG	100MG	Pirimiphos-methyl-oxon	64709-45-1	MET-13064B-250MG	250MG
Phenyl mercuric chloride(Technical)	100-56-1	N-13013-250MG	250MG	2-Pivaloyl-1,3-indandione	83-26-1	N-10512-250MG	250MG
Phenyl mercuric chloride(Technical)	100-56-1	N-13013-250MG	250MG	2-Pivaloyl-1,3-indandione solution	83-26-1	S-10512A1-1ML	1ML
Phenyl mercuric lactate	122-64-5	N-13015-250MG	250MG	2-Pivaloyl-1,3-indandione solution	83-26-1	S-10512U1-1ML	1ML

Description	CAS#	Part #	Size	Description	CAS#	Part #	Size
Pivot Solution	213464-77-8	S-13944B1-1ML	1ML	Propyzamide	23950-58-5	N-13136-250MG	250MG
Plondrel(TM)	5131-24-8	N-11820-100MG	100MG	Propyzamide solution	23950-58-5	S-13136M1-1ML	1ML
Plondrel(TM) Solution	5131-24-8	S-11820A1-1ML	1ML	Propyzamide solution	23950-58-5	S-13136M1-5ML	5ML
Plondrel(TM) Solution	5131-24-8	S-11820J1-1ML	1ML	Proquiquad	189278-12-4	N-14080-50MG	50MG
Plondrel(TM) Solution	5131-24-8	S-11820J1-5ML	5ML	Proslufocarb	52888-80-9	N-13918-10MG	10MG
Plondrel(TM) Solution	5131-24-8	S-11820U1-1ML	1ML	Proslufuron	94125-34-5	N-13137-100MG	100MG
PMA TM	62-38-4	N-13011-250MG	250MG	Proslulfuron solution	94125-34-5	S-13137A1-1ML	1ML
Polaris	2439-99-8	N-12135-100MG	100MG	Protect (TM)	81-84-5	N-10012-500MG	500MG
Polaris Solution	2439-99-8	S-12135F1-1ML	1ML	Protect (TM) Solution	81-84-5	S-10012A1-1ML	1ML
Polaris Solution	2439-99-8	S-12135U1-1ML	1ML	Protect (TM) Solution	81-84-5	S-10012U1-1ML	1ML
Polyram (TM)	9006-42-2	N-13086-250MG	250MG	Prothioconazole	178928-70-6	N-13138-100MG	100MG
Potassium cyanate(Technical)	590-28-3	N-13087-250MG	250MG	Prothioconazole solution	178928-70-6	S-13138B1-1ML	1ML
Potassium isocyanate	590-28-3	N-13087-250MG	250MG	Prothiophos	34643-46-4	N-13577-100MG	100MG
Prallethrin	23031-36-9	N-13088-10MG	10MG	Prothiophos solution	34643-46-4	S-13577A1-1ML	1ML
Prallethrin solution	23031-36-9	S-13088A1-1ML	1ML	Prothiophos solution	34643-46-4	S-13577U1-1ML	1ML
Prallethrin solution	23031-36-9	S-13088K1-1ML	1ML	Prowl	40487-42-1	N-13140-100MG	100MG
Pretilachlor	51218-49-6	N-13089-250MG	250MG	Prowl Solution	40487-42-1	S-13140A1-1ML	1ML
Pretilachlor solution	51218-49-6	S-13089A1-1ML	1ML	Prowl Solution	40487-42-1	S-13140K4-1ML	1ML
Pretilachlor solution	51218-49-6	S-13089U1-1ML	1ML	Pymetrozine	123312-89-0	N-13141-250MG	250MG
Primisulfuron-methyl	86209-51-0	N-13090-100MG	100MG	Pymetrozine solution	123312-89-0	S-13141A1-1ML	1ML
Primisulfuron-methyl solution	86209-51-0	S-13090A1-1ML	1ML	Pymetrozine solution	123312-89-0	S-13141U1-1ML	1ML
Primo(TM)	80844-07-1	N-11941-100MG	100MG	Pyraclafos	89784-60-1	S-13142-10MG	10MG
Primo(TM) Solution	80844-07-1	S-11941A1-1ML	1ML	Pyraclafos solution	89784-60-1	S-13142A1-1ML	1ML
Primo(TM) Solution	80844-07-1	S-11941U1-1ML	1ML	Pyraclostrobin	175013-18-0	N-13143-100MG	100MG
Probenazole	27605-76-1	N-13092-500MG	500MG	Pyraclostrobin solution	175013-18-0	S-13143A1-1ML	1ML
Probenazole solution	27605-76-1	S-13092M1-1ML	1ML	Pyraflufen	129630-17-7	MET-13144A-10MG	10MG
Prochloraz	67747-09-5	N-13093-250MG	250MG	Pyraflufen (acid)	129630-17-7	MET-13144A-10MG	10MG
Prochloraz solution	67747-09-5	S-13093A1-1ML	1ML	Pyraflufen-ethyl	129630-19-9	N-13144-100MG	100MG
Prochloraz solution	67747-09-5	S-13093U1-1ML	1ML	Pyranica(TM)	119168-77-3	N-13502-100MG	100MG
Prociazine	32889-48-8	N-13094-100MG	100MG	Pyrasulfotole	365400-11-9	N-14001-100MG	100MG
Prociazine solution	32889-48-8	S-13094A1-1ML	1ML	Pyrazolynate	58011-68-0	N-13083-50MG	50MG
Prociazine solution	32889-48-8	S-13094T1-1ML	1ML	Pyrazolynate solution	58011-68-0	S-13083A1-1ML	1ML
Procymidone	32809-16-8	N-13095-250MG	250MG	Pyrazon	1698-60-8	N-13145-250MG	250MG
Procymidone solution	32809-16-8	S-13095A1-1ML	1ML	Pyrazon solution	1698-60-8	S-13145M1-1ML	1ML
Procymidone solution	32809-16-8	S-13095U1-1ML	1ML	Pyrazophos	13457-18-6	N-13146-100MG	100MG
Prodiamine	29091-21-2	N-13096-100MG	100MG	Pyrazophos solution	13457-18-6	S-13146A1-1ML	1ML
Prodiamine solution	29091-21-2	S-13096A1-1ML	1ML	Pyrazophos solution	13457-18-6	S-13146U1-1ML	1ML
Prodiamine solution	29091-21-2	S-13096K1-1ML	1ML	Pyrazosulfuron-ethyl	93697-74-6	N-13147-100MG	100MG
Profenofos	41198-08-7	N-13097-250MG	250MG	Pyrazosulfuron-ethyl solution	93697-74-6	S-13147A1-1ML	1ML
Profenofos solution	41198-08-7	S-13097A1-1ML	1ML	Pyrazoxyfen	71561-11-0	N-13148-100MG	100MG
Profenofos solution	41198-08-7	S-13097J1-1ML	1ML	Pyrazoxyfen solution	71561-11-0	S-13148M1-1ML	1ML
Profluralin	26399-36-0	N-13098-250MG	250MG	Pyrethrum	8003-34-7	N-13151-100MG	100MG
Profluralin solution	26399-36-0	S-13098A1-1ML	1ML	Pyrethrum solution	8003-34-7	S-13151A1-1ML	1ML
Profluralin solution	26399-36-0	S-13098J4-1ML	1ML	Pyrethrum solution	8003-34-7	S-13151T1-1ML	1ML
Profluralin solution	26399-36-0	S-13098J4-5ML	5ML	Pyribenzoxim	168088-61-7	N-14070-100MG	100MG
Prohexadione	88805-35-0	N-13099-100MG	100MG	Pyridaben	96489-71-3	N-13152-25MG	25MG
Prohexadione-Calcium	127277-53-6	N-13100-100MG	100MG	Pyridaben solution	96489-71-3	S-13152M1-1ML	1ML
Promecarb	2631-37-0	N-13101-100MG	100MG	Pyridalyl	179101-81-6	N-13153-50MG	50MG
Promecarb solution	2631-37-0	S-13101M4-1ML	1ML	Pyridaphenthion	119-12-0	N-13154-250MG	250MG
Prometon	1610-18-0	N-13102-250MG	250MG	Pyridaphenthion solution	119-12-0	S-13154A1-1ML	1ML
Prometon solution	1610-18-0	S-13102A1-1ML	1ML	Pyridaphenthion solution	119-12-0	S-13154U1-1ML	1ML
Prometon solution	1610-18-0	S-13102T1-1ML	1ML	Pyridate	55512-33-9	N-13155-500MG	500MG
Prometon solution	1610-18-0	S-13102T1-5ML	5ML	Pyridate solution	55512-33-9	S-13155A1-1ML	1ML
Prometryne	7287-19-6	N-13103-250MG	250MG	Pyridate solution	55512-33-9	S-13155T1-1ML	1ML
Prometryne solution	7287-19-6	S-13103A1-1ML	1ML	Pyrifeno	88283-41-4	N-12990-10MG	10MG
Prometryne solution	7287-19-6	S-13103T1-1ML	1ML	Pyrifeno solution	88283-41-4	S-12990M1-1ML	1ML
Prometryne solution	7287-19-6	S-13103T1-5ML	5ML	Pyriquinazon	337458-27-2	N-13158-100MG	100MG
Propachlor	1918-16-7	S-13104-250MG	250MG	Pyrimethanil	53112-28-0	N-13159-100MG	100MG
Propachlor solution	1918-16-7	S-13104A1-1ML	1ML	Pyrimethanil solution	53112-28-0	S-13159A1-1ML	1ML
Propachlor solution	1918-16-7	S-13104T1-1ML	1ML	Pyrimethanil solution	53112-28-0	S-13159U1-1ML	1ML
Propachlor solution	1918-16-7	S-13104T1-5ML	5ML	Pyrimidifen solution	105779-78-0	S-13963A1-1ML	1ML
Propachlor ESA sodium salt	947601-88-9	MET-13104A-25MG	25MG	Pyriminil	53558-25-1	N-13738-100MG	100MG
Propachlor ESA sodium salt Solution	947601-88-9	MET-13104AM1-1ML	1ML	Pyriminil Solution	53558-25-1	S-13738A1-1ML	1ML
Propachlor OA	70628-36-3	MET-13104B-10MG	10MG	Pyriminil Solution	53558-25-1	S-13738U1-1ML	1ML
Propachlor OA Solution	70628-36-3	MET-13104BM1-1ML	1ML	Pyrinuron	53558-25-1	N-13738-100MG	100MG
Propamocarb	24579-73-5	N-13105-250MG	250MG	Pyrinuron solution	53558-25-1	S-13738A1-1ML	1ML
Propamocarb HCl	25606-41-1	N-13106-250MG	250MG	Pyrinuron solution	53558-25-1	S-13738U1-1ML	1ML
Propamocarb HCl solution	25606-41-1	S-13106F1-1ML	1ML	Pyriphenox	88283-41-4	N-12990-10MG	10MG
Propamocarb Hydrochloride	25606-41-1	N-13106F-250MG	250MG	Pyriphenox solution	88283-41-4	S-12990M1-1ML	1ML
Propamocarb Hydrochloride Solution	25606-41-1	S-13106F1-1ML	1ML	Pyriproxyfen	95737-68-1	N-13160-100MG	100MG
Propanil	709-98-8	N-10800-500MG	500MG	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propanil	709-98-8	N-10800-500MG	500MG	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propanil Solution	709-98-8	S-10800A1-1ML	1ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propanil solution	709-98-8	S-10800A1-1ML	1ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propanil solution	709-98-8	S-10800Z4-1ML	1ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propanil solution	709-98-8	S-10800Z4-1ML	1ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propanil solution	709-98-8	S-10800Z4-5ML	5ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propanil solution	709-98-8	S-10800Z4-5ML	5ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propaphos solution	7292-16-2	S-13967A1-1ML	1ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propaquizafop	111479-05-1	N-14007-250MG	250MG	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propaquizafop solution	111479-05-1	S-14007A1-1ML	1ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propargite	2312-35-8	N-12727-100MG	100MG	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propargite solution	2312-35-8	S-12727A1-1ML	1ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propargite solution	2312-35-8	S-12727U1-1ML	1ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
2-Propargylsuccinimidymethyl chrysanthemate	72963-72-5	S-14006A1-1ML	1ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propazine	139-40-2	N-13108-100MG	100MG	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propazine (ring-13C3) Solution	139-40-2	S-FC2236S-1.2ML	1.2ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propazine solution	139-40-2	S-13108A1-1ML	1ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propazine solution	139-40-2	S-13108T1-1ML	1ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propazine solution	139-40-2	S-13108T1-5ML	5ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propazine-2-hydroxy Solution	7374-53-0	MET-13108AAA1-1ML	1ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propetamphos	31218-83-4	N-13197-100MG	100MG	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propetamphos solution	31218-83-4	S-13197A1-1ML	1ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propetamphos solution	31218-83-4	S-13197U1-1ML	1ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propham	122-42-9	N-13109-250MG	250MG	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propham solution	122-42-9	S-13109A4-1ML	1ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Prophos	13194-48-4	N-13110-250MG	250MG	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Prophos solution	13194-48-4	S-13110A1-1ML	1ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Prophos solution	13194-48-4	S-13110U1-1ML	1ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Prophos solution	13194-48-4	S-13110U1-5ML	5ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propiconazole	60207-90-1	N-13576-250MG	250MG	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propiconazole solution	60207-90-1	S-13576A1-1ML	1ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propiconazole solution	60207-90-1	S-13576E1-1ML	1ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propineb(Technical)	9016-72-2	N-13111-250MG	250MG	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propisochlor	86763-47-5	N-13823-50MG	50MG	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propoxur	114-26-1	N-11128-250MG	250MG	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propoxur solution	114-26-1	S-11128A1-1ML	1ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propoxycarbazone-sodium	181274-15-7	N-13119-100MG	100MG	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propoxycarbazone-sodium solution	181274-15-7	S-13119A1-1ML	1ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
S-Propyl butylethylthiocarbamate solution	1114-71-2	S-13192A1-1ML	1ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
S-Propyl butylethylthiocarbamate solution	1114-71-2	S-13192T1-1ML	1ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propylene dichloride	78-87-5	N-10156-1G	1G	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propylene dichloride Solution	78-87-5	S-10156M1-1ML	1ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML
Propylene dichloride Solution	78-87-5	S-10156M1-5ML	5ML	Pyriproxyfen solution	95737-68-1	S-13160M1-1ML	1ML

Description	CAS#	Part #	Size	Description	CAS#	Part #	Size
d-trans-Resmethrin Solution	28434-01-7	S-11206T1-1ML	1ML	Spirotetramat	203313-25-1	N-13225-100MG	100MG
Revus™	374726-62-2	N-12354-100MG	100MG	Spirotetramat solution	203313-25-1	S-13225B1-1ML	1ML
Rifit(TM)	51218-49-6	N-13089-250MG	250MG	Spiroxamine	118134-30-8	N-13226-100MG	100MG
Rifit(TM) Solution	51218-49-6	S-13089A1-1ML	1ML	Spiroxamine solution	118134-30-8	S-13226U1-1ML	1ML
Rifit(TM) Solution	51218-49-6	S-13089U1-1ML	1ML	Spotless(TM)	83657-24-3	N-11782-100MG	100MG
Rimon(TM)	116714-46-6	N-12670-100MG	100MG	Spotless(TM) Solution	83657-24-3	S-11782M1-1ML	1ML
Rimon(TM) Solution	116714-46-6	S-12670A1-1ML	1ML	Stablene	90035-08-8	N-11354-250MG	250MG
Rimon(TM) Solution	116714-46-6	S-12670U1-1ML	1ML	Staple(TM)	123343-16-8	N-13161-500MG	500MG
Risulfuron	122931-48-0	N-13180-100MG	100MG	Staple(TM) Solution	123343-16-8	S-13161F1-1ML	1ML
Risolex(TM)	57018-04-9	N-13578-500MG	500MG	STB	132827-25-9	MET-1138B-500MG	500MG
Risolex(TM) Solution	57018-04-9	S-13578A1-1ML	1ML	Storm(TM)	90035-08-8	N-11975-10MG	10MG
Risolex(TM) Solution	57018-04-9	S-13578U1-1ML	1ML	Storm(TM) Solution	90035-08-8	S-11975A1-1ML	1ML
Ronnel	299-84-3	N-11952-100MG	100MG	Streptomycin sulfate	3810-74-0	N-13228-250MG	250MG
Ronnel Solution	299-84-3	S-11952U1-1ML	1ML	Streptomycin sulfate (2:3) salt	3810-74-0	N-13228-250MG	250MG
Ronnel Solution	299-84-3	S-11952U1-5ML	5ML	Stroban (TM)	8001-50-1	N-13229-100MG	100MG
Rootene	86-86-2	N-10074-1G	1G	Stroban (TM) solution	8001-50-1	S-13229A1-1ML	1ML
Rootene Solution	86-86-2	S-10074M1-1ML	1ML	Stroban (TM) solution	8001-50-1	S-13229K4-1ML	1ML
Rootene Solution	86-86-2	S-10074U1-1ML	1ML	Stroby(TM)	143390-89-0	N-12296-250MG	250MG
Rotenone	83-79-4	N-13184-250MG	250MG	Stroby(TM) Solution	143390-89-0	S-12296A1-1ML	1ML
Rotenone solution	83-79-4	S-13184A1-1ML	1ML	Stroby(TM) Solution	143390-89-0	S-12296U1-1ML	1ML
RPA 203328	142994-06-7	MET-12287B-100MG	100MG	Strongarm(TM)	145701-21-9	N-11679-250MG	250MG
Ruelene (TM)	299-86-5	N-13186-50MG	50MG	Strongarm(TM) Solution	145701-21-9	S-11679A1-1ML	1ML
Ruelene (TM) solution	299-86-5	S-13186A1-1ML	1ML	Succinic acid - cadmium salt	141-00-4	N-11391-250MG	250MG
Ruelene (TM) solution	299-86-5	S-13186U1-1ML	1ML	Succinic acid - dibutyl ester	141-03-7	N-11644-500MG	500MG
Rynaxypyr (TM)	500008-45-7	N-11422-100MG	100MG	Succinic acid - dibutyl ester Solution	141-03-7	S-11644A1-1ML	1ML
S 421	127-90-2	N-13919-10MG	10MG	Succinic acid - dibutyl ester Solution	141-03-7	S-11644U1-1ML	1ML
Saflufenacil	372137-35-4	N-13195-100MG	100MG	Sulcotrione	99105-77-8	N-12890-100MG	100MG
Safrotin(TM)	31218-83-4	N-13197-100MG	100MG	Sulfabenzamide	127-71-9	N-13240-250MG	250MG
Safrotin(TM) Solution	31218-83-4	S-13197A1-1ML	1ML	Sulfachloropyridazine	80-32-0	N-13241-250MG	250MG
Safrotin(TM) Solution	31218-83-4	S-13197U1-1ML	1ML	Sulfachloropyridazine	80-32-0	N-13241-250MG	250MG
Salannin	992-20-1	N-12911-5MG	5MG	Sulfadimethoxine	122-11-2	N-13242-250MG	250MG
Salithion Solution	3811-49-2	S-13964A1-1ML	1ML	Sulfadimethoxine solution	122-11-2	S-13242A1-1ML	1ML
Sanmarton	51630-58-1	N-13201-250MG	250MG	Sulfallate	95-06-7	N-13243-250MG	250MG
Sanmarton Solution	51630-58-1	S-13201A4-1ML	1ML	Sulfallate solution	95-06-7	S-13243M1-1ML	1ML
Sanmarton Solution	51630-58-1	S-13201A4-5ML	5ML	Sulfallate solution	95-06-7	S-13243M1-5ML	5ML
Sanmarton Solution	51630-58-1	S-13201U1-1ML	1ML	Sulfamerazine	127-79-7	N-13244-250MG	250MG
Sanmite*	96489-71-3	N-13152-25MG	25MG	Sulfamethazine	57-68-1	N-13245-250MG	250MG
Schradan	152-16-9	N-12715-50MG	50MG	Sulfamethizole	144-82-1	N-13246-250MG	250MG
Sebuthylazin	7286-69-3	N-12950-250MG	250MG	Sulfamethoxazole	723-46-6	N-13247-250MG	250MG
Sebuthylazine-desethyl	37019-18-4	MET-12950A-100MG	100MG	Sulfanilic acid	121-57-3	N-13249-250MG	250MG
Sedaxane	874967-67-6	N-13068-100MG	100MG	Sulfapyridine	144-83-2	N-13250-250MG	250MG
Sedaxane solution	874967-67-6	S-13068M1-1ML	1ML	Sulfaquinoxaline	59-40-5	N-13251-10MG	10MG
Select*	99129-21-2	N-11488-100MG	100MG	Sulfaquinoxaline solution	59-40-5	S-13251M1-1ML	1ML
Sentinel(TM)	94361-06-5	N-11547-50MG	50MG	Sulfathiazole	72-14-0	N-13252-250MG	250MG
Sentinel(TM) Solution	94361-06-5	S-11547A1-1ML	1ML	Sulfentrazone	122836-35-5	N-13253-250MG	250MG
Sentinel(TM) Solution	94361-06-5	S-11547E1-1ML	1ML	Sulfentrazone solution	122836-35-5	S-13253M1-1ML	1ML
Septiphen	120-32-1	N-12676-100MG	100MG	Sulfometuron methyl	74222-97-2	N-13254-100MG	100MG
Septiphen Solution	120-32-1	S-12676M1-1ML	1ML	Sulfometuron methyl solution	74222-97-2	S-13254U1-1ML	1ML
Di-septon	298-03-3	N-11580-100MG	100MG	Sulfosulfuron	141776-32-1	N-13255-100MG	100MG
Di-septon Solution	298-03-3	S-11580U1-1ML	1ML	Sulfosulfuron solution	141776-32-1	S-13255A1-1ML	1ML
Seradix	133-32-4	N-10728-100MG	100MG	Sulfosulfuron solution	141776-32-1	S-13255U1-1ML	1ML
Seradix Solution	133-32-4	S-10728A1-1ML	1ML	Sulfoxaflo	946578-00-3	N-12883-10MG	10MG
Seradix Solution	133-32-4	S-10728T1-1ML	1ML	Sulfoxaflo solution	946578-00-3	S-12883M1-1ML	1ML
Sethoxydim	74051-80-2	N-13210-100MG	100MG	Sulfoxide	120-62-7	N-13257-100MG	100MG
Sethoxydim solution	74051-80-2	S-13210U1-1ML	1ML	Sulfoxide solution	120-62-7	S-13257M1-1ML	1ML
Sethoxydim-desethoxy	113561-03-8	MET-13210A-50MG	50MG	Sulfoxide solution	120-62-7	S-13257M1-5ML	5ML
Sethoxydim-oxazole	413594-06-6	MET-13210B-50MG	50MG	Sulprofos	35400-43-2	N-13258-100MG	100MG
Siduron	1982-49-6	N-13211-100MG	100MG	Sulprofos solution	35400-43-2	S-13258A1-1ML	1ML
Siduron solution	1982-49-6	S-13211A4-1ML	1ML	Sulprofos solution	35400-43-2	S-13258J1-1ML	1ML
Siduron solution	1982-49-6	S-13211A4-5ML	5ML	Sulprofos solution	35400-43-2	S-13258J1-5ML	5ML
Silafluofen	105024-66-6	N-13914-10MG	10MG	Sulprofos sulfone Solution	92642-35-8	MET-13258AU1-1ML	1ML
Silvex	93-72-1	N-13212-250MG	250MG	Sulprofos sulfone Solution	92642-35-8	MET-13258AU1-5ML	5ML
Silvex solution	93-72-1	S-13212A1-1ML	1ML	Sulprofos-sulfoxide	34643-47-5	MET-13258C-100MG	100MG
Silvex solution	93-72-1	S-13212B1-1ML	1ML	Sumilex(TM)	32809-16-8	N-13095-250MG	250MG
Silvex solution	93-72-1	S-13212B1-5ML	5ML	Sumilex(TM) Solution	32809-16-8	S-13095A1-1ML	1ML
Silvex. 2-ethylhexyl ester	7778-43-0	N-13215-250MG	250MG	Sumilex(TM) Solution	32809-16-8	S-13095U1-1ML	1ML
Silvex. 2-ethylhexyl ester Solution	7778-43-0	S-13215F1-1ML	1ML	Sumithrin	26002-80-2	N-11554-50MG	50MG
Silvex. 2-ethylhexyl ester Solution	7778-43-0	S-13215U1-1ML	1ML	Sumithrin Solution	26002-80-2	S-11554A4-1ML	1ML
Silvex. 2-ethylhexyl ester	32534-95-5	N-13213-100MG	100MG	Sumithrin Solution	26002-80-2	S-11554U1-1ML	1ML
Silvex. 2-ethylhexyl ester Solution	32534-95-5	S-13213A1-1ML	1ML	Super-Rozol	28772-56-7	N-13259-100MG	100MG
Silvex isocetyl ester	32534-95-5	N-13213-100MG	100MG	Super-Rozol solution	28772-56-7	S-13259A1-1ML	1ML
Silvex isocetyl ester solution	32534-95-5	S-13213A1-1ML	1ML	Super-Rozol solution	28772-56-7	S-13259U1-1ML	1ML
Silvex isocetyl ester solution	32534-95-5	S-13213T1-1ML	1ML	Super-tin	76-87-9	N-11966-250MG	250MG
Silvex methyl ester	4841-20-7	N-13214-250MG	250MG	Super-tin Solution	76-87-9	S-11966A1-1ML	1ML
Silvex methyl ester solution	4841-20-7	S-13214M1-1ML	1ML	Supracide	950-37-8	N-12397-250MG	250MG
Silvex methyl ester solution	4841-20-7	S-13214M1-5ML	5ML	Supracide Solution	950-37-8	S-12397A1-1ML	1ML
Simazine	122-34-9	N-13800-100MG	100MG	Supracide Solution	950-37-8	N-12397U1-1ML	1ML
Simazine (ring-13C3) Solution		S-FC2073S-1.2ML	1.2ML	Sweep	1918-18-9	N-13260-250MG	250MG
Simazine solution	122-34-9	S-13800B4-1ML	1ML	Sweep solution	1918-18-9	S-13260A4-1ML	1ML
Simazine solution	122-34-9	S-13800B4-5ML	5ML	Sweep solution	1918-18-9	S-13260A4-5ML	5ML
Simazine solution	122-34-9	S-13800M1-1ML	1ML	Sweep solution	1918-18-9	S-13260U1-1ML	1ML
Simazine solution	122-34-9	S-13800M1-5ML	5ML	Systox	8065-48-3	N-11581-100MG	100MG
Simazine hydroxy	2599-11-300	MET-13800B-25MG	25MG	Systox Solution	8065-48-3	S-11581U1-1ML	1ML
Simconazole	14958-90-7	N-13931-25MG	25MG	2,4,5-T	93-76-5	N-10646-250MG	250MG
Simconazole solution	14958-90-7	S-13931A1-1ML	1ML	2,4,5-T (ring-13C6) Solution		S-FC972S-1.2ML	1.2ML
Simeton	673-04-1	N-12960-50MG	50MG	2,4,5-T solution	93-76-5	S-10646A1-1ML	1ML
Simetryn	1014-70-6	N-13802-50MG	50MG	2,4,5-T solution	93-76-5	S-10646B1-1ML	1ML
Simetryn solution	1014-70-6	S-13802A1-1ML	1ML	2,4,5-T solution	93-76-5	S-10646B1-5ML	5ML
Simetryn solution	1014-70-6	S-13802T1-1ML	1ML	2,4,5-T butoxyethyl ester	2545-59-7	N-10648-100MG	100MG
Simetryn solution	1014-70-6	S-13802T1-5ML	5ML	2,4,5-T butoxyethyl ester solution	2545-59-7	S-10648A4-1ML	1ML
Sodium arsenate	7778-43-0	N-13215-250MG	250MG	2,4,5-T butoxyethyl ester solution	2545-59-7	S-10648A4-5ML	5ML
Sodium arsenate solution	7778-43-0	S-13215F1-1ML	1ML	2,4,5-T butoxyethyl ester solution	2545-59-7	S-10648U1-1ML	1ML
Sodium arsenate solution	7778-43-0	S-13215U1-1ML	1ML	2,4,5-T butoxypolypropylene glycol ester (Technical)		N-10004-1G	1G
Sodium dimethyldithiocarbamate hydrate	207233-95-2	N-12863-100MG	100MG	2,4,5-T n-butyl ester	93-79-8	N-10651-250MG	250MG
Sodium fluoracetate (Technical)	62-74-8	N-13216-250MG	250MG	2,4,5-T n-butyl ester solution	93-79-8	S-10651A1-1ML	1ML
Sodiummethyl-dithiocarbamate-dihydrate	6734-80-1	N-12382-250MG	250MG	2,4,5-T n-butyl ester solution	93-79-8	S-10651J1-1ML	1ML
Solan(TM)	2307-68-8	N-12844-100MG	100MG	2,4,5-T 2-ethylhexyl ester	1928-47-8	N-10647-250MG	250MG
Solan(TM) Solution	2307-68-8	S-12844A1-1ML	1ML	2,4,5-T 2-ethylhexyl ester solution	1928-47-8	S-10647A1-1ML	1ML
Solan(TM) Solution	2307-68-8	S-12844U1-1ML	1ML	2,4,5-T 2-ethylhexyl ester solution	1928-47-8	S-10647T1-1ML	1ML
Solicam™	27314-13-2	N-12668-100MG	100MG	2,4,5-T isocetyl ester	25168-15-4	N-10649-100MG	100MG
Solicam™ Solution	27314-13-2	S-12668A1-1ML	1ML	2,4,5-T isocetyl ester solution	25168-15-4	S-10649A1-1ML	1ML
Solicam™ Solution	27314-13-2	S-12668T1-1ML	1ML	2,4,5-T isocetyl ester solution	25168-15-4	S-10649T1-1ML	1ML
Solicam™ Solution	27314-13-2	S-12668T1-5ML	5ML	2,4,5-T isopropyl ester	93-78-7	N-10650-100MG	100MG
Sonar™	59756-60-4	N-13217-100MG	100MG	2,4,5-T isopropyl ester solution	93-78-7	S-10650A1-1ML	1ML
Sonar™ Solution	59756-60-4	S-13217T1-1ML	1ML	2,4,5-T isopropyl ester solution	93-78-7	S-10650T1-1ML	1ML
Sonar™ Solution	59756-60-4	S-13217T1-5ML	5ML	2,4,5-T-methyl ester	1928-37-6	N-10005-250MG	250MG
Spinetoram	187166-40-1	N-13221-100MG	100MG	2,4,5-T-methyl ester Solution	1928-37-6	S-10005M1-1ML	1ML
Spinosad	168316-95-8	N-13222-100MG	100MG	2,4,5-T-methyl ester Solution	1928-37-6	S-10005M1-5ML	5ML
Spinosad solution	168316-95-8	S-13222A1-1ML	1ML	2,4,5-T PGEE		N-10004-1G	1G
Spinosyn A, 131929-60-7; Spinosyn D, 131929-63-0	168316-95-8	N-13222-100MG	100MG	2,4,5-T-triethylamine salt	57213-69-1	N-10652-250MG	250MG
Spinosyn A, 131929-60-7; Spinosyn D, 131929-63-0	168316-95-8	S-13222A1-1ML	1ML	2,4,5-T-triethylamine salt solution	57213-69-1	S-10652M1-1ML	1ML
Spirodiclofen	148477-71-8	N-13223-100MG	100MG	2,3,6-TBA	50-31-7	N-10605-100MG	100MG
Spiromesifen	283594-90-1	N-13224-100MG	100MG	TCA	76-03-9	N-13652-250MG	250MG
Spiromesifen alcohol	148476-30-6	MET-13224-10MG	10MG	TCA	1582-09-8	N-13689-250MG	250MG

Description	CAS#	Part #	Size	Description	CAS#	Part #	Size
TCA Solution	76-03-9	S-13652T1-1ML	1ML	Thiabendazole	148-79-8	N-13560-250MG	250MG
TCA Solution	76-03-9	S-13652T1-5ML	5ML	Thiabendazole solution	148-79-8	S-13560M4-1ML	1ML
TCA Solution	1582-09-8	S-13689A1-1ML	1ML	Thiacloprid	111988-49-9	N-13561-100MG	100MG
TCA Solution	1582-09-8	S-13689U1-1ML	1ML	Thiacloprid solution	111988-49-9	S-13561A1-1ML	1ML
TCA Solution	1582-09-8	S-13689U1-5ML	5ML	Thiacloprid solution	111988-49-9	S-13561U1-1ML	1ML
TCMTB	21564-17-0	N-13500-10MG	10MG	Thiacloprid-amide	676228-91-4	MET-13561-100MG	100MG
TCMTB solution	21564-17-0	S-13500A1-1ML	1ML	Thiamethoxam	153719-23-4	N-13562-100MG	100MG
TCMTB solution	21564-17-0	S-13500X4-1ML	1ML	Thiamethoxam solution	153719-23-4	S-13562M1-1ML	1ML
TCMTB solution	21564-17-0	S-13500X4-5ML	5ML	2-(Thiazol-4-yl)benzimidazole	148-79-8	N-13560-250MG	250MG
TDE	72-54-8	N-10874-250MG	250MG	2-(Thiazol-4-yl)benzimidazole Solution	148-79-8	S-13560M4-1ML	1ML
TDE Solution	72-54-8	S-10874A1-1ML	1ML	Thiazopyr	117718-60-2	N-13563-500MG	500MG
TDE Solution	72-54-8	S-10874K1-1ML	1ML	Thiazopyr solution	117718-60-2	S-13563A1-1ML	1ML
TDE Solution	72-54-8	S-10874K1-5ML	5ML	Thiazopyr solution	117718-60-2	S-13563U1-1ML	1ML
Tebuconazole	107534-96-3	N-12006-250MG	250MG	Thidiazuron	51707-55-2	N-13564-250MG	250MG
Tebuconazole solution	107534-96-3	S-12006A1-1ML	1ML	Thidiazuron solution	51707-55-2	S-13564A1-1ML	1ML
Tebuconazole solution	107534-96-3	S-12006K1-1ML	1ML	Thidiazuron solution	51707-55-2	S-13564U1-1ML	1ML
Tebufenozide	112410-23-8	N-13501-100MG	100MG	Thien carbazole-methyl	317815-83-1	N-13565-50MG	50MG
Tebufenozide solution	112410-23-8	S-13501A1-1ML	1ML	Thiensiulfuron-methyl	79277-27-3	N-13566-100MG	100MG
Tebufenpyrad	119168-77-3	N-13502-100MG	100MG	Thiensiulfuron-methyl solution	79277-27-3	S-13566A1-1ML	1ML
Tebupirifos	96182-53-5	N-13503-100MG	100MG	Thifluzamide	130000-40-7	N-12884-100MG	100MG
Tebupirifos solution	96182-53-5	S-13503A1-1ML	1ML	Thifluzamide solution	130000-40-7	S-12884M1-1ML	1ML
Tebupirifos solution	96182-53-5	S-13503U1-1ML	1ML	Thiodicarb	59669-26-0	N-13568-250MG	250MG
Tebupirifos	96182-53-5	N-13503-100MG	100MG	Thiodicarb solution	59669-26-0	S-13568A1-1ML	1ML
Tebupirifos Solution	96182-53-5	S-13503A1-1ML	1ML	Thiofanocarb	39196-18-4	N-14058-25MG	25MG
Tebupirifos Solution	96182-53-5	S-13503U1-1ML	1ML	Thiofanocarb Solution	39196-18-4	S-14058A1-1ML	1ML
Tebutam	35256-85-0	N-13504-250MG	250MG	Thiofanox	39196-18-4	N-14058-25MG	25MG
Tebutam solution	35256-85-0	S-13504A1-1ML	1ML	Thiofanox solution	39196-18-4	S-14058A1-1ML	1ML
Tebutam solution	35256-85-0	S-13504U1-1ML	1ML	Thiometon	640-15-3	N-13570-50MG	50MG
Tebuthiuron	34014-18-1	N-13505-250MG	250MG	Thiometon solution	640-15-3	S-13570U1-1ML	1ML
Tebuthiuron solution	34014-18-1	S-13505A1-1ML	1ML	Thionazin	297-97-2	N-13760-50MG	50MG
Tebuthiuron solution	34014-18-1	S-13505T1-1ML	1ML	Thionazin Solution	297-97-2	S-13760A1-1ML	1ML
Tebuthiuron solution	34014-18-1	S-13505T1-5ML	5ML	Thionazin Solution	297-97-2	S-13760J1-1ML	1ML
Tebuthiuron TP 108	39222-73-6	N-13948-100MG	100MG	Thionazin Solution	297-97-2	S-13760J1-5ML	5ML
Tecloftalam	76280-91-6	N-14068-10MG	10MG	Thionazin-O-analog	7359-55-9	MET-13760A-50MG	50MG
Tecnazene	117-18-0	N-10111-250MG	250MG	Thiophanate	23564-06-9	N-13571-100MG	100MG
Tecnazene Solution	117-18-0	S-10111M1-1ML	1ML	Thiophanate solution	23564-06-9	S-13571U1-1ML	1ML
Teflubenzuron	83121-18-0	N-13506-250MG	250MG	Thiophanate-methyl	23564-05-8	N-13572-250MG	250MG
Teflubenzuron solution	83121-18-0	S-13506A1-1ML	1ML	Thiophanate-methyl solution	23564-05-8	S-13572A1-1ML	1ML
Telone	542-75-6	N-10193-500MG	500MG	Thiosultap-disodium	52207-48-4	N-11752-250MG	250MG
Telone Solution	542-75-6	S-10193M1-1ML	1ML	Thiosultap-disodium Solution	52207-48-4	S-11752F1-1ML	1ML
Tembotrione	335104-84-2	N-13508-100MG	100MG	Thiotepe	3689-24-5	N-13540-50MG	50MG
Temephos	3383-96-8	N-10996-100MG	100MG	Thiotepe Solution	3689-24-5	S-13540A1-1ML	1ML
Temephos solution	3383-96-8	S-10996A1-1ML	1ML	Thiotepe Solution	3689-24-5	S-13540J1-1ML	1ML
Temephos solution	3383-96-8	S-10996U1-1ML	1ML	Tiadinil	223580-51-6	N-13282-25MG	25MG
TEPP	107-49-3	N-13543-500MG	500MG	Tiadinil solution	223580-51-6	S-13282A1-1ML	1ML
Tepaloxymid solution	149979-41-9	S-13970A1-1ML	1ML	TIBA	88-82-4	N-10601-250MG	250MG
Terbacil	5902-51-2	N-13509-250MG	250MG	TIBA Solution	88-82-4	S-10601A1-1ML	1ML
Terbacil solution	5902-51-2	S-13509A1-1ML	1ML	TIBA Solution	88-82-4	S-10601T1-1ML	1ML
Terbacil solution	5902-51-2	S-13509T1-1ML	1ML	Til(TM)	60207-90-1	N-13576-250MG	250MG
Terbacil solution	5902-51-2	S-13509T1-5ML	5ML	Til(TM) Solution	60207-90-1	S-13576A1-1ML	1ML
Terbufos	13071-79-9	N-13510-250MG	250MG	Til(TM) Solution	60207-90-1	S-13576E1-1ML	1ML
Terbufos solution	13071-79-9	S-13510M1-1ML	1ML	Tokuthion™	34643-46-4	N-13577-100MG	100MG
Terbufos sulfone	56070-16-7	MET-13510A-50MG	50MG	Tokuthion™ Solution	34643-46-4	S-13577A1-1ML	1ML
Terbufos sulfoxide	10548-10-4	MET-13510B-100MG	100MG	Tokuthion™ Solution	34643-46-4	S-13577U1-1ML	1ML
Terbufos-sulfoxide	10548-10-4	MET-13510B-100MG	100MG	Tolclofos-Methyl	57018-04-9	N-13578-500MG	500MG
Terbumeton	33693-04-8	N-13511-250MG	250MG	Tolclofos-Methyl solution	57018-04-9	S-13578A1-1ML	1ML
Terbumeton solution	33693-04-8	S-13511M1-1ML	1ML	Tolclofos-Methyl solution	57018-04-9	S-13578U1-1ML	1ML
Terbumeton-desethyl	30125-64-5	MET-13511-50MG	50MG	Tolifenpyrad	129558-76-5	N-13579-100MG	100MG
Terbuthylazine	5915-41-3	N-13512-50MG	50MG	Tolylfluand	731-27-1	N-13584-250MG	250MG
Terbuthylazine solution	5915-41-3	S-13512A1-1ML	1ML	N-m-Tolylphthalamic acid	85-72-3	N-12557-100MG	100MG
Terbuthylazine solution	5915-41-3	S-13512T4-1ML	1ML	N-m-Tolylphthalamic acid solution	85-72-3	S-12557A1-1ML	1ML
Terbuthylazine-desethyl	30125-63-4	MET-13512A-50MG	50MG	N-m-Tolylphthalamic acid solution	85-72-3	S-12557T1-1ML	1ML
Terbutol	1918-11-200	N-13513-100MG	100MG	1,3-di-p-tolylurea	621-00-1	N-14078-100MG	100MG
Terbutryne	886-50-0	N-13514-250MG	250MG	Topramezone	210631-68-8	N-13585-100MG	100MG
Terbutryne solution	886-50-0	S-13514A1-1ML	1ML	Toxaphene	8001-35-2	N-13586-250MG	250MG
Terbutryne solution	886-50-0	S-13514T1-1ML	1ML	Toxaphene Solution	8001-35-2	S-13586B6-1ML	1ML
Terbutryne solution	886-50-0	S-13514T1-5ML	5ML	Toxaphene Solution	8001-35-2	S-13586J4-1ML	1ML
Teridox(TM)	50563-36-5	N-11754-250MG	250MG	Toxaphene Solution	8001-35-2	S-13586J4-5ML	5ML
Terpine polychlorinates	8001-50-1	N-13229-100MG	100MG	Toxaphene Solution	8001-35-2	S-13586M1-1ML	1ML
Terpine polychlorinates Solution	8001-50-1	S-13229A1-1ML	1ML	Toxaphene Solution	8001-35-2	S-13586M1-5ML	5ML
Terpine polychlorinates Solution	8001-50-1	S-13229K4-1ML	1ML	Trailoxymid	87820-88-0	N-13579-500MG	500MG
2,3,5,6-Tetrachloroaniline	3481-20-7	N-10530-250MG	250MG	Trailoxymid solution	87820-88-0	S-13587A1-1ML	1ML
2,3,5,6-Tetrachloroaniline solution	3481-20-7	S-10530M1-1ML	1ML	Trailoxymid solution	87820-88-0	S-13587U1-1ML	1ML
2,3,4,5-Tetrachloroaniline solution	634-83-3	S-13990A1-1ML	1ML	Tralometrin	66841-25-6	N-13588-100MG	100MG
2,3,5,6-Tetrachloroanisole solution	6936-40-9	S-13984J1-1ML	1ML	Tralometrin solution	66841-25-6	S-13588A1-1ML	1ML
3,4,5,6-Tetrachloro-N-(2,3-dichlorophenyl)phtha-	76280-91-6	N-14068-10MG	10MG	Tralometrin solution	66841-25-6	S-13588T1-1ML	1ML
Tetrachlorohydroquinone	87-87-6	MET-12159A-250MG	250MG	Transfluthrin	118712-89-3	N-13626-100MG	100MG
1,2,4,5-Tetrachloro-3-nitrobenzene	117-18-0	N-10111-250MG	250MG	Tri(cyclohexyl)-1H-1,2,4-triazol-1-yltin Solution	41083-11-8	S-11111A1-1ML	1ML
1,2,4,5-Tetrachloro-3-nitrobenzene solution	117-18-0	S-10111M1-1ML	1ML	Triadimefon	43121-43-3	N-13636-500MG	500MG
2,3,4,6-Tetrachlorophenol	58-90-2	N-10599-50MG	50MG	Triadimefon solution	43121-43-3	S-13636A1-1ML	1ML
2,3,4,6-Tetrachlorophenol solution	58-90-2	S-10599M1-1ML	1ML	Triadimefon solution	43121-43-3	S-13636T1-1ML	1ML
Tetrachlorvinphos	22248-79-9	N-13537-250MG	250MG	Triadimefon solution	43121-43-3	S-13636T1-5ML	5ML
Tetrachlorvinphos solution	22248-79-9	S-13537A1-1ML	1ML	Triadimenol	55219-65-3	N-11129-250MG	250MG
Tetrachlorvinphos solution	22248-79-9	S-13537X1-1ML	1ML	Triadimenol solution	55219-65-3	S-11129A1-1ML	1ML
Tetrachlorvinphos solution	22248-79-9	S-13537X1-5ML	5ML	Triadimenol solution	55219-65-3	S-11129U1-1ML	1ML
Tetraconazole	112281-77-3	N-13538-100MG	100MG	Tri-allate	2303-17-5	N-13628-250MG	250MG
Tetraconazole solution	112281-77-3	S-13538A1-1ML	1ML	Tri-allate solution	2303-17-5	S-13628A1-1ML	1ML
Tetraconazole solution	112281-77-3	S-13538U1-1ML	1ML	Tri-allate solution	2303-17-5	S-13628K1-1ML	1ML
Tetradifon	116-29-0	N-13539-250MG	250MG	Triasulfuron	82097-50-5	N-13638-100MG	100MG
Tetradifon solution	116-29-0	S-13539A1-1ML	1ML	Triasulfuron solution	82097-50-5	S-13638A1-1ML	1ML
Tetradifon solution	116-29-0	S-13539B3-1ML	1ML	Triazamate	112143-82-5	N-13956-100MG	100MG
Tetradifon solution	116-29-0	S-13539B3-5ML	5ML	1,2,4-Triazole	288-88-0	MET-13576A-250MG	250MG
Tetraethylthiopyrophosphate	3689-24-5	N-13540-50MG	50MG	Triazole Solution	60207-31-0	S-13976A1-1ML	1ML
Tetraethylthiopyrophosphate solution	3689-24-5	S-13540A1-1ML	1ML	N-(2-Triazolyl)sulfanilamide	72-14-0	N-13252-250MG	250MG
Tetraethylthiopyrophosphate solution	3689-24-5	S-13540J1-1ML	1ML	Triazophos	24017-47-8	N-13639-100MG	100MG
Tetraethyl monothiopyrophosphate	645-78-3	N-12496-50MG	50MG	Triazophos solution	24017-47-8	S-13639A1-1ML	1ML
Tetraethyl monothiopyrophosphate Solution	645-78-3	S-12496U1-1ML	1ML	Triazophos solution	24017-47-8	S-13639U1-1ML	1ML
Tetraethylpyrophosphate	107-49-3	N-13543-500MG	500MG	Triazophos solution	24017-47-8	S-13639U1-5ML	5ML
cis-1,2,3,6-Tetrahydrophthalimide	1469-48-3	MET-11399A-250MG	250MG	Tribenuron methyl	101200-48-0	N-13640-250MG	250MG
Tetramethrin	7696-12-0	N-13552-100MG	100MG	Tribenuron methyl solution	101200-48-0	S-13640A1-1ML	1ML
Tetramethrin solution	7696-12-0	S-13552A4-1ML	1ML	2,4,6-Tris-(2,4,6-tribromophenoxy)-1,3,5-triazine	25713-60-4	N-12888-50MG	50MG
Tetramethrin solution	7696-12-0	S-13552A4-5ML	5ML	2,4,6-Tris-(2,4,6-tribromophenoxy)-1,3,5-triazine Sol	25713-60-4	S-12888U1-1ML	1ML
Tetramethrin solution	7696-12-0	S-13552U1-1ML	1ML	Tribunil*	18691-97-9	N-12387-100MG	100MG
d-Tetramethrin	1166-46-7	N-13994-100MG	100MG	Tribunil* Solution	18691-97-9	S-12387A1-1ML	1ML
d-Tetramethrin Solution	1166-46-7	S-13994A1-1ML	1ML	Tribunil* Solution	18691-97-9	S-12387U1-1ML	1ML
(1R)-trans-Tetramethrin	1166-46-7	N-13994-100MG	100MG	S,S,S-Tributyl phosphorotriothioate	78-48-8	N-13194-250MG	250MG
(1R)-trans-Tetramethrin solution	1166-46-7	S-13994A1-1ML	1ML	S,S,S-Tributyl phosphorotriothioate solution	78-48-8	S-13194A1-1ML	1ML
Tetramine solution	80-12-6	S-13995X1-1ML	1ML	S,S,S-Tributyl phosphorotriothioate solution	78-48-8	S-13194T1-1ML	1ML
TFNA (flonicamid metabolite)	158063-66-2	MET-11976A-100MG	100MG	Tributylphosphoro-trithioate	150-50-5	N-13648-250MG	250MG
Thallium sulfate(Technical)	7446-18-6	N-13559-250MG	250MG	1,1,1-Trichloro-2	76-03-9	N-FC92-5-5MG	5MG
Thanite	115-31-1	N-12232-250MG	250MG	Trichloroacetic acid	76-03-9	N-13652-250MG	250MG
Thanite Solution	115-31-1	S-12232A1-1ML	1ML	Trichloroacetic acid solution	76-03-9	S-13652T1-1ML	1ML
Thanite Solution	115-31-1	S-12232U1-1ML	1ML	Trichloroacetic acid solution	76-03-9	S-13652T1-5ML	5ML
Thenylchlor solution	96491-05-3	S-13966A1-1ML	1ML	2,3,6-Trichlorobenzoic acid	50-31-7	N-10605-100MG	100MG


Description	CAS#	Part #	Size	Description	CAS#	Part #	Size
1,1,1-Trichloro-2,2-bis[p-chlorophenyl]ethane Sol	50-29-3	S-10876J1-1ML	1ML	Visor(TM)	117718-60-2	N-13563-500MG	500MG
2,4,4'-Trichloro-2'-hydroxydiphenyl ether Solution	3380-34-5	S-13076M1-1ML	1ML	Visor(TM) Solution	117718-60-2	S-13563A1-1ML	1ML
Trichloromethylsulfenyl chloride Solution	594-42-3	S-12845A1-1ML	1ML	Vista(TM)	81406-37-3	N-12001-100MG	100MG
Trichloromethylsulfenyl chloride Solution	594-42-3	S-12845U1-1ML	1ML	Vista(TM) Solution	81406-37-3	S-12001M1-1ML	1ML
Trichloronate	327-98-0	N-13657-100MG	100MG	Vitamin D3	67-97-0	N-11463-100MG	100MG
Trichloronate solution	327-98-0	S-13657A1-1ML	1ML	Vitamin D3 Solution	67-97-0	S-11463A1-1ML	1ML
Trichloronate solution	327-98-0	S-13657L1-1ML	1ML	Viviful(TM)	127277-53-6	N-13100-100MG	100MG
Trichloronate solution	327-98-0	S-13657L1-5ML	5ML	Vydate	23135-22-0	N-12735-100MG	100MG
Trichloronitromethane	76-06-2	N-11452-100MG	100MG	Vydate Solution	23135-22-0	S-12735A1-1ML	1ML
Trichloronitromethane solution	76-06-2	S-11452B1-1ML	1ML	Warfarin	81-81-2	N-13750-250MG	250MG
Trichloronitromethane solution	76-06-2	S-11452B1-5ML	5ML	Warfarin solution	81-81-2	S-13750A1-1ML	1ML
2,4,6-Trichlorophenol	88-06-2	N-10657-250MG	250MG	Weed Rhop	1928-43-4	N-10356-1G	1G
2,4,6-Trichlorophenol solution	88-06-2	S-10657M1-1ML	1ML	Weed Rhop Solution	1928-43-4	S-10356A1-1ML	1ML
2,4,6-Trichlorophenol solution	88-06-2	S-10657M1-5ML	5ML	Weedtrin-D (tm)	6385-62-2	N-11816-500MG	500MG
2,4,5-Trichlorophenol ethanolamine salt		N-10655-250MG	250MG	Weedtrin-D (tm) Solution	6385-62-2	S-11816F1-1ML	1ML
2,4,5-Trichlorophenol ethanolamine salt solution		S-10655F1-1ML	1ML	Xemium®	907204-31-3	N-12921-100MG	100MG
2,4,5-Trichlorophenol ethanolamine salt solution		S-10655U1-1ML	1ML	XMC	2655-14-3	N-13275-100MG	100MG
(2,4,5-Trichlorophenoxy)acetic acid methyl ester	1928-37-6	N-10005-250MG	250MG	XMC solution	2655-14-3	S-13275A1-1ML	1ML
(2,4,5-Trichlorophenoxy)acetic acid methyl ester Sol	1928-37-6	S-10005M1-1ML	1ML	3,5-Xylyl methylcarbamate	2655-14-3	N-13275-100MG	100MG
(2,4,5-Trichlorophenoxy)acetic acid methyl ester Sol	1928-37-6	S-10005M1-5ML	5ML	3,5-Xylyl methylcarbamate Solution	2655-14-3	S-13275A1-1ML	1ML
3,5,6-Trichloro-2-pyridinol	6515-38-4	MET-11459A-250MG	250MG	Xylycarb	2425-10-700	N-13274-10MG	10MG
Trichlorophenol	52-68-6	N-11843-250MG	250MG	Zardex(TM)	54460-46-7	N-11541-500MG	500MG
Triclopyr	55335-06-3	N-13659-500MG	500MG	Zardex(TM) Solution	54460-46-7	S-11541A1-1ML	1ML
Triclopyr solution	55335-06-3	S-13659A1-1ML	1ML	Z-Dimethomorphe	113210-98-3	N-13752-100MG	100MG
Triclopyr solution	55335-06-3	S-13659T1-1ML	1ML	Zectran	315-18-4	N-13753-500MG	500MG
Triclopyr-butyl	64700-56-7	N-13660-250MG	250MG	Zectran (Mexacarbate)	315-18-4	N-13753-500MG	500MG
Triclopyr 2-butoxyethyl ester	64700-56-7	N-13660-250MG	250MG	Zectran (Mexacarbate) Solution	315-18-4	S-13753A4-1ML	1ML
Triclopyr methyl ester	60825-26-5	N-13661-100MG	100MG	Zectran solution	315-18-4	S-13753A4-1ML	1ML
Triclopyr methyl ester solution	60825-26-5	S-13661A1-1ML	1ML	Zinc phosphate(Technical)	7779-90-0	N-13758-250MG	250MG
Triclopyr methyl ester solution	60825-26-5	S-13661T1-1ML	1ML	Zineb(Technical)	12122-67-7	N-13759-250MG	250MG
Triclosan	3380-34-5	N-13076-100MG	100MG	Zinophos (TM)	297-97-2	N-13760-50MG	50MG
Triclosan solution	3380-34-5	S-13076M1-1ML	1ML	Zinophos (TM) solution	297-97-2	S-13760A1-1ML	1ML
Tricyclazole	41814-78-2	N-11132-50MG	50MG	Ziram	137-30-4	N-13761-250MG	250MG
Tricyclazole solution	41814-78-2	S-11132A1-1ML	1ML	Zoxamide	156052-68-5	N-13762-100MG	100MG
Tricyclazole solution	41814-78-2	S-11132T1-1ML	1ML				
Tricyclazole solution	41814-78-2	S-11132T1-5ML	5ML				
Tricyclohexyltin hydroxide	13121-70-5	N-13662-250MG	250MG				
Tricyclohexyltin hydroxide solution	13121-70-5	S-13662M1-1ML	1ML				
Tridemorph (Mix of Isomers)	81412-43-3	N-13665-100MG	100MG				
Tridemorph (Mix of Isomers) solution	81412-43-3	S-13665A1-1ML	1ML				
Tridemorph (Mix of Isomers) solution	81412-43-3	S-13665U1-1ML	1ML				
Tridiphane solution	58138-08-2	S-13965A1-1ML	1ML				
Trietazine	1912-26-1	N-13666-100MG	100MG				
Trietazine solution	1912-26-1	S-13666A1-1ML	1ML				
Trietazine solution	1912-26-1	S-13666X1-1ML	1ML				
Trietazine-desethyl	38902-68-0	MET-13666A-100MG	100MG				
Trietazine-desethyl Solution	38902-68-0	MET-13666AA1-1ML	1ML				
Trietazine-ethylamine	93567-85-2	N-13989-100MG	100MG				
Trietazine-2-hydroxy	13532-25-7	MET-13666B-10MG	10MG				
Trietazine-2-hydroxy Solution	13532-25-7	MET-13666BA1-1ML	1ML				
Trifloxystrobin	141517-21-7	N-13684-100MG	100MG				
Trifloxystrobin solution	141517-21-7	S-13684A1-1ML	1ML				
Trifloxysulfuron-Sodium	199119-58-9	N-13685-100MG	100MG				
Triflumizole	68694-11-1	N-13686-250MG	250MG				
Triflumuron	64628-44-0	N-13687-100MG	100MG				
Triflumuron solution	64628-44-0	S-13687A1-1ML	1ML				
Triflumuron solution	64628-44-0	S-13687U1-1ML	1ML				
alpha, alpha, alpha-Trifluoro-3'-hydroxy-O-toluanilide		MET-12004C-250MG	250MG				
4-(Trifluoromethyl)nicotinic acid	158063-66-2	MET-11976A-100MG	100MG				
Trifluralin	1582-09-8	N-13689-250MG	250MG				
Trifluralin (di-n-propyl-d14) Solution		S-FD2140S-1.2ML	1.2ML				
Trifluralin solution	1582-09-8	S-13689A1-1ML	1ML				
Trifluralin solution	1582-09-8	S-13689U1-1ML	1ML				
Trifluralin solution	1582-09-8	S-13689U1-5ML	5ML				
Triflusalufuron-methyl	126535-15-7	N-13690-50MG	50MG				
Triflusalufuron-methyl solution	126535-15-7	S-13690A1-1ML	1ML				
Triforine	26644-46-2	N-13691-250MG	250MG				
Triforine solution	26644-46-2	S-13691U1-1ML	1ML				
2,3,5-Triiodobenzoic acid	88-82-4	N-10601-250MG	250MG				
2,3,5-Triiodobenzoic acid solution	88-82-4	S-10601A1-1ML	1ML				
2,3,5-Triiodobenzoic acid solution	88-82-4	S-10601T1-1ML	1ML				
2,3,5-Trimethacarb	2655-15-4	N-10603-50MG	50MG				
3,4,5-Trimethacarb	2686-99-9	N-10769-10MG	10MG				
2,3,5-Trimethacarb Solution	2655-15-4	S-10603A1-1ML	1ML				
2,3,5-Trimethylphenyl methyl carbamate	2655-15-4	N-10603-50MG	50MG				
3,4,5-Trimethylphenyl methyl carbamate	2686-99-9	N-10769-10MG	10MG				
2,3,5-Trimethylphenyl methyl carbamate solution	2655-15-4	S-10603A1-1ML	1ML				
Trinexapac	104273-73-6	MET-13706A-100MG	100MG				
Trinexapac-ethyl	95266-40-3	N-13706-250MG	250MG				
Trinexapac-ethyl solution	95266-40-3	S-13706A1-1ML	1ML				
Trinexapac-ethyl solution	95266-40-3	S-13706T1-1ML	1ML				
Triticonazole	131983-72-7	N-13729-100MG	100MG				
Tycor	64529-56-2	N-13731-100MG	100MG				
Tycor Solution	64529-56-2	S-13731M1-1ML	1ML				
Umbelliferone	93-35-6	MET-13750A-250MG	250MG				
2-Undecanone	112-12-9	N-10520-1G	1G				
2-Undecanone solution	112-12-9	S-10520N1-1ML	1ML				
Uniconazole	83657-22-1	N-14032-50MG	50MG				
Uniconazole solution	83657-22-1	S-14032A1-1ML	1ML				
Uniconazole-P	83657-17-4	N-13735-10MG	10MG				
Uniconazole-P solution	83657-17-4	S-13735M1-1ML	1ML				
Upbeet(TM)	126535-15-7	N-13690-50MG	50MG				
Upbeet(TM) Solution	126535-15-7	S-13690A1-1ML	1ML				
Urea	4300-43-0	MET-11827A-100MG	100MG				
Valifenalate	283159-90-0	N-13281-25MG	25MG				
Valifenalate solution	283159-90-0	S-13281A1-1ML	1ML				
Valiphenal	283159-90-0	N-13281-25MG	25MG				
Valiphenal Solution	283159-90-0	S-13281A1-1ML	1ML				
Valone	83-28-3	N-10384-250MG	250MG				
Valone Solution	83-28-3	S-10384A1-1ML	1ML				
Valone Solution	83-28-3	S-10384U1-1ML	1ML				
Vamidothion	2275-23-2	N-13741-100MG	100MG				
Vamidothion solution	2275-23-2	S-13741A1-1ML	1ML				
Vangard(TM)	121552-61-2	N-11548-250MG	250MG				
Vangard(TM) Solution	121552-61-2	S-11548A1-1ML	1ML				
Vangard(TM) Solution	121552-61-2	S-11548U1-1ML	1ML				
Vanicide-20S (TM)	5902-85-2	N-13742-250MG	250MG				
Vanicide-20S (TM) solution	5902-85-2	S-13742A1-1ML	1ML				
Velpar	51235-04-2	N-13743-250MG	250MG				
Velpar solution	51235-04-2	S-13743A1-1ML	1ML				
Vernolate	1929-77-7	N-13744-250MG	250MG				
Vernolate solution	1929-77-7	S-13744A1-1ML	1ML				
Vinclozolin	50471-44-8	N-13745-250MG	250MG				
Vinclozolin solution	50471-44-8	S-13745A1-1ML	1ML				



Pesticide Standards

Listing by CAS Number with Structure

a	→	Nitenpyram		CAS#150824-47-8	←	k
b	→	Chemical Name: 1,1-Ethenediamine, N-[(6-chloro-3-pyridinyl)methyl]-N-ethyl-N'-methyl-2-nitro-, (1E)-				
c	→	PS: Solid				
d	→	MF: C ₁₁ H ₁₅ ClN ₄ O ₂				
e	→	MW: 270.75				
f	→	MP: 82°C				
g	→	FP: 70°C				
h	→	Alternate Name: Capstar				
i	→	Part Number N-12657-100MG S-12657M1-1ML	Matrix Neat 100ug/ml in Methanol			



← l

← m

- a Product Description
- b Chemical Name
- c Physical State
- d Molecular Formula
- e Molecular Weight
- f Melting Point
- g Flash Point

- h Alternate Name
- i Part Number
- j Matrix
- k Chemical Abstracts Registry Number
- l Structure Image
- m Size

Formaldehyde (in water)
CAS#50-00-0


PS: Liquid
MF: CH₂O
MW: 30.03
MP: -92°C


Alternate Name: Formalin

Part Number	Matrix	Size
N-12012-1G	Neat	1G

4,4'-DDT
CAS#50-29-3
Chemical Name: 1,1,1-Trichloro-2,2-bis[p-chlorophenyl]ethane


PS: Solid
MF: C₁₄H₉Cl₅
MW: 354.49
MP: 108.5°C


Alternate Name: Dichlorodiphenyltrichloroethane

Part Number	Matrix	Size
N-10876-100MG	Neat	100MG
S-10876J1-1ML	100ug/ml in Hexane	1ML
N-FC92-5-5MG	Neat (C13)	5MG
S-FC92S-1.2ML	100ug/ml in n-Nonane (C13)	1.2ML

2,3,6-Trichlorobenzoic acid
CAS#50-31-7


PS: Solid
MF: C₇H₃Cl₃O₂
MW: 225.45
MP: 124.5°C


Alternate Name: 2,3,6-TBA

Part Number	Matrix	Size
N-10605-100MG	Neat	100MG

Niclosamide
CAS#50-65-7
Chemical Name: 2,5-Dichloro-4-nitrosalicylanilide


PS: Solid
MF: C₁₃H₈Cl₂N₂O₄
MW: 327.1
MP: 225°C



Part Number	Matrix	Size
N-12653-250MG	Neat	250MG

Piperonyl butoxide
CAS#51-03-6
Chemical Name: 2-(2-butoxyethoxy)ethyl 6-propyl piperonyl ether


PS: Liquid
MF: C₁₉H₃₀O₅
MW: 338.43



Part Number	Matrix	Size
N-13061-100MG	Neat	100MG
S-13061A1-1ML	100ug/ml in Acetonitrile	1ML
S-13061K1-1ML	100ug/ml in Isooctane	1ML

3,5-Dichlorobenzoic acid
CAS#51-36-5


PS: Solid
MF: C₇H₄Cl₂O₂
MW: 191.01



Part Number	Matrix	Size
N-10770-250MG	Neat	250MG
S-10770A1-1ML	100ug/ml in Acetonitrile	1ML
S-10770T5-1ML	2000ug/ml in t-Butyl methyl ether	1ML
S-10770T5-5ML	2000ug/ml in t-Butyl methyl ether	5ML

Bronopol
CAS#52-51-7
Chemical Name: 2-Bromo-2-nitro-1,3-propanediol


PS: Solid
MF: C₃H₆BrNO₄
MW: 199.99
MP: 131.5°C



Part Number	Matrix	Size
N-11348-250MG	Neat	250MG
S-11348M1-1ML	100ug/ml in Methanol	1ML

Trichlorphon
CAS#52-68-6
Chemical Name: Dimethyl-[2,2,2-trichloro-1-hydroxyethyl]phosphonate

PS: Solid
MF: C₄H₈Cl₃O₄P
MW: 257.45
MP: 77°C



Alternate Name: Dylox (TM)

Part Number	Matrix	Size
N-11843-250MG	Neat	250MG

Famphur CAS#52-85-7

Chemical Name: O,O-Dimethyl-O-p-[dimethylsulfamoyl]phenylphosphorothioate

PS: Solid
MF: C₁₀H₁₆NO₅PS₂
MW: 325.36
MP: 52.5°C




Part Number	Matrix	Size
N-11944-250MG	Neat	250MG
S-11944K1-1ML	100ug/ml in Isooctane	1ML

o,p'-DDD CAS#53-19-0

Chemical Name: 1,1-Dichloro-2-(o-chlorophenyl)-2-(p-chlorophenyl)ethane

PS: Solid
MF: C₁₄H₁₀Cl₄
MW: 320.05
MP: 77°C



Alternate Name: mitotane

Part Number	Matrix	Size
N-12706-250MG	Neat	250MG
S-12706M1-1ML	100ug/ml in Methanol	1ML

(S)-Nicotine CAS#54-11-5

Chemical Name: 1-Methyl-2-[3-pyridyl]pyrrolidine

PS: Liquid
MF: C₁₀H₁₄N₂
MW: 162.23
MP: -79°C




Part Number	Matrix	Size
N-12655-1G	Neat	1G
S-12655M1-1ML	100ug/ml in Methanol	1ML

Fenthion CAS#55-38-9

Chemical Name: O,O-Dimethyl-O-[4-(methylthio)-m-tolyl]phosphorothioate


PS: Liquid
MF: C₁₀H₁₅O₃PS₂
MW: 278.34
MP: 7°C



Part Number	Matrix	Size
N-11964-250MG	Neat	250MG
S-11964S1-1ML	100ug/ml in Propanol	1ML
S-11964S1-5ML	100ug/ml in Propanol	5ML

Phenylmercury nitrate(Technical) CAS#55-68-5


PS: Solid
MF: C₆H₅HgNO₃
MW: 339.71
MP: 176°C


Alternate Name: Merphenyl nitrate

Part Number	Matrix	Size
N-13032-100MG	Neat	100MG

Bis[tri-n-butyltin]oxide(Technical) CAS#56-35-9

PS: Liquid
MF: C₂₄H₅₄OSn₂
MW: 596.11
MP: 45°C




Part Number	Matrix	Size
N-11322-250MG	Neat	250MG

Parathion CAS#56-38-2

Chemical Name: O,O-Diethyl-O-p-nitrophenyl phosphorothioate

PS: Liquid
MF: C₁₀H₁₄NO₅PS
MW: 291.26
MP: 6.1°C




Part Number	Matrix	Size
N-12819-100MG	Neat	100MG
S-12819A1-1ML	100ug/ml in Acetonitrile	1ML
S-12819U1-1ML	100ug/ml in Toluene	1ML
N-FD998-E-0.01G	Neat (Deuterated)	0.01G

Coumaphos CAS#56-72-4

Chemical Name: O,O-Diethyl-O-[3-chloro-4-methyl-2-oxo-2H-1-benzopyran-7-yl]

PS: Solid
MF: C₁₄H₁₆ClO₅PS
MW: 362.77
MP: 95°C




Part Number	Matrix	Size
N-11507-100MG	Neat	100MG
S-11507A1-1ML	100ug/ml in Acetonitrile	1ML
S-11507U1-1ML	100ug/ml in Toluene	1ML
S-11507U1-5ML	100ug/ml in Toluene	5ML

Sulfamethazine CAS#57-68-1

Chemical Name: 4-Amino-N-(4,6-dimethyl-2-pyrimidinyl)benzenesulfonamide

PS: Solid
MF: C₁₂H₁₄N₄O₂S
MW: 278.32
MP: 198°C



Alternate Name: Diazil

Part Number	Matrix	Size
N-13245-250MG	Neat	250MG

Chlordane CAS#57-74-9

Chemical Name: Methyl-2-chloro-9-hydroxyfluorene-9-carboxylate

PS: Liquid
MF: C₁₀H₆Cl₈
MW: 409.8
MP: 104°C



Alternate Name: Octachloro-hexahydro-methano-1H-indene

Part Number	Matrix	Size
N-11425-250MG	Neat	250MG
S-11425J4-1ML	1000ug/ml in Hexane	1ML
S-11425J4-5ML	1000ug/ml in Hexane	5ML
S-11425M1-1ML	100ug/ml in Methanol	1ML
S-FC915-1.2ML	100ug/ml in n-Nonane (C13)	1.2ML

Lindane (BHC gamma isomer) CAS#58-89-9

Chemical Name: 1,2,3,4,5,6-Hexachlorocyclohexane-gamma isomer


PS: Solid
MF: C₆H₆Cl₆
MW: 290.83
MP: 112.5°C



Part Number	Matrix	Size
N-12319-250MG	Neat	250MG
S-12319A1-1ML	100ug/ml in Acetonitrile	1ML
S-FC104S-1.2ML	100ug/ml in n-Nonane (C13)	1.2ML
S-FCD104S-1.2ML	100ug/ml in n-Nonane (C13)	1.2ML

2,3,4,6-Tetrachlorophenol CAS#58-90-2

PS: Solid
MF: C₆H₂Cl₄O
MW: 231.89
MP: 70°C



Alternate Name: Dowicide 6(TM)

Part Number	Matrix	Size
N-10599-50MG	Neat	50MG
S-10599M1-1ML	100ug/ml in Methanol	1ML

Sulfaquinoxaline CAS#59-40-5

Chemical Name: 4-Amino-N-2-quinoxalinybenzenesulfonamide

PS: Solid
MF: C₁₄H₁₂N₄O₂S
MW: 300.36
MP: 247.5°C




Part Number	Matrix	Size
N-13251-10MG	Neat	10MG
S-13251M1-1ML	100ug/ml in Methanol	1ML

Dimethoate CAS#60-51-5

Chemical Name: O,O-Dimethyl-S-[N-methylcarbomoylmethyl]phosphorodithioate

PS: Solid
MF: C₅H₁₂NO₃PS₂
MW: 229.28
MP: 52°C




Part Number	Matrix	Size
N-11758-100MG	Neat	100MG
S-11758J1-1ML	100ug/mL in Hexane	1ML
S-11758J1-5ML	100ug/mL in Hexane	5ML

Dieldrin CAS#60-57-1

Chemical Name: Hexachloro-epoxy-octahydro-endo-exo-dimethanonaphthalene

PS: Solid
MF: C₁₂H₈Cl₆O
MW: 380.91
MP: 175.5°C




Part Number	Matrix	Size
N-11688-100MG	Neat	100MG
S-11688M1-1ML	100ug/ml in Methanol	1ML
S-11688M1-5ML	100ug/ml in Methanol	5ML
S-FC90AS-1.2ML	100ug/mL in Nonane (C13)	1.2ML

Aminotriazole CAS#61-82-5

Chemical Name: 1H-1,2,4-Triazole-3-amino


PS: Solid
MF: C₂H₄N₄
MW: 84.08
MP: 159°C


Alternate Name: Amitrole

Part Number	Matrix	Size
N-11067-250MG	Neat	250MG
S-11067A1-1ML	100ug/ml in Acetonitrile	1ML
S-11067B1-1ML	100ug/ml in Acetone	1ML

Phenyl mercuric acetate(Technical) CAS#62-38-4

PS: Solid
MF: C₈H₈HgO₂
MW: 336.75
MP: 149°C




Alternate Name: PMA TM

Part Number	Matrix	Size
N-13011-250MG	Neat	250MG

Dichlorvos CAS#62-73-7

Chemical Name: 2,2-Dichlorovinyl dimethylphosphate


PS: Liquid
MF: C₄H₇Cl₂O₄P
MW: 220.98



Part Number	Matrix	Size
N-11675-250MG	Neat	250MG
S-11675J4-1ML	1000ug/ml in Hexane	1ML
S-11675U1-1ML	100ug/ml in Toluene	1ML

Sodium fluoroacetate(Technical) CAS#62-74-8

PS: Solid
MF: C₂H₃FO₂Na
MW: 100.02
MP: 200°C




Alternate Name: Na2AsHO4

Part Number	Matrix	Size
N-13216-250MG	Neat	250MG

Carbaryl CAS#63-25-2

Chemical Name: 1-Naphthyl N-methylcarbamate

PS: Solid
MF: C₁₂H₁₁NO₂
MW: 201.24
MP: 142°C




Part Number	Matrix	Size
N-11402-250MG	Neat	250MG
S-11402A4-1ML	1000ug/ml in Acetonitrile	1ML
S-11402A4-5ML	1000ug/ml in Acetonitrile	5ML
S-FC2005S-1.2ML	100ug/ml in n-Nonane (C13)	1.2ML

Cycloheximide CAS#66-81-9

Chemical Name: 3-[2(3,5-Dimethyl-2-oxocyclohexyl)-2-hydroxyethyl]glutarimide

PS: Solid
MF: C₁₅H₂₃NO₄
MW: 281.35
MP: 119.5°C




Part Number	Matrix	Size
N-11534-100MG	Neat	100MG
S-11534A1-1ML	100ug/ml in Acetonitrile	1ML
S-11534U1-1ML	100ug/ml in Toluene	1ML

Cholecalciferol CAS#67-97-0

Chemical Name: (1S,3Z)-3-[(2E)-2-[(1R,3aS,7aR)-1-[(1R)-1,5-dimethylhexyl]octahydro-7a-methyl-4H-inden-4-ylidene]ethylidene]-4-methylenecyclohexanol

PS: Solid
MF: C₂₇H₄₄O
MW: 384.64
MP: 84°C
FP: 214.2°C




Alternate Name: Vitamin D3

Part Number	Matrix	Size
N-11463-100MG	Neat	100MG
S-11463A1-1ML	100ug/ml in Acetonitrile	1ML
S-11463U1-1ML	100ug/ml in Toluene	1ML

Hexachlorophene CAS#70-30-4

Chemical Name: 2,2'-Methylene-bis(3,4,6-trichlorophenol)

PS: Solid
MF: C₁₃H₆Cl₆O₂
MW: 406.91
MP: 164°C




Alternate Name: Hexide

Part Number	Matrix	Size
N-12163-250MG	Neat	250MG
S-12163M1-1ML	100ug/ml in Methanol	1ML
S-12163M1-5ML	100ug/ml in Methanol	5ML

Sulfathiazole CAS#72-14-0

Chemical Name: 4-Amino-N-2-thiazolylbenzenesulfonamide

PS: Solid
MF: C₉H₉N₃O₂S₂
MW: 255.32
MP: 175°C




Alternate Name: N-(2-Triazolyl)sulfanilamide

Part Number	Matrix	Size
N-13252-250MG	Neat	250MG

Endrin CAS#72-20-8
Chemical Name: Hexachloro-epoxy-octahydro-endo-endo-dimethanonaphthalene

PS: Solid

MF: C₁₂H₈Cl₆O


MW: 380.9


Part Number	Matrix	Size
N-11854-100MG	Neat	100MG
S-11854M1-1ML	100ug/ml in Methanol	1ML
S-11854M1-5ML	100ug/ml in Methanol	5ML

Methoxychlor CAS#72-43-5
Chemical Name: 1,1,1-Trichloro-2,2-bis[p-methoxyphenyl]ethane

PS: Solid

MF: C₁₆H₁₅Cl₃O₂
MW: 345.65


MP: 89°C


Part Number	Matrix	Size
N-12403-100MG	Neat	100MG
S-12403M1-1ML	100ug/ml in Methanol	1ML
S-12403M1-5ML	100ug/ml in Methanol	5ML

4,4'-DDD CAS#72-54-8
Chemical Name: 1,1-Dichloro-2,2-bis(p-chlorophenyl)ethane

PS: Solid

MF: C₁₄H₁₀Cl₄
MW: 320.04


MP: 109°C

Alternate Name: TDE

Part Number	Matrix	Size
N-10874-250MG	Neat	250MG
S-10874A1-1ML	100ug/ml in Acetonitrile	1ML
S-10874K1-1ML	100ug/ml in Isooctane	1ML
S-10874K1-5ML	100ug/ml in Isooctane	5ML
S-FD94S-1.2ML	100ug/ml in n-Nonane (Deuterated)	1.2ML

4,4'-DDE CAS#72-55-9
Chemical Name: 1,1-Dichloro-2,2-bis[p-chlorophenyl]ethylene

PS: Solid

MF: C₁₄H₈Cl₄
MW: 318.03


MP: 89°C

Alternate Name: p,p-DDE

Part Number	Matrix	Size
N-10875-100MG	Neat	100MG
S-10875M1-1ML	100ug/ml in Methanol	1ML
S-10875M1-5ML	100ug/ml in Methanol	5ML
N-FC93-5-5MG	Neat (C13)	5MG
S-FC93S-1.2ML	100ug/ml in n-Nonane (C13)	1.2ML

Ethylan CAS#72-56-0
Chemical Name: 1,1-Dichloro-2,2-bis[p-ethylphenyl]ethane


PS: Solid

MF: C₁₈H₂₀Cl₂
MW: 307.26

MP: 56°C


Part Number	Matrix	Size
N-12849-250MG	Neat	250MG
S-12849A1-1ML	100ug/ml in Acetonitrile	1ML
S-12849K4-1ML	1000ug/ml in Isooctane	1ML
S-12849K4-5ML	1000ug/ml in Isooctane	5ML

Dimethylarsinic acid CAS#75-60-5
PS: Solid


MF: C₂H₇AsO₂
MW: 138

Alternate Name: Cacodylic acid

Part Number	Matrix	Size
N-11779-500MG	Neat	500MG
S-11779F4-1ML	1000ug/ml in H2O	1ML

Dalapon CAS#75-99-0
Chemical Name: 2,2-Dichloropropionic acid

PS: Liquid

MF: C₃H₄Cl₂O₂
MW: 142.97

MP: 20°C



Part Number	Matrix	Size
N-11562-250MG	Neat	250MG
S-11562B1-1ML	100ug/ml in Acetone	1ML
S-11562M1-1ML	100ug/ml in Methanol	1ML

Trichloroacetic acid CAS#76-03-9
Chemical Name: S-(2,3,3-trichloro-2-propen-1-yl) N,N-bis(1-methylethyl)carbamothioate

PS: Solid

MF: C₂HCl₃O₂
MW: 163.4


MP: 57.5°C

FP: 113°C

Alternate Name: TCA

Part Number	Matrix	Size
N-13652-250MG	Neat	250MG
S-13652T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-13652T1-5ML	100ug/ml in t-Butyl methyl ether	5ML

Trichloronitromethane
CAS#76-06-2


PS: Liquid
MF: CCl₃NO₂
MW: 164.39
MP: -64°C


Alternate Name: CHLOROPICRIN

Part Number	Matrix	Size
N-11452-100MG	Neat	100MG
S-11452B1-1ML	100ug/ml in Acetone	1ML
S-11452B1-5ML	100ug/ml in Acetone	5ML

Heptachlor
CAS#76-44-8
Chemical Name: 1,4,5,6,7,8,8-heptachloro-3a,4,7,7a-tetrahydro-4,7-methanoindene


PS: Solid
MF: C₁₀H₅Cl₇
MW: 373.32
MP: 95°C


Alternate Name: Heptachlor-tetrahydro-methanoindene

Part Number	Matrix	Size
N-12147-100MG	Neat	100MG
S-12147M1-1ML	100ug/ml in Methanol	1ML
S-12147M1-5ML	100ug/ml in Methanol	5ML
S-FC100S-1.2ML	100ug/ml in n-Nonane (C13)	1.2ML

Fentin hydroxide
CAS#76-87-9
Chemical Name: Triphenyltin hydroxide


PS: Solid
MF: C₁₈H₁₆OSn
MW: 367.03
MP: 119°C


Alternate Name: SUPER-TIN

Part Number	Matrix	Size
N-11966-250MG	Neat	250MG
S-11966A1-1ML	100ug/ml in Acetonitrile	1ML

Gibberellic acid
CAS#77-06-5


PS: Solid
MF: C₁₉H₂₂O₆
MW: 346.37
MP: 233°C


Alternate Name: Gibberellin X

Part Number	Matrix	Size
N-12108-100MG	Neat	100MG
S-12108A1-1ML	100ug/ml in Acetonitrile	1ML
S-12108T1-1ML	100ug/ml in T-Butyl methyl Ether	1ML

Dicyclopentadiene
CAS#77-73-6


PS: Liquid
MF: C₁₀H₁₂
MW: 132.21
MP: 33.6°C


Alternate Name: CYCLOPENTADIENE DIMER

Part Number	Matrix	Size
N-11686-500MG	Neat	500MG
S-11686K1-1ML	100ug/ml in Isooctane	1ML

Dioxathion
CAS#78-34-2
Chemical Name: p-Dioxane-2,3-diyl ethyl phosphorodithioate


PS: Liquid
MF: C₁₂H₂₆O₆P₂S₄
MW: 456.54
MP: -20°C



Part Number	Matrix	Size
N-11790-100MG	Neat	100MG
S-11790B1-1ML	100ug/ml in Acetone	1ML

S,S,S-Tributyl phosphorotrithioate
CAS#78-48-8


PS: Liquid
MF: C₁₂H₂₇OPS₃
MW: 314.54


Alternate Name: DEF

Part Number	Matrix	Size
N-13194-250MG	Neat	250MG
S-13194A1-1ML	100ug/ml in Acetonitrile	1ML
S-13194T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

1,2-Dichloropropane
CAS#78-87-5


PS: Liquid
MF: C₃H₆Cl₂
MW: 112.99
MP: -100.4°C


Alternate Name: Propylene dichloride

Part Number	Matrix	Size
N-10156-1G	Neat	1G
S-10156M1-1ML	100ug/ml in Methanol	1ML
S-10156M1-5ML	100ug/ml in Methanol	5ML


Chloroacetic acid
CAS#79-11-8
Chemical Name: Monochloroacetic acid

PS: Solid

MF: C₂H₃ClO₂
MW: 94.5


Part Number	Matrix	Size
N-11438-250MG	Neat	250MG
S-11438T1-1ML	100ug/mL in tert-Butyl methyl ether	1ML

Tetramine Solution
CAS#80-12-6
Chemical Name: Tetramethylenedisulfotetramine


PS: Liquid


Part Number	Matrix	Size
S-13995X1-1ML	100ug/mL in Methylene chloride	1ML

Sulfachloropyridazine
CAS#80-32-0
Chemical Name: 4-Amino-N-(6-chloro-3-pyridazinyl)benzenesulfonamide

PS: Solid

MF: C₁₀H₉ClN₄O₂S

MW: 284.73

Alternate Name: Sulfachloropyridazine


Part Number	Matrix	Size
N-13241-250MG	Neat	250MG

Ovex
CAS#80-33-1
Chemical Name: p-Chlorophenyl-p-chlorobenzene sulfonate

PS: Solid

MF: C₁₂H₈Cl₂O₃S

MW: 303.17

MP: 86.5°C

Alternate Name: Chlorofenson (TM)

Part Number	Matrix	Size
N-12730-250MG	Neat	250MG
S-12730A1-1ML	100ug/mL in Acetonitrile	1ML
S-12730J1-1ML	100ug/mL in Hexane	1ML


Fenson
CAS#80-38-6
Chemical Name: 4-chlorophenyl benzenesulfonate

PS: Solid

MF: C₁₂H₉ClO₃S

MW: 268.72

MP: 62°C


FP: °C

Alternate Name: Aracid

Part Number	Matrix	Size
N-13812-250MG	Neat	250MG

Warfarin
CAS#81-81-2
Chemical Name: 3-(a-Acetyl benzyl)-4-hydroxy-coumarin

PS: Solid

MF: C₁₉H₁₆O₄
MW: 308.33


MP: 161°C


Part Number	Matrix	Size
N-13750-250MG	Neat	250MG
S-13750A1-1ML	100ug/mL in Acetonitrile	1ML
S-13750U1-1ML	100ug/mL in Toluene	1ML

1-8-Naphthalic anhydride
CAS#81-84-5
Chemical Name: Naphthalene-1,8-dicarboxylic anhydride

PS: Solid


MF: C₁₂H₆O₃
MW: 198.18

MP: 273°C

Alternate Name: Protect (TM)

Part Number	Matrix	Size
N-10012-500MG	Neat	500MG
S-10012A1-1ML	100ug/mL in Acetonitrile	1ML
S-10012U1-1ML	100ug/mL in Toluene	1ML

Diphacinone
CAS#82-66-6
Chemical Name: 2-Diphenylacetyl-1,3-indandione

PS: Solid


MF: C₂₃H₁₆O₃
MW: 340.4


Part Number	Matrix	Size
N-11793-100MG	Neat	100MG
S-11793A1-1ML	100ug/mL in Acetonitrile	1ML

Pentachloronitrobenzene
CAS#82-68-8
Chemical Name: Quintozene

PS: Solid

MF: C₆Cl₅NO₂
MW: 295.34


MP: 144°C

Alternate Name: PCNB

Part Number	Matrix	Size
N-12830-250MG	Neat	250MG
S-12830H5-1ML	2000ug/ml in Ethyl acetate	1ML
S-12830H5-5ML	2000ug/ml in Ethyl acetate	5ML
S-12830K4-1ML	1000ug/ml in Isooctane	1ML
S-12830K4-5ML	1000ug/ml in Isooctane	5ML
S-12830K7-1ML	5000ug/ml in Isooctane	1ML

2-Pivaloyl-1,3-indandione
CAS#83-26-1
Chemical Name: 2-pivaloylindan-1,3-dione

PS: Solid

MF: C₁₄H₁₄O₃
MW: 230.26


MP: 108.5°C

Alternate Name: Pindone

Part Number	Matrix	Size
N-10512-250MG	Neat	250MG
S-10512A1-1ML	100ug/ml in Acetonitrile	1ML
S-10512U1-1ML	100ug/ml in Toluene	1ML

2-Isovaleryl-1,3-indanedione
CAS#83-28-3
Chemical Name: 2-isovaleryl-1,3-dione

PS: Solid

MF: C₁₄H₁₄O₃
MW: 230.26


MP: 68°C

Alternate Name: Valone

Part Number	Matrix	Size
N-10384-250MG	Neat	250MG
S-10384A1-1ML	100ug/ml in Acetonitrile	1ML
S-10384U1-1ML	100ug/ml in Toluene	1ML

Rotenone
CAS#83-79-4
Chemical Name: (2R,6aS,12aS)-1,2,6,6a,12,12a-hexahydro-2-isopropenyl-8,9-dimethoxychromeno[3,4-b]furo[2,3-h]chromen-6-one

PS: Solid


MF: C₂₃H₂₂O₆
MW: 394.41

MP: 165°C



Part Number	Matrix	Size
N-13184-250MG	Neat	250MG
S-13184A1-1ML	100ug/ml in Acetonitrile	1ML

9,10-Anthraquinone
CAS#84-65-1
PS: Solid

MF: C₁₄H₈O₂
MW: 208.2

MP: 286°C

Alternate Name: Corbit

Part Number	Matrix	Size
N-10970-250MG	Neat	250MG
S-10970M1-1ML	100ug/ml in Methanol	1ML


Phenyl mercuric phthalate
CAS#84-70-8
PS: Solid


Part Number	Matrix	Size
N-13016-100MG	Neat	100MG
S-13016A1-1ML	100ug/ml in Acetonitrile	1ML
S-13016U1-1ML	100ug/ml in Toluene	1ML

Fenac
CAS#85-34-7
Chemical Name: 2,3,6-Trichlorophenyl acetic acid

PS: Solid


MF: C₈H₅Cl₃O₂
MW: 239.5

MP: 161°C

Alternate Name: Chlorfenac

Part Number	Matrix	Size
N-11945-250MG	Neat	250MG
S-11945A1-1ML	100ug/ml in Acetonitrile	1ML
S-11945T1-1ML	100ug/ml in T-Butyl methyl	1ML

N-m-Tolylphthalamic acid
CAS#85-72-3
PS: Solid


MF: C₁₅H₁₃NO₃
MW: 255.29

MP: 149°C

Alternate Name: Duraset (TM)

Part Number	Matrix	Size
N-12557-100MG	Neat	100MG
S-12557A1-1ML	100ug/ml in Acetonitrile	1ML
S-12557T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Diphenylacetonitrile
CAS#86-29-3

PS: Solid
MF: C₁₄H₁₁N



Alternate Name: Diphenatril

Part Number	Matrix	Size
N-11800-100MG	Neat	100MG
S-11800M1-1ML	100ug/ml in Methanol	1ML

Azinphos-methyl
CAS#86-50-0

Chemical Name: O,O-dimethyl S-[(4-oxo-1,2,3-benzotriazin-3(4H)-yl)methyl] phosphorodithioate


PS: Solid
MF: C₁₀H₁₂N₃O₃PS₂
MW: 317.3
MP: 72°C


Alternate Name: Guthion

Part Number	Matrix	Size
N-12139-250MG	Neat	250MG
S-12139A1-1ML	100ug/ml in Acetonitrile	1ML
S-12139U1-1ML	100ug/ml in Toluene	1ML
S-12139U1-5ML	100ug/ml in Toluene	5ML

1-Naphthaleneacetamide
CAS#86-86-2


PS: Solid
MF: C₁₂H₁₁NO
MW: 185.24
MP: 182°C


Alternate Name: Rootene

Part Number	Matrix	Size
N-10074-1G	Neat	1G
S-10074M1-1ML	100ug/ml in Methanol	1ML
S-10074U1-1ML	100ug/ml in Toluene	1ML

1-Naphthaleneacetic acid
CAS#86-87-3

PS: Solid
MF: C₁₂H₁₀O₂
MW: 186.21
MP: 134.5°C



Alternate Name: NAA(TM)

Part Number	Matrix	Size
N-10075-1G	Neat	1G
S-10075M1-1ML	100ug/ml in Methanol	1ML
S-10075T1-1ML	100ug/ml in T-Butyl methyl Ether	1ML

a-Naphthylthiourea
CAS#86-88-4

Chemical Name: a-Naphthyl thiocarbamide


PS: Solid
MF: C₁₁H₁₀N₂S
MW: 202.27
MP: 198°C


Alternate Name: ANTU

Part Number	Matrix	Size
N-10983-250MG	Neat	250MG
S-10983A1-1ML	100ug/ml in Acetonitrile	1ML
S-10983U1-1ML	100ug/ml in Toluene	1ML

3-Indoleacetic acid
CAS#87-51-4


PS: Solid
MF: C₁₀H₉NO₂


Alternate Name: Heteroauxin

Part Number	Matrix	Size
N-10727-250MG	Neat	250MG
S-10727A1-1ML	100ug/ml in Acetonitrile	1ML
S-10727T1-1ML	100ug/ml in T-Butyl methyl Ether	1ML

Pentachlorophenol
CAS#87-86-5


PS: Solid
MF: C₆HCl₅O
MW: 266.34
MP: 174°C


Alternate Name: Penta(TM)

Part Number	Matrix	Size
N-12831-500MG	Neat	500MG
S-12831M1-1ML	100ug/ml in Methanol	1ML
S-12831M1-5ML	100ug/ml in Methanol	5ML

2,4,6-Trichlorophenol
CAS#88-06-2


PS: Solid
MF: C₆H₃Cl₃O
MW: 197.45
MP: 69°C


Alternate Name: Dowicide 2S(TM)

Part Number	Matrix	Size
N-10657-250MG	Neat	250MG
S-10657M1-1ML	100ug/mL in Methanol	1ML
S-10657M1-5ML	100ug/mL in Methanol	5ML

2,3,5-Triiodobenzoic acid CAS#88-82-4

PS: Solid
MF: C₇H₃I₃O₂
MP: 225°C




Alternate Name: TIBA

Part Number	Matrix	Size
N-10601-250MG	Neat	250MG
S-10601A1-1ML	100ug/ml in Acetonitrile	1ML
S-10601T1-1ML	100ug/ml in t-Butyl methyl Ether	1ML

Dinoseb CAS#88-85-7

Chemical Name: 2-sec-Butyl-4,6-dinitrophenol


PS: Solid
MF: C₁₀H₁₂N₂O₅
MW: 240.2
MP: 38°C



Part Number	Matrix	Size
N-11786-100MG	Neat	100MG
S-11786A1-1ML	100ug/ml in Acetonitrile	1ML
S-11786B1-1ML	100ug/ml in Acetone	1ML

2,5-Dichloronitrobenzene CAS#89-61-2

PS: Solid
MP: 52°C




Alternate Name: 1,4-Dichloro-2-nitrobenzene

Part Number	Matrix	Size
N-10670-250MG	Neat	250MG
S-10670A1-1ML	100ug/ml in Acetonitrile	1ML
S-10670U1-1ML	100ug/ml in Toluene	1ML

2-Biphenylol CAS#90-43-7

Chemical Name: 2-Phenylphenol

PS: Solid
MF: C₁₂H₁₀O
MW: 170.21
MP: 59°C
FP: 124°C




Alternate Name: Dowcide 1

Part Number	Matrix	Size
N-10280-250MG	Neat	250MG
S-10280A1-1ML	100ug/ml in Acetonitrile	1ML
S-10280M5-1ML	2000ug/ml in Methanol	1ML

Ethoxyquin CAS#91-53-2

Chemical Name: 1,2-Dihydro-6-ethoxy-2,2,4-trimethyl quinoline


PS: Liquid
MF: C₁₄H₁₉NO
MW: 217.3



Part Number	Matrix	Size
N-11877-100MG	Neat	100MG

Biphenyl CAS#92-52-4

PS: Solid
MF: C₁₂H₁₀
MW: 154.2
MP: 69°C



Alternate Name: Phenyl benzene

Part Number	Matrix	Size
N-11207-500MG	Neat	500MG
S-11207J1-1ML	100ug/ml in Hexane	1ML
S-11207J1-5ML	100ug/ml in Hexane	5ML
S-11207M4-1ML	1000ug/ml in Methanol	1ML
S-11207M4-5ML	1000ug/ml in Methanol	5ML

N,N-Diallyl-2-chloroacetamide CAS#93-71-0

PS: Liquid
MF: C₈H₁₂ClNO
MW: 173.64
MP: 145°C




Alternate Name: Allidochlor

Part Number	Matrix	Size
N-12617-250MG	Neat	250MG
S-12617M1-1ML	100ug/ml in Methanol	1ML

Silvex CAS#93-72-1

Chemical Name: 2-[2,4,5-Trichlorophenoxy]propionic acid

PS: Solid
MF: C₉H₇Cl₃O₃
MW: 269.51
MP: 179°C




Alternate Name: Fenoprop (TM)

Part Number	Matrix	Size
N-13212-250MG	Neat	250MG
S-13212A1-1ML	100ug/ml in Acetonitrile	1ML
S-13212B1-1ML	100ug/ml in Acetone	1ML
S-13212B1-5ML	100ug/ml in Acetone	5ML

2,4,5-T CAS#93-76-5
Chemical Name: (2,4,5-trichlorophenoxy)acetic acid

PS: Solid


MF: C₈H₅Cl₃O₃
MW: 255.49

MP: 153°C


Part Number	Matrix	Size
N-10646-250MG	Neat	250MG
S-10646A1-1ML	100ug/ml in Acetonitrile	1ML
S-10646B1-1ML	100ug/ml in Acetone	1ML
S-10646B1-5ML	100ug/ml in Acetone	5ML
S-FC972S-1.2ML	100ug/ml in Methylene chloride (C13)	1.2ML

2,4,5-T isopropyl ester CAS#93-78-7
Chemical Name: 2,4,5-Trichlorophenoxy acetic acid-isopropyl ester

PS: Solid


MF: C₁₁H₁₁Cl₃O₃


Part Number	Matrix	Size
N-10650-100MG	Neat	100MG
S-10650A1-1ML	100ug/ml in Acetonitrile	1ML
S-10650T1-1ML	100ug/ml in t- Butyl methyl Ether	1ML

2,4,5-T n-butyl ester CAS#93-79-8
Chemical Name: 2,4,5-Trichlorophenoxy acetic acid - n-butyl ester

PS: Solid

MF: C₁₂H₁₃Cl₃O₃
MW: 311.6


MP: 46°C


Part Number	Matrix	Size
N-10651-250MG	Neat	250MG
S-10651A1-1ML	100ug/ml in Acetonitrile	1ML
S-10651J1-1ML	100ug/ml in Hexane	1ML

2,4-D isopropyl ester CAS#94-11-1
Chemical Name: Isopropyl 2,4-dichlorophenoxyacetate


PS: Liquid

MF: C₁₁H₁₂Cl₂O₃
MW: 263.12

MP: 5°C

Alternate Name: Esteron 45

Part Number	Matrix	Size
N-10617-250MG	Neat	250MG
S-10617A1-1ML	100ug/ml in Acetonitrile	1ML
S-10617U1-1ML	100ug/ml in Toluene	1ML


Phenyl mercuric benzoate CAS#94-43-9
PS: Solid

MF: C₁₃H₁₀HgO₂
MW: 199.21


Part Number	Matrix	Size
N-13012-100MG	Neat	100MG

4-Chloro-o-tolyloxyacetic acid CAS#94-74-6
PS: Solid

MF: C₉H₇ClO₃
MW: 200.6


MP: 118°C

Alternate Name: MCPA (TM)

Part Number	Matrix	Size
N-10818-250MG	Neat	250MG
S-10818B1-1ML	100ug/mL in Acetone	1ML
S-10818B1-5ML	100ug/mL in Acetone	5ML
S-10818M1-1ML	100ug/ml in Methanol	1ML


2,4-D CAS#94-75-7
Chemical Name: (2,4-dichlorophenoxy)acetic acid

PS: Solid

MF: C₈H₆Cl₂O₃
MW: 221.04

MP: 138°C



Part Number	Matrix	Size
N-10609-250MG	Neat	250MG
S-10609B1-1ML	100ug/ml in Acetone	1ML
S-10609M1-1ML	100ug/ml in Methanol	1ML

2,4-D sec-butyl ester solution CAS#94-75-7
PS: Liquid


Part Number	Matrix	Size
S-10621A4-1ML	1000ug/ml in Acetonitrile	1ML
S-10621A4-5ML	1000ug/ml in Acetonitrile	5ML
S-10621U1-1ML	100ug/ml in Toluene	1ML

2,4-D sec-butyl ester CAS#94-79-1

PS: Liquid
MF: C₁₂H₁₄Cl₂O₃
MW: 277.14
MP: 138°C




Part Number	Matrix	Size
N-10621-250MG	Neat	250MG

2,4-D butyl ester CAS#94-80-4

Chemical Name: Butyl(2,4-dichlorophenoxy)acetate

PS: Liquid
MF: C₁₂H₁₄Cl₂O₃
MW: 277.1
MP: 9°C




Alternate Name: Lirinox

Part Number	Matrix	Size
N-10611-250MG	Neat	250MG
S-10611A1-1ML	100ug/ml in Acetonitrile	1ML
S-10611U1-1ML	100ug/ml in Toluene	1ML

MCPB CAS#94-81-5

Chemical Name: 4-[2-Methyl-4-chlorophenoxy]butyric acid

PS: Solid
MF: C₁₁H₁₃ClO₃
MW: 228.69
MP: 95°C




Part Number	Matrix	Size
N-12357-100MG	Neat	100MG
S-12357A1-1ML	100ug/ml in Acetonitrile	1ML
S-12357T1-1ML	100ug/ml in t-Butyl methyl Ether	1ML

2,4-DB CAS#94-82-6

Chemical Name: 4-(2,4-Dichlorophenoxy)butyric acid

PS: Solid
MF: C₁₀H₁₀Cl₂O₃
MW: 249.09
MP: 118°C




Alternate Name: Embutox

Part Number	Matrix	Size
N-10622-250MG	Neat	250MG
S-10622A1-1ML	100ug/ml in Acetonitrile	1ML
S-10622B1-1ML	100ug/mL in Acetone	1ML
S-10622B1-5ML	100ug/mL in Acetone	5ML

2-Ethyl-1,3-hexanediol CAS#94-96-2

Chemical Name: 2-ethylhexane-1,3-diol

PS: Liquid
MF: C₈H₁₈O₂
MW: 146.22
MP: -40°C




Alternate Name: Ethohexadiol (TM)

Part Number	Matrix	Size
N-10346-500MG	Neat	500MG
S-10346B6-1ML	2500ug/ml in Acetone	1ML
S-10346B6-5ML	2500ug/ml in Acetone	5ML

Sulfallate CAS#95-06-7

Chemical Name: 2-Chloroallyl diethyldithiocarbamate

PS: Liquid
MF: C₈H₁₄ClNS₂
MW: 223.79




Alternate Name: Diethyldithiocarbamic acid 2-Chloroallyl ester

Part Number	Matrix	Size
N-13243-250MG	Neat	250MG
S-13243M1-1ML	100ug/ml in Methanol	1ML
S-13243M1-5ML	100ug/ml in Methanol	5ML

1,2-Dichlorobenzene CAS#95-50-1

Chemical Name: o-dichlorobenzene


PS: Liquid
MF: C₆H₄Cl₂
MW: 147
MP: -16.7°C



Part Number	Matrix	Size
N-10152-1G	Neat	1G
S-10152M1-1ML	100ug/ml in Methanol	1ML
S-10152M1-5ML	100ug/ml in Methanol	5ML

1,2-Dibromo-3-chloropropane CAS#96-12-8

PS: Liquid
MF: C₃H₅Br₂Cl
MW: 236.36
MP: 5°C



Alternate Name: DBCP

Part Number	Matrix	Size
N-10149-250MG	Neat	250MG
S-10149M1-1ML	100ug/ml in Methanol	1ML
S-10149M1-5ML	100ug/ml in Methanol	5ML

2-Imidazolidinethione


CAS#96-45-7

Chemical Name: N,N-Ethylenethiourea

PS: Solid

MF: C₃H₆N₂S

MW: 102.16

MP: 199°C


Part Number	Matrix	Size
N-10381-250MG	Neat	250MG
S-10381A7-1ML	5000ug/ml in Acetonitrile	1ML
S-10381A7-5ML	5000ug/ml in Acetonitrile	5ML
S-10381H1-1ML	100ug/ml in Ethyl acetate containing 0.1% of DTT	1ML
S-10381H1-5ML	100ug/ml in Ethyl acetate containing 0.1% of DTT	5ML


2,4-Dichlorophenyl benzenesulfonate

CAS#97-16-5

PS: Solid

MF: C₁₂H₈Cl₂O₃S

MW: 303.17

MP: 45°C

Alternate Name: Genite (TM)

Part Number	Matrix	Size
N-10628-250MG	Neat	250MG
S-10628A1-1ML	100ug/ml in Acetonitrile	1ML
S-10628U1-1ML	100ug/ml in Toluene	1ML


Dichlofenthion

CAS#97-17-6

Chemical Name: O-2,4-Dichlorophenyl-O,O-diethylphosphorothioate

PS: Liquid

MF: C₁₀H₁₃Cl₂O₃PS

MW: 315.17


Part Number	Matrix	Size
N-11659-250MG	Neat	250MG
S-11659A1-1ML	100ug/ml in Acetonitrile	1ML
S-11659J4-1ML	1000ug/ml in Hexane	1ML
S-11659J4-5ML	1000ug/ml in Hexane	5ML


Dichlorophen

CAS#97-23-4

Chemical Name: 2,2-Methylenebis[4-chlorophenol]

PS: Solid

MF: C₁₃H₁₀Cl₂O₂
MW: 269.13

MP: 177°C



Part Number	Matrix	Size
N-11667-250MG	Neat	250MG
S-11667M4-1ML	1000ug/ml in Methanol	1ML
S-11667M4-5ML	1000ug/ml in Methanol	5ML

4-Allyl-2-methoxyphenol

CAS#97-53-0

PS: Liquid

MF: C₁₀H₁₂O₂
MW: 164.2

MP: -9.2°C

Alternate Name: Eugenol

Part Number	Matrix	Size
N-10803-250MG	Neat	250MG
S-10803M1-1ML	100ug/ml in Methanol	1ML


Dicloran

CAS#99-30-9

Chemical Name: 2,6-Dichloro-4-nitroaniline

PS: Solid

MF: C₆H₄Cl₂N₂O₂
MW: 207.02

MP: 195°C

Alternate Name: DCNA

Part Number	Matrix	Size
N-11678-250MG	Neat	250MG
S-11678A1-1ML	100ug/ml in Acetonitrile	1ML
S-11678U4-1ML	1000ug/ml in Toluene	1ML


3,4-Dichloronitrobenzene

CAS#99-54-7

Chemical Name: 1,2-Dichloro-4-nitrobenzene

PS: Solid

MF: C₆H₃Cl₂NO₂
MW: 192

MP: 43°C


Part Number	Matrix	Size
N-10785-250MG	Neat	250MG
S-10785A1-1ML	100ug/ml in Acetonitrile	1ML
S-10785U1-1ML	100ug/ml in Toluene	1ML


Phenyl mercuric chloride(Technical)

CAS#100-56-1

PS: Solid

MF: C₆H₅ClHg


MW: 313.18

MP: 250°C

Alternate Name: Chlorophenyl mercury

Part Number	Matrix	Size
N-13013-250MG	Neat	250MG

Phenyl mercuric hydroxide(Technical)
CAS#100-57-2


PS: Solid
MF: C₆H₅HgOH
MW: 294.71
MP: 190°C


Alternate Name: Hydroxyphenyl mercury

Part Number	Matrix	Size
N-13014-250MG	Neat	250MG

Anilazine
CAS#101-05-3
Chemical Name: 2,4-Dichloro-6-o-chloroanilino-s-triazine


PS: Solid
MF: C₉H₅Cl₃N₄
MW: 275.52
MP: 159°C



Part Number	Matrix	Size
N-11075-250MG	Neat	250MG
S-11075A1-1ML	100ug/ml in Acetonitrile	1ML
S-11075U1-1ML	100ug/ml in Toluene	1ML
S-11075U1-5ML	100ug/ml in Toluene	5ML

2-(3-Chlorophenoxy)propionic acid
CAS#101-10-0


PS: Solid
MF: C₉H₉ClO₃
MW: 200.63
MP: 113°C


Alternate Name: Cloprop

Part Number	Matrix	Size
N-10254-100MG	Neat	100MG
S-10254A1-1ML	100ug/ml in Acetonitrile	1ML
S-10254T1-1ML	100ug/ml in t-Butyl methyl Ether	1ML

Chlorpropham
CAS#101-21-3
Chemical Name: Isopropyl-N-[m-chlorophenyl]carbamate


PS: Solid
MF: C₁₀H₁₂ClNO₂
MW: 213.68
MP: 40.7°C



Part Number	Matrix	Size
N-11458-250MG	Neat	250MG
S-11458A1-1ML	100ug/ml in Acetonitrile	1ML
S-11458T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-11458T1-5ML	100ug/ml in t-Butyl methyl ether	5ML

Barban
CAS#101-27-9
Chemical Name: 4-Chloro-2-butynyl-m-chlorocarbanilate


PS: Solid
MF: C₁₁H₉Cl₂NO₂
MW: 258.11
MP: 75°C



Part Number	Matrix	Size
N-11123-100MG	Neat	100MG
S-11123M1-1ML	100ug/ml in Methanol	1ML
S-11123M1-5ML	100ug/ml in Methanol	5ML

Fenuron
CAS#101-42-8
Chemical Name: N-Phenyl-N';N'-dimethylurea


PS: Solid
MF: C₉H₁₂N₂O
MW: 164.21
MP: 133°C



Part Number	Matrix	Size
N-11967-250MG	Neat	250MG
S-11967A4-1ML	1000ug/ml in Acetonitrile	1ML
S-11967A4-5ML	1000ug/ml in Acetonitrile	5ML
S-11967U1-1ML	100ug/ml in Toluene	1ML

Phenylmercury borate(Technical)
CAS#102-98-7


PS: Solid
MF: C₆H₅BHgO_{3.2}H
MW: 338.53



Part Number	Matrix	Size
N-13031-100MG	Neat	100MG

Chlorbenside
CAS#103-17-3
Chemical Name: p-Chlorobenzyl-p-chlorophenyl sulfide

PS: Solid
MF: C₁₃H₁₀Cl₂S
MW: 269.19
MP: 75°C



Part Number	Matrix	Size
N-11423-100MG	Neat	100MG
S-11423A1-1ML	100ug/ml in Acetonitrile	1ML
S-11423U1-1ML	100ug/ml in Toluene	1ML


Azobenzene CAS#103-33-3

Chemical Name: 1,2-diphenyldiazene

PS: Solid

MF: C₁₂H₁₀N₂
MW: 182.22

MP: 68°C

FP: 100

Alternate Name: Diphenyldiazene


Part Number	Matrix	Size
N-11110-500MG	Neat	500MG
S-11110M1-1ML	100ug/ml in Methanol	1ML
S-11110M1-5ML	100ug/ml in Methanol	5ML

Phenyl-2-thiourea CAS#103-85-5

PS: Solid

MF: C₇H₇N₂S

MW: 152.22

MP: 154°C

Alternate Name: Phenylthiocarbamide

Part Number	Matrix	Size
N-13020-250MG	Neat	250MG
S-13020A1-1ML	100ug/ml in Acetonitrile	1ML
S-13020U1-1ML	100ug/ml in Toluene	1ML


Cinnamaldehyde CAS#104-55-2

Chemical Name: 3-Phenyl-2-propenal

PS: Liquid

MF: C₉H₈O


MW: 132.15

MP: -7.5°C

Alternate Name: Cinnacure(TM)

Part Number	Matrix	Size
N-11468-500MG	Neat	500MG
S-11468A1-1ML	100ug/ml in Acetonitrile	1ML
S-11468U1-1ML	100ug/ml in Toluene	1ML

Phenylmercuric oleate(Technical) CAS#104-60-9

PS: Liquid

MF: C₂₄H₃₈HgO₂
MW: 559.21



Part Number	Matrix	Size
N-13030-100MG	Neat	100MG
S-13030A1-1ML	100ug/ml in Acetonitrile	1ML
S-13030U1-1ML	100ug/ml in Toluene	1ML

1,4-Dichlorobenzene CAS#106-46-7

Chemical Name: p-dichlorobenzene

PS: Solid

MF: C₆H₄Cl₂
MW: 147


MP: 53.09°C


Part Number	Matrix	Size
N-10216-1G	Neat	1G
S-10216M1-1ML	100ug/ml in Methanol	1ML
S-10216M1-5ML	100ug/ml in Methanol	5ML

1,2-Dibromoethane CAS#106-93-4

PS: Liquid

MF: C₂H₄Br₂
MW: 187.86

MP: 9.97°C

Alternate Name: Ethylene dibromide

Part Number	Matrix	Size
N-10150-1G	Neat	1G
S-10150M1-1ML	100ug/ml in Methanol	1ML

Ethylmercuric chloride(Technical) CAS#107-27-7

PS: Solid

MF: C₂H₅ClHg

MW: 265.11


MP: 193°C

Alternate Name: Chloroethylmercury

Part Number	Matrix	Size
N-11940-100MG	Neat	100MG

Tetraethylpyrophosphate CAS#107-49-3

PS: Liquid

MF: C₈H₂₀O₇P₂
MW: 290.22

Alternate Name: TEPP


Part Number	Matrix	Size
N-13543-500MG	Neat	500MG

2-Undecanone
CAS#112-12-9
PS: Liquid

MF: C₁₁H₂₂O

MW: 170.33


MP: 15°C

FP: 88°C

Alternate Name: Methyl nonyl ketone

Part Number	Matrix	Size
N-10520-1G	Neat	1G
S-10520N1-1ML	100ug/ml in Methanol:Water (90:10)	1ML

2-[2-Butoxyethoxy]ethyl thiocyanate
CAS#112-56-1
PS: Liquid


MF: C₉H₁₇NO₂S

MW: 203.32

Alternate Name: Lethane (TM)

Part Number	Matrix	Size
N-10260-250MG	Neat	250MG
S-10260A1-1ML	100ug/ml in Acetonitrile	1ML
S-10260J4-1ML	1000ug/ml in Hexane	1ML
S-10260J4-5ML	1000ug/ml in Hexane	5ML

MGK 264 (TM)
CAS#113-48-4
Chemical Name: N-(2-ETHYLHEXYL)-5-NORBORNENE-2,3-DICARBOXIMIDE

PS: Liquid


MF: C₁₇H₂₅NO₂
MW: 275.4


Part Number	Matrix	Size
N-12483-250MG	Neat	250MG
S-12483A1-1ML	100ug/ml in Acetonitrile	1ML
S-12483T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-12483T1-5ML	100ug/ml in t-Butyl methyl ether	5ML

Propoxur
CAS#114-26-1
Chemical Name: O-Isopropoxyphenylmethylcarbamate

PS: Solid

MF: C₁₁H₁₅NO₃
MW: 209.25

MP: 91.5°C

Alternate Name: Baygon


Part Number	Matrix	Size
N-11128-250MG	Neat	250MG
S-11128A1-1ML	100ug/ml in Acetonitrile	1ML

a,b-Endosulfan
CAS#115-29-7
Chemical Name: Hexachloro-norbornene-dimethanol cyclic sulfite

PS: Solid

MF: C₉H₆Cl₆O₃S

MW: 406.93


MP: 106°C


Part Number	Matrix	Size
N-10992-500MG	Neat	500MG
S-10992A1-1ML	100ug/ml in Acetonitrile	1ML
S-10992J1-1ML	100ug/ml in Hexane	1ML
S-10992J1-5ML	100ug/ml in Hexane	5ML

Isobornyl thiocyanacetate
CAS#115-31-1
PS: Liquid

MF: C₁₃H₁₉NO₂S

MW: 253.37

MP: -48°C

Alternate Name: Thanite


Part Number	Matrix	Size
N-12232-250MG	Neat	250MG
S-12232A1-1ML	100ug/ml in Acetonitrile	1ML
S-12232U1-1ML	100ug/ml in Toluene	1ML

Dicofol
CAS#115-32-2
Chemical Name: 4,4-Dichloro-a-[trichloromethyl]benzhydrol

PS: Solid

MF: C₁₄H₉Cl₅O

MW: 370.47

MP: 77.5°C

Alternate Name: Kelthane (T.M.)


Part Number	Matrix	Size
N-12290-100MG	Neat	100MG

Phosphan
CAS#115-78-6
Chemical Name: Tributyl[2,4-dichlorobenzyl]phosphonium chloride

PS: Solid

MF: C₁₉H₃₂Cl₂P.Cl


MW: 397.83

MP: 114°C


Part Number	Matrix	Size
N-13040-50MG	Neat	50MG
S-13040M1-1ML	100ug/ml in Methanol	1ML

2-Butyl-2-ethyl-1,3-propanediol
CAS#115-84-4


PS: Solid
MF: C₉H₂₀O₂
MW: 160.29
MP: 43.8°C


Alternate Name: 2-Ethyl-2-butyl-1,3-propanediol

Part Number	Matrix	Size
N-10306-1G	Neat	1G
S-10306M1-1ML	100ug/ml in Methanol	1ML

Fensulfthion
CAS#115-90-2
Chemical Name: O,O-Diethyl O-4-methylsulfinyl phenyl Phosphorothioate


PS: Liquid
MF: C₁₁H₁₇O₄PS₂
MW: 308.36



Part Number	Matrix	Size
N-11963-100MG	Neat	100MG
S-11963M1-1ML	100ug/ml in Methanol	1ML
S-11963M1-5ML	100ug/ml in Methanol	5ML

Cythioate
CAS#115-93-5
Chemical Name: O,O-Dimethyl-O-p-sulfamoylphenylphosphorothioate


PS: Solid
MF: C₈H₁₂NO₅PS₂
MW: 297.3
MP: 73°C



Part Number	Matrix	Size
N-11551-50MG	Neat	50MG
S-11551A1-1ML	100ug/ml in Acetonitrile	1ML
S-11551U1-1ML	100ug/ml in Toluene	1ML

Aldicarb
CAS#116-06-3
Chemical Name: 2-Methyl-2[methylthio]propionaldehyde O-[methylcarbamoyl] oxime


PS: Solid
MF: C₇H₁₄N₂O₂S
MW: 190.27
MP: 99°C



Part Number	Matrix	Size
N-11044-100MG	Neat	100MG
S-11044M1-1ML	100ug/ml in Methanol	1ML

Hexachloroacetone
CAS#116-16-5
Chemical Name: 2-Heptadecyl-2-imidazoline acetate


PS: Liquid
MF: C₃Cl₆O
MW: 264.75
MP: -2°C



Part Number	Matrix	Size
N-12158-500MG	Neat	500MG
S-12158A1-1ML	100ug/ml in Acetonitrile	1ML
S-12158J1-1ML	100ug/ml in Hexane	1ML

Tetradifon
CAS#116-29-0
Chemical Name: p-Chlorophenyl-2,4,5-trichlorophenylsulfone


PS: Solid
MF: C₁₂H₆Cl₄O₂S
MW: 356.06
MP: 146.5°C



Part Number	Matrix	Size
N-13539-250MG	Neat	250MG
S-13539A1-1ML	100ug/ml in Acetonitrile	1ML
S-13539B3-1ML	500ug/ml in Acetone	1ML
S-13539B3-5ML	500ug/ml in Acetone	5ML

1,2,4,5-Tetrachloro-3-nitrobenzene
CAS#117-18-0
Chemical Name: 2,3,5,6-Tetrachloronitrobenzene


PS: Solid
MF: C₆HCl₄NO₂
MW: 260.88
MP: 99.5°C


Alternate Name: Tecnazene

Part Number	Matrix	Size
N-10111-250MG	Neat	250MG
S-10111M1-1ML	100ug/ml in Methanol	1ML

Dichlone
CAS#117-80-6
Chemical Name: 2,3-Dichloro-1,4-naphthalenedione


PS: Solid
MF: C₁₀H₄Cl₂O₂
MW: 227.05
MP: 195°C


Alternate Name: 2,3-Dichloro-1,4-naphthoquinone

Part Number	Matrix	Size
N-11661-250MG	Neat	250MG
S-11661M1-1ML	100ug/ml in Methanol	1ML
S-11661U1-1ML	100ug/ml in Toluene	1ML
S-11661U1-5ML	100ug/ml in Toluene	5ML

Hexachlorobenzene
CAS#118-74-1


PS: Solid
MF: C₆Cl₆
MW: 284.78
MP: 231.8°C
FP: 242°C


Alternate Name: Perchlorobenzene

Part Number	Matrix	Size
N-12159-250MG	Neat	250MG
S-12159K4-1ML	1000ug/ml in Isooctane	1ML
S-12159M1-1ML	100ug/ml in Methanol	1ML
S-12159M1-5ML	100ug/ml in Methanol	5ML

Chloranil
CAS#118-75-2
Chemical Name: Tetrachloro-p-benzoquinone


PS: Solid
MF: C₆Cl₄O₂
MW: 245.88
MP: 292°C



Part Number	Matrix	Size
N-11420-250MG	Neat	250MG
S-11420M1-1ML	100ug/ml in Methanol	1ML
S-11420U1-1ML	100ug/ml in Toluene	1ML

Pyridaphenthion
CAS#119-12-0
Chemical Name: O-(1,6-Dihydro-6-oxo-1-phenylpyridazin-3-yl) O,O-diethyl phosphorothioate


PS: Solid
MF: C₁₄H₁₇N₂O₄PS
MW: 340.36


Alternate Name: Ofunack(TM)

Part Number	Matrix	Size
N-13154-250MG	Neat	250MG
S-13154A1-1ML	100ug/ml in Acetonitrile	1ML
S-13154U1-1ML	100ug/ml in Toluene	1ML

2-Naphthoxyacetic acid
CAS#120-23-0
Chemical Name: (2-naphthoxy)acetic acid


PS: Solid
MF: C₁₂H₁₀O₃
MW: 202.22
MP: 156°C


Alternate Name: Betoxon

Part Number	Matrix	Size
N-10437-250MG	Neat	250MG
S-10437M1-1ML	100ug/ml in Methanol	1ML
S-10437T1-1ML	100ug/ml in t-Butyl methyl Ether	1ML

o-Benzyl-p-chlorophenol
CAS#120-32-1


PS: Solid
MF: C₁₃H₁₁ClO
MW: 218.69
MP: 48.5°C


Alternate Name: Septiphen

Part Number	Matrix	Size
N-12676-100MG	Neat	100MG
S-12676M1-1ML	100ug/ml in Methanol	1ML

Dichlorprop
CAS#120-36-5
Chemical Name: 2-[2,4-Dichlorophenoxy]propionic acid


PS: Solid
MF: C₉H₈Cl₂O₃
MW: 235.07
MP: 117.5°C



Part Number	Matrix	Size
N-11671-250MG	Neat	250MG
S-11671B1-1ML	100ug/ml in Acetone	1ML
S-11671B1-5ML	100ug/ml in Acetone	5ML
S-11671M1-1ML	100ug/ml in Methanol	1ML
S-FC962S-1.2ML	100ug/ml in n-Nonane (C13)	1.2ML

Benzyl benzoate
CAS#120-51-4


PS: Liquid
MF: C₁₄H₁₂O₂
MW: 212.24
MP: 21°C


Alternate Name: Benzoic acid-benzyl ester

Part Number	Matrix	Size
N-11182-100MG	Neat	100MG
S-11182A1-1ML	100ug/ml in Acetonitrile	1ML
S-11182J5-1ML	2000ug/ml in Hexane	1ML
S-11182J7-1ML	5000ug/ml in Hexane	1ML
S-11182J7-5ML	5000ug/ml in Hexane	5ML

Sulfoxide
CAS#120-62-7
Chemical Name: 1-Methyl-2-(3,4-methylene dioxypheyl)ethyl octyl sulfoxide

PS: Liquid
MF: C₁₈H₂₈O₃S
MW: 324.49




Part Number	Matrix	Size
N-13257-100MG	Neat	100MG
S-13257M1-1ML	100ug/ml in Methanol	1ML
S-13257M1-5ML	100ug/ml in Methanol	5ML

Sulfanilic acid
CAS#121-57-3
Chemical Name: 4-Aminobenzenesulfonic acid

PS: Solid

MF: C₆H₇NO₃S

MW: 173.84


MP: 288°C


Part Number	Matrix	Size
N-13249-250MG	Neat	250MG

Malathion
CAS#121-75-5
Chemical Name: S-[1,2-Dicarbethoxyethyl]-O,O-dimethyl dithiophosphate

PS: Liquid

MF: C₁₀H₁₉O₆PS₂
MW: 330.36

MP: 2.9°C



Part Number	Matrix	Size
N-12346-100MG	Neat	100MG
S-12346A1-1ML	100ug/ml in Acetonitrile	1ML
S-12346U1-1ML	100ug/ml in Toluene	1ML
S-12346U1-5ML	100ug/ml in Toluene	5ML
S-FD2118S-1.2ML	100ug/ml in n-Nonane (Deuterated)	1.2ML

Sulfadimethoxine
CAS#122-11-2
Chemical Name: 4-Amino-N-(2,6-dimethoxy-4-pyrimidinyl)benzenesulfonamide

PS: Solid

MF: C₁₂H₁₄N₄O₄S

MW: 310.36

MP: 202°C

Alternate Name: Agribon(TM)


Part Number	Matrix	Size
N-13242-250MG	Neat	250MG
S-13242A1-1ML	100ug/ml in Acetonitrile	1ML

Fenitrothion
CAS#122-14-5
Chemical Name: O,O-Dimethyl-O-[4-nitro-m-tolyl]phosphorothioate

PS: Liquid

MF: C₉H₁₂NO₃PS

MW: 277.24


MP: 0.3°C


Part Number	Matrix	Size
N-11955-250MG	Neat	250MG
S-11955A1-1ML	100ug/ml in Acetonitrile	1ML
S-11955J4-1ML	1000ug/ml in Hexane	1ML

Simazine
CAS#122-34-9
Chemical Name: 2-Chloro-4,6-bis(ethylamino)-S-triazine

PS: Solid

MF: C₇H₁₂ClN₅
MW: 201.66


MP: 225°C


Part Number	Matrix	Size
N-13800-100MG	Neat	100MG
S-13800B4-1ML	1000ug/ml in Acetone	1ML
S-13800B4-5ML	1000ug/ml in Acetone	5ML
S-13800M1-1ML	100ug/ml in Methanol	1ML
S-13800M1-5ML	100ug/ml in Methanol	5ML
S-FC2073S-1.2ML	100ug/ml in Methanol (C13)	1.2ML

Propham
CAS#122-42-9
Chemical Name: Isopropyl-N-phenylcarbamate

PS: Solid


MF: C₁₀H₁₃NO₂
MW: 179.21

MP: 90°C


Part Number	Matrix	Size
N-13109-250MG	Neat	250MG
S-13109A4-1ML	1000ug/ml in Acetonitrile	1ML

Phenoxyacetic acid
CAS#122-59-8
PS: Solid


MF: C₈H₈O₃
MW: 152.16

MP: 98°C

Alternate Name: Phenylum (TM)

Part Number	Matrix	Size
N-13004-250MG	Neat	250MG
S-13004M1-1ML	100ug/ml in Methanol	1ML
S-13004T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Phenyl mercuric lactate
CAS#122-64-5
PS: Solid


MF: C₉H₁₀HgO₃
MW: 366.68

MP: 150°C

Alternate Name: [Lactoyloxy]phenylmercury

Part Number	Matrix	Size
N-13015-250MG	Neat	250MG

p-Chlorophenoxyacetic acid
CAS#122-88-3


PS: Solid
MF: C₈H₇ClO₃
MW: 186.6
MP: 156.5°C


Alternate Name: 4-CPA

Part Number	Matrix	Size
N-12763-250MG	Neat	250MG
S-12763M1-1ML	100ug/ml in Methanol	1ML
S-12763T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Maleic hydrazide
CAS#123-33-1
Chemical Name: 1,2-Dihydro-3,6-pyridazinedione


PS: Solid
MF: C₄H₄N₂O₂
MW: 112.09
MP: 298°C



Part Number	Matrix	Size
N-12350-250MG	Neat	250MG
S-12350M1-1ML	100ug/mL in Methanol	1ML

Methylarsonic acid
CAS#124-58-3


PS: Solid
MF: CH₃AsO₃
MW: 139.97
MP: 160.5°C



Part Number	Matrix	Size
N-12922-100MG	Neat	100MG

Demeton S
CAS#126-75-0
Chemical Name: O,O-Diethyl-S-[2-(ethylthio)ethyl]phosphorothioate


PS: Liquid
MF: C₈H₁₉O₃PS₂
MW: 258.36



Part Number	Matrix	Size
N-11582-100MG	Neat	100MG
S-11582A1-1ML	100ug/ml in Acetonitrile	1ML
S-11582U1-1ML	100ug/ml in Toluene	1ML

Dalapon Sodium
CAS#127-20-8
Chemical Name: 2,2-Dichloropropionic acid, sodium salt


PS: Solid
MF: C₃H₄Cl₂O₂Na
MP: 166.5°C


Alternate Name: Dowpon S(TM)

Part Number	Matrix	Size
N-11564-250MG	Neat	250MG

Sulfabenzamide
CAS#127-71-9
Chemical Name: N-[(4-Aminophenyl)sulfonyl]benzamide


PS: Solid
MF: C₁₃H₁₂N₂O₃S
MW: 276.3
MP: 181°C



Part Number	Matrix	Size
N-13240-250MG	Neat	250MG

Sulfamerazine
CAS#127-79-7
Chemical Name: 4-Amino-N-(4-methyl-2-pyrimidinyl)benzenesulfonamide


PS: Solid
MF: C₁₁H₁₂N₄O₂S
MW: 264.33



Part Number	Matrix	Size
N-13244-250MG	Neat	250MG

S 421
CAS#127-90-2
Chemical Name: BIS(2,3,3,3-TETRACHLOROPROPYL) ETHER

PS: Liquid
MF: C₆H₆Cl₈O




Part Number	Matrix	Size
N-13919-10MG	Neat	10MG

Busan 85 CAS#128-03-0

Chemical Name: Potassium dimethyldithiocarbamate

PS: Solid


MF: C₃H₇NS₂K


Part Number	Matrix	Size
N-11352-100MG	Neat	100MG
S-11352A1-1ML	100ug/ml in Acetonitrile	1ML
S-11352F1-1ML	100ug/ml in Water	1ML
S-11352M1-1ML	100ug/mL in Methanol	1ML

Dibam CAS#128-04-1

PS: Liquid

MF: C₃H₇NS₂Na


MP: 110°C


Part Number	Matrix	Size
N-12955-10MG	Neat	10MG

Dimethyl phthalate CAS#131-11-3

PS: Liquid

MF: C₁₀H₁₀O₄
MW: 194.19

MP: 5.5°C

Alternate Name: Phthalic acid - dimethyl ester


Part Number	Matrix	Size
N-11770-500MG	Neat	500MG
S-11770H1-1ML	100ug/ml in Ethyl acetate	1ML
S-11770H1-5ML	100ug/ml in Ethyl acetate	5ML
S-11770M1-1ML	100ug/ml in Methanol	1ML
S-11770M1-5ML	100ug/ml in Methanol	5ML

Meptyldinocap CAS#131-72-6

Chemical Name: 2-(1-methylheptyl)-4,6-dinitrophenyl (2E)-2-butenate

PS: Liquid

MF: C₁₈H₂₄N₂O₆
MW: 364.4


MP: -22.5°C


Part Number	Matrix	Size
N-13288-25MG	Neat	25MG
S-13288M1-1ML	100ug/mL in Methanol	1ML

N-1-Naphthylphthalamic acid CAS#132-66-1

PS: Solid

MF: C₁₈H₁₃NO₃
MW: 291.31

MP: 185°C

Alternate Name: Naptalam

Part Number	Matrix	Size
N-12507-250MG	Neat	250MG
S-12507A1-1ML	100ug/ml in Acetonitrile	1ML
S-12507T1-1ML	100ug/ml in t-Butyl methyl ether	1ML


Captan CAS#133-06-2

Chemical Name: N-[Trichloromethylthio]-4-cyclohexene-1,2-dicarboximide

PS: Solid

MF: C₉H₈Cl₃NO₂S

MW: 300.59

MP: 178°C


Part Number	Matrix	Size
N-11400-250MG	Neat	250MG
S-11400U1-1ML	100ug/ml in Toluene	1ML
S-11400U1-5ML	100ug/ml in Toluene	5ML


Folpet CAS#133-07-3

Chemical Name: N-[Trichloromethylthio]phthalimide

PS: Solid

MF: C₉H₄Cl₃NO₂S

MW: 296.56

MP: 177°C



Part Number	Matrix	Size
N-12007-250MG	Neat	250MG
S-12007M1-1ML	100ug/ml in Methanol	1ML
S-12007U1-1ML	100ug/ml in Toluene	1ML

3-Indolebutyric acid CAS#133-32-4

Chemical Name: 4-(Indol-3-yl)butyric acid

PS: Solid

MF: C₁₂H₁₃NO₂
MW: 203.24


MP: 124.5°C

Alternate Name: Seradix

Part Number	Matrix	Size
N-10728-100MG	Neat	100MG
S-10728A1-1ML	100ug/ml in Acetonitrile	1ML
S-10728T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Chloramben
CAS#133-90-4
Chemical Name: 3-amino-2,5-dichlorobenzoic acid

PS: Solid

MF: C₇H₅Cl₂NO₂
MW: 206.03


MP: 200°C


Part Number	Matrix	Size
N-11418-100MG	Neat	100MG
S-11418A1-1ML	100ug/ml in Acetonitrile	1ML
S-11418B1-1ML	100ug/mL in Acetone	1ML
S-11418B1-5ML	100ug/mL in Acetone	5ML

Methyl anthranilate
CAS#134-20-3
Chemical Name: 2-Aminobenzoic acid methyl ester

PS: Liquid

MF: C₈H₉NO₂
MW: 151.16

MP: 24.5°C



Part Number	Matrix	Size
N-12414-100MG	Neat	100MG
S-12414M1-1ML	100ug/ml in Methanol	1ML

N,N-Diethyl-m-tolamide
CAS#134-62-3
PS: Liquid

MF: C₁₂H₁₇NO


MW: 191.27

FP: 145°C

Alternate Name: DEET


Part Number	Matrix	Size
N-12618-250MG	Neat	250MG
S-12618A1-1ML	100ug/ml in Acetonitrile	1ML
S-12618B5-1ML	2000ug/ml in Acetone	1ML
S-12618B5-5ML	2000ug/ml in Acetone	5ML

Di-n-propyl isocinchomeronate
CAS#136-45-8
PS: Solid

MF: C₁₃H₁₇NO₄
MW: 251.31

Alternate Name: MGK Repellant 326

Part Number	Matrix	Size
N-12484-250MG	Neat	250MG


Ziram
CAS#137-30-4
Chemical Name: Dimethyl dithiocarbamic acid-zinc salt

PS: Solid

MF: C₆H₁₂N₂S₄Zn

MW: 305.83

MP: 246°C


FP: 215°C


Part Number	Matrix	Size
N-13761-250MG	Neat	250MG

Propazine
CAS#139-40-2
Chemical Name: 2-Chloro-4,6-bis[isopropylamino-s-triazine

PS: Solid

MF: C₉H₁₆ClN₅
MW: 230.09


MP: 229.7°C


Part Number	Matrix	Size
N-13108-100MG	Neat	100MG
S-13108A1-1ML	100ug/ml in Acetonitrile	1ML
S-13108T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-13108T1-5ML	100ug/ml in t-Butyl methyl ether	5ML
S-FC2236S-1.2ML	100ug/ml in n-Nonane (C13)	1.2ML

Monuron TCA
CAS#140-41-0
Chemical Name: 3-(4-chlorophenyl)-1,1-dimethylurea - trichloroacetic acid

PS: Solid


MF: C₉H₁₁ClN₂O₂·C₂HCl₃O₂
MW: 362.05

MP: 78°C


Part Number	Matrix	Size
N-12498-100MG	Neat	100MG

Fenaminosulf
CAS#140-56-7
PS: Solid

MF: C₈H₁₀N₃NaO₃S

MW: 251.26

Alternate Name: Dexon


Part Number	Matrix	Size
N-11586-250MG	Neat	250MG
S-11586M1-1ML	100ug/ml in Methanol	1ML

Aramite CAS#140-57-8
Chemical Name: 2(p-tert-Butylphenoxy)isopropyl-2-chloroethyl sulfite

PS: Liquid


MF: C₁₅H₂₃ClO₄S

MW: 334.87

MP: -31.7°C


Part Number	Matrix	Size
N-11085-10MG	Neat	10MG
S-11085J1-1ML	100ug/ml in Hexane	1ML


Cadmium succinate(Technical) CAS#141-00-4
PS: Solid

MF: C₄H₆O₄.Cd

Alternate Name: Succinic acid - cadmium salt

Part Number	Matrix	Size
N-11391-250MG	Neat	250MG

Dibutyl succinate CAS#141-03-7
PS: Liquid

MF: C₁₂H₂₂O₄
MW: 230.3


MP: -29°C

Alternate Name: Succinic acid - dibutyl ester

Part Number	Matrix	Size
N-11644-500MG	Neat	500MG
S-11644A1-1ML	100ug/ml in Acetonitrile	1ML
S-11644U1-1ML	100ug/ml in Toluene	1ML

Dicrotophos CAS#141-66-2
Chemical Name: O,O-Dimethyl O-(N,N-dimethylcarbamoyl-1-methylvinyl) phosphate

PS: Liquid

MF: C₈H₁₆NO₅P


MW: 237.19


Part Number	Matrix	Size
N-11680-100MG	Neat	100MG
S-11680U1-1ML	100ug/ml in Toluene	1ML
S-11680U1-5ML	100ug/ml in Toluene	5ML

Nabam CAS#142-59-6
Chemical Name: Disodium-ethylene-bis[dithiocarbamate]

PS: Solid

MF: C₄H₆N₂S_{4.2}Na

MW: 256.35



Part Number	Matrix	Size
N-12639-250MG	Neat	250MG
S-12639F1-1ML	100ug/ml in Water	1ML
S-12639U1-1ML	100ug/ml in Toluene	1ML

Chlordecone CAS#143-50-0
Chemical Name: 1,1a,3,3a,4,5,5a,5b,6-decachlorooctahydro-1,3,4-metheno-2H-cyclobuta[cd]pentalen-2-one

PS: Solid

MF: C₁₀Cl₁₀O


MW: 490.68

MP: 350°C


Part Number	Matrix	Size
N-12291-100MG	Neat	100MG
S-12291A1-1ML	100ug/ml in Acetonitrile	1ML
S-12291J1-1ML	100ug/ml in Hexane	1ML
S-12291J1-5ML	100ug/ml in Hexane	5ML
S-12291T4-1ML	1000ug/ml in t-Butyl methyl ether	1ML
S-12291T4-5ML	1000ug/ml in t-Butyl methyl ether	5ML

Disodium methyl arsonate hexahydrate CAS#144-21-8
PS: Solid

MF: CH₃AsNa₂O₃ · 6H₂O

MW: 291.9

Alternate Name: DSMA Hexahydrate


Part Number	Matrix	Size
N-11817-500MG	Neat	500MG
S-11817F1-1ML	100ug/ml in H2O	1ML

Sulfamethizole CAS#144-82-1
Chemical Name: 4-Amino-N-(5-methyl-1,3,4-thiadiazol-2-yl)benzenesulfonamide

PS: Solid

MF: C₉H₁₀N₄O₂S₂
MW: 270.34

MP: 210°C

FP: 400°C



Part Number	Matrix	Size
N-13246-250MG	Neat	250MG

Sulfapyridine
CAS#144-83-2
Chemical Name: 4-Amino-N-2-pyridinylbenzenesulfonamide

PS: Solid

MF: C₁₁H₁₁N₃O₂S

MW: 249.3


MP: 192°C


Part Number	Matrix	Size
N-13250-250MG	Neat	250MG

8-Hydroxyquinoline
CAS#148-24-3
PS: Solid

MF: C₉H₇NO

MW: 145.16

MP: 73°C

Alternate Name: Bioquin


Part Number	Matrix	Size
N-10965-250MG	Neat	250MG
S-10965M1-1ML	100ug/ml in Methanol	1ML

Thiabendazole
CAS#148-79-8
Chemical Name: 2-(4-Thiazolyl)benzimidazole

PS: Solid

MF: C₁₀H₇N₃S

MW: 201.25

MP: 304°C

Alternate Name: 2-(Thiazol-4-yl)benzimidazole

Part Number	Matrix	Size
N-13560-250MG	Neat	250MG
S-13560M4-1ML	1000ug/ml in Methanol	1ML

Tributylphosphoro-trithioite
CAS#150-50-5
PS: Liquid

MF: C₁₂H₂₇P₃S₃
MW: 298.5

Alternate Name: Folex (TM)


Part Number	Matrix	Size
N-13648-250MG	Neat	250MG

Monuron
CAS#150-68-5
Chemical Name: 3-(p-Chlorophenyl)-1,1-dimethylurea

PS: Solid

MF: C₉H₁₁ClN₂O


MW: 198.65

MP: 170.5°C


Part Number	Matrix	Size
N-12497-250MG	Neat	250MG
S-12497A4-1ML	1000ug/ml in Acetonitrile	1ML
S-12497M4-1ML	1000ug/ml in Methanol	1ML

Octamethylpyrophosphoramidate
CAS#152-16-9
PS: Liquid

MF: C₈H₂₄N₄O₃P₂
MW: 286.3

MP: 17°C

Alternate Name: Schradan


Part Number	Matrix	Size
N-12715-50MG	Neat	50MG

Isobenzan
CAS#297-78-9
Chemical Name: 1,3,4,5,6,7,8-octachloro-1,3,3a,4,7,7a-hexahydro-4,7-methanoisobenzofuran

PS: Solid

MF: C₉H₄Cl₈O

MW: 411.79

MP: 120°C



Part Number	Matrix	Size
N-12231-10MG	Neat	10MG
S-12231M1-1ML	100ug/ml in Methanol	1ML

Zinophos (TM)
CAS#297-97-2
Chemical Name: O,O-Diethyl-O-2-pyrazinylphosphorothioate

PS: Liquid

MF: C₈H₁₃N₂O₃PS

MW: 248.26

MP: -1.7°C

Alternate Name: Thionazin


Part Number	Matrix	Size
N-13760-50MG	Neat	50MG
S-13760A1-1ML	100ug/mL in Acetonitrile	1ML
S-13760J1-1ML	100ug/ml in Hexane	1ML
S-13760J1-5ML	100ug/ml in Hexane	5ML

Methyl parathion
CAS#298-00-0
Chemical Name: O,O-Dimethyl-O-p-nitrophenyl phosphorothioate

PS: Solid

MF: C₈H₁₀NO₅PS

MW: 263.21


MP: 35°C


Part Number	Matrix	Size
N-12452-100MG	Neat	100MG
S-12452A1-1ML	100ug/ml in Acetonitrile	1ML
S-12452U1-1ML	100ug/ml in Toluene	1ML

Phorate
CAS#298-02-2
Chemical Name: O,O-Diethyl-S-[(ethylthio)-methyl]phosphorodithioate

PS: Liquid


MF: C₇H₁₇O₂PS₃
MW: 260.34

MP: -15°C


Part Number	Matrix	Size
N-13035-100MG	Neat	100MG
S-13035A1-1ML	100ug/ml in Acetonitrile	1ML
S-13035J1-1ML	100ug/ml in Hexane	1ML

Demeton O
CAS#298-03-3
Chemical Name: O,O-Diethyl O-2-(ethylthio)ethyl phosphorothioate

PS: Liquid


MF: C₈H₁₉O₃PS₂C₈H₁₉O₃PS₂
MW: 258.34

Alternate Name: Di-septon

Part Number	Matrix	Size
N-11580-100MG	Neat	100MG
S-11580U1-1ML	100ug/ml in Toluene	1ML

Disulfoton
CAS#298-04-4
Chemical Name: O,O-Diethyl-S-[2-(ethylthio)ethyl]phosphorodithioate

PS: Liquid

MF: C₈H₁₉O₂PS₃
MW: 274.41

MP: -25°C



Part Number	Matrix	Size
N-11819-250MG	Neat	250MG
S-11819M1-1ML	100ug/ml in Methanol	1ML

Fenchlorphos
CAS#299-84-3
Chemical Name: O,O-Dimethyl-O-[2,4,5-trichlorophenyl]phosphorothioate

PS: Solid

MF: C₈H₆Cl₃O₃PS

MW: 321.57

MP: 41°C

Alternate Name: RONNEL


Part Number	Matrix	Size
N-11952-100MG	Neat	100MG
S-11952U1-1ML	100ug/ml in Toluene	1ML
S-11952U1-5ML	100ug/ml in Toluene	5ML

Ruelene (TM)
CAS#299-86-5
Chemical Name: 4-t-Butyl-chlorophenylmethyl-N-methyl phosphoramidate

PS: Solid

MF: C₁₂H₁₉ClNO₃P

MW: 291.71

MP: 60°C

Alternate Name: Crufomate


Part Number	Matrix	Size
N-13186-50MG	Neat	50MG
S-13186A1-1ML	100ug/ml in Acetonitrile	1ML
S-13186U1-1ML	100ug/ml in Toluene	1ML

Naled
CAS#300-76-5
Chemical Name: 1,2-Dibromo-2,2-dichloroethyl dimethyl phosphate

PS: Liquid

MF: C₄H₇Br₂Cl₂O₄P

MW: 380.79


MP: 26.5°C

Alternate Name: Dibrom

Part Number	Matrix	Size
N-12640-250MG	Neat	250MG
S-12640U1-1ML	100ug/ml in Toluene	1ML

Oxydemeton-methyl
CAS#301-12-2
Chemical Name: O,O-Dimethyl-S-[2-ethylsulfinyethyl]phosphorothioate

PS: Liquid

MF: C₆H₁₅O₄PS₂
MW: 246.3


MP: -20°C


Part Number	Matrix	Size
N-12741-50MG	Neat	50MG

Aldrin CAS#309-00-2

Chemical Name: Hexachlorohexahydro-endo-exo-dimethanonaphthalene

PS: Solid
MF: C₁₂H₈Cl₆
MW: 364.91
MP: 104°C




Part Number	Matrix	Size
N-11049-250MG	Neat	250MG
S-11049M1-1ML	100ug/ml in Methanol	1ML
S-11049M1-5ML	100ug/ml in Methanol	5ML
S-FC89S-1.2ML	100ug/mL in Nonane (C13)	1.2ML

Paraoxon CAS#311-45-5

Chemical Name: Diethyl-p-nitrophenylphosphate

PS: Liquid
MF: C₁₀H₁₄NO₆P
MW: 275.2



Alternate Name: Parathion O-analog

Part Number	Matrix	Size
N-12816-100MG	Neat	100MG
S-12816A1-1ML	100ug/ml in Acetonitrile	1ML
S-12816U1-1ML	100ug/ml in Toluene	1ML

Bromacil CAS#314-40-9

Chemical Name: 5-BROMO-3-SEC-BUTYL-6-METHYLURACIL

PS: Solid
MF: C₉H₁₃BrN₂O₂
MW: 261.15
MP: 157.5°C




Part Number	Matrix	Size
N-11330-250MG	Neat	250MG
S-11330A1-1ML	100ug/ml in Acetonitrile	1ML
S-11330T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Zectran CAS#315-18-4

Chemical Name: 4-(dimethylamino)-3,5-dimethylphenyl N-methylcarbamate

PS: Solid
MF: C₁₂H₁₈N₂O₂
MW: 222.29
MP: 85°C



Alternate Name: Zectran (Mexacarbate)

Part Number	Matrix	Size
N-13753-500MG	Neat	500MG
S-13753A4-1ML	1000ug/ml in Acetonitrile	1ML
S-13753U1-1ML	100ug/ml in Toluene	1ML

BHC (alpha isomer) CAS#319-84-6

Chemical Name: 1,2,3,4,5,6-Hexachlorocyclohexane (alpha isomer)

PS: Solid
MF: C₆H₆Cl₆
MW: 290.85
MP: 158°C




Part Number	Matrix	Size
N-11194-100MG	Neat	100MG
S-11194U1-1ML	100ug/ml in Toluene	1ML

BHC (beta isomer) CAS#319-85-7

Chemical Name: 1,2,3,4,5,6-Hexachlorocyclohexane (beta isomer)

PS: Solid
MF: C₆H₆Cl₆
MW: 290.83
MP: 309°C




Part Number	Matrix	Size
N-11195-100MG	Neat	100MG
S-11195U1-1ML	100ug/ml in Toluene	1ML
S-11195U1-5ML	100ug/ml in Toluene	5ML
S-FC103S-2X1.2ML	50ug/ml in Nonane (C13)	1.2ML

BHC (delta isomer) CAS#319-86-8

Chemical Name: 1,2,3,4,5,6-Hexachlorocyclohexane (delta isomer)

PS: Solid
MF: C₆H₆Cl₆
MW: 290.83
MP: 141.5°C




Part Number	Matrix	Size
N-11196-50MG	Neat	50MG
S-11196U1-1ML	100ug/mL in Toluene	1ML
S-11196U1-5ML	100ug/mL in Toluene	5ML
S-FC105S-1.2ML	100ug/ml in n-Nonane (C13)	1.2ML

Trichloronate CAS#327-98-0

Chemical Name: O-Ethyl O-(2,4,5-trichlorophenyl)ethylphosphonothioate

PS: Liquid
MF: C₁₀H₁₂Cl₃O₂PS
MW: 333.6




Part Number	Matrix	Size
N-13657-100MG	Neat	100MG
S-13657A1-1ML	100ug/ml in Acetonitrile	1ML
S-13657L1-1ML	100ug/ml in Isopropanol	1ML
S-13657L1-5ML	100ug/ml in Isopropanol	5ML

Diuron CAS#330-54-1

Chemical Name: 3-(3,4-Dichlorophenyl)-1,1-dimethylurea

PS: Solid
MF: C₉H₁₀Cl₂N₂O
MW: 233.1
MP: 158°C




Part Number	Matrix	Size
N-11827-250MG	Neat	250MG
S-11827A4-1ML	1000ug/ml in Acetonitrile	1ML
S-11827U1-1ML	100ug/ml in Toluene	1ML

Linuron CAS#330-55-2

Chemical Name: 3-(3,4-Dichlorophenyl)-1-methoxy-1-methylurea

PS: Solid
MF: C₉H₁₀Cl₂N₂O₂
MW: 249.1
MP: 93°C




Part Number	Matrix	Size
N-12322-250MG	Neat	250MG
S-12322A6-1ML	2500ug/ml in Acetonitrile	1ML
S-12322A6-5ML	2500ug/ml in Acetonitrile	5ML

Nicarbazine CAS#330-95-0

Chemical Name: N,N'-Bis(4-nitrophenyl)urea compd, with 4,6-dimethyl-2-pyrimidinol (1:1)

PS: Solid
MF: C₁₃H₁₀N₄O₅·C₆H₈N₂O
MP: 265°C



Alternate Name: Nicrazine

Part Number	Matrix	Size
N-12651-100MG	Neat	100MG

Diazinon CAS#333-41-5

Chemical Name: O,O-Diethyl-O-[2-isopropyl-4-methyl-6-pyridimyl]phosphorothioate

PS: Liquid
MF: C₁₂H₂₁N₂O₃PS
MW: 304.35




Part Number	Matrix	Size
N-11621-250MG	Neat	250MG
S-11621U1-1ML	100ug/ml in Toluene	1ML
S-FD2060S-1.2ML	100ug/ml in n-Nonane (Deuterated)	1.2ML

Z-Mevinphos CAS#338-45-4

Chemical Name: 2-Butenoic acid, 3-[(dimethoxyphosphinyl)oxy]-, methyl ester, (Z)-


PS: Liquid
MF: C₇H₁₃O₆P
MW: 224.16
MP: 6.9°C


Alternate Name: trans-Phosdrin

Part Number	Matrix	Size
N-12906-50MG	Neat	50MG

Capsaicin, Natural CAS#404-86-4

PS: Solid
MF: C₁₈H₂₇NO₃
MW: 305.42
MP: 65°C




Part Number	Matrix	Size
N-11398-100MG	Neat	100MG
S-11398A1-1ML	100ug/ml in Acetonitrile	1ML
S-11398U1-1ML	100ug/ml in Toluene	1ML

Isodrin CAS#465-73-6

Chemical Name: (1R,4S,5R,8S)-1,2,3,4,10,10-hexachloro-1,4,4a,5,8,8a-hexahydro-1,4:5,8-dimethanonaphthalene

PS: Solid
MF: C₁₂H₈Cl₆
MW: 364.9
MP: 240°C




Part Number	Matrix	Size
N-12248-100MG	Neat	100MG
S-12248M1-1ML	100ug/ml in Methanol	1ML
S-12248M1-5ML	100ug/ml in Methanol	5ML

Flurenol CAS#467-69-6

Chemical Name: 9-hydroxy-9H-fluorene-9-carboxylic acid

PS: Solid
MF: C₁₄H₁₀O₃
MW: 226.24
MP: 164°C


Alternate Name: 9-Hydroxyfluorene-9-carboxylic acid


Part Number	Matrix	Size
N-13815-250MG	Neat	250MG

Chlorfenvinphos CAS#470-90-6
Chemical Name: 2-Chloro-1-(2,4-dichlorophenyl)vinyl diethyl phosphate

PS: Liquid

MF: C₁₂H₁₄Cl₃O₄P

MW: 359.56


MP: -23°C


Part Number	Matrix	Size
N-11429-250MG	Neat	250MG
S-11429K1-1ML	100ug/ml in Isooctane	1ML

Binapacryl CAS#485-31-4
Chemical Name: 2-sec-Butyl-4,6-dinitrophenyl-3-methylbut-2-enoate

PS: Solid


MF: C₁₅H₁₈N₂O₆
MW: 322.31

MP: 66°C


Part Number	Matrix	Size
N-11204-100MG	Neat	100MG
S-11204A1-1ML	100ug/ml in Acetonitrile	1ML
S-11204T1-1ML	100ug/ml in t-Butyl methyl ether	1ML


p-Chloromandelic acid CAS#492-86-4
PS: Solid

MF: C₈H₇ClO₃
MW: 186.45

MP: 119°C


Part Number	Matrix	Size
N-12761-250MG	Neat	250MG
S-12761A1-1ML	100ug/ml in Acetonitrile	1ML
S-12761T1-1ML	100ug/ml in t-Butyl methyl ether	1ML


Diethyl dithiobis(thionoformate) CAS#502-55-6
PS: Solid

MF: C₆H₁₀O₂S₄

Alternate Name: EXD

Part Number	Matrix	Size
N-11696-250MG	Neat	250MG
S-11696A1-1ML	100ug/ml in Acetonitrile	1ML
S-11696K1-1ML	100ug/ml in Isooctane	1ML

4-Aminopyridine
CAS#504-24-5
PS: Solid

MF: C₅H₆N₂
MW: 94.12


MP: 158.5°C

Alternate Name: Amitrol

Part Number	Matrix	Size
N-10806-1G	Neat	1G
S-10806M1-1ML	100ug/ml in Methanol	1ML

Chlorobenzilate CAS#510-15-6
Chemical Name: Ethyl-4,4-dichlorobenzilate

PS: Solid

MF: C₁₆H₁₄Cl₂O₃
MW: 325.19

MP: 36°C



Part Number	Matrix	Size
N-11442-100MG	Neat	100MG
S-11442A1-1ML	100ug/ml in Acetonitrile	1ML
S-11442J1-1ML	100ug/ml in Hexane	1ML
S-11442J1-5ML	100ug/ml in Hexane	5ML

Kinetin CAS#525-79-1
Chemical Name: N-(2-furanylmethyl)-1H-purin-6-amine

PS: Solid

MF: C₁₀H₉N₅O

MW: 215.2


MP: 266°C

Alternate Name: 6-Furfurylamino-9H-purine

Part Number	Matrix	Size
N-12294-250MG	Neat	250MG

Butylmesityl oxide oxalate CAS#532-34-3
Chemical Name: Butyl-3,4-dihydro-2,2-dimethyl-4-oxo-1,2-H-pyran-[6-carboxylate]

PS: Liquid


MF: C₁₂H₁₈O₄
MW: 226.27

FP: 112°C


Part Number	Matrix	Size
N-11385-250MG	Neat	250MG

2,4-D ethyl ester
CAS#533-23-3
Chemical Name: Ethyl(2,4-dichlorophenoxy)acetate

PS: Liquid


MF: C₁₀H₁₀Cl₂O₃
MW: 249.1

Alternate Name: (2,4-Dichlorophenoxy)acetic acid - ethyl ester

Part Number	Matrix	Size
N-10614-250MG	Neat	250MG
S-10614A1-1ML	100ug/ml in Acetonitrile	1ML
S-10614J1-1ML	100ug/ml in Hexane	1ML

Dazomet
CAS#533-74-4
Chemical Name: Tetrahydro-3,5-dimethyl-2H-1,3,5-thiadiazine-2-thione

PS: Solid


MF: C₅H₁₀N₂S₂
MW: 162.28

MP: 106°C


Part Number	Matrix	Size
N-11566-250MG	Neat	250MG
S-11566M1-1ML	100ug/ml in Methanol	1ML

4,6-Dinitro-o-cresol (contains ~10% water)
CAS#534-52-1
PS: Solid

MF: C₇H₆N₂O₅
MW: 198.13


MP: 86.5°C

Alternate Name: DNOC

Part Number	Matrix	Size
N-10889-250MG	Neat	250MG
S-10889M4-1ML	1000ug/ml in Methanol	1ML

2-Furanacrylic acid
CAS#539-47-9
Chemical Name: 2-Fural acetic acid

PS: Solid

MF: C₇H₆O₃
MW: 138.13


MP: 142°C


Part Number	Matrix	Size
N-10363-100MG	Neat	100MG
S-10363A1-1ML	100ug/ml in Acetonitrile	1ML
S-10363T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

1,3-Dichloropropene
CAS#542-75-6
Chemical Name: 1,3-Dichloropropylene

PS: Liquid

MF: C₃H₄Cl₂
MW: 110.97

MP: -50°C

Alternate Name: Telone


Part Number	Matrix	Size
N-10193-500MG	Neat	500MG
S-10193M1-1ML	100ug/ml in Methanol	1ML

Neburon
CAS#555-37-3
Chemical Name: 1-n-Butyl-3-[3,4-dichlorophenyl]-1-methylurea

PS: Solid

MF: C₁₂H₁₆Cl₂N₂O

MW: 275.2


MP: 102°C


Part Number	Matrix	Size
N-12647-250MG	Neat	250MG
S-12647A4-1ML	1000ug/ml in Acetonitrile	1ML
S-12647A4-5ML	1000ug/ml in Acetonitrile	5ML

Glyodin(Technical)
CAS#556-22-9
Chemical Name: Ammonium ethylcarbamoylphosphonate

PS: Solid

MF: C₂₀H₄₀N₂.C₂H₄O₂
MW: 368.68

MP: 62°C


Part Number	Matrix	Size
N-12131-250MG	Neat	250MG

Ethanedial dioxime
CAS#557-30-2
PS: Solid

MF: C₂H₄N₂O₂
MW: 88.08

Alternate Name: Glyoxime


Part Number	Matrix	Size
N-11867-250MG	Neat	250MG
S-11867M1-1ML	100ug/ml in Methanol	1ML

Ethion CAS#563-12-2

Chemical Name: S,S-Methylene bis[O,O-diethyl]phosphorodithioate

PS: Liquid

MF: C₉H₂₂O₄P₂S₄
MW: 384.48

MP: -13°C



Part Number	Matrix	Size
N-11873-250MG	Neat	250MG
S-11873M1-1ML	100ug/ml in Methanol	1ML
S-11873M1-5ML	100ug/ml in Methanol	5ML

N-1-Naphthylacetamide CAS#575-36-0

PS: Solid

MF: C₁₂H₁₁NO

MW: 185.23


MP: 159°C


Part Number	Matrix	Size
N-12506-500MG	Neat	500MG
S-12506M1-1ML	100ug/ml in Methanol	1ML

2,6-D Acid CAS#575-90-6

Chemical Name: 2,6-(Dichlorophenoxy)acetic acid

PS: Solid


MF: C₈H₆Cl₂O₃


Part Number	Matrix	Size
N-12996-10MG	Neat	10MG

Allethrin CAS#584-79-2

Chemical Name: Bis(pentafluorophenyl)phenylphosphine

PS: Liquid

MF: C₁₉H₂₆O₃
MW: 302.4


Part Number	Matrix	Size
N-11051-100MG	Neat	100MG
S-11051X1-1ML	0.1 mg/mL in Methylene chloride	1ML
S-11051X1-1ML	0.1 mg/mL in Methylene chloride	1ML
S-11051X1-5ML	100ug/ml in Methylene chloride	5ML

Potassium cyanate(Technical) CAS#590-28-3

PS: Solid

MF: CHNO.K

MW: 81.12

MP: 315°C

Alternate Name: Potassium isocyanate


Part Number	Matrix	Size
N-13087-250MG	Neat	250MG

Perchloromethyl mercaptan CAS#594-42-3

Chemical Name: Trichloromethylsulfenyl chloride

PS: Liquid


MF: CCl₄S

MW: 185.88

Alternate Name: Clairsit(TM)

Part Number	Matrix	Size
N-12845-250MG	Neat	250MG
S-12845A1-1ML	100ug/ml in Acetonitrile	1ML
S-12845U1-1ML	100ug/ml in Toluene	1ML

1,1-Dichloro-1-nitroethane CAS#594-72-9

PS: Liquid

MF: C₂H₃Cl₂NO₂
MW: 143.9

Alternate Name: Ethide (TM)

Part Number	Matrix	Size
N-10122-1G	Neat	1G
S-10122M1-1ML	100ug/ml in Methanol	1ML
S-10122X1-1ML	100ug/ml in Methylene chloride	1ML

Bromoacetyl bromide CAS#598-21-0

PS: Liquid


MF: C₂H₂Br₂O


Part Number	Matrix	Size
N-11333-500MG	Neat	500MG
S-11333A1-1ML	100ug/ml in Acetonitrile	1ML
S-11333U1-1ML	100ug/ml in Toluene	1ML

BHC (mixed isomers)
CAS#608-73-1
Chemical Name: 1,2,3,4,5,6-Hexachlorocyclohexane (Mixed isomers)

PS: Solid


MF: C₆H₆Cl₆
MW: 290.8

MP: 65°C


Part Number	Matrix	Size
N-11198-250MG	Neat	250MG
S-11198A1-1ML	100ug/ml in Acetonitrile	1ML
S-11198U1-1ML	100ug/ml in Toluene	1ML

2,4-Dichloronitrobenzene
CAS#611-06-3
PS: Solid

MF: C₆H₃Cl₂NO₂
MW: 192


MP: 34°C

Alternate Name: asym-Nitro-m-dichlorobenzene

Part Number	Matrix	Size
N-10625-250MG	Neat	250MG
S-10625A1-1ML	100ug/ml in Acetonitrile	1ML
S-10625U1-1ML	100ug/ml in Toluene	1ML

o-Chlorophenoxy acetic acid
CAS#614-61-9
Chemical Name: 3,5-Dinitro-N4,N4-dipropylsulfanilamide

PS: Solid

MF: C₈H₇ClO₃
MW: 186.6


MP: 148.5°C


Part Number	Matrix	Size
N-12679-250MG	Neat	250MG
S-12679M1-1ML	100ug/ml in Methanol	1ML
S-12679T1-1ML	100ug/ml in t-Butyl methyl ether	1ML


3,5-Dichloronitrobenzene
CAS#618-62-2
Chemical Name: 1,3-Dichloro-5-nitrobenzene

PS: Solid


MF: C₆H₃Cl₂NO₂
MW: 192

MP: 65.4°C


Part Number	Matrix	Size
N-10792-250MG	Neat	250MG
S-10792A1-1ML	100ug/ml in Acetonitrile	1ML
S-10792U1-1ML	100ug/ml in Toluene	1ML

1,3-di-p-tolylurea
CAS#621-00-1
Chemical Name: 4,4'-Dimethyldiphenyl urea



Part Number	Matrix	Size
N-14078-100MG	Neat	100MG

2,3,4,5-Tetrachloroaniline Solution
CAS#634-83-3
PS: Liquid


Part Number	Matrix	Size
S-13990A1-1ML	100ug/ml in Acetonitrile	1ML

Thiometon
CAS#640-15-3
Chemical Name: S-(2-(Ethylthio)ethyl) O,O-dimethylphosphorodithionate

PS: Liquid


MF: C₆H₁₅O₂PS₃
MW: 246.35

Alternate Name: Dithiometon

Part Number	Matrix	Size
N-13570-50MG	Neat	50MG
S-13570U1-1ML	100ug/ml in Toluene	1ML

2-Fluoroacetamide
CAS#640-19-7
PS: Solid

MF: C₂H₄FNO

MW: 77.06


MP: 108°C

Alternate Name: Fluoroacetic acid amide

Part Number	Matrix	Size
N-10358-100MG	Neat	100MG
S-10358M1-1ML	100ug/ml in Methanol	1ML

Monothiono TEPP
CAS#645-78-3
Chemical Name: Tetraethyl Thionopyrophosphate

PS: Liquid

MF: C₈H₂₀O₆P₂S


MW: 306.28

Alternate Name: Tetraethyl monothiopyrophosphate

Part Number	Matrix	Size
N-12496-50MG	Neat	50MG
S-12496U1-1ML	100ug/ml in Toluene	1ML

Carbanolate
CAS#671-04-5
Chemical Name: 2-chloro-4,5-dimethylphenyl methylcarbamate

PS: Solid


MF: C₁₀H₁₂ClNO₂
MW: 213.68

MP: 130°C

Alternate Name: Banol

Part Number	Matrix	Size
N-11122-100MG	Neat	100MG
S-11122A1-1ML	100ug/ml in Acetonitrile	1ML
S-11122U1-1ML	100ug/ml in Toluene	1ML

BDMC
CAS#672-99-1
Chemical Name: 4-Bromo-3,5-dimethyl phenyl-N-methyl carbamate

PS: Solid


MF: C₁₀H₁₂BrNO₂
MW: 258.12


Part Number	Matrix	Size
N-11131-100MG	Neat	100MG
S-11131M1-1ML	100ug/ml in Methanol	1ML
S-11131M1-5ML	100ug/ml in Methanol	5ML

Simeton
CAS#673-04-1
Chemical Name: N2,N4-DIETHYL-6-METHOXY-1,3,5-TRIAZINE-2,4-DIAMINE

PS: Solid

MF: C₈H₁₅N₅O

MW: 197.24

Alternate Name: 2,4-BIS(ETHYLAMINO)-6-METHOXY-S-TRIAZINE

Part Number	Matrix	Size
N-12960-50MG	Neat	50MG


Propanil
CAS#709-98-8
Chemical Name: 3,4-Dichloropropionanilide

PS: Solid

MF: C₉H₉Cl₂NO

MW: 218.09

MP: 91°C

FP: 100°C

Alternate Name: Propanil


Part Number	Matrix	Size
N-10800-500MG	Neat	500MG
S-10800A1-1ML	100ug/ml in Acetonitrile	1ML
S-10800Z4-1ML	1000ug/ml in Acetonitrile:Acetone (90:10)	1ML
S-10800Z4-5ML	1000ug/ml in Acetonitrile:Acetone (90:10)	5ML

Sulfamethoxazole
CAS#723-46-6
Chemical Name: 4-Amino-N-(5-methyl-3-isoxazolyl)benzenesulfonamide

PS: Solid

MF: C₁₀H₁₁N₃O₃S

MW: 253.28


MP: 167°C


Part Number	Matrix	Size
N-13247-250MG	Neat	250MG

Tolyfluanid
CAS#731-27-1
Chemical Name: N-Dichlorofluoromethylthio-N',N'-dimethyl-N-p-tolylsulfamide

PS: Solid

MF: C₁₀H₁₃Cl₂FN₂O₂S₂
MW: 347.27


MP: 96°C

Alternate Name: Euparen M(TM)

Part Number	Matrix	Size
N-13584-250MG	Neat	250MG

Phosmet
CAS#732-11-6
Chemical Name: O,O-Dimethyl-S-phthalimido methyl phosphorodithioate

PS: Solid

MF: C₁₁H₁₂NO₄PS₂
MW: 317.32

MP: 72°C

Alternate Name: Imidan (T.M.)


Part Number	Matrix	Size
N-12207-100MG	Neat	100MG
S-12207A1-1ML	100ug/ml in Acetonitrile	1ML
S-12207J1-1ML	100ug/ml in Hexane	1ML
S-12207J1-5ML	100ug/ml in Hexane	5ML
S-FD2129S-1.2ML	100ug/ml in Acetonitrile (Deuterated)	1.2ML

Bensulide CAS#741-58-2

Chemical Name: O,O-bis(1-methylethyl) S-[2-[(phenylsulfonyl)amino]ethyl] phosphorodithioate

PS: Solid

MF: C₁₄H₂₄NO₄PS₃
MW: 397.52


MP: 34.4°C


Part Number	Matrix	Size
N-11140-250MG	Neat	250MG
S-11140M4-1ML	100ug/ml in Methanol	1ML
S-11140M4-5ML	1000ug/ml in Methanol	5ML

S-Ethyl dipropylthiocarbamate CAS#759-94-4

PS: Liquid

MF: C₉H₁₉NOS


MW: 189.3

Alternate Name: EPTC

Part Number	Matrix	Size
N-13188-500MG	Neat	500MG
S-13188A1-1ML	100ug/ml in Acetonitrile	1ML
S-13188T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-13188T1-5ML	100ug/ml in t-Butyl methyl ether	5ML

Carbophenothion CAS#786-19-6

Chemical Name: O,O-Diethyl-S-[p-chlorophenylthio]methyl phosphorodithioate

PS: Liquid


MF: C₁₁H₁₆ClO₂PS₃
MW: 342.87


Part Number	Matrix	Size
N-11408-250MG	Neat	250MG
S-11408M1-1ML	100ug/ml in Methanol	1ML

o,p'-DDT CAS#789-02-6

Chemical Name: 1,1,1-Trichloro-2-[o-chlorophenyl]-2-[p-chlorophenyl]ethane

PS: Solid

MF: C₁₄H₉Cl₅
MW: 354.5



Part Number	Matrix	Size
N-12708-50MG	Neat	50MG
S-12708M1-1ML	100ug/ml in Methanol	1ML
S-12708M1-5ML	100ug/ml in Methanol	5ML
S-FC2307S-1.2ML	100ug/mL in Nonane (C13)	1.2ML

3-Indolepropionic acid CAS#830-96-6

Chemical Name: 3-(3-Indolyl)propanoic acid

PS: Solid

MF: C₁₁H₁₁NO₂
MW: 189.21

MP: 134°C

Alternate Name: b-Indolepropionic acid

Part Number	Matrix	Size
N-10729-250MG	Neat	250MG
S-10729A1-1ML	100ug/ml in Acetonitrile	1ML
S-10729T1-1ML	100ug/ml in t-Butyl methyl ether	1ML


Methoprotrotryne CAS#841-06-5

Chemical Name: N2-Isopropyl-N4-(3-methoxypropyl)-6-methylthio-1,3,5-triazine-2,4-diamine

PS: Solid

MF: C₁₁H₂₁N₅OS

MW: 271.43


MP: 68°C


Part Number	Matrix	Size
N-12401-100MG	Neat	100MG
S-12401A1-1ML	100ug/ml in Acetonitrile	1ML
S-12401T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Demeton-O-methyl CAS#867-27-6

Chemical Name: O,O-Dimethyl O-[2-(ethylthio)ethyl] thiophosphate

PS: Liquid

MF: C₆H₁₅O₃PS₂
MW: 230.3


Part Number	Matrix	Size
N-11583-50MG	Neat	50MG


Terbutryne CAS#886-50-0

Chemical Name: N-[5-(1,1-Dimethylethyl)-1,3,4-thiadiazol-2-yl]-N,N'-dimethylurea

PS: Solid

MF: C₁₀H₁₉N₅S

MW: 241.4

MP: 104°C



Part Number	Matrix	Size
N-13514-250MG	Neat	250MG
S-13514A1-1ML	100ug/ml in Acetonitrile	1ML
S-13514T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-13514T1-5ML	100ug/ml in t-Butyl methyl ether	5ML

Phentinoacetate
CAS#900-95-8
Chemical Name: Triphenyltin acetate

PS: Solid

MF: C₂₀H₁₈O₂Sn

MW: 409


MP: 122°C

Alternate Name: Fentin acetate

Part Number	Matrix	Size
N-13007-250MG	Neat	250MG
S-13007A1-1ML	100ug/ml in Acetonitrile	1ML
S-13007U1-1ML	100ug/ml in Toluene	1ML

Metasystox I (TM)
CAS#919-86-8
Chemical Name: S-2-Ethylthioethyl O,O-dimethyl phosphorothioate


PS: Liquid

MF: C₆H₁₅O₃PS₂
MW: 230.3

MP: -76°C

Alternate Name: Demeton-S-methyl

Part Number	Matrix	Size
N-12384-50MG	Neat	50MG
S-12384A1-1ML	100ug/ml in Acetonitrile	1ML
S-12384U1-1ML	100ug/ml in Toluene	1ML


2-Phenoxypropionic acid
CAS#940-31-8
PS: Solid

MF: C₉H₁₀O₃
MW: 166.2


Part Number	Matrix	Size
N-10507-500MG	Neat	500MG
S-10507M1-1ML	100ug/ml in Methanol	1ML
S-10507T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Fonofos
CAS#944-22-9
Chemical Name: O-Ethyl-S-phenylethylphosphonodithioate

PS: Liquid


MF: C₁₀H₁₅OPS₂
MW: 246.33

Alternate Name: Dyfonate

Part Number	Matrix	Size
N-11842-100MG	Neat	100MG
S-11842A1-1ML	100ug/ml in Acetonitrile	1ML
S-11842J4-1ML	1000ug/ml in Hexane	1ML
S-11842J4-5ML	1000ug/ml in Hexane	5ML

Phosfolan
CAS#947-02-4
Chemical Name: Cyclic ethylene[diethoxyphosphinyl]dithioimidocarbonate

PS: Solid


MF: C₇H₁₄NO₃PS₂
MW: 255.28

MP: 36.5°C


Part Number	Matrix	Size
N-13038-100MG	Neat	100MG
S-13038A1-1ML	100ug/ml in Acetonitrile	1ML
S-13038U1-1ML	100ug/ml in Toluene	1ML

Mephosfolan
CAS#950-10-7
Chemical Name: 2-(Diethoxyphosphinylimino)-4-methyl-1,3-dithiolane

PS: Liquid


MF: C₈H₁₆NO₃PS₂
MW: 269.3


Part Number	Matrix	Size
N-12370-100MG	Neat	100MG
S-12370A1-1ML	100ug/ml in Acetonitrile	1ML
S-12370U1-1ML	100ug/ml in Toluene	1ML

Dimethyl-p-nitrophenylphosphate
CAS#950-35-6
PS: Liquid

MF: C₈H₁₀NO₆P

MW: 247.16


FP: 100°C

Alternate Name: Methyl paraoxon

Part Number	Matrix	Size
N-11775-100MG	Neat	100MG
S-11775A1-1ML	100ug/ml in Acetonitrile	1ML
S-11775T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-11775T1-5ML	100ug/ml in t-Butyl methyl ether	5ML

Methidathion
CAS#950-37-8
Chemical Name: S-[(5-methoxy-2-oxo-1,3,4-thiadiazol-3(2H)-yl)methyl] O,O-dimethyl phosphorodithioate

PS: Solid

MF: C₆H₁₁N₂O₄PS₃
MW: 302.33

MP: 39°C

Alternate Name: Supracide


Part Number	Matrix	Size
N-12397-250MG	Neat	250MG
S-12397A1-1ML	100ug/ml in Acetonitrile	1ML
S-12397U1-1ML	100ug/ml in Toluene	1ML

Diphenamid
CAS#957-51-7
Chemical Name: N,N-Dimethyl-2,2-diphenylacetamide

PS: Solid

MF: C₁₆H₁₇NO

MW: 239.34

MP: 134.5°C



Part Number	Matrix	Size
N-11794-250MG	Neat	250MG
S-11794A1-1ML	100ug/ml in Acetonitrile	1ML
S-11794T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-11794T1-5ML	100ug/ml in t-Butyl methyl ether	5ML

a-Endosulfan
CAS#959-98-8
Chemical Name: (3a,5ab,6a,9a,9ab)-6,7,8,9,10,10-hexachloro-1,5,5a,6,9,9a-hexahydro-6,9-methano-2,4,3-benzodioxathiepin 3-oxide

PS: Solid


MF: C₉H₆Cl₆O₃S

MW: 406.93

MP: 106°C


Part Number	Matrix	Size
N-10979-100MG	Neat	100MG
S-10979A1-1ML	100ug/ml in Acetonitrile	1ML
S-10979J1-1ML	100ug/ml in Hexane	1ML
S-10979J1-5ML	100ug/ml in Hexane	5ML
S-10979T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-10979T1-5ML	100ug/ml in t-Butyl methyl ether	5ML

Salannin
CAS#992-20-1
PS: Solid

MF: C₃₄H₄₄O₉
MW: 596.71


Part Number	Matrix	Size
N-12911-5MG	Neat	5MG

Chlormequat chloride
CAS#999-81-5
Chemical Name: [2-Chloroethyl]trimethylammonium chloride

PS: Solid

MF: C₅H₁₃Cl₂N

MW: 158.07

MP: 239°C



Part Number	Matrix	Size
N-11435-250MG	Neat	250MG
S-11435F1-1ML	100ug/ml in H2O	1ML

Desmetryn
CAS#1014-69-3
Chemical Name: N-methyl-N?-(1-methylethyl)-6-(methylthio)-1,3,5-triazine-2,4-diamine

PS: Solid

MF: C₈H₁₅N₅S

MW: 213.3


MP: 85°C


Part Number	Matrix	Size
N-12973-10MG	Neat	10MG

Simetryn
CAS#1014-70-6
Chemical Name: 2,4-Bis(ethylamino)-6-methylthio-1,3,5-triazine

PS: Solid

MF: C₈H₁₅N₅S

MW: 213.3

Alternate Name: Cymetrin


Part Number	Matrix	Size
N-13802-50MG	Neat	50MG
S-13802A1-1ML	100ug/ml in Acetonitrile	1ML
S-13802T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-13802T1-5ML	100ug/ml in t-Butyl methyl ether	5ML

Heptachlor epoxide (Isomer B)
CAS#1024-57-3
Chemical Name: 1,4,5,6,7,8-Heptachloro-2,3-epoxy-3a,4,7,7a-tetrahydro-4,7-endo-methanoindan

PS: Solid

MF: C₁₀H₅Cl₇O

MW: 389.4

MP: 160°C

Alternate Name: cis-Heptachlor epoxide


Part Number	Matrix	Size
N-12148-50MG	Neat	50MG
S-12148M1-1ML	100ug/ml in Methanol	1ML
S-12148M1-5ML	100ug/ml in Methanol	5ML
S-FC1015-1.2ML	100ug/ml in n-Nonane (C13)	1.2ML

Endosulfan sulfate
CAS#1031-07-8
Chemical Name: Hexachloro-hexahydro-methano-benzodioxathiepin-dioxide

PS: Solid

MF: C₉H₆Cl₆O₄S


MW: 422.95

MP: 181°C


Part Number	Matrix	Size
N-11851-50MG	Neat	50MG
S-11851A1-1ML	100ug/ml in Acetonitrile	1ML
S-11851J1-1ML	100ug/ml in Hexane	1ML
S-11851J1-5ML	100ug/ml in Hexane	5ML

Glyphosate CAS#1071-83-6
Chemical Name: N-(Phosphonomethyl)glycine


PS: Solid
MF: C₃H₈NO₅P
MW: 169.07
MP: 230°C



Part Number	Matrix	Size
N-12133-250MG	Neat	250MG
S-FCN1104S-1.2ML	100ug/ml in Water (C13)	1.2ML

Dichlofluanid CAS#1085-98-9
Chemical Name: 1,1-Dichloro-N-[(dimethylamino)sulfonyl]-1-fluoro-N-phenylmethanesulfenamide


PS: Solid
MF: C₉H₁₁Cl₂FN₂O₂S₂
MW: 333.21
MP: 105°C




Part Number	Matrix	Size
N-11660-250MG	Neat	250MG
S-11660U1-1ML	100ug/ml in Toluene	1ML

Omethoate CAS#1113-02-6
Chemical Name: O,O-Dimethyl-S-(N-methylcarbamoylmethyl)phosphorothioate

PS: Liquid
MF: C₅H₁₂NO₄PS
MW: 213.2
MP: -28°C




Part Number	Matrix	Size
N-12726-100MG	Neat	100MG
S-12726U1-1ML	100ug/ml in Toluene	1ML

Pebulate Solution CAS#1114-71-2
PS: Liquid


Part Number	Matrix	Size
S-13192A1-1ML	100ug/ml in Acetonitrile	1ML
S-13192T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Metolcarb CAS#1129-41-5
Chemical Name: 3-methylphenyl N-methylcarbamate


PS: Solid
MF: C₉H₁₁NO₂
MW: 165.2
MP: 76°C


Alternate Name: MTMC

Part Number	Matrix	Size
N-12479-100MG	Neat	100MG
S-12479A1-1ML	100ug/mL in Acetonitrile	1ML
S-12479U1-1ML	100ug/ml in Toluene	1ML

Cycloate CAS#1134-23-2
Chemical Name: Ethylcyclohexylethyl thiocarbamate


PS: Liquid
MF: C₁₁H₂₁NOS
MW: 215.39
MP: 11.5°C



Part Number	Matrix	Size
N-11522-250MG	Neat	250MG
S-11522A1-1ML	100ug/ml in Acetonitrile	1ML
S-11522T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-11522T1-5ML	100ug/ml in t-Butyl methyl ether	5ML

(1R)-trans-Tetramethrin CAS#1166-46-7
Chemical Name: (1,3-dioxo-4,5,6,7-tetrahydroisindol-2-yl)methyl (1R,3R)-2,2-dimethyl-3-(2-methylprop-1-enyl)cyclopropane-1-carboxylate


PS: Liquid
MF: C₁₉H₂₅NO₄
MW: 331.4


Alternate Name: d-Tetramethrin

Part Number	Matrix	Size
N-13994-100MG	Neat	100MG
S-13994A1-1ML	100ug/mL in Acetonitrile	1ML

2,6-Dichlorobenzonitrile CAS#1194-65-6


PS: Solid
MF: C₇H₃Cl₂N
MW: 172.01
MP: 144°C


Alternate Name: Dichlobenil

Part Number	Matrix	Size
N-10688-250MG	Neat	250MG
S-10688M1-1ML	100ug/ml in Methanol	1ML

6-Benzylaminopurine
CAS#1214-39-7


PS: Solid
MF: C₁₂H₁₁N₅
MW: 225.28
MP: 234°C


Alternate Name: BAP

Part Number	Matrix	Size
N-10956-250MG	Neat	250MG
S-10956A1-1ML	100ug/ml in Acetonitrile	1ML

Flurenol-methyl
CAS#1216-44-0
Chemical Name: Methyl 9-hydroxyfluorene-9-carboxylate


PS: Solid
MP: 161°C



Part Number	Matrix	Size
N-13813-250MG	Neat	250MG

Copper oxychloride(Technical)
CAS#1332-65-6

PS: Solid
MP: 240°C



Part Number	Matrix	Size
N-11502-250MG	Neat	250MG

Copper Naphthenates(Technical)
CAS#1338-02-9


PS: Liquid
MF: UNKNOWN
MW: 405.86

Alternate Name: Naphthenic acid, copper(2+) salt

Part Number	Matrix	Size
N-11501-250MG	Neat	250MG

Neomycin trisulfate hydrate(Technical)
CAS#1405-10-3


PS: Solid
MW: 908.88


Alternate Name: Neomycin trisulfate salt hydrate

Part Number	Matrix	Size
N-12649-250MG	Neat	250MG

Pentachlorophenyl acetate
CAS#1441-02-7


PS: Solid
MF: C₈H₃Cl₅O₂
MW: 308.37


Alternate Name: Rabcon

Part Number	Matrix	Size
N-13824-100MG	Neat	100MG

Carbofuran
CAS#1563-66-2
Chemical Name: 2,3-Dihydro-2,2-dimethyl-benzofuran-7-yl methylcarbamate


PS: Solid
MF: C₁₂H₁₅NO₃
MW: 221.26



Part Number	Matrix	Size
N-11405-250MG	Neat	250MG
S-11405M1-1ML	100ug/ml in Methanol	1ML
S-FC20065-1.2ML	100ug/ml in 1,4-Dioxane (C13)	1.2ML

Trifluralin
CAS#1582-09-8

PS: Solid
MF: C₁₃H₁₆F₃N₃O₄
MW: 335.28
MP: 46°C


Alternate Name: TCA


Part Number	Matrix	Size
N-13689-250MG	Neat	250MG
S-13689A1-1ML	100ug/ml in Acetonitrile	1ML
S-13689U1-1ML	100ug/mL in Toluene	1ML
S-13689U1-5ML	100ug/mL in Toluene	5ML
S-FD2140S-1.2ML	100ug/ml in n-Nonane (Deuterated)	1.2ML

Dodemorph
CAS#1593-77-7
Chemical Name: 4-Cyclododecyl-2,6-dimethylmorpholine

PS: Solid

MF: C₁₈H₃₅NO

MW: 281.5


MP: 71°C


Part Number	Matrix	Size
N-11839-250MG	Neat	250MG
S-11839M1-1ML	100ug/ml in Methanol	1ML

Daminozide
CAS#1596-84-5
Chemical Name: Succinic acid - 2,2-dimethylhydrazide

PS: Solid

MF: C₆H₁₂N₂O₃
MW: 160.17

MP: 154°C



Part Number	Matrix	Size
N-11565-250MG	Neat	250MG
S-11565A1-1ML	100ug/ml in Acetonitrile	1ML
S-11565T1-1ML	100ug/ml t-Butyl methyl ether	1ML

Atraton
CAS#1610-17-9
Chemical Name: 2-[Ethylamino]-4-[isopropylamino]-6-methoxy-s-triazine

PS: Solid

MF: C₉H₁₇N₅O

MW: 211.13

MP: 95°C

Alternate Name: Gestamin


Part Number	Matrix	Size
N-11105-250MG	Neat	250MG
S-11105A1-1ML	100ug/ml in Acetonitrile	1ML
S-11105T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-11105T1-5ML	100ug/ml in t-Butyl methyl ether	5ML

Prometon
CAS#1610-18-0
Chemical Name: 2-Methoxy-4,6-bis[isopropylamino]-s-triazine

PS: Solid

MF: C₁₀H₁₉N₅O

MW: 225.3


MP: 91°C


Part Number	Matrix	Size
N-13102-250MG	Neat	250MG
S-13102A1-1ML	100ug/ml in Acetonitrile	1ML
S-13102T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-13102T1-5ML	100ug/ml in t-Butyl methyl ether	5ML

Aldicarb sulfoxide
CAS#1646-87-3
Chemical Name: 2-Methyl-2-(methylsulfinyl)propanol O-[(methylamino)carbonyl] oxime

PS: Solid

MF: C₇H₁₄N₂O₃S


MW: 206.29


Part Number	Matrix	Size
N-11046-10MG	Neat	10MG
S-11046A1-1ML	100ug/ml in Acetonitrile	1ML

Aldicarb sulfone
CAS#1646-88-4
Chemical Name: 2-methyl-2-(methylsulfonyl)propanal O-[(methylamino)carbonyl] oxime

PS: Solid

MF: C₇H₁₄N₂O₄S

MW: 222.29



Part Number	Matrix	Size
N-11045-10MG	Neat	10MG
S-11045M1-1ML	100ug/ml in Methanol	1ML

3,5-Diiodo-4-hydroxybenzonitrile
CAS#1689-83-4
Chemical Name: 4-Hydroxy-3,5-diiodobenzonitrile

PS: Solid

MF: C₇H₃I₂NO

MW: 370.91

MP: 212°C



Part Number	Matrix	Size
N-10794-100MG	Neat	100MG
S-10794M1-1ML	100ug/ml in Methanol	1ML

Bromoxynil
CAS#1689-84-5
Chemical Name: 3,5-Dibromo-4-hydroxybenzonitrile

PS: Solid

MF: C₇H₃Br₂NO

MW: 276.9


MP: 194°C


Part Number	Matrix	Size
N-11345-250MG	Neat	250MG
S-11345M1-1ML	100ug/ml in Methanol	1ML
S-FC2093S-1.2ML	50ug/ml in n-Nonane (C13)	1.2ML

Bromoxynil octanoic acid ester
CAS#1689-99-2
Chemical Name: 2,6-Dibromo-4-cyanophenyl octanoate

PS: Solid

MF: C₁₅H₁₇Br₂NO₂
MW: 403.15

MP: 45°C



Part Number	Matrix	Size
N-11346-250MG	Neat	250MG
S-11346A1-1ML	100ug/ml in Acetonitrile	1ML
S-11346K4-1ML	1000ug/ml in Isooctane	1ML
S-11346K4-5ML	1000ug/ml in Isooctane	5ML

Pyrazon
CAS#1698-60-8
Chemical Name: 5-Amino-4-chloro-2-phenyl-3(2H)-pyridazone

PS: Solid

MF: C₁₀H₈ClN₃O

MW: 221.6

MP: 206°C



Part Number	Matrix	Size
N-13145-250MG	Neat	250MG
S-13145M1-1ML	100ug/ml in Methanol	1ML

Clopyralid
CAS#1702-17-6
Chemical Name: 3,6-DICHLORO-2-PYRIDINECARBOXYLIC ACID

PS: Solid

MF: C₆H₃Cl₂NO₂
MW: 192


MP: 151°C

FP: 100°C

Alternate Name: LONTREL™

Part Number	Matrix	Size
N-12324-250MG	Neat	250MG
S-12324A1-1ML	100ug/ml in Acetonitrile	1ML
S-12324T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

2,4-D isobutyl ester
CAS#1713-15-1
Chemical Name: (2,4-Dichlorophenoxy)acetic acid isobutyl ester


PS: Liquid

MF: C₁₂H₁₄Cl₂O₃
MW: 277.16

Alternate Name: Isobutyl 2,4-dichlorophenoxyacetate

Part Number	Matrix	Size
N-10532-250MG	Neat	250MG
S-10532A1-1ML	100ug/ml in Acetonitrile	1ML
S-10532U1-1ML	100ug/ml in Toluene	1ML

Fenthion-ethyl
CAS#1716-09-2
Chemical Name: O,O-Diethyl O-(3-methyl-4-(methylthio)phenyl)phosphorothioate

PS: Liquid


MF: C₁₂H₁₉O₃PS₂
MW: 306.4


Part Number	Matrix	Size
N-11965-500MG	Neat	500MG
S-11965A1-1ML	100ug/ml in Acetonitrile	1ML
S-11965U1-1ML	100ug/ml in Toluene	1ML

Monolinuron
CAS#1746-81-2
Chemical Name: 3-(4-Chlorophenyl)-1-methoxy-1-methylurea

PS: Solid

MF: C₉H₁₁ClN₂O₂
MW: 214.6


MP: 80°C

Alternate Name: Aresin(TM)

Part Number	Matrix	Size
N-12494-250MG	Neat	250MG
S-12494A1-1ML	100ug/ml in Acetonitrile	1ML
S-12494U1-1ML	100ug/ml in Toluene	1ML

Akton (TM)
CAS#1757-18-2
Chemical Name: O,O-Diethyl O-[2-chloro-1-[2,5-dichlorophenyl]vinyl]phosphorothioate

PS: Solid

MF: C₁₂H₁₄Cl₃O₃PS


MW: 375.64

Alternate Name: axiom

Part Number	Matrix	Size
N-11042-100MG	Neat	100MG
S-11042A1-1ML	100ug/ml in Acetonitrile	1ML
S-11042U1-1ML	100ug/ml in Toluene	1ML

Mesoprazine
CAS#1824-09-5
Chemical Name: 6-chloro-N-(3-methoxypropyl)-N7-(1-methylethyl)-1,3,5-triazine-2,4-diamine

PS: Solid


MF: C₁₀H₁₈ClN₅O

MW: 259.7


Part Number	Matrix	Size
N-13296-50MG	Neat	50MG
S-13296A1-1ML	100ug/ml in Acetonitrile	1ML

Pentachloroanisole
CAS#1825-21-4
Chemical Name: 2,3,4,5,6-Pentachloroanisole


PS: Solid
MF: C₇H₃Cl₅O
MW: 280.37
MP: 108°C



Part Number	Matrix	Size
N-12826-100MG	Neat	100MG
S-12826M1-1ML	100ug/ml in Methanol	1ML

Nitrofen
CAS#1836-75-5
Chemical Name: 2,4-Dichlorophenyl-p-nitrophenyl ether


PS: Solid
MF: C₁₂H₇Cl₂NO₃
MW: 284.1
MP: 70°C



Part Number	Matrix	Size
N-12663-100MG	Neat	100MG
S-12663M1-1ML	100ug/ml in Methanol	1ML

Chlorthal-dimethyl
CAS#1861-32-1
Chemical Name: Dimethyl-2,3,5,6-tetrachloroterephthalate


PS: Solid
MF: C₁₀H₆Cl₄O₄
MW: 332
MP: 155°C


Alternate Name: DCPA

Part Number	Matrix	Size
N-11462-250MG	Neat	250MG
S-11462B1-1ML	100ug/ml in Acetone	1ML
S-11462B1-5ML	100ug/ml in Acetone	5ML
S-11462M1-1ML	100ug/ml in Methanol	1ML

Benfluralin
CAS#1861-40-1
Chemical Name: N-Butyl-N-ethyl-a,a,a-trifluoro-2,6-dinitro-p-toluidine


PS: Solid
MF: C₁₃H₁₆F₃N₃O₄
MW: 335.32
MP: 65°C



Part Number	Matrix	Size
N-11136-250MG	Neat	250MG
S-11136A1-1ML	100ug/ml in Acetonitrile	1ML
S-11136J4-1ML	1000ug/ml in Hexane	1ML
S-11136J4-5ML	1000ug/ml in Hexane	5ML

p-Bromophenoxy acetic acid
CAS#1878-91-7


PS: Solid
MF: C₈H₇BrO₃



Part Number	Matrix	Size
N-12757-250MG	Neat	250MG
S-12757A1-1ML	100ug/ml in Acetonitrile	1ML
S-12757T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Chloroxynil
CAS#1891-95-8
Chemical Name: 3,5-Dichloro-4-hydroxybenzonitrile


PS: Solid
MF: C₇H₃Cl₂NO
MW: 188.01



Part Number	Matrix	Size
N-11457-500MG	Neat	500MG
S-11457M1-1ML	100ug/ml in Methanol	1ML

Chlorothalonil
CAS#1897-45-6
Chemical Name: Tetrachloroisophthalonitrile


PS: Solid
MF: C₆Cl₄N₂
MW: 265.89
MP: 250°C



Part Number	Matrix	Size
N-11454-250MG	Neat	250MG
S-11454A1-1ML	100ug/ml in Acetonitrile	1ML
S-11454T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Paraquat dichloride hydrate
CAS#1910-42-5
Chemical Name: 1,1'-dimethyl-4,4'-bipyridinium dichloride

PS: Solid
MF: C₁₂H₁₄Cl₂N₂
MW: 257.2


Alternate Name: 1,1-Dimethyl-4,4'-bipyridinium dichloride


Part Number	Matrix	Size
N-12818-500MG	Neat	500MG
S-12818F1-1ML	100ug/ml in Water	1ML
S-12818F17-1ML	1770ug/mL in Water	1ML
S-12818F17-5ML	1770ug/mL in Water	5ML

Atrazine CAS#1912-24-9

Chemical Name: 2-Chloro-4-ethylamino-6-isopropylamino-s-triazine

PS: Solid

MF: C₈H₁₄ClN₅
MW: 215.69

MP: 175°C



Part Number	Matrix	Size
N-11106-250MG	Neat	250MG

Trietazine CAS#1912-26-1

Chemical Name: 6-Chloro-N,N,N'-triethyl-1,3,5-triazine-2,4-diamine

PS: Solid

MF: C₉H₁₆ClN₅
MW: 229.75


MP: 102°C


Part Number	Matrix	Size
N-13666-100MG	Neat	100MG
S-13666A1-1ML	100ug/ml in Acetonitrile	1ML
S-13666X1-1ML	100ug/ml in Methylene chloride	1ML

Dicamba CAS#1918-00-9

Chemical Name: 2-Methoxy-3,6-dichlorobenzoic acid

PS: Solid

MF: C₈H₆Cl₂O₃
MW: 221.04

Alternate Name: Banvel (TM)


Part Number	Matrix	Size
N-11656-250MG	Neat	250MG
S-11656A1-1ML	100ug/ml in Acetonitrile	1ML
S-11656B1-1ML	100ug/ml in Acetone	1ML
S-11656B1-5ML	100ug/ml in Acetone	5ML

Chlorthiamid CAS#1918-13-4

Chemical Name: 2,6-dichlorothiobenzamide

PS: Solid

MF: C₇H₂Cl₂NS

MW: 206.09

Alternate Name: 2,6-DICHLOROBENZENECARBOTHIOAMIDE

Part Number	Matrix	Size
N-12992-10MG	Neat	10MG


Propachlor CAS#1918-16-7

Chemical Name: 2-Chloro-N-isopropylacetanilide

PS: Solid

MF: C₁₁H₁₄ClNO

MW: 211.71

MP: 77°C


Part Number	Matrix	Size
N-13104-250MG	Neat	250MG


S-13104A1-1ML	100ug/ml in Acetonitrile	1ML
S-13104T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-13104T1-5ML	100ug/ml in t-Butyl methyl ether	5ML

SWEP CAS#1918-18-9

Chemical Name: 3,4-Dichlorocarbanilic acid methyl ester

PS: Solid

MF: C₈H₇Cl₂NO₂
MW: 220.06

MP: 113°C



Part Number	Matrix	Size
N-13260-250MG	Neat	250MG

S-13260A4-1ML	1000ug/ml in Acetonitrile	1ML
S-13260A4-5ML	1000ug/ml in Acetonitrile	5ML
S-13260U1-1ML	100ug/ml in Toluene	1ML

(2,4,5-Trichlorophenoxy)acetic acid methyl ester CAS#1928-37-6

Chemical Name: Methyl 2-(2,4,5-trichlorophenoxy)acetate

PS: Solid

MF: C₉H₇Cl₃O₃
MW: 269.51

Alternate Name: 2,4,5-T-methyl ester


Part Number	Matrix	Size
N-10005-250MG	Neat	250MG

S-10005M1-1ML	100ug/ml in Methanol	1ML
S-10005M1-5ML	100ug/ml in Methanol	5ML

2,4-D methyl ester CAS#1928-38-7

Chemical Name: 2,4-Dichlorophenoxy acetic acid-methyl ester

PS: Solid

MF: C₉H₈Cl₂O₃
MW: 235.07



Part Number	Matrix	Size
N-10619-250MG	Neat	250MG

S-10619M1-1ML	100ug/ml in Methanol	1ML
S-10619M1-5ML	100ug/ml in Methanol	5ML

2-Ethylhexyl-2,4-dichlorophenoxy acetate CAS#1928-43-4
Chemical Name: 2-ethylhexyl(2,4-dichlorophenoxy)acetate

PS: Liquid


MF: C₁₆H₂₂Cl₂O₃
MW: 333.25

MP: -37°C

Alternate Name: Weed Rhop

Part Number	Matrix	Size
N-10356-1G	Neat	1G
S-10356A1-1ML	100ug/ml in Acetonitrile	1ML
S-10356J1-1ML	100ug/ml in Hexane	1ML
S-10531A4-1ML	1000ug/ml in Acetonitrile	1ML
S-10531A4-5ML	1000ug/ml in Acetonitrile	5ML


2,4,5-T 2-ethylhexyl ester CAS#1928-47-8
Chemical Name: 2-Ethylhexyl-2,4,5-trichlorophenoxyacetate

PS: Liquid

MF: C₁₆H₂₁Cl₃O₃


Part Number	Matrix	Size
N-10647-250MG	Neat	250MG
S-10647A1-1ML	100ug/ml in Acetonitrile	1ML
S-10647T1-1ML	100ug/ml in t-Butyl methyl ether	1ML


2,3-D methyl ester CAS#1928-56-9
Chemical Name: Methyl-2-(2,3-dichlorophenyl)acetate

MF: C₉H₈Cl₂O₃
MW: 235


Part Number	Matrix	Size
N-12861-100MG	Neat	100MG

2,4-D butoxyethyl ester CAS#1929-73-3
Chemical Name: 2-butoxyethyl(2,4-dichlorophenoxy)acetate

PS: Liquid


MF: C₁₄H₁₈Cl₂O₄
MW: 321.2

Alternate Name: Aqua-Kleen (TM)

Part Number	Matrix	Size
N-10610-250MG	Neat	250MG
S-10610A1-1ML	100ug/ml in Acetonitrile	1ML
S-10610A4-1ML	1000ug/ml in Acetonitrile	1ML
S-10610A4-5ML	1000ug/ml in Acetonitrile	5ML

Vernolate CAS#1929-77-7
Chemical Name: S-PROPYL DIPROPYLTHIOCARBAMATE

PS: Liquid

MF: C₁₀H₂₁NOS

MW: 203.35



Part Number	Matrix	Size
N-13744-250MG	Neat	250MG
S-13744A1-1ML	100ug/ml in Acetonitrile	1ML
S-13744T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-13744T1-5ML	100ug/ml in t-Butyl methyl ether	5ML

Nitrapyrin CAS#1929-82-4
Chemical Name: 2-Chloro-6-(trichloromethyl)pyridine

PS: Solid

MF: C₆H₃Cl₄N

MW: 230.9


MP: 62°C


Part Number	Matrix	Size
N-12659-100MG	Neat	100MG
S-12659A1-1ML	100ug/ml in Acetonitrile	1ML

Chlorbufam CAS#1967-16-4
Chemical Name: 1-Methyl-2-propynyl N-(3-chlorophenyl)carbamate

PS: Solid

MF: C₁₁H₁₀ClNO₂
MW: 223.67


MP: 45°C


Part Number	Matrix	Size
N-12970-10MG	Neat	10MG

Chloroxuron CAS#1982-47-4
Chemical Name: 3-[p-(p-Chlorophenoxy)phenyl]-1,1-dimethylurea

PS: Solid

MF: C₁₅H₁₅ClN₂O₂
MW: 290.75

MP: 151°C



Part Number	Matrix	Size
N-11456-250MG	Neat	250MG
S-11456A1-1ML	100ug/ml in Acetonitrile	1ML

Siduron CAS#1982-49-6
Chemical Name: 1-[2-Methylcyclohexyl]-3-phenylurea

PS: Solid

MF: C₁₄H₂₀N₂O

MW: 232.36

MP: 135°C



Part Number	Matrix	Size
N-13211-100MG	Neat	100MG
S-13211A4-1ML	1000ug/ml in Acetonitrile	1ML
S-13211A4-5ML	1000ug/ml in Acetonitrile	5ML

2,4-D dimethylamine salt CAS#2008-39-1
Chemical Name: Dimethylammonium(2,4-dichlorophenoxy)acetate

PS: Solid

MF: C₈H₆Cl₂O₃C₂H₇N

MW: 266.1

MP: 85°C

Alternate Name: DMA 4

Part Number	Matrix	Size
N-10612-250MG	Neat	250MG
S-10612M1-1ML	100ug/ml in Methanol	1ML

Butylate CAS#2008-41-5
Chemical Name: S-Ethyl N,N-diisobutylthiocarbamate

PS: Liquid

MF: C₁₁H₂₃NOS


MW: 217.41


Part Number	Matrix	Size
N-11383-250MG	Neat	250MG
S-11383A1-1ML	100ug/ml in Acetonitrile	1ML
S-11383T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-11383T1-5ML	100ug/ml in t-Butyl methyl ether	5ML

Aminocarb CAS#2032-59-9
Chemical Name: 4-(Dimethylamino)-m-tolyl methylcarbamate

PS: Solid

MF: C₁₁H₁₆N₂O₂
MW: 208.26

MP: 93°C



Part Number	Matrix	Size
N-11065-250MG	Neat	250MG
S-11065A4-1ML	1000ug/ml in Acetonitrile	1ML
S-11065A4-5ML	1000ug/ml in Acetonitrile	5ML
S-11065U1-1ML	100ug/ml in Toluene	1ML

Methiocarb CAS#2032-65-7
Chemical Name: 3,5-Dimethyl-4-[methylthio]phenylmethylcarbamate

PS: Solid

MF: C₁₁H₁₅NO₂S


MW: 225.3

MP: 119°C


Part Number	Matrix	Size
N-12398-100MG	Neat	100MG
S-12398M1-1ML	100ug/ml in Methanol	1ML

Oxytetracycline hydrochloride CAS#2058-46-0
PS: Solid

MF: C₂₂H₂₄N₂O₉ClH

MW: 496.94

Alternate Name: Aquacycline(TM)


Part Number	Matrix	Size
N-12743-250MG	Neat	250MG

EPN CAS#2104-64-5
Chemical Name: o-Ethyl-o-p-nitrophenyl thiophosphonate

PS: Solid

MF: C₁₄H₁₄NO₄PS

MW: 323.31

MP: 36°C



Part Number	Matrix	Size
N-11858-250MG	Neat	250MG
S-11858A1-1ML	100ug/ml in Acetonitrile	1ML
S-11858J1-1ML	100ug/ml in Hexane	1ML
S-11858J1-5ML	100ug/ml in Hexane	5ML
S-11858J4-1ML	1000ug/ml in Hexane	1ML
S-11858J4-5ML	1000ug/ml in Hexane	5ML

Bromophos methyl CAS#2104-96-3
Chemical Name: 4-Bromo-2,5-dichlorophenyl dimethyl phosphorothioate

PS: Solid

MF: C₈H₆BrCl₂O₃PS


MW: 365.98

MP: 53°C

Alternate Name: Bromophos

Part Number	Matrix	Size
N-11343-50MG	Neat	50MG
S-11343A1-1ML	100ug/ml in Acetonitrile	1ML
S-11343J1-1ML	100ug/ml in Hexane	1ML

Ethyl-1-naphthalene acetate
CAS#2122-70-5
Chemical Name: Ethyl 1-naphthylacetate

PS: Liquid


MF: C₁₄H₁₄O₂
MW: 214.18

Alternate Name: 1-Naphthaleneacetic acid, ethyl ester

Part Number	Matrix	Size
N-11914-250MG	Neat	250MG
S-11914A1-1ML	100ug/ml in Acetonitrile	1ML
S-11914T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Monosodium acid methane arsonate sesquihydrate
CAS#2163-80-6
Chemical Name: sodium hydrogen methylarsonate


PS: Solid

MF: CH₃AsNaO₃ X 1.5 H₂O

MW: 189

Alternate Name: MSMA sesquihydrate

Part Number	Matrix	Size
N-12495-100MG	Neat	100MG
S-12495F1-1ML	100ug/ml in Water	1ML

Chloranocryl
CAS#2164-09-200
Chemical Name: N-(3,4-Dichlorophenyl)methylacrylamide

PS: Solid

Alternate Name: Dicryl


Part Number	Matrix	Size
N-11421-250MG	Neat	250MG

Fluometuron
CAS#2164-17-2
Chemical Name: 1,1-Dimethyl-3-[a,a,a-trifluoro-m-tolyl]urea

PS: Solid

MF: C₁₀H₁₁F₃N₂O

MW: 232.21


MP: 163°C


Part Number	Matrix	Size
N-11992-250MG	Neat	250MG
S-11992A4-1ML	1000ug/ml in Acetonitrile	1ML

Molinate
CAS#2212-67-1
Chemical Name: S-Ethylhexahydro-1H-azepine-1-carbothioate

PS: Liquid

MF: C₉H₁₇NOS


MW: 187.31


Part Number	Matrix	Size
N-12487-250MG	Neat	250MG
S-12487A1-1ML	100ug/ml in Acetonitrile	1ML
S-12487T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-12487T1-5ML	100ug/ml in t-Butyl methyl ether	5ML

Dienochlor
CAS#2227-17-0
Chemical Name: Decachlorobis[2,4-cyclopentadiene-1-yl]


PS: Solid

MF: C₁₀Cl₁₀
MW: 474.6

MP: 122°C


Part Number	Matrix	Size
N-12824-250MG	Neat	250MG
S-12824M1-1ML	100ug/ml in Methanol	1ML

Dimethylvinphos Solution
CAS#2274-67-1
Chemical Name: 2-Chloro-1-(2,4-dichlorophenyl)vinyl dimethyl phosphate

PS: Liquid


Part Number	Matrix	Size
S-14052U1-1ML	100ug/ml in Toluene	1ML

Vamidothion
CAS#2275-23-2
Chemical Name: O,O-Dimethyl S-2-(1-N-methylcarbamoylthio)ethyl thiophosphate

PS: Solid

MF: C₈H₁₈NO₄PS₂
MW: 287.35

MP: 43°C

Alternate Name: Kilval(TM)


Part Number	Matrix	Size
N-13741-100MG	Neat	100MG
S-13741A1-1ML	100ug/ml in Acetonitrile	1ML

Di-allate CAS#2303-16-4
Chemical Name: S-[2,3-Dichloroallyl]-N,N'-diisopropylthiocarbamate

PS: Liquid

MF: C₁₀H₁₇Cl₂NOS

MW: 270.24

MP: 25°C



Part Number	Matrix	Size
N-11587-100MG	Neat	100MG
S-11587M1-1ML	100ug/ml in Methanol	1ML
S-11587M1-5ML	100ug/ml in Methanol	5ML

Tri-allate CAS#2303-17-5
Chemical Name: S-2,3,3-Trichloroallyl diisopropylthiocarbamate

PS: Solid

MF: C₁₀H₁₆Cl₃NOS

MW: 304.7

MP: 29°C

Alternate Name: Eccothal


Part Number	Matrix	Size
N-13628-250MG	Neat	250MG
S-13628A1-1ML	100ug/ml in Acetonitrile	1ML
S-13628K1-1ML	100ug/ml in Isooctane	1ML

Pentanochlor CAS#2307-68-8
Chemical Name: N-(3-Chloro-4-methylphenyl)-2-methylpentanamide

PS: Solid

MF: C₁₃H₁₈ClNO

MW: 239.7


MP: 85°C

Alternate Name: Solan(TM)

Part Number	Matrix	Size
N-12844-100MG	Neat	100MG
S-12844A1-1ML	100ug/ml in Acetonitrile	1ML
S-12844U1-1ML	100ug/ml in Toluene	1ML

Phosalone CAS#2310-17-0
Chemical Name: O,O-Diethyl-S-[-(6-chloro-2-oxobenzoxazolin-3-yl)-methyl]phosphorodithioate

PS: Solid

MF: C₁₂H₁₅ClNO₄PS₂
MW: 367.8


MP: 47.5°C


Part Number	Matrix	Size
N-13036-100MG	Neat	100MG
S-13036M1-1ML	100ug/ml in Methanol	1ML

Propargite CAS#2312-35-8
Chemical Name: 2-[4-(1,1-Dimethylethyl)phenoxy]cyclohexyl-2-propynyl sulfite

PS: Liquid

MF: C₁₉H₂₆O₄S

MW: 350.51

Alternate Name: Omite(TM)

Part Number	Matrix	Size
N-12727-100MG	Neat	100MG
S-12727A1-1ML	100ug/ml in Acetonitrile	1ML
S-12727U1-1ML	100ug/ml in Toluene	1ML

Mirex CAS#2385-85-5
Chemical Name: 1,1a,2,2,3,3a,4,5,5,5a,5b,6-dodecachlorooctahydro-1,3,4-metheno-1H-cyclobuta[cd]pentalene


PS: Solid

MF: C₁₀Cl₁₂
MW: 545.59

MP: 485°C


Part Number	Matrix	Size
N-12486-100MG	Neat	100MG
S-12486M1-1ML	100ug/ml in Methanol	1ML
S-FC2121S-1ML	200ug/ml in Toluene (C13)	1ML

4-Chloro-o-tolyloxyacetic acid methyl ester CAS#2436-73-9
PS: Liquid


MF: C₁₀H₁₁ClO₃
MW: 214.65

Alternate Name: MCPA methyl ester

Part Number	Matrix	Size
N-10819-100MG	Neat	100MG
S-10819M1-1ML	100ug/ml in Methanol	1ML
S-10819M1-5ML	100ug/ml in Methanol	5ML

Glyphosine CAS#2439-99-8
Chemical Name: N-Carboxymethyl-N,N-bis(methylenephosphonic acid)-amine

PS: Solid


MF: C₄H₁₁NO₈P₂
MW: 263.08

MP: 200°C

Alternate Name: Polaris

Part Number	Matrix	Size
N-12135-100MG	Neat	100MG
S-12135F1-1ML	100ug/ml in Water	1ML
S-12135U1-1ML	100ug/ml in Toluene	1ML

Dicaphon CAS#2463-84-5
Chemical Name: O,O-Dimethyl-O-[2-chloro-4-nitrophenyl]phosphorothioate

PS: Solid
MF: C₈H₉ClNO₃PS
MW: 297.66
MP: 53°C



Part Number	Matrix	Size
N-11658-100MG	Neat	100MG
S-11658A1-1ML	100ug/ml in Acetonitrile	1ML
S-11658U1-1ML	100ug/ml in Toluene	1ML

EGT CAS#2514-53-6
Chemical Name: Ethylene bis(trichloroacetate)


PS: Solid
MF: C₆H₄Cl₆O₄
MW: 352.8
MP: 40.3°C



Part Number	Matrix	Size
N-11845-100MG	Neat	100MG

O,O-Dimethyl phosphorochloridothioate CAS#2524-03-0

PS: Liquid
MF: C₂H₆ClO₂PS
MW: 160.56




FP: 105°C
Alternate Name: Dimethyl chlorothiophosphate

Part Number	Matrix	Size
N-12704-250MG	Neat	250MG
S-12704A1-1ML	100ug/ml in Acetonitrile	1ML
S-12704U1-1ML	100ug/ml in Toluene	1ML

Chlorflurecol-methyl ester CAS#2536-31-4
Chemical Name: Methyl-2-chloro-9-hydroxyfluorene-9-carboxylate


PS: Solid
MF: C₁₅H₁₁ClO₃
MW: 274.71
MP: 136°C



Part Number	Matrix	Size
N-11431-100MG	Neat	100MG
S-11431M1-1ML	100ug/ml in Methanol	1ML

Formothion CAS#2540-82-1
Chemical Name: 2-DIMETHOXYPHOSPHINOTHIOYLTHIO-N-FORMYL-N-METHYLAC-ETAMIDE


PS: Liquid
MF: C₆H₁₂NO₄PS₂
MW: 257.27
MP: 25°C


Alternate Name: ANTHIO™

Part Number	Matrix	Size
N-12017-100MG	Neat	100MG
S-12017A1-1ML	100ug/ml in Acetonitrile	1ML
S-12017J1-1ML	100ug/ml in Hexane	1ML

2,4,5-T butoxyethyl ester CAS#2545-59-7
Chemical Name: Butoxyethyl-2,4,5-trichlorophenoxyacetate


PS: Liquid
MF: C₁₄H₁₇Cl₃O₄
MW: 355.66



Part Number	Matrix	Size
N-10648-100MG	Neat	100MG
S-10648A4-1ML	1000ug/ml in Acetonitrile	1ML
S-10648A4-5ML	1000ug/ml in Acetonitrile	5ML
S-10648U1-1ML	100ug/ml in Toluene	1ML

N,N-Dimethylsuccinamic acid CAS#2564-95-6


PS: Solid
MF: C₆H₁₁NO₃
MW: 145.16
MP: 85°C



Part Number	Matrix	Size
N-12630-250MG	Neat	250MG
S-12630A1-1ML	100ug/ml in Acetonitrile	1ML
S-12630T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Etridiazole CAS#2593-15-9
Chemical Name: 5-Ethoxy-3-trichloromethyl-1,2,4-thiadiazole


PS: Liquid
MF: C₅H₅Cl₃N₂OS
MW: 247.5
MP: 19.9°C



Part Number	Matrix	Size
N-13517-250MG	Neat	250MG
S-13517A1-1ML	100ug/ml in Acetonitrile	1ML
S-13517T1-1ML	100ug/mL in t-Butyl methyl ether	1ML

Mecarbam CAS#2595-54-2
Chemical Name: S-(N-Ethoxycarbonyl-N-methylcarbamoylmethyl) O,O-diethyl phosphorodithioate

PS: Liquid

MF: C₁₀H₂₀NO₅PS₂
MW: 329.4

FP: 100°C


Alternate Name: Afos(TM)

Part Number	Matrix	Size
N-12359-50MG	Neat	50MG
S-12359A1-1ML	100ug/ml in Acetonitrile	1ML
S-12359U1-1ML	100ug/ml in Toluene	1ML

Emmi (Technical) CAS#2597-93-5
Chemical Name: Ethylmercurichlorethimide

PS: Solid

MF: C₉H₄Cl₆N₃O₂Hg


MW: 599.37


Part Number	Matrix	Size
N-10224-100MG	Neat	100MG

Promecarb CAS#2631-37-0
Chemical Name: 3-Isopropyl-5-methylphenyl methyl carbamate

PS: Solid


MF: C₁₂H₁₇NO₂
MW: 207.28

MP: 87°C

Alternate Name: Carbamult (TM)

Part Number	Matrix	Size
N-13101-100MG	Neat	100MG
S-13101M4-1ML	1000ug/ml in Methanol	1ML

Isoprocarb CAS#2631-40-5
Chemical Name: 2-Isopropylphenyl methyl carbamate

PS: Solid

MF: C₁₁H₁₅NO₂

Alternate Name: Mipcin*


Part Number	Matrix	Size
N-12264-250MG	Neat	250MG
S-12264A1-1ML	100ug/mL in Acetonitrile	1ML
S-12264U1-1ML	100ug/ml in Toluene	1ML

Cyanophos CAS#2636-26-2
Chemical Name: O-4-Cyanophenyl-O,O-dimethyl phosphorothioate

PS: Liquid

MF: C₉H₁₀NO₃PS

MW: 243.21


MP: 14°C

Alternate Name: O-4-Cyanophenyl

Part Number	Matrix	Size
N-11517-100MG	Neat	100MG
S-11517A1-1ML	100ug/ml in Acetonitrile	1ML
S-11517U1-1ML	100ug/ml in Toluene	1ML

Azinphos-ethyl CAS#2642-71-9
Chemical Name: O,O-diethyl S-[(4-oxo-1,2,3-benzotriazin-3(4H)-yl)methyl] phosphorodithioate

PS: Solid

MF: C₁₂H₁₆N₃O₃PS₂
MW: 345.4


MP: 53°C

Alternate Name: GUTHION ETHYL

Part Number	Matrix	Size
N-12140-100MG	Neat	100MG
S-12140A1-1ML	100ug/mL in Acetonitrile	1ML

XMC CAS#2655-14-3
Chemical Name: 3,5-dimethylphenyl N-methylcarbamate

PS: Solid


MF: C₁₀H₁₃NO₂
MW: 179.2

MP: 99°C

Alternate Name: 3,5-XYLYL METHYLCARBAMATE

Part Number	Matrix	Size
N-13275-100MG	Neat	100MG
S-13275A1-1ML	100ug/mL in Acetonitrile	1ML

2,3,5-Trimethylphenyl methyl carbamate CAS#2655-15-4
PS: Solid

MF: C₁₁H₁₅NO₂
MW: 193.2

MP: 124°C

Alternate Name: 2,3,5-Trimethacarb


Part Number	Matrix	Size
N-10603-50MG	Neat	50MG
S-10603A1-1ML	100ug/ml in Acetonitrile	1ML

Chloroneb
CAS#2675-77-6
Chemical Name: 1,4-Dichloro-2,5-dimethoxybenzene

PS: Solid


MF: C₈H₈Cl₂O₂
MW: 207.06

MP: 133°C

FP: 161°C



Part Number	Matrix	Size
N-11449-250MG	Neat	250MG
S-11449A1-1ML	100ug/ml in Acetonitrile	1ML
S-11449T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-11449T1-5ML	100ug/ml in t-Butyl methyl ether	5ML

3,4,5-Trimethylphenyl methyl carbamate
CAS#2686-99-9
PS: Solid

MF: C₁₁H₁₅NO₂
MW: 193.2

Alternate Name: 3,4,5-Trimethacarb

Part Number	Matrix	Size
N-10769-10MG	Neat	10MG

2,4-DB dimethylamine salt
CAS#2758-42-1
Chemical Name: Dimethylamine-4[2,4-dichlorophenoxy]butyrate


PS: Solid


Part Number	Matrix	Size
N-10623-250MG	Neat	250MG
S-10623M1-1ML	100ug/ml in Methanol	1ML

Quinoclamine
CAS#2797-51-5
Chemical Name: 2-Amino-3-chloro-1,4-naphthalenedione

PS: Solid


MF: C₁₀H₆ClNO₂
MW: 207.62

MP: 202°C

Alternate Name: 2-Amino-3-chloro-1,4-naphthoquinone

Part Number	Matrix	Size
N-13825-100MG	Neat	100MG

Methyl-1-naphthalene acetate
CAS#2876-78-0
Chemical Name: 1-Naphthaleneacetic acid-Methyl ester

PS: Liquid


MF: C₁₃H₁₂O₂
MW: 200.25


Part Number	Matrix	Size
N-12464-500MG	Neat	500MG
S-12464M1-1ML	100ug/ml in Methanol	1ML
S-12464U1-1ML	100ug/ml in Toluene	1ML

3,5-Dichlorobenzoic acid methyl ester
CAS#2905-67-1
Chemical Name: Methyl-3,5-dichlorobenzoate

PS: Solid

MF: C₈H₆Cl₂O₂
MW: 205.04

MP: 58°C



Part Number	Matrix	Size
N-10791-100MG	Neat	100MG
S-10791A1-1ML	100ug/ml in Acetonitrile	1ML
S-10791T5-1ML	2000ug/ml in t-Butyl methyl ether	1ML
S-10791T5-5ML	2000ug/ml in t-Butyl methyl ether	5ML

Chlorpyrifos
CAS#2921-88-2
Chemical Name: O,O-Diethyl-O-[3,5,6-trichloro-2-pyridyl]phosphorothioate

PS: Solid


MF: C₉H₁₁Cl₃NO₃PS

MW: 350.59

MP: 41°C


Part Number	Matrix	Size
N-11459-250MG	Neat	250MG
S-11459A1-1ML	100ug/ml in Acetonitrile	1ML
S-11459K1-1ML	100ug/ml in Isooctane	1ML
S-11459K1-5ML	100ug/ml in Isooctane	5ML

2,3-Dichlorophenoxyacetic acid
CAS#2976-74-1
PS: Solid


MF: C₈H₆Cl₂O₃
MW: 221


Part Number	Matrix	Size
N-12953-10MG	Neat	10MG

Allyl trimethylammonium bromide
CAS#3004-51-1
PS: Solid

MF: C₆H₁₄NBr

MW: 179.9


MP: 170°C


Part Number	Matrix	Size
N-11059-250MG	Neat	250MG
S-11059F1-1ML	100ug/ml in Water	1ML
S-11059U1-1ML	100ug/ml in Toluene	1ML

Metobromuron
CAS#3060-89-7
Chemical Name: 3-[p-Bromophenyl]-1-methyl-1-methoxyurea

PS: Solid

MF: C₉H₁₁BrN₂O₂
MW: 259.1

MP: 95°C


Part Number	Matrix	Size
N-12477-250MG	Neat	250MG
S-12477A1-1ML	100ug/ml in Acetonitrile	1ML
S-12477U1-1ML	100ug/ml in Toluene	1ML

Hexachlorodimethyl sulfone
CAS#3064-70-8
PS: Solid

MF: C₂Cl₆O₂S

MW: 300.78


MP: 36°C


Part Number	Matrix	Size
N-12161-250MG	Neat	250MG
S-12161A1-1ML	100ug/ml in Acetonitrile	1ML
S-12161U1-1ML	100ug/ml in Toluene	1ML

2,3-Dichloronitrobenzene
CAS#3209-22-1
Chemical Name: 1,2-Dichloro-3-nitrobenzene

PS: Solid


MF: C₆H₃Cl₂NO₂
MW: 192

MP: 61.5°C


Part Number	Matrix	Size
N-10584-250MG	Neat	250MG
S-10584A1-1ML	100ug/ml in Acetonitrile	1ML
S-10584U1-1ML	100ug/ml in Toluene	1ML

Aspon (TM)
CAS#3244-90-4
Chemical Name: O,O,O',O'-Tetrapropylthiopyrophosphate

PS: Liquid

MF: C₁₂H₂₈O₅P₂S₂
MW: 378.46



Part Number	Matrix	Size
N-11103-10MG	Neat	10MG
S-11103A1-1ML	100ug/ml in Acetonitrile	1ML
S-11103J4-1ML	1000ug/ml in Hexane	1ML
S-11103J4-5ML	1000ug/ml in Hexane	5ML

Asulam
CAS#3337-71-1
Chemical Name: Methyl sulfanilylcarbamate

PS: Solid

MF: C₈H₁₀N₂O₄S

MW: 230.26

MP: 143°C



Part Number	Matrix	Size
N-11104-250MG	Neat	250MG
S-11104M1-1ML	100ug/ml in Methanol	1ML
S-11104M1-5ML	100ug/ml in Methanol	5ML

Dithianon
CAS#3347-22-6
Chemical Name: 5,10-Dihydro-5,10-dioxonaphtho[2,3-b]1,4-dithi-in-2,3-dicarbonitrile

PS: Solid

MF: C₁₄H₄N₂O₂S₂
MW: 296.33

MP: 220°C


FP: 300°C


Part Number	Matrix	Size
N-11821-250MG	Neat	250MG

Triclosan
CAS#3380-34-5
Chemical Name: 2,4,4'-TRICHLORO-2'-HYDROXYDIPHENYL ETHER

PS: Solid

MF: C₁₂H₇Cl₃O₂
MW: 289.55

MP: 54°C



Part Number	Matrix	Size
N-13076-100MG	Neat	100MG
S-13076M1-1ML	100ug/ml in Methanol	1ML

Temephos CAS#3383-96-8

Chemical Name: O,O,O',O'-Tetramethyl-O,O'-thio-di-p-phenylene phosphorothioate

PS: Solid

MF: C₁₆H₂₀O₆P₂S₃
MW: 466.47


MP: 30°C

Alternate Name: Abate

Part Number	Matrix	Size
N-10996-100MG	Neat	100MG
S-10996A1-1ML	100ug/ml in Acetonitrile	1ML
S-10996U1-1ML	100ug/ml in Toluene	1ML

o,p'-DDE CAS#3424-82-6

Chemical Name: 1,1-Dichloro-2-[p-chlorophenyl]-2-[o-chlorophenyl]ethylene

PS: Solid

MF: C₁₄H₈Cl₄
MW: 318.02



Part Number	Matrix	Size
N-12707-50MG	Neat	50MG
S-12707M1-1ML	100ug/ml in Methanol	1ML
S-12707M1-5ML	100ug/ml in Methanol	5ML
S-FC2306S-1.2ML	100ug/ml in n-Nonane (C13)	1.2ML

Piperalin CAS#3478-94-2

Chemical Name: 3-(2-Methylpiperidino)propyl-3,4-dichlorobenzoate

PS: Liquid

MF: C₁₆H₂₁Cl₂NO₂
MW: 330.3

MP: 208°C

Alternate Name: Pipron


Part Number	Matrix	Size
N-13058-100MG	Neat	100MG
S-13058A1-1ML	100ug/ml in Acetonitrile	1ML
S-13058U1-1ML	100ug/ml in Toluene	1ML

2,3,5,6-Tetrachloroaniline CAS#3481-20-7

PS: Solid

MF: C₆H₃Cl₄N

MW: 230.91

MP: 106°C


Part Number	Matrix	Size
N-10530-250MG	Neat	250MG
S-10530M1-1ML	100ug/ml in Methanol	1ML

2-Hydroxyethyl-n-octyl sulphide CAS#3547-33-9

Chemical Name: 2-(octylthio)ethanol

PS: Liquid

MF: C₁₀H₂₂OS

MW: 190.38



Part Number	Matrix	Size
N-10378-100MG	Neat	100MG
S-10378M1-1ML	100ug/ml in Methanol	1ML

Asomate CAS#3586-60-5

Chemical Name: arsinetriyl tris(dimethyldithiocarbamate)

PS: Solid

MF: C₉H₁₈AsN₃S₆
MW: 435.58


MP: 224°C

Alternate Name: Azomat

Part Number	Matrix	Size
N-13904-50MG	Neat	50MG

2,4-D ethanolamine salt CAS#3599-58-4

Chemical Name: 2-hydroxyethylazanium

PS: Solid


MF: C₈H₆Cl₂O₃·C₂H₇NO


Part Number	Matrix	Size
N-10613-250MG	Neat	250MG
S-10613M1-1ML	100ug/ml in Methanol	1ML

Tetraethyldithiopyrophosphate CAS#3689-24-5

Chemical Name: O,O,O,O-Tetraethyl dithiopyrophosphate

PS: Liquid


MF: C₈H₂₀O₅P₂S₂
MW: 322.32

Alternate Name: Thiotep

Part Number	Matrix	Size
N-13540-50MG	Neat	50MG
S-13540A1-1ML	100ug/ml in Acetonitrile	1ML
S-13540J1-1ML	100ug/ml in Hexane	1ML

Chlorophacinone
CAS#3691-35-8
Chemical Name: 2-[2-(4-Chlorophenyl)-2-phenylacetyl]indan-1,3-dione

PS: Solid


MF: C₂₃H₁₅ClO₃
MW: 374.82

MP: 138°C


Part Number	Matrix	Size
N-11451-250MG	Neat	250MG
S-11451A1-1ML	100ug/ml in Acetonitrile	1ML

Bitrex (TM)
CAS#3734-33-6
Chemical Name: Benzyldiethyl[(2,6-xylylcarbamoyl)methyl]ammonium benzoate

PS: Solid


MF: C₂₈H₃₄N₂O₃
MW: 446.58

Alternate Name: Denatonium benzoate

Part Number	Matrix	Size
N-11324-100MG	Neat	100MG
S-11324F1-1ML	100ug/mL in Water	1ML
S-11324U1-1ML	100ug/ml in Toluene	1ML

Fenchlorim
CAS#3740-92-9
Chemical Name: 4,6-DICHLORO-2-PHENYLPYRIMIDINE

PS: Solid


MF: C₁₀H₆Cl₂N₂
MW: 225.1

MP: 96.9°C


Part Number	Matrix	Size
N-13265-250MG	Neat	250MG

Buturon
CAS#3766-60-7
Chemical Name: N?-(4-chlorophenyl)-N-methyl-N-(1-methyl-2-propyn-1-yl)urea


PS: Liquid

MF: C₁₂H₁₃ClN₂O


Part Number	Matrix	Size
N-12944-10MG	Neat	10MG

Baycarb (TM)
CAS#3766-81-2
Chemical Name: 2-sec-Butylphenyl-N-methylcarbamate


PS: Solid

MF: C₁₂H₁₇NO₂

Alternate Name: BPMC


Part Number	Matrix	Size
N-11126-250MG	Neat	250MG
S-11126M1-1ML	100ug/ml in Methanol	1ML

Streptomycin sulfate
CAS#3810-74-0
PS: Solid

MF: C₂₁H₃₉N₇O_{12.3/2}H₂O₄S

MW: 1457.58

Alternate Name: Streptomycin sulfate (2:3) salt

Part Number	Matrix	Size
N-13228-250MG	Neat	250MG


Dioxabenzofos solution
CAS#3811-49-2
PS: Liquid

Alternate Name: Salithion

Part Number	Matrix	Size
S-13964A1-1ML	100ug/ml in Acetonitrile	1ML

Ioxynil octanoate
CAS#3861-47-0
Chemical Name: 4-Cyano-2,6-di-iodophenyl octanoate

PS: Solid


MF: C₁₅H₁₇I₂NO₂
MW: 497.11

MP: 52°C


Part Number	Matrix	Size
N-12218-100MG	Neat	100MG
S-12218A1-1ML	100ug/ml in Acetonitrile	1ML
S-12218U1-1ML	100ug/ml in Toluene	1ML

Fuberidazole
CAS#3878-19-1
Chemical Name: 2-(2'-Furyl)benzimidazole


PS: Solid
MF: C₁₁H₈N₂O
MW: 184.21
MP: 292°C



Part Number	Matrix	Size
N-12021-100MG	Neat	100MG
S-12021M1-1ML	100ug/ml in Methanol	1ML

Dipropetryn
CAS#4147-51-7
Chemical Name: 2-Ethylthio-4,6-bis(isopropylamino)-1,3,5-triazine


PS: Solid
MF: C₁₁H₂₁N₅S
MW: 255.43
MP: 105°C



Part Number	Matrix	Size
N-11808-100MG	Neat	100MG
S-11808A1-1ML	100ug/ml in Acetonitrile	1ML
S-11808T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Fenuron TCA
CAS#4482-55-7
Chemical Name: 1,1-Dimethyl-3-phenylurea trichloroacetate


PS: Solid
MF: C₉H₁₂N₂O₂Cl₃O₂
MW: 327.61
MP: 65°C


Alternate Name: 1,1-Dimethyl-3-phenylurea trichloroacetate

Part Number	Matrix	Size
N-11968-50MG	Neat	50MG

Nitralin
CAS#4726-14-1
Chemical Name: [Methylsulfonyl]-2,6-dinitro-N,N-dipropylaniline


PS: Solid
MF: C₁₃H₁₉N₃O₆S
MW: 345.4
MP: 151°C



Part Number	Matrix	Size
N-12658-250MG	Neat	250MG
S-12658A1-1ML	100ug/ml in Acetonitrile	1ML
S-12658J1-1ML	100ug/ml in Hexane	1ML

Phthalanillic acid
CAS#4727-29-1
Chemical Name: 2-((Phenylamino)carbonyl)benzoic acid


PS: Solid
MW: 241.26
MP: 169°C


Alternate Name: Nevirol(TM)

Part Number	Matrix	Size
N-13043-250MG	Neat	250MG
S-13043A1-1ML	100ug/ml in Acetonitrile	1ML
S-13043T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Bromophos ethyl
CAS#4824-78-6
Chemical Name: O,O-Diethyl O-2,5-dichloro-4-bromophenylphosphorothioate


PS: Liquid
MF: C₁₀H₁₂BrCl₂O₃PS
MW: 394
FP: 100°C



Part Number	Matrix	Size
N-11342-100MG	Neat	100MG
S-11342J1-1ML	100ug/ml in Hexane	1ML

Silvex methyl ester
CAS#4841-20-7
Chemical Name: 2-[2,4,5-Trichlorophenoxy]propionic acid, methyl ester


PS: Solid
MF: C₈H₅N₅S
MW: 213.34
MP: 83°C



Part Number	Matrix	Size
N-13214-250MG	Neat	250MG
S-13214M1-1ML	100ug/mL in Methanol	1ML
S-13214M1-5ML	100ug/mL in Methanol	5ML

Karbutilate
CAS#4849-32-5
Chemical Name: 3-[[[(dimethylamino)carbonyl]amino]phenyl N-(1,1-dimethylethyl) carbamate


PS: Solid
MF: C₁₄H₂₁N₃O₃



Part Number	Matrix	Size
N-12945-250MG	Neat	250MG

cis-Chlordane
CAS#5103-71-9
PS: Solid


MF: C₁₀H₆Cl₈
MW: 409.8

MP: 104°C

Alternate Name: a-Chlordan

Part Number	Matrix	Size
N-11480-10MG	Neat	10MG
S-11480M1-1ML	100ug/ml in Methanol	1ML
S-11480M1-5ML	100ug/ml in Methanol	5ML

trans-Chlordane
CAS#5103-74-2
PS: Solid

MF: C₁₀H₆Cl₈
MW: 409.8

MP: 104°C

Alternate Name: beta-Chlordan


Part Number	Matrix	Size
N-13615-10MG	Neat	10MG
S-13615M1-1ML	100ug/ml in Methanol	1ML
S-13615M1-5ML	100ug/ml in Methanol	5ML
S-FC908S-1.2ML	100ug/ml in n-Nonane (C13)	1.2ML

Ditalimfos
CAS#5131-24-8
Chemical Name: o,o-Diethyl phthalimidophosphonothioate

PS: Solid

MF: C₁₂H₁₄NO₄PS

MW: 299.3

MP: 83°C

Alternate Name: Plondrel(TM)

Part Number	Matrix	Size
N-11820-100MG	Neat	100MG
S-11820A1-1ML	100ug/ml in Acetonitrile	1ML
S-11820J1-1ML	100ug/ml in Hexane	1ML
S-11820J1-5ML	100ug/ml in Hexane	5ML
S-11820U1-1ML	100ug/ml in Toluene	1ML


Carboxin
CAS#5234-68-4
Chemical Name: 5,6-Dihydro-2-methyl-1,4-oxathiin-3-carboxanilide

PS: Solid

MF: C₁₂H₁₃NO₂S

MW: 235.31

MP: 94°C

FP: 100°C



Part Number	Matrix	Size
N-11410-250MG	Neat	250MG
S-11410A1-1ML	100ug/ml in Acetonitrile	1ML

Oxycarboxin
CAS#5259-88-1
Chemical Name: 2,3-dihydro-6-methyl-5-phenylcarbamoyl-1,4-oxathiine 4,4-dioxide

PS: Solid

MF: C₁₂H₁₃NO₄S

MW: 267.32

MP: 119.5°C


Part Number	Matrix	Size
N-12740-100MG	Neat	100MG
S-12740A1-1ML	100ug/ml in Acetonitrile	1ML
S-12740U1-1ML	100ug/ml in Toluene	1ML


Chlorpyrifos Methyl
CAS#5598-13-0
Chemical Name: O,O-Dimethyl-O-3,5,6-trichloro-2-pyridyl phosphorothioate

PS: Solid

MF: C₇H₇Cl₃NO₃PS

MW: 327.5

MP: 45.5°C


FP: 182°C


Part Number	Matrix	Size
N-11460-250MG	Neat	250MG
S-11460A1-1ML	100ug/ml in Acetonitrile	1ML
S-11460J4-1ML	1000ug/ml in Hexane	1ML

2,4-D isopropylamine salt
CAS#5742-17-6
Chemical Name: Isopropylammonium (2,4-dichlorophenoxy) acetate


PS: Solid

MF: C₁₁H₁₅Cl₂NO₃
MW: 280.15

MP: 121°C


Part Number	Matrix	Size
N-10618-250MG	Neat	250MG
S-10618M1-1ML	100ug/ml in Methanol	1ML

2-[3-Chlorophenoxy]propionamide
CAS#5825-87-6
PS: Solid

MF: C₉H₉ClNO₂
MW: 199.64


Part Number	Matrix	Size
N-10261-1G	Neat	1G
S-10261M1-1ML	100ug/ml in Methanol	1ML

Coumatetralyl
CAS#5836-29-3
Chemical Name: 4-hydroxy-3-(1,2,3,4-tetrahydro-1-naphthalenyl)-2H-1-benzopyran-2-one

PS: Solid


MF: C₁₉H₁₆O₃
MW: 292.3

MP: 172°C

FP: 100°C

Alternate Name: Racumin(TM)

Part Number	Matrix	Size
N-11509-250MG	Neat	250MG
S-11509M1-1ML	100ug/ml in Methanol	1ML



Terbacil
CAS#5902-51-2
Chemical Name: 5-Chloro-3-(1,1-dimethylethyl)-6-methyl-2,4(1H,3H)-pyrimidinedione

PS: Solid

MF: C₉H₁₃ClN₂O₂
MW: 216.67

MP: 175°C

Part Number	Matrix	Size
N-13509-250MG	Neat	250MG
S-13509A1-1ML	100ug/ml in Acetonitrile	1ML
S-13509T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-13509T1-5ML	100ug/ml in t-Butyl methyl ether	5ML


Vanicide-20S (TM)
CAS#5902-85-2
Chemical Name: a,a,a-Trifluoro-2,6-dinitro-N,N-dipropyl-p-toluidine


PS: Solid

MF: C₇H₅NS₂C₂H₇NO

MW: 228.35

Alternate Name: 2-MERCAPTOBENZOTHAZOLE MONOETHANOLAMINE SALT

Part Number	Matrix	Size
N-13742-250MG	Neat	250MG
S-13742A1-1ML	100ug/ml in Acetonitrile	1ML
S-13742U1-1ML	100ug/ml in Toluene	1ML



N-Methyl-N-1-naphthyl acetamide
CAS#5903-13-9
PS: Solid

MF: C₁₃H₁₃NO

MW: 199.25

MP: 96°C

Part Number	Matrix	Size
N-12559-100MG	Neat	100MG
S-12559M1-1ML	100ug/ml in Methanol	1ML


Terbuthylazine
CAS#5915-41-3
Chemical Name: 6-Chloro-N2-(1,1-dimethylethyl)-N4-ethyl-1,3,5-triazine-2,4-diamine


PS: Solid

MF: C₉H₁₆ClN₅
MW: 229.7


MP: 177°C

FP: 100°C

Part Number	Matrix	Size
N-13512-50MG	Neat	50MG
S-13512A1-1ML	100ug/ml in Acetonitrile	1ML
S-13512T4-1ML	1000ug/ml in t-Butyl methyl ether	1ML


Nimbin
CAS#5945-86-8
PS: Liquid

Part Number	Matrix	Size
N-12912-5MG	Neat	5MG



Dinoseb methyl ether
CAS#6099-79-2
Chemical Name: 2-sec-Butyl-4,6-dinitro phenol methyl ether

PS: Liquid

MF: C₁₁H₁₄N₂O₅
MW: 254.24

Alternate Name: DNBP methyl ether

Part Number	Matrix	Size
N-11787-100MG	Neat	100MG
S-11787M1-1ML	100ug/ml in Methanol	1ML
S-11787M1-5ML	100ug/ml in Methanol	5ML



Chlordimeform
CAS#6164-98-3
Chemical Name: N'-(4-chloro-2-methylphenyl)-N,N-dimethylmethanimidamide

PS: Solid

MF: C₁₀H₁₃ClN₂
MW: 196.67


MP: 35°C

Part Number	Matrix	Size
N-11426-250MG	Neat	250MG
S-11426A1-1ML	100ug/ml in Acetonitrile	1ML
S-11426U1-1ML	100ug/ml in Toluene	1ML



2,4-D isopropanol amine salt
CAS#6365-72-6
Chemical Name: Isopropanolamine 2,4-dichlorophenoxyacetate


PS: Solid

MF: C₈H₈Cl₂O₃C₃H₉N


Part Number	Matrix	Size
N-10616-250MG	Neat	250MG
S-10616M1-1ML	100ug/ml in Methanol	1ML


Diquat Dibromide monohydrate
CAS#6385-62-2
PS: Solid

MF: C₁₂H₁₂N₂Br₂
MW: 344.05

MP: 337°C

Alternate Name: Weedtrin-D (tm)

Part Number	Matrix	Size
N-11816-500MG	Neat	500MG
S-11816F1-1ML	100ug/ml in Water	1ML
S-11816F13-1ML	1970ug/ml in Water	1ML
S-11816F13-5ML	1970ug/ml in Water	5ML


2-sec-Butyl-4,6-dinitrophenol triethanolamine salt
CAS#6420-47-9
PS: Solid

MF: C₁₀H₁₂N₂O₅C₆H₁₅NO₃
MW: 389.46

Alternate Name: Elgetol

Part Number	Matrix	Size
N-10517-500MG	Neat	500MG
S-10517F1-1ML	100ug/ml in Water	1ML
S-10517U1-1ML	100ug/ml in Toluene	1ML

Dicamba methyl ester
CAS#6597-78-0
Chemical Name: 2-Methoxy-3,6-dichloromethylbenzoate

PS: Solid

MF: C₉H₈Cl₂O₃
MW: 235.07

Alternate Name: Disugran

Part Number	Matrix	Size
N-11657-100MG	Neat	100MG
S-11657M1-1ML	100ug/ml in Methanol	1ML
S-11657M1-5ML	100ug/ml in Methanol	5ML

Metam sodium
CAS#6734-80-1
PS: Solid

MF: C₂H₄NS₂Na₂H₂O

MW: 165.22

Alternate Name: Sodiummethyl-dithiocarbamate-dihydrate


Part Number	Matrix	Size
N-12382-250MG	Neat	250MG

Monocrotophos
CAS#6923-22-4
Chemical Name: dimethyl (1E)-1-methyl-3-(methylamino)-3-oxo-1-propen-1-yl phosphate

PS: Solid


MF: C₇H₁₄NO₅P

MW: 223.17

MP: 54°C


Part Number	Matrix	Size
N-12493-100MG	Neat	100MG
S-12493A1-1ML	100ug/ml in Acetonitrile	1ML
S-12493U1-1ML	100ug/ml in Toluene	1ML
S-12493U1-5ML	100ug/ml in Toluene	5ML


2,3,5,6-Tetrachloroanisole Solution
CAS#6936-40-9
Chemical Name: 1,2,4,5-TETRACHLORO-3-METHOXYBENZENE

PS: Liquid


Part Number	Matrix	Size
S-13984J1-1ML	100ug/ml in Hexane	1ML

Ormetoprim
CAS#6981-18-6
Chemical Name: 2,4-Diamino-5-(6-methylveratryl)pyrimidine


PS: Solid

MF: C₁₄H₁₈N₄O₂
MW: 274.33


Part Number	Matrix	Size
N-12728-250MG	Neat	250MG
S-12728A1-1ML	100ug/ml in Acetonitrile	1ML
S-12728U1-1ML	100ug/ml in Toluene	1ML

Dioxacarb
CAS#6988-21-2
Chemical Name: 2-(1,3-Dioxolan-2-yl)phenyl methylcarbamate

PS: Solid


MF: C₁₁H₁₃NO₄

Alternate Name: 2-(1,3-Dioxolan-2-yl)phenyl methylcarbamate

Part Number	Matrix	Size
N-11789-250MG	Neat	250MG
S-11789M4-1ML	100ug/ml in Methanol	1ML

Mecoprop
CAS#7085-19-0
Chemical Name: 2-(4-Chloro-2-methylphenoxy)propionic acid

PS: Solid


MF: C₁₀H₁₁ClO₃
MW: 214.6

MP: 94°C

Alternate Name: MCP

Part Number	Matrix	Size
N-12360-100MG	Neat	100MG
S-12360A1-1ML	100ug/ml in Acetonitrile	1ML
S-12360B1-1ML	100ug/ml in Acetone	1ML

Sebuthylazin
CAS#7286-69-3
Chemical Name: 6-chloro-N-ethyl-N'-(1-methylpropyl)-1,3,5-triazine-2,4-diamine


PS: Solid

MF: C₉H₁₆ClN₅
MW: 229.71


Part Number	Matrix	Size
N-12950-250MG	Neat	250MG

Chloramben methyl ester
CAS#7286-84-2
Chemical Name: 3-Amino-2,5-dichloro benzoic acid methyl ester

PS: Solid

MF: C₈H₇Cl₂NO₂
MW: 220.06



Part Number	Matrix	Size
N-11419-100MG	Neat	100MG
S-11419Y1-1ML	100ug/ml in Isooctane:Acetone (90:10)	1ML
S-11419Y1-5ML	100ug/ml in Isooctane:Acetone (90:10)	5ML

Prometryne
CAS#7287-19-6
Chemical Name: 2,4-Bis[isopropylamino]-6-methylthio-s-triazine

PS: Solid

MF: C₁₀H₁₉N₅S

MW: 241.36


MP: 119°C


Part Number	Matrix	Size
N-13103-250MG	Neat	250MG
S-13103A1-1ML	100ug/ml in Acetonitrile	1ML
S-13103T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-13103T1-5ML	100ug/ml in t-Butyl methyl ether	5ML

Monalide
CAS#7287-36-7
Chemical Name: Pentanamide, N-(4-chlorophenyl)-2,2-dimethyl-


PS: Solid

MF: C₁₃H₁₈ClNO

MW: 239.7


Part Number	Matrix	Size
N-13917-10MG	Neat	10MG

Propaphos Solution
CAS#7292-16-2
Chemical Name: 4-(methylthio)phenyl dipropyl phosphate


PS: Liquid


Part Number	Matrix	Size
S-13967A1-1ML	100ug/ml in Acetonitrile	1ML

Endrin aldehyde
CAS#7421-93-4
Chemical Name: 1,2,4-Methenocyclopenta[cd]pentalene-5-carboxaldehyde

PS: Solid

MF: C₁₂H₈Cl₆O

MW: 380.89



Part Number	Matrix	Size
N-11855-10MG	Neat	10MG
S-11855M1-1ML	100ug/ml in Methanol	1ML
S-11855M1-5ML	100ug/ml in Methanol	5ML

Thallium sulfate(Technical)
CAS#7446-18-6
Chemical Name: 2-Ethylamino-4-methylthio-6-tert-butylamino-1,3,5-triazine

PS: Solid

MF: H₂O₄S₂Tl

MW: 504.83


MP: 632°C


Part Number	Matrix	Size
N-13559-250MG	Neat	250MG

Tetramethrin
CAS#7696-12-0
Chemical Name: 3,4,5,6-Tetrahydrophthalimidomethylchrysanthemate

PS: Solid

MF: C₁₉H₂₅NO₄
MW: 331.41


MP: 68°C


Part Number	Matrix	Size
N-13552-100MG	Neat	100MG
S-13552A4-1ML	1000ug/ml in Acetonitrile	1ML
S-13552A4-5ML	1000ug/ml in Acetonitrile	5ML
S-13552U1-1ML	100ug/mL in Toluene	1ML

Crotoxypfos
CAS#7700-17-6
Chemical Name: 1-phenylethyl (2E)-3-[(dimethoxyphosphinyl)oxy]-2-butenate

PS: Liquid

MF: C₁₄H₁₉O₆P

MW: 314.28


Part Number	Matrix	Size
N-11513-50MG	Neat	50MG
S-11513K1-1ML	100ug/ml in Isooctane	1ML
S-11513K1-5ML	100ug/ml in Isooctane	5ML

Ammonium sulfamate(Technical)
CAS#7773-06-0
PS: Solid

MF: H₃NO₃S₂H₃N

MW: 114.13


MP: 131°C

Alternate Name: Ammate (TM)

Part Number	Matrix	Size
N-11070-500MG	Neat	500MG

Sodium arsenate
CAS#7778-43-0
PS: Solid

MF: AsH₃Na₂O₄
MW: 185.91

MP: 57°C


Part Number	Matrix	Size
N-13215-250MG	Neat	250MG
S-13215F1-1ML	100ug/ml in H2O	1ML
S-13215U1-1ML	100ug/ml in Toluene	1ML

Zinc phosphate(Technical)
CAS#7779-90-0
PS: Solid

MF: H₃O₄P₃/2Zn


Part Number	Matrix	Size
N-13758-250MG	Neat	250MG

Lead arsenate(Technical)
CAS#7784-40-9
PS: Solid

MF: AsHO₄Pb

MW: 347.1

MP: 280°C


Alternate Name: Gypsine

Part Number	Matrix	Size
N-12310-100MG	Neat	100MG

Mevinphos
CAS#7786-34-7
Chemical Name: 2-Carbomethoxy-1-methylvinyl dimethyl phosphate

PS: Liquid

MF: C₇H₁₃O₆P


MW: 224.16

Alternate Name: Phosdrin

Part Number	Matrix	Size
N-13037-250MG	Neat	250MG
S-13037A1-1ML	100ug/ml in Acetonitrile	1ML
S-13037U1-1ML	100ug/ml in Toluene	1ML
S-13037U1-5ML	100ug/ml in Toluene	5ML

Camphechlor
CAS#8001-35-2
PS: Liquid

MF: C₁₀H₁₀Cl₈ (Approximate)

MW: 414

MP: 65°C

Alternate Name: Toxaphene

Part Number	Matrix	Size
N-13586-250MG	Neat	250MG
S-13586B6-1ML	2500ug/ml in Acetone	1ML
S-13586J4-1ML	1000ug/ml in Hexane	1ML
S-13586J4-5ML	1000ug/ml in Hexane	5ML
S-13586M1-1ML	100ug/ml in Methanol	1ML
S-13586M1-5ML	100ug/ml in Methanol	5ML

Strobane (TM)
CAS#8001-50-1
Chemical Name: Heptachloro-2,2-dimethyl-3-methylene norbornane

PS: Liquid

Alternate Name: Terpene polychlorinates

Part Number	Matrix	Size
N-13229-100MG	Neat	100MG
S-13229A1-1ML	100ug/ml in Acetonitrile	1ML
S-13229K4-1ML	1000ug/ml in Isooctane	1ML


Pyrethrum
CAS#8003-34-7
PS: Liquid

MW: 316

Part Number	Matrix	Size
N-13151-100MG	Neat	100MG
S-13151A1-1ML	100ug/ml in Acetonitrile	1ML
S-13151T1-1ML	100ug/ml t-Butyl methyl ether	1ML


DDT (mixture p,p' & o,p')
CAS#8017-34-3
Chemical Name: Bis(chlorophenyl)-trichloroethane (mixed isomers)

PS: Solid

MF: C₁₄H₉Cl₅
MW: 354.48



Part Number	Matrix	Size
N-11567-250MG	Neat	250MG
S-11567M1-1ML	100ug/ml in Methanol	1ML

Demeton O&S
CAS#8065-48-3
PS: Liquid

MF: C₈H₁₉O₃PS₂.C₈H₁₉O₃PS₂
MW: 258.34

Alternate Name: Systox

Part Number	Matrix	Size
N-11581-100MG	Neat	100MG
S-11581U1-1ML	100ug/ml in Toluene	1ML


Metaldehyde(Technical)
CAS#9002-91-9
PS: Solid

MF: (C₂H₄O)_x
MP: 246°C

Alternate Name: ACETALDEHYDE POLYMERIZED

Part Number	Matrix	Size
N-12381-1G	Neat	1G

Butoxypolypropylene glycol(Technical)
CAS#9003-13-8
PS: Liquid

MF: (C₃H₆O)_xMULTC₄H₁₀O

MW: 400

Alternate Name: Stabilene

Part Number	Matrix	Size
N-11354-250MG	Neat	250MG

Metiram(Technical)
CAS#9006-42-2
Chemical Name: zinc ammoniate ethylenebis(dithiocarbamate) - poly(ethylenethiuram disulfide)

PS: Solid

MF: C₁₆H₃₃N₁₁S₁₆Zn₃
MW: 1088.7


MP: 156°C

Alternate Name: Polyram (TM)

Part Number	Matrix	Size
N-13086-250MG	Neat	250MG

Propineb(Technical)
CAS#9016-72-2
Chemical Name: Zinc propylene bis(dithiocarbamate)polymer

PS: Solid


MF: (C₅H₈N₂S₄Zn)_x
FP: 70°C


Part Number	Matrix	Size
N-13111-250MG	Neat	250MG

Hymexazole
CAS#10004-44-1
Chemical Name: 3-Hydroxy-5-methylisoxazole

PS: Solid

MF: C₄H₅NO₂
MW: 99

MP: 86°C

Alternate Name: Bucid(TM)


Part Number	Matrix	Size
N-12196-25MG	Neat	25MG
S-12196M1-1ML	100ug/ml in Methanol	1ML

Methamidophos
CAS#10265-92-6
Chemical Name: O,S-Dimethylphosphoramidothioate

PS: Solid

MF: C₂H₈NO₂PS


MW: 141.1

MP: 44.9°C


Part Number	Matrix	Size
N-12393-100MG	Neat	100MG
S-12393A1-1ML	100ug/ml in Acetonitrile	1ML
S-12393U1-1ML	100ug/ml in Toluene	1ML

Copper-8-Quinolinolate
CAS#10380-28-6
PS: Solid

MF: C₁₈H₁₂CuN₂O₂
MW: 351.86


MP: 270°C

Alternate Name: Bis(quinolin-8-olato-O,N)copper

Part Number	Matrix	Size
N-11503-500MG	Neat	500MG

Resmethrin
CAS#10453-86-8
Chemical Name: 5-Benzyl-3-furylmethyl(1RS)-cis-trans-chrysanthemate

PS: Solid

MF: C₂₂H₂₆O₃
MW: 338.4

MP: 56.5°C


Part Number	Matrix	Size
N-13176-250MG	Neat	250MG


S-13176A1-1ML	100ug/ml in Acetonitrile	1ML
S-13176T4-1ML	1000ug/ml in t-Butyl methyl ether	1ML
S-13176T4-5ML	1000ug/ml in t-Butyl methyl ether	5ML

Nitrothal isopropyl
CAS#10552-74-6
Chemical Name: 1,7,3-?Benzenedicarboxylic acid, 5-?nitro-?, 1,7,3-?bis(1-?methylethyl) ester

PS: Solid

MF: C₁₄H₁₇NO₆
MW: 295.3


MP: 65°C

FP: 400°C

Alternate Name: Nitrothale

Part Number	Matrix	Size
N-12665-100MG	Neat	100MG

Carbendazim
CAS#10605-21-7
Chemical Name: 2-Methoxycarbonylamino)-benzimidazole


PS: Solid

MF: C₉H₉N₃O₂
MW: 191.19


Part Number	Matrix	Size
N-11404-100MG	Neat	100MG
S-11404M1-1ML	100ug/ml in Methanol	1ML
S-11404M1-5ML	100ug/ml in Methanol	5ML

Azadirachtin(Technical)
CAS#11141-17-6
PS: Solid

MF: C₃₅H₄₄O₁₆
MW: 720.7

MP: 154°C

Alternate Name: Azadiractin


Part Number	Matrix	Size
N-11107-10MG	Neat	10MG

Zineb(Technical)
CAS#12122-67-7
Chemical Name: Ethylene bis(dithiocarbamic acid)-zinc salt

PS: Solid

MF: C₄H₆N₂S₄Zn


MW: 275.75

MP: 157°C


Part Number	Matrix	Size
N-13759-250MG	Neat	250MG

Maneb(Technical)
CAS#12427-38-2
Chemical Name: Ethylene bis(dithiocarbamic acid)-manganese salt

PS: Solid


MF: C₄H₆MnN₂S₄
MW: 265.31


Part Number	Matrix	Size
N-12355-250MG	Neat	250MG

Ancymidol
CAS#12771-68-5
Chemical Name: a-Cyclopropyl-a-(p-methoxyphenyl)-5-pyrimidine methanol

PS: Solid

MF: C₁₅H₁₆N₂O₂
MW: 256.33


MP: 110°C


Part Number	Matrix	Size
N-11073-100MG	Neat	100MG
S-11073M1-1ML	100ug/ml in Methanol	1ML

Cyanofenphos
CAS#13067-93-1
Chemical Name: O-(4-cyanophenyl) O-ethyl P-phenylphosphonothioate

PS: Liquid

MF: C₁₅H₁₄NO₂PS


MW: 303.32


Part Number	Matrix	Size
N-12881-50MG	Neat	50MG
S-12881M1-1ML	100ug/ml in Methanol	1ML

Terbufos
CAS#13071-79-9
Chemical Name: O,O-Diethyl S-[[[(1,1-dimethylethyl)thio]methyl] phosphorodithioate

PS: Liquid

MF: C₉H₂₁O₂PS₃
MW: 288.4


MP: -29°C


Part Number	Matrix	Size
N-13510-250MG	Neat	250MG
S-13510M1-1ML	100ug/ml in Methanol	1ML

Tricyclohexyltin hydroxide
CAS#13121-70-5
PS: Solid

MF: C₁₈H₃₄OSn

MW: 385.18

MP: 196°C

Alternate Name: Cyhexatin


Part Number	Matrix	Size
N-13662-250MG	Neat	250MG
S-13662M1-1ML	100ug/ml in Methanol	1ML

Phosphamidon
CAS#13171-21-6
Chemical Name: 2-Chlorodiethylcarbamoyl methylvinyl dimethyl phosphate

PS: Liquid

MF: C₁₀H₁₉ClNO₅P


MW: 299.69

MP: -45°C


Part Number	Matrix	Size
N-13039-250MG	Neat	250MG
S-13039A1-1ML	100ug/ml in Acetonitrile	1ML
S-13039U1-1ML	100ug/ml in Toluene	1ML
S-13039U1-5ML	100ug/ml in Toluene	5ML

Bromofenoxim
CAS#13181-17-4
Chemical Name: 3,5-dibromo-4-hydroxybenzaldehyde O-(2,4-dinitrophenyl)oxime

PS: Liquid

MF: C₁₃H₇Br₂N₃O₆



Part Number	Matrix	Size
N-12961-10MG	Neat	10MG

Prophos CAS#13194-48-4

Chemical Name: O-Ethyl-S,S-Dipropylphosphorodithioate

PS: Liquid

MF: C₈H₁₉O₂PS₂
MW: 242.34

MP: -13°C

Alternate Name: Ethoprophos (TM)

Part Number	Matrix	Size
N-13110-250MG	Neat	250MG
S-13110A1-1ML	100ug/ml in Acetonitrile	1ML
S-13110U1-1ML	100ug/ml in Toluene	1ML
S-13110U1-5ML	100ug/ml in Toluene	5ML


Fenbutatin-oxide(Technical) CAS#13356-08-6

Chemical Name: Hexakis-(2-methyl-2-phenylpropyl)distanoxane

PS: Solid

MF: C₆₀H₇₈OSn₂
MW: 1052.7

MP: 138°C

FP: 100°C



Part Number	Matrix	Size
N-11951-250MG	Neat	250MG

Chlorbromuron CAS#13360-45-7

Chemical Name: 3-[4-Bromo-3-chlorophenyl]-1-methoxy-1-methylurea

PS: Solid

MF: C₉H₁₀BrClN₂O₂
MW: 293.57

MP: 95°C



Part Number	Matrix	Size
N-11424-250MG	Neat	250MG
S-11424A1-1ML	100ug/ml in Acetonitrile	1ML
S-11424U1-1ML	100ug/ml in Toluene	1ML

Pyrazophos CAS#13457-18-6

Chemical Name: ethyl 2-diethoxyphosphinothioxy-5-methylpyrazolo[1,5-a]pyrimidine-6-carboxylate

PS: Solid

MF: C₁₄H₂₀N₃O₅PS

MW: 373.4

Alternate Name: Afugan(TM)

Part Number	Matrix	Size
N-13146-100MG	Neat	100MG
S-13146A1-1ML	100ug/ml in Acetonitrile	1ML
S-13146U1-1ML	100ug/ml in Toluene	1ML


Quinalphos CAS#13593-03-8

Chemical Name: O,O-diethyl-O-(2-quinoxalyl)phosphorothioate

PS: Solid

MF: C₁₂H₁₅N₂O₃PS

MW: 289.32

MP: 32°C

Alternate Name: Bayrusil


Part Number	Matrix	Size
N-13168-250MG	Neat	250MG
S-13168A1-1ML	100ug/ml in Acetonitrile	1ML
S-13168U1-1ML	100ug/ml in Toluene	1ML

Desmedipham CAS#13684-56-5

Chemical Name: Ethyl-m-hydroxycarbanilate(ester)

PS: Solid

MF: C₁₆H₁₆N₂O₄
MW: 300.34

MP: 120°C



Part Number	Matrix	Size
N-11584-250MG	Neat	250MG
S-11584M1-1ML	100ug/ml in Methanol	1ML

Phenmedipham CAS#13684-63-4

Chemical Name: 3-[(methoxycarbonyl)amino]phenyl N-(3-methylphenyl)carbamate

PS: Solid

MF: C₁₆H₁₆N₂O₄
MW: 300.3

MP: 143°C



Part Number	Matrix	Size
N-12859-250MG	Neat	250MG
S-12859A1-1ML	100ug/ml in Acetonitrile	1ML
S-12859U1-1ML	100ug/ml in Toluene	1ML

Picloram methyl ester CAS#14143-55-6

Chemical Name: 4-Amino-3,5,6-trichloropicolinic acid methyl ester

PS: Solid


MF: C₇H₅Cl₃N₂O₂
MW: 255.49

MP: 114°C


Part Number	Matrix	Size
N-13051-100MG	Neat	100MG
S-13051M1-1ML	100ug/ml in Methanol	1ML

Difenoxyuron
CAS#14214-32-5
Chemical Name: 3-[4-(4-methoxyphenoxy)phenyl]-1,1-dimethylurea


PS: Liquid

MF: C₁₆H₁₈N₂O₃


Part Number	Matrix	Size
N-12977-10MG	Neat	10MG

Chlorfenprop-methyl
CAS#14437-17-3
Chemical Name: methyl 2,4-dichlorobenzenepropanoate

PS: Liquid


MF: C₁₀H₁₀Cl₂O₂
MW: 233.09


Part Number	Matrix	Size
N-12898-50MG	Neat	50MG

Ferbam
CAS#14484-64-1
Chemical Name: Dimethyl dithiocarbamic acid-ferric salt

PS: Solid

MF: C₉H₁₈FeN₂S₆
MW: 416.5


MP: 180°C

Alternate Name: Fermate

Part Number	Matrix	Size
N-11970-250MG	Neat	250MG
S-11970M1-1ML	100ug/ml in Methanol	1ML
S-11970U1-1ML	100ug/ml in Toluene	1ML

Oxolinic acid
CAS#14698-29-4
Chemical Name: 5-Ethyl-5,8-dihydro-8-oxo[1,3]dioxolo[4,5-g]quinoline-7-carboxylic acid

PS: Solid

MF: C₁₃H₁₁NO₅
MW: 261.24

MP: 314°C

Alternate Name: Dioxacin(TM)


Part Number	Matrix	Size
N-12738-250MG	Neat	250MG
S-12738A1-1ML	100ug/ml in Acetonitrile	1ML

Phoxim
CAS#14816-18-3
Chemical Name: O,O-Diethyl a-cyanobenzylideneamino-oxyphosphonothioate

PS: Liquid

MF: C₁₂H₁₅N₂O₃PS

MW: 298.32


MP: 6.1°C


Part Number	Matrix	Size
N-13041-100MG	Neat	100MG
S-13041A1-1ML	100ug/ml in Acetonitrile	1ML
S-13041U1-1ML	100ug/ml in Toluene	1ML

Simeconazole
CAS#14958-90-7
Chemical Name: (RS)-2-(4-fluorophenyl)-1-(1H-1,2,4-triazol-1-yl)-3-(trimethylsilyl)propan-2-ol

PS: Solid

MF: C₁₄H₂₀FN₃OSi


MW: 293.4


Part Number	Matrix	Size
N-13931-25MG	Neat	25MG

Dichlorprop-P
CAS#15165-67-0
Chemical Name: (+)-2,4-DP/R-2-(2,4-Dichlorophenoxy)propionic acid

PS: Solid

MF: C₉H₈Cl₂O₃
MW: 235.07

MP: 117.5°C

Alternate Name: Duplosan(TM)


Part Number	Matrix	Size
N-11674-250MG	Neat	250MG
S-11674A1-1ML	100ug/ml in Acetonitrile	1ML
S-11674U1-1ML	100ug/ml in Toluene	1ML

Cartap hydrochloride
CAS#15263-52-2
Chemical Name: S,S'-(2-Dimethylaminotrimethylene)bis(thiocarbamate) hydrochloride

PS: Solid

MF: C₇H₁₅N₃O₂S₂ClH

MW: 273.83

MP: 180°C



Part Number	Matrix	Size
N-11412-250MG	Neat	250MG

Napropamide CAS#15299-99-7
Chemical Name: 2-[a-Naphthoxy]-N,N-diethylpropionamide

PS: Solid

MF: C₁₇H₂₁NO₂
MW: 271.36

MP: 75°C


FP: 104°C


Part Number	Matrix	Size
N-11585-250MG	Neat	250MG
S-11585A1-1ML	100ug/ml in Acetonitrile	1ML
S-11585T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-11585T1-5ML	100ug/ml in t-Butyl methyl ether	5ML

Benodanil CAS#15310-01-7
Chemical Name: 2-iodo-N-phenylbenzamide

PS: Solid

MF: C₁₃H₁₀INO


MW: 323.13


Part Number	Matrix	Size
N-12991-10MG	Neat	10MG

Fluorodifen CAS#15457-05-3
Chemical Name: ((4-Amino-3,5-dichloro-6-fluoro-2-pyridinyl)oxy)-,1-methylheptyl ester acetic acid

PS: Solid

MF: C₁₃H₇F₃N₂O₅
MW: 328.22

MP: 94°C



Part Number	Matrix	Size
N-11998-250MG	Neat	250MG
S-11998M1-1ML	100ug/ml in Methanol	1ML

Chlorotoluron CAS#15545-48-9
Chemical Name: 3-(3-Chloro-p-tolyl)-1,1-dimethylurea

PS: Solid

MF: C₁₀H₁₃ClN₂O

MW: 212.7


MP: 147°C

Alternate Name: Dicuran*

Part Number	Matrix	Size
N-11455-250MG	Neat	250MG
S-11455A1-1ML	100ug/ml in Acetonitrile	1ML
S-11455U1-1ML	100ug/ml in Toluene	1ML

a-Chloralose(Technical) CAS#15879-93-3
Chemical Name: 1,2-O-(2,2,2-Trichloroethylidene)-?-D-glucufuranose

PS: Solid

MF: C₈H₁₁Cl₃O₆
MW: 309.54


MP: 187°C

Alternate Name: Alfa-4(TM)

Part Number	Matrix	Size
N-10978-250MG	Neat	250MG

Alachlor CAS#15972-60-8
Chemical Name: 2-Chloro-2,6-diethyl-N-methoxymethylacetanilide

PS: Solid

MF: C₁₄H₂₀ClNO₂
MW: 269.77


MP: 40°C


Part Number	Matrix	Size
N-11043-250MG	Neat	250MG
S-11043M1-1ML	100ug/ml in Methanol	1ML
S-11043M1-5ML	100ug/ml in Methanol	5ML
S-FC2072S-1.2ML	100ug/ml in n-Nonane (C13)	1.2ML

Agribrom CAS#16079-88-2
Chemical Name: 1-Bromo-3-chloro-5,5-dimethyl-2,4-imidazolidinedione

PS: Solid


MF: C₅H₆BrClN₂O₂
MW: 241.49

MP: 163°C


Part Number	Matrix	Size
N-11041-500MG	Neat	500MG

Carbetamide CAS#16118-49-3
Chemical Name: (2R)-N-ethyl-2-[[[(phenylamino)carbonyl]oxy]propanamide


PS: Liquid

MF: C₁₂H₁₆N₂O₃


Part Number	Matrix	Size
N-12971-10MG	Neat	10MG

b-Chloralose(Technical)
CAS#16376-36-6

PS: Solid
MF: C₈H₁₁Cl₃O₆




Part Number	Matrix	Size
N-11117-100MG	Neat	100MG

Mecoprop-P
CAS#16484-77-8

(R)-2-(4-chloro-2-methylphenoxy)propanoic acid

PS: Solid
MF: C₁₀H₁₁ClO₃
MW: 214.6
MP: 94°C




Alternate Name: Astix(TM)

Part Number	Matrix	Size
N-12363-250MG	Neat	250MG
S-12363A1-1ML	100ug/ml in Acetonitrile	1ML
S-12363T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

(2-Chloroethyl)phosphonic acid
CAS#16672-87-0

PS: Solid
MF: C₂H₆ClO₃P
MW: 144.49
MP: 74°C




Alternate Name: Ethephon

Part Number	Matrix	Size
N-10002-250MG	Neat	250MG

Methomyl
CAS#16752-77-5

Chemical Name: S-Methyl-N-[(methylcarbamoyl)oxy]thioacetamide

PS: Solid
MF: C₅H₁₀N₂O₂S
MW: 162.21
MP: 78°C




Part Number	Matrix	Size
N-12399-100MG	Neat	100MG
S-12399A1-1ML	100ug/ml in Acetonitrile	1ML

Ammonium fluosilicate(Technical)
CAS#16919-19-0

Chemical Name: Ammonium hexafluorosilicate

PS: Solid
MF: F₆Si₂H₄N
MW: 178.15




Part Number	Matrix	Size
N-11069-500MG	Neat	500MG

Dioxydemeton-S-methyl
CAS#17040-19-6

Chemical Name: O,O-Dimethyl-S-[2-ethylsulfonylethyl]phosphorothioate

PS: Solid
MF: C₆H₁₅O₅PS₂
MW: 262.3
MP: 48.8°C




Part Number	Matrix	Size
N-11791-100MG	Neat	100MG
S-11791A1-1ML	100ug/ml in Acetonitrile	1ML

Edifenphos
CAS#17109-49-8

Chemical Name: O-Ethyl-S,S-diphenyl phosphorodithioate

PS: Liquid
MF: C₁₄H₁₅O₂PS₂
MW: 310.38
MP: -25°C
FP: 115°C




Part Number	Matrix	Size
N-11844-100MG	Neat	100MG
S-11844A1-1ML	100ug/mL in Acetonitrile	1ML
S-11844U1-1ML	100ug/ml in Toluene	1ML

Bensultap
CAS#17606-31-4

Chemical Name: S,S',-2-Dimethylaminotrimethylenedi(benzenethiosulfonate)

PS: Solid
MF: C₁₇H₂₁NO₄S₄
MW: 431.6



Part Number	Matrix	Size
N-11141-250MG	Neat	250MG
S-11141M1-1ML	100ug/ml in Methanol	1ML

Dalapon methyl ester
CAS#17640-02-7
Chemical Name: 2,2-Dichloropropanoic acid methyl ester


PS: Liquid


Part Number	Matrix	Size
N-11563-100MG	Neat	100MG
S-11563M1-1ML	100ug/ml in Methanol	1ML

Benomyl
CAS#17804-35-2
Chemical Name: Methyl-1-[butylcarbamoyl]-2-benzimidazole carbamate

PS: Solid

MF: C₁₄H₁₈N₄O₃
MW: 290.32

MP: 140°C



Part Number	Matrix	Size
N-11138-100MG	Neat	100MG

Bromopropylate
CAS#18181-80-1
Chemical Name: Isopropyl 4,4-dibromobenzilate

PS: Solid

MF: C₁₇H₁₆Br₂O₃
MW: 428.14

MP: 77°C

FP: 100°C

Alternate Name: Acarol (TM)


Part Number	Matrix	Size
N-10997-250MG	Neat	250MG
S-10997A1-1ML	100ug/ml in Acetonitrile	1ML

Dikegulac acid
CAS#18467-77-1
Chemical Name: 2,3:4,6-di-O-Isopropylidene-a-L-xylo-2-hexulofuranosonic acid

PS: Solid

MF: C₁₂H₁₈O₇
MW: 274.3

MP: 74°C

FP: 100°C



Part Number	Matrix	Size
N-11750-250MG	Neat	250MG
S-11750T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Noruron
CAS#18530-56-8
Chemical Name: 3a,α,4α,5α,7α,7a,α-N,N-Dimethyl-N'-(octahydro-4,7-methano-1H-inden-5-yl)urea/3-(Hexahydro-4,7-methanoindan-5-yl)-1,1-dimethylurea

PS: Solid


MF: C₁₃H₂₂N₂O

MW: 222.3

MP: 38°C


Part Number	Matrix	Size
N-12669-250MG	Neat	250MG
S-12669A1-1ML	100ug/ml in Acetonitrile	1ML
S-12669U1-1ML	100ug/ml in Toluene	1ML

2,4-DB methyl ester
CAS#18625-12-2
Chemical Name: Methyl 4-(2,4-dichlorophenoxy)butyrate

PS: Solid

Alternate Name: 4-(2,4-Dichlorophenoxy)butyric acid methyl ester


Part Number	Matrix	Size
N-10533-100MG	Neat	100MG
S-10533M1-1ML	100ug/ml in Methanol	1ML
S-10533M1-5ML	100ug/ml in Methanol	5ML

Methabenzthiazuron
CAS#18691-97-9
Chemical Name: 1-(1,3-Benzothiazol-2-yl)-1,3-dimethylurea

PS: Solid

MF: C₁₀H₁₁N₃OS

MW: 221.3

MP: 120°C

Alternate Name: Tribunil*


Part Number	Matrix	Size
N-12387-100MG	Neat	100MG
S-12387A1-1ML	100ug/ml in Acetonitrile	1ML
S-12387U1-1ML	100ug/ml Toluene	1ML

Isoxathion
CAS#18854-01-8
Chemical Name: O,O-diethyl O-5-phenylisoxazol-3-yl phosphorothioate

PS: Liquid

MF: C₁₃H₁₆NO₄PS

MW: 313.3


FP: 210°C


Part Number	Matrix	Size
N-13932-25MG	Neat	25MG

Oryzalin CAS#19044-88-3

Chemical Name: 3,5-Dinitro-N,N-dipropylsulfanilamide

PS: Solid
MF: C₁₂H₁₈N₄O₆S
MW: 346.36
MP: 137°C




Part Number	Matrix	Size
N-12729-100MG	Neat	100MG
S-12729A4-1ML	1000ug/ml in Acetonitrile	1ML

Oxadiazon CAS#19666-30-9

Chemical Name: 3-[2,4-dichloro-5-(1-methylethoxy)phenyl]-5-(1,1-dimethylethyl)-1,3,4-oxadiazol-2(3H)-one

PS: Solid
MF: C₁₅H₁₈Cl₂N₂O₃
MW: 345.23
MP: 90°C




Part Number	Matrix	Size
N-12731-250MG	Neat	250MG
S-12731A1-1ML	100ug/ml in Acetonitrile	1ML
S-12731U1-1ML	100ug/ml in Toluene	1ML

Metoxuron CAS#19937-59-8

Chemical Name: N-(3-Chloro-4-methoxyphenyl)-N',N'-dimethylurea

PS: Solid
MF: C₁₀H₁₃ClN₂O₂
MW: 228.7
MP: 124°C



Alternate Name: Deftor(TM)

Part Number	Matrix	Size
N-12480-250MG	Neat	250MG
S-12480A1-1ML	100ug/ml in Acetonitrile	1ML
S-12480U1-1ML	100ug/ml in Toluene	1ML

Metribuzin CAS#21087-64-9

Chemical Name: 4-Amino-6-t-butyl-3-[methylthio]-1,2,4-triazine-5[4H]-one

PS: Solid
MF: C₈H₁₄N₄OS
MW: 214.3
MP: 126.2°C




Part Number	Matrix	Size
N-12481-250MG	Neat	250MG
S-12481A1-1ML	100ug/ml in Acetonitrile	1ML
S-12481T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-12481T1-5ML	100ug/ml in t-Butyl methyl ether	5ML

TCMTB CAS#21564-17-0

Chemical Name: 2-(Thiocyanatomethylthio)benzothiazole

PS: Solid
MF: C₉H₆N₂S₃
MW: 238.36
MP: 36°C



Alternate Name: Busan 30A

Part Number	Matrix	Size
N-13500-10MG	Neat	10MG
S-13500A1-1ML	100ug/ml in Acetonitrile	1ML
S-13500X4-1ML	1000ug/ml in Methylene chloride	1ML
S-13500X4-5ML	1000ug/ml in Methylene chloride	5ML

Leptophos CAS#21609-90-5

Chemical Name: O-[4-Bromo-2,5-dichlorophenyl]-O-methylphenylphosphonothioate

PS: Solid
MF: C₁₃H₁₀BrCl₂O₂PS
MW: 412.07
MP: 70.2°C
FP: 222°C



Alternate Name: Phosvel

Part Number	Matrix	Size
N-12315-100MG	Neat	100MG
S-12315A1-1ML	100ug/ml in Acetonitrile	1ML
S-12315U1-1ML	100ug/ml in Toluene	1ML

Bladex CAS#21725-46-2

Chemical Name: 2-[(4-Chloro-6-ethylamino-s-triazin-2-yl)Amino]-2-methyl propionitrile

PS: Solid
MF: C₉H₁₃ClN₆
MW: 240.7
MP: 167.5°C



Alternate Name: Cyanazine

Part Number	Matrix	Size
N-11327-250MG	Neat	250MG
S-11327A1-1ML	100ug/ml in Acetonitrile	1ML
S-11327T4-1ML	1000ug/ml in t-Butyl methyl ether	1ML
S-11327T4-5ML	1000ug/ml in t-Butyl methyl ether	5ML

Benzoylprop ethyl CAS#22212-55-1

Chemical Name: Ethyl-N-benzoyl-N-(3,4-dichlorophenyl)-DL-alaninate

PS: Solid
MF: C₁₈H₁₇Cl₂NO₃
MW: 366.2




Part Number	Matrix	Size
N-11177-100MG	Neat	100MG
S-11177A4-1ML	1000ug/ml in Acetonitrile	1ML
S-11177A4-5ML	1000ug/ml in Acetonitrile	5ML
S-11177U1-1ML	100ug/ml in Toluene	1ML

Fenamiphos
CAS#22224-92-6
Chemical Name: Ethyl 3-methyl-4-[methylthio]phenyl-[1-methylethyl]phosphoramide

PS: Solid

MF: C₁₃H₂₂NO₃PS

MW: 303.4

MP: 49.2°C

Alternate Name: METHAPHENAMIPHOS


Part Number	Matrix	Size
N-12854-250MG	Neat	250MG
S-12854A1-1ML	100ug/ml in Acetonitrile	1ML
S-12854T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-12854T1-5ML	100ug/ml in t-Butyl methyl ether	5ML

Tetrachlorvinphos
CAS#22248-79-9
Chemical Name: Phosphoric acid, (1Z)-2-chloro-1-(2,4,6-trichlorophenyl)ethenyl dimethyl ester

PS: Solid

MF: C₁₀H₉Cl₄O₄P

MW: 365.95


MP: 97°C


Part Number	Matrix	Size
N-13537-250MG	Neat	250MG
S-13537A1-1ML	100ug/ml in Acetonitrile	1ML
S-13537X1-1ML	100ug/ml in Methylene chloride	1ML
S-13537X1-5ML	100ug/ml in Methylene chloride	5ML

Bendiocarb
CAS#22781-23-3
Chemical Name: 2,2-Dimethyl benzo-1,3-dioxol-4-yl methyl carbamate

PS: Solid


MF: C₁₁H₁₃NO₄
MW: 223.23

MP: 129°C


Part Number	Matrix	Size
N-11135-250MG	Neat	250MG
S-11135U1-1ML	100ug/ml in Toluene	1ML

Dimepax
CAS#22936-75-0
Chemical Name: N-(1,2-dimethylpropyl)-N'-ethyl-6-(methylthio)-1,3,5-triazine-2,4-diamine

PS: Liquid


MF: C₁₁H₂₁N₅S

Alternate Name: DIMETHAMETRYN

Part Number	Matrix	Size
N-12984-10MG	Neat	10MG

Cyprazine
CAS#22936-86-3
Chemical Name: 6-Chloro-N2-cyclopropyl-N4-isopropyl-1,3,5-triazine-2,4-diamine

PS: Solid

MF: C₉H₁₄ClN₅
MW: 227.7


MP: 167°C


Part Number	Matrix	Size
N-11546-50MG	Neat	50MG
S-11546A1-1ML	100ug/ml in Acetonitrile	1ML
S-11546B1-1ML	100ug/ml in Acetone	1ML

Prallethrin
CAS#23031-36-9
Chemical Name: (S)-2-Methyl-4-oxo-3-(2-propynyl)cyclopent-1-yl(1R)-cis-trans-chrysanthemate

PS: Liquid

MF: C₁₉H₂₄O₃
MW: 300.43


FP: 139°C


Part Number	Matrix	Size
N-13088-10MG	Neat	10MG
S-13088A1-1ML	100ug/ml in Acetonitrile	1ML
S-13088K1-1ML	100ug/ml in Isooctane	1ML

Pirimicarb
CAS#23103-98-2
Chemical Name: 2-[Dimethylamino]-5,6-dimethyl-4-pyrimidinyl dimethylcarbamate

PS: Solid

MF: C₁₁H₁₈N₄O₂
MW: 238.29


MP: 90.5°C


Part Number	Matrix	Size
N-13062-250MG	Neat	250MG
S-13062A1-1ML	100ug/ml in Acetonitrile	1ML
S-13062U1-1ML	100ug/ml in Toluene	1ML

Oxyamyl
CAS#23135-22-0
Chemical Name: methyl 2-(dimethylamino)-N-[(methylamino)carbonyl]oxy]-2-oxoethanimidothioate

PS: Solid

MF: C₇H₁₃N₃O₃S


MW: 219.26

Alternate Name: VYDATE

Part Number	Matrix	Size
N-12735-100MG	Neat	100MG
S-12735A1-1ML	100ug/ml in Acetonitrile	1ML

Butachlor
CAS#23184-66-9
Chemical Name: 2-Chloro-2,6-diethyl-N-[butoxymethyl]acetanilide

PS: Liquid

MF: C₁₇H₂₆ClNO₂
MW: 311.9


MP: -5°C


Part Number	Matrix	Size
N-11353-100MG	Neat	100MG
S-11353A1-1ML	100ug/ml in Acetonitrile	1ML
S-11353T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-11353T1-5ML	100ug/ml in t-Butyl methyl ether	5ML

Formetanate hydrochloride
CAS#23422-53-9
Chemical Name: N,N-dimethyl-N'-[3-[[[(methylamino)carbonyl]oxy]phenyl]methanimidamide hydrochloride

PS: Solid

MF: C₁₁H₁₅N₃O₂ClH

MW: 257.8



Part Number	Matrix	Size
N-12015-250MG	Neat	250MG
S-12015A1-1ML	100ug/ml in Acetonitrile	1ML

Pirimiphos-ethyl
CAS#23505-41-1
Chemical Name: O,O-Diethyl-O-[2-(diethylamino)-6-methyl-4-pyrimidinyl]phosphorothioate

PS: Liquid

MF: C₁₃H₂₄N₃O₃PS

MW: 333.38


MP: 15°C


Part Number	Matrix	Size
N-13063-250MG	Neat	250MG
S-13063U1-1ML	100ug/ml in Toluene	1ML

Heptenophos
CAS#23560-59-0
Chemical Name: 7-Chlorobicyclo(3,2,0)hepta-2,6-dien-6-yl dimethyl phosphate

PS: Liquid

MF: C₉H₁₂ClO₄P


MW: 250.6

Alternate Name: Heptenofos

Part Number	Matrix	Size
N-12153-100MG	Neat	100MG
S-12153A1-1ML	100ug/ml in Acetonitrile	1ML
S-12153U1-1ML	100ug/ml in Toluene	1ML

Thiophanate-methyl
CAS#23564-05-8
Chemical Name: dimethyl N,N'-[1,2-phenylenebis(iminocarbonothioyl)]bis[carbamate]

PS: Solid

MF: C₁₂H₁₄N₄O₄S₂
MW: 342.4


MP: 172°C


Part Number	Matrix	Size
N-13572-250MG	Neat	250MG
S-13572A1-1ML	100ug/ml in Acetonitrile	1ML

Thiophanate
CAS#23564-06-9
Chemical Name: Diethyl-4-4-O-Phenylene-bis[3-thioallophanate]

PS: Solid


MF: C₁₄H₁₈N₄O₄S₂
MW: 370.48

MP: 194°C



Part Number	Matrix	Size
N-13571-100MG	Neat	100MG
S-13571U1-1ML	100ug/ml in Toluene	1ML

Mecoprop methyl ester
CAS#23844-56-6
PS: Liquid

MF: C₁₁H₁₃ClO₃
MW: 228.7

FP: 100°C

Alternate Name: MCPP methyl ester

Part Number	Matrix	Size
N-12361-100MG	Neat	100MG
S-12361M1-1ML	100ug/ml in Methanol	1ML
S-12361M1-5ML	100ug/ml in Methanol	5ML


Mecoprop- 2,4,4-trimethylpentylester
CAS#217487-13-3
PS: Liquid

Alternate Name: MCPP 2,4,4-trimethylpentylester

Part Number	Matrix	Size
S-12362B1-1ML	100ug/ml in Acetone	1ML
S-12362B1-5ML	100ug/ml in Acetone	5ML

Ethirimol CAS#23947-60-6

Chemical Name: 5-butyl-2-ethylamino-6-methylpyrimidin-4-ol

PS: Liquid

MF: C₁₁H₁₉N₃O


Part Number	Matrix	Size
N-12986-10MG	Neat	10MG


Propyzamide CAS#23950-58-5

Chemical Name: N-[1,1-Dimethylpropynyl]-3,5-dichlorobenzamide

PS: Solid

MF: C₁₂H₁₁Cl₂NO

MW: 256.13

MP: 155°C

Alternate Name: Kerb (TM)

Part Number	Matrix	Size
N-13136-250MG	Neat	250MG
S-13136M1-1ML	100ug/mL in Methanol	1ML
S-13136M1-5ML	100ug/mL in Methanol	5ML


Triazophos CAS#24017-47-8

Chemical Name: O,O-Diethyl O-1-phenyl-1H-1,2,4-triazol-3-yl phosphorothioate

PS: Liquid

MF: C₁₂H₁₆N₃O₃PS

MW: 313.3


MP: 2°C

Alternate Name: Hostathion(TM)

Part Number	Matrix	Size
N-13639-100MG	Neat	100MG
S-13639A1-1ML	100ug/mL in Acetonitrile	1ML
S-13639U1-1ML	100ug/mL in Toluene	1ML
S-13639U1-5ML	100ug/mL in Toluene	5ML

Dimethachlon CAS#24096-53-5

Chemical Name: N-(3,5-Dichlorophenyl)succinimide

PS: Solid

MF: C₁₀H₇Cl₂NO₂

Alternate Name: Ohric


Part Number	Matrix	Size
N-11753-500MG	Neat	500MG
S-11753M1-1ML	100ug/mL in Methanol	1ML

Piperophos

CAS#24151-93-7

Chemical Name: S-2-methylpiperidinocarbonylmethyl O,O-dipropyl phosphorodithioate


PS: Liquid

MF: C₁₄H₂₈NO₃PS₂


Part Number	Matrix	Size
N-12989-250MG	Neat	250MG

2,6-Diisopropyl naphthalene CAS#24157-81-1

PS: Solid

MF: C₁₆H₂₀
MW: 212.33



Part Number	Matrix	Size
N-10550-250MG	Neat	250MG
S-10550A1-1ML	100ug/mL in Acetonitrile	1ML
S-10550U1-1ML	100ug/mL in Toluene	1ML

Mepiquat chloride CAS#24307-26-4

Chemical Name: N,N-Dimethylpiperidinium chloride

PS: Solid


MF: C₇H₁₆ClN

MW: 149.7


Part Number	Matrix	Size
N-12371-100MG	Neat	100MG
S-12371F1-1ML	100ug/mL in Water	1ML

Isocarbophos Solution CAS#24353-61-5


Chemical Name: O-Methyl O-[2-(isopropoxycarbonyl)phenyl] phosphoramidothioate

PS: Liquid


Part Number	Matrix	Size
S-13933A1-1ML	100ug/mL in Acetonitrile	1ML

Propamocarb
CAS#24579-73-5
Chemical Name: Propyl (3-(dimethylamino)propyl)carbamate


PS: Liquid

MF: C₉H₂₀N₂O₂


Part Number	Matrix	Size
N-13105-250MG	Neat	250MG

Fenfuram
CAS#24691-80-3
Chemical Name: 2-methyl-3-furanilide


PS: Solid

MF: C₁₂H₁₁NO₂
MW: 201.2


Part Number	Matrix	Size
N-14013-250MG	Neat	250MG
S-14013M1-1ML	100ug/mL in Methanol	1ML

Chlormephos
CAS#24934-91-6
Chemical Name: S-(Chloromethyl) O,O-diethyl ester phosphorodithioic acid

PS: Liquid

MF: C₅H₁₂ClO₂PS₂
MW: 234.7

Alternate Name: Dotan(TM)


Part Number	Matrix	Size
N-11434-250MG	Neat	250MG
S-11434A1-1ML	100ug/ml in Acetonitrile	1ML
S-11434U1-1ML	100ug/ml in Toluene	1ML

Bentazon
CAS#25057-89-0
Chemical Name: 3-(1-Methylethyl)-1H-2,1,3-benzothiadiazin-4(3H)-one-2,2-dioxide

PS: Solid

MF: C₁₀H₁₂N₂O₃S


MW: 240.3

MP: 139.4°C


Part Number	Matrix	Size
N-11142-250MG	Neat	250MG
S-11142A4-1ML	1000ug/ml in Acetonitrile	1ML
S-11142A4-5ML	1000ug/ml in Acetonitrile	5ML
S-11142B1-1ML	100ug/ml in Acetone	1ML

Benazolin-ethyl
CAS#25059-80-7
Chemical Name: Ethyl 4-chloro-2-oxo-3(2H)-benzothiazoleacetate


PS: Solid

MF: C₁₁H₁₀ClNO₃S

Alternate Name: Galtak(TM)

Part Number	Matrix	Size
N-11134-250MG	Neat	250MG
S-11134A1-1ML	100ug/ml in Acetonitrile	1ML
S-11134T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

2,4,5-T isooctyl ester
CAS#25168-15-4
Chemical Name: Isooctyl-2,4,5-trichlorophenoxyacetate

PS: Liquid


MF: C₁₆H₂₁Cl₃O₃
MW: 367.72


Part Number	Matrix	Size
N-10649-100MG	Neat	100MG
S-10649A1-1ML	100ug/ml in Acetonitrile	1ML
S-10649T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

2,4-D isooctyl ester
CAS#25168-26-7
Chemical Name: Isooctyl-2,4-dichlorophenoxyacetate

PS: Liquid

MF: C₁₆H₂₂Cl₂O₃
MW: 333.3

MP: 12°C



Part Number	Matrix	Size
N-10615-250MG	Neat	250MG
S-10615A1-1ML	100ug/ml in Acetonitrile	1ML
S-10615T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Isofenphos
CAS#25311-71-1
Chemical Name: 1-methylethyl 2-[[ethoxy[(1-methylethyl)amino]phosphinothioyl]oxy]benzoate

PS: Liquid

MF: C₁₅H₂₄NO₄PS

MW: 345.4

FP: °C



Part Number	Matrix	Size
N-12249-100MG	Neat	100MG
S-12249A1-1ML	100ug/ml in Acetonitrile	1ML
S-12249U1-1ML	100ug/ml in Toluene	1ML

Propamocarb HCl CAS#25606-41-1
Chemical Name: propyl N-[3-(dimethylamino)propyl]carbamate hydrochloride

PS: Solid

MF: C₉H₂₀N₂O₂ClH


MW: 224.77

FP: 118°C

Alternate Name: Propamocarb Hydrochloride


Part Number	Matrix	Size
N-13106-250MG	Neat	250MG
S-13106F1-1ML	100ug/ml in H2O	1ML

2,4,6-Tris-(2,4,6-tribromophenoxy)-1,3,5-triazine CAS#25713-60-4
PS: Solid

MF: C₂₁H₆Br₉N₃O₃
MW: 1067.43

MP: 230°C



Part Number	Matrix	Size
N-12888-50MG	Neat	50MG
S-12888U1-1ML	100ug/ml in Toluene	1ML

Ammonium ethylcarbamoylphosphonate Solution CAS#25954-13-6
PS: Liquid

Alternate Name: Ammonium ethylcarbamoylphosphonate

Part Number	Matrix	Size
S-12018F1-1ML	100 ug/mL in Water	1ML

d-(cis-trans)-Phenothrin CAS#26002-80-2
Chemical Name: 3-Phenoxybenzyl-d-cis,trans-chrysanthemate

PS: Liquid


MF: C₂₃H₂₆O₃
MW: 350.46

Alternate Name: Sumithrin

Part Number	Matrix	Size
N-11554-50MG	Neat	50MG
S-11554A4-1ML	1000ug/ml in Acetonitrile	1ML
S-11554U1-1ML	100ug/ml in Toluene	1ML

Iprobenfos CAS#26087-47-8
Chemical Name: S-Benzyl-O,O-di-isopropyl phosphorothioate

PS: Liquid

MF: C₁₃H₂₁O₃PS

MW: 288.3


Part Number	Matrix	Size
N-12219-250MG	Neat	250MG
S-12219A1-1ML	100ug/ml in Acetonitrile	1ML
S-12219U1-1ML	100ug/ml in Toluene	1ML


Ethofumesate CAS#26225-79-6
Chemical Name: 2-Ethoxy-2,3-dihydro-3,3-dimethyl-5-benzofuranyl methanesulfonate

PS: Solid

MF: C₁₃H₁₈O₃S

MW: 286.3

MP: 70°C


FP: 100°C


Part Number	Matrix	Size
N-11874-250MG	Neat	250MG
S-11874A1-1ML	100ug/ml in Acetonitrile	1ML
S-11874U1-1ML	100ug/ml in Toluene	1ML

Profluralin CAS#26399-36-0
Chemical Name: N-(cyclopropylmethyl)-2,6-dinitro-N-propyl-4-(trifluoromethyl)benzenamine

PS: Solid

MF: C₁₄H₁₆F₃N₃O₄
MW: 347.33


MP: 32.1°C


Part Number	Matrix	Size
N-13098-250MG	Neat	250MG
S-13098A1-1ML	100ug/ml in Acetonitrile	1ML
S-13098J4-1ML	1000ug/ml in Hexane	1ML
S-13098J4-5ML	1000ug/ml in Hexane	5ML

Octhilonone CAS#26530-20-1
Chemical Name: 2-octylisothiazol-3(2H)-one

PS: Liquid

MF: C₁₁H₁₉NOS


MW: 213.3


Part Number	Matrix	Size
N-14012-250MG	Neat	250MG
S-14012X1-1ML	100ug/mL in Methylene chloride	1ML

Triforine CAS#26644-46-2

Chemical Name: 1,4-bis(1-Foramide-2,2,2-trichloroethyl)-piperazine


PS: Solid
MF: C₁₀H₁₄Cl₆N₄O₂
MW: 434.97
MP: 155°C



Part Number	Matrix	Size
N-13691-250MG	Neat	250MG
S-13691U1-1ML	100ug/ml in Toluene	1ML

E-Mevinphos CAS#26718-65-0

PS: Liquid
MW: 224.17
MP: 21°C



Alternate Name: cis-Duraphos(TM)

Part Number	Matrix	Size
N-11481-100MG	Neat	100MG
S-11481A1-1ML	100ug/ml in Acetonitrile	1ML
S-11481T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Norflurazon CAS#27314-13-2

Chemical Name: 4-chloro-5-methylamino-2-(a,a,a-trifluoro-m-tolyl)pyridazin-3(2H)-one

PS: Solid
MF: C₁₂H₉ClF₃N₃O
MW: 303.67



Alternate Name: SOLICAM™

Part Number	Matrix	Size
N-12668-100MG	Neat	100MG
S-12668A1-1ML	100ug/ml in Acetonitrile	1ML
S-12668T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-12668T1-5ML	100ug/ml in t-Butyl methyl ether	5ML

Fthalide CAS#27355-22-2

Chemical Name: 4,5,6,7-Tetrachlorophthalide

PS: Solid
MF: C₈H₂Cl₄O₂
MW: 271.9



Alternate Name: PHTHALIDE

Part Number	Matrix	Size
N-13085-10MG	Neat	10MG

Probenazole CAS#27605-76-1

Chemical Name: 3-(2-Propenyloxy)-1,2-benzisothiazole-1,1-dioxide


PS: Solid
MF: C₁₀H₉NO₃S
MW: 223.26
MP: 138°C
FP: 177°C


Alternate Name: Oryzemat(TM)

Part Number	Matrix	Size
N-13092-500MG	Neat	500MG
S-13092M1-1ML	100ug/ml in Methanol	1ML

trans-Heptachlor epoxide CAS#28044-83-9

PS: Solid
MF: C₁₀H₅Cl₇O
MW: 389.4
FP: 162°C



Alternate Name: trans-Heptachlor epoxide (Isomer A)

Part Number	Matrix	Size
N-13618-50MG	Neat	50MG
S-13618M1-1ML	100ug/ml in Methanol	1ML

Bioallethrin CAS#28057-48-9

(RS)-3-allyl-2-methyl-4-oxocyclopent-2-enyl (1R,3R)-2,2-dimethyl-3-(2-methylprop-1-enyl)cyclopropanecarboxylate


PS: Liquid
MF: C₁₉H₂₆O₃
MW: 302.45
FP: 87°C


Alternate Name: d-trans-Allethrin

Part Number	Matrix	Size
N-11205-250MG	Neat	250MG
S-11205A1-1ML	100ug/ml in Acetonitrile	1ML
S-11205T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Phenylmercury salicylate(Technical) CAS#28086-13-7

PS: Solid
MF: C₁₃H₁₀HgO₃
MW: 414.72



Alternate Name: Mercusol

Part Number	Matrix	Size
N-13033-250MG	Neat	250MG
S-13033A1-1ML	100ug/ml in Acetonitrile	1ML
S-13033U1-1ML	100ug/ml in Toluene	1ML

Cybutryne
CAS#28159-98-0
Chemical Name: N-Cyclopropyl-N'-(1,1-dimethylethyl)-6-(methylthio)-1,3,5-Triazine-2,4-diamine

PS: Solid


MF: C₁₁H₁₉N₃S

MW: 253.37

Alternate Name: Irgarol(TM)

Part Number	Matrix	Size
N-11520-100MG	Neat	100MG
S-11520A1-1ML	100ug/ml in Acetonitrile	1ML
S-11520U1-1ML	100ug/ml in Toluene	1ML

N-Methyl-methamidophos
CAS#28167-49-9
PS: Solid

MF: C₃H₁₀NO₂PS

MW: 155.2

Alternate Name: N-Methyl-Tamaron(TM)


Part Number	Matrix	Size
N-12558-50MG	Neat	50MG
S-12558A1-1ML	100ug/ml in Acetonitrile	1ML
S-12558U1-1ML	100ug/ml in Toluene	1ML

Benthiocarb
CAS#28249-77-6
Chemical Name: S-(4-Chlorophenyl)methyl diethylcarbamothioate

PS: Liquid

MF: C₁₂H₁₆CINOS

MW: 257.8


MP: 3.3°C


Part Number	Matrix	Size
N-11144-250MG	Neat	250MG
S-11144M1-1ML	100ug/ml in Methanol	1ML

s-Bioallethrin
CAS#28434-00-6
Chemical Name: dl-Allylhydroxymethylcyclopentenone ester

PS: Liquid

MF: C₁₉H₂₆O₃
MW: 302.4

FP: 120°C



Part Number	Matrix	Size
N-13187-250MG	Neat	250MG
S-13187A1-1ML	100ug/ml in Acetonitrile	1ML
S-13187T1-1ML	100ug/ml t-Butyl methyl ether	1ML

Bioresmethrin
CAS#28434-01-7
Chemical Name: 5-Benzyl-3-furylmethyl-(+)-trans-chrysanthemate

PS: Liquid

MF: C₂₂H₂₆O₃
MW: 338.48

MP: 30°C

FP: 92°C

Alternate Name: d-trans-Resmethrin


Part Number	Matrix	Size
N-11206-250MG	Neat	250MG
S-11206A1-1ML	100ug/ml in Acetonitrile	1ML
S-11206T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Super-Rozol
CAS#28772-56-7
Chemical Name: 3-[3-(4?-bromo[1,1?-biphenyl]-4-yl)-3-hydroxy-1-phenylpropyl]-4-hydroxy-2H-1-benzopyran-2-one

PS: Solid

MF: C₃₀H₂₃BrO₄
MW: 527.42

MP: 200°C


FP: 218°C

Alternate Name: BROMADIOLON

Part Number	Matrix	Size
N-13259-100MG	Neat	100MG
S-13259A1-1ML	100ug/ml in Acetonitrile	1ML
S-13259U1-1ML	100ug/ml in Toluene	1ML

Dinitramine
CAS#29091-05-2
Chemical Name: N1,N1-diethyl-2,6-dinitro-4-(trifluoromethyl)-1,3-benzenediamine

PS: Liquid


MF: C₁₁H₁₃F₃N₄O₄
MW: 322.2

MP: 98°C


Part Number	Matrix	Size
N-12976-10MG	Neat	10MG

Prodiamine
CAS#29091-21-2
Chemical Name: 5-Dipropylamino-a,a-trifluoro-4,6-dinitro-o-toluidine

PS: Solid


MF: C₁₃H₁₇F₃N₄O₄
MW: 350.3

Alternate Name: Barricade(tm)

Part Number	Matrix	Size
N-13096-100MG	Neat	100MG
S-13096A1-1ML	100ug/ml in Acetonitrile	1ML
S-13096K1-1ML	100ug/ml in Isooctane	1ML

Pirimiphos-methyl
CAS#29232-93-7
Chemical Name: O,O-Dimethyl-O-[2-(diethylamino)-6-methyl-4-pyrimidinyl]phosphorothioate

PS: Liquid

MF: C₁₃H₂₄N₃O₃PS

MW: 305.3



Part Number	Matrix	Size
N-13064-250MG	Neat	250MG
S-13064A1-1ML	100ug/ml in Acetonitrile	1ML
S-13064U1-1ML	100ug/mL in Toluene	1ML

Ethiofencarb
CAS#29973-13-5
Chemical Name: a-Ethylthio-o-tolylmethyl carbamate

PS: Solid

MF: C₁₁H₁₅NO₂S

MW: 225.33


MP: 33.4°C

Alternate Name: Croneton*

Part Number	Matrix	Size
N-11872-100MG	Neat	100MG
S-11872A1-1ML	100ug/ml in Acetonitrile	1ML
S-11872U1-1ML	100ug/ml in Toluene	1ML

Ethidimuron
CAS#30043-49-3
Chemical Name: 1-(5-Ethylsulfonyl-1,3,4-thiadiazol-2-yl)-1,3-dimethylurea

PS: Solid

MF: C₇H₁₄N₄O₃S₂
MW: 264.35


MP: 156°C


Part Number	Matrix	Size
N-11871-100MG	Neat	100MG
S-11871A1-1ML	100ug/ml in Acetonitrile	1ML
S-11871U1-1ML	100ug/ml in Toluene	1ML

Acephate
CAS#30560-19-1
Chemical Name: O,S-Dimethyl acetylphosphoroamidothioate

PS: Solid

MF: C₄H₁₀NO₃PS


MW: 183.16


Part Number	Matrix	Size
N-11002-250MG	Neat	250MG
S-11002A1-1ML	100ug/ml in Acetonitrile	1ML
S-11002B4-1ML	1000ug/ml in Acetone	1ML
S-11002B4-5ML	1000ug/ml in Acetone	5ML

Methacrifos
CAS#30864-28-9
PS: Liquid

MF: C₇H₁₃O₃PS

MW: 240.23


FP: 69°C

Alternate Name: Damfin(TM)

Part Number	Matrix	Size
N-12388-100MG	Neat	100MG
S-12388A1-1ML	100ug/ml in Acetonitrile	1ML
S-12388U1-1ML	100ug/ml in Toluene	1ML

Propetamphos
CAS#31218-83-4
Chemical Name: 1-methylethyl (2E)-3-[[[(ethylamino)methoxyphosphinothioyl]oxy]-2-butenate

PS: Liquid


MF: C₁₀H₂₀NO₄PS

MW: 281.3

Alternate Name: Saftrotin(TM)

Part Number	Matrix	Size
N-13197-100MG	Neat	100MG
S-13197A1-1ML	100ug/ml in Acetonitrile	1ML
S-13197U1-1ML	100ug/ml in Toluene	1ML


Dodemorph acetate
CAS#31717-87-0
Chemical Name: 4-cyclododecyl-2,6-dimethylmorpholine acetate

PS: Liquid

MF: C₁₈H₃₅NO₂C₂H₄O₂


Part Number	Matrix	Size
N-12948-100MG	Neat	100MG


Silvex isooctyl ester
CAS#32534-95-5
PS: Liquid

MF: C₁₇H₂₃Cl₃O₃
MW: 381.75

Alternate Name: Pebulate

Part Number	Matrix	Size
N-13213-100MG	Neat	100MG

Procymidone
CAS#32809-16-8
Chemical Name: N-(3,5-Dichlorophenyl)-1,2-dimethylcyclopropane-1,2-dicarboximide

PS: Solid


MF: C₁₃H₁₁Cl₂NO₂
MW: 284.1

Alternate Name: Sumilex(TM)

Part Number	Matrix	Size
N-13095-250MG	Neat	250MG
S-13095A1-1ML	100ug/ml in Acetonitrile	1ML
S-13095U1-1ML	100ug/ml in Toluene	1ML

Procyazine
CAS#32889-48-8
Chemical Name: 2-[[4-Chloro-6-(cyclopropylamino)-1,3,5-triazin-2-yl]amino]-2-methylpropanenitrile

PS: Solid

MF: C₁₀H₁₃ClN₆
MW: 252.74


MP: 170°C


Part Number	Matrix	Size
N-13094-100MG	Neat	100MG
S-13094A1-1ML	100ug/ml in Acetonitrile	1ML
S-13094T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Amitraz
CAS#33089-61-1
Chemical Name: N-Methyl-N'-2,4-xylyl-N-(N-2,4-xylylformimidoyl)formamidine

PS: Solid

MF: C₁₉H₂₃N₃
MW: 293.45


MP: 86°C


Part Number	Matrix	Size
N-11068-250MG	Neat	250MG
S-11068A1-1ML	100ug/ml in Acetonitrile	1ML

b-Endosulfan
CAS#33213-65-9
PS: Solid

MF: C₉H₆Cl₆O₃S

MW: 406.93


MP: 213.3°C


Part Number	Matrix	Size
N-11118-100MG	Neat	100MG
S-11118A1-1ML	100ug/ml in Acetonitrile	1ML
S-11118J1-1ML	100ug/ml in Hexane	1ML
S-11118J1-5ML	100ug/ml in Hexane	5ML
N-FD203-5-5MG	Neat (Deuterated)	5MG
S-FD203S-1.2ML	100ug/ml in n-Nonane (Deuterated)	1.2ML


Fluchloralin
CAS#33245-39-5
Chemical Name: N-(2-Chloroethyl)-2,6-dinitro-N-propyl-a,a,a-trifluoro-p-toluidine

PS: Solid

MF: C₁₂H₁₃ClF₃N₃O₄
MW: 355.7

MP: 42°C


Part Number	Matrix	Size
N-11983-50MG	Neat	50MG
S-11983A1-1ML	100ug/ml in Acetonitrile	1ML
S-11983J4-1ML	1000ug/ml in Hexane	1ML
S-11983J4-5ML	1000ug/ml in Hexane	5ML

Oxy-chlordene Solution
CAS#33386-84-4
PS: Liquid


Part Number	Matrix	Size
S-14054J1-1ML	100ug/ml in Hexane	1ML

Butralin
CAS#33629-47-9
Chemical Name: 4-(1,1-Dimethylethyl)-N-(1-methylpropyl)-2,6-dinitrobenzenamine

PS: Solid


MF: C₁₄H₂₁N₃O₄
MW: 295.34

MP: 61°C

FP: 36°C


Alternate Name: Amex(TM)

Part Number	Matrix	Size
N-11355-500MG	Neat	500MG
S-11355A1-1ML	100ug/ml in Acetonitrile	1ML
S-11355K1-1ML	100ug/ml in Isooctane	1ML


Terbumeton
CAS#33693-04-8
Chemical Name: 2-Methoxy-4-ethylamino-6-tert-butylamino-s-triazine

PS: Solid

MF: C₁₀H₁₉N₅O

MW: 225.3



Part Number	Matrix	Size
N-13511-250MG	Neat	250MG
S-13511M1-1ML	100ug/ml in Methanol	1ML

Isopropalin
CAS#33820-53-0
Chemical Name: 2,6-Dinitro-N,N-di-n-propyl cumidene

PS: Liquid

MF: C₁₅H₂₃N₃O₄
MW: 309.37


MP: 21°C

FP: 36°C


Part Number	Matrix	Size
N-12265-250MG	Neat	250MG
S-12265A1-1ML	100ug/ml in Acetonitrile	1ML
S-12265J4-1ML	1000ug/ml in Hexane	1ML
S-12265J4-5ML	1000ug/ml in Hexane	5ML

Tebuthiuron
CAS#34014-18-1
PS: Solid

MF: C₉H₁₆N₄O₅
MW: 228.3

MP: 161.5°C

Alternate Name: Fluoroacetic acid-sodium salt


Part Number	Matrix	Size
N-13505-250MG	Neat	250MG
S-13505A1-1ML	100ug/ml in Acetonitrile	1ML
S-13505T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-13505T1-5ML	100ug/ml in t-Butyl methyl ether	5ML

Isoproturon
CAS#34123-59-6
Chemical Name: N,N-dimethyl-N'-[4-(1-methylethyl)phenyl]urea

PS: Solid

MF: C₁₂H₁₈N₂O


MW: 206.3

MP: 158°C

Alternate Name: Alon

Part Number	Matrix	Size
N-12279-100MG	Neat	100MG
S-12279A1-1ML	100ug/ml in Acetonitrile	1ML

Dimefuron
CAS#34205-21-5
Chemical Name: N'-[3-chloro-4-[5-(1,1-dimethylethyl)-2-oxo-1,3,4-oxadiazol-3(2H)-yl]phenyl]-N,N-dimethylurea

PS: Solid


MF: C₁₅H₁₉ClN₄O₃


Part Number	Matrix	Size
N-12972-10MG	Neat	10MG

Acetochlor
CAS#34256-82-1
Chemical Name: 2-Chloro-N-ethoxymethyl-N-(2-ethyl-6-methylphenyl)acetamide

PS: Liquid


MF: C₁₄H₂₀ClNO₂
MW: 269.8

MP: 10.6°C

Alternate Name: Harness(TM)

Part Number	Matrix	Size
N-11013-100MG	Neat	100MG
S-11013M1-1ML	100ug/ml in Methanol	1ML

Methothrin
CAS#34388-29-9
Chemical Name: 4-(methoxymethyl)benzyl (1RS)-cis,trans-2,2-dimethyl-3-(2-methylprop-1-enyl)cyclopropanecarboxylate

PS: Liquid


MF: C₁₉H₂₆O₃
MW: 302.4


Part Number	Matrix	Size
N-14037-100MG	Neat	100MG

Orbencarb
CAS#34622-58-7
Chemical Name: S-(2-Chlorobenzyl) diethylthiocarbamate

PS: Liquid


MF: C₁₂H₁₆ClNOS

MW: 257.8

Alternate Name: Orthobencarb

Part Number	Matrix	Size
N-12951-10MG	Neat	10MG

Prothiophos
CAS#34643-46-4
Chemical Name: O-(2,4-dichlorophenyl) O-ethyl S-propyl phosphorodithioate

PS: Liquid

MF: C₁₁H₁₅Cl₂O₂PS₂
MW: 345.2

Alternate Name: TOKUTHION™


Part Number	Matrix	Size
N-13577-100MG	Neat	100MG
S-13577A1-1ML	100ug/ml in Acetonitrile	1ML
S-13577U1-1ML	100ug/ml in Toluene	1ML

Butocarboxim CAS#34681-10-2
Chemical Name: 3-(methylthio)-2-butanone O-[(methylamino)carbonyl]oxime

PS: Liquid

MF: C₇H₁₄N₂O₂S

MW: 190.3

MP: 37°C



Part Number	Matrix	Size
N-12963-100MG	Neat	100MG
S-12963M1-1ML	100ug/ml in Methanol	1ML

Butoxycarboxim CAS#34681-23-7
Chemical Name: (EZ)-3-mesylbutanone O-methylcarbamoyloxime

PS: Solid

MF: C₇H₁₄N₂O₄S

MW: 222.3


MP: 85°C


Part Number	Matrix	Size
N-13921-100MG	Neat	100MG
S-13921M1-1ML	100ug/mL in Methanol	1ML

Tebutam CAS#35256-85-0
Chemical Name: 2,2-Dimethyl-N-(1-methylethyl)-N-(phenylmethyl)propanamide

PS: Liquid

MF: C₁₅H₂₃NO


MW: 233.39

Alternate Name: Butam(TM)

Part Number	Matrix	Size
N-13504-250MG	Neat	250MG
S-13504A1-1ML	100ug/ml in Acetonitrile	1ML
S-13504U1-1ML	100ug/ml in Toluene	1ML

Diflubenzuron CAS#35367-38-5
Chemical Name: 1-(4-Chlorophenyl)-3-(2,6-difluorobenzoyl)urea

PS: Solid

MF: C₁₄H₉ClF₂N₂O₂
MW: 310.69


MP: 228°C

Alternate Name: 1-(4-Chlorophenyl)-3-(2,6-difluorobenzoyl)urea

Part Number	Matrix	Size
N-11722-250MG	Neat	250MG
S-11722A1-1ML	100ug/ml in Acetonitrile	1ML

Sulprofos CAS#35400-43-2
Chemical Name: O-Ethyl-O-[4-(methylthio)phenyl]-S-propylphosphorothioate

PS: Liquid

MF: C₁₂H₁₉O₂PS₃
MW: 322.45

MP: -15°C



Part Number	Matrix	Size
N-13258-100MG	Neat	100MG
S-13258A1-1ML	100ug/ml in Acetonitrile	1ML
S-13258J1-1ML	100ug/ml in Hexane	1ML
S-13258J1-5ML	100ug/ml in Hexane	5ML

Imazalil CAS#35554-44-0
Chemical Name: 1-[2-(2,4-Dichlorophenyl)-2-(2-propenyloxy)ethyl]-1H-imidazole

PS: Solid

MF: C₁₄H₁₄Cl₂N₂O

MW: 297.2


MP: 52.7°C


Part Number	Matrix	Size
N-12198-100MG	Neat	100MG
S-12198A1-1ML	100ug/ml in Acetonitrile	1ML
S-12198J1-1ML	100ug/ml in Hexane	1ML

Iprodione CAS#36734-19-7
Chemical Name: 3-(3,5-dichlorophenyl)-N-isopropyl-2,4-dioximidazolidine-1-carboxamide

PS: Solid

MF: C₁₃H₁₃Cl₂N₃O₃
MW: 330.17

MP: 136°C



Part Number	Matrix	Size
N-12220-100MG	Neat	100MG
S-12220A1-1ML	100ug/ml in Acetonitrile	1ML

Dichlormid CAS#37764-25-3
Chemical Name: 2,2-Dichloro-N,N-di-2-propenyl acetamide

PS: Liquid

MF: C₈H₁₁Cl₂NO

MW: 208.1

MP: 5°C

Alternate Name: N,N-Diallyl-2,2-dichloroacetamide

Part Number	Matrix	Size
N-11662-100MG	Neat	100MG
S-11662A1-1ML	100ug/ml in Acetonitrile	1ML
S-11662U1-1ML	100ug/ml in Toluene	1ML


Etrifos CAS#38260-54-7

Chemical Name: O-6-ethoxy-2-ethylpyrimidin-4-yl O,O-dimethyl phosphorothioate

PS: Liquid

MF: C₉H₁₇N₂O₄PS

MW: 292.29

MP: -3.35°C



Part Number	Matrix	Size
N-12954-10MG	Neat	10MG

Glyphosate-isopropylammonium CAS#38641-94-0

PS: Solid

MF: C₆H₁₇N₂O₅P

MW: 228.22

MP: 143°C

Alternate Name: Glycel(TM)


Part Number	Matrix	Size
N-12134-250MG	Neat	250MG
S-12134F1-1ML	100ug/ml in Water	1ML
S-12134U1-1ML	100ug/ml in Toluene	1ML

Diethatyl ethyl CAS#38727-55-8

Chemical Name: N-Chloroacetyl-N-(2,6-diethylphenyl)-glycine ethyl ester

PS: Solid

MF: C₁₆H₂₂ClNO₃
MW: 311.84

MP: 50°C



Part Number	Matrix	Size
N-11690-100MG	Neat	100MG
S-11690M1-1ML	100ug/ml in Methanol	1ML

Fosetyl-Aluminum(Technical) CAS#39148-24-8

Chemical Name: Ethyl hydrogen phosphate aluminum salt

PS: Solid

MF: C₆H₁₈AlO₉P₃
MW: 354.11

MP: 215°C

Alternate Name: Fosetyl-Al

Part Number	Matrix	Size
N-12019-100MG	Neat	100MG


Thiofanox CAS#39196-18-4

Chemical Name: 1-(2,2-dimethyl-1-methylthiomethylpropylideneamino-oxy)-N-methylformamide

PS: Solid

MF: C₉H₁₈N₂O₂S

MW: 218.32


MP: 56.5°C

Alternate Name: THIOFANOCARB

Part Number	Matrix	Size
N-14058-25MG	Neat	25MG
S-14058A1-1ML	100ug/ml in Acetonitrile	1ML

Tebuthiuron TP 108 CAS#39222-73-6

PS: Solid

MF: C₆H₁₁N₃S


MW: 157.24


Part Number	Matrix	Size
N-13948-100MG	Neat	100MG

Dinocap CAS#39300-45-3

Chemical Name: 2-Methylheptyl-4,6-dinitrophenyl crotonate


PS: Liquid

MF: C₁₈H₂₄N₂O₆
MW: 364.39


Part Number	Matrix	Size
N-11783-100MG	Neat	100MG
S-11783M1-1ML	100ug/ml in Methanol	1ML
S-11783M1-5ML	100ug/ml in Methanol	5ML

d-trans-Cyphenothrin CAS#39515-40-7


PS: Liquid

MF: C₂₄H₂₅NO₃
MW: 375.5

Alternate Name: Gokilaht(TM)

Part Number	Matrix	Size
N-11561-100MG	Neat	100MG
S-11561A1-1ML	100ug/ml in Acetonitrile	1ML
S-11561T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Fenpropathrin
CAS#39515-41-8
Chemical Name: (RS)-a-cyano-3-phenoxybenzyl 2,2,3,3-tetramethylcyclopropane carboxylate

PS: Solid


MF: C₂₂H₂₃NO₃
MW: 349.4

Alternate Name: Danitol

Part Number	Matrix	Size
N-11960-250MG	Neat	250MG
S-11960A1-1ML	100ug/ml in Acetonitrile	1ML
S-11960T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Oxadiargyl
CAS#39807-15-3
Chemical Name: 5-tert-butyl-3-[2,4-dichloro-5-(prop-2-ynoxy)phenyl]-1,3,4-oxadiazol-2(3H)-one

PS: Solid

MF: C₁₅H₁₄Cl₂N₂O₃
MW: 341.2


MP: 131°C


Part Number	Matrix	Size
N-14011-100MG	Neat	100MG
S-14011A1-1ML	100ug/mL in Acetonitrile	1ML

Pendimethalin
CAS#40487-42-1
Chemical Name: N-[1-Ethylpropyl]-3,4-dimethyl-2,6-dinitrobenzeneamine

PS: Solid


MF: C₁₃H₁₉N₃O₄
MW: 281.31

MP: 56°C

Alternate Name: Prowl

Part Number	Matrix	Size
N-13140-100MG	Neat	100MG
S-13140A1-1ML	100ug/ml in Acetonitrile	1ML
S-13140K4-1ML	1000ug/ml in Isooctane	1ML

Methoprene
CAS#40596-69-8
Chemical Name: Isopropyl-(E,E)-(RS)-11-methoxy-3,7,11-trimethyl-2,4-dodeca-2,4-dienoate

PS: Liquid

MF: C₁₉H₃₄O₃
MW: 310.48



Part Number	Matrix	Size
N-12400-100MG	Neat	100MG
S-12400T4-1ML	1000ug/ml in t-Butyl methyl ether	1ML
S-12400T4-5ML	1000ug/ml in t-Butyl methyl ether	5ML

Azocyclotin (Technical)
CAS#41083-11-8
Chemical Name: Tri(cyclohexyl)-1H-1,2,4-triazol-1-yltin

PS: Solid

MF: C₂₀H₃₅N₃Sn

MW: 436.2


MP: 210°C


Part Number	Matrix	Size
N-11111-250MG	Neat	250MG
S-11111A1-1ML	100ug/ml in Acetonitrile	1ML

Profenofos
CAS#41198-08-7
Chemical Name: O-(4-Bromo-2-chlorophenyl)-O-ethyl-S-propylphosphorothioate

PS: Liquid

MF: C₁₁H₁₅BrClO₃PS

MW: 373.6


Part Number	Matrix	Size
N-13097-250MG	Neat	250MG
S-13097A1-1ML	100ug/ml in Acetonitrile	1ML
S-13097J1-1ML	100ug/ml in Hexane	1ML


Metamitron
CAS#41394-05-2
Chemical Name: 4-Amino-4,5-dihydro-3-methyl-6-phenyl-1,2,4-triazin-5-one

PS: Solid

MF: C₁₀H₁₀N₄O

MW: 202.24

MP: 166.9

FP: 69°C


Part Number	Matrix	Size
N-12383-250MG	Neat	250MG
S-12383M1-1ML	100ug/ml in Methanol	1ML


Bupirimate
CAS#41483-43-6
Chemical Name: 5-Butyl-2-ethylamino-6-methylpyrimidin-4-yl

PS: Solid

MF: C₁₃H₂₄N₄O₃S

MW: 316.47

MP: 50°C


FP: °C

Alternate Name: Nimrod(TM)

Part Number	Matrix	Size
N-11349-100MG	Neat	100MG

Tricyclazole
CAS#41814-78-2
Chemical Name: 5-Methyl-1,2,4-triazolo-(3,4-b)-benzothiazole

PS: Solid

MF: C₉H₇N₃S


MW: 189.2

Alternate Name: BEAM™

Part Number	Matrix	Size
N-11132-50MG	Neat	50MG
S-11132A1-1ML	100ug/mL in Acetonitrile	1ML
S-11132T1-1ML	100ug/ml in Methyl tert-butyl ether	1ML
S-11132T1-5ML	100ug/ml in Methyl tert-butyl ether	5ML

Isazophos
CAS#42509-80-8
Chemical Name: Phosphorothioic acid, O-[5-chloro-1-(1-methylethyl)-1H-1,2,4-triazol-3-yl] O,O-diethyl ester

PS: Liquid

MF: C₉H₁₇ClN₃O₃PS


MW: 313.75

Alternate Name: Brace

Part Number	Matrix	Size
N-12222-100MG	Neat	100MG
S-12222U1-1ML	100ug/ml in Toluene	1ML

Bifenox
CAS#42576-02-3
Chemical Name: Methyl-5(2,4-dichlorophenoxy)-2-nitrobenzoate

PS: Solid

MF: C₁₄H₉Cl₂NO₅
MW: 342.13


MP: 85°C


Part Number	Matrix	Size
N-11202-250MG	Neat	250MG
S-11202M1-1ML	100ug/ml in Methanol	1ML

Kinoprene
CAS#42588-37-4
Chemical Name: Prop-2-ynyl(+/-)-(E,E)-3,7,11-trimethyldodeca-2,4-dienoate

PS: Liquid

MF: C₁₈H₂₈O₂
MW: 276.46

FP: 40.5°C



Part Number	Matrix	Size
N-12295-100MG	Neat	100MG
S-12295M1-1ML	100ug/ml in Methanol	1ML
S-12295T4-1ML	1000ug/ml in t-Butyl methyl ether	1ML
S-12295T4-5ML	1000ug/ml in t-Butyl methyl ether	5ML

Daimuron
CAS#42609-52-9
Chemical Name: N-(4-methylphenyl)-N'-(1-methyl-1-phenylethyl)urea

PS: Solid

MF: C₁₇H₂₀N₂O

MW: 268.4


MP: 203°C

Alternate Name: 1-(1-Methyl-1-phenylethyl)-3-p-tolylurea

Part Number	Matrix	Size
N-13955-100MG	Neat	100MG
S-13955A1-1ML	100ug/mL in Acetonitrile	1ML

Oxyfluorfen
CAS#42874-03-3
Chemical Name: 2-Chloro-1-(3-ethoxy-4-nitrophenoxy)-4-(trifluoromethyl)benzene

PS: Solid

MF: C₁₅H₁₁ClF₃NO₄
MW: 361.72


MP: 85°C


Part Number	Matrix	Size
N-12742-250MG	Neat	250MG
S-12742M1-1ML	100ug/ml in Methanol	1ML

Triadimefon
CAS#43121-43-3
Chemical Name: 1-(4-Chlorophenoxy)-3,3-dimethyl-1-(1H-1,2,4-triazol-1-yl)-2-butanone

PS: Solid

MF: C₁₄H₁₆ClN₃O₂
MW: 293.75

MP: 82.3°C



Part Number	Matrix	Size
N-13636-500MG	Neat	500MG
S-13636A1-1ML	100ug/ml in Acetonitrile	1ML
S-13636T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-13636T1-5ML	100ug/ml in t-Butyl methyl ether	5ML

Difenzoquat methylsulfate
CAS#43222-48-6
Chemical Name: 1,2-Dimethyl-3,5-diphenyl-1H-pyrazolium methyl sulfate

PS: Solid

MF: C₁₇H₁₇N₂CH₃O₄S

MW: 360.11


MP: 155°C

Alternate Name: Difenzoquat

Part Number	Matrix	Size
N-11721-250MG	Neat	250MG
S-11721F1-1ML	100ug/ml in H2O	1ML
S-11721U1-1ML	100ug/ml in Toluene	1ML

Vinclozolin CAS#50471-44-8
Chemical Name: 3-(3,5-dichlorophenyl)-5-ethenyl-2,4-oxazolidinedione

PS: Solid

MF: C₁₂H₉Cl₂NO₃
MW: 286.11


MP: 108°C


Part Number	Matrix	Size
N-13745-250MG	Neat	250MG
S-13745A1-1ML	100ug/ml in Acetonitrile	1ML
S-13745K1-1ML	100ug/ml in Isooctane	1ML

Isoprothiolane CAS#50512-35-1
Chemical Name: Bis(1-methylethyl)-1,3-dithiolan-2-ylidene propanedioate

PS: Solid

MF: C₁₂H₁₈O₄S₂
MW: 290.42


MP: 54.5°C


Part Number	Matrix	Size
N-12278-250MG	Neat	250MG
S-12278A1-1ML	100ug/ml in Acetonitrile	1ML
S-12278T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Dimethachlor CAS#50563-36-5
Chemical Name: 2-Chloro-N-(2-methoxyethyl)aceto-2',6'-xylylidide

PS: Solid


MF: C₁₃H₁₈ClNO₂
MW: 255.77

MP: 45.8°C

Alternate Name: Teridox(TM)

Part Number	Matrix	Size
N-11754-250MG	Neat	250MG

Acifluorfen CAS#50594-66-6
Chemical Name: 5-[2-chloro-4-(trifluoromethyl)phenoxy]-2-nitrobenzoic acid

PS: Solid


MF: C₁₄H₇ClF₃NO₃
MW: 361.67


Part Number	Matrix	Size
N-11027-250MG	Neat	250MG
S-11027A1-1ML	100ug/ml in Acetonitrile	1ML
S-11027B1-1ML	100ug/ml in Acetone	1ML
S-11027B1-5ML	100ug/ml in Acetone	5ML

Acifluorfen methyl ester CAS#50594-67-7
Chemical Name: Methyl 5-[2-chloro-4-(trifluoromethyl)phenoxy]-2-nitrobenzoate

PS: Solid

MF: C₁₅H₉ClF₃NO₃
MW: 375.69


MP: 52.8°C

Alternate Name: Blazer methyl ester

Part Number	Matrix	Size
N-11028-100MG	Neat	100MG
S-11028Y1-1ML	100ug/ml in Isooctane:Acetone (90:10)	1ML
S-11028Y1-5ML	100ug/ml in Isooctane:Acetone (90:10)	5ML

Busan 40 CAS#51026-28-9
Chemical Name: Potassium-N-hydroxymethyl-N-methyl dithiocarbamate

PS: Solid

MF: C₃H₇NOS₂K


MW: 175.32


Part Number	Matrix	Size
N-11351-100MG	Neat	100MG
S-11351A1-1ML	100ug/ml in Acetonitrile	1ML
S-11351F1-1ML	100ug/ml in Water	1ML

Metolachlor CAS#51218-45-2
Chemical Name: 2-chloro-N-(2-ethyl-6-methylphenyl)-N-(2-methoxy-1-methylethyl)acetamide

PS: Liquid

MF: C₁₅H₂₂ClNO₂
MW: 283.8


MP: -62.1°C


Part Number	Matrix	Size
N-12478-100MG	Neat	100MG
S-12478A1-1ML	100ug/ml in Acetonitrile	1ML
S-12478T1-1ML	100ug/mL in t-Butyl methyl ether	1ML
S-12478T1-5ML	100ug/mL in t-Butyl methyl ether	5ML
S-FC2241S-1.2ML	100ug/ml in n-Nonane (C13)	1.2ML

Pretilachlor CAS#51218-49-6
Chemical Name: 2-Chloro-2',6'-diethyl-N-(2-propoxyethyl)acetanilide

PS: Liquid

MF: C₁₇H₂₆ClNO₂
MW: 311.9


FP: 129°C

Alternate Name: Rift(TM)

Part Number	Matrix	Size
N-13089-250MG	Neat	250MG
S-13089A1-1ML	100ug/ml in Acetonitrile	1ML
S-13089U1-1ML	100ug/ml in Toluene	1ML

1-cis-3-Chloroallyl-3,5,7-triaza-1-azonia-adamantane chlorid CAS#51229-78-8
Chemical Name: Hexamethylenetetramine chloroallyl chloride

PS: Solid

MF: C₉H₁₆Cl₂N₄
MW: 251.16


MP: 178°C

Alternate Name: Quaternium 15(TM)

Part Number	Matrix	Size
N-10045-100MG	Neat	100MG
S-10045M1-1ML	100ug/ml in Methanol	1ML

Velpar CAS#51235-04-2
Chemical Name: 3-Cyclohexyl-6-(dimethylamino)-1-methyl-1,3,5-triazine-2,4(1H,3H)-dione

PS: Solid

MF: C₁₂H₂₀N₄O₂
MW: 252.36


MP: 115°C

Alternate Name: Hexazinone

Part Number	Matrix	Size
N-13743-250MG	Neat	250MG
S-13743A1-1ML	100ug/ml in Acetonitrile	1ML
S-13743T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-13743T1-5ML	100ug/ml in t-Butyl methyl ether	5ML

Diclofop methyl CAS#51338-27-3
Chemical Name: Methyl-2-[4-(2,4'-dichlorophenoxyphenoxy)]propanoate

PS: Solid


MF: C₁₆H₁₄Cl₂O₄
MW: 341.19

MP: 39°C


Part Number	Matrix	Size
N-11677-250MG	Neat	250MG
S-11677A1-1ML	100ug/ml in Acetonitrile	1ML
S-11677T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Fenvalerate CAS#51630-58-1
Chemical Name: cyano(3-phenoxyphenyl)methyl 4-chloro-2-(1-methylethyl)benzeneacetate

PS: Solid

MF: C₂₅H₂₂ClNO₃
MW: 419.92

Alternate Name: Sanmarton


Part Number	Matrix	Size
N-13201-250MG	Neat	250MG
S-13201A4-1ML	1000ug/mL in Acetonitrile	1ML
S-13201A4-5ML	1000ug/mL in Acetonitrile	5ML
S-13201U1-1ML	100ug/ml in Toluene	1ML

Thidiazuron CAS#51707-55-2
Chemical Name: N-phenyl-N'-1,2,3-thiadiazol-5-ylurea

PS: Solid

MF: C₉H₈N₄OS


MW: 220.2

MP: 210.5°C

Alternate Name: Dropp(tm)

Part Number	Matrix	Size
N-13564-250MG	Neat	250MG
S-13564A1-1ML	100ug/ml in Acetonitrile	1ML
S-13564U1-1ML	100ug/ml in Toluene	1ML

Dimehypo CAS#52207-48-4
Chemical Name: Disodium S,S'-[2-(dimethylamino)trimethylene]di(thiosulfate)


PS: Solid

MF: C₅H₁₁NO₆S₄Na

Alternate Name: Thiosultap-disodium

Part Number	Matrix	Size
N-11752-250MG	Neat	250MG
S-11752F1-1ML	100ug/ml in H2O	1ML

zeta-Cypermethrin CAS#52315-07-8
Chemical Name: (+)-alpha-Cyano-3-phenoxybenzyl 2,2-dimethyl-3-(2,2-dichlorovinyl)cyclopropane carboxylate


PS: Liquid

MF: C₂₂H₁₉Cl₂NO₃
MW: 416.3

Alternate Name: Fury(TM)

Part Number	Matrix	Size
N-13754-250MG	Neat	250MG
S-13754A1-1ML	100ug/ml in Acetonitrile	1ML
S-13754U1-1ML	100ug/ml in Toluene	1ML

Cypermethrin CAS#52315-07-8
Chemical Name: (RS)-alpha-cyano-3-phenoxybenzyl (1RS,3RS;1RS,3SR)-3-(2,2-dichlorovinyl)-2,2-dimethylcyclopropanecarboxylate

PS: Liquid


MF: C₂₂H₁₉Cl₂NO₃
MW: 416.3


Part Number	Matrix	Size
N-11545-100MG	Neat	100MG
S-11545A1-1ML	100ug/ml in Acetonitrile	1ML
S-11545U1-1ML	100ug/ml in Toluene	1ML

Permethrin CAS#52645-53-1
Chemical Name: (3-phenoxyphenyl)methyl 3-(2,2-dichloroethenyl)-2,2-dimethylcyclopropanecarboxylate

PS: Solid

MF: C₂₁H₂₀Cl₂O₃
MW: 391.29


MP: 34°C


Part Number	Matrix	Size
N-12848-250MG	Neat	250MG
S-12848A1-1ML	100ug/ml in Acetonitrile	1ML
S-12848T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Flamprop isopropyl CAS#52756-22-6
Chemical Name: Isopropyl-N-benzoyl-N-(3-chloro-4-fluorophenyl)-DL-alaninate

PS: Solid

MF: C₁₉H₁₉ClFNO₃
MW: 363.84


MP: 56°C


Part Number	Matrix	Size
N-11972-100MG	Neat	100MG
S-11972A1-1ML	100ug/ml in Acetonitrile	1ML
S-11972T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Flamprop-methyl CAS#52756-25-9
Chemical Name: Methyl-N-benzoyl-N-(3-chloro-4-fluorophenyl)-DL-alaninate

PS: Solid

MF: C₁₇H₁₅ClFNO₃
MW: 335.78


MP: 85°C


Part Number	Matrix	Size
N-11973-100MG	Neat	100MG
S-11973A1-1ML	100ug/ml in Acetonitrile	1ML
S-11973T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Prosulfocarb CAS#52888-80-9
Chemical Name: Benzyl dipropylthiocarbamate

PS: Liquid

MF: C₁₄H₂₁NOS


MW: 251.4


Part Number	Matrix	Size
N-13918-10MG	Neat	10MG

Deltamethrin CAS#52918-63-5
Chemical Name: (S)-a-Cyano-3-phenoxybenzyl (1R,3R)-3-(2,2-dibromovinyl)-2,2-dimethylcyclopropane carboxylate

PS: Solid

MF: C₂₂H₁₉Br₂NO₃
MW: 505.21


MP: 101°C


Part Number	Matrix	Size
N-11579-250MG	Neat	250MG
S-11579A1-1ML	100ug/ml in Acetonitrile	1ML
S-11579T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Pyrimethanil CAS#53112-28-0
Chemical Name: N-(4,6-Dimethylpyrimidin-2-yl)aniline

PS: Solid


MF: C₁₂H₁₃N₃
MW: 199.3

MP: 96.3°C


Part Number	Matrix	Size
N-13159-100MG	Neat	100MG
S-13159A1-1ML	100ug/ml in Acetonitrile	1ML
S-13159U1-1ML	100ug/ml in Toluene	1ML

Dichlorprop-butoxyethyl ester CAS#53404-31-2
Chemical Name: Butoxyethyl 2-(2,4-dichlorophenoxy)propionate

PS: Liquid


MF: C₁₅H₂₀Cl₂O₄


Part Number	Matrix	Size
N-11673-50MG	Neat	50MG

Endrin ketone CAS#53494-70-5
Chemical Name: Hexachloro-decahydro-metheno-3H-cyclopenta(a)pental-ene-3-one

PS: Solid

MF: C₁₂H₈Cl₆O


MW: 380.9


Part Number	Matrix	Size
N-11856-10MG	Neat	10MG
S-11856A1-1ML	100ug/ml in Acetonitrile	1ML
S-11856U1-1ML	100ug/ml in Toluene	1ML

Pyrimuron
CAS#53558-25-1
Chemical Name: 1-(4-nitrophenyl)-3-(3-pyridylmethyl)urea

PS: Solid


MF: C₁₃H₁₂N₄O₃
MW: 272.27

MP: 223°C

Alternate Name: Pyriminil

Part Number	Matrix	Size
N-13738-100MG	Neat	100MG
S-13738A1-1ML	100ug/ml in Acetonitrile	1ML
S-13738U1-1ML	100ug/ml in Toluene	1ML

Mefluidide
CAS#53780-34-0
Chemical Name: 5'-(1,1,1-Trifluoromethanesulfonamido)acet-2',4'-xylylide

PS: Solid


MF: C₁₁H₁₃F₃N₂O₃S


Part Number	Matrix	Size
N-12367-100MG	Neat	100MG
S-12367M1-1ML	100ug/ml in Methanol	1ML

Empenthrin
CAS#54406-48-3
Chemical Name: (2E)-1-ethynyl-2-methyl-2-penten-1-yl 2,2-dimethyl-3-(2-methyl-1-propen-1-yl)cyclopropanecarboxylate

PS: Liquid

MF: C₁₈H₂₆O₂
MW: 274.4


FP: 107°C

Alternate Name: 1-Ethynyl-2-methyl-2-penten-1-yl chrysanthemate

Part Number	Matrix	Size
N-13913-10MG	Neat	10MG

Cycloprate
CAS#54460-46-7
Chemical Name: Hexadecylcyclopropane carboxylate

PS: Liquid

MF: C₂₀H₃₈O₂
MW: 310.58


MP: 19.8°C

Alternate Name: Zardex(TM)

Part Number	Matrix	Size
N-11541-500MG	Neat	500MG
S-11541A1-1ML	100ug/ml in Acetonitrile	1ML
S-11541T4-1ML	1000ug/ml in t-Butyl methyl ether	1ML
S-11541T4-5ML	1000ug/ml in t-Butyl methyl ether	5ML

Chlorethoxyfos
CAS#54593-83-8
Chemical Name: O,O-Diethyl O-(1,2,2,2-tetrachloroethyl) phosphorothioate

PS: Liquid

MF: C₆H₁₁Cl₄O₃PS


MW: 336

Alternate Name: Fortress(TM)

Part Number	Matrix	Size
N-11427-250MG	Neat	250MG
S-11427A1-1ML	100ug/ml in Acetonitrile	1ML
S-11427U1-1ML	100ug/ml in Toluene	1ML

Albendazole
CAS#54965-21-8
Chemical Name: methyl N-[5-(propylthio)-1H-benzimidazol-2-yl]carbamate

PS: Liquid

MF: C₁₂H₁₅N₃O₂S


MP: 208°C


Part Number	Matrix	Size
N-12959-10MG	Neat	10MG

Bitertanol
CAS#55179-31-2
PS: Solid

MF: C₂₀H₂₃N₃O₂
MW: 337.46

MP: 138.6°C


FP: 100°C

Alternate Name: Baycor

Part Number	Matrix	Size
N-11127-100MG	Neat	100MG
S-11127A1-1ML	100ug/ml in Acetonitrile	1ML
S-11127U1-1ML	100ug/ml in Toluene	1ML

Triadimenol
CAS#55219-65-3
PS: Solid

MF: C₁₄H₁₈ClN₃O₂
MW: 295.8

MP: 118°C


FP: 100°C

Alternate Name: Baytan(TM)

Part Number	Matrix	Size
N-11129-250MG	Neat	250MG
S-11129A1-1ML	100ug/ml in Acetonitrile	1ML
S-11129U1-1ML	100ug/ml in Toluene	1ML

Ethalfuralin
CAS#55283-68-6
Chemical Name: N-Ethyl-N-(2-methyl-2-propenyl)-2,6-dinitro-4-(trifluoromethyl)benzenamine

PS: Solid

MF: C₁₃H₁₄F₃N₃O₄
MW: 425.28


MP: 55°C


Part Number	Matrix	Size
N-11865-250MG	Neat	250MG
S-11865A1-1ML	100ug/ml in Acetonitrile	1ML
S-11865J4-1ML	1000ug/ml in Hexane	1ML
S-11865J4-5ML	1000ug/ml in Hexane	5ML

Carbosulfan
CAS#55285-14-8
Chemical Name: 2,3-Dihydro-2,2-dimethylbenzofuran-7-yl(dibutylaminothio)methylcarbamate

PS: Liquid

MF: C₂₀H₃₂N₂O₃S

MW: 380.6



Part Number	Matrix	Size
N-11409-250MG	Neat	250MG
S-11409A1-1ML	100ug/ml in Acetonitrile	1ML

Dimethipin
CAS#55290-64-7
Chemical Name: 2,3-Dihydro-5,6-dimethyl-1,4-dithiin-1,1,4,4-tetraoxide

PS: Solid

MF: C₆H₁₀O₄S₂
MW: 210.3

MP: 167°C


FP: 246°C


Part Number	Matrix	Size
N-11757-100MG	Neat	100MG
S-11757M1-1ML	100ug/ml in Methanol	1ML

Triclopyr
CAS#55335-06-3
Chemical Name: 3,5,6-Trichloro-2-pyridinyloxyacetic acid

PS: Solid

MF: C₇H₄Cl₃NO₃
MW: 256.5

MP: 150.5°C

Alternate Name: Crossbow(TM)

Part Number	Matrix	Size
N-13659-500MG	Neat	500MG
S-13659A1-1ML	100ug/ml in Acetonitrile	1ML
S-13659T1-1ML	100ug/ml in t-Butyl methyl ether	1ML


Pyridate
CAS#55512-33-9
Chemical Name: 6-Chloro-3-phenylpyridazin-4-yl-S-octylthiocarbonate

PS: Liquid

MF: C₁₉H₂₃ClN₂O₂S

MW: 378.95

MP: 27°C


FP: 131°C


Part Number	Matrix	Size
N-13155-500MG	Neat	500MG
S-13155A1-1ML	100ug/ml in Acetonitrile	1ML
S-13155T1-1ML	100ug/ml t-Butyl methyl ether	1ML

Mepronil
CAS#55814-41-0
Chemical Name: 2-Methyl-N-(3-(1-methylethoxy)phenyl)benzamide

PS: Solid


MF: C₁₇H₁₉NO₂
MW: 269.37

MP: 92°C

Alternate Name: Basitac(TM)

Part Number	Matrix	Size
N-12372-100MG	Neat	100MG
S-12372M1-1ML	100ug/mL in Methanol	1ML

Isouron Solution
CAS#55861-78-4
Chemical Name: 3-(5-tert-butylisoxazol-3-yl)-1,1-dimethylurea

PS: Liquid


MF: C₁₀H₁₇N₃O₂
MW: 211.3


Part Number	Matrix	Size
S-13971A1-1ML	100ug/ml in Acetonitrile	1ML

Difenacoum
CAS#56073-07-5
Chemical Name: 3-(3-[1,1'-biphenyl]-4-yl)-1,2,3,4-tetrahydro-1-naphthalenyl)-4-hydroxy-2H-1-benzopyran-2-one

PS: Solid

MF: C₃₁H₂₄O₃
MW: 444.55


MP: 228°C

Alternate Name: Frunax-DS(TM)

Part Number	Matrix	Size
N-11719-100MG	Neat	100MG
S-11719M1-1ML	100ug/ml in Methanol	1ML

Brodifacoum
CAS#56073-10-0
Chemical Name: 3-[3-(4'-Bromobiphenyl-4-yl)-1,2,3,4-tetrahydro-1-naphthyl]-4-hydroxycoumarin

PS: Solid

MF: C₃₁H₂₃BrO₃
MW: 523.44


MP: 201°C


Part Number	Matrix	Size
N-11329-50MG	Neat	50MG
S-11329A1-1ML	100ug/ml in Acetonitrile	1ML

Flurprimidol
CAS#56425-91-3
Chemical Name: a-(1-Methylethyl)-a-[4-trifluoromethoxyphenyl]-5-pyrimidine-methanol

PS: Solid

MF: C₁₅H₁₅F₃N₂O₂
MW: 312.3

MP: 95°C

Alternate Name: Cutless*

Part Number	Matrix	Size
N-12002-100MG	Neat	100MG
S-12002M1-1ML	100ug/ml in Methanol	1ML


Tolclofos-Methyl
CAS#57018-04-9
Chemical Name: O-2,6-Dichloro-p-tolyl O,O-dimethyl phosphorothioate

PS: Solid

MF: C₉H₁₁Cl₂O₃PS

MW: 301.13


MP: 79°C

FP: 210°C

Alternate Name: Risolex(TM)

Part Number	Matrix	Size
N-13578-500MG	Neat	500MG
S-13578A1-1ML	100ug/ml in Acetonitrile	1ML
S-13578U1-1ML	100ug/ml in Toluene	1ML


Dichloroprop methyl ester
CAS#57153-17-0
Chemical Name: 2-(2,4-Dichlorophenoxy)propanoic acid methyl ester

PS: Liquid

MF: C₁₀H₁₀Cl₂O₃
MW: 249.09

Alternate Name: 2,4-DP methyl ester

Part Number	Matrix	Size
N-11672-100MG	Neat	100MG
S-11672M1-1ML	100ug/ml in Methanol	1ML
S-11672M1-5ML	100ug/ml in Methanol	5ML


MCPB methyl ester
CAS#57153-18-1
Chemical Name: methyl 4-(4-chloro-2-methylphenoxy)butanoate

PS: Liquid


Part Number	Matrix	Size
N-12358-250MG	Neat	250MG
S-12358A1-1ML	100ug/ml in Acetonitrile	1ML
S-12358T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

2,4,5-T-triethylamine salt
CAS#57213-69-1
Chemical Name: 2,4,5-Trichlorophenoxyacetic acid triethylamine salt


PS: Solid

MF: C₇H₄Cl₃NO₃·C₆H₁₅N


Part Number	Matrix	Size
N-10652-250MG	Neat	250MG
S-10652M1-1ML	100ug/ml in Methanol	1ML

Pyroquilon
CAS#57369-32-1
Chemical Name: 1,2,5,6-Tetrahydro-4H-pyrrolo[3,2,1-ij]quinolin-4-one


PS: Liquid

MF: C₁₁H₁₁NO


Part Number	Matrix	Size
N-12998-10MG	Neat	10MG

Furalaxyl
CAS#57646-30-7
Chemical Name: Alanine, N-(2,6-dimethylphenyl)-N-(2-furanylcarbonyl)-, methyl ester

PS: Solid


MF: C₁₇H₁₉NO₄
MW: 301.34

Alternate Name: Fongarid

Part Number	Matrix	Size
N-12978-10MG	Neat	10MG

Metalaxyl CAS#57837-19-1
Chemical Name: Methyl N-(methoxyacetyl)-N-(2,6-xylyl)-DL-alaninate

PS: Solid

MF: C₁₅H₂₁NO₄
MW: 279.33


MP: 71°C

Alternate Name: Apron

Part Number	Matrix	Size
N-12380-100MG	Neat	100MG
S-12380A1-1ML	100ug/ml in Acetonitrile	1ML
S-12380U1-1ML	100ug/ml in Toluene	1ML

Acequinocyl CAS#57960-19-7
Chemical Name: 3-Dodecyl-1,4-dihydro-1,4-dioxo-2-naphthyl acetate

PS: Solid

MF: C₂₄H₃₂O₄
MW: 384.5


MP: 59.6°C

Alternate Name: Kanemite(TM)

Part Number	Matrix	Size
N-11003-250MG	Neat	250MG
S-11003A1-1ML	100ug/ml in Acetonitrile	1ML
S-11003U1-1ML	100ug/ml in Toluene	1ML

Cymoxanil CAS#57966-95-7
Chemical Name: 2-Cyano-N-[(ethylamino)carbonyl]-2-methoxyiminoacetamide

PS: Solid


MF: C₇H₁₀N₄O₃
MW: 198.21

MP: 160°C


Part Number	Matrix	Size
N-11544-100MG	Neat	100MG
S-11544A1-1ML	100ug/ml in Acetonitrile	1ML
S-11544U1-1ML	100ug/ml in Toluene	1ML

Flamprop-M-isopropyl CAS#57973-67-8
Chemical Name: isopropyl N-benzoyl-N-(3-chloro-4-fluorophenyl)-D-alaninate

PS: Solid


MF: C₁₉H₁₉ClFNO₃
MW: 363.8


Part Number	Matrix	Size
N-13939-250MG	Neat	250MG
S-13939M1-1ML	100ug/mL in Methanol	1ML

Pyrazolynate CAS#58011-68-0
Chemical Name: (2,4-Dichlorophenyl)[1,3-dimethyl-5-[[[(4-methylphenyl)sulfonyl]oxy]-1H-pyrazol-4-yl]methanone


PS: Solid

MF: C₁₉H₁₆Cl₂N₂O₄S

MW: 439.3


Part Number	Matrix	Size
N-13083-50MG	Neat	50MG
S-13083A1-1ML	100ug/mL in Acetonitrile	1ML

Tridiphane Solution CAS#58138-08-2
Chemical Name: (RS)-2-(3,5-dichlorophenyl)-2-(2,2,2-trichloroethyl)oxirane

PS: Liquid



Part Number	Matrix	Size
S-13965A1-1ML	100ug/ml in Acetonitrile	1ML

Imazalil sulfate CAS#58594-72-2
Chemical Name: (+)-1-[2-(2,4-Dichlorophenyl)-2-(2-propenyloxy)ethyl]-1H-imidazole-3-sulfate

PS: Solid

MF: C₁₄H₁₆Cl₂N₂O₅S

MW: 395.3

MP: 128°C



Part Number	Matrix	Size
N-12199-100MG	Neat	100MG
S-12199F1-1ML	100ug/ml in Water	1ML

Kadethrine CAS#58769-20-3
Chemical Name: Cyclopropanecarboxylic acid, 3-[(E)-(dihydro-2-oxo-3(2H)-thienylidene)methyl]-2,2-dimethyl-, [5-(phenylmethyl)-3-furanyl]methyl ester, (1R,3S)-

PS: Liquid


MF: C₂₃H₂₄O₄S

MW: 396.5

MP: 31°C


Part Number	Matrix	Size
N-12947-10MG	Neat	10MG

Ofurace Solution
CAS#58810-48-3
PS: Liquid


MF: C₁₄H₁₆ClNO₃
MW: 281.7


Part Number	Matrix	Size
S-12988K10-1ML	10ug/ml in Isooctane	1ML

Thiodicarb
CAS#59669-26-0
Chemical Name: Dimethyl N,N'-[thiobis[(methylimino_carbonyloxy)]bis[ethanimidothioate]

PS: Solid

MF: C₁₀H₁₈N₄O₄S₃
MW: 354.48


MP: 173°C


Part Number	Matrix	Size
N-13568-250MG	Neat	250MG
S-13568A1-1ML	100ug/ml in Acetonitrile	1ML

Fluridone
CAS#59756-60-4
Chemical Name: 1-methyl-3-phenyl-5-[3-(trifluoromethyl)phenyl]-4(1H)-pyridinone

PS: Solid

MF: C₁₉H₁₄F₃NO

MW: 329.3

Alternate Name: SONAR™


Part Number	Matrix	Size
N-13217-100MG	Neat	100MG
S-13217T1-1ML	100ug/ml in Methyl tert-butyl ether	1ML
S-13217T1-5ML	100ug/ml in Methyl tert-butyl ether	5ML

Fenarimol
CAS#60168-88-9
Chemical Name: a(2-Chlorophenyl)-a-(4-chlorophenyl)-5-(pyrimidine methanol)

PS: Solid


MF: C₁₇H₁₂Cl₂N₂O

MW: 331.21

MP: 118°C


Part Number	Matrix	Size
N-11948-100MG	Neat	100MG
S-11948A1-1ML	100ug/mL in Acetonitrile	1ML
S-11948T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-11948T1-5ML	100ug/ml in t-Butyl methyl ether	5ML


Azaconazole solution
CAS#60207-31-0
PS: Liquid

MF: C₁₂H₁₁Cl₂N₃O₂
MW: 300.1


Part Number	Matrix	Size
S-13976A1-1ML	100ug/ml in Acetonitrile	1ML

Propiconazole
CAS#60207-90-1
Chemical Name: 1-[[[2-(2,4-dichlorophenyl)-4-propyl-1,3-dioxolan-2-yl]methyl]-1H-1,2,4-triazole


PS: Liquid

MF: C₁₅H₁₇Cl₂N₃O₂
MW: 342.2

Alternate Name: Tilt(TM)

Part Number	Matrix	Size
N-13576-250MG	Neat	250MG
S-13576A1-1ML	100ug/mL in Acetonitrile	1ML
S-13576E1-1ML	100ug/ml in Cyclohexane	1ML

Etaconazole
CAS#60207-93-4
Chemical Name: 1-[[[2-(2,4-dichlorophenyl)-4-ethyl-1,3-dioxolan-2-yl]methyl]-1H-1,2,4-triazole


PS: Solid

MF: C₁₄H₁₅Cl₂N₃O₂
MW: 328.19


Part Number	Matrix	Size
N-12985-100MG	Neat	100MG

Chlorthiophos
CAS#60238-56-4
Chemical Name: O-[dichloro(methylthio)phenyl] O,O-diethyl phosphorothioate

PS: Solid

MF: C₁₁H₁₅Cl₂O₃PS₂
MW: 361.2



Part Number	Matrix	Size
N-12969-100MG	Neat	100MG

Ethyl-3-amino-2,5-dichlorobenzoate CAS#60541-86-8

Chemical Name: 3-Amino-2,5-dichlorobenzoic acid - ethyl ester

PS: Solid

MF: C₉H₉Cl₂NO₂
MW: 234.08

MP: 60°C

Alternate Name: Chloramben ethyl ester


Part Number	Matrix	Size
N-11916-100MG	Neat	100MG
S-11916A1-1ML	100ug/ml in Acetonitrile	1ML
S-11916T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Furmecyclox CAS#60568-05-0

Chemical Name: 3-Furancarboxamide, N-cyclohexyl-N-methoxy-2,5-dimethyl-

PS: Liquid

MF: C₁₄H₂₁NO₃
MW: 251.3

MP: 33°C

Alternate Name: Furmetamid


Part Number	Matrix	Size
N-12979-10MG	Neat	10MG
S-12979M1-1ML	100ug/ml in Methanol	1ML

Triclopyr methyl ester CAS#60825-26-5

Chemical Name: 3,4,6-Trichloro-2-pyridinyloxyacetic acid methyl ester

PS: Solid

MF: C₈H₆Cl₃NO₃
MW: 270.5


MP: 105°C


Part Number	Matrix	Size
N-13661-100MG	Neat	100MG
S-13661A1-1ML	100ug/ml in Acetonitrile	1ML
S-13661T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Flurochloridone CAS#61213-25-0

Chemical Name: 3-chloro-4-(chloromethyl)-1-[3-(trifluoromethyl)phenyl]-2-pyrrolidinone

PS: Solid

MF: C₁₂H₁₀Cl₂F₃NO



Part Number	Matrix	Size
N-12987-100MG	Neat	100MG
S-12987M1-1ML	100ug/ml in Methanol	1ML

Dimepiperate Solution CAS#61432-55-1

Chemical Name: S-1-methyl-1-phenylethyl piperidine-1-carbothioate

PS: Liquid

MF: C₁₅H₂₁NOS


MW: 263.4


Part Number	Matrix	Size
S-13978A1-1ML	100ug/ml in Acetonitrile	1ML

Bentazon methyl derivative CAS#61592-45-8

PS: Solid


MF: C₁₁H₁₄N₂O₃S

MW: 254.3


Part Number	Matrix	Size
N-11143-50MG	Neat	50MG
S-11143Y1-1ML	100ug/ml in Isooctane:Acetone (90:10)	1ML
S-11143Y1-5ML	100ug/ml in Isooctane:Acetone (90:10)	5ML

Cymiazole solution CAS#61676-87-7


PS: Liquid

MF: C₁₂H₁₄N₂S


Part Number	Matrix	Size
S-14051A1-1ML	100ug/ml in Acetonitrile	1ML

cis-Permethrin CAS#61949-76-6


PS: Solid

MF: C₂₁H₂₀Cl₂O₃
MW: 391.29


Part Number	Matrix	Size
N-11483-50MG	Neat	50MG
S-11483A1-1ML	100ug/ml in Acetonitrile	1ML
S-11483T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-FC2216AS-1.2ML	50ug/mL in n-Nonane (C13)	1.2ML

trans-Permethrin
CAS#61949-77-7
Chemical Name: trans-meta-Phenoxybenzyl-3-(2,2-Dichlorovinyl)-2,2-dimethylcyclopropanecarboxylate

PS: Solid

MF: C₂₁H₂₀Cl₂O₃
MW: 391.29



Part Number	Matrix	Size
N-13620-50MG	Neat	50MG
S-13620A1-1ML	100ug/ml in Acetonitrile	1ML
S-13620T1-1ML	100ug/ml in t-Butyl methyl ether	1ML
S-FC2216BS-1.2ML	50ug/mL in n-Nonane (C13)	1.2ML

Endothal monohydrate
CAS#62059-43-2
Chemical Name: 7-Oxabicyclo[2,2,1]heptane-2,3-dicarboxylic acid

PS: Solid

MF: C₈H₁₀O₅H₂O


MW: 204.2

MP: 144°C

Alternate Name: Endothal (unstated stereochemistry)

Part Number	Matrix	Size
N-11852-100MG	Neat	100MG
S-11852A1-1ML	100ug/ml in Acetonitrile	1ML
S-11852A1-5ML	100ug/ml in Acetonitrile	5ML
S-11852F0-1ML	50ug/mL in Water	1ML
S-11852F10-1ML	10 ug/mL in Water	1ML
S-11852F10-5ML	10 ug/mL in Water	5ML

N-[(2,4-Dichlorophenoxy)acetyl]-DL-methionine
CAS#62625-13-2
PS: Solid

MF: C₁₃H₁₅Cl₂NO₄S


MW: 352.23


Part Number	Matrix	Size
N-12505-10MG	Neat	10MG
S-12505A1-1ML	100ug/ml in Acetonitrile	1ML
S-12505T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Fenothiocarb Solution
CAS#62850-32-2
Chemical Name: S-(4-Phenoxybutyl)-N,N-dimethylthiocarbamate

PS: Liquid

MF: C₁₃H₁₉NO₂S


MW: 253.4


Part Number	Matrix	Size
S-13962A1-1ML	100ug/ml in Acetonitrile	1ML

Flumetralin
CAS#62924-70-3
Chemical Name: N-(2-Chloro-6-fluorobenzyl)-N-ethyl-a,a,a-trifluoro-2,6-dinitro-p-toluidine

PS: Solid

MF: C₁₆H₁₂C₁F₄N₃O₄
MW: 421.74


MP: 102°C


Part Number	Matrix	Size
N-11988-250MG	Neat	250MG
S-11988A1-1ML	100ug/ml in Acetonitrile	1ML
S-11988B1-1ML	100ug/ml in Acetone	1ML
S-11988B1-5ML	100ug/ml in Acetone	5ML
S-11988K1-1ML	100ug/ml in Isooctane	1ML

Nuarimol
CAS#63284-71-9
Chemical Name: (±)-2-chloro-4'-fluoro-a-(pyrimidin-5-yl)benzhydryl alcohol


PS: Solid

MF: C₁₇H₁₂ClFN₂O

MW: 314.7


Part Number	Matrix	Size
N-13942-100MG	Neat	100MG
S-13942M1-1ML	100ug/mL in Methanol	1ML


Diofenolan Solution
CAS#63837-33-2
PS: Liquid

MF: C₁₈H₂₀O₄
MW: 300.35


Part Number	Matrix	Size
S-13979A1-1ML	100ug/ml in Acetonitrile	1ML

Anilofos
CAS#64249-01-0
Chemical Name: S-[2-[(4-chlorophenyl)(1-methylethyl)amino]-2-oxoethyl] O,O-dimethyl phosphorodithioate

PS: Solid

MF: C₁₃H₁₉ClNO₃PS₂
MW: 367.8

Alternate Name: ANILOGUARD


Part Number	Matrix	Size
N-13131-100MG	Neat	100MG
S-13131K1-1ML	100ug/mL in Isooctane	1ML

Ethiozin CAS#64529-56-2

Chemical Name: 4-amino-6-(1,1-dimethylethyl)-3-(ethylthio)-1,2,4-triazin-5(4H)-one

PS: Solid

MF: C₉H₁₆N₄OS

MW: 228.3

Alternate Name: Tycor


Part Number	Matrix	Size
N-13731-100MG	Neat	100MG
S-13731M1-1ML	100ug/ml in Methanol	1ML

Triflururon CAS#64628-44-0

Chemical Name: 2-Chloro-N-(((4-(trifluoromethoxy)phenyl)amino)carbonyl)benzamide

PS: Solid

MF: C₁₅H₁₀ClF₃N₂O₃
MW: 358.72


MP: 195°C

Alternate Name: Alsystine(TM)

Part Number	Matrix	Size
N-13687-100MG	Neat	100MG
S-13687A1-1ML	100ug/ml in Acetonitrile	1ML
S-13687U1-1ML	100ug/ml in Toluene	1ML

Triclopyr 2-butoxyethyl ester CAS#64700-56-7

Chemical Name: 2-Butoxyethyl [(3,5,6-trichloropyridin-2-yl)oxy]acetate

PS: Liquid

MF: C₁₃H₁₆Cl₃NO₄
MW: 356.65

Alternate Name: Triclopyr-butotyl

Part Number	Matrix	Size
N-13660-250MG	Neat	250MG


Chlorsulfuron CAS#64902-72-3

Chemical Name: 1-(2-chlorophenylsulfonyl)-3-(4-methoxy-6-methyl-1,3,5-triazin-2-yl)urea

PS: Solid

MF: C₁₂H₁₂ClN₅O₄S

MW: 357.78

MP: 174°C

Alternate Name: Glean

Part Number	Matrix	Size
N-11461-100MG	Neat	100MG


beta-Cypermethrin

CAS#65731-84-2

Chemical Name: cyano(3-phenoxyphenyl)methyl 3-(2,2-dichloroethenyl)-2,2-dimethylcyclopropanecarboxylate

PS: Solid

MF: C₂₂H₁₉Cl₂NO₃
MW: 416.3

MP: 80.5°C

Alternate Name: Agrothrin(TM)


Part Number	Matrix	Size
N-11192-250MG	Neat	250MG
S-11192A1-1ML	100ug/ml in Acetonitrile	1ML
S-11192T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

(S)-Methoprene

CAS#65733-16-6

Chemical Name: Isopropyl-(E,E)-(S)-11-methoxy-3,7,11-trimethyldodeca-2,4-dienoate

PS: Liquid

MF: C₁₉H₃₄O₃
MW: 310.5



Part Number	Matrix	Size
N-13832-100MG	Neat	100MG

S-Hydroprene

CAS#65733-18-8

Chemical Name: ethyl (2E,4E)-3,7,11-trimethyl-2,4-dodecadienoate

PS: Liquid

MF: C₁₇H₃₀O₂
MW: 266.4

FP: 84°C


(7S)-Hydroprene

Part Number	Matrix	Size
N-13189-50MG	Neat	50MG
S-13189A1-1ML	100ug/ml in Acetonitrile	1ML
S-13189J1-1ML	100ug/ml in Hexane	1ML

Alloxydim-Sodium

CAS#66003-55-2

PS: Solid

MF: C₁₇H₂₄NNaO₅
MW: 345.4



Part Number	Matrix	Size
N-11052-50MG	Neat	50MG

Pencycuron
CAS#66063-05-6
Chemical Name: N-[(4-chlorophenyl)methyl]-N-cyclopentyl-N'-phenylurea

PS: Solid

MF: C₁₉H₂₁ClN₂O

MW: 328.87


MP: 125°C

Alternate Name: Monceren

Part Number	Matrix	Size
N-12488-100MG	Neat	100MG
S-12488A1-1ML	100ug/ml in Acetonitrile	1ML
S-12488U1-1ML	100ug/ml in Toluene	1ML

Cyromazine
CAS#66215-27-8
Chemical Name: N-Cyclopropyl-1,3,5-triazine-2,4,6-triamine

PS: Solid

MF: C₆H₁₀N₆
MW: 166.18


MP: 224.9°C


Part Number	Matrix	Size
N-11550-250MG	Neat	250MG
S-11550A1-1ML	100ug/ml in Acetonitrile	1ML
S-11550B1-1ML	100ug/ml in Acetone	1ML

Esfenvalerate
CAS#66230-04-4
Chemical Name: (S)-cyano(3-phenoxyphenyl)methyl (S)-4-chloro-2-(1-methylethyl)benzeneacetate

PS: Solid

MF: C₂₅H₂₂ClNO₃
MW: 419.9


MP: 59°C

Alternate Name: Asana(TM)

Part Number	Matrix	Size
N-11102-100MG	Neat	100MG
S-11102A1-1ML	100ug/ml in Acetonitrile	1ML
S-11102T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Penconazole
CAS#66246-88-6
Chemical Name: 1-(2-(2,4-Dichlorophenyl)pentyl)-1H-1,2,4-triazole

PS: Solid


MF: C₁₃H₁₅Cl₂N₃
MW: 284.21

MP: 61°C

Alternate Name: Ofir(TM)

Part Number	Matrix	Size
N-12822-100MG	Neat	100MG
S-12822A1-1ML	100ug/ml in Acetonitrile	1ML
S-12822U1-1ML	100ug/ml in Toluene	1ML

Flutolanil
CAS#66332-96-5
Chemical Name: a,a,a-Trifluoro-3-isopropoxy-o-toluanilide

PS: Solid


MF: C₁₇H₁₆F₃NO₂
MW: 323.3

Alternate Name: Moncut (TM)

Part Number	Matrix	Size
N-12004-250MG	Neat	250MG
S-12004M1-1ML	100ug/ml in Methanol	1ML

Fenoxaprop ethyl
CAS#66441-23-4
Chemical Name: ethyl 2-[4-[(6-chloro-2-benzoxazolyl)oxy]phenoxy]propanoate

PS: Solid

MF: C₁₈H₁₆ClNO₅
MW: 291.38


MP: 84°C

Alternate Name: Fenoxaprop ethyl (racemate CAS # 82110-72-3)

Part Number	Matrix	Size
N-11957-100MG	Neat	100MG
S-11957M1-1ML	100ug/ml in Methanol	1ML

Tralomethrin
CAS#66841-25-6
Chemical Name: (S)-cyano(3-phenoxyphenyl)methyl (1R,3S)-2,2-dimethyl-3-(1,2,2,2-tetrabromoethyl)cyclopropanecarboxylate

PS: Liquid


MF: C₂₂H₁₉Br₄NO₃
MW: 665

MP: 138°C


Part Number	Matrix	Size
N-13588-100MG	Neat	100MG
S-13588A1-1ML	100ug/ml in Acetonitrile	1ML
S-13588T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Methasulfocarb
CAS#66952-49-6
Chemical Name: Carbamothioic acid, N-methyl-, S-[4-[(methylsulfonyl)oxy]phenyl] ester

PS: Solid


MF: C₉H₁₁NO₄S₂
MW: 261.3


Part Number	Matrix	Size
N-13283-25MG	Neat	25MG
S-13283A1-1ML	100ug/mL in Acetonitrile	1ML

Metazachlor
CAS#67129-08-2
Chemical Name: Acetamide, 2-chloro-N-(2,6-dimethylphenyl)-N-(1H-pyrazol-1-yl-methyl)-

PS: Solid

MF: C₁₄H₁₆ClN₃O


MW: 277.8

Alternate Name: Butisan S(TM)

Part Number	Matrix	Size
N-12385-100MG	Neat	100MG

alpha-Cypermethrin
CAS#67375-30-8
Chemical Name: ?-cyano-3-phenoxybenzyl (1S)-cis-3-(2,2-dichlorovinyl)-2,2-dimethylcyclopropanecarboxylate

PS: Solid

MF: C₂₂H₁₉Cl₂NO₃
MW: 416.3


MP: 78°C


Part Number	Matrix	Size
N-11061-100MG	Neat	100MG
S-11061A1-1ML	100ug/ml in Acetonitrile	1ML
S-11061J1-1ML	100ug/ml in Hexane	1ML

Amdro
CAS#67485-29-4
Chemical Name: tetrahydro-5,5-dimethyl-2(1H)-pyrimidinone [3-[4-(trifluoromethyl)phenyl]-1-[2-[4-(trifluoromethyl)phenyl]ethenyl]-2-propen-1-ylidene]hydrazone

PS: Solid


MF: C₂₅H₂₄F₆N₄
MW: 494.53

MP: 174°C

Alternate Name: Hydramethylnon

Part Number	Matrix	Size
N-11063-100MG	Neat	100MG
S-11063A1-1ML	100ug/ml in Acetonitrile	1ML
S-11063J1-1ML	100ug/ml in Hexane	1ML

Prochloraz
CAS#67747-09-5
Chemical Name: 1-(N-Propyl-N-(2-(2,4,6-trichlorophenoxy)ethyl)carbamoyl)imidazole

PS: Solid


MF: C₁₅H₁₆Cl₃N₃O₂
MW: 376.69

Alternate Name: Ascurit(TM)

Part Number	Matrix	Size
N-13093-250MG	Neat	250MG
S-13093A1-1ML	100ug/ml in Acetonitrile	1ML
S-13093U1-1ML	100ug/ml in Toluene	1ML

Forchlorfenuron
CAS#68157-60-8
Chemical Name: 1-(2-Chloro-4-pyridyl)-3-phenylurea

PS: Solid

MF: C₁₂H₁₀ClN₃O


MW: 247.7


Part Number	Matrix	Size
N-12010-100MG	Neat	100MG
S-12010A1-1ML	100ug/ml in Acetonitrile	1ML

beta-Cyfluthrin
CAS#68359-37-5
Chemical Name: Cyclopropanecarboxylic acid, 3-(2,2-dichloroethenyl)-2,2-dimethyl-, cyano(4-fluoro-3-phenoxyphenyl)methyl ester

PS: Solid


MF: C₂₂H₁₈Cl₂FNO₃
MW: 434.29

MP: 60°C

Alternate Name: beta-Baythroid(TM)

Part Number	Matrix	Size
N-11191-250MG	Neat	250MG
S-11130A4-1ML	1000ug/ml in Acetonitrile	1ML
S-11130A4-5ML	1000ug/ml in Acetonitrile	5ML
S-11130U1-1ML	100ug/mL in Toluene	1ML
S-11191A1-1ML	100ug/ml in Acetonitrile	1ML
S-11191T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Cyfluthrin
CAS#68359-37-5
PS: Solid

MF: C₂₂H₁₈Cl₂FNO₃
MW: 434.29

MP: 60°C

Alternate Name: Cyfoxylate


Part Number	Matrix	Size
N-11130-250MG	Neat	250MG
S-11130A4-1ML	1000ug/ml in Acetonitrile	1ML
S-11130A4-5ML	1000ug/ml in Acetonitrile	5ML
S-11130U1-1ML	100ug/mL in Toluene	1ML

Triflumizole
CAS#68694-11-1
Chemical Name: 4-Chloro-a,a,a-trifluoro-N-(1-imidazol-1-yl-2-propoxyethylidene)-o-toluidine

PS: Solid

MF: C₁₅H₁₅ClF₃N₃O

MW: 345.78

MP: 63.5°C



Part Number	Matrix	Size
N-13686-250MG	Neat	250MG

Buprofezin
CAS#69327-76-0
Chemical Name: 2-tert-Butylimino-3-isopropyl-5-phenyl-1,3,5-thiadiazinan-4-one

PS: Solid

MF: C₁₆H₂₃N₃OS


MW: 305.48

MP: 106.1°C


Part Number	Matrix	Size
N-11350-100MG	Neat	100MG
S-11350A1-1ML	100ug/ml in Acetonitrile	1ML
S-11350U1-1ML	100ug/ml in Toluene	1ML

Fluazifop-methyl
CAS#69335-90-6
Chemical Name: Propanoic acid, 2-[4-[[5-(trifluoromethyl)-2-pyridinyl]oxy]phenoxy]-, methyl ester


PS: Liquid

MF: C₁₆H₁₄F₃NO₄
MW: 341.3


Part Number	Matrix	Size
N-13906-50MG	Neat	50MG

Fluazifop-P-methyl
CAS#69335-90-6
Chemical Name: Propanoic acid, 2-[4-[[5-(trifluoromethyl)-2-pyridinyl]oxy]phenoxy]-, methyl ester


PS: Liquid

MF: C₁₆H₁₄F₃NO₄
MW: 341.3


Part Number	Matrix	Size
N-13908-50MG	Neat	50MG

Fluazifop
CAS#69335-91-7
Chemical Name: 2-[4-[[5-(Trifluoromethyl)-2-pyridinyl]oxy]phenoxy]propanoic acid

PS: Solid


MF: C₁₅H₁₂F₃NO₄
MW: 327.3


Part Number	Matrix	Size
N-11977-10MG	Neat	10MG
S-11977A1-1ML	100ug/ml in Acetonitrile	1ML
S-11977A1-5ML	100ug/ml in Acetonitrile	5ML

Fluroxypyr
CAS#69377-81-7
Chemical Name: 4-Amino-3,5-dichloro-6-fluoro-2-pyridyloxyacetic acid

PS: Solid


MF: C₇H₅Cl₂FN₂O₃
MW: 255.03

MP: 232°C

Alternate Name: Advance(TM)

Part Number	Matrix	Size
N-12000-100MG	Neat	100MG
S-12000M1-1ML	100ug/ml in Methanol	1ML

Flumethrin
CAS#69770-45-2
Chemical Name: cyano(4-fluoro-3-phenoxyphenyl)methyl 3-[2-chloro-2-(4-chlorophenyl)ethenyl]-2,2-dimethylcyclopropanecarboxylate


PS: Liquid

MF: C₂₈H₂₂Cl₂FNO₃
MW: 510.4


Part Number	Matrix	Size
N-13139-100MG	Neat	100MG
S-13139A1-1ML	100ug/mL in Acetonitrile	1ML

Haloxypop (free acid)
CAS#69806-34-4
Chemical Name: (RS)-2-[4-(3-chloro-5-trifluoromethyl-2-pyridyloxy)phenoxy]propionic acid

PS: Solid


MF: C₁₅H₁₁ClF₃NO₄
MW: 361.7


Part Number	Matrix	Size
N-12143-10MG	Neat	10MG

Fluazifop-butyl
CAS#69806-50-4
Chemical Name: Butyl (RS)-2-[4-(5-trifluoromethyl-2-pyridyloxy)phenoxy]propionate

PS: Liquid


MF: C₁₉H₂₀F₃NO₄
MW: 383.4

MP: 13°C

Alternate Name: Fusilade(TM)

Part Number	Matrix	Size
N-11978-250MG	Neat	250MG
S-11978A1-1ML	100ug/ml in Acetonitrile	1ML
S-11978U1-1ML	100ug/ml in Toluene	1ML


Flucythrinate
CAS#70124-77-5
Chemical Name: (+)-Cyano(3-phenoxyphenyl)methyl(+)-4-(difluoromethoxy)-alpha-(1-methylethyl) benzeneacetate

PS: Liquid

MF: C₂₆H₂₃F₂NO₄
MW: 451.48

Alternate Name: Aastar(TM)

Part Number	Matrix	Size
N-11984-100MG	Neat	100MG
S-11984A1-1ML	100ug/ml in Acetonitrile	1ML
S-11984U1-1ML	100ug/ml in Toluene	1ML

Ivermectin
CAS#70288-86-7
Chemical Name: Cardomec(TM)/Eqvalan(TM)/Hyvermectin(TM)/Ivomec(TM)/Mectizan(TM)/Zimecterin(TM)/22,23-Dihydroavermectin 1/22,23-Dihydroabamectin


PS: Solid

Alternate Name: Cardomec(TM)

Part Number	Matrix	Size
N-12289-100MG	Neat	100MG
S-12289A1-1ML	100ug/ml in Acetonitrile	1ML
S-12289U1-1ML	100ug/ml in Toluene	1ML

Mefenoxam
CAS#70630-17-0
Chemical Name: Methyl N-(methoxyacetyl)-N-(2,6-xylyl)-D-alaninate

PS: Liquid


MF: C₁₅H₂₁NO₄
MW: 279.33

MP: -38.7°C

Alternate Name: Metalaxyl-M

Part Number	Matrix	Size
N-12365-100MG	Neat	100MG
S-12365A1-1ML	100ug/ml in Acetonitrile	1ML
S-12365T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Fenoxaprop-P-ethyl
CAS#71283-80-2
Chemical Name: Ethyl(R)-2-[4-[(6-chloro-2-benzoxazolyl)oxy]phenoxy]propanoate


PS: Solid

MF: C₁₈H₁₆ClNO₅
MW: 361.8

Alternate Name: Acclaim(TM)

Part Number	Matrix	Size
N-11958-100MG	Neat	100MG
S-11958A1-1ML	100ug/ml in Acetonitrile	1ML
S-11958U1-1ML	100ug/ml in Toluene	1ML

Chlorfluazuron
CAS#71422-67-8
Chemical Name: 1-[3,5-Dichloro-4-(3-chloro-5-trifluoromethyl-2-pyridyloxy)phenyl]-3-(2,6-difluoro-benzoyl)urea

PS: Solid


MF: C₂₀H₉Cl₃F₃N₃O₃
MW: 540.7

Alternate Name: Atabron(TM)

Part Number	Matrix	Size
N-11430-50MG	Neat	50MG
S-11430A1-1ML	100ug/ml in Acetonitrile	1ML

Pyrazoxyfen
CAS#71561-11-0
Chemical Name: 2-[4-(2,4-Dichlorobenzoyl)-1,3-dimethylpyrazol-5-yloxy]acetophenone

PS: Solid


MF: C₂₀H₁₆Cl₂N₂O₃
MW: 403.28

MP: 112°C


Part Number	Matrix	Size
N-13148-100MG	Neat	100MG
S-13148M1-1ML	100ug/ml in Methanol	1ML


Benalaxyl
CAS#71626-11-4
Chemical Name: Methyl N-phenylacetyl-N-2,6-xylyl-DL-alaninate

PS: Solid

MF: C₂₀H₂₃NO₃
MW: 325.4

Alternate Name: Galben(TM)

Part Number	Matrix	Size
N-11133-250MG	Neat	250MG
S-11133M1-1ML	100ug/ml in Methanol	1ML

Abamectin
CAS#71751-41-2
PS: Solid

MF: C₄₈H₇₂O₁₄, C₄₇H₇₀O₁₄
MP: 150°C

Alternate Name: Avermectin B1


Part Number	Matrix	Size
N-10995-100MG	Neat	100MG
S-10995A1-1ML	100ug/ml in Acetonitrile	1ML

Fomesafen
CAS#72178-02-0
Chemical Name: 5-(2-Chloro-a,a,a-trifluoro-p-tolyloxy)-N-methylsulfonyl-2-nitro-benzamide

PS: Solid

MF: C₁₅H₁₀ClF₃N₂O₆S

MW: 438.8

MP: 220°C



Part Number	Matrix	Size
N-12008-500MG	Neat	500MG
S-12008M1-1ML	100ug/ml in Methanol	1ML

Fenoxycarb
CAS#72490-01-8
Chemical Name: Ethyl (2-(4-phenoxyphenoxy)ethyl)carbamate

PS: Solid

MF: C₁₇H₁₉NO₄
MW: 301.3

MP: 53°C


FP: 224°C

Alternate Name: Comply(TM)

Part Number	Matrix	Size
N-11959-100MG	Neat	100MG
S-11959A1-1ML	100ug/ml in Acetonitrile	1ML

Haloxypop-P methyl
CAS#72619-32-0
Chemical Name: Methyl (R)-2-(4-(3-chloro-5-trifluoromethyl-2-pyridyloxy)phenoxy)propionate

PS: Liquid


MF: C₁₆H₁₃ClF₃NO₄
MW: 375.7

MP: -12.4°C

Alternate Name: Haloxypop-R methyl

Part Number	Matrix	Size
N-12146-50MG	Neat	50MG

Imiprothrin
CAS#72963-72-5
Chemical Name: [2,5-dioxo-3-(2-propynyl)-1-imidazolidinyl]methyl 2,2-dimethyl-3-(2-methyl-1-propenyl)cyclopropanecarboxylate

PS: Liquid

MF: C₁₇H₂₂N₂O₄
MW: 318.4



Part Number	Matrix	Size
N-14006-50MG	Neat	50MG
S-14006A1-1ML	100ug/mL in Acetonitrile	1ML

Mefenacet
CAS#73250-68-7
Chemical Name: 2-(1,3-benzothiazol-2-yloxy)-N-methylacetanilide

PS: Solid

MF: C₁₆H₁₄N₂O₂S

MW: 298.4


MP: 134.8°C

Alternate Name: Hinochloa

Part Number	Matrix	Size
N-12364-250MG	Neat	250MG
S-12364M1-1ML	100ug/ml in Methanol	1ML

Sethoxydim
CAS#74051-80-2
Chemical Name: 2-(1-Ethoxyimino)butyl)-5(2-ethylthio)propyl)-3-hydroxy-2-cyclohexene-1-one

PS: Liquid


MF: C₁₇H₂₉NO₃S

MW: 327.5


Part Number	Matrix	Size
N-13210-100MG	Neat	100MG
S-13210U1-1ML	100ug/ml in Toluene	1ML

Acclonifen
CAS#74070-46-5
Chemical Name: 2-chloro-6-nitro-3-phenoxybenzenamine


PS: Solid

MF: C₁₂H₉ClN₂O₃
MW: 264.7


Part Number	Matrix	Size
N-13991-250MG	Neat	250MG
S-13991A1-1ML	100ug/mL in Acetonitrile	1ML

Clofentezine
CAS#74115-24-5
Chemical Name: 3,6-Bis(2-chlorophenyl)-1,2,4,5-tetrazine

PS: Solid

MF: C₁₄H₈Cl₂N₄
MW: 303.16

Alternate Name: Apollo


Part Number	Matrix	Size
N-11490-100MG	Neat	100MG
S-11490M1-1ML	100ug/ml in Methanol	1ML

Sulfometuron methyl CAS#74222-97-2
Chemical Name: 2-[3-(4,6-Dimethylpyrimidin-2-yl)ureidosulfonyl]benzoic acid

PS: Solid

MF: C₁₅H₁₆N₄O₅S

MW: 364.4

MP: 202°C



Part Number	Matrix	Size
N-13254-100MG	Neat	100MG
S-13254U1-1ML	100ug/ml in Toluene	1ML

Metsulfuron methyl CAS#74223-64-6
Chemical Name: Methyl-2-[[[4-methoxy-6-methyl-1,3,5-triazin-2-yl]amino]carbonyl]amino]sulfonyl]benzoate

PS: Solid


MF: C₁₄H₁₅N₅O₆S

MW: 381.36

MP: 163°C


Part Number	Matrix	Size
N-12482-100MG	Neat	100MG
S-12482A1-1ML	100ug/ml in Acetonitrile	1ML

Tecloftalam CAS#76280-91-6
Chemical Name: 2',3',4',5,6-hexachlorophthalanilic acid

MF: C₁₄H₅Cl₆NO₃
MW: 447.9

Alternate Name: 3,4,5,6-Tetrachloro-N-(2,3-dichlorophenyl)phthalamic acid


Part Number	Matrix	Size
N-14068-10MG	Neat	10MG

Quizalofop ethyl CAS#76578-14-8
Chemical Name: Ethyl-2-[4-[(6-chloro-2-quinoxalinyloxy]-phenoxy]propionate

PS: Solid

MF: C₁₉H₁₇ClN₂O₄
MW: 372.8

MP: 91.7°C

FP: 100°C

Alternate Name: Assure


Part Number	Matrix	Size
N-13174-100MG	Neat	100MG
S-13174A1-1ML	100ug/ml in Acetonitrile	1ML
S-13174U1-1ML	100ug/ml in Toluene	1ML

Flutriafol CAS#76674-21-0
Chemical Name: alpha-(2-Fluorophenyl)-alpha-(4-fluorophenyl)-1H-1,2,4-triazole-1-ethanol

PS: Solid

MF: C₁₆H₁₃F₂N₃O

MW: 301.3


MP: 130°C

Alternate Name: Armour(TM)

Part Number	Matrix	Size
N-12005-100MG	Neat	100MG
S-12005A1-1ML	100ug/ml in Acetonitrile	1ML
S-12005U1-1ML	100ug/ml in Toluene	1ML

gamma-Cyhalothrin CAS#76703-62-3
Chemical Name: (S)-a-cyano-3-phenoxybenzyl (Z)-(1R,3R)-3-(2-chloro-3,3,3-trifluoro-1-propenyl)-2,2-dimethylcyclopropanecarboxylate

PS: Solid

MF: C₂₃H₁₉ClF₃NO₃
MW: 449.9

MP: 55.6°C



Part Number	Matrix	Size
N-12106-100MG	Neat	100MG

Paclobutrazol CAS#76738-62-0
Chemical Name: (2RS,3RS)-1-(4-Chlorophenyl)-4,4-dimethyl-2-(1H-1,2,4-triazol-1-yl)-pentan-3-ol

PS: Solid

MF: C₁₅H₂₀ClN₃O

MW: 293.83

MP: 165°C



Part Number	Matrix	Size
N-12811-500MG	Neat	500MG
S-12811A1-1ML	100ug/ml in Acetonitrile	1ML
S-12811U1-1ML	100ug/ml in Toluene	1ML

Glufosinate-ammonium CAS#77182-82-2
Chemical Name: Ammonium-4-[hydroxy(methyl)phosphinoyl]-DL-homoalaninate

PS: Solid


MF: C₅H₁₅N₂O₄P

MW: 198.2

MP: 215°C

Alternate Name: D,L-Phosphinothricin

Part Number	Matrix	Size
N-12111-100MG	Neat	100MG

(E)-Ferimzone
CAS#77359-18-3
PS: Liquid


MF: C₁₅H₁₈N₄
MW: 254.3


Part Number	Matrix	Size
S-14031A1-1ML	100ug/ml in Acetonitrile	1ML

Lactofen
CAS#77501-63-4
Chemical Name: Ethyl O-[5-(2-chloro-a,a,a-trifluoro-p-tolyloxy)-2-nitrobenzoyl]-DL-lactate

PS: Solid

MF: C₁₉H₁₅ClF₃NO₇
MW: 461.8


MP: 44°C


Part Number	Matrix	Size
N-12304-100MG	Neat	100MG
S-12304A1-1ML	100ug/ml in Acetonitrile	1ML
S-12304U1-1ML	100ug/ml in Toluene	1ML

Fluoroglycofen-ethyl
CAS#77501-90-7
Chemical Name: Benzoic acid, 5-[2-chloro-4-(trifluoromethyl)phenoxy]-2-nitro-, 2-ethoxy-2-oxoethyl ester

PS: Solid

MF: C₁₇H₂₆ClF₃N₂O₇
MW: 311.89

MP: 65°C

Alternate Name: Compete(TM)


Part Number	Matrix	Size
N-11999-100MG	Neat	100MG

Oxadixyl
CAS#77732-09-3
Chemical Name: 2-Methoxy-N-(2-oxo-1,3-oxazolidin-3-yl)aceto-2'-xyli-dide/N-(2,6-Dimethylphenyl)-2-methoxy-N-(2-oxo-3-oxazolidinyl)acetamide

PS: Solid

MF: C₁₄H₁₈N₂O₄
MW: 278.3

MP: 105°C

FP: 100°C

Alternate Name: Anchor (TM)


Part Number	Matrix	Size
N-12732-100MG	Neat	100MG
S-12732A1-1ML	100ug/ml in Acetonitrile	1ML
S-12732U1-1ML	100ug/ml in Toluene	1ML

Hexythiazox
CAS#78587-05-0
Chemical Name: trans-5-(4-Chlorophenyl)-N-cyclohexyl-4-methyl-2-oxo-3-thiaz/idine carboxamide

PS: Solid

MF: C₁₇H₂₁ClN₂O₂S

MW: 352.9


MP: 108°C


Part Number	Matrix	Size
N-12184-100MG	Neat	100MG
S-12184A1-1ML	100ug/ml in Acetonitrile	1ML

Fluazifop-P-butyl
CAS#79241-46-6
Chemical Name: (R)-2-[4-[[5-(Trifluoromethyl)-2-pyridinyl]oxy]phenoxy]propanoic acid

PS: Liquid


MF: C₁₉H₂₀F₃NO₄
MW: 383.4

MP: -20°C


Part Number	Matrix	Size
N-11979-100MG	Neat	100MG
S-11979A1-1ML	100ug/ml in Acetonitrile	1ML
S-11979U1-1ML	100ug/ml in Toluene	1ML

Dichlorprop 2-ethylhexyl ester
CAS#79270-78-3
Chemical Name: 2-ethylhexyl 2-(2,4-dichlorophenoxy)propanoate

PS: Liquid


MF: C₁₇H₂₄Cl₂O₃
MW: 347.28


Part Number	Matrix	Size
N-13821-250MG	Neat	250MG

Thifensulfuron-methyl
CAS#79277-27-3
Chemical Name: methyl 3-(4-methoxy-6-methyl-1,3,5-triazin-2-ylcarbamoylsulfamoyl)thiophene-2-carboxylate

PS: Solid

MF: C₁₂H₁₃N₅O₆S₂
MW: 387.42

MP: 186°C

Alternate Name: Harmony(TM)

Part Number	Matrix	Size
N-13056-100MG	Neat	100MG
S-13056A1-1ML	100ug/ml in Acetonitrile	1ML

Fluazinam CAS#79622-59-6


Chemical Name: 3-Chloro-N-(3-chloro-5-trifluoromethyl-2-pyridyl)-alpha,alpha,alpha-trifluoro-2,6-dinitro-p-toluidine

PS: Solid

MF: C₁₃H₄Cl₂F₆N₄O₄

MW: 465.11

MP: 115°C



Alternate Name: Altima(TM)

Part Number	Matrix	Size
N-11980-100MG	Neat	100MG
S-11980A1-1ML	100ug/ml in Acetonitrile	1ML
S-11980U1-1ML	100ug/ml in Toluene	1ML


Hexaconazole CAS#79983-71-4

Chemical Name: (RS)-2-(2,4-Dichlorophenyl)-1-(1H-1,2,4-triazol-1-yl)hexan-2-ol

PS: Solid

MF: C₁₄H₁₇Cl₂N₃O

MW: 314.2



Alternate Name: Anvil(TM)

Part Number	Matrix	Size
N-12165-100MG	Neat	100MG
S-12165A1-1ML	100ug/ml in Acetonitrile	1ML
S-12165U1-1ML	100ug/ml in Toluene	1ML

Diafenthiuron CAS#80060-09-9


Chemical Name: N-(2,6-Bis(1-methylethyl)-4-phenoxyphenyl)-N'-(1,1-dimethylethyl)thiourea

PS: Solid

MF: C₂₃H₃₂N₂O₂

MW: 384.63

FP: 149°C



Alternate Name: Pegasus(TM)

Part Number	Matrix	Size
N-11612-250MG	Neat	250MG
S-11612A1-1ML	100ug/ml in Acetonitrile	1ML
S-11612U1-1ML	100ug/ml in Toluene	1ML

Etofenprox CAS#80844-07-1


Chemical Name: 2-(4-Ethoxyphenyl)-2-methylpropyl 3-phenoxybenzyl ether

PS: Solid

MF: C₂₅H₂₈O₃

MW: 376.5

MP: 37.4°C



Alternate Name: Primo(TM)

Part Number	Matrix	Size
N-11941-100MG	Neat	100MG
S-11941A1-1ML	100ug/ml in Acetonitrile	1ML
S-11941U1-1ML	100ug/ml in Toluene	1ML


Imazapyr CAS#81334-34-1

Chemical Name: 2-[4,5-Dihydro-4-methyl-4-(1-methylethyl)-5-oxo-1H-imidazol-2-yl]-3-pyridinecarboxylic acid

PS: Solid

MF: C₁₃H₁₅N₃O₃

MW: 261.28



Alternate Name: Arsenal(TM)

Part Number	Matrix	Size
N-12203-100MG	Neat	100MG
S-12203A1-1ML	100ug/ml in Acetonitrile	1ML

Imazaquin CAS#81335-37-7


Chemical Name: (RS)-2-(4-Isopropyl-4-methyl-5-oxo-2-imidazolin-2-yl)quinoline-3-carboxylic acid

PS: Solid

MF: C₁₇H₁₇N₃O₃

MW: 311.34

MP: 219°C



Part Number	Matrix	Size
N-12204-100MG	Neat	100MG
S-12204A1-1ML	100ug/ml in Acetonitrile	1ML
S-12204T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Imazethapyr CAS#81335-77-5


Chemical Name: (RS)-5-Ethyl-2-(4-isopropyl-4-methyl-5-oxo-2-imidazolin-2-yl)nicotinic acid

PS: Solid

MF: C₁₅H₁₉N₃O₃

MW: 289.37

MP: 172°C



Part Number	Matrix	Size
N-12205-100MG	Neat	100MG
S-12205M1-1ML	100ug/ml in Methanol	1ML

Imazamethabenz-methyl CAS#81405-85-8


Chemical Name: Methyl-2-(4,5-dihydro-4-methyl-4-(1-methylethyl)-5-oxo-1H-imidazol-2-yl)-4(or 5)-methyl ester benzoic acid

PS: Solid

MF: C₁₆H₂₀N₂O₃

MW: 288.3

MP: 113°C




Alternate Name: Assert(TM)

Part Number	Matrix	Size
N-12200-100MG	Neat	100MG
S-12200A1-1ML	100ug/ml in Acetonitrile	1ML
S-12200U1-1ML	100ug/ml in Toluene	1ML

Fluroxypyr-meptyl
CAS#81406-37-3
Chemical Name: 1-methylheptyl 2-[(4-amino-3,5-dichloro-6-fluoro-2-pyridinyl)oxy] acetate

PS: Solid

MF: C₁₅H₂₁Cl₂FN₂O₃
MW: 367.2

Alternate Name: Vista(TM)


Part Number	Matrix	Size
N-12001-100MG	Neat	100MG
S-12001M1-1ML	100ug/ml in Methanol	1ML

Tridemorph (Mix of Isomers)
CAS#81412-43-3
Chemical Name: Dimethyl-2,6 tridecyl-4 morpholine

PS: Liquid

MF: C₁₉H₃₉NO

MW: 297.59

FP: 142°C


n = 10, 11, 12 or 13

Alternate Name: Calixin(TM)


Part Number	Matrix	Size
N-13665-100MG	Neat	100MG
S-13665A1-1ML	100ug/ml in Acetonitrile	1ML
S-13665U1-1ML	100ug/ml in Toluene	1ML

Command
CAS#81777-89-1
Chemical Name: 2-(2-Chlorobenzyl)-4,4-dimethylisoxazolidin-3-one

PS: Liquid

MF: C₁₂H₁₄ClNO₂
MW: 239.72

MP: 25°C

FP: 157°C

Alternate Name: Clomazone


Part Number	Matrix	Size
N-11495-100MG	Neat	100MG
S-11495M1-1ML	100ug/ml in Methanol	1ML

Triasulfuron
CAS#82097-50-5
Chemical Name: 1-[2-(2-chloroethoxy)phenylsulfonyl]-3-(4-methoxy-6-methyl-1,3,5-triazin-2-yl)urea

PS: Solid

MF: C₁₄H₁₆C₁N₅O₅S

MW: 401.86


MP: 186°C


Part Number	Matrix	Size
N-13638-100MG	Neat	100MG
S-13638A1-1ML	100ug/ml in Acetonitrile	1ML

Isoxaben
CAS#82558-50-7
Chemical Name: N-[3-(1-Ethyl-1-methylpropyl)isoxazol-5-yl]-2,6-dimethoxybenzamide

PS: Solid

MF: C₁₈H₂₄N₂O₄
MW: 332.4


MP: 176°C


Part Number	Matrix	Size
N-12285-250MG	Neat	250MG
S-12285M1-1ML	100ug/ml in Methanol	1ML

Benfuracarb
CAS#82560-54-1
Chemical Name: Ethyl N-[2,3-dihydro-2,2-dimethylbenzofuran-7-yloxy]carbamoyl(methyl)aminothio]-N-isopropyl-beta-alaninate

PS: Liquid

MF: C₂₀H₃₀N₂O₅S


MW: 410.5

Alternate Name: Furacon(TM)

Part Number	Matrix	Size
N-11137-100MG	Neat	100MG

Bifenthrin
CAS#82657-04-3
Chemical Name: (2-methyl[1,1'-biphenyl]-3-yl)methyl (1R,3R)-rel-3-[(1Z)-2-chloro-3,3,3-trifluoro-1-propen-1-yl]-2,2-dimethylcyclopropanecarboxylate

PS: Solid

MF: C₂₃H₂₂ClF₃O₂
MW: 422.87


MP: 69°C

Alternate Name: Capture

Part Number	Matrix	Size
N-11203-100MG	Neat	100MG
S-11203A1-1ML	100ug/ml in Acetonitrile	1ML
S-11203T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Bensulfuron-methyl
CAS#83055-99-6
Chemical Name: Methyl 2-[[[[(4,6-dimethoxy-2-pyrimidinyl)amino]carbonyl]amino]sulfonyl]methyl]benzoate

PS: Solid


MF: C₁₆H₁₈N₄O₇S

MW: 410.44

Alternate Name: Escuri(TM)

Part Number	Matrix	Size
N-11139-100MG	Neat	100MG
S-11139A1-1ML	100ug/ml in Acetonitrile	1ML

Fluazifop-P
CAS#83066-88-0
Chemical Name: (2R)-2-[4-[[5-(Trifluoromethyl)-2-pyridinyl]oxy]phenoxy]propanoic acid

PS: Solid


MF: C₁₅H₁₂F₃NO₄
MW: 327.3


Part Number	Matrix	Size
N-13907-50MG	Neat	50MG

Teflubenzuron
CAS#83121-18-0
Chemical Name: N-[[[(3,5-Dichloro-2,4-difluorophenyl)amino]carbonyl]-2,6-difluorobenzamide

PS: Solid

MF: C₁₄H₆Cl₂F₄N₂O₂
MW: 381.12


MP: 220°C

Alternate Name: Nomolt(TM)

Part Number	Matrix	Size
N-13506-250MG	Neat	250MG
S-13506A1-1ML	100ug/ml in Acetonitrile	1ML

Diflufenican
CAS#83164-33-4
Chemical Name: N-(2,4-Difluorophenyl)-2-(3-(trifluoromethyl)phenoxy)-3-pyridine-carboxamide

PS: Solid

MF: C₁₉H₁₁F₅N₂O₂
MW: 394.32

MP: 160°C

Alternate Name: Legacy(TM)

Part Number	Matrix	Size
N-11723-100MG	Neat	100MG


Uniconazole-P
CAS#83657-17-4
Chemical Name: (E)-(S)-1-(4-chlorophenyl)-4,4-dimethyl-2-(1H-1,2,4-triazol-1-yl)pent-1-en-3-ol

PS: Solid

MF: C₁₅H₁₆ClN₃O

MW: 291.81

MP: 152.1°C


FP: 195°C


Part Number	Matrix	Size
N-13735-10MG	Neat	10MG
S-13735M1-1ML	100ug/ml in Methanol	1ML

Uniconazole
CAS#83657-22-1
Chemical Name: (E)-(RS)-1-(4-chlorophenyl)-4,4-dimethyl-2-(1H-1,2,4-triazol-1-yl)pent-1-en-3-ol

PS: Solid

MF: C₁₅H₁₆ClN₃O

MW: 291.78



Part Number	Matrix	Size
N-14032-50MG	Neat	50MG
S-14032A1-1ML	100ug/ml in Acetonitrile	1ML

Diniconazole
CAS#83657-24-3
Chemical Name: (E)-(RS)-1-(2,4-Dichlorophenyl)-4,4-dimethyl-2-(1H-1,2,4-triazol-1-yl)pent-1-en-3-ol

PS: Solid

MF: C₁₅H₁₇Cl₂N₃O


MW: 326.25

MP: 134°C

Alternate Name: Spotless(TM)

Part Number	Matrix	Size
N-11782-100MG	Neat	100MG
S-11782M1-1ML	100ug/ml in Methanol	1ML

Esbiothrin
CAS#84030-86-4
Chemical Name: (R,S)-3-Allyl-2-methyl-4-oxocyclopent-2-enyl-(1R)-trans-2,2-dimethyl-3-(2-methylprop-1-enyl)cyclopropanecarboxylate

PS: Liquid


MF: C₁₉H₂₆O₃
MW: 302.4


Part Number	Matrix	Size
N-11864-100MG	Neat	100MG

Quinclorac
CAS#84087-01-4
Chemical Name: 3,7-Dichloroquinoline-8-carboxylic acid

PS: Solid

MF: C₁₀H₅Cl₂NO₂
MW: 242.1

MP: 274°C


Part Number	Matrix	Size
N-13169-500MG	Neat	500MG
S-13169A1-1ML	100ug/ml in Acetonitrile	1ML
S-13169T1-1ML	100ug/ml t-Butyl methyl ether	1ML


Itraconazole
CAS#84625-61-6

Chemical Name: (±)-1-[(RS)-sec-butyl]-4-[p-[4-[p-[[[(2R,4S)-rel-2-(2,4-dichlorophenyl)-2-(1H-1,2,4-triazol-1-yl)methyl]-1,3-dioxolan-4-yl]methoxy]phenyl]-1-piperazinyl]phenyl]-7-2,1,2,4-triazolin-5-one

PS: Solid

MF: C₃₅H₃₈Cl₂N₈O₄

MW: 705.6



Part Number	Matrix	Size
N-13958-100MG	Neat	100MG

Flusilazole
CAS#85509-19-9


Chemical Name: 1-[[Bis(4-fluorophenyl)methylsilyl]methyl]-1H-1,2,4-triazole

PS: Solid

MF: C₁₆H₁₅F₂N₃Si

MW: 315.43

MP: 55°C



Alternate Name: Nustar(tm)

Part Number	Matrix	Size
N-12003-100MG	Neat	100MG
S-12003A1-1ML	100ug/ml in Acetonitrile	1ML
S-12003U1-1ML	100ug/ml in Toluene	1ML

Primisulfuron-methyl
CAS#86209-51-0


Chemical Name: 2-Chloro-2',6'-diethyl-N-(2-propoxyethyl)acetanilide

PS: Solid

MF: C₁₅H₁₂F₄N₄O₇S

MW: 468.34

MP: 194.8°C



Alternate Name: Beacon(TM)

Part Number	Matrix	Size
N-13090-100MG	Neat	100MG
S-13090A1-1ML	100ug/ml in Acetonitrile	1ML

Hexaflumuron
CAS#86479-06-3


Chemical Name: 1-[3,5-Dichloro-4-(1,1,2,2-tetrafluoroethoxy)phenyl]3-(2,6-difluorobenzoyl)urea

PS: Solid

MF: C₁₆H₈Cl₂F₆N₂O₃

MW: 461.15

MP: 202°C



Alternate Name: Consult(TM)


Part Number	Matrix	Size
N-12167-250MG	Neat	250MG
S-12167A1-1ML	100ug/ml in Acetonitrile	1ML
S-12167U1-1ML	100ug/ml in Toluene	1ML

Imibenconazole
CAS#86598-92-7

PS: Liquid

MF: C₁₇H₁₃Cl₃N₄S

MW: 411.7



Part Number	Matrix	Size
S-13968A1-1ML	100ug/ml in Acetonitrile	1ML

Propisochlor
CAS#86763-47-5

Chemical Name: 2-chloro-N-(2-ethyl-6-methylphenyl)-N-[(1-methylethoxy)methyl]acetamide


PS: Solid

MF: C₁₂H₂₀N₄O₂

MW: 252.36

MP: 116°C

FP: 175°C



Part Number	Matrix	Size
N-13823-50MG	Neat	50MG

Diethofencarb
CAS#87130-20-9

Chemical Name: 1-Methylethyl (3,4-diethoxyphenyl) carbamate


PS: Solid

MF: C₁₄H₂₁NO₄

MW: 267.3

MP: 100.3°C

FP: 140°C



Alternate Name: Isopropyl 3,4-diethoxycarbanilate

Part Number	Matrix	Size
N-11691-250MG	Neat	250MG
S-11691M1-1ML	100ug/mL in Methanol	1ML

S-Metolachlor
CAS#87392-12-9

Chemical Name: 2-chloro-N-(2-ethyl-6-methylphenyl)-N-[(1S)-2-methoxy-1-methylethyl]acetamide


PS: Liquid

MF: C₁₅H₂₂ClNO₂

MW: 283.8

MP: -61.1°C

FP: 190°C




Part Number	Matrix	Size
N-13191-100MG	Neat	100MG
S-13191A1-1ML	100ug/mL in Acetonitrile	1ML
S-13191T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Flumiclorac-pentyl CAS#87546-18-7

Chemical Name: Pentyl[2-chloro-5-(cyclohex-1-ene-1,2-dicarboximido)-4-fluorophenoxy]acetate

PS: Solid

MF: C₂₁H₂₃ClFNO₅
MW: 423.9

FP: 68°C



Part Number	Matrix	Size
N-11990-500MG	Neat	500MG
S-11990A1-1ML	100ug/ml in Acetonitrile	1ML
S-11990T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Dimethenamid CAS#87674-68-8

Chemical Name: (RS)-2-Chloro-N-(2,4-dimethyl-3-thienyl)-N-(2-methoxy-1-methylethyl)acetamide

PS: Liquid


MF: C₁₂H₁₈ClNO₂S

MW: 275.8

Alternate Name: Frontier(TM)

Part Number	Matrix	Size
N-11755-100MG	Neat	100MG
S-11755A1-1ML	100ug/ml in Acetonitrile	1ML
S-11755U1-1ML	100ug/ml in Toluene	1ML

Cinmethylin solution CAS#87818-31-3

PS: Liquid

MF: C₁₈H₂₆O₂
MW: 274.4



Part Number	Matrix	Size
S-14050M1-1ML	100ug/ml in Methanol	1ML

Tralkoxydim CAS#87820-88-0

Chemical Name: 2-[1-(Ethoxyimino) propyl]-3-hydroxy-5-mesitylcyclohex-2-enone

PS: Solid

MF: C₂₀H₂₇NO₃
MW: 329.4

MP: 106°C

Alternate Name: Grasp(TM)

Part Number	Matrix	Size
N-13587-500MG	Neat	500MG
S-13587A1-1ML	100ug/ml in Acetonitrile	1ML
S-13587U1-1ML	100ug/ml in Toluene	1ML


Pyriphenox CAS#88283-41-4

Chemical Name: 1-(2,4-Dichlorophenyl)-2-(3-pyridinyl)ethanone O-methylloxime

PS: Liquid

MF: C₁₄H₁₂Cl₂N₂O

MW: 295.16


FP: 145°C

Alternate Name: PYRIFENOX

Part Number	Matrix	Size
N-12990-10MG	Neat	10MG
S-12990M1-1ML	100ug/mL in Methanol	1ML

Cloquintocet CAS#88349-88-6

Chemical Name: 1-methylhexyl (5-chloroquinolin-8-yloxy)acetate

PS: Solid


MF: C₁₁H₈ClNO₃
MW: 237.64


Part Number	Matrix	Size
N-13820-50MG	Neat	50MG

Fluxofenim CAS#88485-37-4

Chemical Name: 1-(4-chlorophenyl)-2,2,2-trifluoroethanone O-(1,3-dioxolan-2-yl)methylloxime

PS: Liquid

MF: C₁₂H₁₁ClF₃NO₃
MW: 309.7



Part Number	Matrix	Size
N-13940-250MG	Neat	250MG
S-13940M1-1ML	100ug/mL in Methanol	1ML

Myclobutanil CAS#88671-89-0

Chemical Name: (RS)-2-(4-chlorophenyl)-2-(1H-1,2,4-triazol-1-ylmethyl)hexane-nitrile

PS: Solid

MF: C₁₅H₁₇ClN₄
MW: 288.78


MP: 70.9°C


Part Number	Matrix	Size
N-13261-100MG	Neat	100MG
S-13261A1-1ML	100ug/ml in Acetonitrile	1ML
S-13261U1-1ML	100ug/ml in Toluene	1ML

Prohexadione
CAS#88805-35-0
Chemical Name: 3,5-Dioxo-4-propionylcyclohexanecarboxylic acid

PS: Solid


MF: C₁₀H₁₂O₅
MW: 212.2

MP: 123°C


Part Number	Matrix	Size
N-13099-100MG	Neat	100MG

Ferimzone
CAS#89269-64-7
Chemical Name: (Z)-2'-methylacetophenone 4,6-dimethylpyrimidin-2-ylhydrazone

PS: Solid


MF: C₁₅H₁₈N₄
MW: 254.33


Part Number	Matrix	Size
N-14030-100MG	Neat	100MG
S-14030A1-1ML	100ug/ml in Acetonitrile	1ML

Pyraclofos
CAS#89784-60-1
Chemical Name: (RS)-O-1-(4-CHLOROPHENYL)PYRAZOL-4-YL O-ETHYL S-PROPYL PHOSPHOROTHIOATE

PS: Liquid

MF: C₁₄H₁₈ClN₂O₃PS

MW: 360.82

Alternate Name: BOLTAGE


Part Number	Matrix	Size
N-13142-10MG	Neat	10MG
S-13142A1-1ML	100ug/ml in Acetonitrile	1ML

Flocoumafen
CAS#90035-08-8
Chemical Name: 4-Hydroxy-3-[1,2,3,4-tetrahydro-3-[4-(4-trifluoromethylbenzyl-oxy)phenyl]-1-naphthyl] coumarin

PS: Solid

MF: C₃₃H₂₅F₃O₄
MW: 542.54

MP: 181°C


FP: 23°C

Alternate Name: Storm(TM)

Part Number	Matrix	Size
N-11975-10MG	Neat	10MG
S-11975A1-1ML	100ug/ml in Acetonitrile	1ML

Quinmerac
CAS#90717-03-6
Chemical Name: 7-Chloro-3-methyl-8-quinolinecarboxylic acid

PS: Solid

MF: C₁₁H₈ClNO₂
MW: 221.64

MP: 244°C



Part Number	Matrix	Size
N-13925-10MG	Neat	10MG

Chlorimuron ethyl
CAS#90982-32-4
Chemical Name: Ethyl-2-[[[(4-chloro-6-methoxy-2-pyrimidinyl)amino]carbonyl]amino]sulfonyl]benzoate

PS: Solid


MF: C₁₅H₁₅ClN₄O₆S

MW: 414.82

MP: 181°C



Part Number	Matrix	Size
N-11432-100MG	Neat	100MG
S-11432A1-1ML	100ug/ml in Acetonitrile	1ML

lambda-Cyhalothrin
CAS#91465-08-6
PS: Solid

MF: C₂₃H₁₉ClF₃NO₃
MW: 449.86

Alternate Name: Karate (TM)

Part Number	Matrix	Size
N-12307-100MG	Neat	100MG
S-12307A1-1ML	100ug/ml in Acetonitrile	1ML
S-12307J1-1ML	100ug/ml in Hexane	1ML

Trietazine-ethylamine
CAS#93567-85-2
PS: Liquid

MF: C₁₁H₂₂N₆
MW: 238.3



Part Number	Matrix	Size
N-13989-100MG	Neat	100MG

Pyrazosulfuron-ethyl
CAS#93697-74-6
Chemical Name: Ethyl 5-(4,6-dimethoxypyrimidin-2-ylcarbamoylsulfamoyl)-1-methylpyrazole-4-carboxylate

PS: Solid

MF: C₁₄H₁₈N₆O₇S

MW: 414.44

MP: 181°C



Part Number	Matrix	Size
N-13147-100MG	Neat	100MG
S-13147A1-1ML	100ug/ml in Acetonitrile	1ML

Prosulfuron
CAS#94125-34-5
Chemical Name: 1-(4-Methoxy-6-methyl-1,3,5-triazin-2-yl)-3-[2-(3,3,3-trifluoropropyl)phenylsulfonyl]urea

PS: Solid

MF: C₁₅H₁₆F₃N₅O₄S

MW: 419.38

MP: 155°C



Part Number	Matrix	Size
N-13137-100MG	Neat	100MG
S-13137A1-1ML	100ug/ml in Acetonitrile	1ML

Cyproconazole
CAS#94361-06-5
Chemical Name: (2RS,3RS;2RS,3SR)-2-(4-Chlorophenyl)-3-cyclopropyl-1-(1H-1,2,4-triazol-1-yl)butan-2-ol

PS: Solid

MF: C₁₅H₁₈ClN₃O

MW: 291.8

MP: 106.2°C

Alternate Name: Sentinel(TM)


Part Number	Matrix	Size
N-11547-50MG	Neat	50MG
S-11547A1-1ML	100ug/ml in Acetonitrile	1ML
S-11547E1-1ML	100ug/ml in Cyclohexane	1ML

Cinosulfuron
CAS#94593-91-6
Chemical Name: 1-(4,6-dimethoxy-1,3,5-triazin-2-yl)-3-[2-(2-methoxyethoxy)phenylsulfonyl]urea

PS: Solid

MF: C₁₅H₁₉N₅O₇S

MW: 413.4

MP: 127°C



Part Number	Matrix	Size
N-11470-100MG	Neat	100MG

Trinexapac-ethyl
CAS#95266-40-3
Chemical Name: Ethyl 4-cyclopropyl(hydroxy)methylene-3,5-dioxocyclohexanecarboxylate

PS: Solid

MF: C₁₃H₁₆O₅
MW: 252.3


MP: 36°C

FP: 133°C


Part Number	Matrix	Size
N-13706-250MG	Neat	250MG
S-13706A1-1ML	100ug/ml in Acetonitrile	1ML
S-13706T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Cadusafos
CAS#95465-99-9
Chemical Name: O-Ethyl-S,S-bis(1-methylpropyl)phosphorodithioate


PS: Liquid

MF: C₁₀H₂₃O₂PS₂
MW: 270.4

Alternate Name: Apache(TM)

Part Number	Matrix	Size
N-11392-100MG	Neat	100MG
S-11392A1-1ML	100ug/ml in Acetonitrile	1ML
S-11392U1-1ML	100ug/ml in Toluene	1ML

Fenoxaprop
CAS#95617-09-7
Chemical Name: 2-[4-[(6-chloro-2-benzoxazolyl)oxy]phenoxy]propanoic acid

PS: Solid


MP: 166°C


Part Number	Matrix	Size
N-13830-50MG	Neat	50MG

Pyriproxyfen
CAS#95737-68-1
Chemical Name: 2-(1-Methyl-2-(4-phenoxyphenoxy)ethoxy)pyridine

PS: Solid

MF: C₂₀H₁₉NO₃
MW: 321.37

MP: 45°C

Alternate Name: Admiral(TM)

Part Number	Matrix	Size
N-13160-100MG	Neat	100MG
S-13160M1-1ML	100ug/ml in Methanol	1ML

Haloxypop-P (free acid) CAS#95977-29-0


Chemical Name: (R)-2-[4-[3-chloro-5-(trifluoromethyl)-2-pyridyloxy]phenoxy]propanoic acid

PS: Solid

MF: C₁₅H₁₁ClF₃NO₄

MW: 361.7

MP: 107°C



Alternate Name: Haloxypop-R (free acid)

Part Number	Matrix	Size
N-12145-50MG	Neat	50MG

Tebupirimfos CAS#96182-53-5


Chemical Name: Phosphorothioic acid, O-[2-(1,1-dimethylethyl)-5-pyrimidinyl] O-ethyl O-(1-methylethyl) ester

PS: Liquid

MF: C₁₃H₂₃N₂O₃PS

MW: 318.4

FP: 93°C



Alternate Name: Tebupirimphos

Part Number	Matrix	Size
N-13503-100MG	Neat	100MG
S-13503A1-1ML	100ug/ml in Acetonitrile	1ML
S-13503U1-1ML	100ug/ml in Toluene	1ML

Pyridaben CAS#96489-71-3


Chemical Name: 2-tert-Butyl-5-(4-tert-butylbenzylthio)-4-chloropyridazin-3(2H)-one

PS: Solid

MF: C₁₉H₂₅ClN₂O₂S

MW: 364.94

MP: 111°C



Alternate Name: Sanmite*

Part Number	Matrix	Size
N-13152-25MG	Neat	25MG
S-13152M1-1ML	100ug/ml in Methanol	1ML


Thenylchlor Solution CAS#96491-05-3

Chemical Name: 2-chloro-N-(3-methoxy-2-thenyl)-2',6'-dimethylacetanilide

PS: Liquid

MF: C₁₆H₁₈ClNO₂S

MW: 323.8



Part Number	Matrix	Size
S-13966A1-1ML	100ug/ml in Acetonitrile	1ML


Flurtamone Solution CAS#96525-23-4

Chemical Name: (RS)-5-methylamino-2-phenyl-4-(a,a,a-trifluoro-m-tolyl)furan-3(2H)-one

PS: Liquid

MF: C₁₈H₁₄F₃NO₂

MW: 333.3



Part Number	Matrix	Size
S-13981A1-1ML	100ug/ml in Acetonitrile	1ML

Ethametsulfuron-methyl CAS#97780-06-8


Chemical Name: Methyl 2-((4-ethoxy-6-methylamino-1,3,5-triazin-2-yl)carbamoyl-sulfamoyl)benzoate

PS: Solid

MF: C₁₅H₁₈N₆O₆S

MW: 410.4

MP: 194°C



Alternate Name: Muster(TM)

Part Number	Matrix	Size
N-11866-100MG	Neat	100MG
S-11866A1-1ML	100ug/ml in Acetonitrile	1ML

Dithiopyr CAS#97886-45-8


Chemical Name: S,S'-Dimethyl 2-difluoromethyl-4-isobutyl-6-trifluoromethylpyridine-3,5-dicarbothioate

PS: Solid

MF: C₁₅H₁₆F₅NO₂S₂

MW: 401.4

MP: 65°C



Alternate Name: Dimension(TM)

Part Number	Matrix	Size
N-11822-100MG	Neat	100MG
S-11822A1-1ML	100ug/ml in Acetonitrile	1ML
S-11822U1-1ML	100ug/ml in Toluene	1ML


Benoxacor CAS#98730-04-2

Chemical Name: 4-(2,2-dichloroacetyl)-3,4-dihydro-3-methyl-2H-1,4-benzoxazine

PS: Liquid

MF: C₁₁H₁₁Cl₂NO₂


MW: 260.1



Part Number	Matrix	Size
N-12957-10MG	Neat	10MG

Fosthiazate
CAS#98886-44-3
Chemical Name: S-sec-Butyl O-ethyl (2-oxo-1,3-thiazolidine-3-yl)phosphonothioate

PS: Liquid


MF: C₉H₁₈NO₃PS₂
MW: 283.3

Alternate Name: Cierto(TM)

Part Number	Matrix	Size
N-12020-30MG	Neat	30MG
S-12020A1-1ML	100ug/ml in Acetonitrile	1ML
S-12020U1-1ML	100ug/ml in Toluene	1ML

Flumetsulam
CAS#98967-40-9
Chemical Name: 2,6'-Difluoro-5-methyl[1,2,4]triazolo[1,5-a]pyrimidine-2-sulfonanilide

PS: Solid

MF: C₁₂H₉F₂N₅O₂S

MW: 325.3



Part Number	Matrix	Size
N-11989-100MG	Neat	100MG
S-11989A1-1ML	100ug/ml in Acetonitrile	1ML
S-11989U1-1ML	100ug/ml in Toluene	1ML

Sulcotrione
CAS#99105-77-8
Chemical Name: 2-[2-chloro-4-(methylsulfonyl)benzoyl]-1,3-cyclohexanedione

PS: Solid

MF: C₁₄H₁₃ClO₅S

MW: 328.8


MP: 139°C


Part Number	Matrix	Size
N-12890-100MG	Neat	100MG

Clethodim
CAS#99129-21-2
Chemical Name: 2-[1-[[[(2E)-3-chloro-2-propen-1-yl]oxy]imino]propyl]-5-[2-(ethylthio)propyl]-3-hydroxy-2-cyclohexen-1-one

PS: Liquid

MF: C₁₇H₂₆ClNO₃S


MW: 359.95

Alternate Name: Select*

Part Number	Matrix	Size
N-11488-100MG	Neat	100MG

Cloquintocet-mexyl
CAS#99607-70-2
Chemical Name: 1-Methylhexyl (5-chloroquinolin-8-yloxy)acetate

PS: Solid

MF: C₁₈H₂₂ClNO₃
MW: 335.8


MP: 69.4°C


Part Number	Matrix	Size
N-11491-100MG	Neat	100MG

Isofenphos-Methyl
CAS#99675-03-3
Chemical Name: 1-methylethyl 2-[[methoxy[(1-methylethyl)amino]phosphinothioyl]oxy]benzoate

PS: Liquid

MW: 331.37


FP: 194°C


Part Number	Matrix	Size
N-12250-50MG	Neat	50MG

Quizalofop-P-ethyl
CAS#100646-51-3
Chemical Name: Ethyl (R)-2-[4-(6-chloroquinoxalin-2-yloxy)phenoxy]propionate

PS: Solid

MF: C₁₉H₁₇ClN₂O₄
MW: 372.83


MP: 76°C

Alternate Name: Assure II(TM)

Part Number	Matrix	Size
N-13175-100MG	Neat	100MG
S-13175A1-1ML	100ug/ml in Acetonitrile	1ML
S-13175U1-1ML	100ug/ml in Toluene	1ML

Halosulfuron-methyl
CAS#100784-20-1
Chemical Name: Methyl-3-chloro-5(4,6-dimethoxypyrimidin-2-ylcarbamoysulfamoyl)-1-methylpyrazole-4-carboxylate

PS: Solid


MF: C₁₃H₁₅ClN₆O₇S

MW: 434.8


Part Number	Matrix	Size
N-12142-100MG	Neat	100MG
S-12142A1-1ML	100ug/ml in Acetonitrile	1ML

Acrinathrin
CAS#101007-06-1
Chemical Name: (S)-a-cyano-3-phenoxybenzyl (Z)-(1R,3S)-2,2-dimethyl-3-[2-(2,2,2-trifluoro-1-trifluoromethylethoxycarbonyl)vinyl]cyclopropanecarboxylate

PS: Solid


MF: C₂₆H₂₁F₆NO₅
MW: 541.4


Part Number	Matrix	Size
N-13928-100MG	Neat	100MG
S-13928A1-1ML	100ug/ml in Acetonitrile	1ML

Tribenuron methyl
CAS#101200-48-0
Chemical Name: Methyl-2-[4-methoxy-6-methyl-1,3,5-triazin-2-yl(methyl)carbamoylsulfamoyl]benzoate

PS: Solid

MF: C₁₅H₁₇N₅O₆S

MW: 395.4

Alternate Name: METOMETURON


Part Number	Matrix	Size
N-13640-250MG	Neat	250MG
S-13640A1-1ML	100ug/ml in Acetonitrile	1ML

Cycloxydime
CAS#101205-02-1
Chemical Name: (±)-2-[1-(ethoxyimino)butyl]-3-hydroxy-5-thian-3-ylcyclohex-2-enone

PS: Solid

MF: C₁₇H₂₇NO₃S

MW: 325.47


MP: 41°C


Part Number	Matrix	Size
N-13916-10MG	Neat	10MG

Flufenoxuron
CAS#101463-69-8
Chemical Name: N-((4-(2-Chloro-4-(trifluoromethyl)phenoxy)-2-fluorophenyl)carbamoyl)-2,6-difluorobenzamide

PS: Solid


MF: C₂₁H₁₁ClF₆N₂O₃
MW: 488.79

MP: 169°C

Alternate Name: Cascade(TM)

Part Number	Matrix	Size
N-11987-50MG	Neat	50MG
S-11987A1-1ML	100ug/ml in Acetonitrile	1ML

tau-Fluvalinate
CAS#102851-06-9
Chemical Name: (RS)-a-cyano-3-phenoxybenzyl N-(2-chloro-a,a,a-trifluoro-p-tolyl)-D-valinate


PS: Liquid

MF: C₂₆H₂₂ClF₃N₂O₂
MW: 502.93

Alternate Name: Mavrik(TM)

Part Number	Matrix	Size
N-13263-100MG	Neat	100MG
S-13263A1-1ML	100ug/ml in Acetonitrile	1ML
S-13263K1-1ML	100ug/ml in Isooctane	1ML

Lufenuron
CAS#103055-07-8
Chemical Name: (RS)-1-[2,5-Dichloro-4-(1,1,2,3,3,3-hexafluoropropoxy)phenyl]-3-(2,6-difluorobenzoyl)urea

PS: Solid


MF: C₁₇H₈Cl₂F₈N₂O₃
MW: 511.2


Part Number	Matrix	Size
N-12325-100MG	Neat	100MG
S-12325A1-1ML	100ug/ml in Acetonitrile	1ML

Flumioxazin
CAS#103361-09-7
Chemical Name: N-(7-Fluoro-3,4-dihydro-3-oxo-4-prop-2-ynyl-2H-1,4-benzoxazin-6-yl)cyclohex-1-ene-1,2-dicarboxamide

PS: Solid


MF: C₁₉H₁₅FN₂O₄
MW: 354.3

MP: 201°C

Alternate Name: Flumyazin(TM)

Part Number	Matrix	Size
N-11991-100MG	Neat	100MG
S-11991A1-1ML	100ug/ml in Acetonitrile	1ML
S-11991U1-1ML	100ug/ml in Toluene	1ML

2,4-D propylene glycol ester
CAS#103564-78-9
Chemical Name: 2,4-Dichlorophenoxy acetic acid, propylene glycol ester

PS: Liquid


MF: C₁₁H₁₂Cl₂O₄
MW: 279.12


Part Number	Matrix	Size
N-10620-250MG	Neat	250MG
S-10620A1-1ML	100ug/ml in Acetonitrile	1ML
S-10620J1-1ML	100ug/ml in Hexane	1ML

Carproamid Soluton
CAS#104030-54-8
Chemical Name: 2,2-Dichloro-N-[1-(4-chlorophenyl)ethyl]-1-ethyl-3-methylcyclopropanecarboxamide

PS: Liquid

MF: C₁₃H₁₈Cl₃NO

MW: 334.7



Part Number	Matrix	Size
S-13974A1-1ML	100ug/ml in Acetonitrile	1ML

Flazasulfuron
CAS#104040-78-0
Chemical Name: 1-(4,6-Dimethoxypyrimidin-2-yl)-3-(3-trifluoromethyl-2-pyridylsulfonyl)urea

PS: Solid

MF: C₁₃H₁₂F₃N₅O₅S


MW: 407.3

MP: 180°C


Part Number	Matrix	Size
N-11974-50MG	Neat	50MG

Imazapic
CAS#104098-48-8
Chemical Name: (RS)-2-(4,5-Dihydro-4-isopropyl-4-methyl-5-oxoimidazol-2-yl)-5-methylnicotinic acid

PS: Solid

MF: C₁₄H₁₇N₃O₃
MW: 275.3

Alternate Name: AC 263, 222(TM)


Part Number	Matrix	Size
N-12202-250MG	Neat	250MG
S-12202A1-1ML	100ug/ml in Acetonitrile	1ML

Mesotrione
CAS#104206-82-8
Chemical Name: 2-(4-Mesyl-2-nitrobenzoyl)cyclohexane-1,3-dione

PS: Solid

MF: C₁₄H₁₃NO₇S

MW: 339.3


MP: 165°C


Part Number	Matrix	Size
N-12377-100MG	Neat	100MG

Silafluofen
CAS#105024-66-6
Chemical Name: 4-[3-[(4-Ethoxyphenyl)dimethylsilyl]propyl]-1-fluoro-2-phenoxybenzene

PS: Liquid

MF: C₂₅H₂₉FO₃Si


MW: 408.6


Part Number	Matrix	Size
N-13914-10MG	Neat	10MG

Clodinafop-propargyl
CAS#105512-06-9
Chemical Name: Prop-2-ynyl (R)-2-[4-(5-chloro-3-fluoropyridin-2-yloxy)phenoxy]propionate

PS: Solid


MF: C₁₅H₁₂F₂N₄O₃
MW: 349.8

MP: 59.5°C


Part Number	Matrix	Size
N-11489-100MG	Neat	100MG
S-11489A1-1ML	100ug/mL in Acetonitrile	1ML
S-11489U1-1ML	100ug/ml in Toluene	1ML

Pyrimidifen Solution
CAS#105779-78-0
Chemical Name: 5-Chloro-N-[2-[4-(2-ethoxyethyl)-2,3-dimethylphenoxy]ethyl]-6-ethylpyrimidin-4-amine

PS: Liquid


MF: C₂₀H₂₈ClN₃O₂
MW: 377.9


Part Number	Matrix	Size
S-13963A1-1ML	100ug/ml in Acetonitrile	1ML

Flusulfamide
CAS#106917-52-6
Chemical Name: 2',4-dichloro-a,a,a-trifluoro-4'-nitro-m-toluenesulfonamide

PS: Solid

MF: C₁₃H₇Cl₂F₃N₂O₄S

MW: 415.2



Part Number	Matrix	Size
N-14029-100MG	Neat	100MG
S-14029A1-1ML	100ug/mL in Acetonitrile	1ML

Tebuconazole
CAS#107534-96-3
Chemical Name: a-(2-(4-Chlorophenyl)ethyl)-a-(1,1-dimethylethyl)-1H-1,2,4-triazole-1-ethanol

PS: Solid

MF: C₁₆H₂₂ClN₃O

MW: 307.82


MP: 105°C

Alternate Name: Folicur(TM)

Part Number	Matrix	Size
N-12006-250MG	Neat	250MG
S-12006A1-1ML	100ug/ml in Acetonitrile	1ML
S-12006K1-1ML	100ug/ml in Isooctane	1ML

Diflufenzopyr
CAS#109293-97-2
Chemical Name: 2-[1-[4-(3,5-difluorophenyl)semicarbazono]ethyl]nicotinic acid

PS: Solid

MF: C₁₅H₁₂F₂N₄O₃
MW: 334.3


MP: 135.5°C


Part Number	Matrix	Size
N-11724-100MG	Neat	100MG

Diflufenzopyr sodium salt
CAS#109293-98-3
Chemical Name: sodium 2-[(1E)-1-[2-[[[3,5-difluorophenyl]amino]carbonyl]hydrazinylidene]ethyl]-3-pyridinecarboxylate

PS: Solid

MF: C₁₅H₁₁F₂N₄NaO₃
MW: 356.26


MP: 150°C


Part Number	Matrix	Size
N-13822-100MG	Neat	100MG

Mepanipyrim
CAS#110235-47-7
Chemical Name: 2-Pyrimidinamine, 4-methyl-N-phenyl-6-(1-propyn-1-yl)-

PS: Solid

MF: C₁₄H₁₃N₃
MW: 223.3


MP: 132.8°C

Alternate Name: N-(4-methyl-6-prop-1-ynylpyrimidin-2-yl)aniline

Part Number	Matrix	Size
N-12369-100MG	Neat	100MG

Dimethomorph
CAS#110488-70-5
Chemical Name: (E,Z)-4-[3-(4-chlorophenyl)-3-(3,4-dimethoxyphenyl)acryloyl]-morpholine

PS: Solid

MF: C₂₁H₂₂ClNO₄
MW: 387.9


MP: 127°C

Alternate Name: Acrobat(TM)

Part Number	Matrix	Size
N-11759-100MG	Neat	100MG
S-11759A1-1ML	100ug/ml in Acetonitrile	1ML

Propaquizafop
CAS#111479-05-1
Chemical Name: 2-isopropylideneamino-oxyethyl (R)-2-[4-(6-chloroquinoxalin-2-yl)oxy]phenoxy]propionate

PS: Solid


MP: 66.3°C

MF: C₂₂H₂₂ClN₃O₅
MW: 443.9


Part Number	Matrix	Size
N-14007-250MG	Neat	250MG
S-14007A1-1ML	100ug/mL in Acetonitrile	1ML

Halfenprox
CAS#111872-58-3
Chemical Name: 2-(4-bromodifluoromethoxyphenyl)-2-methylpropyl 3-phenoxybenzyl ether

PS: Liquid

MF: C₂₄H₂₃BrF₂O₃
MW: 477.3

Alternate Name: BROFENPROX


Part Number	Matrix	Size
N-13079-10MG	Neat	10MG
S-13079M1-1ML	100ug/mL in Methanol	1ML

Thiacloprid
CAS#111988-49-9
Chemical Name: (3-((6-CHLORO-3-PYRIDINYL)METHYL)-2-THIAZOLIDINYLIDENE) CYANAMIDE

PS: Solid

MF: C₁₀H₉ClN₄S

MW: 252.7

MP: 136°C



Part Number	Matrix	Size
N-13561-100MG	Neat	100MG
S-13561A1-1ML	100ug/mL in Acetonitrile	1ML
S-13561U1-1ML	100ug/ml in Toluene	1ML

Nicosulfuron
CAS#111991-09-4
Chemical Name: 2-[[4,6-dimethoxypyrimidin-2-ylcarbamoyl]sulfamoyl]-N,N-dimethylnicotinamide

PS: Solid

MF: C₁₅H₁₈N₆O₆S

MW: 410.4

MP: 169-172°C

Alternate Name: Accent (TM)

Part Number	Matrix	Size
N-10998-100MG	Neat	100MG


Triazamate
CAS#112143-82-5
Chemical Name: ethyl 2-[[1-[(dimethylamino)carbonyl]-3-(1,1-dimethylethyl)-1H-1,2,4-triazol-5-yl]thio]acetate

PS: Solid

MF: C₁₃H₂₂N₄O₃S

MW: 314.4


MP: 52.1-53.3°C

FP: 189°C


Part Number	Matrix	Size
N-13956-100MG	Neat	100MG

Halofenozide
CAS#112226-61-6
Chemical Name: N-tert-butyl-N-(4-Chlorobenzoyl)benzohydrazide

PS: Solid

MF: C₁₈H₁₉ClN₂O₂
MW: 330.8

Alternate Name: Mach 2(TM)


Part Number	Matrix	Size
N-12141-500MG	Neat	500MG
S-12141M1-1ML	100ug/ml in Methanol	1ML

Tetraconazole
CAS#112281-77-3
Chemical Name: 1-[2-(2,4-Dichlorophenyl)-3-(1,1,2,2-tetrafluoroethoxy)propyl]-1H-1,2,4-triazole

PS: Liquid

MF: C₁₃H₁₁Cl₂F₄N₃O

MW: 372.1


MP: 6°C

Alternate Name: Arpege(TM)

Part Number	Matrix	Size
N-13538-100MG	Neat	100MG
S-13538A1-1ML	100ug/ml in Acetonitrile	1ML
S-13538U1-1ML	100ug/ml in Toluene	1ML

Tebufenozide
CAS#112410-23-8
Chemical Name: 1-tert-Butyl-1-(3,5-dimethylbenzoyl)-2-(4-ethylbenzoyl)hydrazine

PS: Solid


MF: C₂₂H₂₈N₂O₂
MW: 352.47

MP: 191°C

Alternate Name: Confirm(TM)

Part Number	Matrix	Size
N-13501-100MG	Neat	100MG
S-13501A1-1ML	100ug/ml in Acetonitrile	1ML

Cyclanilide
CAS#113136-77-9
PS: Solid


MF: C₁₁H₉Cl₂NO₃
MW: 274.1

MP: 195.5°C


Part Number	Matrix	Size
N-11521-100MG	Neat	100MG
S-11521A1-1ML	100ug/ml in Acetonitrile	1ML
S-11521U1-1ML	100ug/ml in Toluene	1ML

Z-Dimethomorph
CAS#113210-98-3
PS: Solid

MF: C₂₁H₂₂ClNO₄
MW: 387.9


MP: 127°C

Alternate Name: Z-Dimethomorphe

Part Number	Matrix	Size
N-13752-100MG	Neat	100MG
S-13752M1-1ML	100ug/ml in Methanol	1ML

Imazamox
CAS#114311-32-9
Chemical Name: (RS)-2-(4-Isopropyl-4-methyl-5-oxo-2-imidazolin-2-yl)-5-methoxy-methylnicotinic acid

PS: Solid


MF: C₁₅H₁₉N₃O₃
MW: 305.3

MP: 166°C

Alternate Name: Raptor(TM)

Part Number	Matrix	Size
N-12201-250MG	Neat	250MG
S-12201A1-1ML	100ug/ml in Acetonitrile	1ML
S-12201T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

Fenbuconazole
CAS#114369-43-6
Chemical Name: 4-(4-Chlorophenyl)-2-phenyl-2-(1H-1,2,4-triazol-1-ylmethyl)butyronitrile

PS: Solid


MF: C₁₉H₁₇ClN₄
MW: 336.8

Alternate Name: FENETHANIL

Part Number	Matrix	Size
N-11950-100MG	Neat	100MG
S-11950M1-1ML	100ug/ml in Methanol	1ML

Fenoxanil
CAS#115852-48-7
Chemical Name: Propanamide, N-(1-cyano-1,2-dimethylpropyl)-2-(2,4-dichlorophenoxy)-

PS: Solid

MF: C₁₅H₁₈Cl₂N₂O₂
MW: 329.2

MP: 70°C



Part Number	Matrix	Size
N-11956-100MG	Neat	100MG

Bromuconazole
CAS#116255-48-2
Chemical Name: 1-((4-Bromo-2-(2,4-dichlorophenyl)tetrahydro-2-furanyl)methyl)-1H-1,2,4-triazole

PS: Solid

MF: C₁₃H₁₂BrCl₂N₃O

MW: 377.09


MP: 84°C

Alternate Name: Bromoconazole

Part Number	Matrix	Size
N-11347-100MG	Neat	100MG

Novaluron
CAS#116714-46-6
Chemical Name: (+/-)-N-[[[3-Chloro-4-[1,1,2-trifluoro-2-(trifluoromethoxy)ethoxy]-phenyl]amino]carbonyl]-2,6-difluorobenzamide

PS: Solid

MF: C₁₇H₉ClF₈N₂O₄
MW: 492.7


MP: 176°C

Alternate Name: Rimon(TM)

Part Number	Matrix	Size
N-12670-100MG	Neat	100MG
S-12670A1-1ML	100ug/mL in Acetonitrile	1ML
S-12670U1-1ML	100ug/ml in Toluene	1ML

Fluthiacet-methyl
CAS#117337-19-6
Chemical Name: methyl 2-[[[2-chloro-4-fluoro-5-[(tetrahydro-3-oxo-1H,3H-[1,3,4]thiadiazolo[3,4-a]pyridazin-1-ylidene)amino]phenyl]thio]acetate

PS: Solid

MF: C₁₅H₁₅ClF₃N₃O₃S₂
MW: 403.9


MP: 104.6°C


Part Number	Matrix	Size
N-13129-25MG	Neat	25MG
S-13129A1-1ML	100ug/mL in Acetonitrile	1ML

Picoxystrobin
CAS#117428-22-5
Chemical Name: methyl (E)-3-methoxy-2-[2-(6-trifluoromethyl-2-pyridyloxymethyl)phenyl]acrylate

PS: Solid

MF: C₁₈H₁₆F₃N₂O₄
MW: 367.3

MP: 75°C



Part Number	Matrix	Size
N-13927-100MG	Neat	100MG
S-13927A1-1ML	100ug/mL in Acetonitrile	1ML

Thiazopyr
CAS#117718-60-2
Chemical Name: Methyl 2-difluoromethyl-5-(4,5-dihydro-1,3-thiazol-2-yl)-4-isobutyl-6-trifluoromethylnicotinate

PS: Solid

MF: C₁₆H₁₇F₅N₂O₂S

MW: 396.4


MP: 77.3°C

Alternate Name: Visor(TM)

Part Number	Matrix	Size
N-13563-500MG	Neat	500MG
S-13563A1-1ML	100ug/ml in Acetonitrile	1ML
S-13563U1-1ML	100ug/ml in Toluene	1ML

Spiroxamine
CAS#118134-30-8
Chemical Name: 8-tert-Butyl-1,4-dioxaspiro[4,5]decan-2-ylmethyl(ethyl)(propyl)amine

PS: Liquid

MF: C₁₈H₃₅N₂O₂
MW: 297.5

FP: 147°C


Part Number	Matrix	Size
N-13226-100MG	Neat	100MG
S-13226U1-1ML	100ug/ml in Toluene	1ML

Transfluthrin CAS#118712-89-3


Chemical Name: 2,3,5,6-Tetrafluorobenzyl (+)-(1R-trans)-2,2-dimethyl-3-(2,2-dichlorovinyl)cyclopropanecarboxylate

PS: Solid

MF: C₁₅H₁₂Cl₂F₄O₂

MW: 371.2

MP: 32°C



Alternate Name: Bayothrin

Part Number	Matrix	Size
N-13626-100MG	Neat	100MG

Tebufenpyrad CAS#119168-77-3


Chemical Name: 4-Chloro-N-[[4-(1,1-dimethylethyl)phenyl]methyl]-3-ethyl-1-methyl-1H-pyrazole-5-carboxamide

PS: Solid

MF: C₁₈H₂₄ClN₃O

MW: 333.9

MP: 61°C



Alternate Name: Pyranica(TM)

Part Number	Matrix	Size
N-13502-100MG	Neat	100MG

Difenoconazole CAS#119446-68-3

Chemical Name: cis,trans-3-Chloro-4-[4-methyl-2-(1H-1,2,4-triazol-1-ylmethyl)-1,3-dioxolan-2-yl]phenyl 4-chlorophenyl ether


PS: Solid

MF: C₁₉H₁₇Cl₂N₃O₃

MW: 406.3

MP: 78.6°C

FP: 63°C



Part Number	Matrix	Size
N-11720-250MG	Neat	250MG
S-11720M1-1ML	100ug/ml in Methanol	1ML

Icaridin CAS#119515-38-7

Chemical Name: 1-Methylpropyl 2-(2-hydroxyethyl)-1-piperidinecarboxylate


PS: Liquid

MF: C₁₂H₂₃NO₃

MW: 229.32

MP: -170°C

FP: 142°C



Alternate Name: Picaridin


Part Number	Matrix	Size
N-12197-100MG	Neat	100MG
S-12197A1-1ML	100ug/mL in Acetonitrile	1ML

Emamectin CAS#119791-41-2

PS: Solid

MF: C₄₉H₇₅NO₁₃ (B₁a) C₄₈H₇₃NO₁₃ (B₁b)

MW: 886.1



Part Number	Matrix	Size
N-11847-100MG	Neat	100MG
S-11847A1-1ML	100ug/ml in Acetonitrile	1ML

Fipronil CAS#120068-37-3


Chemical Name: 5-Amino-[2,6-dichloro-4-(trifluoromethyl)phenyl]-4-[(1R,S)-(trifluoromethyl)sulfinyl]-1H-pyrazole-3-carbonitrile

PS: Solid

MF: C₁₂H₄Cl₂F₆N₄OS

MW: 437.15

MP: 200°C



Part Number	Matrix	Size
N-11971-100MG	Neat	100MG
S-11971A1-1ML	100ug/ml in Acetonitrile	1ML
S-11971U1-1ML	100ug/ml in Toluene	1ML

Cyazofamid CAS#120116-88-3


Chemical Name: 4-Chloro-2-cyano-N,N-dimethyl-5-p-tolylimidazole-1-sulfonamide

PS: Solid

MF: C₁₃H₁₃C₁N₄O₂S

MW: 324.79

MP: 152.7°C



Part Number	Matrix	Size
N-11519-100MG	Neat	100MG


Amidosulfuron CAS#120923-37-7

Chemical Name: 1-(4,6-dimethoxypyrimidin-2-yl)-3-mesyl(methyl)sulfamoylurea

PS: Solid

MF: C₉H₁₅N₅O₇S₂

MW: 369.4



Alternate Name: Gratil


Part Number	Matrix	Size
N-13819-100MG	Neat	100MG

Fenazaquin
CAS#120928-09-8

Chemical Name: 4-tert-Butylphenethyl quinazolin-4-yl ether

PS: Solid

MF: C₂₀H₂₂N₂O

MW: 306.4



Part Number	Matrix	Size
N-11949-100MG	Neat	100MG
S-11949A1-1ML	100ug/ml in Acetonitrile	1ML

Cyprodinil
CAS#121552-61-2

Chemical Name: 4-Cyclopropyl-6-methyl-N-phenylpyrimidin-2-amine

PS: Solid

MF: C₁₄H₁₅N₃
MW: 225.3

MP: 75.9°C

Alternate Name: Vangard(TM)


Part Number	Matrix	Size
N-11548-250MG	Neat	250MG
S-11548A1-1ML	100ug/mL in Acetonitrile	1ML
S-11548U1-1ML	100ug/ml in Toluene	1ML

Furilazole
CAS#121776-33-8

Chemical Name: Ethanone, 2,2-dichloro-1-[5-(2-furanyl)-2,2-dimethyl-3-oxazolidinyl]-

PS: Solid

MF: C₁₁H₁₃Cl₂NO₃
MW: 278.1

MP: 96°C



Part Number	Matrix	Size
N-13817-10MG	Neat	10MG

Cyhalofop
CAS#122008-78-0

Chemical Name: (R)-2-[4-(4-cyano-2-fluorophenoxy)phenoxy]propionic acid

PS: Solid

MF: C₁₆H₁₂FNO₄
MW: 301.27

MP: 122°C



Part Number	Matrix	Size
N-14038-25MG	Neat	25MG
S-14038A1-1ML	100ug/ml in Acetonitrile	1ML

Cyhalofop-butyl
CAS#122008-85-9

Chemical Name: Butyl (R)-2-[4-(4-cyano-2-fluorophenoxy)phenoxy]propionate

PS: Solid

MF: C₂₀H₂₀FNO₄
MW: 357.4

MP: 49.5°C


Part Number	Matrix	Size
N-11543-50MG	Neat	50MG
S-11543A1-1ML	100ug/ml in Acetonitrile	1ML
S-11543U1-1ML	100ug/ml in Acetonitrile	1ML


Chlorfenapyr
CAS#122453-73-0

Chemical Name: 4-Bromo-2-(4-chlorophenyl)-1-ethoxymethyl-5-trifluoromethylpyrrole-3-carbonitrile

PS: Solid

MF: C₁₅H₁₁BrClF₃N₂O

MW: 407.6

MP: 100°C


Part Number	Matrix	Size
N-11428-100MG	Neat	100MG
S-11428A1-1ML	100ug/ml in Acetonitrile	1ML
S-11428U1-1ML	100ug/ml in Toluene	1ML


Imazosulfuron
CAS#122548-33-8

Chemical Name: 1-(2-chloroimidazo[1,2-a]pyridin-3-ylsulfonyl)-3-(4,6-dimethoxypyrimidin-2-yl)urea

PS: Solid

MF: C₁₄H₁₃ClN₆O₅S

MW: 412.8

MP: 178.6°C


Part Number	Matrix	Size
N-13818-100MG	Neat	100MG
S-13818A1-1ML	100ug/ml in Acetonitrile	1ML


Sulfentrazone
CAS#122836-35-5

Chemical Name: N-(2,4-Dichloro-5-(4-difluoromethyl)-4,5-dihydro-3-methyl-5-oxo-1H-1,2,4-triazol-1-yl)phenyl)-Methanesulfonamide

PS: Solid

MF: C₁₁H₁₀Cl₂F₂N₄O₃S

MW: 387.2


MP: 121°C

Alternate Name: Authority*

Part Number	Matrix	Size
N-13253-250MG	Neat	250MG
S-13253M1-1ML	100ug/ml in Methanol	1ML

Rimsulfuron
CAS#122931-48-0
Chemical Name: 1-(4,6-Dimethoxypyrimidin-2-yl)-3-(3-ethylsulfonyl-2-pyridylsulfonyl)urea

PS: Solid

MF: C₁₄H₁₇N₅O₇S₂
MW: 431.48

MP: 173°C



Part Number	Matrix	Size
N-13180-100MG	Neat	100MG

Pymetrozine
CAS#123312-89-0
Chemical Name: (E)-4,5-DIHYDRO-6-METHYL-4-(3-PYRIDYLMETHYLENEAMINO)-1,2,4-TRIAZIN-3(2H)-ONE

PS: Solid

MF: C₁₀H₁₁N₅O

MW: 217.23


MP: 217°C


Part Number	Matrix	Size
N-13141-250MG	Neat	250MG
S-13141A1-1ML	100ug/ml in Acetonitrile	1ML
S-13141U1-1ML	100ug/ml in Toluene	1ML

Pyriithiobac-Sodium
CAS#123343-16-8
Chemical Name: Sodium 2-Chloro-6-(4,6-dimethoxypyrimidin-2-ylthio)benzoate

PS: Solid


MW: 348.73

MP: 233°C

Alternate Name: Staple(TM)

Part Number	Matrix	Size
N-13161-500MG	Neat	500MG
S-13161F1-1ML	100ug/ml in H2O	1ML

4-(2-Phenylethoxy)-quinazoline
CAS#124427-60-7
PS: Solid

MF: C₁₆H₁₄N₂O

MW: 250.3



Part Number	Matrix	Size
N-12908-10MG	Neat	10MG

Quinoxifen
CAS#124495-18-7
Chemical Name: 5,7-Dichloro-4-quinolyl 4-fluorophenyl ether

PS: Solid

MF: C₁₅H₈C₁₂FN₂O

MW: 308.14


MP: 106°C

Alternate Name: 5,7-Dichloro-4-(4-fluorophenoxy)quinoline

Part Number	Matrix	Size
N-13173-100MG	Neat	100MG
S-13173A1-1ML	100ug/ml in Acetonitrile	1ML
S-13173U1-1ML	100ug/ml in Toluene	1ML

Metconazole
CAS#125116-23-6
Chemical Name: 5-[(4-Chlorophenyl)methyl]-2,2-dimethyl-1-(1H-1,2,4-triazol-1-ylmethyl)cyclopentanol

PS: Solid

MF: C₁₇H₁₂ClN₃O

MW: 319.8



Part Number	Matrix	Size
N-12386-100MG	Neat	100MG

Bispyribac Sodium
CAS#125401-92-5
Chemical Name: Sodium 2,6-bis(4,6-dimethoxypyrimidin-2-yloxy)benzoate

PS: Solid

MF: C₁₉H₁₇N₄O₈Na

MW: 452.39


MP: 224°C

Alternate Name: Nominee(TM)

Part Number	Matrix	Size
N-11323-100MG	Neat	100MG

Triflusulfuron-methyl
CAS#126535-15-7
Chemical Name: 2-[4-Dimethylamino-6-(2,2,2-trifluoroethoxy)-1,3,5-triazin-2-ylcarbamoylsulfamoyl]-m-toluic acid

PS: Solid

MF: C₁₇H₁₉F₃N₆O₆S


MW: 492.4

Alternate Name: Upbeet(TM)

Part Number	Matrix	Size
N-13690-50MG	Neat	50MG
S-13690A1-1ML	100ug/ml in Acetonitrile	1ML

Ethoxysulfuron
CAS#126801-58-9
Chemical Name: 1-(4,6-dimethoxypyrimidin-2-yl)-3-(2-ethoxyphenoxysulfonyl)urea

PS: Solid

MF: C₁₅H₁₈N₄O₇S


MW: 398.4


Part Number	Matrix	Size
N-11878-100MG	Neat	100MG

Fenhexamid
CAS#126833-17-8
Chemical Name: N-(2,3-dichloro-4-hydroxyphenyl)-1-methyl-cyclohexanecarboxamide

PS: Solid


MF: C₁₄H₁₇Cl₂NO₂
MW: 302.2

MP: 153°C

Alternate Name: Elevate(TM)

Part Number	Matrix	Size
N-11954-250MG	Neat	250MG
S-11954A1-1ML	100ug/ml in Acetonitrile	1ML

Prohexadione-Calcium
CAS#127277-53-6
Chemical Name: Calcium 3-oxido-5-oxo-4-propionylcyclohex-3-enecarboxylate

PS: Solid


MF: C₁₀H₁₀CaO₅
MW: 250.3

Alternate Name: Viviful(TM)

Part Number	Matrix	Size
N-13100-100MG	Neat	100MG

Carfentrazone-ethyl
CAS#128639-02-1
Chemical Name: Ethyl (RS)-2-chloro-3-[2-chloro-5-(4-difluoromethyl-4,5-dihydro-3-methyl-5-oxo-1H-1,2,4-triazol-1-yl)-4-fluorophenyl]propionate

PS: Liquid


MF: C₁₅H₁₄Cl₂F₃N₃O₃
MW: 412.2

MP: -22.1°C

Alternate Name: Aim(TM)

Part Number	Matrix	Size
N-11411-10MG	Neat	10MG
S-11411A1-1ML	100ug/ml in Acetonitrile	1ML
S-11411U1-1ML	100ug/ml in Toluene	1ML

Tolfenpyrad
CAS#129558-76-5
Chemical Name: 4-Chloro-3-ethyl-1-methyl-N-[[4-(4-methylphenoxy)phenyl]methyl]-1H-pyrazole-5-carboxamide

PS: Solid


MF: C₂₁H₂₂ClN₃O₂
MW: 383.9


Part Number	Matrix	Size
N-13579-100MG	Neat	100MG

Pyraflufen-ethyl
CAS#129630-19-9
Chemical Name: Ethyl [2-chloro-5-(4-chloro-5-difluoromethoxy-1-methylpyrazol-3-yl)-4-fluorophenoxy]acetate

PS: Solid

MF: C₁₅H₁₃Cl₂F₃N₂O₄
MW: 413.2


MP: 126°C


Part Number	Matrix	Size
N-13144-100MG	Neat	100MG

Amicarbazone
CAS#129909-90-6
Chemical Name: 4-amino-N-(1,1-dimethylethyl)-4,5-dihydro-3-(1-methylethyl)-5-oxo-1H-1,2,4-triazole-1-carboxamide

PS: Solid

MF: C₁₀H₁₉N₅O₂
MW: 241.3

MP: 137.5°C


Part Number	Matrix	Size
N-13075-100MG	Neat	100MG
S-13075A1-1ML	100ug/mL in Acetonitrile	1ML

Thifluzamide
CAS#130000-40-7
Chemical Name: 2',6'-dibromo-2-methyl-4'-trifluoromethoxy-4-trifluoromethyl-1,3-thiazole-5-carboxanilide

PS: Solid

MF: C₁₃H₆Br₂F₆N₂O₂S

MW: 528.1


MP: 177.9-178.6°C


Part Number	Matrix	Size
N-12884-100MG	Neat	100MG
S-12884M1-1ML	100ug/ml in Methanol	1ML

Fludioxonil CAS#131341-86-1
Chemical Name: 4-(2,2-Difluoro-1,3-benzodioxol-4-yl)-1H-pyrrole-3-carbonitrile

PS: Solid

MF: C₁₂H₆F₂N₂O₂
MW: 248.2


MP: 199.8°C

Alternate Name: Celest(TM)

Part Number	Matrix	Size
N-11985-100MG	Neat	100MG
S-11985A1-1ML	100ug/ml in Acetonitrile	1ML
S-11985U1-1ML	100ug/ml in Toluene	1ML

Famoxadone CAS#131807-57-3
Chemical Name: 3-Anilino-5-methyl-5-(4-phenoxyphenyl)-1,3-oxazolidine-2,4-dione

PS: Solid


MF: C₂₂H₁₈N₂O₄
MW: 374.4

MP: 140.3°C


Part Number	Matrix	Size
N-11943-100MG	Neat	100MG
S-11943A1-1ML	100ug/ml in Acetonitrile	1ML
S-11943U1-1ML	100ug/ml in Toluene	1ML

Azoxystrobin CAS#131860-33-8
Chemical Name: (E)-2-[6-(2-cyanophenoxy)pyrimidin-4-yloxy]phenyl]-3-methoxyacrylate

PS: Solid

MF: C₂₂H₁₇N₃O₅
MW: 403.4

Alternate Name: Heritage(TM)


Part Number	Matrix	Size
N-11113-100MG	Neat	100MG
S-11113A1-1ML	100ug/ml in Acetonitrile	1ML
S-11113U1-1ML	100ug/ml in Toluene	1ML

Triticonazole CAS#131983-72-7
Chemical Name: (±)-(E)-5-(4-Chlorobenzylidene)-2,2-dimethyl-1-(1H-1,2,4-triazol-1-ylmethyl)cyclopentanol

PS: Solid

MF: C₁₇H₂₀ClN₃O

MW: 317.82


MP: 140°C


Part Number	Matrix	Size
N-13729-100MG	Neat	100MG

Indanofan CAS#133220-30-1
Chemical Name: 1H-Indene-1,3(2H)-dione, 2-[[2-(3-chlorophenyl)-2-oxiranyl]methyl]-2-ethyl-

PS: Solid

MF: C₂₀H₁₇ClO₃
MW: 340.8

MP: 60°C

Alternate Name: Epoprodan


Part Number	Matrix	Size
N-12209-100MG	Neat	100MG

Epoxiconazole CAS#133855-98-8
Chemical Name: (2R,3R)-1-[3-(2-Chlorophenyl)-2,3-epoxy-2-(4-fluorophenyl)propyl]-1H-1,2,4-triazole

PS: Solid

MF: C₁₇H₁₃ClFN₃O


MW: 329.8

MP: 136-137°C


Part Number	Matrix	Size
N-11859-100MG	Neat	100MG

Fenproxiimate CAS#134098-61-6
Chemical Name: 1,1-dimethylethyl 4-[[[(E)-[(1,3-dimethyl-5-phenoxy-1H-pyrazol-4-yl)methylene]amino]oxy]methyl]benzoate

PS: Solid


MF: C₂₄H₂₇N₃O₄
MW: 421.5

Alternate Name: Acaban(TM)

Part Number	Matrix	Size
N-11962-100MG	Neat	100MG
S-11962A1-1ML	100ug/ml in Acetonitrile	1ML
S-11962U1-1ML	100ug/ml in Toluene	1ML

Acibenzolar-S-methyl CAS#135158-54-2
Chemical Name: Benzo(1,2,3)thiadiazole-7-carbothioic acid-S-methyl ester

PS: Solid

MF: C₈H₆N₂OS₂
MW: 210.28


MP: 132.9°C


Part Number	Matrix	Size
N-11026-250MG	Neat	250MG
S-11026A1-1ML	100ug/ml in Acetonitrile	1ML
S-11026U1-1ML	100ug/ml in Toluene	1ML

Acetamiprid
CAS#135410-20-7
Chemical Name: (E)-N-[(6-Chloro-3-pyridinyl)methyl]-N'-cyano-N-methylethanamide

PS: Solid

MF: C₁₀H₁₁ClN₄
MW: 222.7


MP: 98.9°C

Alternate Name: Mospilan(TM)

Part Number	Matrix	Size
N-11009-100MG	Neat	100MG
S-11009A1-1ML	100ug/ml in Acetonitrile	1ML
S-11009U1-1ML	100ug/mL in Toluene	1ML

Mefenpyr-diethyl
CAS#135590-91-9
Chemical Name: Diethyl (RS)-1-(2,4-dichlorophenyl)-5-methyl-2-pyrazoline-3,5-dicarboxylate

PS: Solid

MF: C₁₆H₁₈Cl₂N₂O₄
MW: 373.2

MP: 50°C



Part Number	Matrix	Size
N-12366-100MG	Neat	100MG
S-12366A1-1ML	100ug/ml in Acetonitrile	1ML
S-12366U1-1ML	100ug/ml in Toluene	1ML

Cyclosulfamuron
CAS#136849-15-5
Chemical Name: 1-[2-(Cyclopropylcarbonyl)phenylsulfamoyl]-3-(4,6-dimethoxyphenyl)urea

PS: Solid


MF: C₁₇H₁₉N₅O₆S

MW: 421.4

MP: 160.9°C


Part Number	Matrix	Size
N-11542-250MG	Neat	250MG
S-11542A1-1ML	100ug/ml in Acetonitrile	1ML
S-11542U1-1ML	100ug/ml in Toluene	1ML

Picolinafen solution
CAS#137641-05-5
PS: Liquid


MF: C₁₉H₁₂F₄N₂O₂


Part Number	Matrix	Size
S-13982M1-1ML	100ug/ml in Methanol	1ML

Imidacloprid
CAS#138261-41-3
Chemical Name: 1-(6-Chloro-3-pyridylmethyl)-N-nitroimidazolidin-2-ylideneamine

PS: Solid

MF: C₉H₁₀ClN₅O₂
MW: 255.7


MP: 144°C


Part Number	Matrix	Size
N-12206-100MG	Neat	100MG
S-12206M1-1ML	100ug/ml in Methanol	1ML

Metaflumizone
CAS#139968-49-3
Chemical Name: 2-(2-(4-Cyanophenyl)-1-(3-(trifluoromethyl)phenyl)ethylidene)-N-(4-(trifluoromethoxy)phenyl)-hydrazinecarboxamide

PS: Solid


MF: C₂₄H₁₆F₆N₄O₂
MW: 506.4

MP: 197°C


Part Number	Matrix	Size
N-12379-100MG	Neat	100MG

Iprovalicarb
CAS#140923-17-7
Chemical Name: isopropyl 2-methyl-1-[(1-p-tolylethyl)carbamoyl]-(S)-propylcarbamate

PS: Liquid

MF: C₁₈H₂₈N₂O₃
MW: 320.4



Part Number	Matrix	Size
N-13924-10MG	Neat	10MG

Isoxaflutole
CAS#141112-29-0
Chemical Name: 5-Cyclopropyl-1,2-oxazol-4-yl a,a,a-trifluoro-2-mesyl-p-tolyl ketone

PS: Solid

MF: C₁₅H₁₂F₃NO₄S

MW: 359.3

MP: 138°C

Alternate Name: Balance

Part Number	Matrix	Size
N-12287-100MG	Neat	100MG
S-12287U1-1ML	100ug/ml in Toluene	1ML

Trifloxystrobin CAS#141517-21-7


Chemical Name: N-[[[4,6-dimethoxy-2-pyrimidinyl]amino]carbonyl]-3-(2,2,2-trifluoroethoxy)-2-pyridinesulfonamide

PS: Solid

MF: C₂₀H₁₉F₃N₂O₄

MW: 408.4

MP: 72.9°C



Alternate Name: Flint(TM)

Part Number	Matrix	Size
N-13684-100MG	Neat	100MG
S-13684A1-1ML	100ug/ml in Acetonitrile	1ML

Sulfosulfuron CAS#141776-32-1


Chemical Name: 1-(2-Ethylsulfonylimidazo[1,2-a]pyridin-3-ylsulfonyl)-3-(4,6-dimethoxy-2-pyrimidin-2-yl)urea

PS: Solid

MF: C₁₆H₁₈N₆O₇S₂

MW: 470.48

MP: 201°C



Alternate Name: Maverick(TM)

Part Number	Matrix	Size
N-13255-100MG	Neat	100MG
S-13255A1-1ML	100ug/ml in Acetonitrile	1ML
S-13255U1-1ML	100ug/ml in Toluene	1ML

Flufenacet CAS#142459-58-3


Chemical Name: N-(4-Fluorophenyl)-N-(1-methylethyl)-2-((5-(trifluoromethyl)-1,3,4-thiadiazol-2-yl)oxy)acetamide

PS: Solid

MF: C₁₄H₁₃F₄N₃O₂S

MW: 364.34

MP: 75°C



Alternate Name: Fluthiamide

Part Number	Matrix	Size
N-11986-100MG	Neat	100MG
S-11986M1-1ML	100ug/ml in Methanol	1ML

Kresoxim-methyl CAS#143390-89-0


Chemical Name: Methyl(E)-2-methoxyimino-[2-(o-tolylloxymethyl)phenyl]acetate

PS: Solid

MF: C₁₈H₁₉NO₄

MW: 313.4

MP: 97.2°C



Alternate Name: Strobry(TM)


Part Number	Matrix	Size
N-12296-250MG	Neat	250MG
S-12296A1-1ML	100ug/ml in Acetonitrile	1ML
S-12296U1-1ML	100ug/ml in Toluene	1ML

Chromafenozide solution CAS#143807-66-3

PS: Liquid

MF: C₂₄H₃₀N₂O₃

MW: 394.5



Part Number	Matrix	Size
S-13973A1-1ML	100ug/ml in Acetonitrile	1ML


Iodosulfuron-methyl-sodium CAS#144550-36-7

Chemical Name: methyl 4-iodo-2-[3-(4-methoxy-6-methyl-1,3,5-triazin-2-yl)ureido-sulfonyl]benzoate, sodium salt

PS: Solid

MF: C₁₄H₁₃IN₃NaO₆S

MW: 529.2



Part Number	Matrix	Size
N-13923-10MG	Neat	10MG

Diclosulam CAS#145701-21-9


Chemical Name: N-(2,6-Dichlorophenyl)-5-ethoxy-7-fluoro-[1,2,4]-triazolo-[1,5-c]pyrimidine-2-sulfonamide

PS: Solid

MF: C₁₃H₁₀C₁₂FN₅O₃S

MW: 406.2

MP: 218°C



Alternate Name: Strongarm(TM)

Part Number	Matrix	Size
N-11679-250MG	Neat	250MG
S-11679A1-1ML	100ug/ml in Acetonitrile	1ML


Florasulam CAS#145701-23-1

Chemical Name: 2',6'-trifluoro-5-methoxy[1,2,4]triazolo[1,5-c]pyrimidine-2-sulfonamide

PS: Liquid

MF: C₁₂H₈F₃N₅O₃S

MW: 359.3




Part Number	Matrix	Size
N-13922-10MG	Neat	10MG

Cloransulam-methyl
CAS#147150-35-4
Chemical Name: Methyl 3-chloro-2-(5-ethoxy-7-fluoro[1,2,4]triazolo[1,5-c]-pyrimidin-2-ylsulfonamido)benzoate

PS: Solid

MF: C₁₅H₁₃O₅N₅SClF


MW: 429.8

MP: 216°C

Alternate Name: FirstRate

Part Number	Matrix	Size
N-11492-100MG	Neat	100MG
S-11492A1-1ML	100ug/ml in Acetonitrile	1ML
S-11492U1-1ML	100ug/ml in Toluene	1ML

Spirodiclofen
CAS#148477-71-8
Chemical Name: 3-(2,4-Dichlorophenyl)-2-oxo-1-oxaspiro[4,5]dec-3-en-4-yl 2,2-dimethylbutanoate


PS: Solid

MF: C₂₁H₂₄C₁₂O₄
MW: 411.3

Alternate Name: Daniemon(TM)

Part Number	Matrix	Size
N-13223-100MG	Neat	100MG

Simeconazole solution
CAS#149508-90-7
PS: Liquid

MF: C₁₄H₂₀FN₃OSi


MW: 293.4


Part Number	Matrix	Size
S-13931A1-1ML	100ug/ml in Acetonitrile	1ML

Bifenazate
CAS#149877-41-8
Chemical Name: 1-Methylethyl 2-(4-methoxy[1,1'-biphenyl]-3-yl)hydrazinecarboxylate

PS: Solid


MF: C₁₇H₂₀N₂O₃
MW: 300.35

MP: 120°C

Alternate Name: Floramite(TM)

Part Number	Matrix	Size
N-11201-50MG	Neat	50MG
S-11201A1-1ML	100ug/mL in Acetonitrile	1ML
S-11201U1-1ML	100ug/ml in Toluene	1ML


Dimoxystrobin Solution
CAS#149961-52-4
Chemical Name: (E)-2-(2,5-Dimethylphenoxyethyl)-?methoxyimino-N-methyl-phenylacetamide

PS: Liquid

MF: C₁₉H₂₂N₂O₃
MW: 326.4


Part Number	Matrix	Size
S-13926A1-1ML	100ug/ml in Acetonitrile	1ML


Caloxydim Solution
CAS#149979-41-9
PS: Liquid

MF: C₁₇H₂₄ClNO₄
MW: 341.8


Part Number	Matrix	Size
S-13970A1-1ML	100ug/ml in Acetonitrile	1ML

Aminopyralid
CAS#150114-71-9
Chemical Name: 4-Amino-3,6-dichloropyridine-2-carboxylic acid

PS: Liquid

MF: C₆H₄Cl₂N₂O₂
MW: 207



Part Number	Matrix	Size
N-11066-100MG	Neat	100MG

Nitenpyram
CAS#150824-47-8
Chemical Name: 1,1-Ethenediamine, N-[(6-chloro-3-pyridinyl)methyl]-N-ethyl-N'-methyl-2-nitro-, (1E)-

PS: Solid

MF: C₁₁H₁₅ClN₄O₂
MW: 270.75

MP: 82°C

FP: 70°C

Alternate Name: Capstar

Part Number	Matrix	Size
N-12657-100MG	Neat	100MG
S-12657M1-1ML	100ug/ml in Methanol	1ML

sec-Acetochlor OA Solution
CAS#152019-74-4
Chemical Name: S-Metolachlor Metabolite CGA 50720


PS: Liquid

Part Number	Matrix	Size
S-13946M0-1ML	50ug/mL in Methanol	1ML

Etiozazole
CAS#153233-91-1
Chemical Name: (RS)-5-tert-Butyl-2-[2-(2,6-difluorophenyl)-4,5-dihydro-1,3-oxazol-4-yl]phenetole

PS: Solid

MF: C₂₁H₂₃F₂NO₂
MW: 359.4

MP: 101°C



Part Number	Matrix	Size
N-11942-50MG	Neat	50MG
S-11942A1-1ML	100ug/ml in Acetonitrile	1ML
S-11942U1-1ML	100ug/ml in Toluene	1ML

Thiamethoxam
CAS#153719-23-4
Chemical Name: 3-(2-Chloro-1,3-thiazol-5-ylmethyl)-1,3,5-oxadiazinan-4-ylidene(nitro)amine

PS: Solid


MF: C₈H₁₀ClN₅O₃S

MW: 291.7

MP: 139.1°C

Alternate Name: Actara(TM)

Part Number	Matrix	Size
N-13562-100MG	Neat	100MG
S-13562M1-1ML	100ug/ml in Methanol	1ML

Emamectin benzoate
CAS#155569-91-8
PS: Solid


MP: 141-146°C


Part Number	Matrix	Size
N-12920-100MG	Neat	100MG

Zoxamide
CAS#156052-68-5
Chemical Name: (RS)-3,5-Dichloro-N-(3-chloro-1-ethyl-1-methyl-2-oxopropyl)-p-toluidide/3,5-Dichloro-N-(3-chloro-1-ethyl-1-methyl-2-oxopropyl)-4-methyl-benzamide

PS: Solid

MF: C₁₄H₁₆Cl₃NO₂
MW: 336.6

MP: 159.5°C



Part Number	Matrix	Size
N-13762-100MG	Neat	100MG

Flonicamid
CAS#158062-67-0
Chemical Name: 4-Trifluoromethylnicotinamide

PS: Solid

MF: C₉H₆F₃N₃O

MW: 229.2


MP: 156°C


Part Number	Matrix	Size
N-11976-50MG	Neat	50MG
S-11976A1-1ML	100ug/ml in Acetonitrile	1ML
S-11976U1-1ML	100ug/ml in Toluene	1ML

Fentrazamide
CAS#158237-07-1
Chemical Name: 4-(2-chlorophenyl)-N-cyclohexyl-N-ethyl-4,5-dihydro-5-oxo-1H-tetrazole-1-carboxamide

PS: Solid


MF: C₁₆H₂₀ClN₅O₂
MW: 349.8

MP: 79°C


Part Number	Matrix	Size
N-14059-100MG	Neat	100MG
S-14059A1-1ML	100ug/ml in Acetonitrile	1ML

Methoxyfenozide
CAS#161050-58-4
Chemical Name: 3-Methoxy-2-methylbenzoic acid 2-(3,5-dimethylbenzoyl)-2-(1,1-dimethylethyl)hydrazide

PS: Solid


MF: C₂₂H₂₈N₂O₃
MW: 368.5


Part Number	Matrix	Size
N-12406-100MG	Neat	100MG
S-12406A1-1ML	100ug/ml in Acetonitrile	1ML

Fenamidone
CAS#161326-34-7
Chemical Name: (S)-1-Anilino-4-methyl-2-methylthio-4-phenylimidazolin-5-one

PS: Solid

MF: C₁₇H₁₇N₃OS


MW: 311.4


Part Number	Matrix	Size
N-11946-100MG	Neat	100MG
S-11946A1-1ML	100ug/ml in Acetonitrile	1ML
S-11946U1-1ML	100ug/ml in Toluene	1ML


Dimethenamid-P
CAS#163515-14-8
Chemical Name: (S)-2-Chloro-N-(2,4-dimethyl-3-thienyl)-N-(2-methoxy-1-methyl-ethyl)acetamide

PS: Liquid

MF: C₁₂H₁₈ClNO₂S


MW: 275.8

FP: 79°C

Alternate Name: Outlook(TM)

Part Number	Matrix	Size
N-11756-100MG	Neat	100MG
S-11756A1-1ML	100ug/ml in Acetonitrile	1ML
S-11756U1-1ML	100ug/ml in Toluene	1ML

Isoxadifen-ethyl
CAS#163520-33-0
Chemical Name: Ethyl 4,5-dihydro-5,5-diphenyl-1,2-oxazole-3-carboxylate

PS: Liquid


MF: C₁₈H₁₇NO₃
MW: 295.3


Part Number	Matrix	Size
N-12286-100MG	Neat	100MG

Dinotefuran
CAS#165252-70-0
Chemical Name: (RS)-1-methyl-2-nitro-3-(tetrahydro-3-furylmethyl)guanidine

PS: Solid


MF: C₇H₁₄N₄O₃
MW: 202.2

MP: 107.5°C


Part Number	Matrix	Size
N-11788-50MG	Neat	50MG
S-11788M1-1ML	100ug/ml in Methanol	1ML
S-11788U1-1ML	100ug/ml in Toluene	1ML


Pyribenzoxim
CAS#168088-61-7
Chemical Name: Benzophenone O-[2,6-bis(4,6-dimethoxypyrimidin-2-yl)oxy]benzoyl]oxime

PS: Solid

MF: C₃₂H₂₇N₅O₈
MW: 609.6

Alternate Name: Diphenylmethanone O-[2,6-bis[(4,6-dimethoxy-2-pyrimidinyl)oxy]benzoyl]oxime

Part Number	Matrix	Size
N-14070-100MG	Neat	100MG

Spinosad
CAS#168316-95-8
PS: Solid


MF: C₄₁H₆₅NO₁₀
MW: 732

Alternate Name: Spinosyn A, 131929-60-7; Spinosyn D, 131929-63-0

Part Number	Matrix	Size
N-13222-100MG	Neat	100MG
S-13222A1-1ML	100ug/ml in Acetonitrile	1ML

Foramsulfuron
CAS#173159-57-4
Chemical Name: 2-[[[(4,6-dimethoxy-2-pyrimidinyl)amino]carbonyl]amino]sulfonyl]-4-(formylamino)-N,N-dimethylbenzamide

PS: Solid

MF: C₁₇H₂₀N₆O₇S


MP: 126°C


Part Number	Matrix	Size
N-12009-100MG	Neat	100MG
S-12009A1-1ML	100ug/ml in Acetonitrile	1ML

S-Indoxacarb
CAS#173584-44-6
Chemical Name: Methyl (S)-N-[7-chloro-2,3,4a,5-tetrahydro-4a-(methoxycarbonyl)indeno[1,2-e]-[1,3,4]oxadiazin-2-ylcarbonyl]-4'-(trifluoromethoxy)carbanilate

PS: Solid

MF: C₂₂H₁₇ClF₃N₃O₇
MW: 527.8


MP: 88.1°C


Part Number	Matrix	Size
N-12213-100MG	Neat	100MG
S-12213A1-1ML	100ug/ml in Acetonitrile	1ML
S-12213U1-1ML	100ug/ml in Toluene	1ML

Pyraclostrobin
CAS#175013-18-0
Chemical Name: Methyl N-[2-[1-(4-chlorophenyl)-1H-pyrazol-3-yloxy]methyl]phenyl (N-methoxy)carbamate

PS: Solid

MF: C₁₉H₁₈ClN₃O₄
MW: 387.8

MP: 63.7°C

Alternate Name: Cabrio EG(TM)


Part Number	Matrix	Size
N-13143-100MG	Neat	100MG
S-13143A1-1ML	100ug/ml in Acetonitrile	1ML

Benthiavalicarb-isopropyl
CAS#177406-68-7
Chemical Name: Isopropyl [(S)-1-[[[(1R)-1-(6-fluoro-1,3-benzothiazol-2-yl)ethyl]carbamoyl]-2-methylpropyl]carbamate

PS: Solid

MF: C₁₈H₂₄FN₃O₃S

MW: 381.5


MP: 153.1°C


Part Number	Matrix	Size
N-14076-50MG	Neat	50MG

Prothioconazole
CAS#178928-70-6
Chemical Name: 2-[(2RS)-2-(1-Chlorocyclopropyl)-3-(2-chlorophenyl)-2-hydroxypropyl]-2H-1,2,4-triazole-3(4H)-thione

PS: Solid


MF: C₁₄H₁₅Cl₂N₃OS

MW: 344.3


Part Number	Matrix	Size
N-13138-100MG	Neat	100MG
S-13138B1-1ML	100ug/mL in Acetone	1ML

Pyridalyl
CAS#179101-81-6
Chemical Name: 2,6-Dichloro-4-(3,3-dichloroallyloxy)phenyl 3-[5-(trifluoromethyl)-2-pyridyloxy]propyl ether

PS: Liquid


MF: C₁₈H₁₄Cl₄F₃NO₃
MW: 491.1

Alternate Name: Overture(TM)

Part Number	Matrix	Size
N-13153-50MG	Neat	50MG

Cyflufenamid
CAS#180409-60-3
Chemical Name: (Z)-N-[2-(cyclopropylmethoxyimino)-2,3-difluoro-6-(trifluoromethyl)benzyl]-2-phenylacetamide

PS: Solid

MF: C₂₀H₁₇F₅N₂O₂
MW: 412.4


MP: 61.5-62.5°C


Part Number	Matrix	Size
N-13132-25MG	Neat	25MG

Propoxycarbazone-sodium
CAS#181274-15-7
Chemical Name: sodium {[2-(methoxycarbonyl)phenyl]sulfonyl}[(4,5-dihydro-4-methyl-5-oxo-3-propoxy-1H-1,2,4-triazol-1-yl)carbonyl]azanide

PS: Solid

MF: C₁₅H₁₇N₄NaO₇S

MW: 420.4



Part Number	Matrix	Size
N-13119-100MG	Neat	100MG
S-13119A1-1ML	100ug/ml in Acetonitrile	1ML

Flucarbazone-sodium
CAS#181274-17-9
Chemical Name: 4,5-Dihydro-3-methoxy-4-methyl-5-oxo-N-(2-trifluoromethoxyphenylsulfonyl)-1H-1,2,4-triazole-1-carboxamide sodium salt

PS: Solid

MF: C₁₂H₁₀F₃N₄NaO₆S

MW: 418.3

MP: 200°C



Part Number	Matrix	Size
N-11982-250MG	Neat	250MG
S-11982F1-1ML	100ug/ml in H2O	1ML

Ethiprole CAS#181587-01-9

Chemical Name: 5-amino-1-[2,6-dichloro-4-(trifluoromethyl)phenyl]-4-(ethyl-sulfinyl)-1H-pyrazole-3-carbonitrile

PS: Solid

MF: C₁₃H₉Cl₂F₃N₄OS

MW: 397.2



Part Number	Matrix	Size
N-12870-100MG	Neat	100MG

Penthiopyrad CAS#183675-82-3

Chemical Name: N-[2-(1,3-dimethylbutyl)-3-thienyl]-1-methyl-3-(trifluoromethyl)-1H-pyrazole-4-carboxamide

PS: Solid

MF: C₁₆H₂₀F₃N₃OS


MW: 359.4


Part Number	Matrix	Size
N-17917-50MG	Neat	50MG

Spinetoram CAS#187166-40-1

PS: Solid

MF: C₄₂H₆₉NO₁₀
MW: 748

MP: 143.4°C


Part Number	Matrix	Size
N-13221-100MG	Neat	100MG


Boscalid CAS#188425-85-6

Chemical Name: 2-Chloro-N-(4'-chlorobiphenyl-2-yl)nicotinamide

PS: Solid

MF: C₁₈H₁₂Cl₂N₂O

MW: 343.2


MP: 142.8°C

Alternate Name: Nicobifen

Part Number	Matrix	Size
N-11328-100MG	Neat	100MG
S-11328A1-1ML	100ug/mL in Acetonitrile	1ML

Proquinazid CAS#189278-12-4

Chemical Name: 6-Iodo-2-propoxy-3-propyl-4(3H)quinazolinone

PS: Solid


MF: C₁₄H₁₇IN₂O₂
MW: 372.2


Part Number	Matrix	Size
N-14080-50MG	Neat	50MG

Picloram CAS#1918-02-1

Chemical Name: 4-Amino-3,5,6-trichloropropionic acid

PS: Solid

MF: C₆H₃Cl₃N₂O₂
MW: 241.48



Part Number	Matrix	Size
N-13050-250MG	Neat	250MG
S-13050A1-1ML	100ug/ml in Acetonitrile	1ML
S-13050A4-1ML	1000ug/ml in Acetonitrile	1ML
S-13050A4-5ML	1000ug/ml in Acetonitrile	5ML
S-13050B1-1ML	100ug/ml in Acetone	1ML
S-13050B1-5ML	100ug/ml in Acetone	5ML

Terbutol CAS#1918-11-2

Chemical Name: 2,6-bis(1,1-dimethylethyl)-4-methylphenyl N-methylcarbamate

PS: Solid

MF: C₁₇H₂₇NO₂
MW: 277.45

MP: 200°C

Alternate Name: Azak


Part Number	Matrix	Size
N-13513-100MG	Neat	100MG

Trifloxysulfuron-Sodium CAS#199119-58-9

Chemical Name: sodium (4,6-dimethoxypyrimidin-2-yl)[[3-(2,2,2-trifluoroethoxy)pyridin-2-yl]sulfonyl]carbamoyl]azanide

PS: Solid

MF: C₁₄H₁₃F₃N₅NaO₆S


MP: 170°C

Alternate Name: Envoke(TM)

Part Number	Matrix	Size
N-13685-100MG	Neat	100MG

Spirotetramat
CAS#203313-25-1
Chemical Name: cis-4-(Ethoxycarbonyloxy)-8-methoxy-3-(2,5-xylyl)-1-azaspiro[4,5]dec-3-en-2-one


PS: Solid

MF: C₂₁H₂₇NO₅
MW: 373.45

MP: 142°C



Part Number	Matrix	Size
N-13225-100MG	Neat	100MG
S-13225B1-1ML	100ug/mL in Acetone	1ML


Fipronil desulfinyl Solution
CAS#205650-65-3
Chemical Name: 5-Amino-1-[2,6-dichloro-4-(trifluoromethyl)phenyl]-4-(trifluoromethyl)-1H-pyrazole-3-carbonitrile

PS: Liquid


Part Number	Matrix	Size
S-13920B1-1ML	100ug/mL in Acetone	1ML

Sodium dimethyldithiocarbamate hydrate
CAS#207233-95-2
Chemical Name: Carbamodithioic acid, N,N-dimethyl-, sodium salt, hydrate


PS: Solid


Part Number	Matrix	Size
N-12863-100MG	Neat	100MG

Mesosulfuron-methyl
CAS#208465-21-8
Chemical Name: 2-[[4,6-Dimethoxypyrimidin-2-ylcarbamoyl]sulfamoyl]-a-(methanesulfonamido)-p-toluate

PS: Solid

MF: C₁₇H₂₁N₅O₉S₂
MW: 503.5


MP: 195.4°C


Part Number	Matrix	Size
N-12376-100MG	Neat	100MG
S-12376A1-1ML	100ug/ml in Acetonitrile	1ML

Topramezone
CAS#210631-68-8
Chemical Name: [3-(4,5-Dihydro-1,2-oxazol-3-yl)-4-mesyl-o-tolyl](5-hydroxy-1-methyl-1H-pyrazol-4-yl)methanone

PS: Solid

MF: C₁₆H₁₇N₃O₅S

MW: 363.4

Alternate Name: Impact(TM)


Part Number	Matrix	Size
N-13585-100MG	Neat	100MG

Clothianidin
CAS#210880-92-5
Chemical Name: (E)-1-(2-Chloro-1,3-thiazol-5-ylmethyl)-3-methyl-2-nitroguanidine

PS: Solid

MF: C₆H₈ClN₅O₂S


MW: 249.68

MP: 176.8°C


Part Number	Matrix	Size
N-11493-100MG	Neat	100MG
S-11493A1-1ML	100ug/ml in Acetonitrile	1ML

Flumorph
CAS#211867-47-9
Chemical Name: 3-(3,4-Dimethoxyphenyl)-3-(4-fluorophenyl)-1-(4-morpholinyl)-2-propen-1-one


PS: Solid

MF: C₂₁H₂₂FNO₄
MW: 371.4


Part Number	Matrix	Size
N-13280-25MG	Neat	25MG
S-13280A1-1ML	100ug/mL in Acetonitrile	1ML

Orthosulfamuron solution
CAS#213464-77-8
PS: Liquid

MF: C₁₆H₂₀N₆O₆S


MW: 424.4


Part Number	Matrix	Size
S-13944B1-1ML	100ug/ml in Benzene	1ML

Lenacil
CAS#2164-08-100
Chemical Name: 3-cyclohexyl-1,5,6,7-tetrahydrocyclopentapyrimidine-2,4(3H)-di-one

PS: Solid

MF: C₁₃H₁₈N₂O₂
MW: 234.33

MP: 315.6°C



Part Number	Matrix	Size
N-12314-100MG	Neat	100MG
S-12314M1-1ML	100ug/ml in Methanol	1ML

Penoxsulam
CAS#219714-96-2
Chemical Name: 3-(2,2-Difluoroethoxy)-N-(5,8-dimethoxy[1,2,4]triazolo[1,5-c]pyrimidin-2-yl)-a,a,a-trifluorotoluene-2-sulfonamide

PS: Solid

MF: C₁₆H₁₄F₅N₅O₅S


MW: 483.4

MP: 212°C

Alternate Name: Clipper

Part Number	Matrix	Size
N-12823-50MG	Neat	50MG

Metrafenone
CAS#220899-03-6
Chemical Name: (3-bromo-6-methoxy-2-methylphenyl)(2,3,4-trimethoxy-6-methylphenyl)methanone

PS: Solid


MF: C₁₉H₂₁BrO₅
MW: 409.3


Part Number	Matrix	Size
N-12897-100MG	Neat	100MG

Cyprosulfamide
CAS#221667-31-8
Chemical Name: N-[4-(Cyclopropylcarbamoyl)phenylsulfonyl]-o-anisamide

PS: Solid

MF: C₁₈H₁₈N₂O₅S

MW: 374.41


Part Number	Matrix	Size
N-11549-100MG	Neat	100MG

10% OFF


10 or more neat standards of the same item

Tiadinil
CAS#223580-51-6
Chemical Name: N-(3-Chloro-4-methylphenyl)-4-methyl-1,2,3-thiadiazole-5-carboxamide

PS: Solid


MF: C₁₁H₁₀ClN₃OS

MW: 267.7

MP: 112.2°C


Part Number	Matrix	Size
N-13282-25MG	Neat	25MG
S-13282A1-1ML	100ug/mL in Acetonitrile	1ML

Fluacrypyrim solution
CAS#229977-93-9
PS: Liquid

MF: C₂₀H₂₁F₃N₂O₅



Part Number	Matrix	Size
S-13980M1-1ML	100ug/ml in Methanol	1ML

Flurenol-butyl CAS#2314-09-2

Chemical Name: butyl 9-hydroxy-9H-fluorene-9-carboxylate

PS: Solid

MF: C₁₈H₁₈O₃
MW: 282.3

MP: 71°C

Alternate Name: Butyl flurenol

Part Number	Matrix	Size
N-13814-250MG	Neat	250MG


Fluopicolide CAS#239110-15-7

Chemical Name: 2,6-Dichloro-N-[[3-chloro-5-(trifluoromethyl)-2-pyridinyl]methyl]benzamide

PS: Solid

MF: C₁₄H₈C₁₃F₃N₂O

MW: 383.35

MP: 149°C


Part Number	Matrix	Size
N-11993-100MG	Neat	100MG


Captafol CAS#2425-06-1

Chemical Name: N-[1,1,2,2-Tetrachloroethylthio]-4-cyclohexene-1,2-dicarboximide

PS: Solid

MF: C₁₀H₉Cl₄NO₂S

MW: 349.1


MP: 160°C


Part Number	Matrix	Size
N-11399-250MG	Neat	250MG

MPMC CAS#2425-10-7

Chemical Name: 3,4-dimethylphenyl N-methylcarbamate

PS: Solid

MF: C₁₀H₁₃NO₂
MW: 179.2

Alternate Name: XYLYLCARB


Part Number	Matrix	Size
N-13274-10MG	Neat	10MG

Quinomethionate CAS#2439-01-2

Chemical Name: S,S-(6-Methylquinoxaline-2,3-diyl) dithiocarbonate

PS: Solid

MF: C₁₀H₆N₂OS₂
MW: 234.29

MP: 171°C

Alternate Name: Chinomethionate


Part Number	Matrix	Size
N-13172-100MG	Neat	100MG

Dodine (Technical) CAS#2439-10-3

Chemical Name: n-Dodecylguanidin acetate

PS: Solid

MF: C₁₃H₂₉N₃C₂H₄O₂
MW: 287.45

MP: 136°C



Part Number	Matrix	Size
N-11840-250MG	Neat	250MG

Pinoxaden CAS#243973-20-8

Chemical Name: 8-(2,6-Diethyl-p-tolyl)-1,2,4,5-tetrahydro-7-oxo-7H-pyrazolo[1,2-d][1,4,5]oxadiazepin-9-yl 2,2-dimethylpropionate

PS: Solid

MF: C₂₃H₃₂N₂O₄
MW: 400.5

MP: 120°C

Alternate Name: Axial(TM)


Part Number	Matrix	Size
N-13057-100MG	Neat	100MG

O,O-Diethyl phosphorochloridothioate CAS#2524-04-1

PS: Liquid

MF: C₄H₁₀ClO₂PS


MW: 188.62

FP: 112°C


Part Number	Matrix	Size
N-12703-250MG	Neat	250MG
S-12703A1-1ML	100ug/ml in Acetonitrile	1ML
S-12703U1-1ML	100ug/ml Toluene	1ML

Metamifop
CAS#256412-89-2
Chemical Name: (R)-2-[4-(6-chloro-1,3-benzoxazol-2-yloxy)phenoxy]-2-fluoro-N-methylpropionanilide

PS: Solid


MF: C₂₃H₁₈ClFN₂O₄
MW: 440.9


Part Number	Matrix	Size
N-13084-25MG	Neat	25MG

Phenthoate
CAS#2597-03-7
Chemical Name: Ethyl-2-dimethoxyphosphinothioyl(phenyl)acetate

PS: Liquid

MF: C₁₂H₁₇O₄PS₂
MW: 320.38


MP: 17°C


Part Number	Matrix	Size
N-13006-100MG	Neat	100MG
S-13006A1-1ML	100ug/ml in Acetonitrile	1ML
S-13006U1-1ML	100ug/ml in Toluene	1ML

Flubendiamide
CAS#272451-65-7
Chemical Name: 3-iodo-N'-(2-Mesyl-1,1-dimethylethyl)-N-[4-[1,2,2,2-tetrafluoro-1-(trifluoromethyl)-ethyl]-o-tolyl]phthalamide

PS: Solid

MF: C₂₃H₂₂F₇IN₂O₄S

MW: 682.4


Part Number	Matrix	Size
N-11981-100MG	Neat	100MG

(2-Bromoethyl)trimethyl ammonium bromide CAS#2758-06-7
PS: Solid

MF: C₅H₁₃Br₂N

MW: 246.97

MP: 244°C


Alternate Name: BETAB

Part Number	Matrix	Size
N-10001-1G	Neat	1G

Valifenalate
CAS#283159-90-0
Chemical Name: -b-Alanine, N-[(1-methylethoxy)carbonyl]-L-valyl-3-(4-chlorophenyl)-, methyl ester

PS: Solid

MF: C₁₉H₂₇ClN₂O₅
MW: 398.9


MP: 147°C

Alternate Name: VALIPHENAL

Part Number	Matrix	Size
N-13281-25MG	Neat	25MG
S-13281A1-1ML	100ug/mL in Acetonitrile	1ML

Spiromesifen
CAS#283594-90-1
Chemical Name: 3-Mesityl-2-oxo-1-oxaspiro[4,4]non-3-en-4-yl 3,3-dimethylbutyrate

PS: Solid

MF: C₂₃H₃₀O₄
MW: 370.5

MP: 100°C

Alternate Name: Forbid(TM)

Part Number	Matrix	Size
N-13224-100MG	Neat	100MG


20% OFF

5 or more solutions or mixtures of the same item

Thiencarbazone-methyl
CAS#317815-83-1
Chemical Name: Methyl 4-[(4,5-dihydro-3-methoxy-4-methyl-5-oxo-1H-1,2,4-triazol-1-yl)carbonylsulfamoyl]-5-methylthiophene-3-carboxylate

PS: Solid


MF: C₁₂H₁₄N₄O₇S₂
MW: 390.4

MP: 206°C


Part Number	Matrix	Size
N-13565-50MG	Neat	50MG

Fluensulfone CAS#318290-98-1
Chemical Name: 5-chloro-1,3-thiazol-2-yl 3,4,4-trifluorobut-3-en-1-yl sulfone

PS: Solid


MF: C₇H₅ClF₃NO₂S₂
MW: 291.69

Alternate Name: 5-chloro-2-(3,4,4-trifluorobut-3-en-1-ylsulfonyl)-1,3-thiazole

Part Number	Matrix	Size
N-14084-100MG	Neat	100MG


Tembotrione CAS#335104-84-2
Chemical Name: 2-[2-Chloro-4-(methylsulfonyl)-3-[(2,2,2-trifluoroethoxy)methyl]benzoyl]-1,3-cyclohexanedione

PS: Solid


MF: C₁₇H₁₆ClF₃O₆S

MP: 123°C


Part Number	Matrix	Size
N-13508-100MG	Neat	100MG

Pyrifluquinazon CAS#337458-27-2
Chemical Name: 1-Acetyl-3,4-dihydro-3-[(3-pyridinylmethyl)amino]-6-[1,2,2,2-tetrafluoro-1-(trifluoromethyl)ethyl]-2(1H)-quinazolinone


PS: Solid

MF: C₁₉H₁₅F₇N₄O₂
MW: 464.3


Part Number	Matrix	Size
N-13158-100MG	Neat	100MG

Amisulbrom Solution CAS#348635-87-0
Chemical Name: 3-(3-Bromo-6-fluoro-2-methylindol-1-ylsulfonyl)-N,N-dimethyl-1H-1,2,4-triazole-1-sulfonamide


PS: Liquid

MF: C₁₃H₁₃BrFN₅O₄S₂
MW: 466.3


Part Number	Matrix	Size
S-13960M1-1ML	100ug/mL in Methanol	1ML

Amobam CAS#3566-10-7
Chemical Name: Ethylene bis[dithiocarbamic acid],diammonium salt


PS: Solid

MF: C₄H₈N₂S₄
MW: 246.46

Alternate Name: DIAMMONIUM ETHYLENEBIS(DITHIOCARBAMATE)

Part Number	Matrix	Size
N-11071-100MG	Neat	100MG
S-11071M1-1ML	100ug/ml in Methanol	1ML

Fluoxastrobin CAS#361377-29-9
Chemical Name: (E)-[2-[6-(2-chlorophenoxy)-5-fluoropyrimidin-4-yloxy]phenyl](5,6-dihydro-1,4,2-dioxazin-3-yl)methanone O-methyloxime

PS: Solid

MF: C₂₁H₁₆ClFN₄O₅
MW: 458.8



Part Number	Matrix	Size
N-13816-100MG	Neat	100MG

Pyrasulfotole CAS#365400-11-9
Chemical Name: (5-hydroxy-1,3-dimethyl-1H-pyrazol-4-yl)[2-(methylsulfonyl)-4-(trifluoromethyl)phenyl]methanone

PS: Solid

MF: C₁₄H₁₃F₃N₂O₄S

MW: 362.3

MP: 201°C



Part Number	Matrix	Size
N-14001-100MG	Neat	100MG

Saflufenacil
CAS#372137-35-4
Chemical Name: N'-[2-Chloro-4-fluoro-5-[1,2,3,6-tetrahydro-3-methyl-2,6-dioxo-4-(trifluoromethyl)pyrimidin-1-yl]benzoyl]-N-isopropyl-N-methylsulfamide

PS: Solid

MF: C₁₇H₁₇ClF₄N₄O₅S

MW: 500.9


MP: 190°C


Part Number	Matrix	Size
N-13195-100MG	Neat	100MG

Mandipropamid
CAS#374726-62-2
Chemical Name: Benzeneacetamide, 4-chloro-N-[2-[3-methoxy-4-(2-propyn-1-yloxy)phenyl]ethyl-a-(2-propyn-1-yloxy)-

PS: Solid

MF: C₂₃H₂₂ClNO₄
MW: 411.9

MP: 96.4-97.3°C

Alternate Name: REVUS™


Part Number	Matrix	Size
N-12354-100MG	Neat	100MG

Benazolin
CAS#3813-05-6
Chemical Name: 4-chloro-2-oxo-3(2H)-benzothiazoleacetic acid

PS: Solid

MF: C₉H₆ClNO₃S

MW: 243.67


MP: 193°C

Alternate Name: Herbazolin

Part Number	Matrix	Size
N-12901-100MG	Neat	100MG

Cyflumetofen
CAS#400882-07-7
Chemical Name: 2-methoxyethyl a-cyano-a-[4-(1,1-dimethylethyl)phenyl]-7-oxo-2-(trifluoromethyl)benzenepropanoate

PS: Solid

MF: C₂₄H₂₄F₃NO₄
MW: 447.5

MP: 77.9-81.7 °C

Alternate Name: DANISARABA


Part Number	Matrix	Size
N-13972-50MG	Neat	50MG

Pyroxsulam
CAS#422556-08-9
Chemical Name: N-(5,7-Dimethoxy[1,2,4]triazolo[1,5-a]pyrimidin-2-yl)-2-methoxy-4-(trifluoromethyl)-3-pyridinesulfonamide

PS: Solid

MF: C₁₄H₁₃F₃N₆O₅S

MW: 434.4


MP: 208°C


Part Number	Matrix	Size
N-13164-100MG	Neat	100MG

Pyroxa sulfone
CAS#447399-55-5
Chemical Name: 5-(difluoromethoxy)-1-methyl-3-(trifluoromethyl)pyrazol-4-yl-methyl 4,5-dihydro-5,5-dimethyl-1,2-oxazol-3-yl sulfone

PS: Solid

MF: C₁₂H₁₄F₅N₃O₄S

MW: 391.3



Part Number	Matrix	Size
N-13070-100MG	Neat	100MG
S-13070M1-1ML	100ug/ml in Methanol	1ML

Penflufen
CAS#494793-67-8
Chemical Name: 2'-[(RS)-1,3-Dimethylbutyl]-5-fluoro-1,3-dimethylpyrazole-4-carboxanilide

PS: Solid

MF: C₁₈H₂₄FN₃O


MW: 317.41

MP: 111.1°C

Alternate Name: EVERGOL™

Part Number	Matrix	Size
N-14035-100MG	Neat	100MG
S-14035A1-1ML	100ug/ml in Acetonitrile	1ML

Chlorantranilprole
CAS#500008-45-7
Chemical Name: 3-Bromo-4'-chloro-1-(3-chloro-2-pyridyl)-2'-methyl-6'-(methylcarbamoyl)pyrazole-5-carboxanilide

PS: Solid


MF: C₁₈H₁₄BrCl₂N₅O₂
MW: 483.2

Alternate Name: Rynaxypyr (TM)

Part Number	Matrix	Size
N-11422-100MG	Neat	100MG

Chloropropylate
CAS#5836-10-2
Chemical Name: 1-methylethyl 4-chloro-a-(4-chlorophenyl)-a-hydroxybenzeneacetate

PS: Solid

MF: C₁₇H₁₆Cl₂O₃
MW: 339.23


MP: 73°C


Part Number	Matrix	Size
N-11453-50MG	Neat	50MG

Fluopyram
CAS#658066-35-4
Chemical Name: N-[2-[3-chloro-5-(trifluoromethyl)-2-pyridinyl]ethyl]-2-(trifluoromethyl)benzamide

PS: Solid


MF: C₁₆H₁₁ClF₆N₂O

MW: 396.71


Part Number	Matrix	Size
N-13988-50MG	Neat	50MG
S-13988A1-1ML	100ug/mL in Acetonitrile	1ML

Cyantraniliprole
CAS#736994-63-1
Chemical Name: 3-bromo-1-(3-chloro-2-pyridyl)-4'-cyano-2'-methyl-6'-(methylcarbamoyl)pyrazole-5-carboxanilide

PS: Liquid


MF: C₁₉H₁₄BrClN₆O₂
MW: 473.71


Part Number	Matrix	Size
N-12886-10MG	Neat	10MG

Mancozeb (Technical)
CAS#8018-01-7
Chemical Name: [[1,2-Ethanedithiolbis(carbamodithioato)](2-)]manganese/zinc complex

PS: Solid


MF: C₄H₆MnN₂S₄·C₄H₆N₂S₄Zn

MP: 192°C


Part Number	Matrix	Size
N-12353-250MG	Neat	250MG

Aminocyclopyrachlor
CAS#858956-08-8
Chemical Name: 6-amino-5-chloro-2-cyclopropylpyrimidine-4-carboxylic acid

PS: Solid

MF: C₈H₈ClN₃O₂
MW: 213.6

Alternate Name: IMPRELIS™


Part Number	Matrix	Size
N-13826-250MG	Neat	250MG

Sedaxane
CAS#874967-67-6
Chemical Name: N-[2-[1,1'-bicyclopropyl]-2-ylphenyl]-3-(difluoromethyl)-1-methyl-1H-pyrazole-4-carboxamide

PS: Solid

MF: C₁₈H₁₉F₂N₃O

MW: 331.4

MP: 121.4°C



Part Number	Matrix	Size
N-13068-100MG	Neat	100MG
S-13068M1-1ML	100 ug/mL in Methanol	1ML

Isfetamid
CAS#875915-78-9
Chemical Name: N-[1,1-dimethyl-2-(4-isopropoxy-o-tolyl)-2-oxoethyl]-3-methylthiophene-2-carboxamide

PS: Solid

MF: C₂₀H₂₅NO₃S

MW: 359.48


MP: 103.5°C


Part Number	Matrix	Size
N-14079-100MG	Neat	100MG

Isopyrazam Solution
CAS#881685-58-1
Chemical Name: 3-(Difluoromethyl)-1-methyl-N-[1,2,3,4-tetrahydro-9-(1-methylethyl)-1,4-methanonaphthalen-5-yl]-1H-pyrazole-4-carboxamide

PS: Liquid

MF: C₂₀H₂₃F₂N₃O

MW: 259.4



Part Number	Matrix	Size
S-13996A1-1ML	100ug/ml in Acetonitrile	1ML

Fluxapyroxad
CAS#907204-31-3
Chemical Name: 3-(difluoromethyl)-1-methyl-N-(3,4,5'-trifluorobiphenyl-2-yl)pyrazole-4-carboxamide

PS: Solid

MF: C₁₈H₁₂F₅N₃O

MW: 381.3


MP: 157°C

Alternate Name: XEMIU®

Part Number	Matrix	Size
N-12921-100MG	Neat	100MG

Sulfoxaflor
CAS#946578-00-3
Chemical Name: N-[methyloxido[1-[6-(trifluoromethyl)-3-pyridinyl]ethyl]-g-4-sulfa-nylidene]cyanamide

PS: Solid


MF: C₁₀H₁₀F₃N₃OS

MW: 227.26


Part Number	Matrix	Size
N-12883-10MG	Neat	10MG
S-12883M1-1ML	100ug/ml in Methanol	1ML

Indaziflam
CAS#950782-86-2
Chemical Name: N-[(1R,2S)-2,3-dihydro-2,6-dimethyl-1H-inden-1-yl]-6-(1-fluoroethyl)-1,3,5-triazine-2,4-diamine

PS: Solid


MF: C₁₆H₂₀FN₅
MW: 301.36


Part Number	Matrix	Size
N-13074-100MG	Neat	100MG
S-13074A1-1ML	100ug/mL in Acetonitrile	1ML

Flupyradifurone
CAS#951659-40-8
Chemical Name: 4-[(6-chloro-3-pyridylmethyl)(2,2-difluoroethyl)amino]furan-2(5H)-one


PS: Solid

MF: C₁₂H₁₁ClF₂N₂O₂
MW: 288.68

MP: 69°C


Part Number	Matrix	Size
N-13941-100MG	Neat	100MG
S-13941A1-1ML	100ug/mL in Acetonitrile	1ML

2,4-D Choline salt
CAS#1048373-72-3
Chemical Name: 2,4-dichlorophenoxyacetic acid


PS: Solid


Part Number	Matrix	Size
N-13985-100MG	Neat	100MG
S-13985F1-1ML	100ug/mL in Water	1ML

Benzovindiflupyr Solution
CAS#1072957-71-1
Chemical Name: N-[(1R,4SR)-9-(dichloromethylene)-1,2,3,4-tetrahydro-1,4-methanonaphthalen-5-yl]-3-(difluoromethyl)-1-methylpyrazole-4-carboxamide

PS: Liquid

MF: C₁₈H₁₅Cl₂F₂N₃O


MW: 398.23


Part Number	Matrix	Size
S-14010X1-1ML	100ug/mL in Methylene chloride	1ML

Fluazaindolizine Solution
CAS#1254304-22-7
Chemical Name: 8-chloro-N-[(2-chloro-5-methoxyphenyl)sulfonyl]-6-(trifluoromethyl)imidazo[1,2-a]pyridine-2-carboxamide


PS: Liquid

MF: C₁₆H₁₀Cl₂F₃N₃O₄S

MW: 468.23


Part Number	Matrix	Size
S-14042A1-1ML	100ug/ml in Acetonitrile	1ML

DBP
CAS#84-74-2
Chemical Name: Di-n-butyl phthalate


PS: Liquid

Alternate Name: DBP

Part Number	Matrix	Size
S-11589A1-1ML	100ug/ml in Acetonitrile	1ML
S-11589J1-1ML	100ug/ml in Hexane	1ML
S-11589J1-5ML	100ug/ml in Hexane	5ML
N-11589-1G	Neat	1G

2,4,5-T butoxypolypropylene glycol ester (Technical)

(2,4,5-Trichlorophenoxy)acetic acid butoxypolypropylene glycol ester


PS: Liquid


Alternate Name: 2,4,5-T PGBEE

Part Number	Matrix	Size
N-10004-1G	Neat	1G

2,4,5-Trichlorophenol ethanolamine salt solution


PS: Liquid



Part Number	Matrix	Size
S-10655F1-1ML	100ug/ml in Water	1ML

2,4,5-Trichlorophenol ethanolamine salt

PS: Solid



Part Number	Matrix	Size
N-10655-250MG	Neat	250MG

Endothal-PFPH Solution

PS: Liquid


Part Number	Matrix	Size
S-11853T1-1ML	100ug/ml in t-Butyl methyl ether	1ML

4,6-Dinitro-2-sec-butylphenol ethanolamine salt
Chemical Name: 2-(sec-Butyl)-4,6-dinitrophenol, ethanolamine salt

PS: Solid

MF: C₁₂H₁₈N₃O₆

MW: 300.29


Alternate Name: 2-(sec-Butyl)-4,6-dinitrophenol, ethanolamine salt

Part Number	Matrix	Size
N-10888-1G	Neat	1G


o,p'-Dicofol Solution

CAS#10606-46-9

PS: Liquid

MF: C₁₄H₉Cl₅O

MW: 370.5




Part Number	Matrix	Size
S-13929K4-1ML	1000ug/ml in Isooctane	1ML

Analytical Standards and Certified Reference Material

WWW.CHEMSERVICE.COM
CHEM SERVICE INC




Metabolite Standards Listing by Parent Compound with Structure

a	→	N-12249		Isofenphos	← f
b	→	Isofenphos-des-N-isopropyl			
c	→	CAS# 25205-08-7			
			e ↓		← g
d	→	Part Number	Matrix	Size	← h
		MET-12249C-100MG	Neat	100MG	
		MET-12249CJ1-1ML	100 ug/ml in Hexane	1ML	

a Parent Compound Part Number

e Matrix

b Metabolite Description


f Parent Compound Description

c Metabolite Chemical Abstracts Registry Number


g Structure Image

h Size


d Metabolite Part Number

N-10688 2,6-Dichlorobenzonitrile
2,6-Dichlorobenzoic acid methyl ester
CAS# 14920-87-7



Part Number	Matrix	Size
MET-10688C-250MG	Neat	250MG

N-10688 2,6-Dichlorobenzonitrile
2,6-Dichlorobenzamide
CAS# 2008-58-4



Part Number	Matrix	Size
MET-10688B-1G	Neat	1G

N-10688 2,6-Dichlorobenzonitrile
2,6-Dichlorobenzoic acid
CAS# 50-30-6



Part Number	Matrix	Size
MET-10688A-250MG	Neat	250MG

N-10794 3,5-Diiodo-4-hydroxybenzonitrile
Ioxynil-methyl
CAS# 3336-40-1



Part Number	Matrix	Size
MET-10794A-100MG	Neat	100MG

N-10818 4-Chloro-o-tolyloxyacetic acid
4-Chloro-2-methylphenol
CAS# 1570-64-5



Part Number	Matrix	Size
MET-10818A-1G	Neat	1G

N-10876 4,4'-DDT
p,p'-DDMU
CAS# 1022-22-6



Part Number	Matrix	Size
MET-12810B-10MG	Neat	10MG

N-10992 a,b-Endosulfan
Endosulfan alcohol
CAS# 2157-19-9



Part Number	Matrix	Size
MET-11848A-50MG	Neat	50MG

N-10992 a,b-Endosulfan
Endosulfan ether
CAS# 3369-52-6



Part Number	Matrix	Size
MET-11848B-50MG	Neat	50MG

N-10992 a,b-Endosulfan
Endosulfan lactone
CAS# 3868-61-9



Part Number	Matrix	Size
MET-11848C-50MG	Neat	50MG

N-11009 Acetamiprid
Acetamiprid-N-desmethyl
CAS# 190604-92-3



Part Number	Matrix	Size
MET-11009A-10MG	Neat	10MG

N-11013 Acetochlor
Acetochlor ESA sodium salt
CAS# 187022-11-3



Part Number	Matrix	Size
MET-11013A-25MG	Neat	25MG
MET-11013AM1-1ML	100 ug/ml in Methanol	1ML

N-11013 Acetochlor
Acetochlor OA
CAS# 194992-44-4


Part Number	Matrix	Size
MET-11013B-10MG	Neat	10MG
MET-11013BM1-1ML	100 ug/ml in Methanol	1ML

N-11013 Acetochlor
2-Chloro-N-(2-ethyl-6-methylphenyl)
CAS# 32428-71-0



Part Number	Matrix	Size
MET-11013-50MG	Neat	50MG

N-11027 Acifluorfen
Acifluorfen, amino


Part Number	Matrix	Size
MET-11027A-10MG	Neat	10MG

N-11128 **Propoxur**


2-Isopropoxyphenol
CAS# 4812-20-8



Part Number	Matrix	Size
MET-11128A-250MG	Neat	250MG

N-11135 **Bendiocarb**


2,2-Dimethyl-1,3-benzodioxole-4-ol



Part Number	Matrix	Size
MET-11135A-500MG	Neat	500MG

N-11138 **Benomyl**


STB
CAS# 132827-25-9



Part Number	Matrix	Size
MET-11138B-500MG	Neat	500MG

N-11138 **Benomyl**


2-Aminobenzimidazole
CAS# 934-32-7



Part Number	Matrix	Size
MET-11138A-250MG	Neat	250MG

N-11142 **Bentazon**


2-Amino-N-isopropyl benzamide
CAS# 30391-89-0



Part Number	Matrix	Size
MET-11142C-50MG	Neat	50MG

N-11345 **Bromoxynil**


Bromoxynil-methyl ether
CAS# 3336-39-8



Part Number	Matrix	Size
MET-11345B-100MG	Neat	100MG

N-11345 **Bromoxynil**


Bromoxynil-heptanoate
CAS# 56634-95-8



Part Number	Matrix	Size
MET-11345A-50MG	Neat	50MG

N-11353 **Butachlor**


Butachlor ESA sodium salt
CAS# 1173022-75-7



Part Number	Matrix	Size
MET-11353A-25MG	Neat	25MG
MET-11353AM1-1ML	100ug/mL in Methanol	1ML

N-11399 **Captafol**


cis-1,2,3,6-Tetrahydrophthalimide
CAS# 1469-48-3



Part Number	Matrix	Size
MET-11399A-250MG	Neat	250MG

N-11405 **Carbofuran**


2,3-Dihydro-2,2-dimethylbenzofuran-7-ol
CAS# 1563-38-8



Part Number	Matrix	Size
MET-11405B-250MG	Neat	250MG

N-11405 **Carbofuran**


Carbofuran-3-keto
CAS# 16709-30-1



Part Number	Matrix	Size
MET-11405E-100MG	Neat	100MG

N-11405 **Carbofuran**


2,3-Dihydro-2,2-dimethyl benzofuran-3,7-diol
CAS# 17781-15-6



Part Number	Matrix	Size
MET-11405DA1-1ML	100 ug/ml in Acetonitrile	1ML

N-11405 **Carbofuran**


Carbofuranphenol-3-keto
CAS# 17781-16-7




Part Number	Matrix	Size
MET-11405F-50MG	Neat	50MG
MET-11405FM1-1ML	100ug/mL in Methanol	1ML

N-11408 **Carbophenothion**


Carbophenothion sulfone
CAS# 16662-85-4




Part Number	Matrix	Size
MET-11408D-25MG	Neat	25MG

N-11408 Carbophenothion
Carbophenothion sulfoxide
CAS# 17297-40-4



Part Number	Matrix	Size
MET-11408C-100MG	Neat	100MG

N-11408 Carbophenothion
Carbophenothion methyl-O-analog
CAS# 7173-84-4



Part Number	Matrix	Size
MET-11408BH10-1ML	10 ug/ml in Ethyl Acetate	1ML

N-11408 Carbophenothion
Carbophenothion methyl-O-analog
CAS# 7332-32-3



Part Number	Matrix	Size
MET-11408B-50MG	Neat	50MG

N-11408 Carbophenothion
Carbophenothion methyl
CAS# 953-17-3



Part Number	Matrix	Size
MET-11408AM1-1ML	100 ug/ml in Methanol	1ML

N-11423 Chlorbenside
Chlorbenside sulfoxide
CAS# 7047-28-1



Part Number	Matrix	Size
MET-11423B-100MG	Neat	100MG
MET-11423BX1-1ML	100 ug/ml in Methylene Chloride	1ML

N-11423 Chlorbenside
Chlorbenside sulfone
CAS# 7082-99-7



Part Number	Matrix	Size
MET-11423A-25MG	Neat	25MG
MET-11423AJ1-1ML	100ug/mL in Hexane	1ML

N-11425 Chlordane
Oxychlordane
CAS# 27304-13-8



Part Number	Matrix	Size
MET-11425AM1-1ML	100 ug/ml in Methanol	1ML

N-11425 Chlordane
Chlordene
CAS# 3734-48-3



Part Number	Matrix	Size
MET-11425D-250MG	Neat	250MG

N-11426 Chlordimeform
4-Chloro-2-methylaniline
CAS# 95-69-2



Part Number	Matrix	Size
MET-11426A-250MG	Neat	250MG
MET-11426AM1-1ML	100ug/mL in Methanol	1ML

N-11442 Chlorobenzilate
4,4'-Dichlorobenzhydrol
CAS# 90-97-1



Part Number	Matrix	Size
MET-11442A-250MG	Neat	250MG

N-11442 Chlorobenzilate
4,4'-Dichlorobenzophenone
CAS# 90-98-2



Part Number	Matrix	Size
MET-11442B-100MG	Neat	100MG

N-11455 Chlorotoluron
1-(3-Chloro-4-methylphenyl)-3-methylurea
CAS# 22175-22-0



Part Number	Matrix	Size
MET-11455A-100MG	Neat	100MG


N-11459 Chlorpyrifos
Chlorpyrifos Oxon
CAS# 5598-15-2



Part Number	Matrix	Size
MET-11459B-100MG	Neat	100MG


N-11459 Chlorpyrifos
3,5,6-Trichloro-2-pyridinol
CAS# 6515-38-4



Part Number	Matrix	Size
MET-11459A-250MG	Neat	250MG


N-11460		Chlorpyrifos Methyl
Chlorpyrifos methyl-O-analog CAS# 5598-52-7		
		
Part Number	Matrix	Size
MET-11460B-50MG	Neat	50MG


N-11462		Chlorthal-dimethyl
Monomethyl tetrachloroterephthalate CAS# 887-54-7		
		
Part Number	Matrix	Size
MET-11462A-100MG	Neat	100MG


N-11490		Clofentezine
2-Chlorobenzonitrile CAS# 873-32-5		
		
Part Number	Matrix	Size
MET-11490A-250MG	Neat	250MG


N-11507		Coumaphos
Coumaphos-O-analog CAS# 321-54-0		
		
Part Number	Matrix	Size
MET-11507D-50MG	Neat	50MG


N-11521		Cyclanilide
2,4-Dichloroaniline CAS# 554-00-7		
		
Part Number	Matrix	Size
MET-11521A-1G	Neat	1G


N-11522		Cycloate
N-Ethylcyclohexylamine CAS# 5459-93-8		
		
Part Number	Matrix	Size
MET-11522A-250MG	Neat	250MG


N-11543		Cyhalofop-butyl
Cyhalofop-amide		
		
Part Number	Matrix	Size
MET-11543A-50MG	Neat	50MG


N-11579		Deltamethrin
cis-Decamethrinic acid methyl ester CAS# 113830-50-5		
		
Part Number	Matrix	Size
MET-11579C-10MG	Neat	10MG


N-11579		Deltamethrin
3-(2,2-Dibromovinyl)-2,2-dimethyl-(1-cyclopropyl)-2-hydroxypropanoic acid CAS# 53179-78-5		
		
Part Number	Matrix	Size
MET-11579A-10MG	Neat	10MG
MET-11579BM1-1ML	100ug/mL in Methanol	1ML


N-11621		Diazinon
2-Isopropyl-6-methyl-4-pyrimidinol CAS# 2814-20-2		
		
Part Number	Matrix	Size
MET-11621B-250MG	Neat	250MG

N-11621		Diazinon
Diethyl phosphate CAS# 598-02-7		
		
Part Number	Matrix	Size
MET-11621C-100MG	Neat	100MG


N-11621		Diazinon
4-(Ethylthio)-6-methyl-2-(1-methylethyl)-pyrimidine CAS# 77738-92-2		
		
Part Number	Matrix	Size
MET-11621D-10MG	Neat	10MG

N-11621		Diazinon-O-analog
Diazinon-O-analog CAS# 962-58-3		
		
Part Number	Matrix	Size
MET-11621A-100MG	Neat	100MG


N-11656		Dicamba
3,6-Dichloro-2-hydroxy benzoic acid CAS# 3401-80-7		
		
Part Number	Matrix	Size
MET-11656B-50MG	Neat	50MG

N-11658 Dicapthon
2-Chloro-4-nitrophenol
CAS# 619-08-9



Part Number	Matrix	Size
MET-11658A-250MG	Neat	250MG

N-11677 Diclofop methyl
Diclofop acid
CAS# 40843-25-2


Part Number	Matrix	Size
MET-11677B-100MG	Neat	100MG

N-11721 Difenzoquat methylsulfate
Difenzoquat-desmethyl



Part Number	Matrix	Size
MET-11721C-100MG	Neat	100MG

N-11722 Diflubenazuron
p-Chlorophenylurea
CAS# 140-38-5



Part Number	Matrix	Size
MET-11722A-100MG	Neat	100MG

N-11754 Dimethachlor
Dimethachlor ESA sodium salt


Part Number	Matrix	Size
MET-11754A-25MG	Neat	25MG
MET-11754AM1-1ML	100ug/mL in Methanol	1ML

N-11755 Dimethenamid
Dimethenamid ESA sodium salt
CAS# 205939-58-8 (acid)



Part Number	Matrix	Size
MET-11755A-25MG	Neat	25MG
MET-11755AM1-1ML	100 ug/ml in Methanol	1ML

N-11755 Dimethenamid
Dimethenamid OA
CAS# 87674-68-8



Part Number	Matrix	Size
MET-11755B-10MG	Neat	10MG
MET-11755BM1-1ML	100ug/mL in Methanol	1ML

N-11819 Disulfoton
Disulfoton sulfone
CAS# 2497-06-500



Part Number	Matrix	Size
MET-11819A-100MG	Neat	100MG

N-11819 Disulfoton
Disulfoton sulfoxide
CAS# 2497-07-600



Part Number	Matrix	Size
MET-11819B-100MG	Neat	100MG

N-11827 Diuron
1-(3,4-Dichlorophenyl)-3-methylurea
CAS# 3567-62-2



Part Number	Matrix	Size
MET-11827C-250MG	Neat	250MG

N-11827 Diuron
N,N'-Bis(3,4-dichlorophenyl)urea
CAS# 4300-43-0



Part Number	Matrix	Size
MET-11827A-100MG	Neat	100MG

N-11827 Diuron
3-(4-Chlorophenyl)methyl urea
CAS# 5352-88-5



Part Number	Matrix	Size
MET-11827B-100MG	Neat	100MG


N-11840 Diodine (Technical)
Creatine
CAS# 57-00-1



Part Number	Matrix	Size
MET-11840A-250MG	Neat	250MG


N-11872 Ethiofencarb
Ethiofencarb sulfoxide
CAS# 53380-22-6



Part Number	Matrix	Size
MET-11872A-10MG	Neat	10MG


N-11872		Ethiofencarb
Ethiofencarb sulfone CAS# 53380-23-7		
		
Part Number	Matrix	Size
MET-11872B-10MG	Neat	10MG
MET-11872BA1-1ML	100 ug/ml in Acetonitrile	1ML


N-11941		Etofenprox
Etofenprox-desethyl CAS# 80854-21-3		
		
Part Number	Matrix	Size
MET-11941AM1-1ML	100ug/ml in Methanol	1ML


N-11944		Famphur
p-Hydroxy-N,N-dimethylbenzene sulfonamide CAS# 15020-57-2		
		
Part Number	Matrix	Size
MET-11944B-10MG	Neat	10MG


N-11944		Famphur
Famphur-O-analog CAS# 960-25-8		
		
Part Number	Matrix	Size
MET-11944A-50MG	Neat	50MG


N-11946		Fenamidon
Fenamidon Metabolite CAS# 332855-88-6		
		
Part Number	Matrix	Size
MET-11946-100MG	Neat	100MG


N-11948		Fenarimol
2,4-Dichlorobenzhydrol CAS# 43171-49-9		
		
Part Number	Matrix	Size
MET-11948A-100MG	Neat	100MG


N-11948		Fenarimol
2,4'-Dichlorobenzophenone CAS# 85-29-0		
		
Part Number	Matrix	Size
MET-11948-100MG	Neat	100MG


N-11952		Fenchlorphos
Fenchlorphos-oxon CAS# 3983-45-7		
		
Part Number	Matrix	Size
MET-11952A-100MG	Neat	100MG


N-11955		Fenitrothion
Fenitrothion-O-analog CAS# 2255-17-6		
		
Part Number	Matrix	Size
MET-11955AJ1-1ML	100 ug/ml in Hexane	1ML


N-11955		Fenitrothion
3-Methyl-4-nitrophenol CAS# 2581-34-2		
		
Part Number	Matrix	Size
MET-11955B-100MG	Neat	100MG

N-11957		Fenoxaprop ethyl
6-Chloro-2,3-dihydrobenzoxazol-2-one CAS# 19932-84-4		
		
Part Number	Matrix	Size
MET-11957A-100MG	Neat	100MG


N-11963		Fensulfothion
Fensulfothion-sulfone CAS# 14255-72-2		
		
Part Number	Matrix	Size
MET-11963A-50MG	Neat	50MG

N-11963		Fensulfothion
Fensulfothion-sulfide CAS# 3070-15-3		
		
Part Number	Matrix	Size
MET-11963B-100MG	Neat	100MG


N-11963		Fensulfothion
Fensulfothion oxon sulfone CAS# 6132-17-8		
		
Part Number	Matrix	Size
MET-11963CA1-1ML	100ug/ml in Acetonitrile	1ML

N-11963
Fensulfothion
Fensulfothion Oxon
CAS# 6552-21-2



Part Number	Matrix	Size
MET-11963DA1-1ML	100ug/mL in Acetonitrile	1ML

N-11964
Fenthion
Fenthion-oxon-sulfone
CAS# 14086-35-2



Part Number	Matrix	Size
MET-11964E-100MG	Neat	100MG

N-11964
Fenthion
4-Methylsulfonylphenol
CAS# 14763-60-1



Part Number	Matrix	Size
MET-11964G-50MG	Neat	50MG

N-11964
Fenthion
Fenthion sulfoxide
CAS# 3761-41-9



Part Number	Matrix	Size
MET-11964D-100MG	Neat	100MG
MET-11964BX1-1ML	100 ug/ml in Methylene Chloride	1ML

N-11964
Fenthion
Fenthion-sulfone
CAS# 3761-42-0



Part Number	Matrix	Size
MET-11964A-100MG	Neat	100MG

N-11964
Fenthion
Fenthion-O-analog
CAS# 6552-12-100



Part Number	Matrix	Size
MET-11964CK1-1ML	100 ug/ml in Isooctane	1ML

N-11964
Fenthion
Fenthion oxon sulfoxide
CAS# 6552-13-2



Part Number	Matrix	Size
MET-11964FA1-1ML	100ug/mL in Acetonitrile	1ML

N-11971
Fipronil
Fipronil sulfide
CAS# 120067-83-6



Part Number	Matrix	Size
MET-11971B-25MG	Neat	25MG
MET-11971BA1-1ML	100 ug/ml in Acetonitrile	1ML

N-11971
Fipronil
Fipronil sulfone
CAS# 120068-36-2



Part Number	Matrix	Size
MET-11971A-50MG	Neat	50MG
MET-11971AA1-1ML	100 ug/ml in Acetonitrile	1ML

N-11976
Flonicamid
4-(Trifluoromethyl)nicotinic acid
CAS# 158063-66-2



Part Number	Matrix	Size
MET-11976A-100MG	Neat	100MG

N-11986
Flufenacet
Flufenacet OA
CAS# 201668-31-7



Part Number	Matrix	Size
MET-11986B-10MG	Neat	10MG
MET-11986BM1-1ML	100 ug/ml in Methanol	1ML

N-11986
Flufenacet
Flufenacet ESA sodium salt
CAS# 947601-87-8



Part Number	Matrix	Size
MET-11986A-25MG	Neat	25MG
MET-11986AM1-1ML	100 ug/ml in Methanol	1ML


N-12004
Flutolanil
Methyl o-trifluoromethylbenzoate
CAS# 344-96-7



Part Number	Matrix	Size
MET-12004A-100MG	Neat	100MG


N-12004
Flutolanil
Methyl o-trifluoromethylbenzoate
CAS# 344-96-7



Part Number	Matrix	Size
MET-12004A-100MG	Neat	100MG


N-12007		Folpet
Phthalamic acid CAS# 88-97-1		
		
Part Number	Matrix	Size
MET-12007A-100MG	Neat	100MG


N-12017		Formothion
(Dimethoxythiophosphinoylthio)acetic acid CAS# 1113-01-5		
		
Part Number	Matrix	Size
MET-12017A-100MG	Neat	100MG


N-12133		Glyphosate
Aminomethyl phosphonic acid CAS# 1066-51-9		
		
Part Number	Matrix	Size
MET-12133A-100MG	Neat	100MG


N-12135		Glyphosine
Iminobis (methylphosphonic acid) CAS# 6711-48-4		
		
Part Number	Matrix	Size
MET-12135A-500MG	Neat	500MG


N-12139		Azinphos-methyl
Azinphos-methyl oxon CAS# 961-22-8		
		
Part Number	Matrix	Size
MET-12139AU1-1ML	100 ug/ml in Toluene	1ML


N-12159		Hexachlorobenzene
Pentachlorothiophenol CAS# 133-49-3		
		
Part Number	Matrix	Size
MET-12159B-100MG	Neat	100MG


N-12159		Hexachlorobenzene
Tetrachlorohydroquinone CAS# 87-87-6		
		
Part Number	Matrix	Size
MET-12159A-250MG	Neat	250MG


N-12200		Imazamethabenz-methyl
Imazamethabenz (free acid) CAS# 100728-84-5		
		
Part Number	Matrix	Size
MET-12200AA0-1ML	50ug/ml in Acetonitrile	1ML


N-12206		Imidacloprid
Imidacloprid olefin CAS# 115086-54-9		
		
Part Number	Matrix	Size
MET-12206-10MG	Neat	10MG


N-12206		Imidacloprid
Imidacloprid urea CAS# 120868-66-8		
		
Part Number	Matrix	Size
MET-12206A-50MG	Neat	50MG

N-12207		Phosmet
Phosmet-oxon CAS# 3735-33-9		
		
Part Number	Matrix	Size
MET-12207AU1-1ML	100 ug/ml in Toluene	1ML


N-12220		Iprodione
3,5-Dichloroaniline CAS# 626-43-7		
		
Part Number	Matrix	Size
MET-12220B-250MG	Neat	250MG

N-12249		Isofenphos
Isofenphos-des-N-isopropyl CAS# 25205-08-7		
		
Part Number	Matrix	Size
MET-12249C-100MG	Neat	100MG
MET-12249CJ1-1ML	100 ug/ml in Hexane	1ML


N-12249		Isofenphos
Isofenphos-des-N-isopropyl-O-analogue CAS# 31120-83-9		
		
Part Number	Matrix	Size
MET-12249BJ1-1ML	100 ug/ml in Hexane	1ML

N-12279 Isoproturon
Monodesmethyl Isoproturon
CAS# 34123-57-4



Part Number	Matrix	Size
MET-12279A-250MG	Neat	250MG

N-12287 Isoxaflutole
2-Methylsulfonyl-4-trifluoromethyl
CAS# 142994-06-7



Part Number	Matrix	Size
MET-12287B-100MG	Neat	100MG

N-12287 Isoxaflutole
1-(Methylsulfonyl-4-trifluoromethyl
CAS# 143701-75-1



Part Number	Matrix	Size
MET-12287A-100MG	Neat	100MG

N-12315 Leptophos
Leptophos oxon
CAS# 25006-32-0



Part Number	Matrix	Size
MET-12315A-50MG	Neat	50MG

N-12346 Malathion
Malathion-O-analog
CAS# 1634-78-2



Part Number	Matrix	Size
MET-12346C-100MG	Neat	100MG

N-12346 Malathion
Isomalathion
CAS# 3344-12-500



Part Number	Matrix	Size
MET-12346AJ1-1ML	100 ug/ml in Hexane	1ML

N-12346 Malathion
Malathion monocarboxylic acid
CAS# 35884-76-5


Part Number	Matrix	Size
MET-12346D-100MG	Neat	100MG
MET-12346BT1-1ML	100 ug/ml in T-Butylmethyl ether	1ML

N-12348 MCPA
MCPA sodium salt
CAS# 3653-48-3



Part Number	Matrix	Size
MET-12348-250MG	Neat	250MG

N-12371 Mepiquat chloride
N-Methylpiperidine
CAS# 626-67-5



Part Number	Matrix	Size
MET-12371A-250MG	Neat	250MG

N-12382 Metam sodium
Methyl isothiocyanate
CAS# 556-61-6



Part Number	Matrix	Size
MET-12392A-250MG	Neat	250MG

N-12383 Metamitron
Pentachloroaniline
CAS# 527-20-8



Part Number	Matrix	Size
MET-12383A-100MG	Neat	100MG

N-12398 Methiocarb
Methiocarb sulfone
CAS# 2179-25-1



Part Number	Matrix	Size
MET-12398A-50MG	Neat	50MG


N-12398 Methiocarb
Methiocarb sulfoxide
CAS# 2635-10-100



Part Number	Matrix	Size
MET-12398BU1-1ML	100 ug/ml in Toluene	1ML


N-12403 Methoxychlor
p,p'-Methoxychlor-olefin
CAS# 2132-70-9



Part Number	Matrix	Size
MET-12403B-50MG	Neat	50MG
MET-12403BM1-1ML	100 ug/ml in Methanol	1ML


N-12403		Methoxychlor
o,p'-Methoxychlor CAS# 30667-99-3		
		
Part Number	Matrix	Size
MET-12403A-100MG	Neat	100MG


N-12403		Methoxychlor
o,p'-Methoxychlor CAS# 72-43-5		
		
Part Number	Matrix	Size
MET-12403AM1-1ML	100 ug/ml in Methanol	1ML


N-12478		Metolachlor
Metolachlor OA CAS# 152019-73-3		
		
Part Number	Matrix	Size
MET-12478D-10MG	Neat	10MG
MET-12478DM1-1ML	100 ug/ml in Methanol	1ML


N-12478		Metolachlor
2-[(2-Ethyl-6-methylphenyl)amino]-1-propanol CAS# 61520-53-4		
		
Part Number	Matrix	Size
MET-12478A-50MG	Neat	50MG


N-12478		Metolachlor
Metolachlor ESA sodium salt CAS# 947601-85-6		
		
Part Number	Matrix	Size
MET-12478C-25MG	Neat	25MG
MET-12478CM1-1ML	100 ug/ml in Methanol	1ML


N-12497		Monuron
N,N'-Bis(p-chlorophenyl)urea CAS# 1219-99-4		
		
Part Number	Matrix	Size
MET-12497A-500MG	Neat	500MG


N-12507		N-1-Naphthylphthalamic acid
N-(1-Naphthyl)phthalimide CAS# 5333-99-3		
		
Part Number	Matrix	Size
MET-12507A-250MG	Neat	250MG


N-12659		Nitrapyrin
6-Chloro-2-picolinic acid methyl ester CAS# 6636-55-1		
		
Part Number	Matrix	Size
MET-12659BJ1-1ML	100 ug/ml in Hexane	1ML


N-12706		o,p'-DDD
o,p'-DDD olefin CAS# 3424-82-6		
		
Part Number	Matrix	Size
MET-12706A-10MG	Neat	10MG


N-12708		o,p'-DDT
p,p'-DDA CAS# 83-05-6		
		
Part Number	Matrix	Size
MET-12708A-100MG	Neat	100MG

N-12727		Propargite
2-(4-tert-Butylphenoxy)cyclohexanol CAS# 1942-71-8		
		
Part Number	Matrix	Size
MET-12727A-100MG	Neat	100MG


N-12727		Propargite
1,2-Cyclohexanediol (cis & trans mixture) CAS# 931-17-9		
		
Part Number	Matrix	Size
MET-12727B-250MG	Neat	250MG

N-12731		Oxadiazon
Oxadiazon-hydroxy CAS# 39807-19-7		
		
Part Number	Matrix	Size
MET-12731A-50MG	Neat	50MG
MET-12731AU10-1ML	10 ug/ml in Toluene	1ML


N-12731		Oxadiazon
Pivaloyl 2-(2,4-dichloro-5-hydroxyphenyl)hydrazide CAS# 39807-20-0		
		
Part Number	Matrix	Size
MET-12731BM1-1ML	100ug/mL in Methanol	1ML

N-12735 Oxyamyl
Oxamyl oxime
CAS# 30558-43-1



Part Number	Matrix	Size
MET-12735-10MG	Neat	10MG

N-12848 Permethrin
m-Phenoxybenzyl alcohol
CAS# 13826-35-2



Part Number	Matrix	Size
MET-12848B-100MG	Neat	100MG

N-12848 Permethrin
3-Phenoxybenzoic acid
CAS# 3739-38-6



Part Number	Matrix	Size
MET-12848C-100MG	Neat	100MG

N-12848 Permethrin
3-(2,2-Dichlorovinyl)-2,2-dimethyl-(1-cycloprop.)carb. acid
CAS# 55701-05-8



Part Number	Matrix	Size
MET-12848A-10MG	Neat	10MG
MET-12848AM1-1ML	100 ug/ml in Methanol	1ML

N-12854 Fenamiphos
Fenamiphos sulfoxide
CAS# 31972-43-7



Part Number	Matrix	Size
MET-11947A-50MG	Neat	50MG

N-12854 Fenamiphos
Fenamiphos sulfone
CAS# 31972-44-8



Part Number	Matrix	Size
MET-11947B-100MG	Neat	100MG

N-12890 Sulcotrione
2-Chloro-4-(methylsulfonyl)benzoic acid
CAS# 53250-83-2



Part Number	Matrix	Size
MET-12890A-100MG	Neat	100MG

N-12950 Sebuthylazin
Sebuthylazine-desethyl
CAS# 37019-18-4



Part Number	Matrix	Size
MET-12950A-100MG	Neat	100MG

N-13035 Phorate
Phorate sulfoxide
CAS# 2588-03-600



Part Number	Matrix	Size
MET-13035C-100MG	Neat	100MG

N-13035 Phorate
Phorate sulfone
CAS# 2588-04-700



Part Number	Matrix	Size
MET-13035B-100MG	Neat	100MG
MET-13035BJ1-1ML	100 ug/ml in Hexane	1ML

N-13035 Phorate
Phorate-oxon
CAS# 2600-69-3



Part Number	Matrix	Size
MET-13035AU1-1ML	100 ug/ml in Toluene	1ML

N-13063 Pirimiphos-ethyl
Pirimiphos-ethyl-oxon
CAS# 36378-61-7



Part Number	Matrix	Size
MET-13063A-250MG	Neat	250MG


N-13064 Pirimiphos-methyl
N,N-diethyl-4-methoxy-6-methyl-2-Pyrimidinamine
CAS# 111697-06-4



Part Number	Matrix	Size
MET-13064E-100MG	Neat	100MG


N-13064 Pirimiphos-methyl
2-Diethylamino-6-methyl pyrimidin-4-ol
CAS# 42487-72-9



Part Number	Matrix	Size
MET-13064A-100MG	Neat	100MG


N-13064		Pirimiphos-methyl
2-Ethylamino-6-methyl-4(3H)-Pyrimidinone		
CAS# 5734-69-0		
		
Part Number	Matrix	Size
MET-13064C-100MG	Neat	100MG


N-13064		Pirimiphos-methyl
Pirimiphos-methyl-oxon		
CAS# 64709-45-1		
		
Part Number	Matrix	Size
MET-13064B-250MG	Neat	250MG


N-13064		Pirimiphos-methyl
Pirimiphos-methyl-N-desethyl		
CAS# 67018-59-1		
		
Part Number	Matrix	Size
MET-13064D-100MG	Neat	100MG


N-13068		Sedaxane
3-(Difluoromethyl)-1-methyl-1H-pyrazole-4-carboxylic acid		
CAS# 176969-34-9		
		
Part Number	Matrix	Size
MET-13068BA1-1ML	100 ug/ml in Acetonitrile	1ML


N-13079		Halfenprox
De-Bromo-halfenprox		
		
Part Number	Matrix	Size
MET-13079A-10MG	Neat	10MG


N-13097		Profenofos
4-Bromo-2-chlorophenol		
CAS# 3964-56-5		
		
Part Number	Matrix	Size
MET-13097A-500MG	Neat	500MG


N-13104		Propachlor
Propachlor OA		
CAS# 70628-36-3		
		
Part Number	Matrix	Size
MET-13104B-10MG	Neat	10MG
MET-13104BM1-1ML	100 ug/ml in Methanol	1ML


N-13104		Propachlor
Propachlor ESA sodium salt		
CAS# 947601-88-9		
		
Part Number	Matrix	Size
MET-13104A-25MG	Neat	25MG
MET-13104AM1-1ML	100 ug/ml in Methanol	1ML


N-13108		Propazine
Propazine-2-hydroxy		
CAS# 7374-53-0		
		
Part Number	Matrix	Size
MET-13108AAA1-1ML	100 ug/ml in Acetonitrile:water(75:25)	1ML


N-13140		Pendimethalin
N-Nitrosopendimethalin		
CAS# 68897-50-7		
		
Part Number	Matrix	Size
MET-13140A-100MG	Neat	100MG

N-13142		Pyraclufos
1-(4-Chlorophenyl)-1H-pyrazol-4-ol		
CAS# 77458-30-1		
		
Part Number	Matrix	Size
MET-13142A-100MG	Neat	100MG


N-13144		Pyraflufen-ethyl
Pyraflufen		
CAS# 129630-17-7		
		
Part Number	Matrix	Size
MET-13144A-10MG	Neat	10MG

N-13174		Quizalofop ethyl
Quizalofop		
CAS# 76578-12-6		
		
Part Number	Matrix	Size
MET-13174A-50MG	Neat	50MG


N-13174		Quizalofop ethyl
Quizalofop methyl		
CAS# 76578-13-7		
		
Part Number	Matrix	Size
MET-13174B-50MG	Neat	50MG

N-13176
Resmethrin
5-(Phenylmethyl)-3-furanmethanol
CAS# 20416-09-5



Part Number	Matrix	Size
MET-13176A-100MG	Neat	100MG

N-13186
Ruelene (TM)
4-tert-Butyl-2-chlorophenol
CAS# 98-28-2



Part Number	Matrix	Size
MET-13186A-10MG	Neat	10MG

N-13201
Fenvalerate
2-(p-Chlorophenyl)-3-methyl butyric acid
CAS# 2012-74-0



Part Number	Matrix	Size
MET-13201A-100MG	Neat	100MG

N-13210
Sethoxydim
Sethoxydim-desethoxy
CAS# 113561-03-8



Part Number	Matrix	Size
MET-13210A-50MG	Neat	50MG

N-13210
Sethoxydim
Sethoxydim-oxazole
CAS# 413594-06-6



Part Number	Matrix	Size
MET-13210B-50MG	Neat	50MG

N-13213
Silvex isooctyl ester
Atrazine desisopropyl
CAS# 1007-28-9



Part Number	Matrix	Size
MET-13213A-50MG	Neat	50MG

N-13213
Silvex isooctyl ester
2-Chloro-4,6-diamino-1,3,5-triazine
CAS# 3397-62-4



Part Number	Matrix	Size
MET-13213C-100MG	Neat	100MG

N-13224
Spiromesifen
Spiromesifen alcohol
CAS# 148476-30-6



Part Number	Matrix	Size
MET-13224-10MG	Neat	10MG

N-13258
Sulprofos
Sulprofos-sulfoxide
CAS# 34643-47-5



Part Number	Matrix	Size
MET-13258C-100MG	Neat	100MG

N-13258
Sulprofos
Sulprofos sulfone Solution
CAS# 92642-35-8



Part Number	Matrix	Size
MET-13258AU1-1ML	100 ug/ml in Toluene	1ML
MET-13258AU1-5ML	100 ug/ml in Toluene	5ML

N-13510
Terbufos
Terbufos sulfoxide
CAS# 10548-10-4



Part Number	Matrix	Size
MET-13510B-100MG	Neat	100MG

N-13510
Terbufos
Terbufos sulfone
CAS# 56070-16-7



Part Number	Matrix	Size
MET-13510A-50MG	Neat	50MG


N-13511
Terbumeton
Terbumeton-desethyl
CAS# 30125-64-5



Part Number	Matrix	Size
MET-13511-50MG	Neat	50MG


N-13512
Terbuthylazine
Terbuthylazine-desethyl
CAS# 30125-63-4



Part Number	Matrix	Size
MET-13512A-50MG	Neat	50MG


N-13561		Thiacloprid
Thiacloprid-amide CAS# 676228-91-4		
		
Part Number	Matrix	Size
MET-13561-100MG	Neat	100MG


N-13562		Thiamethoxam
N-Desmethylthiamethoxam CAS# 171103-04-1		
		
Part Number	Matrix	Size
MET-13562A-50MG	Neat	50MG
MET-13562AA1-1ML	100ug/mL in Acetonitrile	1ML


N-13576		Propiconazole
1,2,4-Triazole CAS# 288-88-0		
		
Part Number	Matrix	Size
MET-13576A-250MG	Neat	250MG


N-13576		Propiconazole
1-(2,4-Dichlorophenyl)-2-(1H-1,2,4-triazole-1-yl)-ethanone CAS# 58905-16-1		
		
Part Number	Matrix	Size
MET-13576B-100MG	Neat	100MG


N-13584		Tolyfluanid
DMST CAS# 66840-71-9		
		
Part Number	Matrix	Size
MET-13584AM1-1ML	100ug/ml in Methanol	1ML


N-13666		Trietazine
Trietazine-2-hydroxy CAS# 13532-25-7		
		
Part Number	Matrix	Size
MET-13666B-10MG	Neat	10MG
MET-13666BA1-1ML	100 ug/ml in Acetonitrile	1ML


N-13666		Trietazine
Trietazine-desethyl CAS# 38902-68-0		
		
Part Number	Matrix	Size
MET-13666A-100MG	Neat	100MG
MET-13666AA1-1ML	100ug/mL in Acetonitrile	1ML


N-13706		Trinexapac-ethyl
Trinexapac CAS# 104273-73-6		
		
Part Number	Matrix	Size
MET-13706A-100MG	Neat	100MG


N-13750		Warfarin
4-Hydroxycoumarin CAS# 1076-38-6		
		
Part Number	Matrix	Size
MET-13750B-250MG	Neat	250MG


N-13750		Warfarin
7-Hydroxycoumarin CAS# 93-35-6		
		
Part Number	Matrix	Size
MET-13750A-250MG	Neat	250MG

N-13760		Zinophos (TM)
2-Hydroxypyrazine CAS# 6270-63-9		
		
Part Number	Matrix	Size
MET-13760BM1-1ML	100 ug/ml in Methanol	1ML


N-13760		Zinophos (TM)
Thionazin-O-analog CAS# 7359-55-9		
		
Part Number	Matrix	Size
MET-13760A-50MG	Neat	50MG

N-13800		Simazine
Simazine hydroxy CAS# 2599-11-300		
		
Part Number	Matrix	Size
MET-13800B-25MG	Neat	25MG


N-13062		Pirimicarb
Pirimicarb-desmethyl CAS# 30614-22-3		
		
Part Number	Matrix	Size
MET-13801A-50MG	Neat	50MG

N-13830 Fenoxaprop
Fenoxaprop-P
CAS# 113158-40-0



Part Number	Matrix	Size
MET-13802-100MG	Neat	100MG

N-13932 Isoxathion
3-Hydroxy-5-phenylisoxazole
CAS# 939-05-9



Part Number	Matrix	Size
MET-13932A-10MG	Neat	10MG

N-14030 Fermizone
2-Hydrazino-4,6-dimethylpyrimidine
CAS# 23906-13-0



Part Number	Matrix	Size
MET-14030A-100MG	Neat	100MG
MET-14030AM1-1ML	100 ug/ml in Methanol	1ML

N-14013 Fenfuram
2-Methyl-3-furic acid
CAS# 6947-94-0



Part Number	Matrix	Size
MET-14013A-100MG	Neat	100MG

N-14035 Penflufen
2-(1,3-dimethylbutyl)-Benzenamine
CAS# 203448-76-4



Part Number	Matrix	Size
MET-14035A-10MG	Neat	10MG

N-14038 Cyhalofop
Cyhalofop-diacid
CAS# 252564-94-6



Part Number	Matrix	Size
MET-14038A-100MG	Neat	100MG

N-13972 Cyflumetofen
2-Methoxyethyl 2-[4-(tert-butyl)phenyl]
CAS# 400882-57-7


Part Number	Matrix	Size
MET-13987A-50MG	Neat	50MG

N-12206 Imidicloprid
6-Chloronicotinic acid
CAS# 5326-23-8


Part Number	Matrix	Size
MET-12066B-100MG	Neat	100MG
MET-12066BM1-1ML	100ug/mL in Methanol	1ML

N-11013 S-Metolachlor
sec-Aceto chlor OA
CAS# 152019-74-4


Part Number	Matrix	Size
S-13946M0-1ML	50ug/mL in Methanol	1ML

CHEM SERVICE INC Custom Solutions and Mixture Form

ACCOUNT NUMBER _____ DATE _____

CUSTOMER NAME _____ COMPANY NAME _____

ADDRESS _____

CITY/STATE/ZIP _____

PHONE _____ FAX _____ EMAIL _____

If reordering a custom standard please provide Part# or Quote# _____

EPA METHOD

Chem Service, Inc. • Attn: Custom Standards • P.O. Box 599 • West Chester, PA 19381-0599

Fax#: (610) 692-8729 or Email us: info@chemservice.com

	ANALYTE	CATALOG #	CAS#	CONCENTRATION
1.				
2.				
3.				
4.				
5.				
6.				
7.				
8.				
9.				
10.				
11.				
12.				
13.				
14.				
15.				
16.				
17.				
18.				
19.				
20.				
21.				
22.				

TOTAL # OF ANALYTES: _____

DESIRED CONCENTRATION AND SOLVENT: _____

VOLUME (check one): 1mL _____ 2mL _____ 5mL _____ 10mL _____ Other _____ # of Units: _____

MSDS provided with all standards.

Amber Ampules used unless specified otherwise: _____

DOCUMENTATION (check one):

☐ **Gravimetric Certification:** Guarantees the analytes in the custom standard to be made within $\pm 5\%$ of specifications.

☐ **Quantitative Certification:** Includes verification of the concentration and peak identification.

All mixtures are prepared in accordance with our ISO 9001 & ISO 17025 registration. All analytes are certified $\geq 98\%$ unless noted. Corrections are made to 100% purity.

Government/State/International Agency Methods

The Government/State/International Agency chemical standard section provides the necessary standards to run various published methods. These solutions are tested and certified according to a stringent quality criterion. They follow the guidelines for standard preparation for the indicated methods and conform to the high standards of our ISO certification program. We retain samples of the test solutions for on going stability studies.

First the purity and identification of each compound is established before the compound is utilized for a solution or mixture. For stock and custom solutions, corrections are made for any compound having a starting purity of less than 99%. Certificate of Analyses are provided on request. All mixtures and solutions are shipped with Safety Data Sheets. All of the solutions are traceable by a lot number.

Provided with each solution is a list of each of the individual compounds in the solution. Each compound is listed by their part number and the individual compounds can be purchased, as individual materials at 98+% purity, if so desired.

Standards

- EPA 500 Series Methods for Drinking Water (504 through 555).
- EPA 600 Series Methods for Wastewater Pollutants (608 through 680).
- EPA 1600 Series Methods for Wastewater Pollutants (1618 through 1660).
- EPA 8000 and Series Methods for Monitoring Organic Pollutants in Groundwater, Wastewater, and Solid Waste.
- EPA Contract Laboratory Program Standards for Target Compound List in Ground Water, Sediment, and Soil.
- International Standards for Canada (MISA), German DIN Standards, Europe Regulations, and ISO.

EPA 500 Series Methods for Drinking Water

Method #	Method Description	Page
504, 504.1	EDB, DBCP, and 1,2,3-Trichloropropane (123TCP) (Microextraction & GC)	181
505	Organohalide Pesticides and Commercial (PCB) Products (Microextraction & GC)	181
507	Nitrogen and Phosphorus-Containing Pesticides (GC)	181-182
508, 508A, 508.1	Chlorinated Pesticides (GC/ECD)	182-183
509	Ethylene Thiourea (ETU) in Water (GC)	183
515.1, 515.2,		
515.3, 515.4	Chlorinated Acids in Water (GC/ECD)	184-185
525.2	Organic Compounds in Drinking Water (LSE, GC/MS)	185-186
526	Semi-Volatile Organic Compounds in Drinking Water (GC/MS)	187
527	Pesticide and Flame Retardants in Drinking Water (GC/MS)	187
531.1	N-Methylcarbamoyloximes and N-Methylcarbamates in Water (HPLC). 152	187
532	Phenylurea Pesticides (HPLC)	187
535	Chloroacetanilide and Other Acetamide Herbicide Degradates in Drinking Water	188
547	Glyphosate in Drinking Water (HPLC)	188
548, 548.1	Endothall in Drinking Water (LSE, GC/ECD)	188
549	Diquat And Paraquat in Drinking Water (LSE, & HPLC/UV)	188
555	Chlorinated Acids in Drinking	189
608, 608.1, 608.2	Organochlorine Pesticides & PCB's (GC/ECD)	190
610	Polynuclear Aromatic Hydrocarbons (GC/FID), (HPLC)	191
614, 614.1	Organophosphorous Pesticides in Water (GC/MS)	191
615	Chlorinated Herbicides in Wastewaters (GC/ECD)	191-192
617	Organohalide Pesticides & PCB's in Wastewaters (GC/ECD)	192
619	Triazine Pesticides/Herbicides in Wastewaters (GC)	192
622	Organophosphorous Pesticides in Wastewaters (GC)	193
625	Base/Neutrals & Acids, Semivolatiles, Pesticides, Aroclors (GC/MS)	193-194
627	Dinitroaniline Pesticides in Wastewaters (GC)	194
629	Cyanazine (HPLC UV) with optional Florisil column cleanup procedure	194
630, 630.1	Dithiocarbamate Pesticides in Wastewaters (UV/VIS)	194
631	Benomyl & Carbendazim in Wastewater (HPLC)	194
632	Carbamate & Urea Pesticides in Wastewater (HPLC)	195
634	Thiocarbamate Pesticides in Wastewaters (GC)	195
635	Determination of Rotenone in Wastewaters (HPLC/UV)	195
636	Determination of Bensulide in Wastewaters (HPLC/UV)	195
638	Oryzal in Wastewaters (HPLC/UV)	196
639	Bendiocarb in Wastewaters (HPLC/UV)	196
640	Mercaptobenzothiazole in Wastewaters (HPLC/UV)	196
641	Thiabendazole in Wastewaters (HPLC)	196
643	Bentazon in Wastewaters (HPLC)	196
644	Picloram in Wastewater (HPLC/UV)	196
645	Amine Pesticides & Lethane in Wastewater (GC/MS)	197
680	Pesticides & PCB Congeners (GC/MS)	197
1618	Organohalide Pesticides, Organophosphorus Pesticides and Phenoxy-acid Herbicides by GC	197-198
1656	Organohalide Pesticides by GC/HSD	198
1657	Organophosphorus Pesticides by GC/FPD	199
1660	Pyrethrins and Pyrethroids by HPLC/UV	199
8080	Organochlorine Pesticides & PCB's (GC/ECD)	199-200
8140	Organophosphorous Pesticides (GC)	200
8141A	Organophosphorous Compounds (GC)	200-201
8150B	Chlorinated Herbicides (GC)	201-202
8151	Chlorinated Herbicides by GC using Methylation or Pentafluorobenzoylation Derivatization	202
8250A	Semivolatile Organic Compounds (GC/MS)	203-204
8270B	Semivolatile Organic Compounds (GC/MS)	203-204
8318	N-methylcarbamates (HPLC)	204
8321	Solvent Extractable Non-volatile Compounds (HPLC/TSP/MS)	205
CLP	Contract Lab Program Standards	205-207
USP 561	Chlorinated and Phosphorus Pesticides	207
ISO 6468	Organochlorine Insecticides	207
ISO 10695	Nitrogen and Phosphorous Compounds	208
ISO 11369	Plant Treatment Agents	208
ISO 15913	Phenoxyalkanoic Herbicides	208
Canadian	Canadian Drinking Water Standards	208
DIN 38407	German Standard Methods	209
European	European Regulation Standards	209-211
Cannabis	Pesticide Residue Standards	211

METHOD 504, 504.1 - EDB & DBCP

Method 504 is applicable for the determination of 1,2-Dibromoethane (EDB) and 1,2-Dibromo-3-chloropropane (DBCP) in finished drinking water, raw source water, or drinking water in any treatment stage by microextraction and gas chromatography.

Method 504.1 is applicable to the determination of 1,2-Dibromoethane (EDB), 1,2-Dibromo-3-chloropropane and 1,2,3-Trichloropropane (123TCP) in water by microextraction and gas chromatography in finished drinking water and groundwater:

EDB/DBCP MIXTURE**EPA METHOD 504, 8010B, 8011**

1,2-Dibromo-3-chloropropane (DBCP)
1,2-Dibromoethane (EDB)

2000ug/mL in Methanol (2 components)
M-CSHC16M5-1ML 1 mL Ampule

EDB/DBCP/123TCA STANDARDS MIXTURE**EPA METHOD 504.1**

1,2-Dibromo-3-chloropropane (DBCP)
1,2-Dibromoethane (EDB)
1,2,3-Trichloropropane

2000ug/mL in Methanol (3 components)
M-CSHC16AM5-1ML 1 mL Ampule

METHOD 505 - ORGANOHALIDE PESTICIDES & COMMERCIAL PCBs

Method 505 is applicable to the determination of organohalide pesticides and commercial polychlorinated biphenyl (PCB) products in finished drinking water, during intermediate stages of treatment, and raw source water, by microextraction and gas chromatography.

ORGANOHALIDE PESTICIDES MIXTURE #1**EPA METHOD 505**

Alachlor	Aldrin
Dieldrin	Endrin
Heptachlor	Heptachlor epoxide
Hexachlorobenzene	Lindane (BHC gamma isomer)
Methoxychlor	

100ug/mL in Toluene (9 components)
M-OHP5051U1-1ML 1 mL Ampule

ORGANOHALIDE PESTICIDES MIXTURE #2**EPA METHOD 505**

Hexachlorocyclopentadiene	Heptachlor
cis-Chlordane	cis-Nonachlor
trans-Nonachlor	trans-Chlordane
Atrazine	Simazine

100ug/mL in Methanol (8 components)
M-OHP5052M1-1ML 1 mL Ampule

TOXAPHENE STANDARDS**EPA METHOD 505**

N-13586-250MG	Toxaphene
S-13586M1-1ML	100ug/mL in Methanol
S-13586M1-5ML	100ug/mL in Methanol
S-13586J4-1ML	1000ug/mL in Hexane
S-13586J4-5ML	1000ug/mL in Hexane
S-13586B6-1ML	2500ug/mL in Aceone

INSTRUMENT PERFORMANCE CHECK STANDARDS**EPA METHOD 505**

N-11854-100MG	Endrin
S-11854M1-1ML	100ug/mL in Methanol
S-11854M1-5ML	100ug/mL in Methanol
N-11855-10MG	Endrin aldehyde
S-11855M1-1ML	100ug/mL in Methanol
S-11855M1-5ML	100ug/mL in Methanol
N-11856-10MG	Endrin ketone
S-11856A1-1ML	100ug/mL in Acetonitrile
S-11856U1-5ML	100ug/mL in Toluene

METHOD 507 - NITROGEN/PHOSPHOROUS PESTICIDES

Method 507 is applicable to the determination of Nitrogen and Phosphorus containing pesticides in water by gas chromatography in ground water and finished drinking water.

NITROGEN/PHOSPHOROUS PESTICIDES MIXTURE**#1 - EPA METHOD 507, 525.2**

Ametryne	Cycloate
Disulfoton	Tributylphosphorotriothioate
Prometon	Fenamiphos
*Inhibited with Hydrquinone	

1000ug/mL in tert-Butyl methyl ether (6 components)
M-PM5071T4-1ML 1 mL Ampule

NITROGEN/PHOSPHOROUS PESTICIDES MIXTURE**#2 - EPA METHOD 507, 525.2**

Atrazine	Diphenamid
Mevinphos	Prometryne
Propazine	Prophos
S-Ethyl dipropylthiocarbamate	Terbutryne
Triadimefon	

1000ug/mL in tert-Butyl methyl ether (9 components)
M-PM5072T4-1ML 1 mL Ampule

NITROGEN/PHOSPHOROUS PESTICIDES MIXTURE**#3 - EPA METHOD 507, 525.2**

Butachlor	Carboxin
Diazinon	Metolachlor
Metribuzin	MKG 264 (TM)
Norflurazon	Terbufos
Vernolate	

1000ug/mL in tert-Butyl methyl ether (9 components)
M-PM5073T4-1ML 1 mL Ampule

NITROGEN/PHOSPHOROUS PESTICIDES MIXTURE**#4 - EPA METHOD 507, 525.2**

Alachlor	Atraton
Bromacil	Butylate
Chlorpropham	Molinate
Propyzamide	Tetrachlorvinphos
Tricyclazole	Velpar

1000ug/mL in Acetone (10 components)
M-PM5074B4-1ML 1 mL Ampule

NITROGEN/PHOSPHOROUS PESTICIDES MIXTURE

#5 - EPA METHOD 507, 525.2

Dichlorvos	Fenarimol
Fluridone	Napropamide
s-Propyl butylethylthiocarbamate	Simetryn
Tebuthiuron	Terbacil
1000ug/mL in tert-Butylmethyl ether	(8 components)
M-PM5075T4-1ML	1 mL Ampule

SIMAZINE STANDARDS

EPA METHOD 507, 8141A

N-13800-500MG	Simazine
S-13800M1-1ML	100ug/mL in Methanol
S-13800M1-5ML	100ug/mL in Methanol
S-13800B4-1ML	1000ug/mL in Acetone
S-13800B4-5ML	1000ug/mL in Acetone

SURROGATE STANDARDS

EPA METHOD 507

N-10114-1G	1,3-Dimethyl-2-nitrobenzene
S-10114T5-1ML	2000ug/mL in tert-Butylmethyl ether
S-10114T5-5ML	2000ug/mL in tert-Butylmethyl ether
S-10114T6-1ML	2500ug/mL in tert-Butylmethyl ether
S-10114T6-5ML	2500ug/mL in tert-Butylmethyl ether

LABORATORY PERFORMANCE CHECK MIXTURE

EPA METHOD 507

15ug/mL	Atrazine
500ug/mL	Bromacil
250ug/mL	1,3-Dimethyl-2-nitrobenzene
30ug/mL	Prometon
250ug/mL	Triphenyl phosphate
5ug/mL	Vernolate

Varied concentration in tert-Butylmethyl ether (6 components)
M-LPC507T99-1ML 1 mL Ampule

INTERNAL STANDARDS

EPA METHOD 507

N-13708-1G	Triphenyl phosphate
S-13708T3-1ML	500ug/mL in tert-Butylmethyl ether
S-13708T3-5ML	500ug/mL in tert-Butylmethyl ether
N-12690-1G	o-Nitrotoluene
S-12690A4-1ML	1000ug/mL in Acetonitrile
S-12690A4-5ML	1000ug/mL in Acetonitrile
S-12690T6-1ML	2500ug/mL in tert-Butylmethyl ether
S-12690T6-5ML	2500ug/mL in tert-Butylmethyl ether

METHOD 508 - CHLORINATED PESTICIDES AND PCBs

Method 508 is applicable to the determination of chlorinated pesticides in ground water and finished drinking water by gas chromatography with an electron capture detector.

Method 508.1 is applicable for the determination of twenty-nine chlorinated pesticides, three herbicides and four organohalides in drinking water, ground water, and drinking water in any treatment stage using liquid-solid extraction and electron capture gas chromatography.

Method 508A is applicable for screening of Polychlorinated Biphenyls by gas chromatography and may be used for screening finished drinking water, raw source water, or drinking water in any treatment stage for polychlorinated biphenyls (PCBs).

ORGANOCHLORINE PESTICIDES MIXTURE #1

EPA METHOD 508, 508.1, 525.2, 608, 617, 625, 8080A/8081, 8250A/8270B, CLP

Aldrin (TM)	b-Endosulfan	BHC (alpha isomer)	BHC (beta isomer)
BHC (delta isomer)	cis-Chlordane	4,4'-DDD	4,4'-DDE
4,4'-DDT	Dieldrin	α-Endosulfan	Endosulfan sulfate
Endrin	Endrin aldehyde	Endrin ketone	Heptachlor
Heptachlor epoxide (Isomer B)	Lindane (BHC gamma isomer)	Methoxychlor	trans-Chlordane
1000ug/mL in tert-Butyl methyl ether	(20 components)		
M-OCP5081T4-1ML	1 mL Ampule		

ORGANOCHLORINE PESTICIDES MIXTURE #2

EPA METHOD 508/508.1, 525.2

Chlorobenzilate	Chloroneb
Chlorothalonil	Chlorthal-dimethyl
Etridiazole	Hexachlorobenzene
Permethrin	Propachlor
Trifluralin	
1000ug/mL in tert-Butyl methyl ether	(9 components)
M-OCP5082T4-1ML	1 mL Ampule

ORGANOCHLORINE PESTICIDES MIXTURE #3

EPA METHOD 508.1

Alachlor	Atrazine
Bladex	Hexachlorocyclopentadiene
Metolachlor	Metribuzin
Simazine	
500 ug/ml in Ethyl acetate	(7 components)
M-OCP5083H3-1ML	1 mL Ampule

ORGANOCHLORINE PESTICIDES MIXTURES

EPA METHOD 508/508.1, 608, 617, 625, 8080A/8081, 8250A/8270B, CLP

Aldrin (TM)	b-Endosulfan	BHC (alpha isomer)	BHC (beta isomer)
BHC (delta isomer)	4,4'-DDD	4,4'-DDE	4,4'-DDT
Dieldrin	α-Endosulfan	Endosulfan sulfate	Endrin
Endrin aldehyde	Heptachlor	Heptachlor epoxide (Isomer B)	
Lindane (BHC gamma isomer)	Methoxychlor		
100 ug/ml in Toluene:Hexane (50:50)		(17 components)	
M-OCP8080AC1-1ML		1 mL Ampule	
100 ug/ml in Toluene:Hexane (50:50), less Methoxychlor		(16 components)	
M-PPO8AC1-1ML		1 mL Ampule	
2000 ug/ml in Toluene:Hexane (50:50), less Methoxychlor		(16 components)	
M-PPHC5AC5-1ML		1 mL Ampule	

LABORATORY PERFORMANCE CHECK MIXTURE

EPA METHOD 508/508.1

40ug/mL	BHC (delta isomer)
50ug/mL	Chlorothalonil
2ug/mL	Chlorpyrifos
50ug/mL	Chlorthal-dimethyl

Varied concentration in tert-Butylmethyl ether (4 components)
M-LPC508T99-1ML 1mL Ampule

2-CHLOROBIPHENYL STANDARDS

EPA METHOD 508A

BZ-1-50MG	2-Chlorobiphenyl
BZ-1J1-2ML	100ug/mL in Hexane

SURROGATE STANDARDS

EPA METHOD 508/508.1

N-10877-500MG	4,4'-Dibromobiphenyl
S-10877H5-1ML	2000ug/mL in Ethyl acetate
S-10877H5-5ML	2000ug/mL in Ethyl acetate
S-10877J1-1ML	100ug/mL in Hexane
S-10877J1-5ML	100ug/mL in Hexane
S-10877X4-1ML	1000ug/mL in Methylene chloride
S-10877X4-5ML	1000ug/mL in Methylene chloride

INTERNAL STANDARDS

EPA METHOD 508/508.1

N-12830-250MG	Pentachloronitrobenzene
S-12830M1-1ML	100ug/mL in Methanol
S-12830M1-5ML	100ug/mL in Methanol
S-12830K4-1ML	1000ug/mL in Isooctane
S-12830K4-5ML	1000ug/mL in Isooctane
S-12830H5-1ML	2000ug/mL in Ethyl acetate
S-12830H5-5ML	2000ug/mL in Ethyl acetate
S-12830K7-1ML	5000ug/mL in Isooctane
S-12830K7-5ML	5000ug/mL in Isooctane

DEGRADATION CALIBRATION MIXTURE

EPA METHOD 508/508.1, 525.2, 608, 617, 625, 1618,
1656, 8080A/8081, 8250A/8270B, CLP

200ug/mL	4,4'-DDT
100ug/mL	Endrin

Varied concentration in tert-Butylmethyl ether (2 components)
M-DC508T99-1ML 1mL Ampule

DEGRADATION PRODUCTS MIXTURE

EPA METHOD 508/508.1, 525.2, 608, 617, 625, 1618,
1656, 8080A/8081, 8250A/8270B, CLP

200ug/mL	4,4'-DDD
200ug/mL	4,4'-DDE
100ug/mL	Endrin aldehyde
100ug/mL	Endrin ketone

Varied concentration in tert-Butylmethyl ether (4 components)
M-DP508T99-1ML 1mL Ampule

DECACHLOROBIPHENYL STANDARDS

EPA METHOD 508A

BZ-209-10MG	Decachlorobiphenyl
BZ-209B3-5ML	500ug/mL in Acetone
BZ-209B3-1ML	500ug/mL in Acetone
BZ-209J1-1ML	100ug/mL in Hexane
BZ-209J4-1ML	1000ug/mL in Hexane
BZ-209U4-1ML	1000ug/mL in Toluene
BZ-209U4-5ML	1000ug/mL in Toluene

INTERNAL STANDARDS

EPA METHOD 508A

N-11092-50MG	Arochlor 1260
S-11092J4-1ML	500ug/mL in Hexane
S-11092J4-5ML	500ug/mL in Hexane
S-11092K4-1ML	500ug/mL in Isooctane
S-11092K4-5ML	500ug/mL in Isooctane
S-11092M1-1ML	500ug/mL in Methanol
S-11092M1-5ML	500ug/mL in Methanol
S-11092M5-1ML	500ug/mL in Methanol
S-11092M5-5ML	500ug/mL in Methanol
S-11092V1-1ML	500ug/mL in Transformer oil
S-11092V1-5ML	500ug/mL in Transformer oil

METHOD 509 - ETHYLENE THIOUREA

Method 509 is applicable for the determination of Ethylene Thiourea (ETU) in water using gas chromatography with a Nitrogen-Phosphorus Detector.

ETHYLENE THIOUREA STANDARDS

EPA METHOD 509

N-10381-250MG	2-Imidazolidinethione
S-10381A7-1ML	5000ug/mL in Acetonitrile
S-10381A7-5ML	5000ug/mL in Acetonitrile
S-10381H1-1ML	100ug/mL in Ethyl acetate containing 0.1% DTT
S-10381H1-5ML	100ug/mL in Ethyl acetate containing 0.1% DTT

INTERNAL STANDARDS

EPA METHOD 509

N-10790-1G	3,4,5,6-Tetrahydro-2-pyrimidinethiol
S-10790H1-1ML	100ug/mL in Ethyl acetate containing 0.1% DTT
S-10790H1-5ML	100ug/mL in Ethyl acetate containing 0.1% DTT

FREE RADICAL SCAVENGER STANDARDS

EPA METHOD 509

N-11823-100MG	Dithiothreitol
S-11823H4-1ML	1000ug/mL in Ethyl acetate

SURROGATE STANDARDS

EPA METHOD 509

N-13135-250MG	Propylene thiourea
S-13135H1-1ML	100ug/mL in Ethyl acetate containing 0.1% DTT
S-13135H1-5ML	100ug/mL in Ethyl acetate containing 0.1% DTT

METHOD 515.1/515.2/513.3 - CHLORINATED ACIDS

Method 515.1 is applicable to the determination of certain chlorinated acids in ground water and finished drinking water.

Method 515.2 is a gas chromatographic (GC) method applicable to the determination of chlorinated acids in ground water, and finished drinking water, using liquid-solid extraction, derivation and GC with an electron capture detector.

Method 515.3 is a gas chromatographic method applicable to the determination of chlorinated acids in drinking water, ground water, raw source water, and water at any intermediate stage by liquid-liquid extraction, derivatization and gas chromatography with electron capture detection.

CHLORINATED ACIDS MIXTURE**EPA METHOD 515.1**

100ug/mL	Acifluorfen	200ug/mL	Bentazon	100ug/mL	Chloramben
200ug/mL	2,4-D	800ug/mL	2,4-DB	1300ug/mL	Dalapon
100ug/mL	Dicamba	100ug/mL	3,5-Dichlorobenzoic acid	300ug/mL	Dichlorprop
200ug/mL	Dinoseb	100ug/mL	4-Nitrophenol	100ug/mL	Pentachlorophenol
100ug/mL	Picloram	100ug/mL	Silvex	100ug/mL	2,4,5-T (TM)
20ug/mL	Tetrachloroterephthalic acid				

Varied concentration in Acetone (16 components)
M-CA515B99-1ML 1mL Ampule

CHLORINATED ACIDS MIXTURE #2**EPA METHOD 515.1**

Acifluorfen	Bentazon	Chloramben	2,4-D
2,4-DB	Dalapon	Dicamba	3,5-Dichlorobenzoic acid
Dichlorprop	Dinoseb	4-Nitrophenol	Pentachlorophenol
Picloram	Silvex	2,4,5-T (TM)	Tetrachloroterephthalic acid

100ug/mL in Acetone (16 components)
M-CA515AB1-1ML 1mL Ampule

CHLORINATED ACIDS MIXTURE #3**EPA METHOD 515.1**

500ug/mL	Acifluorfen	1000ug/mL	Bentazon	500ug/mL	Chloramben
1000ug/mL	2,4-D	1000ug/mL	2,4-DB	1000ug/mL	Dalapon
500ug/mL	Dicamba	500ug/mL	3,5-Dichlorobenzoic acid	1000ug/mL	Dichlorprop
1000ug/mL	Dinoseb	1000ug/mL	4-Nitrophenol	100ug/mL	Pentachlorophenol
1000ug/mL	Picloram	250ug/mL	Silvex	250ug/mL	2,4,5-T (TM)
500ug/mL	Tetrachloroterephthalic acid				

Varied concentration in Acetone (16 components)
M-CA5153B99-1ML 1mL Ampule

CHLORINATED HERBICIDES MIXTURE - A**EPA METHOD 515.1, 8150B**

2,4-D	
Silvex	
2,4,5-T	
2000ug/mL in Acetone (3 components)	
M-CSHC14AB5-1ML 1mL Ampule	

CHLORINATED HERBICIDES ESTERS MIXTURE**EPA METHOD 515.1, 8150B**

2,4-D methyl ester	
Silvex methyl ester	
2,4,5-T methyl ester	
2000ug/mL in Hexane (3 components)	
M-CSHC14JB5-1ML 1mL Ampule	

METHYLATED CHLORINATED ACIDS MIXTURE**EPA METHOD 515.1**

Acifluorfen methyl ester	Bentazon methyl derivative	Chloramben methyl ester
Chlorthal-dimethyl	2,4-D methyl ester	2,4-DB methyl ester
Dalapon methyl ester	Dicamba methyl ester	3,5-Dichlorobenzoic acid methyl ester
Dichlorprop methyl ester	Dinoseb methyl ether	Pentachloroanisole
Picloram methyl ester	p-Nitroanisole	Silvex methyl ester
(2,4,5-Trichlorophenoxy)acetic acid methyl ester		

1000ug/mL in tert-Butylmethyl ether (16 components)
M-MCA515T4-1ML 1mL Ampule

SDS's

In order to meet the GHS requirements, Chem Service has upgraded MSDS's to SDS's. We ship the English language SDS with every order. We offer SDS's in many other language on our website www.chemservice.com.

METHYLATED CHLORINATED ACIDS MIXTURE #2

EPA METHOD 515.1

100ug/mL	Acifluorfen methyl ester	200ug/mL	Bentazon methyl derivative	100ug/mL	Chloramben methyl ester
20ug/mL	Chlorthal-dimethyl	200ug/mL	2,4-D methyl ester	800ug/mL	2,4-DB methyl ester
1300ug/mL	Dalapon methyl ester	100ug/mL	Dicamba methyl ester		
100ug/mL	3,5-Dichlorobenzoic acid methyl ester			300ug/mL	Dichlorprop methyl ester
200ug/mL	Dinoseb methyl ether	100ug/mL	Pentachloroanisole	100ug/mL	Picloram methyl ester
100ug/mL	p-Nitroanisole	100ug/mL	Silvex methyl ester		
100ug/mL	(2,4,5-Trichlorophenoxy)acetic acid methyl ester				

Varied concentration in tert-Butylmethyl ether

(16 components)

M-MCA515AT99-1ML

1mL Ampule

METHYLATED CHLORINATED ACIDS MIXTURE #3

EPA METHOD 515.4

50ug/mL	Acifluorfen methyl ester	100ug/mL	Bentazon methyl derivative	50ug/mL	Chloramben methyl ester
100ug/mL	2,4-D methyl ester	100ug/mL	2,4-DB methyl ester	100ug/mL	DCPA (Dacthal)
100ug/mL	Dalapon methyl ester	50ug/mL	Dicamba methyl ester		
50ug/mL	3,5-Dichlorobenzoic acid methyl ester			100ug/mL	Dichlorprop methyl ester
100ug/mL	Dinoseb methyl ester	10ug/mL	Pentachloroanisole	50ug/mL	Picloram methyl ester
50ug/mL	Quinclorac	25ug/mL	Silvex methyl ester		
25ug/mL	(2,4,5-Trichlorophenoxy)acetic acid methyl ester				

Varied concentration in tert-Butylmethyl ether

(16 components)

M-CS515AT99-1ML

1mL Ampule

METHYL DERIVATIVES MIXTURE A

EPA METHOD 515.2

100ug/mL	Chlorthal-dimethyl		
500ug/mL	3,5-Dichlorobenzoic acid methyl ester		
100ug/mL	Dichlorprop methyl ester		
200ug/mL	Dinoseb methyl ether		
100ug/mL	Pentachloroanisole		
100ug/mL	(2,4,5-Trichlorophenoxy)acetic acid methyl ester		

Varied concentration in Methanol

(6 components)

M-EPA5152AM99-1ML

1mL Ampule

METHYL DERIVATIVES MIXTURE B

EPA METHOD 515.2

200ug/mL	Acifluorfen methyl ester		
1000ug/mL	Bentazon methyl derivative		
100ug/mL	2,4-D methyl ester		
1000ug/mL	2,4-DB methyl ester		
300ug/mL	Dicamba methyl ester		
300ug/mL	Picloram methyl ester		
100ug/mL	Silvex methyl ester		

Varied concentration in Methanol

(7 components)

M-EPA5152BAM99-1ML

1mL Ampule

SURROGATE STANDARDS

EPA METHOD 515.1/515.2/515.3/515.4

N-10536-100MG	2,4'-Dichlorophenylacetic acid
S-10536B5-1ML	2000ug/mL in Acetone
S-10536B5-5ML	2000ug/mL in Acetone
S-10536T1-1ML	100ug/mL in tert-Butylmethyl ether
S-10536T1-5ML	100ug/mL in tert-Butylmethyl ether
S-10536T7-1ML	5000ug/mL in tert-Butylmethyl ether
S-10536T7-5ML	5000ug/mL in tert-Butylmethyl ether

N-10537-100MG 2,4'-Dichlorophenylacetic acid methyl ester

S-10537B1-1ML	100ug/mL in Acetone
S-10537B5-1ML	2000ug/mL in Acetone
S-10537M7-1ML	5000ug/mL in Methanol
S-10537M7-5ML	5000ug/mL in Methanol

INTERNAL STANDARDS

EPA METHOD 515.1/515.2/515.3/515.4

N-10867-100MG	4,4'-Dibromooctafluorobiphenyl
S-10867M7-1ML	5000ug/mL in Methanol
S-10867M7-5ML	5000ug/mL in Methanol
S-10867T1-1ML	100ug/mL in tert-Butylmethyl ether
S-10867T1-5ML	100ug/mL in tert-Butylmethyl ether
S-10867T5-1ML	2000ug/mL in tert-Butylmethyl ether
S-10867X5-1ML	2000ug/mL in Methylene chloride
S-10867X5-5ML	2000ug/mL in Methylene chloride

METHOD 525.2/525.3 (UCMR) - ORGANIC COMPOUNDS

Method 525.2 is a general purpose method applicable for the identification and measurement of organic compounds in source water, drinking water in any stage of treatment as well as drinking water. This method can determine a wide range of organic compounds through using a liquid phase extraction and capillary column gas chromatography and mass spectrometry method.

NITROGEN/PHOSPHOROUS PESTICIDES MIXTURE #1 - EPA METHOD 525.2, 507

Ametryne	Cyloate
Disulfoton	Tributylphosphorotrithioate
Prometon	Fenamiphos

*Inhibited with Hydroquinone

1000ug/mL in tert-Butyl methyl ether (6 components)

M-PM5071T4-1ML

1mL Ampule

NITROGEN/PHOSPHOROUS PESTICIDES MIXTURE #2 - EPA METHOD 525.2, 507

Atrazine	Diphenamid
Mevinphos	Prometryne
Propazine	Prophos
S-Ethyl dipropylthiocarbamate	Terbutryne
Triadimefon	

1000ug/mL in tert-Butyl methyl ether (9 components)

M-PM5072T4-1ML

1mL Ampule

NITROGEN/PHOSPHOROUS PESTICIDES MIXTURE #3 - EPA METHOD 525.2, 507

Butachlor	Carboxin
Diazinon	Metolachlor
Metribuzin	MGK 264 (TM)
Norflurazon	Terbufos
Vernolate	

1000ug/mL in tert-Butyl methyl ether (9 components)

M-PM5073T4-1ML

1mL Ampule

NITROGEN/PHOSPHOROUS PESTICIDES MIXTURE #4 - EPA METHOD 525.2, 507

Alachlor	Atraton
Bromacil	Butylate
Chlorpropham	Molinate
Propyzamide	Tetrachlorvinphos
Tricyclazole	Velpar

1000ug/mL in Acetone (10 components)

M-PM5074B4-1ML

1mL Ampule

ORGANOCHLORINE PESTICIDES MIXTURE #1

EPA METHOD 525.2, 508, 508.1, 608, 617, 625, 8080A/8081, 8250A/8270B, CLP

Aldrin (TM)	b-Endosulfan	BHC (alpha isomer)	BHC (beta isomer)
BHC (delta isomer)	cis-Chlordane	4,4'-DDD	4,4'-DDE
4,4'-DDT	Dieldrin	α-Endosulfan	Endosulfan sulfate
Endrin	Endrin aldehyde	Endrin ketone	Heptachlor
Heptachlor epoxide (Isomer B)	Lindane (BHC gamma isomer)	Methoxychlor	trans-Chlordane
1000ug/mL in tert-Butyl methyl ether	(20 components)		
M-OCP5081T4-1ML	1mL Ampule		

ORGANOCHLORINE PESTICIDES MIXTURE #2

EPA METHOD 525.2, 508/508.1

Chlorobenzilate	Chloroneb
Chlorothalonil	Chlorthal-dimethyl
Etridiazole	Hexachlorobenzene
Permethrin	Propachlor
Trifluralin	
1000ug/mL in tert-Butyl methyl ether	(9 components)
M-OCP5082T4-1ML	1mL Ampule

ORGANOCHLORINE PESTICIDES MIXTURE #2

EPA METHOD 525.2

Alachlor	Aldrin (TM)
Atrazine	cis-Chlordane
trans-Chlordane	Endrin
Heptachlor	Simazine
Heptachlor epoxide (Isomer B)	trans-Nonachlor
Lindane (BHC gamma isomer)	Methoxychlor
500ug/mL in Acetone	(12 components)
M-OCPh525B3-1ML	1mL Ampule

DEGRADATION CALIBRATION MIXTURE

EPA METHOD 525.2, 508/508.1, 608, 617, 625, 1618, 1656, 8080A/8081, 8250A/8270B, CLP

200ug/mL	4,4'-DDT
100ug/mL	Endrin
Varied concentration in tert-Butylmethyl ether	(2 components)
M-DC508T99-1ML	1mL Ampule

DEGRADATION PRODUCTS MIXTURE

EPA METHOD 525.2, 508/508.1, 608, 617, 625, 1618, 1656, 8080A/8081, 8250A/8270B, CLP

200ug/mL	4,4'-DDD
200ug/mL	4,4'-DDE
100ug/mL	Endrin aldehyde
100ug/mL	Endrin ketone
Varied concentration in tert-Butylmethyl ether	(4 components)
M-DP5081T99-1ML	1mL Ampule

PERFORMANCE CHECK MIXTURE

EPA METHOD 525.2

4,4'-DDT	
Decafluorotriphenylphosphine	
Endrin	
500ug/mL in Toluene -	(3 components)
M-PCM5251U3-1ML	1mL Ampule
M-PCM5251U3-5ML	5mL Ampule

DINITROTOLUENE STANDARDS

EPA METHOD 525.2

N-10643-1G	2,4-Dinitrotoluene
S-10643A4-1ML	1000ug/mL in Acetonitrile
S-10643A4-5ML	1000ug/mL in Acetonitrile
S-10643M1-1ML	100ug/mL in Methanol
S-10643M1-5ML	100ug/mL in Methanol
N-10697-1G	2,4-Dinitrotoluene
S-10697A4-1ML	1000ug/mL in Acetonitrile
S-10697M1-1ML	100ug/mL in Methanol

NITROGEN/PHOSPHOROUS PESTICIDES MIXTURE #5 - EPA METHOD 507, 525.2

Dichlorvos	Fenarimol
Fluridone	Napropamide
s-Propyl butylethylthiocarbamate	Simetryn
Tebuthiuron	Terbacil
1000ug/mL in tert-Butylmethyl ether	(8 components)
M-PM5075T4-1ML	1mL Ampule

INTERNAL STANDARDS

EPA METHOD 525.2

N-13150-10MG	Pyrene-d ₁₀
S-13150B3-1ML	500ug/mL in Acetone

SIMAZINE STANDARDS

EPA METHOD 525.2

N-13800-100MG	Simazine
S-13800M1-1ML	100ug/mL in Methanol
S-13800M1-5ML	100ug/mL in Methanol
S-13800B4-1ML	1000ug/mL in Acetone
S-13800B4-5ML	1000ug/mL in Acetone

DISULFOTON SULFONE & DISULFOTON SULFOXIDE STANDARDS

EPA METHOD 525.2

MET-11819A-100MG	Disulfoton sulfone
MET-11819B-100MG	Disulfoton sulfoxide

UCMR-4 PESTICIDE MIXTURE 1

100ug/mL	BHC (alpha isomer)
300ug/mL	Chlorpyrifos
2000ug/mL	Dimethipin
300ug/mL	Propos
500ug/mL	Oxyfluorfen
3000ug/mL	Profenofos
2000ug/mL	Tebuconazole
400ug/mL	Permethrin
700ug/mL	S,S,S-Tributyl phosphotriothioate
Varied Concentration in Methanol	(9 components)
M-UCMRPESTMIX1-1ML	1mL Ampule

METHOD 526 - SEMIVOLATILE ORGANICS

Method 526 is a general purpose method applicable for the identification and measurement of select Semi-Volatile Organic Compounds in drinking water by solid phase extraction and capillary column gas chromatography coupled with mass spectrometry (GC/MS)

SEMIVOLATILE MIXTURE**EPA METHOD 526**

Acetochlor 1,2-Diphenylhydrazine Prometon	Bladex Disulfoton Terbufos	Diazinon Fonofos 2,4,6-Trichlorophenol	2,4-Dichlorophenol Nitrobenzene
1000ug/mL in Acetone M-CSEPA526PDSB4-1ML	(11 components) 1mL Ampule		

METHOD 527 - PESTICIDES AND FLAME RETARDANTS

Method 527 is a general purpose method applicable for the identification and measurement of select pesticides and flame retardants in drinking water by solid phase extraction and capillary column gas chromatography coupled with mass spectrometry (GC/MS)

PESTICIDES MIXTURE**EPA METHOD 527**

Atrazine Mirex Bifenthrin Nitrofen 500ug/mL in Ethyl acetate M-CS5274H3-1ML	Esbjöl Norflurazon Bromacil Oxychlorthane (16 components) 1mL Ampule	Esfenvalarate Prometryne Fenvalarate Propazine	Hexazinone Thiobencarb Kepone Vinclozolin
--	---	---	--

METHOD 531 - CARBAMATE PESTICIDES

Method 531 is a general purpose method applicable for the identification and measurement of select N-methylcarbamoyloximes and N-methylcarbamates in finished drinking water by direct aqueous injection high performance liquid chromatography (HPLC) with post-column derivatization.

LABORATORY PERFORMANCE CHECK MIXTURE**EPA METHOD 531.1**

Acetonitrile Aldicarb sulfoxide BDMC 3-Hydroxycarbofuran Methiocarb	
Varied concentration in Acetonitrile M-LPC531A99-1ML	(5 components) 1mL Ampule

BDMC STANDARDS**EPA METHOD 531**

N-11131-100MG	BDMC
S-11131M1-1ML	100ug/mL in Methanol
S-11131M1-5ML	100ug/mL in Methanol

CARBAMATE PESTICIDE CALIBRATION MIXTURE**EPA METHOD 531.2**

Aldicarb Oxamyl Methomyl 3-Hydroxyfuran	Carbaryl Aldicarb sulfone Propoxur 1-Naphthol	Methiocarb Carbofuran Aldicarb sulfoxide
100ug/mL in Acetonitrile M-CS5312A1-1ML	(11 components) 1mL Ampule	

METHOD 532 - PHENYLUREA PESTICIDES

Method 532 is a general purpose method applicable for the identification and measurement of select phenylurea compounds in drinking water by solid phase extraction and high performance liquid chromatography (HPLC) with UV detection.

PHENYLUREA PESTICIDES MIXTURE**EPA METHOD 532**

DiFlubenzuron Linuron Tebuthiuron	Diuron Propanil Thidiazuron	Fluometuron Siduron
200ug/mL in Methanol:Acetonitrile M-ESTDAK2-1ML	(8 components) 1mL Ampule	

METHOD 535 - DEGRADATES OF CHLOROACETANILIDE

Method 535 is applicable to the determination of ethanesulfonic acid (ESA) and oxanilic acid (OA) degradates of chloroacetanilide and other acetamide herbicides in drinking water by solid phase extraction and liquid chromatography/tandem mass spectrometry (LC/MS/MS).

ETHANESULFONIC ACID STANDARDS**EPA METHOD 535**

MET-11013A-25MG	Acetochlor ESA sodium salt
MET-11013AM1-1ML	100ug/mL in Methanol
MET-11043D-25MG	Alachlor ESA sodium salt
MET-11043DM1-1ML	100ug/mL in Methanol
MET-11353A-25MG	Butachlor ESA sodium salt
MET-11353AM1-1ML	100ug/mL in Methanol
MET-11754A-25MG	Dimethachlor ESA sodium salt
MET-11754AM1-1ML	100ug/mL in Methanol
MET-11755A-25MG	Dimethenamid ESA sodium salt
MET-11755AM1-1ML	100ug/mL in Methanol
MET-11986A-25MG	Flufenacet ESA sodium salt
MET-11986AM1-1ML	100ug/mL in Methanol
MET-12478C-25MG	Metolachlor ESA sodium salt
MET-12478CM1-1ML	100ug/mL in Methanol
MET-13104A-25MG	Propachlor ESA sodium salt
MET-13104AM1-1ML	100ug/mL in Methanol

OXANILIC ACID STANDARDS**EPA METHOD 535**

MET-11013B-10MG	Acetochlor OA
MET-11013BM1-1ML	100ug/mL in Methanol
MET-11043E-10MG	Alachlor OA
MET-11043EM1-1ML	100ug/mL in Methanol
MET-11755B-10MG	Dimethenamid OA
MET-11755BM1-1ML	100ug/mL in Methanol
MET-11986B-10MG	Flufenacet OA
MET-11986BM1-1ML	100ug/mL in Methanol
MET-12478D-10MG	Metolachlor OA
MET-12478DM1-1ML	100ug/mL in Methanol
MET-13104B-10MG	Propachlor OA
MET-13104BM1-1ML	100ug/mL in Methanol

METHOD 547 - GLYPHOSATE STANDARDS

Method 547 is applicable for the identification and measurement of N-methylcarbamoyloximes and N-methylcarbamates in ground water and finished drinking water, by direct aqueous injection HPLC, post-column derivatization and fluorescence detection.

GLYPHOSATE STOCK STANDARDS**EPA METHOD 547**

S-12133F4-1ML	Glyphosate
	1000ug/mL in Water

METHOD 548, 548.1 - ENDOTHAL STANDARDS

Method 548 is applicable in the N-Methylcarbamoyloximes and N-Methylcarbamates in ground water and finished drinking water, by aqueous derivatization, liquid-solid extraction, and gas chromatography with electron-capture detection. Method 548.1 is applicable for the identification and simultaneous measurement of Endothal in drinking water sources and finished drinking water by ion exchange extraction, acidic methanol methylation, and gas chromatography/mass spectrometry.

ENDOTHAL STOCK STANDARDS**EPA METHOD 548/548.1**

N-11852-100MG	Endothal monohydrate
S-11852A1-1ML	100ug/mL in Acetonitrile
S-11852A1-5ML	100ug/mL in Acetonitrile
S-11852F0-1ML	50ug/mL in Water
S-11852F10-1ML	10ug/mL in Water
S-11852F10-5ML	10ug/mL in Water

ENDOTHAL-PFPH DERIVATIVE STANDARDS**EPA METHOD 548**

N-11853-100MG	Endothal-PFPH
S-11853T1-1ML	100ug/mL in tert-Butylmethyl ether
S-11853T1-5ML	100ug/mL in tert-Butylmethyl ether

ENDOTHAL METHYL DERIVATIVE STANDARDS**EPA METHOD 548.1**

N-11764-10MG	Dimethyl Endothal
S-11853U1-1ML	100ug/mL in Toluene
S-11853U1-5ML	100ug/mL in Toluene

METHOD 549 - DIQUAT & PARAQUAT

Methods 549 is an HPLC method applicable to the determination of Diquat (1,1'-Ethylene-2,2'-bipyridilium dibromide salt) and Paraquat (1,1'-dimethyl-4,4'-bipyridilium dichloride salt) in drinking water sources and finished drinking water by liquid-solid extraction and HPLC with ultraviolet detection.

DIQUAT STOCK STANDARDS**EPA METHOD 549.2**

N-11816-500MG	Diquat dibromide monohydrate
S-11816F1-1ML	100ug/mL in Water
S-11816F13-1ML	1970ug/mL in Water
S-11816F13-5ML	1970ug/mL in Water

PARAQUAT STOCK STANDARDS**EPA METHOD 549.2**

N-12818-500MG	Paraquat dichloride hydrate
S-12818F1-1ML	100ug/mL in Water
S-12818F17-1ML	1770ug/mL in Water
S-12818F17-5ML	1770ug/mL in Water

METHOD 555 - CHLORINATED ACIDS

Method 555 is a general purpose method applicable for the identification and measurement of select chlorinated acids in ground water and finished drinking water by high performance liquid chromatography (HPLC) coupled with a photodiode array ultraviolet detector.

CHLORINATED ACIDS HPLC, MIXTURE A**EPA METHOD 555**

Acifluorfen
Bentazon
Chloramben
2,4-D

Dicamba
Dichlorprop
Picloram
2,4,5-TP

1000ug/ml in Acetonitrile
M-CS5551A4-1ML

(8 components)
1 mL Ampule

CHLORINATED ACIDS HERBICIDE MIXTURE**EPA METHOD 555**

2,4-Dichlorophenoxyacetic acid

2,4,5-TP

1000ug/ml in Acetonitrile
M-CS5553A4-1ML

(2 components)
1 mL Ampule

CHLORINATED ACIDS HPLC, MIXTURE B**EPA METHOD 555**

2,4-DB
3,5-Dichlorobenzoic acid
Dinoseb
MCPP

MCPP
4-Nitrophenol
Pentachlorophenol
2,4,5-T

1000ug/ml in Acetonitrile
M-CS5552A4-1ML

(8 components)
1 mL Ampule

CHEM SERVICE_{INC.}

Certified Analytical Reference Standards
and Chemicals in Small Quantities

Custom
solutions and mixtures

*RELIABLE**PROFESSIONAL***Certificate of Analysis with Each Custom Product****GRAVAMETRIC CERTIFICATION**

By weight - Guarantees the analytes in the custom standard to be made within +/-5% of specifications.

or**QUANTITATIVE CERTIFICATION**

With a chromatogram - Includes verification of the concentration and peak identification.

Chem Service provides custom solutions and mixtures to meet your special analytical needs. We provide accurate, precise and economical custom standards mix. Certified neat reference standards are used to prepare mixtures. Balances used to prepare stock and custom standards are calibrated to NIST traceable weights.

METHOD 608 - ORGANOCHLORINE PESTICIDES

Method 608 is applicable to the determination of organochlorine pesticides and PCBs by gas chromatography with measurement by an electron capture detector.

Method 608.1 is applicable to the determination of organochlorine pesticides in industrial and municipal wastewater by gas chromatography with measurement by an electron capture detector.

Method 608.2 is applicable to the determination of organochlorine pesticides in industrial and municipal wastewater by extraction, analysis by gas chromatography, and measurement by an electron capture detector.

ORGANOCHLORINE PESTICIDES MIXTURES**EPA METHOD 608, 508/508.1, 617, 625, 8080A/8081, 8250A/8270B, CLP**

Aldrin (TM)	b-Endosulfan	BHC (alpha isomer)	BHC (beta isomer)
BHC (delta isomer)	4,4'-DDD	4,4'-DDE	4,4'-DDT
Dieldrin	α-Endosulfan	Endosulfan sulfate	Endrin
Endrin aldehyde	Heptachlor	Heptachlor epoxide (Isomer B)	
Lindane (BHC gamma isomer)	Methoxychlor		

100 ug/ml in Toluene:Hexane (50:50) (17 components)

M-OCP8080AC1-1ML**1mL Ampule**

100 ug/ml in Toluene:Hexane (50:50), less Methoxychlor

(16 components)

M-PPO8AC1-1ML**1mL Ampule**

2000 ug/ml in Toluene:Hexane (50:50), less Methoxychlor

(16 components)

M-PPHC5AC5-1ML**1mL Ampule****DEGRADATION CALIBRATION MIXTURE****EPA METHOD 608, 508/508.1, 525.2, 617, 625, 1618, 1656, 8080A/8081, 8250A/8270B, CLP**

200ug/mL	4,4'-DDT
100ug/mL	Endrin

Varied concentration in tert-Butylmethyl ether (2 components)

M-DC508T99-1ML **1mL Ampule****TUNING STANDARDS****EPA METHOD 608**

N-11571-10MG	Decafluorotriphenylphosphine
S-11571B0-1ML	50ug/mL in Acetone
S-11571B0-5ML	50ug/mL in Acetone
S-11571X0-1ML	50ug/mL in Methylene chloride
S-11571X5-1ML	2000ug/mL in Methylene chloride
S-11571X5-5ML	2000ug/mL in Methylene chloride

DEGRADATION PRODUCTS MIXTURE**EPA METHOD 608, 508/508.1, 525.2, 617, 625, 1618, 1656, 8080A/8081, 8250A/8270B, CLP**

200ug/mL	4,4'-DDD
200ug/mL	4,4'-DDE
100ug/mL	Endrin aldehyde
100ug/mL	Endrin ketone

Varied concentration in tert-Butylmethyl ether (4 components)

M-DP5081T99-1ML **1mL Ampule****PESTICIDES - CONTROL SAMPLE MIXTURE****EPA METHOD 608, 8080A/8081**

20ug/mL	Aldrin	20ug/mL	BHC (alpha isomer)
20ug/mL	BHC (beta isomer)	100ug/mL	4,4'-DDD
20ug/mL	4,4'-DDE	100ug/mL	4,4'-DDT
20ug/mL	Dieldrin	20ug/mL	α-Endosulfan
100ug/mL	b-Endosulfan	100ug/mL	Endosulfan sulfate
100ug/mL	Endrin	20ug/mL	Heptachlor
20ug/mL	Heptachlor epoxide (Isomer B)		

Varied concentration in Toluene (13 components)

M-CSM8080U99-1ML **1mL Ampule****ORGANOCHLORINE PESTICIDES MIXTURE #1****EPA METHOD 608.1**

200ug/mL	Chlorobenzilate
40ug/mL	Chloroneb
200ug/mL	Chloropropylate
40ug/mL	1,2-Dibromo-3-chloropropane
40ug/mL	Etridiazole
60ug/mL	PCNB
1000ug/mL	Propachlor

Varied concentration in Isooctane (7 components)

M-CS6081K99-1ML **1mL Ampule****ORGANOCHLORINE PESTICIDES MIXTURE****EPA METHOD 608.2**

Chlorthal-dimethyl	Chlorothalonil
cis-Permethrin	trans-Permethrin
Methoxychlor	Dicloran

100ug/mL in Isooctane (6 components)

M-EPA6082OCPK1-1ML **1mL Ampule****ORGANOCHLORINE PESTICIDES MIXTURE #2****EPA METHOD 608.2**

1ug/mL	Chlorothalonil
3ug/mL	DCPA
2ug/mL	Dichloran
40ug/mL	Methoxychlor
400ug/mL	Permethrin

Varied concentration in Isooctane (5 components)

M-CS6082K99-1ML **1mL Ampule**

METHOD 610 - POLYNUCLEAR AROMATIC HYDROCARBONS

Method 610 is applicable for determination of certain polynuclear aromatic hydrocarbons (PAH) in municipal and industrial discharges. This method provides for both high performance liquid chromatographic (HPLC) and gas chromatographic (GC) approaches for the determination of PAHs.

**POLYNUCLEAR AROMATIC HYDROCARBONS (PAH)
EPA METHOD 610, 550.1, 8100, 8270B, 8310, CLP**

Acenaphthene Benzo(b)fluoranthene Chrysene Indeno(1,2,3-C,D)pyrene	Acenaphthylene Benzo(k)fluoranthene 1,2:5,6-Dibenzanthracene Naphthalene	Anthracene 1,12-Benzoperylene Fluoranthene Phenanthrene	1,2-Benzanthracene Benzo(a)pyrene Fluorene Pyrene
100ug/mL in Methanol - M-PPH10M1-1ML M-PPH10M1-5ML	(16) components 1mL Ampule 5mL Ampule		
200ug/mL in Acetonitrile - M-PNA550A2-1ML	(16) components 1mL Ampule		
2000 ug/mL in CH ₂ Cl ₂ :Benzene (50:50) - M-PPHC6AD5-1ML M-PPHC6AD5-5ML	(16) components 1mL Ampule 5mL Ampule		

METHOD 614, 614.1 - ORGANOPHOSPHOROUS PESTICIDES

Method 614 is applicable for determination of organophosphorus pesticides in industrial and municipal wastewater by gas chromatography. Method 614 represents a revision of a previously promulgated U.S. EPA method for organophosphorus pesticides. This method has been in widespread use since its promulgation, and represents the state of the art for the analysis of organophosphorus pesticides.

Method 614.1 is applicable for determination of organophosphorus pesticides in municipal and industrial wastewater by gas chromatography with measurement accomplished by a nitrogen/ phosphorous specific detector.

ORGANOPHOSPHOROUS PESTICIDES MIXTURE #1**EPA METHOD 614**

Demeton S Diazinon	Disulfoton Ethion	Guthion® Malathion	Methyl parathion Parathion®
1000ug/mL Hexane - M-OPP6141J4-1ML	(8 components) 1mL Ampule		

ORGANOPHOSPHOROUS PESTICIDES MIXTURE #2**EPA METHOD 614**

Dioxathion	EPN	Ethion	Terbufos
1000ug/mL Hexane - M-OPP6142J4-1ML	(4 components) 1mL Ampule		

METHOD 615 - CHLORINATED HERBICIDES

Method 615 is applicable for determination of chlorinated herbicides, as well as their salts and esters, through extraction and by electron capture gas chromatography (GC/ECD).

CHLORINATED HERBICIDES MIXTURE**EPA METHOD 615, 8151**

Acifluorfen Dicamba Picloram 2,4-D Mecoprop	2,4-DB Pentachlorophenol 4-Chloro-o-toloxoacetic acid Dinoseb Tetrachloroterephthalic acid	4-Nitrophenol Chloramben Dichlorprop 2,4,5-T®	Bentazon 3,5-Dichlorobenzoic acid Silvex Dalapon
1000ug/mL in Acetone M-CH8151B4-1ML	(18 components) 1mL Ampule		

METHYLATED CHLORINATED HERBICIDES MIXTURE**EPA METHOD 615, 8151**

Acifluorfen methyl ester Dichlorprop methyl ester 4-Chloro-o-toloxoacetic acid methyl ester p-Nitroanisole Dalapon methyl ester Silvex methyl ester	3,5-Dichlorobenzoic acid methyl ester Chloramben methyl ester Mecoprop methyl ester 2,4-D methyl ester Picloram methyl ester Dicamba methyl ester	Bentazon methyl derivative Dinoseb methyl ether Chlorthal Pentachloroanisole 2,4-DB methyl ester (2,4,5-Trichlorophenoxy)acetic acid methyl ester
1000ug/mL in Isooctane:Acetone (90:10) M-MCH8151Y4-1ML	(18 components) 1mL Ampule	

CHLORINATED HERBICIDES - CONTROL SAMPLE MIXTURE**EPA METHOD 615, 8150B**

10,000ug/mL	4-Chloro-o-tolxyloxyacetic acid	100ug/mL	2,4-D
250ug/mL	Dalapon	100ug/mL	2,4-DB
10ug/mL	Dicamba	100ug/mL	Dichlorprop
50ug/mL	Dinoseb	10,000ug/mL	Mecoprop
10ug/mL	Silvex	10ug/mL	2,4,5-T®

Varied concentration in Acetone -

(10 components)

M-CSM81501B99-1ML**1mL Ampule****METHYLATED CHLORINATED HERBICIDES - CONTROL SAMPLE MIXTURE****EPA METHOD 615, 8150B**

10000ug/mL	4-Chloro-o-tolxyloxyacetic acid methyl ester	100ug/mL	2,4-D methyl ester
250ug/mL	Dalapon methyl ester	100ug/mL	2,4-DB methyl ester
10ug/mL	Dicamba methyl ester	100ug/mL	Dichlorprop methyl ester
50ug/mL	Dinoseb methyl ether	10000ug/mL	Mecoprop methyl ester
10ug/mL	Silvex methyl ester		
10ug/mL	(2,4,5-Trichlorophenoxy)acetic acid methyl ester		

Varied concentration in Hexane -

(10 components)

M-MCSM81501J99-1ML**1mL Ampule****METHOD 617 - ORGANOHALIDE PESTICIDES**

Method 617 is applicable to the separation and measurement of organohalide pesticides and PCBs in industrial and municipal discharges by electron capture gas chromatography.

ORGANOCHLORINE PESTICIDES MIXTURES**EPA METHOD 617, 608, 508/508.1, 625, 8080A/8081, 8250A/8270B, CLP**

Aldrin (TM)	b-Endosulfan	BHC (alpha isomer)	BHC (beta isomer)
BHC (delta isomer)	4,4'-DDD	4,4'-DDE	4,4'-DDT
Dieldrin	a-Endosulfan	Endosulfan sulfate	Endrin
Endrin aldehyde	Heptachlor	Heptachlor epoxide (Isomer B)	
Lindane (BHC gamma isomer)	Methoxychlor		

100 ug/mL in Toluene:Hexane (50:50)

(17 components)

M-OCF8080AC1-1ML**1mL Ampule****DEGRADATION CALIBRATION MIXTURE****EPA METHOD 617, 608, 508/508.1, 525.2, 625, 1618, 1656, 8080A/8081, 8250A/8270B, CLP**

200ug/mL	4,4'-DDT
100ug/mL	Endrin

Varied concentration in tert-Butylmethyl ether

(2 components)

M-DC508T99-1ML**1mL Ampule****DEGRADATION PRODUCTS MIXTURE****EPA METHOD 617, 608, 508/508.1, 525.2, 625, 1618, 1656, 8080A/8081, 8250A/8270B, CLP**

200ug/mL	4,4'-DDD
200ug/mL	4,4'-DDE
100ug/mL	Endrin aldehyde
100ug/mL	Endrin ketone

Varied concentration in tert-Butylmethyl ether

(4 components)

M-DP5081T99-1ML**1mL Ampule****METHOD 619 - TRIAZINE PESTICIDES**

Method 619 is applicable to triazine pesticides' extraction, separation and measurement in industrial and municipal wastewaters by gas chromatography with a thermionic bead detector in the nitrogen mode.

TRIAZINE PESTICIDES MIXTURE**EPA METHOD 619**

Ametryne	Atraton	Atrazine	Prometon	Prometryne
Propazine	Simetryn	Simazine	Terbutylazine	Terbutryne

500ug/mL Acetone -

(10 components)

M-TP6191B3-1ML**1mL Ampule****M-TP6191B3-5ML****5mL Ampule****ISO Certified**

The certification of the Chem Service products is through the ISO Guide 34, ISO 9001 and ISO 17025, Quality Management System. It is applicable to the design, development, production, distribution and servicing of organic neat and synthetic reference materials.

METHOD 622 - ORGANOPHOSPHOROUS PESTICIDES

Method 622 is applicable to the determination of organophosphorus pesticides in industrial and municipal wastewaters by extraction, separation and measurement by gas chromatography with a thermionic bead, or flame photometric detector, in the phosphorus mode.

Method 622.1 is applicable to the determination, separation, and measurement of thiophosphate pesticides in industrial and municipal wastewater discharge, as provided under 40 CFR 136.1, by alkali flame detector gas chromatography (GC/AFD).

ORGANOPHOSPHOROUS PESTICIDES MIXTURE #2**EPA METHOD 622**

Dichlorvos
Mevinphos
Tetrachlorvinphos

1000ug/mL in Acetone (3 components)
M-OPP6221B4-1ML 1mL Ampule

ORGANOPHOSPHOROUS PESTICIDES MIXTURE #4**EPA METHOD 622**

Naled
Dichlorvos
Tetrachlorvinphos

1000ug/mL in Hexane (3 components)
M-OPP6223J4-1ML 1mL Ampule

ORGANOPHOSPHOROUS PESTICIDES MIXTURE #3**EPA METHOD 622**

Chlorpyrifos Chlorpyrifos Methyl
Diazinon Fenchlorphos
Methyl parathion Prophos
Tributylphosphorotrithioate

1000ug/mL in Hexane (7 components)
M-OPP6222J4-1ML 1mL Ampule

MERPPOS CALIBRATION STANDARDS**EPA METHOD 622, 507, 1657, 8140/8141A, 8321**

N-13194-250MG S,S,S-Tributyl phosphorotrithioate
S-13194A1-1ML 100ug/mL in Acetonitrile
S-13194T1-1ML 100ug/mL in tert-Butylmethylether

METHOD 625/1625 - SEMI-VOLATILE TOXIC ORGANICS

Method 625 is used to determine the semivolatile toxic organic pollutants associated with the Clean Water Act (as amended 1987), the Resource Conservation and Recovery Act (as amended 1986), the Comprehensive Environmental Response, Compensation and Liability Act (as amended 1986), and other (base/neutrals and acids) compounds amenable to being partitioned into an organic solvent and analyzed by GC/MS.

Method 1625 (Revision B 1991) is designed to determine the semivolatile toxic organic pollutants associated with the 1976 Consent Decree, and additional compounds amenable to extraction and analysis by capillary column gas chromatography/mass spectrometry (GC/MS). The method is designed to meet the survey requirements of Effluent Guidelines Division (EGD) and the National Pollutants Discharge Elimination System (NPDES) under 40 CFR 138.1

DEGRADATION CALIBRATION MIXTURE**EPA METHOD 625, 608, 508/508.1, 525.2, 617, 1618, 1656, 8080A/8081, 8250A/8270B, CLP**

200ug/mL 4,4'-DDT
100ug/mL Endrin

Varied concentration in tert-Butylmethyl ether (2 components)
M-DC508T99-1ML 1mL Ampule

TUNING STANDARDS MIXTURE**EPA METHOD 625, 8250A/8270B, CLP**

Benzidine
4,4'-DDT
Decafluorotriphenylphosphine
Pentachlorophenol

1000ug/mL in Methylene chloride (4 components)
M-CLPTS1X4-1ML 1mL Ampule
M-CLPTS1X4-5ML 5mL Ampule

DEGRADATION PRODUCTS MIXTURE**EPA METHOD 625, 608, 508/508.1, 525.2, 617, 1618, 1656, 8080A/8081, 8250A/8270B, CLP**

200ug/mL 4,4'-DDD
200ug/mL 4,4'-DDE
100ug/mL Endrin aldehyde
100ug/mL Endrin ketone

Varied concentration in tert-Butylmethyl ether (4 components)
M-DP5081T99-1ML 1mL Ampule

ORGANOCHLORINE PESTICIDES MIXTURES**EPA METHOD 625, 608, 508/508.1, 617, 8080A/8081, 8250A/8270B, CLP**

Aldrin (TM)	b-Endosulfan	BHC (alpha isomer)	BHC (beta isomer)
BHC (delta isomer)	4,4'-DDD	4,4'-DDE	4,4'-DDT
Dieldrin	α-Endosulfan	Endosulfan sulfate	Endrin
Endrin aldehyde	Heptachlor	Heptachlor epoxide (Isomer B)	
Lindane (BHC gamma isomer)	Methoxychlor		

100 ug/ml in Toluene:Hexane (50:50) (17 components)
M-OCF8080AC1-1ML 1mL Ampule
100 ug/ml in Toluene:Hexane (50:50), less Methoxychlor (16 components)
M-PPO8AC1-1ML 1mL Ampule
2000 ug/ml in Toluene:Hexane (50:50), less Methoxychlor (16 components)
M-PPHC5AC5-1ML 1mL Ampule

PERFORMANCE STANDARDS

EPA METHOD 625

N-11158-100MG	Benzidine
S-11158G1-1ML	100ug/mL in Ethanol
S-11158G1-5ML	100ug/mL in Ethanol
N-10876-100MG	4,4'-DDT
S-10876J1-1ML	100 ug/ml in Hexane

N-11571-10MG	Decafluorotriphenylphosphine
S-11571BO-1ML	50ug/ml in Acetone
S-11571BO-5ML	50ug/ml in Acetone
S-11571X0-1ML	50ug/ml in Methylene Chloride
S-11571X5-1ML	2000ug/ml in Methylene chloride

N-12831-500MG	Pentachlorophenol
S-12831M1-1ML	100 ug/ml in Methanol
S-12831M1-5ML	100 ug/ml in Methanol

METHOD 627 - DINITROANILINE PESTICIDES

Method 627 is applicable to the determination of dinitroaniline pesticides in industrial and municipal discharges as provided under 40 CFR 136.1 by gas chromatography. However, it fails to distinguish between Benfluralin, Ethalfluralin, and Trifluralin. When more than one of these materials may be present in a sample, the results are reported as Trifluralin.

CALIBRATION STANDARDS

EPA METHOD 627

N-13689-250MG	Trifluralin
S-13689A1-1ML	100ug/mL in Acetonitrile
S-13689U1-1ML	100ug/mL in Toluene
S-13689U1-5ML	100ug/mL in Toluene

METHOD 629 - BLADEx (CYANAZINE)

Method 629 is applicable to the determination of Cyanazine in municipal and industrial wastewater. It is a high-performance liquid chromatographic (HPLC) method. HPLC conditions are described which allow for separation and measurement of cyanazine in the extract by HPLC with a UV detector. The method provides for an optional Florisil column cleanup procedure to aid in the elimination or reduction of interferences.

BLADEx (CYANAZINE) STANDARDS

EPA METHOD 629

N-11327-250MG	Bladex
S-11327A1-1ML	100ug/ml in Acetonitrile
S-11327T4-5ML	1000ug/ml in tert-butylmethyl ether

METHOD 630/630.1 - DITHIOCARBAMATE PESTICIDES

Method 630 is applicable to the determination of dithiocarbamate pesticides. It is the colorimetric method for the determination of the dithiocarbamate pesticides compounds in industrial and municipal discharges as provided under 40 CFR 136.1. (Note that this method fails to distinguish between the individual dithiocarbamates, for all results are reported as Ziram). Carbon disulfide is an interferant.

Method 630.1 is applicable the determination of the dithiocarbamate pesticide compounds listed in municipal and industrial discharges as provided under 40 CFR 136.1, by total residue gas chromatography. This method includes procedures to purge CS₂ from the wastewater prior to acid hydrolysis of the sample. A vortex evaporator is used for CS₂ removal.

DITHIOCARBAMATE PESTICIDES CALIBRATION STANDARDS

EPA METHOD 630/630.1

N-13761-250MG	Ziram
----------------------	-------

INTERFERENCE CHECK STANDARDS

EPA METHOD 630.1

N-11406-1G	Carbon disulfide
S-11406J1-1ML	100ug/ml in Hexane
S-11406M1-1ML	100ug/ml in Methanol
S-11406M1-5ML	100ug/ml in Methanol

METHOD 631 - BENOMYL AND CARBENDAZIM

Method 631 is applicable to the determination of Benomyl and Carbendazim by high performance liquid chromatography (HPLC) in industrial and municipal discharges as provided under 40 CFR 136.1. Note that Benomyl is hydrolyzed to carbendazim, and both compounds are measured and reported as carbendazim.

BENOMYL STANDARDS

EPA METHOD 631

N-11138-100MG	Benomyl
----------------------	---------

CARBENDAZIM STANDARDS

EPA METHOD 631

N-11404-100MG	Carbendazim
S-11404M1-1ML	100ug/ml in Methanol
S-11404M1-5ML	100ug/ml in Methanol

METHOD 632 - CARBAMATE UREA PESTICIDES

Method 632 is applicable to the determination of carbamate and urea pesticides in industrial and municipal wastewater by high performance liquid chromatography (HPLC). Note that this method provides for an optional Florisil column cleanup procedure, to aid in the elimination or reduction of interferences, which may be encountered (caused by contaminants in solvents, reagents, glassware and other sample processing apparatus that lead to discrete artifacts or elevated baselines in liquid chromatograms).

CARBAMATE & UREA PESTICIDES MIXTURE**EPA METHOD 632**

Aminocarb Chlorpropham Carbofuran	Fenuron Linuron Methomyl	Neburon Propoxur Diuron	Barban Carbaryl Monuron	Fluometuron Methiocarb Swep	Oxamyl Siduron Propham
---	--------------------------------	-------------------------------	-------------------------------	-----------------------------------	------------------------------

1000ug/mL in Acetonitrile (18 components)
M-CUP6321A4-1ML 1mL Ampule

CARBAMATE & UREA PESTICIDES MIXTURE**EPA METHOD 632**

Aminocarb Monuron Carbofuran Propham	Diuron Siduron Linuron	Methomyl Carbaryl Oxamyl	Propoxur Neburon Zectran	Barban Swep Chlorpropham	Fenuron Methiocarb
---	------------------------------	--------------------------------	--------------------------------	--------------------------------	-----------------------

1000ug/mL in Acetonitrile (19 components)
M-CUP6322A4-1ML 1mL Ampule

ADDITIONAL ANALYTES**EPA METHOD 632**

N-11968-50MG Fenuron TCA
N-12498-100MG Monuron TCA

METHOD 634 - THIOCARBAMATE PESTICIDES

Method 634 is applicable to the determination of thiocarbamate pesticides in industrial and municipal wastewaters. Sample extraction and concentration steps in this method are similar to other 600 Series Methods. The separation and measurement of the compounds in the extract is by alkali flame detector (AFD) gas chromatography. To aid in the elimination of interferences an optional silica gel column cleanup procedure is provided. This method describes analytical conditions for a second gas chromatographic column that can be used to confirm measurements made with the primary column, as well as the gas chromatograph/mass spectrometer (GC/MS) criteria appropriate for the qualitative confirmation of these compound.

THIOCARBAMATE PESTICIDES MIXTURE**EPA METHOD 634**

Butylate
 Cycloate
 S-Ethyl dipropylthiocarbamate
 Molinate
 S-Propyl butylethylthiocarbamate
 Vernolate

1000ug/mL in tert-Butylmethyl ether (6 components)
M-TCP507T4-1ML 1mL Ampule

INTERNAL STANDARDS**EPA METHOD 634**

N-13708-1G Triphenyl phosphate
S-13708T3-1ML 500ug/ml in T-butylmethyl ether
S-13708T3-5ML 500ug/ml in T-butylmethyl ether

SURROGATE STANDARDS**EPA METHOD 634**

N-10114-1G 1,3-Dimethyl-2-nitrobenzene
S-10114T6-1ML 2500ug/ml in tert-butylmethyl ether
S-10114T6-5ML 2500ug/ml in tert-butylmethyl ether
S-10114T5-1ML 2000ug/ml in tert-butylmethyl ether
S-10114T5-5ML 2000ug/ml in tert-butylmethyl ether

METHOD 635 - ROTENONE

Method 635 is applicable to the determination of rotenone (as well as separation and measurement of the compounds in the extract by HPLC/UV) in industrial and municipal wastewaters.

ROTENONE STANDARDS**EPA METHOD 635**

N-13184-250MG Rotenone
S-13184A1-1ML 100ug/ml in Acetonitrile

METHOD 636 - BENSULIDE

Method 636 is applicable to the determination of Bensulide pesticide as well as separation and measurement of the compounds in the extract by HPLC-UV in municipal and industrial discharges as provided under 40 CFR 136.1.

BENSULIDE STANDARDS**EPA METHOD 636**

N-11140-250MG Bensulide
S-11140M4-1ML 1000ug/ml in Methanol
S-11140M4-5ML 1000ug/ml in Methanol

METHOD 638 - ORYZALIN

Method 638 is applicable to the determination of Oryzalin pesticide (as well as separation and measurement of the compounds in the extract by HPLC/UV in municipal and industrial discharges as provided under 40 CFR 136.1.

ORYZALIN STANDARDS**EPA METHOD 638**

N-12729-100MG Oryzalin
S-12729A4-1ML 1000ug/ml in Acetonitrile

METHOD 639 - BENDIOCARB

Method 639 is applicable to the determination of bendiocarb pesticide (as well as separation and measurement of the compounds in the extract by HPLC/UV) in municipal and industrial discharges as provided under 40 CFR 136.1.

BENDIOCARB STANDARDS**EPA METHOD 639**

N-11135-250MG Bendiocarb
S-11135U1-1ML 1000ug/ml in Toluene

INTERNAL STANDARDS**EPA METHOD 639**

N-12416-1G Methyl benzoate
S-12416M3-1ML 500ug/ml in Methanol

METHOD 640 - 2-MERCAPTOBENZOTHIAZOLE

Method 640 is applicable to the determination of Mercaptobenzothiazole (as well as separation and measurement of the compounds in the extract by HPLC/UV) in municipal and industrial discharges as provided under 40 CFR 136.1. This method provides a silica gel column cleanup procedure to aid in the elimination of interferences.

2-MERCAPTOBENZOTHIAZOLE STANDARDS**EPA METHOD 640**

N-10386-1G 2-Mercaptobenzothiazole
S-10386M4-1ML 1000ug/ml in Methanol
S-10386M4-5ML 1000ug/ml in Methanol

METHOD 641 - THIABENDAZOLE

Method 641 is applicable to the determination of thiabendazole (as well as the separation and accurate measurement of thiabendazole by direct aqueous injection HPLC with fluorescence detection) in municipal and industrial wastewater.

THIABENDAZOLE STANDARDS**EPA METHOD 641**

N-13560-250MG Thiabendazole
S-13560M4-1ML 1000ug/ml in Methanol

METHOD 643 - BENTAZON

Method 643 is applicable to the determination of Bensulide pesticide as well as separation and measurement of the compounds in the extract by HPLC-UV in municipal and industrial discharges as provided under 40 CFR 136.1.

BENTAZON STANDARDS**EPA METHOD 643**

N-11142-250MG Bentazon
S-11142A4-1ML 1000ug/ml in Acetonitrile
S-11142A4-5ML 1000ug/ml in Acetonitrile
S-11142B1-1ML 100ug/ml in Acetone

METHOD 644 - PICLORAM

Method 644 is applicable to the determination of Picloram in municipal and industrial wastewater by extraction, followed by analysis by HPLC with ultraviolet (UV) detection. An alkaline back-extraction is used as necessary to eliminate interferences.

PICLORAM STANDARDS**EPA METHOD 644**

N-13050-250MG Picloram
S-13050A4-1ML 1000ug/ml in Acetonitrile
S-13050A4-5ML 1000ug/ml in Acetonitrile
S-13050B1-1ML 100ug/ml in Acetone
S-13050B1-5ML 100ug/ml in Acetone

INTERNAL STANDARDS**EPA METHOD 644**

N-13200-1G Salicylic acid
S-13200L5-1ML 2000ug/ml in Isopropanol
S-13200L5-5ML 2000ug/ml in Isopropanol

METHOD 645 - AMINE PESTICIDES & LETHANE

Method 645 is applicable to the determination of amine pesticides and Lethane in municipal and industrial wastewater. Analysis is by gas chromatography. Measurement of the pesticides is accomplished with a nitrogen/phosphorous specific detector. This method provides identification by gas chromatography/mass spectrometer (GC/MS).

INTERNAL STANDARDS**EPA METHOD 645**

N-13708-1G	Triphenyl phosphate
S-13708T3-1ML	500ug/ml in T-butylmethyl ether
S-13708T3-5ML	500ug/ml in T-butylmethyl ether

INTERNAL STANDARDS**EPA METHOD 644**

N-13200-1G	Salicylic acid
S-13200L5-1ML	2000ug/ml in Isopropanol
S-13200L5-5ML	2000ug/ml in Isopropanol

METHOD 680 - PESTICIDE & PCB CONGENERS

Method 680 is applicable to the determination of pesticides and PCB congeners by gas chromatography and mass spectrometry.

ORGANOCHLORINE PESTICIDES MIXTURE #1**EPA METHOD 680, 508, 508.1, 525.2, 608, 617, 625, 8080A/8081, 8250A/8270B, CLP**

Aldrin (TM)	b-Endosulfan	BHC (alpha isomer)	BHC (beta isomer)
BHC (delta isomer)	cis-Chlordane	4,4'-DDD	4,4'-DDE
4,4'-DDT	Dieldrin	α-Endosulfan	Endosulfan sulfate
Endrin	Endrin aldehyde	Endrin ketone	Heptachlor
Heptachlor epoxide (Isomer B)	Lindane (BHC gamma isomer)	Methoxychlor	trans-Chlordane
1000ug/mL in tert-Butyl methyl ether	(20 components)		
M-OCP5081T4-1ML	1mL Ampule		

METHOD 1618 - ORGANOHALIDE & ORGANOPHOSPHOROUS PESTICIDES

Method 1618 is applicable to the determination of organohalide pesticides, organophosphorus pesticides, and phenoxy-acid herbicides. This method applies a wide bore capillary column gas chromatography procedure for the analysis with halogen-specific and organophosphorus detectors.

ORGANOPHOSPHOROUS PESTICIDES MIXTURE #1**EPA METHOD 1618**

50ug/mL	Coumaphos	100ug/mL	Diazinon
50ug/mL	Dichlorvos	100ug/mL	Dimethoate
100ug/mL	Dylox®	100ug/mL	EPN
100ug/mL	Fenchlorphos	200ug/mL	Fensulfothion
100ug/mL	Azinphos-methyl	100ug/mL	Guthion Ethyl
100ug/mL	Malathion	100ug/mL	Methyl parathion
100ug/mL	Parathion®	100ug/mL	Phorate
100ug/mL	Prophos	50ug/mL	Sulprofos
100ug/mL	Terbufos	100ug/mL	Tetrachlorvinphos
200ug/mL	S.S.S-Tributylphosphorotriothioate		

Varied concentration in Isooctane (19 components)
M-OPP16181K99-1ML **1mL Ampule**

ORGANOPHOSPHOROUS PESTICIDES MIXTURE #2**EPA METHOD 1618**

50ug/mL	Chlorfenvinphos	50ug/mL	Chlorpyrifos
100ug/mL	Chlorpyrifos Methyl	200ug/mL	Crotoxyphos
200ug/mL	Demeton S	100ug/mL	Naled
100ug/mL	Dichlofenthion	600ug/mL	Dioxathion
100ug/mL	Disulfoton	100ug/mL	Ethion
200ug/mL	Famphur	100ug/mL	Fenthion
100ug/mL	Mevinphos	200ug/mL	Phosmet
100ug/mL	Leptophos	50ug/mL	Tetraethylthiopyrophosphate
100ug/mL	Tokuthion®	100ug/mL	Trichloronate

Varied concentration in Isooctane - (18 components)
M-OPP16182K99-1ML **1mL Ampule**

CHLORINATED HERBICIDES MIXTURE #1**EPA METHOD 1618, 8150B, 8321**

2,4-D	Dalapon	Dichlorprop	4-Chloro-o-tolyloxyacetic acid
Silvex	2,4-DB	Dicamba	Dinoseb
Mecoprop	2,4,5-T®		

100ug/mL in Acetone - (10 components)
M-CH8150B1-1ML **1mL Ampule**

METHYLATED CHLORINATED HERBICIDES MIXTURE #1**EPA METHOD 1618, 8150B**

2,4-D methyl ester	Dinoseb methyl ether	2,4-DB methyl ester	4-Chloro-o-toloxoy acetic acid methyl ester
Dalapon methyl ester	Mecoprop methyl ester	Dicamba methyl ester	Silvex, methyl ester
Dichlorprop methyl ester	2,4,5-(Trichlorophenoxy) acetic acid methyl ester		

100ug/mL in Methanol -
M-MCH8150M1-1ML

(10 components)
1mL Ampule

DIBUTYL CHLORENDATE STANDARDS**EPA METHOD 1618**

N-11638-250MG	Dibutyl chlorendate
S-11638B5-1ML	2000ug/mL in Acetone
S-11638B5-5ML	2000ug/mL in Acetone
S-11638X4-1ML	1000ug/mL in Methylene Chloride
S-11638X4-5ML	1000ug/mL in Methylene Chloride

SPE CARTRIDGE CALIBRATION STANDARDS**EPA METHOD 1618**

N-10657-100MG	2,4-Dichlorophenylacetic acid
S-10657M1-1ML	100ug/mL in Methanol
S-10657M1-5ML	100ug/mL in Methanol

SURROGATE STANDARDS FOR PHENOXY ACID HERBICIDES**EPA METHOD 1618**

N-10536-100MG	2,4-Dichlorophenylacetic acid
S-10536B5-1ML	2000ug/mL in Acetone
S-10536B5-5ML	2000ug/mL in Acetone
S-10536T1-1ML	100ug/mL in tert-Butylmethyl ether
S-10536T1-5ML	100ug/mL in tert-Butylmethyl ether
S-10536T7-1ML	5000ug/mL in tert-Butylmethyl ether
S-10536T7-5ML	5000ug/mL in tert-Butylmethyl ether

NON-ITD ORGANOPHOSPHATE COMPOUND**EPA METHOD 1618**

MET-11408AM1-1ML	Carbophenothion methyl
	1000ug/mL in Methanol

DEGRADATION CALIBRATION MIXTURE**EPA METHOD 1618, 608, 508/508.1, 525.2,617, 625, 1656, 8080A/8081, 8250A/8270B, CLP**

200ug/mL	4,4'-DDT
100ug/mL	Endrin

Varied concentration in tert-Butylmethyl ether (2 components)
M-DC508T99-1ML **1mL Ampule**

GPC CALIBRATION MIXTURE**EPA METHOD 1618, 1656, 1657**

15,000ug/mL	Bis(2-ethylhexyl)phthalate
300,000ug/mL	Corn oil
1,400ug/mL	Pentachlorophenol
100ug/mL	Perylene
500ug/mL	Sulfur-sublimed

Varied concentration in Isooctane (5 components)
M-GPC16571K99-1ML **1mL Ampule**

DEGRADATION PRODUCTS MIXTURE**EPA METHOD 1618, 608, 508/508.1, 525.2,617, 625, 1656, 8080A/8081, 8250A/8270B, CLP**

200ug/mL	4,4'-DDD
200ug/mL	4,4'-DDE
100ug/mL	Endrin aldehyde
100ug/mL	Endrin ketone

Varied concentration in tert-Butylmethyl ether (4 components)
M-DP5081T99-1ML **1mL Ampule**

METHOD 1656 - ORGANOHALIDE PESTICIDES

Method 1656 is applicable to determination of the organo-halide pesticides and polychlorinated biphenyls (PCBs) and other compounds amenable to extraction and analysis by wide-bore capillary column gas chromatography (GC) with halogen-specific detectors (HSD).

DEGRADATION CALIBRATION MIXTURE**EPA METHOD 1656, 608, 508/508.1, 525.2,617, 625, 1618, 8080A/8081, 8250A/8270B, CLP**

200ug/mL	4,4'-DDT
100ug/mL	Endrin

Varied concentration in tert-Butylmethyl ether (2 components)
M-DC508T99-1ML **1mL Ampule**

GPC CALIBRATION MIXTURE**EPA METHOD 1656, 1618, 1657**

15,000ug/mL	Bis(2-ethylhexyl)phthalate
300,000ug/mL	Corn oil
1,400ug/mL	Pentachlorophenol
100ug/mL	Perylene
500ug/mL	Sulfur-sublimed

Varied concentration in Isooctane (5 components)
M-GPC16571K99-1ML **1mL Ampule**

DEGRADATION PRODUCTS MIXTURE**EPA METHOD 1656, 608, 508/508.1, 525.2,617, 625, 1656, 8080A/8081, 8250A/8270B, CLP**

200ug/mL	4,4'-DDD
200ug/mL	4,4'-DDE
100ug/mL	Endrin aldehyde
100ug/mL	Endrin ketone

Varied concentration in tert-Butylmethyl ether (4 components)
M-DP5081T99-1ML **1mL Ampule**

DIBUTYL CHLORENDATE STANDARDS**EPA METHOD 1656**

N-11638-250MG	Dibutyl chlorendate
S-11638B5-1ML	2000ug/mL in Acetone
S-11638B5-5ML	2000ug/mL in Acetone
S-11638X4-1ML	1000ug/mL in Methylene Chloride
S-11638X4-5ML	1000ug/mL in Methylene Chloride

METHOD 1657 - ORGANOPHOSPHOROUS PESTICIDES

Method 1657 is applicable to the determination in waters, soils, sediments, and sludges, organophosphorus pesticides and additional compounds amenable to extraction and analysis by wide-bore capillary column gas chromatography (GC) with a flame photometric detector (FPD).

SURROGATE SPIKING MIXTURE**EPA METHOD 1657**

Tributyl phosphate
Triphenyl phosphate

2000ug/mL in Acetone (2 components)
M-SS1657B5-1ML 1mL Ampule

SURROGATE STANDARDS**EPA METHOD 1657**

N-13194-250MG S,S,S-Tributyl phosphorotrithioate
S-13194A1-1ML 100ug/mL in Acetonitrile
S-13194T1-1ML 100ug/mL in tert-Butylmethyl ether

SURROGATE STANDARDS**EPA METHOD 1657**

N-13645-1G Tributyl phosphate
S-13645B5-1ML 2000ug/mL in Acetone
S-13645B5-5ML 2000ug/mL in Acetone

N-13708-1G Triphenyl phosphate
S-13708T3-1ML 500ug/mL in tert-Butylmethyl ether
S-13708T3-5ML 500ug/mL in tert-Butylmethyl ether

ADDITIONAL ANALYTE**EPA METHOD 1657**

MET-11408AM1-1ML Carbophenothion methyl
1000ug/mL in Methanol

METHOD 1660 - PYRETHRINS & PYRETHROIDS

Method 1660 is applicable to the determination of pyrethrins and pyrethroids in municipal and industrial wastewater by extraction and high-performance liquid chromatography (HPLC) with an ultra-violet detector (UV). Separation and measurement of the pyrethrins and pyrethroids is performed by reverse-phase C18 column HPLC with a multiple-wavelength UV detector.

PYRETHRINS & PYRETHROIDS MIXTURE #1**EPA METHOD 1660**

400ug/mL Cyfluthrin
400ug/mL d-(cis-trans)-Phenothrin
200ug/mL Fenvalerate
400ug/mL Tetramethrin

Varied concentration in Acetonitrile (4 components)
M-PP16601A99-1ML 1mL Ampule

PYRETHRINS & PYRETHROIDS MIXTURE #2**EPA METHOD 1660**

Allethrin Resmethrin Permethrin

400ug/L in Acetonitrile (3 components)
M-PP16602A15-1ML 1mL Ampule

PYRETHRUM STANDARDS**EPA METHOD 1660**

N-13151-100MG Pyrethrum
S-13151A1-1ML 100ug/mL in Acetonitrile
S-13151T1-1ML 100ug/mL in tert-Butylmethyl ether

METHOD 8080A/8081/8081A - ORGANOCHLORINE PESTICIDES

Method 8080 is applicable for the determination of organochlorine pesticides and PCBs utilizing a GC with detection by an electron capture detector (ECD) or a halogen-specific detector (HSD). The sensitivity of Method 8080 usually depends on the level of interferences rather than on instrumental limitations.

Method 8080A is applicable for the determination of organochlorine pesticides and polychlorinated biphenyls by gas chromatography.

Method 8081 is applicable for determination of organochlorine pesticides and PCBs as arochlors by gas chromatography utilizing capillary column technique.

ORGANOCHLORINE PESTICIDES MIXTURES**EPA METHOD 8080A/8081, 608, 508/508.1, 617, 625, 8250A/8270B, CLP**

Aldrin (TM)	b-Endosulfan	BHC (alpha isomer)	BHC (beta isomer)
BHC (delta isomer)	4,4'-DDD	4,4'-DDE	4,4'-DDT
Dieldrin	α-Endosulfan	Endosulfan sulfate	Endrin
Endrin aldehyde	Heptachlor	Heptachlor epoxide (Isomer B)	
Lindane (BHC gamma isomer)	Methoxychlor		

100 ug/mL in Toluene:Hexane (50:50) (17 components)
M-OCp8080AC1-1ML 1mL Ampule
100 ug/mL in Toluene:Hexane (50:50), less Methoxychlor (16 components)
M-PPO8AC1-1ML 1mL Ampule
2000 ug/mL in Toluene:Hexane (50:50), less Methoxychlor (16 components)
M-PPHC5AC5-1ML 1mL Ampule

PESTICIDES - CONTROL SAMPLE MIXTURE**EPA METHOD 8080A/8081, 608**

20ug/mL Aldrin	20ug/mL BHC (alpha isomer)
20ug/mL BHC (beta isomer)	100ug/mL 4,4'-DDD
20ug/mL 4,4'-DDE	100ug/mL 4,4'-DDT
20ug/mL Dieldrin	20ug/mL α-Endosulfan
100ug/mL b-Endosulfan	100ug/mL Endosulfan sulfate
100ug/mL Endrin	20ug/mL Heptachlor
20ug/mL Heptachlor epoxide (Isomer B)	

Varied concentration in Toluene (13 components)
M-CSM8080U99-1ML 1mL Ampule

DEGRADATION CALIBRATION MIXTURE**EPA METHOD 8080A/8081,608,508/508.1, 525.2, 617, 625, 1618, 1656, 8250A/8270B, CLP**200ug/mL 4,4'-DDT
100ug/mL EndrinVaried concentration in tert-Butylmethyl ether (2 components)
M-DC508T99-1ML 1mL Ampule**SURROGATE STANDARDS MIXTURE #1****EPA METHOD 8080A/8081/8081A, CLP**Decachlorobiphenyl
2,4,5,6-Tetrachloro-m-xylene2000ug/mL in Toluene (2 components)
M-SM8080U5-1ML 1mL Ampule**DEGRADATION PRODUCTS MIXTURE****EPA METHOD 8080A/8081 608, 508/508.1, 525.2, 617, 625, 1618, 1656,8250A/8270B, CLP**200ug/mL 4,4'-DDD
200ug/mL 4,4'-DDE
100ug/mL Endrin aldehyde
100ug/mL Endrin ketoneVaried concentration in tert-Butylmethyl ether (4 components)
M-DP5081T99-1ML 1mL Ampule**SURROGATE STANDARDS MIXTURE #2****EPA METHOD 8080A/8081/8081A, CLP**Dibutyl chlorendate
2,4,5,6-Tetrachloro-m-xylene2000ug/mL in Acetone (2 components)
M-SM8080AB5-1ML 1mL Ampule**METHOD 8140/8141A/B - ORGANOPHOSPHOROUS PESTICIDES**

Method 8141 is applicable for determination of organophosphorous compounds by gas chromatography utilizing a capillary column technique.

ORGANOPHOSPHOROUS PESTICIDES**EPA METHOD 8140/8141A/8141B**

Chlorpyrifos	Dichlorvos	Guthion®	Sulprofos
Coumaphos	Disulfoton	Methyl parathion	Tetrachlorvinphos
Demeton-S	Fenchlorophos	Phorate	Tokuthion®
Diazinon	Fensulfothion	Phosdrin	Tributylphosphorotrithioite
Dibrom	Fenthion	Prophos	Trichloronate

100ug/ml in Methylene chloride (20) components
M-OPP8140X1-1ML 1mL Ampule**ORGANOPHOSPHOROUS COMPOUNDS MIXTURE****EPA METHOD 8141A/8141B**

Aspon	Dichlofenthion	Ethion	Imidan
Carbophenothion	Dicrotophos	Famphur	Phosphamidon
Chlorfenvinphos	Dioxathion	Fenitrothion	Terbufos
Chlorpyrifos methyl	Dyfonate®	Guthion Ethyl	Zinophos
Crotoxyphos	Dylox®	Leptophos	

100ug/ml in Toluene (19) components
M-OPC8141A1U1-1ML 1mL Ampule**SPECIAL ORGANOPHOSPHOROUS MIXTURE****EPA METHOD 8141A/8141B**Diazinon Methyl parathion
Ethion Parathion
Malathion100ug/ml in Cyclohexane (5 components)
M-CR81411E1-1ML 1mL Ampule**TRIAZINE HERBICIDES MIXTURES****EPA METHOD 8141A/8141B**Atrazine
Simazine200ug/ml in Acetone (2 components)
M-TH81411B2-1ML 1mL Ampule**ORGANOPHOSPHOROUS COMPOUNDS****EPA METHOD 8141A/8141B (SUPPLEMENT TO EPA METHOD 8140)**

Dimethoate	Malathion	Parathion
EPN	Monocrotophos	Tetraethyldithiopyrophosphate

100ug/ml in Toluene - (6 components)
M-OPP8141U1-1ML 1mL Ampule**TEPP STANDARDS****EPA METHOD 8141A/8141B****N-13543-500MG** Tetraethylpyrophosphate**EPN STANDARDS****EPA METHOD 8141A/8141B**

N-11858-250MG EPN
S-11858A1-1ML 100ug/ml in Acetonitrile
S-11858J4-1ML 1000ug/ml in Hexane
S-11858J4-5ML 1000ug/ml in Hexane
S-11858J1-1ML 100ug/ml in Hexane
S-11858J1-5ML 100ug/ml in Hexane

ATRAZINE STANDARDS**EPA METHOD 8141A/8141B**

N-11106-250MG Atrazine
S-11106T1-1ML 100ug/ml in T-butylmethyl ether
S-11106T1-5ML 100ug/ml in T-butylmethyl ether

DIBROM DEGRADATION STANDARDS**EPA METHOD 8141A/8141B**

N-11675-250MG Dichlorvos
S-11675J4-1ML 1000ug/ml in Hexane
S-11675J4-5ML 1000ug/ml in Hexane
S-11675U1-1ML 100ug/ml in Toluene
S-11675U1-5ML 100ug/ml in Toluene

TEPP DEGRADATION STANDARDS**EPA METHOD 8141A/8141B**

N-13670-1G	Triethyl phosphate
S-13670J1-1ML	100ug/ml in Hexane
S-13670J1-5ML	100ug/ml in Hexane

SURROGATE STANDARDS**EPA METHOD 8140/8141A/8141B**

N-10816-1G	4-Chloro-3-nitrobenzotrifluoride
S-10816K4-1ML	1000ug/ml in Isooctane
S-10816K4-5ML	1000ug/ml in Isooctane
N-13645-1G	Tributyl phosphate
S-13645B5-1ML	2000ug/ml in Acetone
S-13645B5-5ML	2000ug/ml in Acetone
N-13708-1G	Triphenyl phosphate
S-13708T3-1ML	500ug/ml in tert-Butylmethyl ether
S-13708T3-5ML	500ug/ml in tert-Butylmethyl ether

SIMAZINE STANDARDS**EPA METHOD 8141A/8141B**

N-13800-500MG	Simazine
S-13800B4-1ML	1000ug/ml in Acetone
S-13800B4-5ML	1000ug/ml in Acetone
S-13800M1-1ML	100ug/ml in Methanol

MERPHOS CALIBRATION STANDARDS**EPA METHOD 8140/8141A, 622, 507, 1657, 8321**

N-13194-250MG	S,S,S-Tributyl phosphorotrithioate
S-13194A1-1ML	100ug/mL in Acetonitrile
S-13194T1-1ML	100ug/mL in tert-Butylmethyl ether

METHOD 8150B - CHLORINATED HERBICIDES

Method 8150 is applicable for the determination of chlorinated herbicides. Method 8150 provides extraction, esterification, and GC conditions for the analysis. When Method 8150 is used to analyze unfamiliar samples, compound identifications should be supported by at least one additional qualitative technique. This method describes analytical conditions for a second gas chromatographic column that can be used to confirm measurements made with the primary column.

Method 8150A is applicable for the determination of chlorinated herbicides by gas chromatography. It provides GC/MS criteria for qualitative confirmation of compound identifications.

Method 8150B is applicable for the determination of chlorinated herbicides by gas chromatography. The sensitivity of the 8150 methods usually depends on the level of interferences rather than on the instrumental limitations

CHLORINATED HERBICIDES MIXTURE #1**EPA METHOD 8150B, 1618, 8321**

2,4-D	Dalapon	Dichlorprop	4-Chloro-o-tolylloxyacetic acid
Silvex	2,4-DB	Dicamba	Dinoseb
Mecoprop	2,4,5-T®		

100ug/mL in Acetone - (10 components)
M-CH8150B1-1ML 1mL Ampule

METHYLATED CHLORINATED HERBICIDES MIXTURE #1**EPA METHOD 8150B, 1618**

2,4-D methyl ester	Dinoseb methyl ether	2,4-DB methyl ester	4-Chloro-o-tolylloxy acetic acid methyl ester
Dalapon methyl ester	Mecoprop methyl ester	Dicamba methyl ester	Silvex, methyl ester
Dichlorprop methyl ester	2,4,5-(Trichlorophenoxy) acetic acid methyl ester		

100ug/mL in Methanol (10 components)
M-MCH8150M1-1ML 1mL Ampule

CHLORINATED HERBICIDES MIXTURE - A**EPA METHOD 8150B, 515.1**

2,4-D
 Silvex
 2,4,5-T

2000ug/mL in Acetone (3 components)
M-CSHC14AB5-1ML 1mL Ampule

CHLORINATED HERBICIDES ESTERS MIXTURE**EPA METHOD 8150B, 515.1**

2,4-D methyl ester
 Silvex methyl ester
 2,4,5-T methyl ester

2000ug/mL in Hexane (3 components)
M-CSHC14JB5-1ML 1mL Ampule

METHYLATED CHLORINATED HERBICIDES - CONTROL SAMPLE MIXTURE**EPA METHOD 8150B, 615**

10000ug/mL	4-Chloro-o-tolylloxyacetic acid methyl ester
250ug/mL	Dalapon methyl ester
10ug/mL	Dicamba methyl ester
50ug/mL	Dinoseb methyl ether
10ug/mL	Silvex methyl ester

Varied concentration in Hexane - (10 components)
M-MCSM81501J99-1ML 1mL Ampule

100ug/mL	2,4-D methyl ester
100ug/mL	2,4-DB methyl ester
100ug/mL	Dichlorprop methyl ester
10000ug/mL	Mecoprop methyl ester
10ug/mL	(2,4,5-Trichlorophenoxy)acetic acid methyl ester

CHLORINATED HERBICIDES - CONTROL SAMPLE MIXTURE**EPA METHOD 8150B, 615**

10,000ug/mL	4-Chloro-o-toloxoyacetic acid	100ug/mL	2,4-D
250ug/mL	Dalapon	100ug/mL	2,4-DB
10ug/mL	Dicamba	100ug/mL	Dichlorprop
50ug/mL	Dinoseb	10,000ug/mL	Mecoprop
10ug/mL	Silvex	10ug/mL	2,4,5-T®

Varied concentration in Acetone

(10 components)

M-CSM81501B99-1ML**1mL Ampule****SURROGATE STANDARDS****EPA METHOD 8150B**

N-10536-100MG	2,4'-Dichlorophenylacetic acid	N-10537-100MG	2,4'-Dichlorophenylacetic acid methyl ester
S-10536B5-1ML	2000ug/mL in Acetone	S-10537B1-1ML	100ug/mL in Acetone
S-10536B5-5ML	2000ug/mL in Acetone	S-10537B1-5ML	100ug/mL in Acetone
S-10536T1-1ML	100ug/mL in tert-Butylmethyl ether	S-10537B5-1ML	2000ug/mL in Acetone
S-10536T1-5ML	100ug/mL in tert-Butylmethyl ether	S-10537M7-1ML	5000ug/mL in Methanol
S-10536T7-1ML	5000ug/mL in tert-Butylmethyl ether	S-10537M7-5ML	5000ug/mL in Methanol
S-10536T7-5ML	5000ug/mL in tert-Butylmethyl ether		

METHOD 8151/8151A - CHLORINATED HERBICIDES

Method 8151 is applicable for the determination of chlorinated herbicides by GC, using methylation or pentafluorobenzoylation derivitization, by capillary column technique.

Method 8151A is applicable for determination of chlorinated herbicides by GC using methylation or pentafluorobenzoylation derivitization. This method utilizes a capillary column technique with the GC

CHLORINATED HERBICIDES MIXTURE**EPA METHOD**

Acifluorfen	2,4-DB	4-Nitrophenol	Bentazon
Dicamba	Pentachlorophenol	Chloramben	3,5-Dichlorobenzoic acid
Picloram	4-Chloro-o-toloxoyacetic acid	Dichlorprop	Silvex
2,4-D	Dinoseb	2,4,5-T®	Dalapon
Mecoprop	Tetrachloroterephthalic acid		

1000ug/mL in Acetone -

(18 components)

M-CH8151B4-1ML**1mL Ampule****METHYLATED CHLORINATED HERBICIDES MIXTURE****EPA METHOD 8151, 615**

Acifluorfen methyl ester	3,5-Dichlorobenzoic acid methyl ester	Bentazon methyl derivative
Dichlorprop methyl ester	Chloramben methyl ester	Dinoseb methyl ether
4-Chloro-o-toloxoyacetic acid methyl ester	Mecoprop methyl ester	Chlorthal
p-Nitroanisole	2,4-D methyl ester	Pentachloroanisole
Dalapon methyl ester	Picloram methyl ester	2,4-DB methyl ester
Silvex methyl ester	Dicamba methyl ester	(2,4,5-Trichlorophenoxy)acetic acid methyl ester

1000ug/mL in Isooctane:Acetone (90:10)

(18 components)

M-MCH8151Y4-1ML**1mL Ampule****SURROGATE STANDARDS****EPA METHOD 8151**

N-10536-100MG	2,4'-Dichlorophenylacetic acid
S-10536B5-1ML	2000ug/mL in Acetone
S-10536B5-5ML	2000ug/mL in Acetone
S-10536T1-1ML	100ug/mL in tert-Butylmethyl ether
S-10536T1-5ML	100ug/mL in tert-Butylmethyl ether
S-10536T7-1ML	5000ug/mL in tert-Butylmethyl ether
S-10536T7-5ML	5000ug/mL in tert-Butylmethyl ether

N-10216-1G

1,4-Dichlorobenzene

S-10216M1-1ML

100ug/mL in Methanol

S-10216M1-5ML

100ug/mL in Methanol

SPIKING STANDARDS**EPA METHOD 8151**

N-11345-250MG	Bromoxynil
S-11345M1-1ML	100ug/mL in Methanol

INTERNAL STANDARDS**EPA METHOD 8151**

N-10867-100MG	4,4'-Dibromooctafluorobiphenyl
S-10867M7-1ML	5000ug/mL in Methanol
S-10867M7-5ML	5000ug/mL in Methanol
S-10867T1-1ML	100ug/mL in tert-Butylmethyl ether
S-10867T1-5ML	100ug/mL in tert-Butylmethyl ether
S-10867T5-1ML	2000ug/mL in tert-Butylmethyl ether
S-10867X5-1ML	2000ug/mL in Methylene chloride
S-10867X5-5ML	2000ug/mL in Methylene chloride

METHOD 8250A/8270B - SEMIVOLATILE ORGANICS

Method 8250A is applicable for determination of semivolatile organic compounds in extracts prepared from all types of solid waste matrices, soils, and ground water. Direct injection of a sample may be used in limited applications. Method 8250A can be used to quantify most neutral, acidic, and basic organic compounds that are soluble in methylene chloride, and capable of being eluted without derivatization as sharp peaks, from a gas chromatographic packed column by GC/MS.

Method 8270B is applicable for determination of semivolatile organic compounds by gas chromatography/mass spectrometry (GC/MS) by capillary column technique. This method is used to determine the concentration of semivolatile organic compounds in extracts prepared from all types of solid waste matrices, soils, and ground water. Direct injection of a sample may be used in limited applications. Prior to using this method, the samples should be prepared for chromatography using the appropriate sample preparation and cleanup methods. This method describes chromatographic conditions that will allow for the separation of the compounds in the extract, and for their qualitative and quantitative analysis, by mass spectrometry.

SEMI-VOLATILE MIXTURE #1**EPA METHOD 8270B**

Chlorobenzilate
Famphur
Tetraethyldithio pyrophosphate

Disulfoton
Phorate
Dinoseb

Parathion®
Dimethoate
Methyl parathion

Di-allate
Kepone®
Zinophos®

1000ug/mL in Methylene chloride
M-SV82701X4-1ML

(12 components)
1mL Ampule

ORGANOCHLORINE PESTICIDES MIXTURE #1**EPA METHOD 8250A/8270B, 508, 508.1, 525.2, 608, 617, 625, 8080A/8081, CLP**

Aldrin
BHC (delta isomer)
4,4'-DDT
Endrin
Heptachlor epoxide (Isomer B)

b-Endosulfan
cis-Chlordane
Dieldrin
Endrin aldehyde
Lindane (BHC gamma isomer)

BHC (alpha isomer)
4,4'-DDD
α-Endosulfan
Endrin ketone
Methoxychlor

BHC (beta isomer)
4,4'-DDE
Endosulfan sulfate
Heptachlor
trans-Chlordane

1000ug/mL in tert-Butyl methyl ether
M-OC5081T4-1ML

(20 components)
1mL Ampule

ORGANOCHLORINE PESTICIDES MIXTURES**EPA METHOD 8250A/8270B, 508/508.1, 608, 617, 625, 8080A/8081, CLP**

Aldrin (TM)
BHC (delta isomer)
Dieldrin
Endrin aldehyde
Lindane (BHC gamma isomer)

b-Endosulfan
4,4'-DDD
α-Endosulfan
Heptachlor
Methoxychlor

BHC (alpha isomer)
4,4'-DDE
Endosulfan sulfate
Heptachlor epoxide (Isomer B)

BHC (beta isomer)
4,4'-DDT
Endrin

100 ug/ml in Toluene:Hexane (50:50)

M-OC8080AC1-1ML

(17 components)
1mL Ampule

100 ug/ml in Toluene:Hexane (50:50), less Methoxychlor

M-PP08AC1-1ML

(16 components)
1mL Ampule

2000 ug/ml in Toluene:Hexane (50:50), less Methoxychlor

M-PPHC5AC5-1ML

(16 components)
1mL Ampule

ORGANOPHOSPHOROUS PESTICIDES HIGH CONCENTRATION MIXTURE**EPA METHOD 8270B**

Dimethoate
Disulfoton
Famphur
Methyl parathion
Phorate

Parathion
Tetraethyldithiopyrophosphate
O,O,O-Triethylphosphorothioate
Thionazin

2000ug/ml in Hexane:Acetone (80:20)

M-CSHC13AG5-1ML

(9 components)
1mL Ampule

DEGRADATION PRODUCTS MIXTURE**EPA METHOD 8250A/8270B, 508/508.1, 525.2, 608, 617, 625, 1618, 1656, 8080A/8081, CLP**

200ug/mL
200ug/mL
100ug/mL
100ug/mL

4,4'-DDD
4,4'-DDE
Endrin aldehyde
Endrin ketone

Varied concentration in tert-Butylmethyl ether
M-DP5081T99-1ML

(4 components)
1mL Ampule

ORGANOCHLORINE PESTICIDE MATRIX SPIKING MIXTURE**EPA METHOD 8250A/8270B**

200ug/mL Aldrin
500ug/mL 4,4'-DDT
500ug/mL Dieldrin
500ug/mL Endrin
200ug/mL Heptachlor
200ug/mL Lindane

Varied Concentration in Acetone

M-OCPSB99-1ML

(6 components)
1mL Ampule

DEGRADATION CALIBRATION MIXTURE**EPA METHOD 8250A/8270B, 508/508.1, 525.2, 608, 617, 625, 1618, 1656, 8080A/8081, CLP**

200ug/mL
100ug/mL

4,4'-DDT
Endrin

Varied concentration in tert-Butylmethyl ether
M-DC508T99-1ML

(2 components)
1mL Ampule

POLYNUCLEAR AROMATIC HYDROCARBONS (PAH) EPA METHOD 8270B, 550.1, 610, 8100, 8310, CLP

Acenaphthene Benzo(b)fluoranthene Chrysene Indeno(1,2,3-C,D)pyrene	Acenaphthylene Benzo(k)fluoranthene 1,2:5,6-Dibenzanthracene Naphthalene	Anthracene 1,12-Benzoperylene Fluoranthene Phenanthrene	1,2-Benzanthracene Benzo(a)pyrene Fluorene Pyrene
100ug/mL in Methanol M-PPH10M1-1ML M-PPH10M1-5ML	(16) components 1mL Ampule 5mL Ampule		
200ug/mL in Acetonitrile M-PNA550A2-1ML	(16) components 1mL Ampule		
2000 ug/mL in CH ₂ Cl ₂ :Benzene (50:50) M-PPHC6AD5-1ML M-PPHC6AD5-5ML	(16) components 1mL Ampule 5mL Ampule		

TUNING STANDARDS MIXTURE

EPA METHOD 8250A/8270B, 625, CLP

Benzidine
4,4'-DDT
Decafluorotriphenylphosphine
Pentachlorophenol

1000ug/mL in Methylene chloride M-CLPTS1X4-1ML M-CLPTS1X4-5ML	(4 components) 1mL Ampule 5mL Ampule
---	--

METHOD 8318 - n-METHYLCARBAMATES

Method 8318 is applicable for determination of n-Methylcarbamates by high performance liquid chromatography (HPLC). N-Methylcarbamates are extracted from aqueous samples with Methylene chloride, from soils, oily solid waste, and oils with acetonitrile. After separation, the target analytes are hydrolyzed and derivatized post-column, then quantified fluorometrically.

CARBAMATE PESTICIDES EPA METHOD 8318, 531.1

Aldicarb Aldicarb sulfone Oxamyl	Carbaryl Carbofuran	Methiocarb Methomyl	Propoxur Aldicarb sulfoxide	3-Hydroxycarbofuran
100ug/mL in Acetonitrile- M-CP8318A1-1ML	(10 components) 1mL Ampule			
100ug/mL in Acetonitrile- M-CP8318A1-5ML	(10 components) 5mL Ampule			

CARBAMATE PESTICIDES EPA METHOD, 531.1

Aldicarb Aldicarb sulfone Dioxacarb	Carbaryl Carbofuran Promecarb	Methiocarb Methomyl	Propoxur 3-Hydroxycarbofuran
1000ug/mL in Acetonitrile- M-CP8318A4-1ML	(10 components) 1mL Ampule		

LABORATORY PERFORMANCE CHECK MIXTURE

EPA METHOD 8318

Acetonitrile
Aldicarb sulfoxide
BDMC
3-Hydroxycarbofuran
Methiocarb

Varied concentration in Acetonitrile M-LPC531A99-1M	(5 components) 1mL Ampule
---	-------------------------------------

BDMC STANDARDS

EPA METHOD 8318

N-11131-100MG	BDMC
S-11131M1-1ML	100ug/mL in Methanol
S-11131M1-5ML	100ug/mL in Methanol

METHOD 8321/8321A/8321B - ORGANOPHOSPHOROUS

Method 8321 is applicable to the determination, by high performance liquid chromatography (HPLC), coupled with either thermospray-mass spectrometry (TSP-MS), and/or ultraviolet (UV), of disperse azo dyes, organophosphorus compounds, and Tris-(2,3-dibromopropyl) phosphate, in wastewater, ground water, sludge, soil/sediment matrices, and chlorinated phenoxyacid compounds, and their esters. Quantitative analysis is performed by TSP/MS, using an external standard approach. Sample extracts can be analyzed by direct injection into the thermospray or onto a liquid chromatographic-thermospray interface. Note that detection is achieved both by negative ionization (discharge electrode) and positive ionization, with a single quadrupole mass spectrometer. Since this method is based on an HPLC technique, the use of ultraviolet (UV) detection is optional for routine samples

CHLORINATED HERBICIDES MIXTURE #1**EPA METHOD 8321, 1618, 8150B,**

2,4-D	Dalapon	Dichlorprop	4-Chloro-o-tolxyloxyacetic acid
Silvex	2,4-DB	Dicamba	Dinoseb
Mecoprop	2,4,5-T®		

100ug/mL in Acetone - (10 components)
M-CH8150B1-1ML 1mL Ampule

MERPPOS CALIBRATION STANDARDS**EPA METHOD 8321, 622, 507, 1657, 8140/8141A,**

N-13194-250MG	S,S,S-Tributyl phosphorotriithioate
S-13194A1-1ML	100ug/mL in Acetonitrile
S-13194T1-1ML	100ug/mL in tert-Butylmethylether

ALKALOIDS**EPA METHOD 8321A/B**

N-11393-1G	Caffeine
S-11393A7-1ML	5000ug/mL in Acetonitrile
S-11393A7-5ML	5000ug/mL in Acetonitrile
N-13231-100MG	Strychnine
S-13231M1-1ML	100ug/mL in methanol
S-13231M1-5ML	100ug/mL in methanol

INDIVIDUAL ESTER STANDARDS**EPA METHOD 8321**

N-10610-250MG	2,4-D butoxyethyl ester
S-10610A1-1ML	100ug/mL in Acetonitrile
S-10610A4-1ML	1000ug/mL in Acetonitrile
S-10610A4-5ML	1000ug/mL in Acetonitrile

N-10621-250MG	2,4-sec-butyl ester
S-10621A4-1ML	1000ug/mL in Acetonitrile
S-10621A4-5ML	1000ug/mL in Acetonitrile
S-10621U1-1ML	100ug/mL in Toluene

N-10648-100MG	2,4,5-T butoxyethyl ester
S-10648A4-1ML	1000ug/mL in Acetonitrile
S-10648A4-5ML	1000ug/mL in Acetonitrile
S-10648U1-1ML	100ug/mL in Toluene

N-10531-100MG	2,4-D 2-ethylhexyl ester
S-10531A4-1ML	1000ug/mL in Acetonitrile
S-10531A4-5ML	1000ug/mL in Acetonitrile

METHOD CLP - SEMIVOLATILES

The U.S. EPA Contract Laboratory Program (CLP) was developed under the authority of the Comprehensive Environmental Response, Compensation, and Liability Act (CERCLA) and under the 1986 Superfund Amendments and Reauthorization Act (SARA). The analytical methods were designed for identifying and quantifying the organic compounds on the EPA's Target Compound List (TLC) in ground water, sediment, and soil from abandoned hazardous waste sites.

The following Chem Service reference standards are prepared for analyzing the volatile compounds, semi-volatile compounds, as well as pesticide compounds included on the EPA Target Compound List. Chem Service is committed to providing you with the highest quality reference materials for the CLP methods, as well as to furnishing you with high quality standards documentation

**MIXTURE #5 - PESTICIDES
EPA METHOD CLP SEMIVOLATILES**

Aldrin®	BHC (alpha isomer)	BHC (beta isomer)	BHC (delta isomer)
4,4'-DDD	4,4'-DDE	4,4'-DDT	Dieldrin
α-Endosulfan	b-Endosulfan	Endosulfan sulfate	Endrin
Endrin aldehyde	Heptachlor	Heptachlor epoxide (Isomer B)	Lindane

2000ug/mL in Toluene:Hexane (50:50) - (16) components)
M-PPHC5AC5-1ML 1mL Ampule

TUNING STANDARDS MIXTURE**EPA METHOD CLP, 625, 8250A/8270B**

Benzidine
4,4'-DDT
Decafluorotriphenylphosphine
Pentachlorophenol

1000ug/mL in Methylene chloride (4 components)
M-CLPTS1X4-1ML 1mL Ampule
M-CLPTS1X4-5ML 5mL Ampule

METHOD CLP - PESTICIDES**AROCHLOR 1016 & AROCHLOR 1260 MIXTURE****CLP PESTICIDES**

Arochlor 1016 Arochlor 1260

200ug/ml in Hexane (2 components)
M-CSMCLPPJ2-1ML **1mL Ampule****DEGRADATION PRODUCTS MIXTURE****EPA METHOD CLP 508/508.1, 525.2, 608, 617, 625, 1618, 1656, 8080A/8081, 8250A/8270B**200ug/mL 4,4'-DDD
200ug/mL 4,4'-DDE
100ug/mL Endrin aldehyde
100ug/mL Endrin ketoneVaried concentration in tert-Butylmethyl ether - (4 components)
M-DP5081T99-1ML **1mL Ampule****PESTICIDES RESOLUTION CHECK MIXTURE****CLP PESTICIDES**10ug/mL Chlordane
20ug/mL 4,4'-DDE
20ug/mL Dieldrin
10ug/mL Endosulfan I
20ug/mL Endosulfan sulfate
20ug/mL Endrin ketone
100ug/mL Methoxychlor
20ug/mL Decachlorobiphenyl
20ug/mL 2,4,5,6-Tetrachloro-m-xyleneVaried Concentration in Hexane (9 components)
M-CSCLP1J99-1ML **1mL Ampule****PESTICIDES PERFORMANCE EVALUATION MIXTURE****CLP PESTICIDES**50ug/mL Aldrin
50ug/mL Lindane
100ug/mL Dieldrin
100ug/mL 4,4-DDT
100ug/mL Endrin
50ug/mL HeptachlorVaried Concentration in Methanol (6 components)
M-CSCLP3M99-1ML **1mL Ampule****SURROGATE STANDARDS MIXTURE #1****EPA METHOD CLP, 8080A/8081/8081A**Decachlorobiphenyl
2,4,5,6-Tetrachloro-m-xylene2000ug/mL in Toluene (2 components)
M-SM8080U5-1ML **1mL Ampule****SURROGATE STANDARDS MIXTURE #2****EPA METHOD CLP, 8080A/8081/8081A**Dibutyl chlorendate
2,4,5,6-Tetrachloro-m-xylene2000ug/mL in Acetone (2 components)
M-SM8080AB5-1ML **1mL Ampule****DEGRADATION CALIBRATION MIXTURE****EPA METHOD CLP 508/508.1, 525.2, 608, 617, 625, 1618, 1656, 8080A/8081, 8250A/8270B**200ug/mL 4,4'-DDT
100ug/mL EndrinVaried concentration in tert-Butylmethyl ether - (2 components)
M-DC508T99-1ML **1mL Ampule****GPC CALIBRATION CHECK MIXTURE****CLP SEMIVOLATILES**500ug/mL Bis(2-ethylhexyl)phthalate
25,000ug/mL Corn oil
100ug/mL Methoxychlor
20ug/mL Perylene
80ug/mL Sulfur-SublimedVaried Concentration in Methylene chloride (5 components)
M-CLP14X99-1ML **1mL Ampule****PESTICIDES PERFORMANCE EVALUATION MIXTURE****CLP PESTICIDES**10ug/mL α-BHC
10ug/mL β-BHC
10ug/mL Lindane
100ug/mL 4,4-DDT
50ug/mL Endrin
250ug/mL Methoxychlor
20ug/mL Decachlorobiphenyl
20ug/mL 2,4,5,6-Tetrachloro-m-xyleneVaried Concentration in Hexane (8 components)
M-CSCLP2J99-1ML **1mL Ampule****PESTICIDES PERFORMANCE EVALUATION MIXTURE****CLP PESTICIDES**1ug/mL Aldrin
2ug/mL 4,4-DDT
2ug/mL Endrin
1ug/mL Dibutyl chlorendateVaried Concentration in Isooctane (4 components)
M-CSCLP4K99-1ML **1mL Ampule****FLORISIL CARTRIDGE CHECK STANDARDS****CLP PESTICIDES****N-10654-250MG** 2,4,5-Trichlorophenol
S-10654J1-1ML 100ug/mL in Hexane
S-10654J1-5ML 100ug/mL in Hexane
S-10654M1-1ML 100ug/mL in Methanol
S-10654M1-5ML 100ug/mL in Methanol**PESTICIDES STANDARD MIXTURE A****CLP PESTICIDES**100ug/mL BHC (alpha isomer)
200ug/mL 4,4'-DDD
200ug/mL 4,4'-DDT
200ug/mL Decachlorobiphenyl
200ug/mL Dieldrin
100ug/mL α-Endosulfan
200ug/mL Endrin
100ug/mL Heptachlor
100ug/mL Lindane (BHC gamma isomer)
1000ug/mL Methoxychlor
100ug/mL 2,4,5,6-Tetrachloro-m-xyleneVaried Concentration in Hexane (11 components)
M-CLP18J99-1ML **1mL Ampule**

ORGANOCHLORINE PESTICIDES MIXTURES**EPA METHOD CLP, 608, 508/508.1, 617, 625, 8080A/8081, 8250A/8270B**

Aldrin (TM)	b-Endosulfan	BHC (alpha isomer)	BHC (beta isomer)
BHC (delta isomer)	4,4'-DDD	4,4'-DDE	4,4'-DDT
Dieldrin	α-Endosulfan	Endosulfan sulfate	Endrin
Endrin aldehyde	Heptachlor	Heptachlor epoxide (Isomer B)	
Lindane (BHC gamma isomer)	Methoxychlor		

100 ug/ml in Toluene:Hexane (50:50) (17 components)

M-OCp8080AC1-1ML 1mL Ampule

100 ug/ml in Toluene:Hexane (50:50), less Methoxychlor (16 components)

M-PP08AC1-1ML 1mL Ampule

2000 ug/ml in Toluene:Hexane (50:50), less Methoxychlor (16 components)

M-PPHC5AC5-1ML 1mL Ampule**SURROGATE STANDARDS****CLP PESTICIDES**

BZ-209-10MG	Decachlorobiphenyl	N-11635-250MG	Dibutyl chlorendate
BZ-209J1-1ML	100ug/ml in Hexane	S-11638B5-1ML	2000ug/ml in Acetone
BZ-209J1-5ML	100ug/ml in Hexane	S-11638B5-5ML	2000ug/ml in Acetone
BZ-209B3-1ML	500ug/ml in Acetone	S-11638X4-1ML	1000ug/ml in Methylene Chloride
BZ-209B3-5ML	500ug/ml in Acetone	S-11638X4-5ML	1000ug/ml in Methylene chloride
BZ-209U4-1ML	1000ug/ml in Toluene		
BZ-209U4-5ML	1000ug/ml in Toluene	N-10542-100MG	2,4,5,6-Tetrachloro-m-xylene
BZ-209J4-1ML	1000ug/ml in Hexane	S-10542B3-1ML	500ug/ml in Acetone
BZ-209J4-5ML	1000ug/ml in Hexane	S-10542B3-5ML	500ug/ml in Acetone
		S-10542M2-1ML	200ug/ml in Methanol
		S-10542M2-5ML	200ug/ml in Methanol

USP 561 STANDARDS**CHLORINATED PESTICIDES MIXTURE****ISO 561**

Alachlor	Bromopropylate	Chlordane	Chlorfenvinphos
Chlorpyrifos	Chlorpyrifos Methyl	Chlorthal-dimethyl	DDT (mixture p,p' & o,p')
Dichlorofluanid	Dicofol	Endrin	
Heptachlor	Hexachlorobenzene	BHC (mixed isomers)	Lindane (BHC gamma isomer)
Methoxychlor	Mirex	Pendimethalin	Pentachloroanisole
Piperonyl butoxide	Procymidone	Prothiophos	Pentachloronitrobenzene
S 421	1,2,4,5-Tetrachloro-3-nitrobenzene	Tetradifon	Vinclozolin

100ug/mL in Toluene - (28 components)

M-USP561CU-1ML 1mL Ampule**PHOSPHOROUS PESTICIDES MIXTURE****ISO 561**

Acephate	Azinphos-ethyl	Azinphos-methyl	Bromophos ethyl
Bromophos methyl	Diazinon	Dichlorvos	Dimethoate
Omethoate	Ethion	Etrinfos	Fenchlorphos
Fenitrothion	Fensulfothion	Fenthion	Fonofos
Malathion	Mecarbam	Methacrifos	Methamidophos
Methidathion	Monocrotophos	Parathion	Methyl parathion
Phosalone	Phosmet	Pirimiphos-ethyl	Pirimiphos-methyl
Profenofos	Quinalphos		

100ug/mL in Toluene - (28 components)

M-USP561PU-1ML 1mL Ampule**ISO 6468 - ORGANOCHLORINE INSECTICIDES, PCB'S AND CHLOROBENZENE**

ISO 6468 is for determining certain organochlorine insecticides, polychlorinated biphenyls (PCBs) and chlorobenzenes (except the mono- and dichlorobenzenes) in drinking water, ground water, surface waters and waste waters. The method is applicable to samples containing up to 0,05 g/l of suspended solids.

ORGANOCHLORINE INSECTICIDE MIXTURE**ISO 6468**

Aldrin	BHC (alpha isomer)	BHC (beta isomer)	BHC (delta isomer)
4,4'-DDD	o,p'-DDD	o,p'-DDE	4,4'-DDE
o,p'-DDT	4,4'-DDT	Dieldrin	α-Endosulfan
b-Endosulfan	Endrin	Heptachlor	Heptachlor epoxide (Isomer B)
Methoxychlor	Lindane (BHC gamma isomer)	trans-Heptachlor epoxide	

10ug/mL in Hexane - (19 components)

M-ISO6468PSTJ10-1ML 1mL Ampule

ISO 10695 - NITROGEN AND PHOSPHOROUS COMPOUNDS

ISO 10695 is a general purpose method applicable for the identification and measurement of select organic nitrogen and phosphorus compounds in drinking water, ground water, surface waters, and waste water by liquid-liquid extraction or liquid-solid extraction phase using gas chromatography (GC).

ORGANIC NITROGEN/PHOSPHORUS PESTICIDE MIXTURE

ISO 10695

Simazine Bladex Parathion	Propazine Methyl parathion Terbutylazine	Trifluralin Pendimethalin Metazachlor	Atrazine Vinclozolin Sebutylazin
---------------------------------	--	---	--

10ug/mL in Acetone (12 components)
M-ISO10695PST1B10-1ML **1mL Ampule**

ISO 11369 - PLANT TREATMENT AGENTS

ISO 11369 is a general purpose method applicable for the identification and measurement of select plant treatment agents by solid-liquid extraction using high performance liquid chromatography (HPLC).

HERBICIDE MIXTURE

ISO 11369

Simazine Atrazine Diuron Velpar Sebutylazin	Chlorotoluron Metoxuron Linuron Terbutylazine	Monolinuron Bladex Isoproturon Atrazine desethyl	Methabenzthiazuron Metobromuron Metolachlor Metazachlor
---	--	---	--

10ug/mL in Ethyl acetate (17 components)
M-ISO11369PS1H10-1ML **1mL Ampule**

ISO 15913 - PHENOXYALKANOIC HERBICIDES

ISO 15913 is a general purpose method applicable for the identification and measurement of select phenoxyalkanoic herbicides, including bentazones and hydroxybenzonitriles in ground and drinking water by solid phase extraction and derivatization using gas chromatography coupled with mass spectrometry (GC/MS).

PHENOXYALKANOIC HERBICIDE ACID MIXTURE

ISO 15913

Clopyralid Picloram Triclopyr Fluroxypyr Pentachlorophenol	Dicamba 2,3,6-Trichlorobenzoic acid Diclofop acid Benazolin Fenoxaprop	Tetrachloroterephthalic acid 4,6-Dinitro-o-cresol Fluazifop Imazapyr 3,5-Diiodo-4-hydroxybenzonitrile
--	--	---

10ug/mL in Acetone (15 components)
M-ISO15913FAB10-1ML **1mL Ampule**

INTERNATIONAL STANDARDS - CANADA

CANADIAN DRINKING WATER TRIAZINE MIXTURE

Aalachlor Atrazine Cyanazine Metolachlor	Metribuzin Prometryne Simazine
---	--------------------------------------

500ug/mL in Acetone (7 components)
M-CSCAN1B3-1ML **1mL Ampule**

CANADIAN DRINKING WATER CARBAMATE MIXTURE

Aldicarb Carbaryl Triallate	Bendiocarb Carbofuran
-----------------------------------	--------------------------

100ug/mL in Acetonitrile (5 components)
M-CSCAN3A1-1ML **1mL Ampule**

CANADIAN DRINKING WATER PHENOXYACID HERBICIDE MIXTURE

Bromoxynil 2,4-Dichlorophenol	Pentachlorophenol 2,3,4,6-Tetrachlorophenol	2,4-D Diclofop methyl	Picloram 2,4,6-Trichlorophenol	Dicamba	2,4,5-T Dinoseb
----------------------------------	--	--------------------------	-----------------------------------	---------	--------------------

1000ug/mL in Acetone (11 components)
M-CSCAN2B4-1ML **1mL Ampule**

CANADIAN DRINKING WATER CHLORINATED PESTICIDES MIXTURE

Aldrin Methoxychlor 4,4'-DDT	Heptachlor 4,4'-DDE cis-Chlordane	Lindane Oxychlordane Dieldrin	Heptachlor epoxide 2,4'-DDT trans-Chlordane	2,4'-DDE Trifluralin
------------------------------------	---	-------------------------------------	---	-------------------------

200ug/mL in Hexane:Toluene (50:50) (14 components)
M-CSCAN4AC2-1ML **1mL Ampule**

CANADIAN DRINKING WATER ORGANOPHOSPHOROUS PESTICIDES MIXTURE

Abate Phorate	Malathion Dimethoate	Chlorpyrifos Terbufos	Parathion Guthion	Diazinon
------------------	-------------------------	--------------------------	----------------------	----------

1000ug/mL in Acetonitrile (9 components)
M-CSCAN5B4-1ML **1mL Ampule**

DIN 38407 - GERMAN STANDARD METHODS

DIN 38407 Glyphosate and AMPA MIX (Part 22) is a German standard method applicable for the identification and measurement of glyphosate and aminomethyl phosphonic acid (AMPA) in water, waste water, and sludge by post-column derivation using high performance liquid chromatography (HPLC) with fluorescence detection.

DIN 38407 Pesticide Free Acid Mixture & Methyl Esters Mixture (Part 14) is a German standard method applicable for the identification and measurement of phenoxyalkancarboxylic acids in water, waste water and sludge by solid-liquid extraction and derivation using gas chromatography coupled with mass spectrometry (GC/MS).

DIN 38407 Pesticide Mixture (Part 2) is a German standard method application for the identification and measurement of low-volatility halogenated hydrocarbons in water, waste water and sludge using gas chromatography (GC).

DIN 38407 Pesticide Mixture (Part 11) is a German standard method application for the identification and measurement of select organic plant protecting agents in water, waste water and sludge.

GLYPHOSATE AND AMPA MIXTURE (DIN 38407-22)

Aminomethyl phosphonic acid
Glyphosate

100ug/mL in Water (2 components)
M-DIN38407GAF1-1ML 1mL Ampule

PESTICIDE FREE ACID MIXTURE (DIN 38407-14)

Dichlorprop MCPA	Mecoprop 2,4-D	Silvex MCPB	2,4,5-T (TM) 2,4-DB
---------------------	-------------------	----------------	------------------------

500ug/mL in Acetonitrile (8 components)
M-DIN38407FAM3-1ML 1mL Ampule

METHYL ESTERS MIXTURE

2,4-DB methyl ester	(2,4,5-Trichlorophenoxy)acetic acid methyl ester	2,4-D methyl ester
4-Chloro-o-tolyloxyacetic acid methyl ester	Mecoprop methyl ester	Silvex methyl ester
Dichlorprop methyl ester	MCPB methyl ester	

500ug/mL in Hexane (8 components)
M-DIN38407ME1J3-1ML 1mL Ampule

PESTICIDE MIXTURE (DIN 38407-2)

trans-Heptachlor epoxide	Aldrin	BHC (alpha isomer)	BHC (beta isomer)
b-Endosulfan	o,p'-DDE	4,4'-DDT	Lindane (BHC gamma isomer)
Dieldrin	Endrin	Methoxychlor	4,4'-DDD
4,4'-DDE	Heptachlor	o,p'-DDT	α-Endosulfan
Heptachlor epoxide (Isomer B)			

10ug/mL in Hexane (17 components)
M-DIN38407PS1J10-1ML 1mL Ampule

PESTICIDE MIXTURE (DIN 38407-11)

Simazine	Propazine	Monuron	Chlorotoluron
Trifluralin	Alachlor	Atrazine	Metoxuron
Bladex	Metobromuron	Pendimethalin	Chlorfenvinphos
Vinclozolin	Metolachlor	Parathion	Terbutylazine
Metazachlor	Sebutylazine	2,4,5-T	2,4,5-T-amine
2,4-D	2,4-D sec-butyl ester	4-Chloro-o-tolyloxyacetic acid	

10ug/mL in Acetonitrile (23 components)
M-DIN3840711A10-1ML 1mL Ampule

PESTICIDE MIXTURE (DIN 38407-11)

Simazine	Propazine	Monuron	Chlorotoluron
Trifluralin	Alachlor	Atrazine	Metoxuron
Bladex	Metobromuron	Pendimethalin	Chlorfenvinphos
Vinclozolin	Metolachlor	Parathion	Terbutylazine
Metazachlor	Sebutylazine	2,4,5-T	4-Chloro-o-tolyloxyacetic acid
2,4-D			

10ug/mL in Acetonitrile (21 components)
M-DIN38407PS2A10-1ML 1mL Ampule

EUROPEAN REGULATION STANDARDS

EUROPEAN ORGANOCHLORINE PESTICIDE MIXTURE #1

Heptachlor epoxide (Isomer B)	1,2,4,5-Tetrachloro-3-nitrobenzene	Hexachlorobenzene
Pyrazon	trans-Heptachlor epoxide	Aldrin
BHC (alpha isomer)	BHC (beta isomer)	BHC (delta isomer)
o,p'-DDE	4,4'-DDT	o,p'-DDD
Lindane (BHC gamma isomer)	Endrin	4,4'-DDD
4,4'-DDE	Heptachlor	o,p'-DDT
Tetrachlorvinphos		

10ug/mL in Toluene (19 components)
M-EUPESTMIX1U10-1ML 1mL Ampule

EUROPEAN ORGANOPHOSPHOROUS PESTICIDE MIXTURE #2

100 µg/mL	Methamidophos	150 µg/mL	Phorate sulfone
150 µg/mL	Phorate	150 µg/mL	Disulfoton
100 µg/mL	Fenthion	100 µg/mL	Thiometon
100 µg/mL	Carbophenothion	10 µg/mL	Tetrachlorvinphos
100 µg/mL	Phorate Sulfone		

Varied concentration in Toluene

(9 components)

M-EUPESTMIX2U99-1ML**1mL Ampule****EUROPEAN ORGANOPHOSPHOROUS PESTICIDE MIXTURE #3**

100ug/mL	Quinalphos	100ug/mL	Heptenophos
100ug/mL	Triazophos	100ug/mL	Chlorpyrifos
100ug/mL	Chlorfenvinphos	100ug/mL	Dimethoate
100ug/mL	Dichlorvos	10ug/mL	Tetrachlorvinphos

Varied concentration in Toluene

(9 components)

M-EUPESTMIX3U99-1ML**1mL Ampule****EUROPEAN PESTICIDE MIXTURE #4**

100 µg/mL	Pendimethalin	500 µg/mL	Fenvalerate
500 µg/mL	Cypermethrin	500 µg/mL	Permethrin
500 µg/mL	Deltamethrin	500 µg/mL	Cyfluthrin
100 µg/mL	Tefluthrin	500 µg/mL	lambda-Cyhalothrin
10 µg/mL	Tetrachlorvinphos	100 µg/mL	Dicloran

Varied concentration in Toluene

(10 components)

M-EUPESTMIX4U99-1ML**1mL Ampule****EUROPEAN PESTICIDE MIXTURE #5**

Endosulfan sulfate	Alachlor	Bromopropylate	Pentachloronitrobenzene
trans-Permethrin	Zectran	Cycloprate	Pentachloroaniline
Deltamethrin	cis-Permethrin	Carbophenothion	Pyrethrum

10ug/mL in Toluene

(12 components)

M-EUPESTMIX5U10-1ML**1mL Ampule****EUROPEAN ORGANOPHOSPHOROUS PESTICIDE MIXTURE #6**

Malathion	Fenitrothion	Phosalone	
Chlorpyrifos	Pirimiphos-methyl	Methyl parathion	Diazinon
Chlorfenvinphos	Chlorpyrifos Methyl	Ethion	Parathion
Dichlorvos	Carbophenothion	Azinphos-methyl	Fonofos
Methidathion			

10ug/mL in Toluene

(16 components)

M-EUPESTMIX6U10-1ML**1mL Ampule****EUROPEAN PESTICIDE MIXTURE #7**

Endosulfan sulfate	1,2,4,5-Tetrachloro-3-nitrobenzene	Fenitrothion	Phosalone
Methyl parathion	4,4'-DDT	Deltamethrin	Chlorpyrifos Methyl
Ethion	Lindane (BHC gamma isomer)	Dieldrin	Dichlorvos
Pentachloronitrobenzene			

10ug/mL in Toluene

(13 components)

M-EUPESTMIX7U10-1ML**1mL Ampule****EUROPEAN PESTICIDE MIXTURE #8**

Aldicarb	Dichlorprop	Simazine	Chlorotoluron
Atrazine	Metoxuron	Bladex	Bentazon
Metobromuron	b-Endosulfan	Isoproturon	1,1-Dichloropropene
Lindane (BHC gamma isomer)	Terbutylazine	Metazachlor	Mecoprop
Sebutylazin	4-Chloro-o-tolxyloxyacetic acid	2,4-D	2,4-D sec-butyl ester
a-Endosulfan			

10ug/mL in Ethyl acetate

(21 components)

M-EUPESTMIX8H10-1ML**1mL Ampule****EUROPEAN ORGANOCHLORINE PESTICIDE MIXTURE #9**

Heptachlor epoxide (Isomer B)	Hexachlorobenzene	Aldrin (TM)	BHC (alpha isomer)
BHC (beta isomer)	b-Endosulfan	4,4'-DDT	Lindane (BHC gamma isomer)
Dieldrin	Endrin	Methoxychlor	4,4'-DDD
4,4'-DDE	Heptachlor	o,p'-DDT	a-Endosulfan

10ug/mL in Cyclohexane

(16 components)

M-EUPESTMIX9E10-1ML**1mL Ampule****EUROPEAN NITROGEN HERBICIDE MIXTURE #10**

Simazine	Propazine	Pyrazon	Atrazine
Metoxuron	Bromacil	Terbutylazine	Atrazine desethyl

10ug/mL in Acetonitrile

(8 components)

M-EUPESTMIX10A11-1ML**1mL Ampule**

EUROPEAN ORGANOCHLORINE PESTICIDE MIXTURE #11

Heptachlor epoxide (Isomer B)
BHC (delta isomer)
Dieldrin
4,4'-DDE

Mirex
o,p'-DDE
Endrin

BHC (alpha isomer)
o,p'-DDD
Methoxychlor

BHC (beta isomer)
Lindane (BHC gamma isomer)
4,4'-DDD

10ug/mL in Isooctane
M-EUPESTMIX11K10-1ML

(13 components)
1mL Ampule

EUROPEAN PHENOXYCARBOXYLIC ACID HERBICIDE MIXTURE #12

Dichlorprop
MCPB

Mecoprop
2,4-DB

2,4-D
4-Chloro-o-tolxyloxyacetic acid

10ug/mL in Acetonitrile
M-EUPESTMIX12A10-1ML

(6 components)
1mL Ampule

EUROPEAN BHC ISOMER PESTICIDE MIXTURE #14

BHC (alpha isomer)
BHC (epsilon isomer) Solution

BHC (beta isomer)

BHC (delta isomer)

Lindane (BHC gamma isomer)

1ug/mL in Cyclohexane
M-EUPESTMIX14E10-1ML

(5 components)
1mL Ampule

CANNABIS STANDARDS

Cannabis Standards are a general purpose method applicable for the identification and measurement of select pesticide residue in cannabis products in Colorado, Nevada, and Oregon, respectively. Note: Piperonyl butoxide and pyrethrin may be used in cannabis. Due to misuse, they have been incorporated into these mixtures.

COLORADO PESTICIDE STANDARD MIXTURE 1

Abamectin
Imazalil
Permethrin
Tebuconazole

Azoxystrobin
Imidacloprid
Spinosad

Bifenazate
Malathion
Spiromesifen

Etoazole
Myclobutanil
Spirotetramat

100ug/mL in Acetonitrile
M-COPESTMIX1A1-1ML

(13 components)
1mL Ampule

NEVADA PESTICIDE STANDARD MIXTURE 1

Abamectin
Captan
Dimethomorph
Fludioxonil
Pentachloronitrobenzene
Spinosad

Acequinocyl
beta-Cyfluthrin
Etoazole
Imidacloprid
Piperonyl butoxide
Spirotetramat

Bifenazate
Cypermethrin
Fenhexamid
Myclobutanil
Pyrethrum
Thiamethoxam

Bitenthrin
Daminozide
Flonicamid
Paclobutrazol
Spinetoram
Trifloxystrobin

100ug/mL in Acetonitrile
M-NVPESTMIX1A1-1ML

(24 components)
1mL Ampule

OREGON PESTICIDE STANDARD MIXTURE 1

Abamectin
Aldicarb
Bifenthrin
Chlorantraniliprole
Cypermethrin

Acephate
Azoxystrobin
Boscalid
Chlorfenapyr
Daminozide

Acequinocyl
beta-Cyfluthrin
Carbaryl
Chlorpyrifos
Diazinon

Acetamiprid
Bifenazate
Carbofuran
Clofentezine
Dichlorvos

100ug/mL in Acetonitrile
M-ORPESTMIX1A1-1ML

(20 components)
1mL Ampule

OREGON PESTICIDE STANDARD MIXTURE 2

Dimethoate
Fenoxycarb
Fludioxonil
Kresoxim-methyl
Methomyl

Etofenprox
Fenpyroximate
Hexythiazox
Malathion
Methyl parathion

Propos
Fipronil
Imazalil
Metalaxyl
MGK 264 (TM)

Etoazole
Flonicamid
Imidacloprid
Methiocarb
Myclobutanil

100ug/mL in Acetonitrile
M-ORPESTMIX2A1-1ML

(20 components)
1mL Ampule

OREGON PESTICIDE STANDARD MIXTURE 3

Naled
Phosmet
Propoxur
Spiromesifen
Thiacloprid

Oxamyl
Piperonyl butoxide
Pyrethrum
Spirotetramat
Thiamethoxam

Paclobutrazol
Prallethrin
Pyridaben
Spiroxamine
Trifloxystrobin

Permethrin
Propiconazole
Spinosad
Tebuconazole

100ug/mL in Acetonitrile
M-ORPESTMIX3A1-1ML

(19 components)
1mL Ampule



Pesticide & Metabolite Standards Listing by Part Number



MET-10688A-250MG



2,6-Dichlorobenzoic acid



50-30-6



250MG

a Part Number

b Description

c Chemical Abstracts Registry
Number

d Size

Part #	Description	CAS #	Size	Part #	Description	CAS #	Size
MET-10688A-250MG	2,6-Dichlorobenzoic acid	50-30-6	250MG	MET-11963DA1-1ML	Fensulfiothion Oxon Solution	6552-21-2	1ML
MET-10688B-1G	2,6-Dichlorobenzamide	2008-58-4	1G	MET-11964A-100MG	Fenthion-sulfone	3761-42-0	100MG
MET-10688C-250MG	2,6-Dichlorobenzoic acid methyl ester	14920-87-7	250MG	MET-11964BX1-1ML	Fenthion sulfoxide Solution	3761-41-9	1ML
MET-10794A-100MG	loxynil-methyl	3336-40-1	100MG	MET-11964CK1-1ML	Fenthion-O-analog Solution	6552-12-100	1ML
MET-10818A-1G	4-Chloro-2-methylphenol	1570-64-5	1G	MET-11964D-100MG	Fenthion sulfoxide	3761-41-9	100MG
MET-11009A-10MG	Acetamidiprid-N-desmethyl	190604-92-3	10MG	MET-11964E-100MG	Fenthion-oxon-sulfone	14086-35-2	100MG
MET-11013-50MG	2-Chloro-N-(2-ethyl-6-methylphenyl) acetamide	32428-71-0	50MG	MET-11964FA1-1ML	Fenthion oxon sulfoxide Solution	6552-13-2	1ML
MET-11013A-25MG	Acetochlor ESA sodium salt	187022-11-3	25MG	MET-11964G-50MG	4-Methylsulfonylphenol	14763-60-1	50MG
MET-11013AM1-1ML	Acetochlor ESA sodium salt Solution	187022-11-3	1ML	MET-11971A-50MG	Fipronil sulfone	120068-36-2	50MG
MET-11013B-10MG	Acetochlor OA	194992-44-4	10MG	MET-11971AA1-1ML	Fipronil sulfone Solution	120068-36-2	1ML
MET-11013BM1-1ML	Acetochlor OA Solution	194992-44-4	1ML	MET-11971B-25MG	Fipronil sulfide	120067-83-6	25MG
MET-11027A-10MG				MET-11971BA1-1ML	Fipronil sulfide Solution	120067-83-6	1ML
MET-11043A-100MG	2-Chloro-2',6'-diethylacetanilide	6967-29-9	100MG	MET-11976A-100MG	4-(Trifluoromethyl)nicotinic acid	158063-66-2	100MG
MET-11043B-100MG	Alachlor acetate	166407-15-4	100MG	MET-11986A-25MG	Flufenacet ESA sodium salt	947601-87-8	25MG
MET-11043C-500MG	2,6-Diethylaniline	579-66-8	500MG	MET-11986AM1-1ML	Flufenacet ESA sodium salt Solution	947601-87-8	1ML
MET-11043D-25MG	Alachlor ESA sodium salt	140939-15-7	25MG	MET-11986B-10MG	Flufenacet OA	201668-31-7	10MG
MET-11043DM1-1ML	Alachlor ESA sodium salt Solution	140939-15-7	1ML	MET-11986BM1-1ML	Flufenacet OA Solution	201668-31-7	1ML
MET-11043E-10MG	Alachlor OA	171262-17-2	10MG	MET-12004A-100MG	Methyl o-trifluoromethylbenzoate	344-96-7	100MG
MET-11043EM1-1ML	Alachlor OA Solution	171262-17-2	1ML	MET-12004C-250MG	alpha, alpha, alpha-Trifluoro-3'-hydroxy-O-toluanilide		250MG
MET-11044A-10MG	Aldicarb oxime	1646-75-9	10MG	MET-12007A-100MG	Phthalamic acid	88-97-2	100MG
MET-11065-100MG	Aminocarb-phenol	14143-25-0	100MG	MET-12017A-100MG	(Dimethoxythiophosphorylthio)acetic acid	1113-01-5	100MG
MET-11068A-50MG	N-(2,4-Dimethylphenyl)formamide	60397-77-5	50MG	MET-12066B-100MG	6-Chloronicotinic acid	5326-23-8	100MG
MET-11068B-10MG	Amiratz Metabolite BTS 27271	33089-74-6	10MG	MET-12066BM1-1ML	6-Chloronicotinic acid Solution	5326-23-8	1ML
MET-11085A-100MG	1-(4-tert-Butylphenoxy)propan-2-ol	2416-30-0	100MG	MET-12133A-100MG	Aminomethyl phosphonic acid	1066-51-9	100MG
MET-11106B-25MG	Atrazine desethyl	6190-65-4	25MG	MET-12135A-500MG	Iminobis (methylphosphonic acid)	6711-48-4	500MG
MET-11106D-50MG	Atrazine desethyl-2-hydroxy	19988-24-0	50MG	MET-12139AU1-1ML	Azinphos-methyl oxon Solution	961-22-8	1ML
MET-11106DAM1-1ML	Atrazine-desethyl-2-hydroxy Solution	19988-24-0	1ML	MET-12159A-250MG	Tetrachlorohydroquinone	87-87-6	250MG
MET-11106E-50MG	Atrazine desethyl desisopropyl-2-hydroxy	645-92-1	50MG	MET-12159B-100MG	Pentachlorothiophenol	133-49-3	100MG
MET-11106F-100MG	Atrazine-2-hydroxy	2163-68-0	100MG	MET-12200AA0-1ML	Imazathiazabenz (free acid) Solution	100728-84-5	1ML
MET-11128A-250MG	2-Isopropoxyphenol	4812-20-8	250MG	MET-12206-10MG	Imidacloprid olefin	115086-54-9	10MG
MET-11135A-500MG	2,2-Dimethyl-1,3-benzodioxole-4-ol		500MG	MET-12206A-50MG	Imidacloprid urea	120868-66-8	50MG
MET-11138A-250MG	2-Aminobenzimidazole	934-32-7	250MG	MET-12207AU1-1ML	Phosmet-oxon Solution	3735-33-9	1ML
MET-11138B-50MG	STB	132827-25-9	500MG	MET-12220B-250MG	3,5-Dichloroaniline	626-43-7	250MG
MET-11142C-50MG	2-Amino-N-isopropyl benzamide	30391-89-0	50MG	MET-12249AJ1-1ML	Isofenphos-des-N-isopropyl Solution	15205-08-7	1ML
MET-11345A-50MG	Bromoxynil-heptanoate	56634-95-8	50MG	MET-12249BJ1-1ML	Isofenphos-des-N-isopropyl-O-analogue Solution	31120-83-9	1ML
MET-11345B-100MG	Bromoxynil-methyl ether	3336-39-8	100MG	MET-12249C-100MG	Isofenphos-des-N-isopropyl	15205-08-7	100MG
MET-11353A-25MG	Butachlor ESA sodium salt	1173022-75-7	25MG	MET-12279A-250MG	Monodesmethyl isoproturon	34123-57-4	250MG
MET-11353AM1-1ML	Butachlor ESA sodium salt Solution		1ML	MET-12287A-100MG	1-(Methylsulfonyl-4-trifluoromethyl phenyl)-2-cyano-3-cyclop	143701-75-1	100MG
MET-11399A-250MG	cis-1,2,3,6-Tetrahydrophthalimide	1469-48-3	250MG	MET-12287B-100MG	2-Methylsulfonyl-4-trifluoromethyl benzoic acid	142994-06-7	100MG
MET-11405B-250MG	2,3-Dihydro-2,2-dimethylbenzofuran-7-ol	1563-38-8	250MG	MET-12315A-50MG	Leptophos oxon	25006-32-0	50MG
MET-11405DA1-1ML	2,3-Dihydro-2,2-dimethyl benzofuran-3,7-diol Solution	17781-15-6	1ML	MET-12346AJ1-1ML	Isomalathion Solution	3344-12-500	1ML
MET-11405E-100MG	Carbofuran-3-keto	16709-30-1	100MG	MET-12346BT1-1ML	Malathion monocarboxylic acid Solution	35884-76-5	1ML
MET-11405F-50MG	Carbofuranphenol-3-keto	17781-16-7	50MG	MET-12346C-100MG	Malathion-O-analog	1634-78-2	100MG
MET-11405FM1-1ML	Carbofuranphenol-3-keto Solution	17781-16-7	1ML	MET-12346D-100MG	Malathion monocarboxylic acid	35884-76-5	100MG
MET-11408AM1-1ML	Carbophenothion methyl Solution	953-17-3	1ML	MET-12348-250MG	MCPA sodium salt	3653-48-3	250MG
MET-11408B-50MG	Carbophenothion methyl-O-analog	7332-32-3	50MG	MET-12371A-250MG	N-Methylpiperidine	626-67-5	250MG
MET-11408BH10-1ML	Carbophenothion methyl-O-analog Solution	7173-84-4	1ML	MET-12383A-100MG	Pentachloroaniline	527-20-8	100MG
MET-11408C-100MG	Carbophenothion sulfoxide	17297-40-4	100MG	MET-12392A-250MG	Methyl isothiocyante	556-61-6	250MG
MET-11408D-25MG	Carbophenothion sulfone	16662-85-4	25MG	MET-12398A-50MG	Methiocarb sulfone	2179-25-1	50MG
MET-11423A-25MG	Chlorbenside sulfone	7082-99-7	25MG	MET-12398BU1-1ML	Methiocarb sulfoxide Solution	2599-11-300	1ML
MET-11423AJ1-1ML	Chlorbenside sulfone Solution	7082-99-7	1ML	MET-12403A-100MG	o,p'-Methoxychlor	30667-99-3	100MG
MET-11423B-100MG	Chlorbenside sulfoxide	7047-28-1	100MG	MET-12403AM1-1ML	o,p'-Methoxychlor Solution	72-43-5	1ML
MET-11423BX1-1ML	Chlorbenside sulfoxide Solution	7047-28-1	1ML	MET-12403B-50MG	p,p'-Methoxychlor-olefin	2132-70-9	50MG
MET-11425AM1-1ML	Oxychlorane Solution	27304-13-8	1ML	MET-12403BM1-1ML	p,p'-Methoxychlor-olefin Solution	2132-70-9	1ML
MET-11425D-250MG	Chlordane	3734-48-3	250MG	MET-12478A-50MG	2-[(2-Ethyl-6-methylphenyl)amino]-1-propanol	61520-53-4	50MG
MET-11426A-250MG	4-Chloro-2-methylaniline	95-69-2	250MG	MET-12478C-25MG	Metolachlor ESA sodium salt	947601-85-6	25MG
MET-11426AM1-1ML	4-Chloro-2-methylaniline Solution	95-69-2	1ML	MET-12478CM1-1ML	Metolachlor ESA sodium salt Solution	947601-85-6	1ML
MET-11442A-250MG	4,4'-Dichlorobenzhydrol	90-97-1	250MG	MET-12478D-10MG	Metolachlor OA	152019-73-3	10MG
MET-11442B-100MG	4,4'-Dichlorobenzophenone	90-98-2	100MG	MET-12478DM1-1ML	Metolachlor OA Solution	152019-73-3	1ML
MET-11455A-100MG	1-(3-Chloro-4-methylphenyl)-3-methylurea	22175-22-0	100MG	MET-12497A-500MG	N,N'-Bis(p-chlorophenyl)urea	1219-99-4	500MG
MET-11459A-250MG	3,5,6-Trichloro-2-pyridinol	6515-38-4	250MG	MET-12507A-250MG	N-(1-Naphthyl)phthalimide	5333-99-3	250MG
MET-11459B-100MG	Chlorpyrifos Oxon	5598-15-2	100MG	MET-12659BJ1-1ML	6-Chloro-2-picolinic acid methyl ester Solution	6636-55-1	1ML
MET-11460B-50MG	Chlorpyrifos methyl-O-analog	5598-52-7	50MG	MET-12706A-10MG	o,p'-DDD olefin	3424-82-6	10MG
MET-11462A-100MG	Monomethyl tetrachloroterephthalate	887-54-7	100MG	MET-12708A-100MG	p,p'-DDA	83-05-6	100MG
MET-11490A-250MG	2-Chlorobenzonitrile	873-32-5	250MG	MET-12727A-100MG	2-(4-tert-Butylphenoxy)cyclohexanol	1942-71-8	100MG
MET-11507D-50MG	Coumaphos-O-analog	321-54-0	50MG	MET-12727B-250MG	1,2-Cyclohexanediol (cis & trans mixture)	931-17-9	250MG
MET-11521A-1G	2,4-Dichloroaniline	554-00-7	1G	MET-12731A-50MG	Oxadiazon-hydroxy	39807-19-7	50MG
MET-11522A-250MG	N-Ethylcyclohexylamine	5459-93-8	250MG	MET-12731AU10-1ML	Oxadiazon-hydroxy Solution	39807-19-7	1ML
MET-11543A-50MG	Cyhalofop-amide		50MG	MET-12731BM1-1ML	Pivaloyl 2-(2,4-dichloro-5-hydroxyphenyl)hydrazide	39807-20-0	1ML
MET-11579A-10MG	3-(2,2-Dibromovinyl)-2,2-dimethyl-(1-cyclopropyl)carb. acid(cis)	53179-78-5	10MG	MET-12735-10MG	Oxamyl oxime	30558-43-1	10MG
MET-11579BM1-1ML	3-(2,2-Dibromovinyl)-2,2-dimethyl-(1-cyclopropyl)carb. acid(cis)	53179-78-5	1ML	MET-12810B-10MG	p,p'-DDMU	1022-22-6	10MG
MET-11579C-10MG	cis-Decamethrinic acid methyl ester	113830-50-5	10MG	MET-12848A-10MG	3-(2,2-Dichlorovinyl)-2,2-dimethyl-(1-cyclopropyl)carb. acid	55701-05-8	10MG
MET-11621A-100MG	Diazinon-O-analog	962-58-3	100MG	MET-12848AM1-1ML	3-(2,2-Dichlorovinyl)-2,2-dimethyl-(1-cyclopropyl)carb. acid	55701-05-8	1ML
MET-11621B-250MG	2-Isopropyl-6-methyl-4-pyrimidinol	2814-20-2	250MG	MET-12848B-100MG	m-Phenoxybenzyl alcohol	13826-35-2	100MG
MET-11621C-100MG	Diethyl phosphate	598-02-7	100MG	MET-12848C-100MG	3-Phenoxybenzoic acid	3739-38-6	100MG
MET-11621D-10MG	4-(Ethylthio)-6-methyl-2-(1-methylethyl)-pyrimidine	77738-92-2	10MG	MET-12890A-100MG	2-Chloro-4-(methylsulfonyl)benzoic acid	53250-83-2	100MG
MET-11656B-50MG	3,6-Dichloro-2-hydroxy benzoic acid	3401-80-7	50MG	MET-12950A-100MG	Sebuthyazine-desethyl	37019-18-4	100MG
MET-11658A-250MG	2-Chloro-4-nitrophenol	619-08-9	250MG	MET-13035AU1-1ML	Phorate-oxon Solution	2600-69-3	1ML
MET-11677B-100MG	Diclofop acid	40843-25-2	100MG	MET-13035B-100MG	Phorate sulfone	2588-04-700	100MG
MET-11721C-100MG	Difenzoquat-desmethyl		100MG	MET-13035BJ1-1ML	Phorate sulfone Solution	2588-04-700	1ML
MET-11722A-100MG	p-Chlorophenylurea	140-38-5	100MG	MET-13035C-100MG	Phorate sulfoxide	2588-03-600	100MG
MET-11754A-25MG	Dimethachlor ESA sodium salt		25MG	MET-13063A-250MG	Pirimiphos-ethyl-oxon	36378-61-7	250MG
MET-11754AM1-1ML	Dimethachlor ESA sodium salt Solution		1ML	MET-13064A-100MG	2-Diethylamino-6-methyl pyrimidin-4-ol	42487-72-9	100MG
MET-11755A-25MG	Dimethenamid ESA sodium salt	205939-58-8 (acid)	25MG	MET-13064B-250MG	Pirimiphos-methyl-oxon	64709-45-1	250MG
MET-11755AM1-1ML	Dimethenamid ESA sodium salt Solution	205939-58-8 (acid)	1ML	MET-13064C-100MG	2-Ethylamino-6-methyl-4(3H)-Pyrimidinone	5734-69-0	100MG
MET-11755B-10MG	Dimethenamid OA		10MG	MET-13064D-100MG	Pirimiphos-methyl-N-desethyl	67018-59-1	100MG
MET-11755BM1-1ML	Dimethenamid OA Solution		1ML	MET-13064E-100MG	N,N-diethyl-4-methoxy-6-methyl-2-Pyrimidinamine	111697-06-4	100MG
MET-11819A-100MG	Disulfoton sulfone	2497-06-500	100MG	MET-13068AA1-1ML	3-(difluoromethyl)-1-methyl-1H-Pyrazole-4-carboxamide	925689-10-7	1ML
MET-11819B-100MG	Disulfoton sulfoxide	2497-07-600	100MG	MET-13068BA1-1ML	3-(Difluoromethyl)-1-methyl-1H-pyrazole-4-carboxylic acid	176969-34-9	1ML
MET-11827A-100MG	N,N'-Bis(3,4-dichlorophenyl)urea	4300-43-0	100MG	MET-13079A-100MG	De-Bromo-halphenprol		10MG
MET-11827B-100MG	3-(4-Chlorophenyl)methyl urea	5352-88-5	100MG	MET-13097A-500MG	4-Bromo-2-chlorophenol	3964-56-5	500MG
MET-11827C-250MG	1-(3,4-Dichlorophenyl)-3-methylurea	3567-62-2	250MG	MET-13104A-25MG	Propachlor ESA sodium salt	947601-88-9	25MG
MET-11840A-250MG	Creatine	57-00-1	250MG	MET-13104AM1-1ML	Propachlor ESA sodium salt Solution	947601-88-9	1ML
MET-11848A-50MG	Endosulfan alcohol	2157-19-9	50MG	MET-13104B-10MG	Propachlor OA	70628-36-3	10MG
MET-11848B-50MG	Endosulfan ether	3369-52-6	50MG	MET-13104BM1-1ML	Propachlor OA Solution	70628-36-3	1ML
MET-11848C-50MG	Endosulfan lactone	3868-61-9	50MG	MET-13108AAA1-1ML	Propazine-2-hydroxy Solution	7374-53-0	1ML
MET-11872A-10MG	Ethiofencarb sulfoxide	53380-22-6	10MG	MET-13140A-100MG	N-Nitrosopendimethalin	68897-50-7	100MG
MET-11872B-10MG	Ethiofencarb sulfone	53380-23-7	10MG	MET-13142A-200MG	1-(4-Chlorophenyl)-1H-pyrazol-4-ol	17458-30-1	100MG
MET-11872BA1-1ML	Ethiofencarb sulfone Solution	53380-23-7	1ML	MET-13144A-10MG	Pyraflufen	129630-17-7	10MG
MET-11941AM1-1ML	Etofenprox-desethyl Solution	80854-21-3	1ML	MET-13174A-50MG	Quizalofop	76578-12-6	50MG
MET-11944A-50MG	Famphur-O-analog	960-25-8	50MG	MET-13174B-50MG	Quizalofop methyl	76578-13-7	50MG
MET-11944B-10MG	p-Hydroxy-N,N-dimethylbenzene sulfonamide	15020-57-2	10MG	MET-13176A-100MG	5-(Phenylmethyl)-3-furanmethanol	20416-09-5	100MG
MET-11946-100MG	Fenamidone Metabolite	332855-88-6	100MG	MET-13186A-10MG	4-tert-Butyl-2-chlorophenol	98-28-2	10MG
MET-11947A-50MG	Fenamiphos sulfoxide	31972-43-7	50MG	MET-13201A-100MG	2-(p-Chlorophenyl)-3-methyl butyric acid	2012-74-0	100MG
MET-11947B-100MG	Fenamiphos sulfone	31972-44-8	100MG	MET-13210A-50MG	Sethoxydim-desethoxy	113561-03-8	50MG
MET-11948-100MG	2,4'-Dichlorobenzophenone	85-29-0	100MG	MET-13210B-50MG	Sethoxydim-oxazole	413594-06-6	50MG
MET-11948A-100MG	2,4-Dichlorobenzhydrol	43171-49-9	100MG	MET-13213A-50MG	Atrazine desisopropyl	1007-28-9	50MG
MET-11952A-100MG	Fenchlorphos-oxon</						

Part #	Description	CAS #	Size	Part #	Description	CAS #	Size
MET-13511-50MG	Terbumeton-desethyl	30125-64-5	50MG	N-10876-100MG	4,4'-DDT	50-29-3	100MG
MET-13512A-50MG	Terbutylazine-desethyl	30125-63-4	50MG	N-10888-1G	4,6-Dinitro-2-sec-butylphenol ethanolamine salt	534-52-1	250MG
MET-13561-100MG	Thiacloprid-amide	676228-91-4	100MG	N-10889-250MG	4,6-Dinitro-o-cresol (contains ~10% water)	1214-39-7	250MG
MET-13562A-50MG	N-Desmethylthiamethoxam	171103-04-1	50MG	N-10956-250MG	6-Benzylaminopurine	148-24-3	250MG
MET-13562AA1-1ML	N-Desmethylthiamethoxam Solution	171103-04-1	1ML	N-10965-250MG	8-Hydroxyquinoline	84-65-1	250MG
MET-13576A-250MG	1,2,4-Triazole	288-88-0	250MG	N-10970-250MG	9,10-Anthraquinone	15879-93-3	250MG
MET-13576B-100MG	1-(2,4-Dichlorophenyl)-2-(1H-1,2,4-triazole-1-yl)-ethanone	58905-16-1	100MG	N-10978-250MG	a-Chloralose(Technical)	959-98-8	100MG
MET-13584AM1-1ML	DMST Solution	66840-71-9	1ML	N-10979-100MG	a-Endosulfan	86-88-4	250MG
MET-13666A-100MG	Trietazine-desethyl	38902-68-0	100MG	N-10983-250MG	a-Naphthylthiourea	115-29-7	500MG
MET-13666AA1-1ML	Trietazine-desethyl Solution	38902-68-0	1ML	N-10992-500MG	a,b-Endosulfan	71751-41-2	100MG
MET-13666B-10MG	Trietazine-2-hydroxy	13532-25-7	10MG	N-10995-100MG	Abamectin	3383-96-8	100MG
MET-13666BA1-1ML	Trietazine-2-hydroxy Solution	13532-25-7	1ML	N-10996-100MG	Temephos	18181-80-1	250MG
MET-13706A-100MG	Trinezapac	104273-73-6	100MG	N-10997-250MG	Bromopropylate	111991-09-4	100MG
MET-13750A-250MG	7-Hydroxycoumarin	93-35-6	250MG	N-10998-100MG	Nicosulfuron	30560-19-1	250MG
MET-13750B-250MG	4-Hydroxycoumarin	1076-38-6	250MG	N-11002-250MG	Acephate	57960-19-7	250MG
MET-13760A-50MG	Thionazin-O-analog	7359-55-9	50MG	N-11003-250MG	Acequinocyl	135410-20-7	100MG
MET-13760BM1-1ML	2-Hydroxypyrazine Solution	6270-63-9	1ML	N-11009-100MG	Acetamidiprid	34256-82-1	100MG
MET-13800B-25MG	Simazine hydroxy	2599-11-300	25MG	N-11013-100MG	Acetochlor	131558-54-2	250MG
MET-13801A-50MG	Pirimicarb-desmethyl	30614-22-3	50MG	N-11026-250MG	Acibenzolar-S-methyl	50594-66-6	250MG
MET-13802-100MG	Fenoxaprop-P	113158-40-0	100MG	N-11027-250MG	Acifluorfen	50594-67-7	100MG
MET-13932A-10MG	3-Hydroxy-5-phenylisoxazole	939-05-9	10MG	N-11028-100MG	Acifluorfen methyl ester	16079-88-2	500MG
MET-13987A-50MG	2-Methoxyethyl 2-[4-(tert-butyl)phenyl]cyanoacetate	400882-57-7	50MG	N-11041-500MG	Agribrom	1757-18-2	100MG
MET-14013A-100MG	2-Methyl-3-furic acid	6947-94-0	100MG	N-11042-100MG	Akton (TM)	15972-60-8	250MG
MET-14030A-100MG	2-Hydrazino-4,6-dimethylpyrimidine	23906-13-0	100MG	N-11043-250MG	Alachlor	116-06-3	100MG
MET-14030AM1-1ML	2-Hydrazino-4,6-dimethylpyrimidine Solution	23906-13-0	1ML	N-11044-100MG	Aldicarb	1646-88-4	10MG
MET-14035A-10MG	2-(1,3-dimethylbutyl)-Benzenamine	203448-76-4	10MG	N-11045-10MG	Aldicarb sulfone	1646-87-3	10MG
MET-14038A-100MG	Cyhalofop-diacid	252564-94-6	100MG	N-11046-10MG	Aldicarb sulfoxide	309-00-2	250MG
N-10001-1G	(2-Bromoethyl)trimethyl ammonium bromide	2758-06-700	1G	N-11049-250MG	Aldrin	584-79-2	100MG
N-10002-250MG	(2-Chloroethyl)phosphonic acid	16672-87-0	250MG	N-11051-100MG	Allethrin	66003-55-2	50MG
N-10004-1G	2,4,5-T butoxypropylene glycol ester (Technical)	1928-37-6	1G	N-11052-50MG	Alloxydim-Sodium	3004-51-1	250MG
N-10005-250MG	(2,4,5-Trichlorophenoxy)acetic acid methyl ester	81-84-5	250MG	N-11059-250MG	Allyl trimethylammonium bromide	67375-30-8	100MG
N-10012-500MG	1-8-Naphthalic anhydride	51229-78-8	100MG	N-11061-100MG	alpha-Cypermethrin	67485-29-4	100MG
N-10045-100MG	1-cis-3-Chloroallyl-3,5,7-triazol-1-azonia-adamantane chlorid	86-86-2	1G	N-11063-100MG	Amdro	2032-59-9	250MG
N-10074-1G	1-Naphthaleneacetamide	86-87-3	1G	N-11065-250MG	Aminocarb	150114-71-9	100MG
N-10075-1G	1-Naphthaleneacetic acid	117-18-0	250MG	N-11066-100MG	Aminopyralid	81-82-5	250MG
N-10111-250MG	1,2,4,5-Tetrachloro-3-nitrobenzene	594-72-9	1G	N-11067-250MG	Aminotriazole	33089-61-1	250MG
N-10122-1G	1,1-Dichloro-1-nitroethane	96-12-8	250MG	N-11068-250MG	Amitraz	16919-19-0	500MG
N-10149-250MG	1,2-Dibromo-3-chloropropane	106-93-4	1G	N-11069-500MG	Ammonium fluosilicate(Technical)	7773-06-0	500MG
N-10150-1G	1,2-Dibromoethane	95-50-1	1G	N-11070-500MG	Ammonium sulfate(Technical)	3566-10-700	100MG
N-10152-1G	1,2-Dichlorobenzene	78-87-5	1G	N-11071-100MG	Amobam	12771-68-5	100MG
N-10156-1G	1,2-Dichloropropane	542-75-6	500MG	N-11073-100MG	Ancymidol	101-05-3	250MG
N-10193-500MG	1,3-Dichloropropene	106-46-7	1G	N-11075-250MG	Anilazine	140-57-8	10MG
N-10216-1G	1,4-Dichlorobenzene	2597-93-5	100MG	N-11085-10MG	Aramite	66230-04-4	100MG
N-10224-100MG	Emmi (Technical)	101-10-0	100MG	N-11102-100MG	Esfenvalerate	3244-90-4	10MG
N-10254-100MG	2-(3-Chlorophenoxy)propionic acid	112-56-1	250MG	N-11103-10MG	Aspon (TM)	3337-71-1	250MG
N-10260-250MG	[2-Butoxyethoxy]ethyl thiocyanate	5825-87-6	1G	N-11104-250MG	Asulam	1610-17-9	250MG
N-10261-1G	2-[3-Chlorophenoxy]propionamide	90-43-7	250MG	N-11105-250MG	Atraton	1912-24-9	250MG
N-10280-250MG	2-Biphenylol	115-84-4	1G	N-11106-250MG	Atrazine	11141-17-6	10MG
N-10306-1G	2-Butyl-2-ethyl-1,3-propanediol	94-96-2	500MG	N-11107-10MG	Azadirachtin(Technical)	103-33-3	500MG
N-10346-500MG	2-Ethyl-1,3-hexanediol	1928-43-4	1G	N-11110-500MG	Azobenzene	41083-11-8	250MG
N-10356-1G	2-Ethylhexyl-2,4-dichlorophenoxy acetate	640-19-7	100MG	N-11111-250MG	Azocyclotin (Technical)	131860-33-8	100MG
N-10358-100MG	2-Fluoroacetamide	539-47-9	100MG	N-11113-100MG	Azoxystrobin	16376-36-6	100MG
N-10363-100MG	2-Furanacrylic acid	3547-33-9	100MG	N-11117-100MG	b-Chloralose(Technical)	32123-65-9	100MG
N-10378-100MG	2-Hydroxyethyl-n-octyl sulphide	96-45-7	250MG	N-11118-100MG	b-Endosulfan	671-04-5	100MG
N-10381-250MG	2-Imidazolidinethione	83-28-3	250MG	N-11122-100MG	Carbanolate	101-27-9	100MG
N-10384-250MG	2-Isovaleryl-1,3-indandione	120-23-0	250MG	N-11123-100MG	Barban	3766-81-2	250MG
N-10437-250MG	2-Naphthoxyacetic acid	940-31-8	500MG	N-11126-250MG	Baycarb (TM)	55179-31-2	100MG
N-10507-500MG	2-Phenoxypropionic acid	83-26-1	250MG	N-11127-100MG	Bifentanol	114-26-1	250MG
N-10512-250MG	2-Pivaloyl-1,3-indandione	6420-47-9	500MG	N-11128-250MG	Propoxur	55219-65-3	250MG
N-10517-500MG	2-sec-Butyl-4,6-dinitrophenol triethanolamine salt	112-12-9	1G	N-11129-250MG	Triadenol	68359-37-5	250MG
N-10520-1G	2-Undecanone	3481-20-7	250MG	N-11130-250MG	Cyfluthrin	672-99-1	100MG
N-10530-250MG	2,3,5,6-Tetrachloroaniline	1713-15-1	250MG	N-11131-100MG	BDMC	41814-78-2	50MG
N-10532-250MG	2,4-D isobutyl ester	18625-12-2	100MG	N-11132-50MG	Tricyclazole	71626-11-4	250MG
N-10533-100MG	2,4-DB methyl ester	24157-81-1	250MG	N-11133-250MG	Benalaxyl	25059-80-7	250MG
N-10550-250MG	2,6-Diisopropyl-naphthalene	3209-22-1	250MG	N-11134-250MG	Benazolin-ethyl	22781-23-3	250MG
N-10584-250MG	2,3-Dichloronitrobenzene	58-90-2	50MG	N-11135-250MG	Bendiocarb	1861-40-1	250MG
N-10599-50MG	2,3,4,6-Tetrachlorophenol	88-82-4	250MG	N-11136-250MG	Benfluralin	82560-54-1	100MG
N-10601-250MG	2,3,5-Triiodobenzoic acid	2655-15-4	50MG	N-11137-100MG	Benfuracarb	17804-35-2	100MG
N-10603-50MG	2,3,5-Trimethylphenyl methyl carbamate	50-31-7	100MG	N-11138-100MG	Benomyl	83055-99-6	100MG
N-10605-100MG	2,3,6-Trichlorobenzoic acid	94-75-7	250MG	N-11139-100MG	Bensulfuron-methyl	741-58-2	250MG
N-10609-250MG	2,4-D	1929-73-3	250MG	N-11140-250MG	Bensulide	17606-31-4	250MG
N-10610-250MG	2,4-D butoxyethyl ester	94-80-4	250MG	N-11141-250MG	Bensultap	25057-89-0	250MG
N-10611-250MG	2,4-D butyl ester	2008-39-1	250MG	N-11142-250MG	Bentazon	61592-45-8	50MG
N-10612-250MG	2,4-D dimethylamine salt	3599-58-4	250MG	N-11143-50MG	Bentazon methyl derivative	28249-77-6	250MG
N-10613-250MG	2,4-D ethanalamine salt	533-23-3	250MG	N-11144-250MG	Benthicarb	22212-55-1	100MG
N-10614-250MG	2,4-D ethyl ester	25168-26-7	250MG	N-11177-100MG	Benzoilprop ethyl	120-51-4	100MG
N-10615-250MG	2,4-D isooctyl ester	6365-72-6	250MG	N-11182-100MG	Benzy benzoate	68359-37-5	250MG
N-10616-250MG	2,4-D isopropanol amine salt	94-11-1	250MG	N-11191-250MG	beta-Cyfluthrin	6731-84-2	250MG
N-10617-250MG	2,4-D isopropyl ester	5742-17-6	250MG	N-11192-250MG	beta-Cypermethrin	319-84-6	100MG
N-10618-250MG	2,4-D isopropylamine salt	1928-38-7	250MG	N-11194-100MG	BHC (alpha isomer)	319-85-7	100MG
N-10619-250MG	2,4-D methyl ester	103564-78-9	250MG	N-11195-100MG	BHC (beta isomer)	319-86-8	50MG
N-10620-250MG	2,4-D propylene glycol ester	94-79-1	250MG	N-11196-50MG	BHC (delta isomer)	608-73-1	250MG
N-10621-250MG	2,4-D sec-butyl ester	94-82-6	250MG	N-11198-250MG	BHC (mixed isomers)	149877-41-8	50MG
N-10622-250MG	2,4-DB	2758-42-1	250MG	N-11201-50MG	Bifenox	42576-02-3	250MG
N-10623-250MG	2,4-DB dimethylamine salt	611-06-3	250MG	N-11202-250MG	Bifenthrin	82657-04-3	100MG
N-10625-250MG	2,4-Dichloronitrobenzene	97-16-5	250MG	N-11203-100MG	Binapacryl	485-31-4	100MG
N-10628-250MG	2,4-Dichlorophenyl benzenesulfonate	93-76-5	250MG	N-11204-100MG	Bioallethrin	28057-48-9	250MG
N-10646-250MG	2,4,5-T	1928-47-8	250MG	N-11205-250MG	Bioresmethrin	28434-01-7	250MG
N-10647-250MG	2,4,5-T 2-ethylhexyl ester	2545-59-7	100MG	N-11206-250MG	Biphenyl	92-52-4	500MG
N-10648-100MG	2,4,5-T butoxyethyl ester	25168-15-4	100MG	N-11207-500MG	Bis[tri-n-butyltin]oxide(Technical)	125401-92-5	100MG
N-10649-100MG	2,4,5-T isooctyl ester	93-78-7	100MG	N-11322-250MG	Bispyribac Sodium	3734-33-6	100MG
N-10650-100MG	2,4,5-T isopropyl ester	93-79-8	250MG	N-11323-100MG	Bitrex (TM)	21725-46-2	250MG
N-10651-250MG	2,4,5-T n-butyl ester	57213-69-1	250MG	N-11324-100MG	Boscalid	188425-85-6	100MG
N-10652-250MG	2,4,5-T triethylamine salt	88-06-2	250MG	N-11327-250MG	Brodifacoum	56073-10-0	50MG
N-10655-250MG	2,4,5-Trichlorophenol ethanolamine salt	89-61-2	250MG	N-11328-100MG	Bromacil	314-40-9	250MG
N-10657-250MG	2,4,6-Trichlorophenol	1194-65-6	250MG	N-11329-50MG	Bromoacetyl bromide	598-21-0	500MG
N-10670-250MG	2,6-Dichlorobenzonitrile	87-51-4	250MG	N-11330-250MG	Bromophos ethyl	4824-78-6	100MG
N-10688-250MG	3-Indoleacetic acid	133-32-4	100MG	N-11333-500MG	Bromophos methyl	2104-96-3	50MG
N-10727-250MG	3-Indolebutyric acid	830-96-6	250MG	N-11342-100MG	Bromoxynil	1689-84-5	250MG
N-10728-100MG	3-Indolepropionic acid	2686-99-9	10MG	N-11343-50MG	Bromoxynil octanoic acid ester	1689-99-2	250MG
N-10729-250MG	3,4,5-Trimethylphenyl methyl carbamate	51-36-5	250MG	N-11345-250MG	Bromuconazole	116255-48-2	100MG
N-10769-10MG	3,5-Dichlorobenzoic acid	99-54-7	250MG	N-11347-100MG	Bronopol	52-51-7	250MG
N-10770-250MG	3,5-Dichloronitrobenzene	2905-67-1	100MG	N-11348-250MG	Bupirimate	41483-43-6	100MG
N-10785-250MG	3,5-Dichloronitrobenzene methyl ester	618-62-2	250MG	N-11349-100MG	Buprofezin	69327-76-0	100MG
N-10791-100MG	3,5-Dichloronitrobenzene	1689-83-4	100MG	N-11350-100MG	Busan 40	51026-28-9	100MG
N-10792-250MG	3,5-Diiodo-4-hydroxybenzonitrile	709-98-8	500MG	N-11351-100MG	Busan 85	128-03-0	100MG
N-10794-100MG	Propanil	97-53-0	250MG	N-11352-100MG	Butachlor	23184-66-9	100MG
N-10800-500MG	4-Allyl-2-methoxyphenol	504-24-5	1G	N-11353-100MG	Butoxypolpropylene glycol(Technical)	9003-13-8	250MG
N-10803-250MG	4-Aminopyridine	94-74-6	250MG	N-11354-250MG	Butralin	33629-47-9	500MG
N-10806-1G	4-Chloro-o-tolxyloxyacetic acid	2436-73-9	100MG	N-11355-500MG	Butylate	2008-41-5	250MG
N-10818-250MG	4-Chloro-o-tolxyloxyacetic acid methyl ester	72-54-8	250MG	N-11356-500MG	Butylmesityl oxide oxalate	532-34-3	250MG
N-10819-100MG	4,4'-DDD	72-55-9	100MG	N-11383-250MG	Cadmium succinate(Technical)	141-00-4	250MG
N-10874-250MG	4,4'-DDE			N-11391-250MG			

Part #	Description	CAS #	Size	Part #	Description	CAS #	Size
N-11392-100MG	Cadusafos	95465-99-9	100MG	N-11678-250MG	Dicloran	99-30-9	250MG
N-11398-100MG	Capsaicin, Natural	404-86-4	100MG	N-11679-250MG	Diclosulam	145701-21-9	250MG
N-11399-250MG	Captafol	2425-06-100	250MG	N-11680-100MG	Dicrotophos	141-66-2	100MG
N-11400-250MG	Captan	133-06-2	250MG	N-11686-500MG	Dicyclopentadiene	77-73-6	500MG
N-11402-250MG	Carbaryl	63-25-2	250MG	N-11688-100MG	Dieldrin	60-57-1	100MG
N-11404-100MG	Carbendazim	10605-21-7	100MG	N-11690-100MG	Diethatyl ethyl	38727-55-8	100MG
N-11405-250MG	Carbofuran	1563-66-2	250MG	N-11691-250MG	Diethofencarb	87130-20-9	250MG
N-11408-250MG	Carbophenothion	786-19-6	250MG	N-11696-250MG	Diethyl dithiobis(thionoformate)	502-55-6	250MG
N-11409-250MG	Carbosulfan	55285-14-8	250MG	N-11719-100MG	Difenacoum	56073-07-5	100MG
N-11410-250MG	Carboxin	5234-68-4	250MG	N-11720-250MG	Difenoconazole	119446-68-3	250MG
N-11411-10MG	Carfentrazone-ethyl	128639-02-1	10MG	N-11721-250MG	Difenzoquat methylsulfate	43222-48-6	250MG
N-11412-250MG	Cartap hydrochloride	15263-52-2	250MG	N-11722-250MG	Diffubenzuron	35367-38-5	250MG
N-11418-100MG	Chloramben	133-90-4	100MG	N-11723-100MG	Diffufenican	83164-33-4	100MG
N-11419-100MG	Chloramben methyl ester	7286-84-2	100MG	N-11724-100MG	Diffufenzopyr	109293-97-2	100MG
N-11420-250MG	Chloranil	118-75-2	250MG	N-11750-250MG	Dikegulac acid	18467-77-1	250MG
N-11421-250MG	Chloranocryl		250MG	N-11752-250MG	Dimehypo	52207-48-4	250MG
N-11422-100MG	Chlorantraniliprole	500008-45-7	100MG	N-11753-500MG	Dimethachlon	24096-53-5	500MG
N-11423-100MG	Chlorbenside	103-17-3	100MG	N-11754-250MG	Dimethachlor	50563-36-5	250MG
N-11424-250MG	Chlorbromuron	13360-45-7	250MG	N-11755-100MG	Dimethenamid	87674-68-8	100MG
N-11425-250MG	Chlordane	57-74-9	250MG	N-11756-100MG	Dimethenamid-P	163515-14-8	100MG
N-11426-250MG	Chlordimeform	6164-98-3	250MG	N-11757-100MG	Dimethipin	55290-64-7	100MG
N-11427-250MG	Chlorethoxyfos	54593-83-8	250MG	N-11758-100MG	Dimethoate	60-51-5	100MG
N-11428-100MG	Chlorfenapyr	122453-73-0	100MG	N-11759-100MG	Dimethomorph	110488-70-5	100MG
N-11429-250MG	Chlorfenvinphos	470-90-6	250MG	N-11770-500MG	Dimethyl phthalate	131-11-3	500MG
N-11430-50MG	Chlorflazuron	71422-67-8	50MG	N-11775-100MG	Dimethyl-p-nitrophenylphosphate	950-35-6	100MG
N-11431-100MG	Chlorflurecol-methyl ester	2536-31-4	100MG	N-11779-500MG	Dimethylarsinic acid	75-60-5	500MG
N-11432-100MG	Chlorimuron ethyl	90982-32-4	100MG	N-11782-100MG	Diniconazole	83657-24-3	100MG
N-11434-250MG	Chlormephos	24934-91-6	250MG	N-11783-100MG	Dinocap	39300-45-3	100MG
N-11435-250MG	Chloromequat chloride	999-81-5	250MG	N-11786-100MG	Dinoseb	88-85-7	100MG
N-11438-250MG	Chloroacetic acid	79-11-8	250MG	N-11787-100MG	Dinoseb methyl ether	6099-79-2	100MG
N-11442-100MG	Chlorobenzilate	510-15-6	100MG	N-11788-50MG	Dinotefuran	165252-70-0	50MG
N-11449-250MG	Chloroneb	2675-77-6	250MG	N-11789-250MG	Dioxacar	6986-21-2	250MG
N-11451-250MG	Chlorophacinone	3691-35-8	250MG	N-11790-100MG	Dioxathion	78-34-2	100MG
N-11452-100MG	Trichloronitromethane	76-06-2	100MG	N-11791-100MG	Dioxydemeton-S-methyl	17040-19-6	100MG
N-11453-500MG	Chloropropylate	5836-10-200	500MG	N-11793-100MG	Diphacinone	82-66-6	100MG
N-11454-250MG	Chlorothalonil	1897-45-6	250MG	N-11794-250MG	Diphenamid	957-51-7	250MG
N-11455-250MG	Chlorotoluron	15545-48-9	250MG	N-11800-100MG	Diphenylacetoneitrile	86-29-3	100MG
N-11456-250MG	Chloroxuron	1982-47-4	250MG	N-11808-100MG	Dipropetryn	4147-51-7	100MG
N-11457-500MG	Chloroxylin	1891-95-8	500MG	N-11816-500MG	Diquat Dibromide monohydrate	6385-62-2	500MG
N-11458-250MG	Chlorpropham	101-21-3	250MG	N-11817-500MG	Sodium methyl arsonate hexahydrate	144-21-8	500MG
N-11459-250MG	Chlorpyrifos	2921-88-2	250MG	N-11819-250MG	Disulfoton	298-04-4	250MG
N-11460-250MG	Chlorpyrifos Methyl	5598-13-0	250MG	N-11820-100MG	Ditalimfos	5131-24-8	100MG
N-11461-100MG	Chlorsulfuron	64902-72-3	100MG	N-11821-250MG	Dithianon	3347-22-6	250MG
N-11462-250MG	Chlorthal-dimethyl	1861-32-1	250MG	N-11822-100MG	Dithiopyr	97886-45-8	100MG
N-11463-100MG	Cholecalciferol	67-97-0	100MG	N-11827-250MG	Diuron	330-54-1	250MG
N-11468-500MG	Cinnamaldehyde	104-55-2	500MG	N-11839-250MG	Dodemorph	1593-77-7	250MG
N-11470-100MG	Cinosulfuron	94593-91-6	100MG	N-11840-250MG	Dodine (Technical)	2439-10-300	250MG
N-11480-100MG	cis-Chlordane	5103-71-9	10MG	N-11842-100MG	Fonofos	944-22-9	100MG
N-11481-100MG	E-Mevinphos	26718-65-0	100MG	N-11843-250MG	Trichlorphon	52-68-6	250MG
N-11483-50MG	cis-Permethrin	61949-76-6	50MG	N-11844-100MG	Edifenphos	17109-49-8	100MG
N-11488-100MG	Clethodim	99129-21-2	100MG	N-11845-100MG	EGT	2514-53-6	100MG
N-11489-100MG	Clodinafop-propargyl	105512-06-9	100MG	N-11847-100MG	Emamectin	119791-41-2	100MG
N-11490-100MG	Clofentezine	74115-24-5	100MG	N-11851-500MG	Endosulfan sulfate	1031-07-8	500MG
N-11491-100MG	Cloquintocet-mexyl	99607-70-2	100MG	N-11852-100MG	Endothal monohydrate	62059-43-2	100MG
N-11492-100MG	Cloransulam-methyl	147150-35-4	100MG	N-11854-100MG	Endrin	72-20-8	100MG
N-11493-100MG	Clothianidin	210880-92-5	100MG	N-11855-10MG	Endrin aldehyde	7421-93-4	10MG
N-11495-100MG	Command	81777-89-1	100MG	N-11856-10MG	Endrin ketone	53494-70-5	10MG
N-11501-250MG	Copper Naphthenates(Technical)	1338-02-9	250MG	N-11858-250MG	EPN	2104-64-5	250MG
N-11502-250MG	Copper oxychloride(Technical)	1332-65-6	250MG	N-11859-100MG	Epoxiconazole	133855-98-8	100MG
N-11503-500MG	Copper-8-Quinolinolate	10380-28-6	500MG	N-11864-100MG	Esbiothrin	84030-86-4	100MG
N-11507-100MG	Coumaphos	56-72-4	100MG	N-11865-250MG	Ethalfuralin	55283-68-6	250MG
N-11509-250MG	Coumatetralyl	5836-29-3	250MG	N-11866-100MG	Ethametsulfuron-methyl	97780-06-8	100MG
N-11513-50MG	Crotoxyphos	7700-17-6	50MG	N-11867-100MG	Ethanedial dioxime	557-30-2	100MG
N-11517-100MG	Cyanophos	2636-26-2	100MG	N-11867-250MG	Ethanedial dioxime	557-30-2	250MG
N-11519-100MG	Cyazofamid	120116-88-3	100MG	N-11871-100MG	Ethidimuron	30043-49-3	100MG
N-11520-100MG	Cybutryne	28159-98-0	100MG	N-11872-100MG	Ethiofencarb	29973-13-5	100MG
N-11521-100MG	Cyclanilide	113136-77-9	100MG	N-11873-250MG	Ethion	563-12-2	250MG
N-11522-250MG	Cycloate	1134-23-2	250MG	N-11874-250MG	Ethofumesate	26225-79-6	250MG
N-11534-100MG	Cycloheximide	66-81-9	100MG	N-11877-100MG	Ethoxyquin	91-53-2	100MG
N-11541-500MG	Cycloprate	54460-46-7	500MG	N-11878-100MG	Ethoxysulfuron	126801-58-9	100MG
N-11542-250MG	Cyclosulfamuron	136849-15-5	250MG	N-11914-250MG	Ethyl-1-naphthalene acetate	2122-70-5	250MG
N-11543-500MG	Cyhalofop-butyl	122008-85-9	500MG	N-11916-100MG	Ethyl-3-amino-2,5-dichlorobenzoate	60541-86-8	100MG
N-11544-100MG	Cymoxanil	57966-95-7	100MG	N-11940-100MG	Ethylmercuric chloride(Technical)	107-27-7	100MG
N-11545-100MG	Cypermethrin	52315-07-8	100MG	N-11941-100MG	Etofenprox	80844-07-1	100MG
N-11546-500MG	Cyprazine	22936-86-3	500MG	N-11942-500MG	Etoxadone	153233-91-1	500MG
N-11547-500MG	Cyproconazole	94361-06-5	500MG	N-11943-100MG	Famoxadone	131807-57-3	100MG
N-11548-250MG	Cyprodinil	121552-61-2	250MG	N-11944-250MG	Famphur	52-85-7	250MG
N-11549-100MG	Cyprosulfamide	221667-31-8	100MG	N-11945-250MG	Fenac	85-34-7	250MG
N-11550-250MG	Cyromazine	66215-27-8	250MG	N-11946-100MG	Fenamidone	161326-34-7	100MG
N-11551-500MG	Cythioate	115-93-5	500MG	N-11948-100MG	Fenarimol	60168-88-9	100MG
N-11554-500MG	d-(cis-trans)-Phenothrin	26002-80-2	500MG	N-11949-100MG	Fenazaquin	120928-09-8	100MG
N-11561-100MG	d-trans-Cyphenothrin	39515-40-7	100MG	N-11950-100MG	Fenbuconazole	114369-43-6	100MG
N-11562-250MG	Dalapon	75-99-0	250MG	N-11951-250MG	Fenbutatin-oxide(Technical)	13536-08-6	250MG
N-11563-100MG	Dalapon methyl ester	17640-02-7	100MG	N-11952-100MG	Fenchlorphos	299-84-3	100MG
N-11564-250MG	Dalapon Sodium	127-20-8	250MG	N-11954-250MG	Fenhexamid	126833-17-8	250MG
N-11565-250MG	Daminozide	1596-84-5	250MG	N-11955-250MG	Fenitrothion	122-14-5	250MG
N-11566-250MG	Dazomet	533-74-4	250MG	N-11956-100MG	Fenoxanil	115852-48-7	100MG
N-11567-250MG	DDT (mixture p,p' & o,p')	8017-34-3	250MG	N-11957-100MG	Fenoxaprop ethyl	66441-23-4	100MG
N-11579-250MG	Deltamethrin	52918-63-5	250MG	N-11958-100MG	Fenoxaprop-P-ethyl	71283-80-2	100MG
N-11580-100MG	Demeton O	298-03-3	100MG	N-11959-100MG	Fenoxycarb	72490-01-8	100MG
N-11581-100MG	Demeton O&S	8065-48-3	100MG	N-11960-250MG	Fenpropathrin	39515-41-8	250MG
N-11582-100MG	Demeton S	126-75-0	100MG	N-11962-100MG	Fenpyroximate	134098-61-6	100MG
N-11583-500MG	Demeton-O-methyl	867-27-6	500MG	N-11963-100MG	Fensulfthion	115-90-2	100MG
N-11584-250MG	Desmedipham	13684-56-5	250MG	N-11964-250MG	Fenthion	55-38-9	250MG
N-11585-250MG	Napropamide	15299-99-7	250MG	N-11965-500MG	Fenthion-ethyl	1716-09-2	500MG
N-11586-250MG	Fenaminosulf	140-56-7	250MG	N-11966-250MG	Fentin hydroxide	76-87-9	250MG
N-11587-100MG	Di-allyl	2303-16-4	100MG	N-11967-250MG	Fenuron	101-42-8	250MG
N-11612-250MG	Diafenthiuron	80060-09-9	250MG	N-11968-50MG	Fenuron TCA	4482-55-7	50MG
N-11621-250MG	Diazinon	333-41-5	250MG	N-11970-250MG	Ferbam	14484-64-1	250MG
N-11644-500MG	Dibutyl succinate	141-03-7	500MG	N-11971-100MG	Fipronil	120068-37-3	100MG
N-11656-250MG	Dicamba	1918-00-9	250MG	N-11972-100MG	Flamprop isopropyl	52756-22-6	100MG
N-11657-100MG	Dicamba methyl ester	6597-78-0	100MG	N-11973-100MG	Flamprop-methyl	52756-25-9	100MG
N-11658-100MG	Dicapthion	2463-84-5	100MG	N-11974-500MG	Flazasulfuron	104040-78-0	500MG
N-11659-250MG	Dichlofenthiol	97-17-6	250MG	N-11975-10MG	Floucomafen	90035-08-8	10MG
N-11660-250MG	Dichlofluanid	1085-98-9	250MG	N-11976-50MG	Flonicamid	158062-67-0	500MG
N-11661-250MG	Dichlone	117-80-6	250MG	N-11977-10MG	Fluazifop	69335-91-7	100MG
N-11662-100MG	Dichlormid	37764-25-3	100MG	N-11978-250MG	Fluazifop-butyl	69806-50-4	250MG
N-11667-250MG	Dichlorophen	97-23-4	250MG	N-11979-100MG	Fluazifop-P-butyl	79241-46-6	100MG
N-11671-250MG	Dichlorprop	120-36-5	250MG	N-11980-100MG	Fluazinam	79622-59-6	100MG
N-11672-100MG	Dichlorprop methyl ester	57153-17-0	100MG	N-11981-100MG	Flubendiamide	272451-65-7	100MG
N-11673-500MG	Dichlorprop-butoxyethyl ester	53404-31-2	500MG	N-11982-250MG	Flucarbazone-sodium	181274-17-9	250MG
N-11674-250MG	Dichlorprop-P	15165-67-0	250MG	N-11983-50MG	Fluchloralin	33245-39-5	50MG
N-11675-250MG	Dichlorvos	62-73-7	250MG	N-11984-100MG	Flucythrinate	70124-77-5	100MG
N-11677-250MG	Diclofop methyl	51338-27-3	250MG	N-11985-100MG	Fludioxonil	131341-86-1	100MG

Part #	Description	CAS #	Size	Part #	Description	CAS #	Size
N-11986-100MG	Flufenacet	142459-58-3	100MG	N-12372-100MG	Mepronil	55814-41-0	100MG
N-11987-50MG	Flufenoxuron	101463-69-8	50MG	N-12376-100MG	Mesosulfuron-methyl	208465-21-8	100MG
N-11988-250MG	Flumetralin	62942-70-3	250MG	N-12377-100MG	Mesotrione	104206-82-8	100MG
N-11989-100MG	Flumetsulam	98967-40-9	100MG	N-12379-100MG	Metaflumizone	139968-49-3	100MG
N-11990-500MG	Flumiclorac-pentyl	87546-18-7	500MG	N-12380-100MG	Metalaxyl	57837-19-1	100MG
N-11991-100MG	Flumioxazin	103361-09-7	100MG	N-12381-1G	Metaldehyde(Technical)	9002-91-9	1G
N-11992-250MG	Fluometuron	2164-17-2	250MG	N-12382-250MG	Metam sodium	6734-80-1	250MG
N-11993-100MG	Fluopicolide	239110-15-7	100MG	N-12383-250MG	Metamitron	41394-05-2	250MG
N-11998-250MG	Fluorodifen	15457-05-3	250MG	N-12384-50MG	Metasystox I (TM)	919-86-8	50MG
N-11999-100MG	Fluoroglycofen-ethyl	77501-90-7	100MG	N-12385-100MG	Metazachlor	67129-08-2	100MG
N-12000-100MG	Fluroxypyr	69377-81-7	100MG	N-12386-100MG	Metconazole	125116-23-6	100MG
N-12001-100MG	Fluroxypyr-meptyl	81406-37-3	100MG	N-12387-100MG	Methabenzthiazuron	18691-97-9	100MG
N-12002-100MG	Flurprimidol	56425-91-3	100MG	N-12388-100MG	Methacrifos	30864-28-9	100MG
N-12003-100MG	Flusilazole	85509-19-9	100MG	N-12393-100MG	Methamidophos	10265-92-6	100MG
N-12004-250MG	Flutolanil	66332-96-5	250MG	N-12397-250MG	Methidathion	950-37-8	250MG
N-12005-100MG	Flutriafol	76674-21-0	100MG	N-12398-100MG	Methiocarb	2032-65-7	100MG
N-12006-250MG	Tebuconazole	107534-96-3	250MG	N-12399-100MG	Methomyl	16752-77-5	100MG
N-12007-250MG	Folpet	133-07-3	250MG	N-12400-100MG	Methoprene	40596-69-8	100MG
N-12008-500MG	Fomesafen	72178-02-0	500MG	N-12401-100MG	Methoprotryne	841-06-5	100MG
N-12009-100MG	Foramsulfuron	173159-57-4	100MG	N-12403-100MG	Methoxychlor	72-43-5	100MG
N-12010-100MG	Forchlorfenuron	68157-60-8	100MG	N-12406-100MG	Methoxyfenozide	161050-58-4	100MG
N-12012-1G	Formaldehyde (in water)	50-00-0	1G	N-12414-100MG	Methyl anthranilate	134-20-3	100MG
N-12015-250MG	Formetanate hydrochloride	23422-53-9	250MG	N-12452-100MG	Methyl parathion	298-00-0	100MG
N-12017-100MG	Formothion	2540-82-1	100MG	N-12464-500MG	Methyl-1-naphthalene acetate	2876-78-0	500MG
N-12019-100MG	Fosetyl-Aluminum(Technical)	39148-24-8	100MG	N-12477-250MG	Metobromuron	3060-89-7	250MG
N-12020-30MG	Fosfiazate	98886-44-3	30MG	N-12478-100MG	Metolachlor	51218-45-2	100MG
N-12021-100MG	Fuberidazole	3878-19-1	100MG	N-12479-100MG	Metolcarb	1129-41-5	100MG
N-12106-100MG	gamma-Cyhalothrin	76703-62-3	100MG	N-12480-250MG	Metoxuron	19937-59-8	250MG
N-12108-100MG	Gibberellic acid	77-06-5	100MG	N-12481-250MG	Metribuzin	21087-64-9	250MG
N-12111-100MG	Glufosinate-ammonium	77182-82-2	100MG	N-12482-100MG	Metsulfuron methyl	74223-64-6	100MG
N-12131-250MG	Glyodin(Technical)	556-22-9	250MG	N-12483-250MG	MGK 264 (TM)	113-48-4	250MG
N-12133-250MG	Glyphosate	1071-83-6	250MG	N-12484-250MG	Di-n-propyl isocinchomeronate	136-45-8	250MG
N-12134-250MG	Glyphosate-isopropylammonium	38641-94-0	250MG	N-12486-100MG	Mirex	2385-85-5	100MG
N-12135-100MG	Glyphosine	2439-99-8	100MG	N-12487-250MG	Molinate	2212-67-1	250MG
N-12139-250MG	Azinphos-methyl	86-50-0	250MG	N-12488-100MG	Pencycuron	66063-05-6	100MG
N-12140-100MG	Azinphos-ethyl	2642-71-9	100MG	N-12493-100MG	Monocrotophos	6923-22-4	100MG
N-12141-500MG	Halofenozide	112226-61-6	500MG	N-12494-250MG	Monolinuron	1746-81-2	250MG
N-12142-100MG	Halosulfuron-methyl	100784-20-1	100MG	N-12495-100MG	Monosodium acid methane arsonate sesquihydrate	2163-80-6	100MG
N-12143-10MG	Haloxyfop (free acid)	69806-34-4	10MG	N-12496-50MG	Monothiono TEPP	645-78-3	50MG
N-12145-50MG	Haloxyfop-P (free acid)	95977-29-0	50MG	N-12497-250MG	Monuron	150-68-5	250MG
N-12146-50MG	Haloxyfop-P methyl	72619-32-0	50MG	N-12498-100MG	Monuron TCA	140-41-0	100MG
N-12147-100MG	Heptachlor	76-44-8	100MG	N-12505-10MG	N-[(2,4-Dichlorophenoxy)acetyl]-DL-methionine	62625-13-2	10MG
N-12148-50MG	Heptachlor epoxide (Isomer B)	1024-57-3	50MG	N-12506-500MG	N-1-Naphthylacetamide	575-36-0	500MG
N-12153-100MG	Heptenophos	23560-59-0	100MG	N-12507-250MG	N-1-Naphthylphthalamic acid	132-66-1	250MG
N-12158-500MG	Hexachloroacetone	116-16-5	500MG	N-12557-100MG	N-m-Tolylphthalamic acid	85-72-3	100MG
N-12159-250MG	Hexachlorobenzene	118-74-1	250MG	N-12558-50MG	N-Methyl-methamidophos	28167-49-9	50MG
N-12161-250MG	Hexachlorodimethyl sulfone	3064-70-8	250MG	N-12559-100MG	N-Methyl-N-1-naphthyl acetamide	5903-13-9	100MG
N-12163-250MG	Hexachlorophene	70-30-4	250MG	N-12617-250MG	N,N-Diallyl-2-chloroacetamide	93-71-0	250MG
N-12165-100MG	Hexaconazole	79983-71-4	100MG	N-12618-250MG	N,N-Diethyl-m-toluamide	134-62-3	250MG
N-12167-250MG	Hexaflumuron	86479-06-3	250MG	N-12630-250MG	N,N-Dimethylsuccinamic acid	2564-95-6	250MG
N-12184-100MG	Hexythiazox	78587-05-0	100MG	N-12639-250MG	Nabam	142-59-6	250MG
N-12196-25MG	Hymexazole	100044-41-1	25MG	N-12640-250MG	Naled	300-76-5	250MG
N-12197-100MG	Icaridin	119515-38-7	100MG	N-12647-250MG	Neburon	555-37-3	250MG
N-12198-100MG	Imazalil	35554-44-0	100MG	N-12649-250MG	Neomycin trisulfate hydrate(Technical)	1405-10-3	250MG
N-12199-100MG	Imazalil sulfate	58594-72-2	100MG	N-12651-100MG	Nicarbazin	330-95-0	100MG
N-12200-100MG	Imazamethabenz-methyl	81405-85-8	100MG	N-12653-250MG	Niclosamide	50-65-7	250MG
N-12201-250MG	Imazamox	114311-32-9	250MG	N-12655-1G	(S)-Nicotine	54-11-5	1G
N-12202-250MG	Imazapic	104098-48-8	250MG	N-12657-100MG	Nitenpyram	150824-47-8	100MG
N-12203-100MG	Imazapyr	81334-34-1	100MG	N-12658-250MG	Nitralin	4726-14-1	250MG
N-12204-100MG	Imazaquin	81335-37-7	100MG	N-12659-100MG	Nitrapyrin	1929-82-4	100MG
N-12205-100MG	Imazethapyr	81335-77-5	100MG	N-12663-100MG	Nitrofen	1836-75-5	100MG
N-12206-100MG	Imidacloprid	138261-41-3	100MG	N-12665-100MG	Nitrothal isopropyl	10552-74-6	100MG
N-12207-100MG	Phosmet	732-11-6	100MG	N-12668-100MG	Norflurazon	27314-13-2	100MG
N-12209-100MG	Indanofan	133220-30-1	100MG	N-12669-250MG	Noruron	18530-56-8	250MG
N-12213-100MG	S-Indoxacarb	173584-44-6	100MG	N-12670-100MG	Novaluron	116714-46-6	100MG
N-12218-100MG	Ioxynil octanoate	3861-47-0	100MG	N-12676-100MG	o-Benzyl-p-chlorophenol	120-32-1	100MG
N-12219-250MG	Iprobenfos	26087-47-8	250MG	N-12679-250MG	o-Chlorophenoxy acetic acid	614-61-9	250MG
N-12220-100MG	Iprodione	36734-19-7	100MG	N-12703-250MG	O,O-Diethyl phosphorochloridothioate	2524-04-100	250MG
N-12222-100MG	Isazophos	42509-80-8	100MG	N-12704-250MG	O,O-Dimethyl phosphorochloridothioate	2524-03-0	250MG
N-12231-10MG	Isobenzan	297-78-9	10MG	N-12706-250MG	o,p'-DDD	53-19-0	250MG
N-12232-250MG	Isobornyl thiocyanacetate	115-31-1	250MG	N-12707-50MG	o,p'-DDE	3424-82-6	50MG
N-12248-100MG	Isodrin	465-73-6	100MG	N-12708-50MG	o,p'-DDT	789-02-6	50MG
N-12249-100MG	Isofenphos	25311-71-1	100MG	N-12715-50MG	Octamethylpyrophosphoramide	152-16-9	50MG
N-12250-50MG	Isofenphos-Methyl	99675-03-3	50MG	N-12726-100MG	Omethoate	1113-02-6	100MG
N-12264-250MG	Isoprocarb	2631-40-5	250MG	N-12727-100MG	Propargite	2312-35-8	100MG
N-12265-250MG	Isopropalin	33820-53-0	250MG	N-12728-250MG	Ormetoprim	6981-18-6	250MG
N-12278-250MG	Isothiolane	50512-35-1	250MG	N-12729-100MG	Oryzalin	19044-88-3	100MG
N-12279-100MG	Isopturon	34123-59-6	100MG	N-12730-250MG	Ovex	80-33-1	250MG
N-12285-250MG	Isoxaben	82558-50-7	250MG	N-12731-250MG	Oxadiazon	19666-30-9	250MG
N-12286-100MG	Isxadifen-ethyl	163520-33-0	100MG	N-12732-100MG	Oxadixyl	77732-09-3	100MG
N-12287-100MG	Isxafflute	141112-29-0	100MG	N-12735-100MG	Oxyamyl	23135-22-0	100MG
N-12289-100MG	Ivermectin	70288-86-7	100MG	N-12738-250MG	Oxolinic acid	14698-29-4	250MG
N-12290-100MG	Dicofol	115-32-2	100MG	N-12740-100MG	Oxycarboxin	5259-88-1	100MG
N-12291-100MG	Chlordecone	143-50-0	100MG	N-12741-50MG	Oxydemeton-methyl	301-12-2	50MG
N-12294-250MG	Kinetin	525-79-1	250MG	N-12742-250MG	Oxyfluorfen	42874-03-3	250MG
N-12295-100MG	Kinoprene	42588-37-4	100MG	N-12743-250MG	Oxytetracycline hydrochloride	2058-46-0	250MG
N-12296-250MG	Kresoxim-methyl	143390-89-0	250MG	N-12757-250MG	p-Bromophenoxy acetic acid	1878-91-7	250MG
N-12304-100MG	Lactofen	77501-63-4	100MG	N-12761-250MG	p-Chloromandelic acid	492-86-4	250MG
N-12307-100MG	lambda-Cyhalothrin	91465-08-6	100MG	N-12763-250MG	p-Chlorophenoxyacetic acid	122-88-3	250MG
N-12310-100MG	Lead arsenate(Technical)	7784-40-9	100MG	N-12811-500MG	Paclobutrazol	76738-62-0	500MG
N-12314-100MG	Lenacil	2164-08-100	100MG	N-12816-100MG	Paraoxon	311-45-5	100MG
N-12315-100MG	Leptophos	21609-90-5	100MG	N-12818-500MG	Paraquat dichloride hydrate	1910-42-5	500MG
N-12319-250MG	Lindane (BHC gamma isomer)	58-89-9	250MG	N-12819-100MG	Parathion	56-38-2	100MG
N-12322-250MG	Linuron	330-55-2	250MG	N-12822-100MG	Penconazole	66246-88-6	100MG
N-12324-250MG	Clopyralid	1702-17-6	250MG	N-12823-50MG	Penoxsulam	219714-96-2	50MG
N-12325-100MG	Lufenuron	103055-07-8	100MG	N-12824-250MG	Dienochlor	2227-17-0	250MG
N-12346-100MG	Malathion	121-75-5	100MG	N-12826-100MG	Pentachloroanisole	1825-21-4	100MG
N-12350-250MG	Maleic hydrazide	123-33-1	250MG	N-12830-250MG	Pentachloronitrobenzene	82-68-8	250MG
N-12353-250MG	Mancozeb (Technical)	8018-01-700	250MG	N-12831-500MG	Pentachlorophenol	87-86-5	500MG
N-12354-100MG	Mandipropamid	374726-62-2	100MG	N-12844-100MG	Pentanochlor	2307-68-8	100MG
N-12355-250MG	Maneb(Technical)		250MG	N-12845-250MG	Perchloromethyl mercaptan	594-42-3	250MG
N-12357-100MG	MCBP	94-81-5	100MG	N-12848-250MG	Permethrin	52645-53-1	250MG
N-12358-250MG	MCPB methyl ester	57153-18-1	250MG	N-12849-250MG	Ethylan	72-56-0	250MG
N-12359-50MG	Mecarbam	2595-54-2	50MG	N-12854-250MG	Fenamiphos	22224-92-6	250MG
N-12360-100MG	Mecoprop	7085-19-0	100MG	N-12859-250MG	Phenmedipham	13684-63-4	250MG
N-12361-100MG	Mecoprop methyl ester	23844-56-6	100MG	N-12861-100MG	2,3-D methyl ester	103-26-4	100MG
N-12363-250MG	Mecoprop-P	16484-77-8	250MG	N-12863-100MG	Sodium dimethyldithiocarbamate hydrate	207233-95-2	100MG
N-12364-250MG	Mefenacet	73250-68-7	250MG	N-12870-100MG	Ethiprole	181587-01-9	100MG
N-12365-100MG	Mefenoxam	70630-17-0	100MG	N-12881-50MG	Cyanofenphos	13067-93-1	50MG
N-12366-100MG	Mefenpyr-diethyl	135590-91-9	100MG	N-12883-10MG	Sulfoxaflo	946578-00-3	10MG
N-12367-100MG	Mefluidide	53780-34-0	100MG	N-12884-100MG	Thifluzamide	130000-40-7	100MG
N-12369-100MG	Mepanipyrim	110235-47-7	100MG	N-12886-10MG	Cyantraniliprole	736994-63-1	10MG
N-12370-100MG	Mephosfolan	950-10-7	100MG	N-12888-50MG	2,4,6-Tris-(2,4,6-tribromophenoxy)-1,3,5-triazine	25713-60-4	50MG
N-12371-100MG	Mepiquat chloride	24307-26-4	100MG	N-12890-100MG	Sulcotrione	99105-77-8	100MG

Part #	Description	CAS #	Size	Part #	Description	CAS #	Size
N-12897-100MG	Metrafenone	220899-03-6	100MG	N-13139-100MG	Flumethrin	69770-45-2	100MG
N-12898-50MG	Chlorfenprop-methyl	14437-17-3	50MG	N-13140-100MG	Pendimethalin	40487-42-1	100MG
N-12901-100MG	Benazolin	3813-05-600	100MG	N-13141-250MG	Pyrimethalin	123312-89-0	250MG
N-12906-50MG	Z-Mevinphos	338-45-4	50MG	N-13142-100MG	Pyraclafos	89784-60-1	10MG
N-12908-10MG	4-(2-phenylethoxy)-quinazoline	124427-60-7	10MG	N-13143-100MG	Pyraclastrobin	175013-18-0	100MG
N-12911-5MG	Salannin	992-20-1	5MG	N-13144-100MG	Pyraflufen-ethyl	129630-19-9	100MG
N-12912-5MG	Nimbin	5945-86-8	5MG	N-13145-250MG	Pyrazon	1698-60-8	250MG
N-12920-100MG	Emamectin benzoate	155569-91-8	100MG	N-13146-100MG	Pyrazophos	13457-18-6	100MG
N-12921-100MG	Fluxapyroxad	907204-31-3	100MG	N-13147-100MG	Pyrazosulfuron-ethyl	93697-74-6	100MG
N-12922-100MG	Methylarsonic acid	124-58-3	100MG	N-13148-100MG	Pyrazoxyfen	71561-11-0	100MG
N-12944-10MG	Buturon	3766-60-7	10MG	N-13151-100MG	Pyrethrum	8003-34-7	100MG
N-12945-250MG	Karbutilate	4849-32-5	250MG	N-13152-25MG	Pyridaben	96489-71-3	25MG
N-12947-10MG	Kadethrine	58769-20-3	10MG	N-13153-50MG	Pyridalyl	179101-81-6	50MG
N-12948-100MG	Dodemorph acetate	31717-87-0	100MG	N-13154-250MG	Pyridaphenthion	119-12-0	250MG
N-12950-250MG	Sebutylazin	7286-69-3	250MG	N-13155-500MG	Pyridate	55512-33-9	500MG
N-12951-10MG	Orbencarb	34622-58-7	10MG	N-13158-100MG	Pyrifluquinazon	337458-27-2	100MG
N-12953-10MG	2,3-Dichlorophenoxyacetic acid	2976-74-1	10MG	N-13159-100MG	Pyrimethanil	53112-28-0	100MG
N-12954-10MG	Etrinfos	38260-54-7	10MG	N-13160-100MG	Pyriproxyfen	95737-68-1	100MG
N-12955-10MG	Dibam	128-04-1	10MG	N-13161-500MG	Pyriothiobac-Sodium	123343-16-8	500MG
N-12957-10MG	Benoxacor	98730-04-2	10MG	N-13164-100MG	Pyroxsulam	422556-08-9	100MG
N-12959-10MG	Albendazole	54965-21-8	10MG	N-13168-250MG	Quinalphos	13593-03-8	250MG
N-12960-50MG	Simeton	673-04-1	50MG	N-13169-500MG	Quinciorac	84087-01-4	500MG
N-12961-10MG	Bromofenoxim	13181-17-4	10MG	N-13172-100MG	Quinomethionate	2439-01-200	100MG
N-12963-100MG	Butocarboxim	34681-10-2	100MG	N-13173-100MG	Quinoxifen	124495-18-7	100MG
N-12969-100MG	Chlorthiophos	60238-56-4	100MG	N-13174-100MG	Quizalofop ethyl	76578-14-8	100MG
N-12970-10MG	Chlorbufam	1967-16-4	10MG	N-13175-100MG	Quizalofop-P-ethyl	100646-51-3	100MG
N-12971-10MG	Carbetamide	16118-49-3	10MG	N-13176-250MG	Resmethrin	10453-86-8	250MG
N-12972-10MG	Dimefuron	34205-21-5	10MG	N-13180-100MG	Rimsulfuron	122931-48-0	100MG
N-12973-10MG	Desmetryn	1014-69-3	10MG	N-13184-250MG	Rotenone	83-79-4	250MG
N-12976-10MG	Dinitramine	29091-05-2	10MG	N-13186-50MG	Rutene (TM)	299-86-5	50MG
N-12977-10MG	Difenoxuron	14214-32-5	10MG	N-13187-250MG	s-Bioallethrin	28434-00-6	250MG
N-12978-10MG	Furalaxyl	57646-30-7	10MG	N-13188-500MG	S-Ethyl dipropylthiocarbamate	759-94-4	500MG
N-12979-10MG	Furmecyclox	60568-05-0	10MG	N-13189-50MG	S-Hydroprene	65733-18-8	50MG
N-12984-10MG	Dimepax	22936-75-0	10MG	N-13191-100MG	S-Metolachlor	87392-12-9	100MG
N-12985-100MG	Etaconazole	60207-93-4	100MG	N-13194-250MG	S,S,S-Tributyl phosphorotrithioate	78-48-8	250MG
N-12986-10MG	Ethirimol	23947-60-6	10MG	N-13195-100MG	Saflufenacil	372137-35-4	100MG
N-12987-100MG	Flurochloridone	61213-25-0	100MG	N-13197-100MG	Propetamphos	31218-83-4	100MG
N-12989-250MG	Piperophos	24151-93-7	250MG	N-13201-250MG	Fenvalerate	51630-58-1	250MG
N-12990-10MG	Pyriphenox	88283-41-4	10MG	N-13210-100MG	Sethoxydim	74051-80-2	100MG
N-12991-10MG	Benodanil	15310-01-7	10MG	N-13211-100MG	Siduron	1982-49-6	100MG
N-12992-10MG	Chlorthiamid	1918-13-4	10MG	N-13212-250MG	Silvex	93-72-1	250MG
N-12996-10MG	2,6-D Acid	575-90-6	10MG	N-13213-100MG	Silvex isooctyl ester	32534-95-5	100MG
N-12998-10MG	Pyroquilon	57369-32-1	10MG	N-13214-250MG	Silvex methyl ester	4841-20-7	250MG
N-13004-250MG	Phenoxyacetic acid	122-59-8	250MG	N-13215-250MG	Sodium arsenate	7778-43-0	250MG
N-13006-100MG	Phenthoate	2597-03-700	100MG	N-13216-250MG	Sodium fluoroacetate(Technical)	62-74-8	250MG
N-13007-250MG	Phentioacetate	900-95-8	250MG	N-13217-100MG	Fluridone	59756-60-4	100MG
N-13011-250MG	Phenyl mercuric acetate(Technical)	62-38-4	250MG	N-13221-100MG	Spinetoram	187166-40-1	100MG
N-13012-100MG	Phenyl mercuric benzoate	94-43-9	100MG	N-13222-100MG	Spinodas	168316-95-8	100MG
N-13013-250MG	Phenyl mercuric chloride(Technical)	100-56-1	250MG	N-13223-100MG	Spirodiclofen	148477-71-8	100MG
N-13014-250MG	Phenyl mercuric hydroxide(Technical)	100-57-2	250MG	N-13224-100MG	Spiromesifen	283594-90-1	100MG
N-13015-250MG	Phenyl mercuric lactate	122-64-5	250MG	N-13225-100MG	Spirotetramat	203313-25-1	100MG
N-13016-100MG	Phenyl mercuric phthalate	84-70-8	100MG	N-13226-100MG	Spiroxamine	118134-30-8	100MG
N-13020-250MG	Phenyl-2-thiourea	103-85-5	250MG	N-13228-250MG	Streptomycin sulfate	3810-74-0	250MG
N-13030-100MG	Phenylmercuric oleate(Technical)	104-60-9	100MG	N-13229-100MG	Strobane (TM)	8001-50-1	100MG
N-13031-100MG	Phenylmercury borate(Technical)	102-98-7	100MG	N-13240-250MG	Sulfabenzamide	127-71-9	250MG
N-13032-100MG	Phenylmercury nitrate(Technical)	55-68-5	100MG	N-13241-250MG	Sulfachloropyridazine	80-32-0	250MG
N-13033-250MG	Phenylmercury salicylate(Technical)	28086-13-7	250MG	N-13242-250MG	Sulfadimethoxine	122-11-2	250MG
N-13035-100MG	Phorate	298-02-2	100MG	N-13243-250MG	Sulfallate	95-06-7	250MG
N-13036-100MG	Phosalone	2310-17-0	100MG	N-13244-250MG	Sulfamerazine	127-79-7	250MG
N-13037-250MG	Mevinphos	7786-34-7	250MG	N-13245-250MG	Sulfamethazine	57-68-1	250MG
N-13038-100MG	Phosfolan	947-02-4	100MG	N-13246-250MG	Sulfamethizole	144-82-1	250MG
N-13039-250MG	Phosphamidon	13171-21-6	250MG	N-13247-250MG	Sulfamethoxazole	723-46-6	250MG
N-13040-50MG	Phosphan	115-78-6	50MG	N-13249-250MG	Sulfanilic acid	121-57-3	250MG
N-13041-100MG	Phoxim	14816-18-3	100MG	N-13250-250MG	Sulfapyridine	144-83-2	250MG
N-13043-250MG	Phthalanillic acid	4727-29-1	250MG	N-13251-10MG	Sulfaquinoxaline	59-40-5	10MG
N-13050-250MG	Picloram	1918-02-100	250MG	N-13252-250MG	Sulfathiazole	72-14-0	250MG
N-13051-100MG	Picloram methyl ester	14143-55-6	100MG	N-13253-250MG	Sulfentrazone	122836-35-5	250MG
N-13056-100MG	Thifensulfuron-methyl	79277-27-3	100MG	N-13254-100MG	Sulfometuron methyl	74222-97-2	100MG
N-13057-100MG	Pinoxaden	243973-20-8	100MG	N-13255-100MG	Sulfosulfuron	141776-32-1	100MG
N-13058-100MG	Piperalin	3478-94-2	100MG	N-13257-100MG	Sulfoxide	120-62-7	100MG
N-13061-100MG	Piperonyl butoxide	51-03-6	100MG	N-13258-100MG	Sulprofos	35400-43-2	100MG
N-13062-250MG	Pririmicarb	23103-98-2	250MG	N-13259-100MG	Super-Rozol	28772-56-7	100MG
N-13063-250MG	Pririmiphos-ethyl	23505-41-1	250MG	N-13260-250MG	SWEP	1918-18-9	250MG
N-13064-250MG	Pririmiphos-methyl	29232-93-7	250MG	N-13261-100MG	Myclobutanil	88671-89-0	100MG
N-13068-100MG	Sedaxane	874967-67-6	100MG	N-13263-100MG	tau-Fluvalinate	102851-06-9	100MG
N-13070-100MG	Pyroxasulfone	447399-55-5	100MG	N-13265-250MG	Fenchlorim	3740-92-9	250MG
N-13074-100MG	Indaziflam	950782-86-2	100MG	N-13274-10MG	MPMC	2425-10-700	10MG
N-13075-100MG	Amicarbazone	129909-90-6	100MG	N-13275-100MG	XMC	2655-14-3	100MG
N-13076-100MG	Triclosan	3380-34-5	100MG	N-13280-25MG	Flumorph	211867-47-9	25MG
N-13079-10MG	Halfenprox	111872-58-3	10MG	N-13281-25MG	Valifenalate	283159-90-0	25MG
N-13083-50MG	Pyrazolynate	58011-68-0	50MG	N-13282-25MG	Tiadinil	223580-51-6	25MG
N-13084-25MG	Metamifop	256412-89-2	25MG	N-13283-25MG	Methasulfocarb	66952-49-6	25MG
N-13085-10MG	Fthalide	27355-22-2	10MG	N-13288-25MG	Meptyldinocap	131-72-6	25MG
N-13086-250MG	Metiram(Technical)	9006-42-2	250MG	N-13296-50MG	Mesoprazine	1824-09-5	50MG
N-13087-250MG	Potassium cyanate(Technical)	590-28-3	250MG	N-13500-10MG	TCMTB	21564-17-0	10MG
N-13088-10MG	Prallethrin	23031-36-9	10MG	N-13501-100MG	Tebufenozide	112410-23-8	100MG
N-13089-250MG	Pretilachlor	51218-49-6	250MG	N-13502-100MG	Tebufenpyrad	191968-77-3	100MG
N-13090-100MG	Primisulfuron-methyl	86209-51-0	100MG	N-13503-100MG	Tebupirifos	96182-53-5	100MG
N-13092-500MG	Probenazole	27605-76-1	500MG	N-13504-250MG	Tebutam	35256-85-0	250MG
N-13093-250MG	Prochloraz	67747-09-5	250MG	N-13505-250MG	Tebuthiuron	34014-18-1	250MG
N-13094-100MG	Procymazine	32889-48-8	100MG	N-13506-250MG	Teflubenzuron	83121-18-0	250MG
N-13095-250MG	Procymidone	32809-16-8	250MG	N-13508-100MG	Tembotrione	335104-84-2	100MG
N-13096-100MG	Prodiame	29091-21-2	100MG	N-13509-250MG	Terbacil	5902-51-2	250MG
N-13097-250MG	Profenofos	41198-08-7	250MG	N-13510-250MG	Terbufos	13071-79-9	250MG
N-13098-250MG	Profluralin	26399-36-0	250MG	N-13511-250MG	Terbumeton	33693-04-8	250MG
N-13099-100MG	Prohexadione	88805-35-0	100MG	N-13512-50MG	Terbutylazine	5915-41-3	50MG
N-13100-100MG	Prohexadione-Calcium	127277-53-6	100MG	N-13513-100MG	Terbutol	1918-11-200	100MG
N-13101-100MG	Promecarb	2631-37-0	100MG	N-13514-250MG	Terbutryne	886-50-0	250MG
N-13102-250MG	Prometon	1610-18-0	250MG	N-13517-250MG	Etridiazole	2593-15-9	250MG
N-13103-250MG	Prometryne	7287-19-6	250MG	N-13537-250MG	Tetrachlorvinphos	22248-79-9	250MG
N-13104-250MG	Propachlor	1918-16-7	250MG	N-13538-100MG	Tetraconazole	112281-77-3	100MG
N-13105-100MG	Propamocarb	24579-73-5	250MG	N-13539-250MG	Tetradifon	116-29-0	250MG
N-13106-250MG	Propamocarb HCl	25606-41-1	250MG	N-13540-50MG	Tetraethylthiopyrophosphate	3689-24-5	50MG
N-13108-100MG	Propazine	139-40-2	100MG	N-13543-500MG	Tetraethylpyrophosphate	107-49-3	500MG
N-13109-250MG	Propham	122-42-9	250MG	N-13552-100MG	Tetramethrin	7696-12-0	100MG
N-13110-250MG	Prophos	13194-48-4	250MG	N-13559-250MG	Thallium sulfate(Technical)	7446-18-6	250MG
N-13111-250MG	Propineb(Technical)	9016-72-2	250MG	N-13560-250MG	Thiabendazole	148-79-8	250MG
N-13119-100MG	Propoxyacarbazine-sodium	181274-15-7	100MG	N-13561-100MG	Thiacloprid	111988-49-9	100MG
N-13129-25MG	Fluthiacet-methyl	117337-19-6	25MG	N-13562-100MG	Thiamethoxam	153719-23-4	100MG
N-13131-100MG	Anilofos	64249-01-0	100MG	N-13563-500MG	Thiazopry	117718-60-2	500MG
N-13132-25MG	Cyflufenamid	180409-60-3	25MG	N-13564-250MG	Thidiazuron	51707-55-2	250MG
N-13136-250MG	Propyzamide	23950-58-5	250MG	N-13565-50MG	Thiencarbazine-methyl	317815-83-1	50MG
N-13137-100MG	Prosulfuron	94125-34-5	100MG	N-13568-250MG	Thiodicarb	59669-26-0	250MG
N-13138-100MG	Prothioconazole	178928-70-6	100MG	N-13570-50MG	Thiometon	640-15-3	50MG

Part #	Description	CAS #	Size	Part #	Description	CAS #	Size
N-13571-100MG	Thiophanate	23564-06-9	100MG	N-14032-50MG	Uniconazole	83657-22-1	50MG
N-13572-250MG	Thiophanate-methyl	23564-05-8	250MG	N-14035-100MG	Penflufen	494793-67-8	100MG
N-13576-250MG	Propiconazole	60207-90-1	250MG	N-14037-100MG	Methothrin	34388-29-9	100MG
N-13577-100MG	Prothiophos	34643-46-4	100MG	N-14038-25MG	Cyhalofop	122008-78-0	25MG
N-13578-500MG	Tolclofos-Methyl	57018-04-9	500MG	N-14058-25MG	Thiofanox	39196-18-4	25MG
N-13579-100MG	Tolfenpyrad	129558-76-5	100MG	N-14059-100MG	Fentrazamide	158237-07-1	100MG
N-13584-250MG	Tolyfluanid	731-27-1	250MG	N-14068-10MG	Tecloftalam	76280-91-6	10MG
N-13585-100MG	Topramezone	210631-68-8	100MG	N-14070-100MG	Pyribenzoxim	168088-61-7	100MG
N-13586-250MG	Camphochlor	8001-35-2	250MG	N-14076-50MG	Benthiavalicarb-isopropyl	177406-68-7	50MG
N-13587-500MG	Trailoxoydim	87820-88-0	500MG	N-14078-100MG	1,3-di-p-tolylurea	621-00-1	100MG
N-13588-100MG	Talomehrin	66841-25-6	100MG	N-14079-100MG	Isofetamid	875915-78-9	100MG
N-13615-10MG	trans-Chlordane	5103-74-2	10MG	N-14080-50MG	Proquinazid	189278-12-4	50MG
N-13618-50MG	trans-Heptachlor epoxide	28044-83-9	50MG	N-14084-100MG	Fluensulfone	318290-98-1	100MG
N-13620-50MG	trans-Permethrin	61949-77-7	50MG	N-17917-50MG	Penthiopyrad	183675-82-3	50MG
N-13626-100MG	Transfluthrin	118712-89-3	100MG	N-FC92-5-5MG	4,4'-DDT (13C12)	104215-84-1	5MG
N-13628-250MG	Tri-alleate	2303-17-5	250MG	N-FC93-5-5MG	4,4'-DDE (13C12)		5MG
N-13636-500MG	Triadimefon	43121-43-3	500MG	N-FD202-E-0.01G	a-Endosulfan (d4)	203645-57-2	0.01G
N-13638-100MG	Triasulfuron	82097-50-5	100MG	N-FD203-5-5MG	b-Endosulfan (d4)	203716-99-8	5MG
N-13639-100MG	Triazophos	24017-47-8	100MG	N-FD2061-E-0.01G	Dichlorvos (dimethyl-d6)	203645-53-8	0.01G
N-13640-250MG	Tribenuron methyl	101200-48-0	250MG	N-FD2208-5-5MG	Atrazine (ethylamine-d5)	163165-75-1	5MG
N-13648-250MG	Tributylphosphoro-trithioite	150-50-5	250MG	N-FD998-E-0.01G	Parathion (TM) (diethyl-d10)	350820-04-1	0.01G
N-13652-250MG	Trichloroacetic acid	76-03-9	250MG	S-10005M1-1ML	(2,4,5-Trichlorophenoxy)acetic acid methyl ester Solution	1928-37-6	1ML
N-13657-100MG	Trichloronate	327-98-0	100MG	S-10005M1-5ML	(2,4,5-Trichlorophenoxy)acetic acid methyl ester Solution	1928-37-6	5ML
N-13659-500MG	Triclopyr	55335-06-3	500MG	S-10012A1-1ML	1-8-Naphthalic anhydride Solution	81-84-5	1ML
N-13660-250MG	Triclopyr 2-butoxyethyl ester	64700-56-7	250MG	S-10012U1-1ML	1-8-Naphthalic anhydride Solution	81-84-5	1ML
N-13661-100MG	Triclopyr methyl ester	60825-26-5	100MG	S-10045M1-1ML	1-cis-3-Chloroallyl-357-triaza-1-azonia-adamantane chloride	51229-78-8	1ML
N-13662-250MG	Tricyclohexylin hydroxide	13121-70-5	250MG	S-10074U1-1ML	1-Naphthaleneacetamide Solution	86-86-2	1ML
N-13665-100MG	Tridemorph (Mix of isomers)	81412-43-3	100MG	S-10074U1-5ML	1-Naphthaleneacetamide Solution	86-86-2	5ML
N-13666-100MG	Trietazine	1912-26-1	100MG	S-10075T1-1ML	1-Naphthaleneacetic acid Solution	86-87-3	1ML
N-13684-100MG	Trifloxystrobin	141517-21-7	100MG	S-10075T1-5ML	1-Naphthaleneacetic acid Solution	86-87-3	5ML
N-13685-100MG	Trifloxysulfuron-Sodium	199119-58-9	100MG	S-10111M1-1ML	1,2,4,5-Tetrachloro-3-nitrobenzene Solution	117-18-0	1ML
N-13686-250MG	Triflumizole	68694-11-1	250MG	S-10122M1-1ML	1,1-Dichloro-1-nitroethane Solution	594-72-9	1ML
N-13687-100MG	Triflumuron	64628-44-0	100MG	S-10122X1-1ML	1,1-Dichloro-1-nitroethane Solution	594-72-9	1ML
N-13689-250MG	Trifluralin	1582-09-8	250MG	S-10149M1-1ML	1,2-Dibromo-3-chloropropane Solution	96-12-8	1ML
N-13690-50MG	Triflurosulfuron-methyl	126535-15-7	50MG	S-10149M1-5ML	1,2-Dibromo-3-chloropropane Solution	96-12-8	5ML
N-13691-250MG	Triforine	26644-46-2	250MG	S-10150M1-1ML	1,2-Dibromomethane Solution	106-93-4	1ML
N-13706-250MG	Trinepaxac-ethyl	95266-40-3	250MG	S-10152M1-1ML	1,2-Dichlorobenzene Solution	95-50-1	1ML
N-13729-100MG	Triticonazole	131983-72-7	100MG	S-10152M1-5ML	1,2-Dichlorobenzene Solution	95-50-1	5ML
N-13731-100MG	Ethiozin	64529-56-2	100MG	S-10156M1-1ML	1,2-Dichloropropane Solution	78-87-5	1ML
N-13735-10MG	Uniconazole-P	83657-17-4	10MG	S-10156M1-5ML	1,2-Dichloropropane Solution	78-87-5	5ML
N-13738-100MG	Pyriminuron	53558-25-1	100MG	S-10193M1-1ML	1,3-Dichloropropene Solution	542-75-6	1ML
N-13741-100MG	Vamidothion	2275-23-2	100MG	S-10216M1-1ML	1,4-Dichlorobenzene Solution	106-46-7	1ML
N-13742-250MG	Vanicide-20S (TM)	5902-85-2	250MG	S-10216M1-5ML	1,4-Dichlorobenzene Solution	106-46-7	5ML
N-13743-250MG	Velpar	51235-04-2	250MG	S-10254A1-1ML	2-(3-Chlorophenoxy)propionic acid Solution	101-10-0	1ML
N-13744-250MG	Vernolate	1929-77-7	250MG	S-10254T1-1ML	2-(3-Chlorophenoxy)propionic acid Solution	101-10-0	1ML
N-13745-250MG	Vinclozolin	50471-44-8	250MG	S-10260A1-1ML	2-(2-Butoxyethoxy)ethyl thiocyanate Solution	112-56-1	1ML
N-13750-250MG	Warfarin	81-81-2	250MG	S-10260J4-1ML	2-(2-Butoxyethoxy)ethyl thiocyanate Solution	112-56-1	1ML
N-13752-100MG	Z-Dimethomorph	113210-98-3	100MG	S-10260J4-5ML	2-(2-Butoxyethoxy)ethyl thiocyanate Solution	112-56-1	5ML
N-13753-500MG	Zectran	315-18-4	500MG	S-10261M1-1ML	2-(3-Chlorophenoxy)propionamide Solution	5825-87-6	1ML
N-13754-250MG	zeta-Cypermethrin	52315-07-8	250MG	S-10280A1-1ML	2-Biphenylol Solution	90-43-7	1ML
N-13758-250MG	Zinc phosphate(Technical)	7779-90-0	250MG	S-10280M5-1ML	2-Biphenylol Solution	90-43-7	1ML
N-13759-250MG	Zineb(Technical)	112122-67-7	250MG	S-10306M1-1ML	2-Butyl-2-ethyl-1,3-propanediol Solution	115-84-4	1ML
N-13760-50MG	Zinophos (TM)	297-97-2	50MG	S-10346B6-1ML	2-Ethyl-1,3-hexanediol Solution	94-96-2	1ML
N-13761-250MG	Ziram	137-30-4	250MG	S-10346B6-5ML	2-Ethyl-1,3-hexanediol Solution	94-96-2	5ML
N-13762-100MG	Zoxamide	156052-68-5	100MG	S-10356A1-1ML	2-Ethylhexyl-2,4-dichlorophenoxy acetate Solution	1928-43-4	1ML
N-13800-100MG	Simazine	122-34-9	100MG	S-10356J1-1ML	2-Ethylhexyl-2,4-dichlorophenoxy acetate Solution	1928-43-4	1ML
N-13802-50MG	Simetryn	1014-70-6	50MG	S-10358M1-1ML	2-Fluoroacetamide Solution	640-19-7	1ML
N-13812-250MG	Fenson	80-38-6	250MG	S-10363A1-1ML	2-Furanacrylic acid Solution	539-47-9	1ML
N-13813-250MG	Flurenol-methyl	1216-44-0	250MG	S-10363T1-1ML	2-Furanacrylic acid Solution	539-47-9	1ML
N-13815-250MG	Flurenol	467-69-6	250MG	S-10378M1-1ML	2-Hydroxyethyl-n-octyl sulphide Solution	3547-33-9	1ML
N-13816-100MG	Fluoxastrobin	361377-29-9	100MG	S-10381A7-1ML	2-Imidazolidinethione Solution	96-45-7	1ML
N-13817-10MG	Furilazole	121776-33-8	10MG	S-10381A7-5ML	2-Imidazolidinethione Solution	96-45-7	5ML
N-13818-100MG	Imazosulfuron	122548-33-8	100MG	S-10381H1-1ML	2-Imidazolidinethione Solution	96-45-7	1ML
N-13819-100MG	Aquidosulfuron	120923-37-7	100MG	S-10381H1-5ML	2-Imidazolidinethione Solution	96-45-7	5ML
N-13820-50MG	Clomintocet	88349-88-6	50MG	S-10384A1-1ML	2-Isovaleryl-1,3-indandione Solution	83-28-3	1ML
N-13821-250MG	Dichlorprop 2-ethylhexyl ester	79270-78-3	250MG	S-10384U1-1ML	2-Isovaleryl-1,3-indandione Solution	83-28-3	1ML
N-13822-100MG	Diffufenzopyr sodium salt	109293-98-3	100MG	S-10437M1-1ML	2-Naphthoxyacetic acid Solution	120-23-0	1ML
N-13823-50MG	Propisochlor	86763-47-5	50MG	S-10437T1-1ML	2-Naphthoxyacetic acid Solution	120-23-0	1ML
N-13824-100MG	Pentachlorophenyl acetate	1441-02-7	100MG	S-10507M1-1ML	2-Phenoxypropionic acid Solution	940-31-8	1ML
N-13825-100MG	Quinoclamine	2797-51-5	100MG	S-10507T1-1ML	2-Phenoxypropionic acid Solution	940-31-8	1ML
N-13826-250MG	Aminocyclopyrachlor	858956-08-8	250MG	S-10512A1-1ML	2-Pivaloyl-1,3-indandione Solution	83-26-1	1ML
N-13830-50MG	Fenoxaprop	95617-09-7	50MG	S-10512U1-1ML	2-Pivaloyl-1,3-indandione Solution	83-26-1	1ML
N-13832-100MG	(S)-Methoprene	65733-16-6	100MG	S-10517F1-1ML	2-sec-Butyl-4,6-dinitrophenol triethanolamine salt Solution	6420-47-9	1ML
N-13904-50MG	Asomate	3586-60-5	50MG	S-10517U1-1ML	2-sec-Butyl-4,6-dinitrophenol triethanolamine salt Solution	6420-47-9	1ML
N-13906-50MG	Fluazifop-methyl	69335-90-6	50MG	S-10520N1-1ML	2-Undecanone Solution	112-12-9	1ML
N-13907-50MG	Fluazifop-P	83066-88-0	50MG	S-10530M1-1ML	2,3,5,6-Tetrachloroaniline Solution	3481-20-7	1ML
N-13908-50MG	Fluazifop-P-methyl	69335-90-6	50MG	S-10531A4-1ML	2,4-D 2-ethylhexyl ester Solution	1928-43-4	1ML
N-13913-10MG	Empenthrin	54406-48-3	10MG	S-10531A4-5ML	2,4-D 2-ethylhexyl ester Solution	1928-43-4	5ML
N-13914-100MG	Silafluofen	105024-66-6	100MG	S-10532A1-1ML	2,4-D isobutyl ester Solution	1713-15-1	1ML
N-13916-10MG	Cycloxydime	101205-02-1	10MG	S-10532U1-1ML	2,4-D isobutyl ester Solution	1713-15-1	1ML
N-13917-10MG	Monalide	7287-36-7	10MG	S-10533M1-1ML	2,4-DB methyl ester Solution	18625-12-2	1ML
N-13918-10MG	Prosulfofcarb	52888-80-9	10MG	S-10533M1-5ML	2,4-DB methyl ester Solution	18625-12-2	5ML
N-13919-10MG	S 421	127-90-2	10MG	S-10550A1-1ML	2,6-Diisopropyl-1-naphthalene Solution	24157-81-1	1ML
N-13921-100MG	Butoxycarboxim	34681-23-7	100MG	S-10550U1-1ML	2,6-Diisopropyl-1-naphthalene Solution	24157-81-1	1ML
N-13922-10MG	Florasulam	145701-23-1	10MG	S-10584A1-1ML	2,3-Dichloronitrobenzene Solution	3209-22-1	1ML
N-13923-10MG	Iodosulfuron-methyl-sodium	144550-36-7	10MG	S-10584U1-1ML	2,3-Dichloronitrobenzene Solution	3209-22-1	1ML
N-13924-10MG	Iprovalicarb	140923-17-7	10MG	S-10599M1-1ML	2,3,4,6-Tetrachlorophenol Solution	58-90-2	1ML
N-13925-10MG	Quinmerac	90717-03-6	10MG	S-10601A1-1ML	2,3,5-Triiodobenzoic acid Solution	88-82-4	1ML
N-13927-100MG	Picoxystrobin	117428-22-5	100MG	S-10601T1-1ML	2,3,5-Triiodobenzoic acid Solution	88-82-4	1ML
N-13928-100MG	Acrinathrin	101007-06-1	100MG	S-10603A1-1ML	2,3,5-Trimethylphenyl methyl carbamate Solution	2655-15-4	1ML
N-13931-25MG	Simeconazole	14958-90-7	25MG	S-10609B1-1ML	2,4-D Solution	94-75-7	1ML
N-13932-25MG	Isoxathion	18854-01-8	25MG	S-10609M1-1ML	2,4-D Solution	94-75-7	1ML
N-13939-250MG	Flamprop-M-isopropyl	57973-67-8	250MG	S-10610A1-1ML	2,4-D butoxyethyl ester Solution	1929-73-3	1ML
N-13940-250MG	Fluxofenim	88485-37-4	250MG	S-10610A4-1ML	2,4-D butoxyethyl ester Solution	1929-73-3	1ML
N-13941-100MG	Flupyradifurone	951659-40-8	100MG	S-10610A4-5ML	2,4-D butoxyethyl ester Solution	1929-73-3	5ML
N-13942-100MG	Nuarimol	63284-71-9	100MG	S-10611A1-1ML	2,4-D butyl ester Solution	94-80-4	1ML
N-13948-100MG	Tebuthiuron TP 108	39222-73-6	100MG	S-10611U1-1ML	2,4-D butyl ester Solution	94-80-4	1ML
N-13955-100MG	Daimuron	42609-52-9	100MG	S-10612M1-1ML	2,4-D dimethylamine salt Solution	2008-39-1	1ML
N-13956-100MG	Triazamate	112143-82-5	100MG	S-10613M1-1ML	2,4-D ethanamine salt Solution	3599-58-4	1ML
N-13958-100MG	Itraconazole	84625-61-6	100MG	S-10614A1-1ML	2,4-D ethyl ester Solution	533-23-3	1ML
N-13972-50MG	Cyflumetofen	400882-07-7	50MG	S-10614J1-1ML	2,4-D ethyl ester Solution	533-23-3	1ML
N-13985-100MG	2,4-D Choline salt	1048373-72-3	100MG	S-10615A1-1ML	2,4-D isooctyl ester Solution	25168-26-7	1ML
N-13988-50MG	Fuopyram	658066-35-4	50MG	S-10615T1-1ML	2,4-D isooctyl ester Solution	25168-26-7	1ML
N-13989-100MG	Trietazine-ethylamine	93567-85-2	100MG	S-10616M1-1ML	2,4-D isopropanol amine salt Solution	6365-72-6	1ML
N-13991-250MG	Aclofenfen	74070-46-5	250MG	S-10617A1-1ML	2,4-D isopropyl ester Solution	94-11-1	1ML
N-13994-100MG	(1R)-trans-Tetramethrin	1166-46-7	100MG	S-10617U1-1ML	2,4-D isopropyl ester Solution	94-11-1	1ML
N-14001-100MG	Pyrasulfotole	365400-11-9	100MG	S-10618M1-1ML	2,4-D isopropylamine salt Solution	5742-17-6	1ML
N-14006-50MG	Imiprophrin	72963-72-5	50MG	S-10619M1-1ML	2,4-D methyl ester Solution	1928-38-7	1ML
N-14007-250MG	Propaquizafop	111479-05-1	250MG	S-10619M1-5ML	2,4-D methyl ester Solution	1928-38-7	5ML
N-14011-100MG	Oxadiazyl	39807-15-3	100MG	S-10620A1-1ML	2,4-D propylene glycol ester Solution	103564-78-9	1ML
N-14012-250MG	Octhillinone	26530-20-1	250MG	S-10620J1-1ML	2,4-D propylene glycol ester Solution	103564-78-9	1ML
N-14013-250MG	Fenfuram	24691-80-3	250MG	S-10621A4-1ML	2,4-D sec-butyl ester Solution	94-75-7	1ML
N-14029-100MG	Flusulfamide	106917-52-6	100MG	S-10621A4-5ML	2,4-D sec-butyl ester Solution	94-75-7	5ML
N-14030-100MG	Ferimzone	89269-64-7	100MG	S-10621U1-1ML	2,4-D sec-butyl ester Solution	94-75-7	1

Part #	Description	CAS #	Size	Part #	Description	CAS #	Size
S-10622A1-1ML	2,4-DB Solution	94-82-6	1ML	S-11063J1-1ML	Amdro Solution	67485-29-4	1ML
S-10622B1-1ML	2,4-DB Solution	94-82-6	1ML	S-11065A4-1ML	Aminocarb Solution	2032-59-9	1ML
S-10622B1-5ML	2,4-DB Solution	94-82-6	1ML	S-11065A4-5ML	Aminocarb Solution	2032-59-9	5ML
S-10623M1-1ML	2,4-DB dimethylamine salt Solution	2758-42-1	1ML	S-11065U1-1ML	Aminocarb Solution	2032-59-9	1ML
S-10625A1-1ML	2,4-Dichloronitrobenzene Solution	611-06-3	1ML	S-11067A1-1ML	Aminotriazole Solution	61-82-5	1ML
S-10625U1-1ML	2,4-Dichloronitrobenzene Solution	611-06-3	1ML	S-11067B1-1ML	Aminotriazole Solution	61-82-5	1ML
S-10628A1-1ML	2,4-Dichlorophenyl benzenesulfonate Solution	97-16-5	1ML	S-11068A1-1ML	Amiltraz Solution	33089-61-1	1ML
S-10628U1-1ML	2,4-Dichlorophenyl benzenesulfonate Solution	97-16-5	1ML	S-11071M1-1ML	Amoban Solution	3566-10-700	1ML
S-10646A1-1ML	2,4,5-T Solution	93-76-5	1ML	S-11073M1-1ML	Ancymidol Solution	12771-68-5	1ML
S-10646B1-1ML	2,4,5-T Solution	93-76-5	1ML	S-11075A1-1ML	Anilazine Solution	101-05-3	1ML
S-10646B1-5ML	2,4,5-T Solution	93-76-5	5ML	S-11075U1-1ML	Anilazine Solution	101-05-3	1ML
S-10647A1-1ML	2,4,5-T 2-ethylhexyl ester Solution	1928-47-8	1ML	S-11075U1-5ML	Anilazine Solution	101-05-3	5ML
S-10647T1-1ML	2,4,5-T 2-ethylhexyl ester Solution	1928-47-8	1ML	S-11085J1-1ML	Aramite Solution	140-57-8	1ML
S-10648A4-1ML	2,4,5-T butoxyethyl ester Solution	2545-59-7	1ML	S-11102A1-1ML	Esfenvalerate Solution	66230-04-4	1ML
S-10648A4-5ML	2,4,5-T butoxyethyl ester Solution	2545-59-7	5ML	S-11102T1-1ML	Esfenvalerate Solution	66230-04-4	1ML
S-10648U1-1ML	2,4,5-T butoxyethyl ester Solution	2545-59-7	1ML	S-11103A1-1ML	Aspon (TM) Solution	3244-90-4	1ML
S-10649A1-1ML	2,4,5-T isooctyl ester Solution	25168-15-4	1ML	S-11103J4-1ML	Aspon (TM) Solution	3244-90-4	1ML
S-10649T1-1ML	2,4,5-T isooctyl ester Solution	25168-15-4	1ML	S-11103J4-5ML	Aspon (TM) Solution	3244-90-4	5ML
S-10650A1-1ML	2,4,5-T isopropyl ester Solution	93-78-7	1ML	S-11104M1-1ML	Asulam Solution	3337-71-1	1ML
S-10650T1-1ML	2,4,5-T isopropyl ester Solution	93-78-7	1ML	S-11104M1-5ML	Asulam Solution	3337-71-1	5ML
S-10651A1-1ML	2,4,5-T n-butyl ester Solution	93-79-8	1ML	S-11105A1-1ML	Atraton Solution	1610-17-9	1ML
S-10651J1-1ML	2,4,5-T n-butyl ester Solution	93-79-8	1ML	S-11105T1-1ML	Atraton Solution	1610-17-9	1ML
S-10652M1-1ML	2,4,5-T triethylamine salt Solution	57213-69-1	1ML	S-11105T1-5ML	Atraton Solution	1610-17-9	5ML
S-10655F1-1ML	2,4,5-Trichlorophenol ethanolamine salt Solution		1ML	S-11110M1-1ML	Azobenzene Solution	103-33-3	1ML
S-10655U1-1ML	2,4,5-Trichlorophenol ethanolamine salt Solution		1ML	S-11110M1-5ML	Azobenzene Solution	103-33-3	5ML
S-10657M1-1ML	2,4,6-Trichlorophenol Solution	88-06-2	1ML	S-11111A1-1ML	Azocyclotin (Technical) Solution	41083-11-8	1ML
S-10657M1-5ML	2,4,6-Trichlorophenol Solution	88-06-2	5ML	S-11113A1-1ML	Azoxystrobin Solution	13860-33-8	1ML
S-10670A1-1ML	2,5-Dichloronitrobenzene Solution	89-61-2	1ML	S-11113U1-1ML	Azoxystrobin Solution	13860-33-8	1ML
S-10670U1-1ML	2,5-Dichloronitrobenzene Solution	89-61-2	1ML	S-11118A1-1ML	b-Endosulfan Solution	33213-65-9	1ML
S-10688M1-1ML	2,6-Dichlorobenzonitrile Solution	1194-65-6	1ML	S-11118J1-1ML	b-Endosulfan Solution	33213-65-9	1ML
S-10727A1-1ML	3-Indoleacetic acid Solution	87-51-4	1ML	S-11118J1-5ML	b-Endosulfan Solution	33213-65-9	5ML
S-10727T1-1ML	3-Indoleacetic acid Solution	87-51-4	1ML	S-11122A1-1ML	Carbanolate Solution	671-04-5	1ML
S-10728A1-1ML	3-Indolebutyric acid Solution	133-32-4	1ML	S-11122U1-1ML	Carbanolate Solution	671-04-5	1ML
S-10728T1-1ML	3-Indolebutyric acid Solution	133-32-4	1ML	S-11123M1-1ML	Barban Solution	101-27-9	1ML
S-10729A1-1ML	3-Indolepropionic acid Solution	830-96-6	1ML	S-11123M1-5ML	Barban Solution	101-27-9	5ML
S-10729T1-1ML	3-Indolepropionic acid Solution	830-96-6	1ML	S-11126M1-1ML	Baycarb (TM) Solution	3766-81-2	1ML
S-10770A1-1ML	3,5-Dichlorobenzoic acid Solution	51-36-5	1ML	S-11127A1-1ML	Bitertanol Solution	55179-31-2	1ML
S-10770T1-1ML	3,5-Dichlorobenzoic acid Solution	51-36-5	1ML	S-11127U1-1ML	Bitertanol Solution	55179-31-2	1ML
S-10770T1-5ML	3,5-Dichlorobenzoic acid Solution	51-36-5	5ML	S-11128A1-1ML	Propoxur Solution	114-26-1	1ML
S-10785A1-1ML	3,4-Dichloronitrobenzene Solution	99-54-7	1ML	S-11129A1-1ML	Triadimenol Solution	55219-65-3	1ML
S-10785U1-1ML	3,4-Dichloronitrobenzene Solution	99-54-7	1ML	S-11129U1-1ML	Triadimenol Solution	55219-65-3	1ML
S-10791A1-1ML	3,5-Dichlorobenzoic acid methyl ester Solution	2905-67-1	1ML	S-11130A4-1ML	Cyfluthrin Solution	68359-37-5	1ML
S-10791T1-1ML	3,5-Dichlorobenzoic acid methyl ester Solution	2905-67-1	1ML	S-11130A4-5ML	Cyfluthrin Solution	68359-37-5	5ML
S-10791T1-5ML	3,5-Dichlorobenzoic acid methyl ester Solution	2905-67-1	5ML	S-11130U1-1ML	Cyfluthrin Solution	68359-37-5	1ML
S-10792A1-1ML	3,5-Dichloronitrobenzene Solution	618-62-2	1ML	S-11131M1-1ML	BDMC Solution	672-99-1	1ML
S-10792U1-1ML	3,5-Dichloronitrobenzene Solution	618-62-2	1ML	S-11131M1-5ML	BDMC Solution	672-99-1	5ML
S-10794M1-1ML	3,5-Diiodo-4-hydroxybenzonitrile Solution	1689-83-4	1ML	S-11132A1-1ML	Tricyclazole Solution	41814-78-2	1ML
S-10800A1-1ML	Propanil Solution	709-98-8	1ML	S-11132T1-1ML	Tricyclazole Solution	41814-78-2	1ML
S-10800Z4-1ML	Propanil Solution	709-98-8	1ML	S-11132T1-5ML	Tricyclazole Solution	41814-78-2	5ML
S-10800Z4-5ML	Propanil Solution	709-98-8	5ML	S-11133M1-1ML	Benalaxyl Solution	71626-11-4	1ML
S-10803M1-1ML	4-Allyl-2-methoxyphenol Solution	97-53-0	1ML	S-11134A1-1ML	Benazolin-ethyl Solution	25059-80-7	1ML
S-10806M1-1ML	4-Aminopyridine Solution	504-24-5	1ML	S-11134T1-1ML	Benazolin-ethyl Solution	25059-80-7	1ML
S-10818B1-1ML	4-Chloro-o-tolylxyacetic acid Solution	94-74-6	1ML	S-11135U1-1ML	Bendiocarb Solution	22781-23-3	1ML
S-10818B1-5ML	4-Chloro-o-tolylxyacetic acid Solution	94-74-6	5ML	S-11136A1-1ML	Benfluralin Solution	1861-40-1	1ML
S-10818M1-1ML	4-Chloro-o-tolylxyacetic acid Solution	94-74-6	1ML	S-11136J4-1ML	Benfluralin Solution	1861-40-1	1ML
S-10819M1-1ML	4-Chloro-o-tolylxyacetic acid methyl ester Solution	2436-73-9	1ML	S-11136J4-5ML	Benfluralin Solution	1861-40-1	5ML
S-10819M1-5ML	4-Chloro-o-tolylxyacetic acid methyl ester Solution	2436-73-9	5ML	S-11139A1-1ML	Bensulfuron-methyl Solution	83055-99-6	1ML
S-10874A1-1ML	4,4'-DDD Solution	72-54-8	1ML	S-11140M4-1ML	Bensulfuron-methyl Solution	83055-99-6	1ML
S-10874K1-1ML	4,4'-DDD Solution	72-54-8	1ML	S-11140M4-5ML	Bensulfuron-methyl Solution	83055-99-6	5ML
S-10874K1-5ML	4,4'-DDD Solution	72-54-8	5ML	S-11141M1-1ML	Bensulfuron-methyl Solution	83055-99-6	1ML
S-10875M1-1ML	4,4'-DDE Solution	72-55-9	1ML	S-11142A4-1ML	Bentazon Solution	25057-89-0	1ML
S-10875M1-5ML	4,4'-DDE Solution	72-55-9	5ML	S-11142A4-5ML	Bentazon Solution	25057-89-0	5ML
S-10876J1-1ML	4,4'-DDT Solution	50-29-3	1ML	S-11142B1-1ML	Bentazon Solution	25057-89-0	1ML
S-10888F1-1ML	4,6-Dinitro-2-sec-butylphenol ethanolamine salt Solution		1ML	S-11143Y1-1ML	Bentazon methyl derivative Solution	61592-45-8	1ML
S-10888U1-1ML	4,6-Dinitro-2-sec-butylphenol ethanolamine salt Solution		1ML	S-11143Y1-5ML	Bentazon methyl derivative Solution	61592-45-8	5ML
S-10889M4-1ML	4,6-Dinitro-o-cresol Solution	534-52-1	1ML	S-11144M1-1ML	Benthiocarb Solution	28249-77-6	1ML
S-10956A1-1ML	6-Benzylaminopurine Solution	1214-39-7	1ML	S-11177A4-1ML	Benzoilprop ethyl Solution	22212-55-1	1ML
S-10965M1-1ML	8-Hydroxyquinoline Solution	148-24-3	1ML	S-11177A4-5ML	Benzoilprop ethyl Solution	22212-55-1	5ML
S-10970M1-1ML	9,10-Anthraquinone Solution	84-65-1	1ML	S-11177U1-1ML	Benzoilprop ethyl Solution	22212-55-1	1ML
S-10979A1-1ML	a-Endosulfan Solution	959-98-8	1ML	S-11182A1-1ML	Benzoilprop ethyl Solution	22212-55-1	1ML
S-10979J1-1ML	a-Endosulfan Solution	959-98-8	1ML	S-11182J5-1ML	Benzoilprop ethyl Solution	22212-55-1	1ML
S-10979J1-5ML	a-Endosulfan Solution	959-98-8	5ML	S-11182J7-1ML	Benzoilprop ethyl Solution	22212-55-1	1ML
S-10979T1-1ML	a-Endosulfan Solution	959-98-8	1ML	S-11182J7-5ML	Benzoilprop ethyl Solution	22212-55-1	5ML
S-10979T1-5ML	a-Endosulfan Solution	959-98-8	5ML	S-11191A1-1ML	Beta-Cyfluthrin Solution	68359-37-5	1ML
S-10983A1-1ML	a-Naphthylthiourea Solution	86-88-4	1ML	S-11191T1-1ML	Beta-Cyfluthrin Solution	68359-37-5	1ML
S-10983U1-1ML	a-Naphthylthiourea Solution	86-88-4	1ML	S-11192A1-1ML	Beta-Cypermethrin Solution	65731-84-2	1ML
S-10992A1-1ML	a,b-Endosulfan Solution	115-29-7	1ML	S-11192T1-1ML	Beta-Cypermethrin Solution	65731-84-2	1ML
S-10992J1-1ML	a,b-Endosulfan Solution	115-29-7	1ML	S-11194U1-1ML	BHC (alpha isomer) Solution	319-84-6	1ML
S-10992J1-5ML	a,b-Endosulfan Solution	115-29-7	5ML	S-11195U1-1ML	BHC (beta isomer) Solution	319-85-7	1ML
S-10995A1-1ML	Abamectin Solution	71751-41-2	1ML	S-11195U1-5ML	BHC (beta isomer) Solution	319-85-7	5ML
S-10996A1-1ML	Temephos Solution	3383-96-8	1ML	S-11196U1-1ML	BHC (delta isomer) Solution	319-86-8	1ML
S-10996U1-1ML	Temephos Solution	3383-96-8	1ML	S-11196U1-5ML	BHC (delta isomer) Solution	319-86-8	5ML
S-10997A1-1ML	Bromopropylate Solution	18181-80-1	1ML	S-11198A1-1ML	BHC (mixed isomers) Solution	608-73-1	1ML
S-11002A1-1ML	Acephate Solution	30560-19-1	1ML	S-11198U1-1ML	BHC (mixed isomers) Solution	608-73-1	1ML
S-11002B4-1ML	Acephate Solution	30560-19-1	1ML	S-11201A1-1ML	Bifenazate Solution	149877-41-8	1ML
S-11002B4-5ML	Acephate Solution	30560-19-1	5ML	S-11201U1-1ML	Bifenazate Solution	149877-41-8	1ML
S-11003A1-1ML	Acequinocyl Solution	57960-19-7	1ML	S-11202M1-1ML	Bifenox Solution	42576-02-3	1ML
S-11003U1-1ML	Acequinocyl Solution	57960-19-7	1ML	S-11203A1-1ML	Bifenthrin Solution	82657-04-3	1ML
S-11009A1-1ML	Acetamiprid Solution	135410-20-7	1ML	S-11203T1-1ML	Bifenthrin Solution	82657-04-3	1ML
S-11009U1-1ML	Acetamiprid Solution	135410-20-7	1ML	S-11204A1-1ML	Binapacryl Solution	485-31-4	1ML
S-11013M1-1ML	Acetochlor Solution	34256-82-1	1ML	S-11204T1-1ML	Binapacryl Solution	485-31-4	1ML
S-11026A1-1ML	Acibenzolar-S-methyl Solution	135158-54-2	1ML	S-11205A1-1ML	Bioallethrin Solution	28057-48-9	1ML
S-11026U1-1ML	Acibenzolar-S-methyl Solution	135158-54-2	1ML	S-11205T1-1ML	Bioallethrin Solution	28057-48-9	1ML
S-11027A1-1ML	Acifluorfen Solution	50594-66-6	1ML	S-11206A1-1ML	Bioresmethrin Solution	28434-01-7	1ML
S-11027B1-1ML	Acifluorfen Solution	50594-66-6	1ML	S-11206T1-1ML	Bioresmethrin Solution	28434-01-7	1ML
S-11027B1-5ML	Acifluorfen Solution	50594-66-6	5ML	S-11207J1-1ML	Biphenyl Solution	92-52-4	1ML
S-11028Y1-1ML	Acifluorfen methyl ester Solution	50594-67-7	1ML	S-11207J1-5ML	Biphenyl Solution	92-52-4	5ML
S-11028Y1-5ML	Acifluorfen methyl ester Solution	50594-67-7	5ML	S-11207M4-1ML	Biphenyl Solution	92-52-4	1ML
S-11042A1-1ML	Akton (TM) Solution	1757-18-2	1ML	S-11207M4-5ML	Biphenyl Solution	92-52-4	5ML
S-11042U1-1ML	Akton (TM) Solution	1757-18-2	1ML	S-11324F1-1ML	Bitrex (TM) Solution	3734-33-6	1ML
S-11043M1-1ML	Alachlor Solution	15972-60-8	1ML	S-11324U1-1ML	Bitrex (TM) Solution	3734-33-6	1ML
S-11043M1-5ML	Alachlor Solution	15972-60-8	5ML	S-11327A1-1ML	Bladex Solution	21725-46-2	1ML
S-11044M1-1ML	Aldicarb Solution	116-06-3	1ML	S-11327T4-1ML	Bladex Solution	21725-46-2	1ML
S-11045M1-1ML	Aldicarb sulfone Solution	1646-88-4	1ML	S-11327T4-5ML	Bladex Solution	21725-46-2	5ML
S-11046A1-1ML	Aldicarb sulfoxide Solution	1646-87-3	1ML	S-11328A1-1ML	Boscalid Solution	188425-85-6	1ML
S-11049M1-1ML	Aldrin Solution	309-00-2	1ML	S-11329A1-1ML	Brodifacoum Solution	56073-10-0	1ML
S-11049M1-5ML	Aldrin Solution	309-00-2	5ML	S-11330A1-1ML	Bromacil Solution	314-40-9	1ML
S-11051X1-1ML	Allethrin Solution	584-79-2	1ML	S-11330T1-1ML	Bromacil Solution	314-40-9	1ML
S-11051X1-5ML	Allethrin Solution	584-79-2	5ML	S-11333A1-1ML	Bromoacetyl bromide Solution	598-21-0	1ML
S-11059F1-1ML	Allyl trimethylammonium bromide Solution	3004-51-1	1ML	S-11333U1-1ML	Bromoacetyl bromide Solution	598-21-0	1ML
S-11059U1-1ML	Allyl trimethylammonium bromide Solution	3004-51-1	1ML	S-11342J1-1ML	Bromophos ethyl Solution	4824-78-6	1ML
S-11061A1-1ML	alpha-Cypermethrin Solution	67375-30-8	1ML	S-11343A1-1ML	Bromophos methyl Solution	2104-96-3	1ML
S-11061J1-1ML	alpha-Cypermethrin Solution	67375-30-8	1ML	S-11343J1-1ML	Bromophos methyl Solution	2104-96-3	1ML
S-11063A1-1ML	Amdro Solution	67485-29-4	1ML	S-11345M1-1ML	Bromoxynil Solution	1689-84-5	1ML

Part #	Description	CAS #	Size	Part #	Description	CAS #	Size
S-11346A1-1ML	Bromoxynil octanoic acid ester Solution	1689-99-2	1ML	S-11520A1-1ML	Cybutryne Solution	28159-98-0	1ML
S-11346K4-1ML	Bromoxynil octanoic acid ester Solution	1689-99-2	1ML	S-11520U1-1ML	Cybutryne Solution	28159-98-0	1ML
S-11346K4-5ML	Bromoxynil octanoic acid ester Solution	1689-99-2	5ML	S-11521A1-1ML	Cyclanilide Solution	113136-77-9	1ML
S-11348M1-1ML	Bronopol Solution	52-51-7	1ML	S-11521U1-1ML	Cyclanilide Solution	113136-77-9	1ML
S-11350A1-1ML	Bupropfen Solution	69327-76-0	1ML	S-11522A1-1ML	Cycloate Solution	1134-23-2	1ML
S-11350U1-1ML	Bupropfen Solution	69327-76-0	1ML	S-11522T1-1ML	Cycloate Solution	1134-23-2	1ML
S-11351A1-1ML	Busan 40 Solution	51026-28-9	1ML	S-11522T1-5ML	Cycloate Solution	1134-23-2	5ML
S-11351F1-1ML	Busan 40 Solution	51026-28-9	1ML	S-11534A1-1ML	Cycloheximide Solution	66-81-9	1ML
S-11352A1-1ML	Busan 85 Solution	128-03-0	1ML	S-11534U1-1ML	Cycloheximide Solution	66-81-9	1ML
S-11352F1-1ML	Busan 85 Solution	128-03-0	1ML	S-11541A1-1ML	Cycloproprate Solution	54460-46-7	1ML
S-11352M1-1ML	Busan 85 Solution	128-03-0	1ML	S-11541T4-1ML	Cycloproprate Solution	54460-46-7	1ML
S-11353A1-1ML	Butachlor Solution	23184-66-9	1ML	S-11541T4-5ML	Cycloproprate Solution	54460-46-7	5ML
S-11353T1-1ML	Butachlor Solution	23184-66-9	1ML	S-11542A1-1ML	Cyclosulfamuron Solution	136849-15-5	1ML
S-11353T1-5ML	Butachlor Solution	23184-66-9	5ML	S-11542U1-1ML	Cyclosulfamuron Solution	136849-15-5	1ML
S-11355A1-1ML	Butralin Solution	33629-47-9	1ML	S-11543A1-1ML	Cyhalofop-butyl Solution	122008-85-9	1ML
S-11355K1-1ML	Butralin Solution	33629-47-9	1ML	S-11543U1-1ML	Cyhalofop-butyl Solution	122008-85-9	1ML
S-11383A1-1ML	Butylate Solution	2008-41-5	1ML	S-11544A1-1ML	Cymoxanil Solution	57966-95-7	1ML
S-11383T1-1ML	Butylate Solution	2008-41-5	1ML	S-11544U1-1ML	Cymoxanil Solution	57966-95-7	1ML
S-11383T1-5ML	Butylate Solution	2008-41-5	5ML	S-11545A1-1ML	Cypermethrin Solution	52315-07-8	1ML
S-11392A1-1ML	Cadusafos Solution	95465-99-9	1ML	S-11545U1-1ML	Cypermethrin Solution	52315-07-8	1ML
S-11392U1-1ML	Cadusafos Solution	95465-99-9	1ML	S-11546A1-1ML	Cyprazine Solution	22936-86-3	1ML
S-11398A1-1ML	Capsaicin, Natural Solution	404-86-4	1ML	S-11546B1-1ML	Cyprazine Solution	22936-86-3	1ML
S-11398U1-1ML	Capsaicin, Natural Solution	404-86-4	1ML	S-11547A1-1ML	Cyproconazole Solution	94361-06-5	1ML
S-11400U1-1ML	Captan Solution	133-06-2	1ML	S-11547E1-1ML	Cyproconazole Solution	94361-06-5	1ML
S-11400U1-5ML	Captan Solution	133-06-2	5ML	S-11548A1-1ML	Cyprodinil Solution	121552-61-2	1ML
S-11402A4-1ML	Carbaryl Solution	63-25-2	1ML	S-11548U1-1ML	Cyprodinil Solution	121552-61-2	1ML
S-11402A4-5ML	Carbaryl Solution	63-25-2	5ML	S-11550A1-1ML	Cyromazine Solution	66215-27-8	1ML
S-11404M1-1ML	Carbendazim Solution	10605-21-7	1ML	S-11550B1-1ML	Cyromazine Solution	66215-27-8	1ML
S-11404M1-5ML	Carbendazim Solution	10605-21-7	5ML	S-11551A1-1ML	Cythioate Solution	115-93-5	1ML
S-11405M1-1ML	Carbofuran Solution	1563-66-2	1ML	S-11551U1-1ML	Cythioate Solution	115-93-5	1ML
S-11408M1-1ML	Carbophenothion Solution	786-19-6	1ML	S-11554A4-1ML	d-(cis-trans)-Phenothrin Solution	26002-80-2	1ML
S-11409A1-1ML	Carbosulfan Solution	55285-14-8	1ML	S-11554U1-1ML	d-(cis-trans)-Phenothrin Solution	26002-80-2	1ML
S-11410A1-1ML	Carboxin Solution	5234-68-4	1ML	S-11561A1-1ML	d-trans-Cyphenothrin Solution	39515-40-7	1ML
S-11411A1-1ML	Carfentrazone-ethyl Solution	128639-02-1	1ML	S-11561T1-1ML	d-trans-Cyphenothrin Solution	39515-40-7	1ML
S-11411U1-1ML	Carfentrazone-ethyl Solution	128639-02-1	1ML	S-11562B1-1ML	Dalapon Solution	75-99-0	1ML
S-11418A1-1ML	Chloramben Solution	133-90-4	1ML	S-11562M1-1ML	Dalapon Solution	75-99-0	1ML
S-11418B1-1ML	Chloramben Solution	133-90-4	1ML	S-11563M1-1ML	Dalapon methyl ester Solution	17640-02-7	1ML
S-11418B1-5ML	Chloramben Solution	133-90-4	5ML	S-11565A1-1ML	Daminozide Solution	1596-84-5	1ML
S-11419Y1-1ML	Chloramben methyl ester Solution	7286-84-2	1ML	S-11565T1-1ML	Daminozide Solution	1596-84-5	1ML
S-11419Y1-5ML	Chloramben methyl ester Solution	7286-84-2	5ML	S-11566M1-1ML	Dazomet Solution	533-74-4	1ML
S-11420M1-1ML	Chloranil Solution	118-75-2	1ML	S-11567M1-1ML	DDT (mixture p,p' & o,p') Solution	8017-34-3	1ML
S-11420U1-1ML	Chloranil Solution	118-75-2	1ML	S-11579A1-1ML	Deltamethrin Solution	52918-63-5	1ML
S-11423A1-1ML	Chlorbenseid Solution	103-17-3	1ML	S-11579T1-1ML	Deltamethrin Solution	52918-63-5	1ML
S-11423U1-1ML	Chlorbenseid Solution	103-17-3	1ML	S-11580U1-1ML	Demeton O Solution	298-03-3	1ML
S-11424A1-1ML	Chlorbromuron Solution	13360-45-7	1ML	S-11581U1-1ML	Demeton O&S Solution	8065-48-3	1ML
S-11424U1-1ML	Chlorbromuron Solution	13360-45-7	1ML	S-11582A1-1ML	Demeton S Solution	126-75-0	1ML
S-11425J4-1ML	Chlordane Solution	57-74-9	1ML	S-11582U1-1ML	Demeton S Solution	126-75-0	1ML
S-11425J4-5ML	Chlordane Solution	57-74-9	5ML	S-11584M1-1ML	Desmedipham Solution	13684-56-5	1ML
S-11425M1-1ML	Chlordane Solution	57-74-9	1ML	S-11585A1-1ML	Napropamide Solution	15299-99-7	1ML
S-11426A1-1ML	Chlordimeform Solution	6164-98-3	1ML	S-11585T1-1ML	Napropamide Solution	15299-99-7	1ML
S-11426U1-1ML	Chlordimeform Solution	6164-98-3	1ML	S-11585T1-5ML	Napropamide Solution	15299-99-7	5ML
S-11427A1-1ML	Chlorethoxyfos Solution	54593-83-8	1ML	S-11586M1-1ML	Fenaminosulf Solution	140-56-7	1ML
S-11427U1-1ML	Chlorethoxyfos Solution	54593-83-8	1ML	S-11587M1-1ML	Di-allate Solution	2303-16-4	1ML
S-11428A1-1ML	Chlorfenapyr Solution	122453-73-0	1ML	S-11587M1-5ML	Di-allate Solution	2303-16-4	5ML
S-11428U1-1ML	Chlorfenapyr Solution	122453-73-0	1ML	S-11589A1-1ML	Di-n-butyl phthalate Solution	84-74-2	1ML
S-11429K1-1ML	Chlorfenvinphos Solution	470-90-6	1ML	S-11589J1-1ML	Di-n-butyl phthalate Solution	84-74-2	1ML
S-11430A1-1ML	Chlorflazuron Solution	71422-67-8	1ML	S-11589J1-5ML	Di-n-butyl phthalate Solution	84-74-2	5ML
S-11431M1-1ML	Chlorfluorecol-methyl ester Solution	2536-31-4	1ML	S-11612A1-1ML	Diafenthion Solution	80060-09-9	1ML
S-11432A1-1ML	Chlorimuron ethyl Solution	90982-32-4	1ML	S-11612U1-1ML	Diafenthion Solution	80060-09-9	1ML
S-11434A1-1ML	Chlorimphos Solution	24934-91-6	1ML	S-11621U1-1ML	Diazinon Solution	333-41-5	1ML
S-11434U1-1ML	Chlorimphos Solution	24934-91-6	1ML	S-11644A1-1ML	Dibutyl succinate Solution	141-03-7	1ML
S-11435F1-1ML	Chlorimequat chloride Solution	999-81-5	1ML	S-11644U1-1ML	Dibutyl succinate Solution	141-03-7	1ML
S-11438T1-1ML	Chloroacetic acid Solution	79-11-8	1ML	S-11656A1-1ML	Dicamba Solution	1918-00-9	1ML
S-11442A1-1ML	Chlorobenzilate Solution	510-15-6	1ML	S-11656B1-1ML	Dicamba Solution	1918-00-9	1ML
S-11442J1-1ML	Chlorobenzilate Solution	510-15-6	1ML	S-11656B1-5ML	Dicamba Solution	1918-00-9	5ML
S-11442J1-5ML	Chlorobenzilate Solution	510-15-6	5ML	S-11657M1-1ML	Dicamba methyl ester Solution	6597-78-0	1ML
S-11449A1-1ML	Chloroneb Solution	2675-77-6	1ML	S-11657M1-5ML	Dicamba methyl ester Solution	6597-78-0	1ML
S-11449T1-1ML	Chloroneb Solution	2675-77-6	1ML	S-11658A1-1ML	Dicaphthos Solution	2463-84-5	1ML
S-11449T1-5ML	Chloroneb Solution	2675-77-6	5ML	S-11658U1-1ML	Dicaphthos Solution	2463-84-5	1ML
S-11451A1-1ML	Chlorophacinone Solution	3691-35-8	1ML	S-11659A1-1ML	Dichlofenthion Solution	97-17-6	1ML
S-11452B1-1ML	Trichloronitromethane Solution	76-06-2	1ML	S-11659J4-1ML	Dichlofenthion Solution	97-17-6	1ML
S-11452B1-5ML	Trichloronitromethane Solution	76-06-2	5ML	S-11659J4-5ML	Dichlofenthion Solution	97-17-6	5ML
S-11454A1-1ML	Chlorothalonil Solution	1897-45-6	1ML	S-11660U1-1ML	Dichlofluandil Solution	1085-98-9	1ML
S-11454T1-1ML	Chlorothalonil Solution	1897-45-6	1ML	S-11661M1-1ML	Dichlone Solution	117-80-6	1ML
S-11455A1-1ML	Chlorotoluron Solution	15545-48-9	1ML	S-11661U1-1ML	Dichlone Solution	117-80-6	1ML
S-11455U1-1ML	Chlorotoluron Solution	15545-48-9	1ML	S-11661U1-5ML	Dichlone Solution	117-80-6	5ML
S-11456A1-1ML	Chloroxuron Solution	1982-47-4	1ML	S-11662A1-1ML	Dichlormid Solution	37764-25-3	1ML
S-11457M1-1ML	Chloroxynil Solution	1891-95-8	1ML	S-11662U1-1ML	Dichlormid Solution	37764-25-3	1ML
S-11458A1-1ML	Chlorpropham Solution	101-21-3	1ML	S-11667M4-1ML	Dichlorophen Solution	97-23-4	1ML
S-11458T1-1ML	Chlorpropham Solution	101-21-3	1ML	S-11667M4-5ML	Dichlorophen Solution	97-23-4	5ML
S-11458T1-5ML	Chlorpropham Solution	101-21-3	5ML	S-11671B1-1ML	Dichlorprop Solution	120-36-5	1ML
S-11459A1-1ML	Chlorpyrifos Solution	2921-88-2	1ML	S-11671B1-5ML	Dichlorprop Solution	120-36-5	5ML
S-11459K1-1ML	Chlorpyrifos Solution	2921-88-2	1ML	S-11671M1-1ML	Dichlorprop Solution	120-36-5	1ML
S-11459K1-5ML	Chlorpyrifos Solution	2921-88-2	5ML	S-11672M1-1ML	Dichlorprop methyl ester Solution	57153-17-0	1ML
S-11460A1-1ML	Chlorpyrifos Methyl Solution	5598-13-0	1ML	S-11672M1-5ML	Dichlorprop methyl ester Solution	57153-17-0	5ML
S-11460J4-1ML	Chlorpyrifos Methyl Solution	5598-13-0	1ML	S-11674A1-1ML	Dichlorprop-P Solution	15165-67-0	1ML
S-11462B1-1ML	Chlorthal-dimethyl Solution	1861-32-1	1ML	S-11674U1-1ML	Dichlorprop-P Solution	15165-67-0	1ML
S-11462B1-5ML	Chlorthal-dimethyl Solution	1861-32-1	5ML	S-11675J4-1ML	Dichlorvos Solution	62-73-7	1ML
S-11462M1-1ML	Chlorthal-dimethyl Solution	1861-32-1	1ML	S-11675U1-1ML	Dichlorvos Solution	62-73-7	1ML
S-11463A1-1ML	Cholecalciferol Solution	67-97-0	1ML	S-11677A1-1ML	Diclofop methyl Solution	51338-27-3	1ML
S-11463U1-1ML	Cholecalciferol Solution	67-97-0	1ML	S-11677T1-1ML	Diclofop methyl Solution	51338-27-3	1ML
S-11468A1-1ML	Cinnamaldehyde Solution	104-55-2	1ML	S-11678A1-1ML	Dicloran Solution	99-30-9	1ML
S-11468U1-1ML	Cinnamaldehyde Solution	104-55-2	1ML	S-11678U4-1ML	Dicloran Solution	99-30-9	1ML
S-11480M1-1ML	cis-Chlordane Solution	5103-71-9	1ML	S-11679A1-1ML	Diclosulam Solution	145701-21-9	1ML
S-11480M1-5ML	cis-Chlordane Solution	5103-71-9	5ML	S-11680U1-1ML	Dicrotophos Solution	141-66-2	1ML
S-11481A1-1ML	E-Mevinphos Solution	26718-65-0	1ML	S-11680U1-5ML	Dicrotophos Solution	141-66-2	5ML
S-11481T1-1ML	E-Mevinphos Solution	26718-65-0	1ML	S-11686K1-1ML	Dicyclopentadiene Solution	77-73-6	1ML
S-11483A1-1ML	cis-Permethrin Solution	61949-76-6	1ML	S-11688M1-1ML	Dieldrin Solution	60-57-1	1ML
S-11483T1-1ML	cis-Permethrin Solution	61949-76-6	1ML	S-11688M1-5ML	Dieldrin Solution	60-57-1	5ML
S-11489A1-1ML	Clofinafop-propargyl Solution	105512-06-9	1ML	S-11690M1-1ML	Diethyl ethyl Solution	38727-55-8	1ML
S-11489U1-1ML	Clofinafop-propargyl Solution	105512-06-9	1ML	S-11691M1-1ML	Diethofencarb Solution	87130-20-9	1ML
S-11490M1-1ML	Clofentazine Solution	74115-24-5	1ML	S-11696A1-1ML	Diethyl dithiobis(thionoformate) Solution	502-55-6	1ML
S-11492A1-1ML	Cloransulam-methyl Solution	147150-35-4	1ML	S-11696K1-1ML	Diethyl dithiobis(thionoformate) Solution	502-55-6	1ML
S-11492U1-1ML	Cloransulam-methyl Solution	147150-35-4	1ML	S-11719M1-1ML	Difenacoum Solution	56073-07-5	1ML
S-11493A1-1ML	Clothianidin Solution	210880-92-5	1ML	S-11720M1-1ML	Difenoconazole Solution	119446-68-3	1ML
S-11495M1-1ML	Command Solution	81777-89-1	1ML	S-11721F1-1ML	Difenoquat methylsulfate Solution	43222-48-6	1ML
S-11507A1-1ML	Coumaphos Solution	56-72-4	1ML	S-11721U1-1ML	Difenoquat methylsulfate Solution	43222-48-6	1ML
S-11507U1-1ML	Coumaphos Solution	56-72-4	1ML	S-11722A1-1ML	Diffubenzuron Solution	35367-38-5	1ML
S-11509M1-1ML	Coumatetralyl Solution	5836-29-3	5ML	S-11750T1-1ML	Dikegulac acid Solution	18467-77-1	1ML
S-11513K1-1ML	Crotoxyphos Solution	7700-17-6	1ML	S-11752F1-1ML	Dimehypo Solution	52207-48-4	1ML
S-11513K1-5ML	Crotoxyphos Solution	7700-17-6	5ML	S-11753M1-1ML	Dimethachlon Solution	24096-53-5	1ML
S-11517A1-1ML	Cyanophos Solution	2636-26-2	1ML	S-11755A1-1ML	Dimethenamid Solution	87674-68-8	1ML
S-11517U1-1ML	Cyanophos Solution	2636-26-2	1ML	S-11755U1-1ML	Dimethenamid Solution	87674-68-8	1ML
				S-11756A1-1ML	Dimethenamid-P Solution	163515-14-8	1ML

Part #	Description	CAS #	Size	Part #	Description	CAS #	Size
S-11756U1-1ML	Dimethenamid-P Solution	163515-14-8	1ML	S-11958U1-1ML	Fenoxaprop-P-ethyl Solution	71283-80-2	1ML
S-11757M1-1ML	Dimethipin Solution	55290-64-7	1ML	S-11959A1-1ML	Fenoxycarb Solution	72490-01-8	1ML
S-11758J1-1ML	Dimethoate Solution	60-51-5	1ML	S-11960A1-1ML	Fenpropathrin Solution	39515-41-8	1ML
S-11758J1-5ML	Dimethoate Solution	60-51-5	5ML	S-11960T1-1ML	Fenpropathrin Solution	39515-41-8	1ML
S-11759A1-1ML	Dimethomorph Solution	110488-70-5	1ML	S-11962A1-1ML	Fenpyroximate Solution	134098-61-6	1ML
S-11770H1-1ML	Dimethyl phthalate Solution	131-11-3	1ML	S-11962U1-1ML	Fenpyroximate Solution	134098-61-6	1ML
S-11770H1-5ML	Dimethyl phthalate Solution	131-11-3	5ML	S-11963M1-1ML	Fensulfothion Solution	115-90-2	1ML
S-11770M1-1ML	Dimethyl phthalate Solution	131-11-3	1ML	S-11963M1-5ML	Fensulfothion Solution	115-90-2	5ML
S-11770M1-5ML	Dimethyl phthalate Solution	131-11-3	5ML	S-11964S1-1ML	Fenthion Solution	55-38-9	1ML
S-11775A1-1ML	Dimethyl-p-nitrophenylphosphate Solution	950-35-6	1ML	S-11964S1-5ML	Fenthion Solution	55-38-9	5ML
S-11775T1-1ML	Dimethyl-p-nitrophenylphosphate Solution	950-35-6	1ML	S-11965A1-1ML	Fenthion-ethyl Solution	1716-09-2	1ML
S-11775T1-5ML	Dimethyl-p-nitrophenylphosphate Solution	950-35-6	5ML	S-11965U1-1ML	Fenthion-ethyl Solution	1716-09-2	1ML
S-11779F4-1ML	Dimethylarsinic acid Solution	75-60-5	1ML	S-11966A1-1ML	Fentin hydroxide Solution	76-87-9	1ML
S-11782M1-1ML	Diniconazole Solution	83657-24-3	1ML	S-11967A4-1ML	Fenuron Solution	101-42-8	1ML
S-11783M1-1ML	Dinocap Solution	39300-45-3	1ML	S-11967A4-5ML	Fenuron Solution	101-42-8	5ML
S-11783M1-5ML	Dinocap Solution	39300-45-3	5ML	S-11967U1-1ML	Fenuron Solution	101-42-8	1ML
S-11786A1-1ML	Dinoseb Solution	88-85-7	1ML	S-11970M1-1ML	Ferbam Solution	14484-64-1	1ML
S-11786B1-1ML	Dinoseb Solution	88-85-7	1ML	S-11970U1-1ML	Ferbam Solution	14484-64-1	1ML
S-11787M1-1ML	Dinoseb methyl ether Solution	6099-79-2	1ML	S-11971A1-1ML	Fipronil Solution	120068-37-3	1ML
S-11787M1-5ML	Dinoseb methyl ether Solution	6099-79-2	5ML	S-11971U1-1ML	Fipronil Solution	120068-37-3	1ML
S-11788M1-1ML	Dinotefuran Solution	165252-70-0	1ML	S-11972A1-1ML	Flamprop isopropyl Solution	52756-22-6	1ML
S-11788U1-1ML	Dinotefuran Solution	165252-70-0	1ML	S-11972T1-1ML	Flamprop isopropyl Solution	52756-22-6	1ML
S-11789M4-1ML	Dioxacarb Solution	6988-21-2	1ML	S-11973A1-1ML	Flamprop-methyl Solution	52756-25-9	1ML
S-11790B1-1ML	Dioxathion Solution	78-34-2	1ML	S-11973T1-1ML	Flamprop-methyl Solution	52756-25-9	1ML
S-11791A1-1ML	Dioxydemeton-S-methyl Solution	17040-19-6	1ML	S-11975A1-1ML	Floucoumafen Solution	90035-08-8	1ML
S-11793A1-1ML	Diphacinone Solution	82-66-6	1ML	S-11976A1-1ML	Flonicamid Solution	158062-67-0	1ML
S-11794A1-1ML	Diphenamid Solution	957-51-7	1ML	S-11976U1-1ML	Flonicamid Solution	158062-67-0	1ML
S-11794T1-1ML	Diphenamid Solution	957-51-7	1ML	S-11977A1-1ML	Fluazifop Solution	69335-91-7	1ML
S-11794T1-5ML	Diphenamid Solution	957-51-7	5ML	S-11977A1-5ML	Fluazifop Solution	69335-91-7	5ML
S-11800M1-1ML	Diphenylacetoneitrile Solution	86-29-3	1ML	S-11978A1-1ML	Fluazifop-butyl Solution	69806-50-4	1ML
S-11808A1-1ML	Dipropetryn Solution	4147-51-7	1ML	S-11978U1-1ML	Fluazifop-butyl Solution	69806-50-4	1ML
S-11808T1-1ML	Dipropetryn Solution	4147-51-7	1ML	S-11979A1-1ML	Fluazifop-P-butyl Solution	79241-46-6	1ML
S-11816F1-1ML	Diquat Dibromide monohydrate Solution	6385-62-2	1ML	S-11979U1-1ML	Fluazifop-P-butyl Solution	79241-46-6	1ML
S-11816F13-1ML	Diquat Dibromide monohydrate Solution	6385-62-2	1ML	S-11980A1-1ML	Fluazinam Solution	79622-59-6	1ML
S-11816F13-5ML	Diquat Dibromide monohydrate Solution	6385-62-2	5ML	S-11980U1-1ML	Fluazinam Solution	79622-59-6	1ML
S-11817F1-1ML	Sodium methyl arsonate hexahydrate Solution	144-21-8	1ML	S-11982F1-1ML	Flucarbazone-sodium Solution	181274-17-9	1ML
S-11819M1-1ML	Disulfoton Solution	298-04-4	1ML	S-11983A1-1ML	Fluchloralin Solution	33245-39-5	1ML
S-11820A1-1ML	Ditalimfos Solution	5131-24-8	1ML	S-11983J4-1ML	Fluchloralin Solution	33245-39-5	1ML
S-11820J1-1ML	Ditalimfos Solution	5131-24-8	1ML	S-11983J4-5ML	Fluchloralin Solution	33245-39-5	5ML
S-11820J1-5ML	Ditalimfos Solution	5131-24-8	5ML	S-11984A1-1ML	Flucythrinate Solution	70124-77-5	1ML
S-11820U1-1ML	Ditalimfos Solution	5131-24-8	1ML	S-11984U1-1ML	Flucythrinate Solution	70124-77-5	1ML
S-11822A1-1ML	Dithiopyr Solution	97886-45-8	1ML	S-11985A1-1ML	Fludioxonil Solution	131341-86-1	1ML
S-11822U1-1ML	Dithiopyr Solution	97886-45-8	1ML	S-11985U1-1ML	Fludioxonil Solution	131341-86-1	1ML
S-11827A4-1ML	Diuron Solution	330-54-1	1ML	S-11986M1-1ML	Flufenacet Solution	142459-58-3	1ML
S-11827U1-1ML	Diuron Solution	330-54-1	1ML	S-11987A1-1ML	Flufenoxuron Solution	101463-69-8	1ML
S-11839M1-1ML	Dodermorph Solution	1593-77-7	1ML	S-11988A1-1ML	Flumetralin Solution	62924-70-3	1ML
S-11842A1-1ML	Fonofos Solution	944-22-9	1ML	S-11988B1-1ML	Flumetralin Solution	62924-70-3	1ML
S-11842J4-1ML	Fonofos Solution	944-22-9	1ML	S-11988B1-5ML	Flumetralin Solution	62924-70-3	5ML
S-11842J4-5ML	Fonofos Solution	944-22-9	5ML	S-11988K1-1ML	Flumetralin Solution	62924-70-3	1ML
S-11844A1-1ML	Edifenphos Solution	17109-49-8	1ML	S-11989A1-1ML	Flumetsulam Solution	98967-40-9	1ML
S-11844U1-1ML	Edifenphos Solution	17109-49-8	1ML	S-11989U1-1ML	Flumetsulam Solution	98967-40-9	1ML
S-11847A1-1ML	Emamectin Solution	119791-41-2	1ML	S-11990A1-1ML	Flumiclorac-pentyl Solution	87546-18-7	1ML
S-11851A1-1ML	Endosulfan sulfate Solution	1031-07-8	1ML	S-11990T1-1ML	Flumiclorac-pentyl Solution	87546-18-7	1ML
S-11851J1-1ML	Endosulfan sulfate Solution	1031-07-8	1ML	S-11991A1-1ML	Flumioxazin Solution	103361-09-7	1ML
S-11851J1-5ML	Endosulfan sulfate Solution	1031-07-8	5ML	S-11991U1-1ML	Flumioxazin Solution	103361-09-7	1ML
S-11852A1-1ML	Endothal monohydrate Solution	62059-43-2	1ML	S-11992A4-1ML	Fluometuron Solution	2164-17-2	1ML
S-11852A1-5ML	Endothal monohydrate Solution	62059-43-2	5ML	S-11998M1-1ML	Fluorodifen Solution	15457-05-3	1ML
S-11852F0-1ML	Endothal monohydrate Solution	62059-43-2	1ML	S-12000M1-1ML	Fluroxypyr Solution	69377-81-7	1ML
S-11852F10-1ML	Endothal monohydrate Solution	62059-43-2	1ML	S-12001M1-1ML	Fluroxypyr-methyl Solution	81406-37-3	1ML
S-11852F10-5ML	Endothal monohydrate Solution	62059-43-2	5ML	S-12002M1-1ML	Flurprimidol Solution	56425-91-3	1ML
S-11852U1-1ML	Endothal monohydrate Solution	62059-43-2	1ML	S-12003A1-1ML	Flusilazole Solution	85509-19-9	1ML
S-11853T1-1ML	Endothal-PFPH Solution		1ML	S-12003U1-1ML	Flusilazole Solution	85509-19-9	1ML
S-11853T1-5ML	Endothal-PFPH Solution		5ML	S-12004M1-1ML	Flutolanil Solution	66332-96-5	1ML
S-11854M1-1ML	Endrin Solution	72-20-8	1ML	S-12005A1-1ML	Flutriafol Solution	76674-21-0	1ML
S-11854M1-5ML	Endrin Solution	72-20-8	5ML	S-12005U1-1ML	Flutriafol Solution	76674-21-0	1ML
S-11855M1-1ML	Endrin aldehyde Solution	7421-93-4	1ML	S-12006A1-1ML	Tebuconazole Solution	107534-96-3	1ML
S-11855M1-5ML	Endrin aldehyde Solution	7421-93-4	5ML	S-12006K1-1ML	Tebuconazole Solution	107534-96-3	1ML
S-11856A1-1ML	Endrin ketone Solution	53494-70-5	1ML	S-12007M1-1ML	Folpet Solution	133-07-3	1ML
S-11856U1-1ML	Endrin ketone Solution	53494-70-5	1ML	S-12007U1-1ML	Folpet Solution	133-07-3	1ML
S-11858A1-1ML	EPN Solution	2104-64-5	1ML	S-12008M1-1ML	Fomesafen Solution	72178-02-0	1ML
S-11858J1-1ML	EPN Solution	2104-64-5	1ML	S-12009A1-1ML	Foramsulfuron Solution	173159-57-4	1ML
S-11858J1-5ML	EPN Solution	2104-64-5	5ML	S-12010A1-1ML	Forchlorfenuron Solution	68157-60-8	1ML
S-11858J4-1ML	EPN Solution	2104-64-5	1ML	S-12015A1-1ML	Formetanate hydrochloride Solution	23422-53-9	1ML
S-11858J4-5ML	EPN Solution	2104-64-5	5ML	S-12017A1-1ML	Formothion Solution	2540-82-1	1ML
S-11865A1-1ML	Ethalfuralin Solution	55283-68-6	1ML	S-12017J1-1ML	Formothion Solution	2540-82-1	1ML
S-11865J4-1ML	Ethalfuralin Solution	55283-68-6	1ML	S-12018F1-1ML	Fosamine Ammonium Solution	25954-13-6	1ML
S-11865J4-5ML	Ethalfuralin Solution	55283-68-6	5ML	S-12020A1-1ML	Fosthiazate Solution	98886-44-3	1ML
S-11866A1-1ML	Ethametsulfuron-methyl Solution	97780-06-8	1ML	S-12020U1-1ML	Fosthiazate Solution	98886-44-3	1ML
S-11867M1-1ML	Ethanedial dioxime Solution	557-30-2	1ML	S-12021M1-1ML	Fuberidazole Solution	3878-19-1	1ML
S-11871A1-1ML	Ethidimuron Solution	30043-49-3	1ML	S-12108A1-1ML	Gibberellic acid Solution	77-06-5	1ML
S-11871U1-1ML	Ethidimuron Solution	30043-49-3	1ML	S-12108T1-1ML	Gibberellic acid Solution	77-06-5	1ML
S-11872A1-1ML	Ethiofencarb Solution	29973-13-5	1ML	S-12134F1-1ML	Glyphosate-isopropylammonium Solution	38641-94-0	1ML
S-11872U1-1ML	Ethiofencarb Solution	29973-13-5	1ML	S-12134U1-1ML	Glyphosate-isopropylammonium Solution	38641-94-0	1ML
S-11873M1-1ML	Ethion Solution	563-12-2	1ML	S-12135F1-1ML	Glyphosine Solution	2439-99-8	1ML
S-11873M1-5ML	Ethion Solution	563-12-2	5ML	S-12135U1-1ML	Glyphosine Solution	2439-99-8	1ML
S-11874A1-1ML	Ethofumesate Solution	26225-79-6	1ML	S-12139A1-1ML	Azinphos-methyl Solution	86-50-0	1ML
S-11874U1-1ML	Ethofumesate Solution	26225-79-6	1ML	S-12139U1-1ML	Azinphos-methyl Solution	86-50-0	1ML
S-11914A1-1ML	Ethyl-1-naphthalene acetate Solution	2122-70-5	1ML	S-12139U1-5ML	Azinphos-methyl Solution	86-50-0	5ML
S-11914T1-1ML	Ethyl-1-naphthalene acetate Solution	2122-70-5	1ML	S-12140A1-1ML	Guthion Ethyl Solution	2642-71-9	1ML
S-11916A1-1ML	Ethyl-3-amino-2,5-dichlorobenzoate Solution	60541-86-8	1ML	S-12141M1-1ML	Halofenozide Solution	112226-61-6	1ML
S-11916T1-1ML	Ethyl-3-amino-2,5-dichlorobenzoate Solution	60541-86-8	1ML	S-12142A1-1ML	Halosulfuron-methyl Solution	100784-20-1	1ML
S-11941A1-1ML	Etofenprox Solution	80844-07-1	1ML	S-12147M1-1ML	Heptachlor Solution	76-44-8	1ML
S-11941U1-1ML	Etofenprox Solution	80844-07-1	1ML	S-12147M1-5ML	Heptachlor Solution	76-44-8	5ML
S-11942A1-1ML	Etoxadone Solution	153233-91-1	1ML	S-12148M1-1ML	Heptachlor epoxide (Isomer B) Solution	1024-57-3	1ML
S-11942U1-1ML	Etoxadone Solution	153233-91-1	1ML	S-12148M1-5ML	Heptachlor epoxide (Isomer B) Solution	1024-57-3	5ML
S-11943A1-1ML	Famoxadone Solution	131807-57-3	1ML	S-12153A1-1ML	Heptenophos Solution	23560-59-0	1ML
S-11943U1-1ML	Famoxadone Solution	131807-57-3	1ML	S-12153U1-1ML	Heptenophos Solution	23560-59-0	1ML
S-11944K1-1ML	Famphur Solution	52-85-7	1ML	S-12158A1-1ML	Hexachloroacetone Solution	116-16-5	1ML
S-11945A1-1ML	Fenac Solution	85-34-7	1ML	S-12158J1-1ML	Hexachloroacetone Solution	116-16-5	1ML
S-11945T1-1ML	Fenac Solution	85-34-7	1ML	S-12159K4-1ML	Hexachlorobenzene Solution	118-74-1	1ML
S-11946A1-1ML	Fenamidone Solution	161326-34-7	1ML	S-12159M1-1ML	Hexachlorobenzene Solution	118-74-1	1ML
S-11946U1-1ML	Fenamidone Solution	161326-34-7	1ML	S-12159M1-5ML	Hexachlorobenzene Solution	118-74-1	5ML
S-11948A1-1ML	Fenarimol Solution	60168-88-9	1ML	S-12161A1-1ML	Hexachlorodimethyl sulfone Solution	3064-70-8	1ML
S-11948T1-1ML	Fenarimol Solution	60168-88-9	1ML	S-12161U1-1ML	Hexachlorodimethyl sulfone Solution	3064-70-8	1ML
S-11948T1-5ML	Fenarimol Solution	60168-88-9	5ML	S-12163M1-1ML	Hexachlorophene Solution	70-30-4	1ML
S-11949A1-1ML	Fenazaquin Solution	120928-09-8	1ML	S-12163M1-5ML	Hexachlorophene Solution	70-30-4	5ML
S-11950M1-1ML	Fenbuconazole Solution	114369-43-6	1ML	S-12165A1-1ML	Hexaconazole Solution	79983-71-4	1ML
S-11952U1-1ML	Fenchlorphos Solution	299-84-3	1ML	S-12165U1-1ML	Hexaconazole Solution	79983-71-4	1ML
S-11952U1-5ML	Fenchlorphos Solution	299-84-3	5ML	S-12167A1-1ML	Hexaflumuron Solution	86479-06-3	1ML
S-11954A1-1ML	Fenhexamid Solution	126833-17-8	1ML	S-12167U1-1ML	Hexaflumuron Solution	86479-06-3	1ML
S-11955A1-1ML	Fenitrothion Solution	122-14-5	1ML	S-12184A1-1ML	Hexythiazox Solution	78587-05-0	1ML
S-11955J4-1ML	Fenitrothion Solution	122-14-5	1ML	S-12196M1-1ML	Hymexazole Solution	10004-44-1	1ML
S-11957M1-1ML	Fenoxaprop ethyl Solution	66441-23-4	1ML	S-12197A1-1ML	Icaridin Solution	119515-38-7	1ML
S-11958A1-1ML	Fenoxaprop-P-ethyl Solution	71283-80-2	1ML	S-12198A1-1ML	Imazail Solution	35554-44-0	1ML

Part #	Description	CAS #	Size	Part #	Description	CAS #	Size
S-12198J1-1ML	Imazalil Solution	35554-44-0	1ML	S-12403M1-5ML	Methoxychlor Solution	72-43-5	5ML
S-12199F1-1ML	Imazalil sulfate Solution	58594-72-2	1ML	S-12406A1-1ML	Methoxyfenozide Solution	161050-58-4	1ML
S-12200A1-1ML	Imazamethabenz-methyl Solution	81405-85-8	1ML	S-12414M1-1ML	Methyl anthranilate Solution	134-20-3	1ML
S-12200U1-1ML	Imazamethabenz-methyl Solution	81405-85-8	1ML	S-12452A1-1ML	Methyl parathion Solution	298-00-0	1ML
S-12201A1-1ML	Imazamox Solution	114311-32-9	1ML	S-12452U1-1ML	Methyl parathion Solution	298-00-0	1ML
S-12201T1-1ML	Imazamox Solution	114311-32-9	1ML	S-12464M1-1ML	Methyl-1-naphthalene acetate Solution	2876-78-0	1ML
S-12202A1-1ML	Imazapic Solution	104098-48-8	1ML	S-12464U1-1ML	Methyl-1-naphthalene acetate Solution	2876-78-0	1ML
S-12203A1-1ML	Imazapyr Solution	81334-34-1	1ML	S-12477A1-1ML	Metobromuron Solution	3060-89-7	1ML
S-12204A1-1ML	Imazaquin Solution	81335-37-7	1ML	S-12477U1-1ML	Metobromuron Solution	3060-89-7	1ML
S-12204T1-1ML	Imazaquin Solution	81335-37-7	1ML	S-12478A1-1ML	Metolachlor Solution	51218-45-2	1ML
S-12205M1-1ML	Imazethapyr Solution	81335-77-5	1ML	S-12478T1-1ML	Metolachlor Solution	51218-45-2	1ML
S-12206M1-1ML	Imidacloprid Solution	138261-41-3	1ML	S-12478T1-5ML	Metolachlor Solution	51218-45-2	5ML
S-12207A1-1ML	Phosmet Solution	732-11-6	1ML	S-12479A1-1ML	Metolcarb Solution	1129-41-5	1ML
S-12207J1-1ML	Phosmet Solution	732-11-6	1ML	S-12479U1-1ML	Metolcarb Solution	1129-41-5	1ML
S-12207J1-5ML	Phosmet Solution	732-11-6	5ML	S-12480A1-1ML	Metoxuron Solution	19937-59-8	1ML
S-12213A1-1ML	S-Indoxacarb Solution	173584-44-6	1ML	S-12480U1-1ML	Metoxuron Solution	19937-59-8	1ML
S-12213U1-1ML	S-Indoxacarb Solution	173584-44-6	1ML	S-12481A1-1ML	Metribuzin Solution	21087-64-9	1ML
S-12218A1-1ML	Ioxynil octanoate Solution	3861-47-0	1ML	S-12481T1-1ML	Metribuzin Solution	21087-64-9	1ML
S-12218U1-1ML	Ioxynil octanoate Solution	3861-47-0	1ML	S-12481T1-5ML	Metribuzin Solution	21087-64-9	5ML
S-12219A1-1ML	Iprobenfos Solution	26087-47-8	1ML	S-12482A1-1ML	Metsulfuron methyl Solution	74223-64-6	1ML
S-12219U1-1ML	Iprobenfos Solution	26087-47-8	1ML	S-12483A1-1ML	MGK 264 (TM) Solution	113-48-4	1ML
S-12220A1-1ML	Iprodione Solution	36734-19-7	1ML	S-12483T1-1ML	MGK 264 (TM) Solution	113-48-4	1ML
S-12222U1-1ML	Isazophos Solution	42509-80-8	1ML	S-12483T1-5ML	MGK 264 (TM) Solution	113-48-4	5ML
S-12231M1-1ML	Isobenzan Solution	297-78-9	1ML	S-12486M1-1ML	Mirex Solution	2385-85-5	1ML
S-12232A1-1ML	Isobornyl thiocynoacetate Solution	115-31-1	1ML	S-12487A1-1ML	Molinate Solution	2212-67-1	1ML
S-12232U1-1ML	Isobornyl thiocynoacetate Solution	115-31-1	1ML	S-12487T1-1ML	Molinate Solution	2212-67-1	1ML
S-12248M1-1ML	Isodrin Solution	465-73-6	1ML	S-12487T1-5ML	Molinate Solution	2212-67-1	5ML
S-12248M1-5ML	Isodrin Solution	465-73-6	5ML	S-12488A1-1ML	Pencycuron Solution	66063-05-6	1ML
S-12249A1-1ML	Isofenphos Solution	25311-71-1	1ML	S-12488U1-1ML	Pencycuron Solution	66063-05-6	1ML
S-12249U1-1ML	Isofenphos Solution	25311-71-1	1ML	S-12493A1-1ML	Monocrotophos Solution	6923-22-4	1ML
S-12264A1-1ML	Isoprocarb Solution	2631-40-5	1ML	S-12493U1-1ML	Monocrotophos Solution	6923-22-4	1ML
S-12264U1-1ML	Isoprocarb Solution	2631-40-5	1ML	S-12493U1-5ML	Monocrotophos Solution	6923-22-4	5ML
S-12265A1-1ML	Isopropalin Solution	33820-53-0	1ML	S-12494A1-1ML	Monolinuron Solution	1746-81-2	1ML
S-12265J4-1ML	Isopropalin Solution	33820-53-0	1ML	S-12494U1-1ML	Monolinuron Solution	1746-81-2	1ML
S-12265J4-5ML	Isopropalin Solution	33820-53-0	5ML	S-12495F1-1ML	Monosodium acid methane arsonate sesquihydrate	2163-80-6	1ML
S-12278A1-1ML	Isothiolane Solution	50512-35-1	1ML	S-12496U1-1ML	Monothiono TEPP Solution	645-78-3	1ML
S-12278T1-1ML	Isothiolane Solution	50512-35-1	1ML	S-12497A4-1ML	Monuron Solution	150-68-5	1ML
S-12279A1-1ML	Isopturon Solution	34123-59-6	1ML	S-12497M4-1ML	Monuron Solution	150-68-5	1ML
S-12285M1-1ML	Isoxaben Solution	82558-50-7	1ML	S-12505A1-1ML	N-[(2,4-Dichlorophenoxy)acetyl]-DL-methionine Solution	62625-13-2	1ML
S-12287U1-1ML	Isoxalutole Solution	141112-29-0	1ML	S-12505T1-1ML	N-[(2,4-Dichlorophenoxy)acetyl]-DL-methionine Solution	62625-13-2	1ML
S-12289A1-1ML	Ivermectin Solution	70288-86-7	1ML	S-12506M1-1ML	N-1-Naphthylacetamide Solution	575-36-0	1ML
S-12289U1-1ML	Ivermectin Solution	70288-86-7	1ML	S-12507A1-1ML	N-1-Naphthylphthalamic acid Solution	132-66-1	1ML
S-12291A1-1ML	Chlordecone Solution	143-50-0	1ML	S-12507T1-1ML	N-1-Naphthylphthalamic acid Solution	132-66-1	1ML
S-12291J1-1ML	Chlordecone Solution	143-50-0	1ML	S-12557A1-1ML	N-m-Tolylphthalamic acid Solution	85-72-3	1ML
S-12291J1-5ML	Chlordecone Solution	143-50-0	5ML	S-12557T1-1ML	N-m-Tolylphthalamic acid Solution	85-72-3	1ML
S-12291T4-1ML	Chlordecone Solution	143-50-0	1ML	S-12558A1-1ML	N-Methyl-methamidophos Solution	28167-49-9	1ML
S-12291T4-5ML	Chlordecone Solution	143-50-0	5ML	S-12558U1-1ML	N-Methyl-methamidophos Solution	28167-49-9	1ML
S-12295M1-1ML	Kinoprene Solution	42588-37-4	1ML	S-12559M1-1ML	N-Methyl-N-1-naphthyl acetamide Solution	5903-13-9	1ML
S-12295T4-1ML	Kinoprene Solution	42588-37-4	1ML	S-12617M1-1ML	N,N-Diallyl-2-chloroacetamide Solution	93-71-0	1ML
S-12295T4-5ML	Kinoprene Solution	42588-37-4	5ML	S-12618A1-1ML	N,N-Diethyl-m-toluamide Solution	134-62-3	1ML
S-12296A1-1ML	Kresoxim-methyl Solution	143390-89-0	1ML	S-12618B5-1ML	N,N-Diethyl-m-toluamide Solution	134-62-3	1ML
S-12296U1-1ML	Kresoxim-methyl Solution	143390-89-0	1ML	S-12618B5-5ML	N,N-Diethyl-m-toluamide Solution	134-62-3	5ML
S-12304A1-1ML	Lactofen Solution	77501-63-4	1ML	S-12630A1-1ML	N,N-Dimethylsuccinamic acid Solution	2564-95-6	1ML
S-12304U1-1ML	Lactofen Solution	77501-63-4	1ML	S-12630T1-1ML	N,N-Dimethylsuccinamic acid Solution	2564-95-6	1ML
S-12307A1-1ML	lambda-Cyhalothrin Solution	91465-08-6	1ML	S-12639F1-1ML	Nabam Solution	142-59-6	1ML
S-12307J1-1ML	lambda-Cyhalothrin Solution	91465-08-6	1ML	S-12639U1-1ML	Nabam Solution	142-59-6	1ML
S-12314M1-1ML	Lenacil Solution	2164-08-100	1ML	S-12640U1-1ML	Naled Solution	300-76-5	1ML
S-12315A1-1ML	Leptophos Solution	21609-90-5	1ML	S-12647A4-1ML	Neburon Solution	555-37-3	1ML
S-12315U1-1ML	Leptophos Solution	21609-90-5	1ML	S-12647A4-5ML	Neburon Solution	555-37-3	5ML
S-12319A1-1ML	Lindane (BHC gamma isomer) Solution	58-89-9	1ML	S-12655M1-1ML	Nicotine Solution	54-11-5	1ML
S-12322A6-1ML	Linuron Solution	330-55-2	1ML	S-12657M1-1ML	Nitenpyram Solution	150824-47-8	1ML
S-12322A6-5ML	Linuron Solution	330-55-2	5ML	S-12658A1-1ML	Nitrin Solution	4726-14-1	1ML
S-12324A1-1ML	Clpyralid Solution	1702-17-6	1ML	S-12658J1-1ML	Nitrin Solution	4726-14-1	1ML
S-12324T1-1ML	Clpyralid Solution	1702-17-6	1ML	S-12659A1-1ML	Nitrapyrin Solution	1929-82-4	1ML
S-12325A1-1ML	Lufenuron Solution	103055-07-8	1ML	S-12663M1-1ML	Nitrofen Solution	1836-75-5	1ML
S-12346A1-1ML	Malathion Solution	121-75-5	1ML	S-12668A1-1ML	Norflurazon Solution	27314-13-2	1ML
S-12346U1-1ML	Malathion Solution	121-75-5	1ML	S-12668T1-1ML	Norflurazon Solution	27314-13-2	1ML
S-12346U1-5ML	Malathion Solution	121-75-5	5ML	S-12668T1-5ML	Norflurazon Solution	27314-13-2	5ML
S-12350M1-1ML	Maleic hydrazide Solution	123-33-1	1ML	S-12669A1-1ML	Noruron Solution	18530-56-8	1ML
S-12357A1-1ML	MCPB Solution	94-81-5	1ML	S-12669U1-1ML	Noruron Solution	18530-56-8	1ML
S-12357T1-1ML	MCPB Solution	94-81-5	1ML	S-12670A1-1ML	Novaluron Solution	116714-46-6	1ML
S-12358A1-1ML	MCPB methyl ester Solution	57153-18-1	1ML	S-12670U1-1ML	Novaluron Solution	116714-46-6	1ML
S-12358T1-1ML	MCPB methyl ester Solution	57153-18-1	1ML	S-12676M1-1ML	o-Benzyl-p-chlorophenol Solution	120-32-1	1ML
S-12359A1-1ML	Mecarbam Solution	2595-54-2	1ML	S-12679M1-1ML	o-Chlorophenoxy acetic acid Solution	614-61-9	1ML
S-12359U1-1ML	Mecarbam Solution	2595-54-2	1ML	S-12679T1-1ML	o-Chlorophenoxy acetic acid Solution	614-61-9	1ML
S-12360A1-1ML	Mecoprop Solution	7085-19-0	1ML	S-12703A1-1ML	O,O-Diethyl phosphorochloridothioate Solution	2524-04-100	1ML
S-12360B1-1ML	Mecoprop Solution	7085-19-0	1ML	S-12703U1-1ML	O,O-Diethyl phosphorochloridothioate Solution	2524-04-100	1ML
S-12361M1-1ML	Mecoprop methyl ester Solution	23844-56-6	1ML	S-12704A1-1ML	O,O-Dimethyl phosphorochloridothioate Solution	2524-03-0	1ML
S-12361M1-5ML	Mecoprop methyl ester Solution	23844-56-6	5ML	S-12704U1-1ML	O,O-Dimethyl phosphorochloridothioate Solution	2524-03-0	1ML
S-12362B1-1ML	Mecoprop-2,4,4-trimethylpentyl ester Solution	23844-56-6	1ML	S-12706M1-1ML	o,p'-DDD Solution	53-19-0	1ML
S-12362B1-5ML	Mecoprop-2,4,4-trimethylpentyl ester Solution	23844-56-6	5ML	S-12707M1-1ML	o,p'-DDD Solution	3424-82-6	1ML
S-12363A1-1ML	Mecoprop-P Solution	16484-77-8	1ML	S-12707M1-5ML	o,p'-DDE Solution	3424-82-6	5ML
S-12363T1-1ML	Mecoprop-P Solution	16484-77-8	1ML	S-12708M1-1ML	o,p'-DDT Solution	789-02-6	1ML
S-12364M1-1ML	Mefenacet Solution	73250-68-7	1ML	S-12708M1-5ML	o,p'-DDT Solution	789-02-6	5ML
S-12365A1-1ML	Mefenoxam Solution	70630-17-0	1ML	S-12720A1-1ML	Octyl aldehyde (DNPH Derivative) Solution		1ML
S-12365T1-1ML	Mefenoxam Solution	70630-17-0	1ML	S-12720A1-5ML	Octyl aldehyde (DNPH Derivative) Solution		5ML
S-12366A1-1ML	Mefenpyr-diethyl Solution	135590-91-9	1ML	S-12720W4-1ML	Octyl aldehyde (DNPH Derivative) Solution		1ML
S-12366U1-1ML	Mefenpyr-diethyl Solution	135590-91-9	1ML	S-12720W4-5ML	Octyl aldehyde (DNPH Derivative) Solution		5ML
S-12367M1-1ML	Meftuicide Solution	53780-34-0	1ML	S-12726U1-1ML	Omethoate Solution	1113-02-6	1ML
S-12370A1-1ML	Mephosfolan Solution	950-10-7	1ML	S-12727A1-1ML	Propargite Solution	2312-35-8	1ML
S-12370U1-1ML	Mephosfolan Solution	950-10-7	1ML	S-12727U1-1ML	Propargite Solution	2312-35-8	1ML
S-12371F1-1ML	Mepiquat chloride Solution	24307-26-4	1ML	S-12728A1-1ML	Ormetoprim Solution	6981-18-6	1ML
S-12372M1-1ML	Mepionil Solution	55814-41-0	1ML	S-12728U1-1ML	Ormetoprim Solution	6981-18-6	1ML
S-12376A1-1ML	Mesosulfuron-methyl Solution	208465-21-8	1ML	S-12729A4-1ML	Oryzalin Solution	19044-88-3	1ML
S-12380A1-1ML	Metalaxyl Solution	57837-19-1	1ML	S-12730A1-1ML	Oxex Solution	80-33-1	1ML
S-12380U1-1ML	Metalaxyl Solution	57837-19-1	1ML	S-12730J1-1ML	Oxex Solution	80-33-1	1ML
S-12383M1-1ML	Metamitron Solution	41394-05-2	1ML	S-12731A1-1ML	Oxadiazon Solution	19666-30-9	1ML
S-12384A1-1ML	Metasystox I (TM) Solution	919-86-8	1ML	S-12731U1-1ML	Oxadiazon Solution	19666-30-9	1ML
S-12384U1-1ML	Metasystox I (TM) Solution	919-86-8	1ML	S-12732A1-1ML	Oxadixyl Solution	77732-09-3	1ML
S-12387A1-1ML	Methabenzthiazuron Solution	18691-97-9	1ML	S-12732U1-1ML	Oxadixyl Solution	77732-09-3	1ML
S-12387U1-1ML	Methabenzthiazuron Solution	18691-97-9	1ML	S-12735A1-1ML	Oxamyl Solution	23135-22-0	1ML
S-12388A1-1ML	Methacrifos Solution	30864-28-9	1ML	S-12738A1-1ML	Oxolinic acid Solution	14698-29-4	1ML
S-12388U1-1ML	Methacrifos Solution	30864-28-9	1ML	S-12740A1-1ML	Oxycarboxin Solution	5259-88-1	1ML
S-12393A1-1ML	Methamidophos Solution	10265-92-6	1ML	S-12740U1-1ML	Oxycarboxin Solution	5259-88-1	1ML
S-12393U1-1ML	Methamidophos Solution	10265-92-6	1ML	S-12742M1-1ML	Oxyfluorfen Solution	42874-03-3	1ML
S-12397A1-1ML	Methidathion Solution	950-37-8	1ML	S-12757A1-1ML	p-Bromophenoxy acetic acid Solution	1878-91-7	1ML
S-12397U1-1ML	Methidathion Solution	950-37-8	1ML	S-12757T1-1ML	p-Bromophenoxy acetic acid Solution	1878-91-7	1ML
S-12398M1-1ML	Methiocarb Solution	2032-65-7	1ML	S-12761A1-1ML	p-Chloromandelic acid Solution	492-86-4	1ML
S-12399A1-1ML	Methomyl Solution	16752-77-5	1ML	S-12761T1-1ML	p-Chloromandelic acid Solution	492-86-4	1ML
S-12400T4-1ML	Methoprene Solution	40596-69-8	1ML	S-12763M1-1ML	p-Chlorophenoxyacetic acid Solution	122-88-3	1ML
S-12400T4-5ML	Methoprene Solution	40596-69-8	5ML	S-12763T1-1ML	p-Chlorophenoxyacetic acid Solution	122-88-3	1ML
S-12401A1-1ML	Methoprotryne Solution	841-06-5	1ML	S-12811A1-1ML	Paclobutrazol Solution	76738-62-0	1ML
S-12401T1-1ML	Methoprotryne Solution	841-06-5	1ML	S-12811U1-1ML	Paclobutrazol Solution	76738-62-0	1ML
S-12403M1-1ML	Methoxychlor Solution	72-43-5	1ML	S-12816A1-1ML	Paraaxon Solution	311-45-5	1ML

Part #	Description	CAS #	Size	Part #	Description	CAS #	Size
S-12816U1-1ML	Paraoxon Solution	311-45-5	1ML	S-13098J4-5ML	Profluralin Solution	26399-36-0	5ML
S-12818F1-1ML	Paraquat dichloride Solution	1910-42-5	1ML	S-13101M4-1ML	Promecarb Solution	2631-37-0	1ML
S-12818F17-1ML	Paraquat dichloride hydrate Solution	1910-42-5	1ML	S-13102A1-1ML	Prometon Solution	1610-18-0	1ML
S-12818F17-5ML	Paraquat dichloride hydrate Solution	1910-42-5	5ML	S-13102T1-1ML	Prometon Solution	1610-18-0	1ML
S-12819A1-1ML	Parathion Solution	56-38-2	1ML	S-13102T1-5ML	Prometon Solution	1610-18-0	5ML
S-12819U1-1ML	Parathion Solution	56-38-2	1ML	S-13103A1-1ML	Prometryne Solution	7287-19-6	1ML
S-12822A1-1ML	Penconazole Solution	66246-88-6	1ML	S-13103T1-1ML	Prometryne Solution	7287-19-6	1ML
S-12822U1-1ML	Penconazole Solution	66246-88-6	1ML	S-13103T1-5ML	Prometryne Solution	7287-19-6	5ML
S-12824M1-1ML	Dienochlor Solution	2227-17-0	1ML	S-13104A1-1ML	Propachlor Solution	1918-16-7	1ML
S-12826M1-1ML	Pentachloroanisole Solution	1825-21-4	1ML	S-13104T1-1ML	Propachlor Solution	1918-16-7	1ML
S-12830H5-1ML	Pentachloronitrobenzene Solution	82-68-8	1ML	S-13104T1-5ML	Propachlor Solution	1918-16-7	5ML
S-12830H5-5ML	Pentachloronitrobenzene Solution	82-68-8	5ML	S-13106F1-1ML	Propamocarb HCl Solution	25606-41-1	1ML
S-12830K4-1ML	Pentachloronitrobenzene Solution	82-68-8	1ML	S-13108A1-1ML	Propazine Solution	139-40-2	1ML
S-12830K4-5ML	Pentachloronitrobenzene Solution	82-68-8	5ML	S-13108T1-1ML	Propazine Solution	139-40-2	1ML
S-12830K7-1ML	Pentachloronitrobenzene Solution	82-68-8	1ML	S-13108T1-5ML	Propazine Solution	139-40-2	5ML
S-12830K7-5ML	Pentachloronitrobenzene Solution	82-68-8	5ML	S-13109A4-1ML	Propham Solution	122-42-9	1ML
S-12830M1-1ML	Pentachloronitrobenzene Solution	82-68-8	1ML	S-13110A1-1ML	Prophos Solution	13194-48-4	1ML
S-12830M1-5ML	Pentachloronitrobenzene Solution	82-68-8	5ML	S-13110U1-1ML	Prophos Solution	13194-48-4	1ML
S-12831M1-1ML	Pentachlorophenol Solution	87-86-5	1ML	S-13110U1-5ML	Prophos Solution	13194-48-4	5ML
S-12831M1-5ML	Pentachlorophenol Solution	87-86-5	5ML	S-13119A1-1ML	Propoxycarbazono-sodium Solution	181274-15-7	1ML
S-12844A1-1ML	Pentachlorophenol Solution	2307-68-8	1ML	S-13129A1-1ML	Fluthiacil-methyl Solution	117337-19-6	1ML
S-12844U1-1ML	Pentachlorophenol Solution	2307-68-8	1ML	S-13131K1-1ML	Anilofos Solution	64249-01-0	1ML
S-12845A1-1ML	Perchloromethyl mercaptan Solution	594-42-3	1ML	S-13136M1-1ML	Propyzamide Solution	23950-58-5	1ML
S-12845U1-1ML	Perchloromethyl mercaptan Solution	594-42-3	1ML	S-13136M1-5ML	Propyzamide Solution	23950-58-5	5ML
S-12848A1-1ML	Permethrin Solution	52645-53-1	1ML	S-13137A1-1ML	Proslufuron Solution	94125-34-5	1ML
S-12848T1-1ML	Permethrin Solution	52645-53-1	1ML	S-13138B1-1ML	Prothioconazole Solution	178928-70-6	1ML
S-12849A1-1ML	Ethylan Solution	72-56-0	1ML	S-13139A1-1ML	Flumethrin Solution	69770-45-2	1ML
S-12849K4-1ML	Ethylan Solution	72-56-0	1ML	S-13140A1-1ML	Pendimethalin Solution	40487-42-1	1ML
S-12849K4-5ML	Ethylan Solution	72-56-0	5ML	S-13140K4-1ML	Pendimethalin Solution	40487-42-1	1ML
S-12854A1-1ML	Fenamiphos Solution	22224-92-6	1ML	S-13141A1-1ML	Pymetrozine Solution	123312-89-0	1ML
S-12854T1-1ML	Fenamiphos Solution	22224-92-6	1ML	S-13141U1-1ML	Pymetrozine Solution	123312-89-0	1ML
S-12854T1-5ML	Fenamiphos Solution	22224-92-6	5ML	S-13142A1-1ML	Pyraclafos Solution	89784-60-1	1ML
S-12859A1-1ML	Phenmedipham Solution	13684-63-4	1ML	S-13143A1-1ML	Pyraclafos Solution	89784-60-1	1ML
S-12859U1-1ML	Phenmedipham Solution	13684-63-4	1ML	S-13145M1-1ML	Pyrazon Solution	175013-18-0	1ML
S-12881M1-1ML	Cyanofenphos Solution	13067-93-1	1ML	S-13146A1-1ML	Pyrazon Solution	1698-60-8	1ML
S-12883M1-1ML	Sulfoxaflo Solution	946578-00-3	1ML	S-13146A1-5ML	Pyrazon Solution	1698-60-8	5ML
S-12884M1-1ML	Thifluzamide Solution	130000-40-7	1ML	S-13147A1-1ML	Pyrazophos Solution	13457-18-6	1ML
S-12888U1-1ML	2,4,6-Tris-(2,4,6-tribromophenoxy)-1,3,5-triazine	25713-60-4	1ML	S-13148M1-1ML	Pyrazophos Solution	13457-18-6	1ML
S-12963M1-1ML	Butocarbonyl Solution	34681-10-2	1ML	S-13151A1-1ML	Pyrazosulfuron-ethyl Solution	93697-74-6	1ML
S-12979M1-1ML	Furmecyclo Solution	60568-05-0	1ML	S-13151T1-1ML	Pyrazoxifen Solution	71561-11-0	1ML
S-12987M1-1ML	Flurochloridone Solution	61213-25-0	1ML	S-13152M1-1ML	Pyrethrum Solution	8003-34-7	1ML
S-12988K10-1ML	Ofurace Solution	58810-48-3	1ML	S-13152M1-5ML	Pyrethrum Solution	8003-34-7	5ML
S-12990M1-1ML	Pyriphenox Solution	88283-41-4	1ML	S-13154A1-1ML	Pyridaben Solution	96489-71-3	1ML
S-13004A1-1ML	Phenoxyacetic acid Solution	122-59-8	1ML	S-13154U1-1ML	Pyridaphenthion Solution	119-12-0	1ML
S-13004T1-1ML	Phenoxyacetic acid Solution	122-59-8	1ML	S-13155A1-1ML	Pyridate Solution	55512-33-9	1ML
S-13006U1-1ML	Phenthoate Solution	2597-03-700	1ML	S-13155T1-1ML	Pyridate Solution	55512-33-9	1ML
S-13007A1-1ML	Phenthoacetate Solution	900-95-8	1ML	S-13159A1-1ML	Pyrimethanil Solution	53112-28-0	1ML
S-13007U1-1ML	Phenthoacetate Solution	900-95-8	1ML	S-13159U1-1ML	Pyrimethanil Solution	53112-28-0	1ML
S-13016A1-1ML	Phenyl mercuric phthalate Solution	84-70-8	1ML	S-13160M1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13016U1-1ML	Phenyl mercuric phthalate Solution	84-70-8	1ML	S-13161F1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13020A1-1ML	Phenyl-2-thiourea	103-85-5	1ML	S-13168A1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13020U1-1ML	Phenyl-2-thiourea	103-85-5	1ML	S-13168U1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13030A1-1ML	Phenylmercuric oleate Solution	104-60-9	1ML	S-13169A1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13030U1-1ML	Phenylmercuric oleate Solution	104-60-9	1ML	S-13169T1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13033A1-1ML	Phenylmercury salicylate Solution	28086-13-7	1ML	S-13173A1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13033U1-1ML	Phenylmercury salicylate Solution	28086-13-7	1ML	S-13173U1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13035A1-1ML	Phorate Solution	298-02-2	1ML	S-13174A1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13035J1-1ML	Phorate Solution	298-02-2	1ML	S-13174U1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13036M1-1ML	Phosalone Solution	2310-17-0	1ML	S-13175A1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13037A1-1ML	Mevinphos Solution	7786-34-7	1ML	S-13175U1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13037U1-1ML	Mevinphos Solution	7786-34-7	1ML	S-13176A1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13037U1-5ML	Mevinphos Solution	7786-34-7	5ML	S-13176T4-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13038A1-1ML	Phosfolan Solution	947-02-4	1ML	S-13176T4-5ML	Pyriproxyfen Solution	95737-68-1	5ML
S-13038U1-1ML	Phosfolan Solution	947-02-4	1ML	S-13184A1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13039A1-1ML	Phosphamidon Solution	13171-21-6	1ML	S-13186A1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13039U1-1ML	Phosphamidon Solution	13171-21-6	1ML	S-13186U1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13039U1-5ML	Phosphamidon Solution	13171-21-6	5ML	S-13187A1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13040M1-1ML	Phosphan Solution	115-78-6	1ML	S-13187T1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13041A1-1ML	Phoxim Solution	14816-18-3	1ML	S-13188A1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13041U1-1ML	Phoxim Solution	14816-18-3	1ML	S-13188T1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13043A1-1ML	Phthalanilic acid Solution	4727-29-1	1ML	S-13188T1-5ML	Pyriproxyfen Solution	95737-68-1	5ML
S-13043T1-1ML	Phthalanilic acid Solution	4727-29-1	1ML	S-13189A1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13050A1-1ML	Picloram Solution	1918-02-100	1ML	S-13189J1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13050A4-1ML	Picloram Solution	1918-02-100	1ML	S-13191A1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13050A4-5ML	Picloram Solution	1918-02-100	5ML	S-13191T1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13050B1-1ML	Picloram Solution	1918-02-100	1ML	S-13192A1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13050B1-5ML	Picloram Solution	1918-02-100	5ML	S-13192T1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13051M1-1ML	Picloram methyl ester Solution	14143-55-6	1ML	S-13194A1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13056A1-1ML	Thifensulfuron-methyl Solution	79277-27-3	1ML	S-13194T1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13058A1-1ML	Piperalin Solution	3478-94-2	1ML	S-13197A1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13058U1-1ML	Piperalin Solution	3478-94-2	1ML	S-13197U1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13061A1-1ML	Piperonyl butoxide Solution	51-03-6	1ML	S-13201A4-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13061K1-1ML	Piperonyl butoxide Solution	51-03-6	1ML	S-13201A4-5ML	Pyriproxyfen Solution	95737-68-1	5ML
S-13062A1-1ML	Pirimicarb Solution	23103-98-2	1ML	S-13201U1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13062U1-1ML	Pirimicarb Solution	23103-98-2	1ML	S-13211A4-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13063U1-1ML	Pirimiphos-ethyl Solution	23505-41-1	1ML	S-13211A4-5ML	Pyriproxyfen Solution	95737-68-1	5ML
S-13064A1-1ML	Pirimiphos-methyl Solution	29232-93-7	1ML	S-13212A1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13064U1-1ML	Pirimiphos-methyl Solution	29232-93-7	1ML	S-13212B1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13068M1-1ML	Sedaxane Solution	874967-67-6	1ML	S-13212B1-5ML	Pyriproxyfen Solution	95737-68-1	5ML
S-13070M1-1ML	Pyroxa-sulfone Solution	447399-55-5	1ML	S-13213A1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13074A1-1ML	Indazifam Solution	950782-86-2	1ML	S-13213T1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13075A1-1ML	Amicarbazone Solution	129909-90-6	1ML	S-13214M1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13076M1-1ML	Triclosan Solution	3380-34-5	1ML	S-13214M1-5ML	Pyriproxyfen Solution	95737-68-1	5ML
S-13079M1-1ML	Halifenox Solution	111872-58-3	1ML	S-13215F1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13083A1-1ML	Pyrazolynate Solution	58011-68-0	1ML	S-13215U1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13088A1-1ML	Prallethrin Solution	23031-36-9	1ML	S-13217T1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13088K1-1ML	Prallethrin Solution	23031-36-9	1ML	S-13217T1-5ML	Pyriproxyfen Solution	95737-68-1	5ML
S-13089A1-1ML	Pretilachlor Solution	51218-49-6	1ML	S-13222A1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13089U1-1ML	Pretilachlor Solution	51218-49-6	1ML	S-13225B1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13090A1-1ML	Primisulfuron-methyl Solution	86209-51-0	1ML	S-13226U1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13092M1-1ML	Probenazole Solution	27605-76-1	1ML	S-13229A1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13093A1-1ML	Prochloraz Solution	67747-09-5	1ML	S-13229K4-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13093U1-1ML	Prochloraz Solution	67747-09-5	1ML	S-13229K4-5ML	Pyriproxyfen Solution	95737-68-1	5ML
S-13094A1-1ML	Procyazine Solution	32889-48-8	1ML	S-13242A1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13094T1-1ML	Procyazine Solution	32889-48-8	1ML	S-13243M1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13095A1-1ML	Procyimidine Solution	32809-16-8	1ML	S-13243M1-5ML	Pyriproxyfen Solution	95737-68-1	5ML
S-13095U1-1ML	Procyimidine Solution	32809-16-8	1ML	S-13251M1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13096A1-1ML	Prodiamine Solution	29091-21-2	1ML	S-13253M1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13096K1-1ML	Prodiamine Solution	29091-21-2	1ML	S-13254U1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13097A1-1ML	Profenofos Solution	41198-08-7	1ML	S-13255A1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13097J1-1ML	Profenofos Solution	41198-08-7	1ML	S-13257M1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
S-13098A1-1ML	Profluralin Solution	26399-36-0	1ML	S-13257M1-5ML	Pyriproxyfen Solution	95737-68-1	5ML
S-13098J4-1ML	Profluralin Solution	26399-36-0	1ML	S-13258A1-1ML	Pyriproxyfen Solution	95737-68-1	1ML
				S-13258J1-1ML	Pyriproxyfen Solution	95737-68-1	1ML

Part #	Description	CAS #	Size	Part #	Description	CAS #	Size
S-13258J1-5ML	Sulprofos Solution	35400-43-2	5ML	S-13689U1-1ML	Trifluralin Solution	1582-09-8	1ML
S-13259A1-1ML	Super-Rozol Solution	28772-56-7	1ML	S-13689U1-5ML	Trifluralin Solution	1582-09-8	5ML
S-13259U1-1ML	Super-Rozol Solution	28772-56-7	1ML	S-13690A1-1ML	Triflurosulfuron-methyl Solution	125635-15-7	1ML
S-13260A4-1ML	SWEP Solution	1918-18-9	1ML	S-13691U1-1ML	Triflorine Solution	26644-46-2	1ML
S-13260A4-5ML	SWEP Solution	1918-18-9	5ML	S-13706A1-1ML	Trinexapac-ethyl Solution	95266-40-3	1ML
S-13260U1-1ML	SWEP Solution	1918-18-9	1ML	S-13706T1-1ML	Trinexapac-ethyl Solution	95266-40-3	1ML
S-13261A1-1ML	Myclobutanil Solution	88671-89-0	1ML	S-13731M1-1ML	Ethiozin Solution	64529-56-2	1ML
S-13261U1-1ML	Myclobutanil Solution	88671-89-0	1ML	S-13735M1-1ML	Uniconazole-P Solution	83657-17-4	1ML
S-13263A1-1ML	tau-Fluvalinate Solution	102851-06-9	1ML	S-13738A1-1ML	Pyriminuron Solution	53558-25-1	1ML
S-13263K1-1ML	tau-Fluvalinate Solution	102851-06-9	1ML	S-13738U1-1ML	Pyriminuron Solution	53558-25-1	1ML
S-13275A1-1ML	XMC Solution	2655-14-3	1ML	S-13741A1-1ML	Vamidithion Solution	2275-23-2	1ML
S-13280A1-1ML	Flumorph Solution	211867-47-9	1ML	S-13742A1-1ML	Vanicide-20S (TM) Solution	5902-85-2	1ML
S-13281A1-1ML	Valifenalate Solution	283159-90-0	1ML	S-13742U1-1ML	Vanicide-20S (TM) Solution	5902-85-2	1ML
S-13282A1-1ML	Tiadinil Solution	223580-51-6	1ML	S-13743A1-1ML	Velpar Solution	51235-04-2	1ML
S-13283A1-1ML	Methasulfocarb Solution	66952-49-6	1ML	S-13743T1-1ML	Velpar Solution	51235-04-2	1ML
S-13288M1-1ML	Meptyldinocap Solution	131-72-6	1ML	S-13743T1-5ML	Velpar Solution	51235-04-2	5ML
S-13296A1-1ML	Mesoprazine Solution	1824-09-5	1ML	S-13744A1-1ML	Vernolate Solution	1929-77-7	1ML
S-13500A1-1ML	TCMTB Solution	21564-17-0	1ML	S-13744T1-1ML	Vernolate Solution	1929-77-7	1ML
S-13500X4-1ML	TCMTB Solution	21564-17-0	1ML	S-13744T1-5ML	Vernolate Solution	1929-77-7	5ML
S-13500X4-5ML	TCMTB Solution	21564-17-0	5ML	S-13745A1-1ML	Vinclozolin Solution	50471-44-8	1ML
S-13501A1-1ML	Tebufenozide Solution	112410-23-8	1ML	S-13745K1-1ML	Vinclozolin Solution	50471-44-8	1ML
S-13503A1-1ML	Tebupirimfos Solution	96182-53-5	1ML	S-13750A1-1ML	Warfarin Solution	81-81-2	1ML
S-13503U1-1ML	Tebupirimfos Solution	96182-53-5	1ML	S-13750U1-1ML	Warfarin Solution	81-81-2	1ML
S-13504A1-1ML	Tebutam Solution	35256-85-0	1ML	S-13752M1-1ML	Z-Dimethomorph Solution	113210-98-3	1ML
S-13504U1-1ML	Tebutam Solution	35256-85-0	1ML	S-13753A4-1ML	Zectran Solution	315-18-4	1ML
S-13505A1-1ML	Tebuthiuron Solution	34014-18-1	1ML	S-13753U1-1ML	Zectran Solution	315-18-4	1ML
S-13505T1-1ML	Tebuthiuron Solution	34014-18-1	1ML	S-13754A1-1ML	zeta-Cypermethrin Solution	52315-07-8	1ML
S-13505T1-5ML	Tebuthiuron Solution	34014-18-1	5ML	S-13754U1-1ML	zeta-Cypermethrin Solution	52315-07-8	1ML
S-13506A1-1ML	Teflubenzuron Solution	83121-18-0	1ML	S-13760A1-1ML	Zinophos (TM) Solution	297-97-2	1ML
S-13509A1-1ML	Terbacil Solution	5902-51-2	1ML	S-13760J1-1ML	Zinophos (TM) Solution	297-97-2	1ML
S-13509T1-1ML	Terbacil Solution	5902-51-2	1ML	S-13760J1-5ML	Zinophos (TM) Solution	297-97-2	5ML
S-13509T1-5ML	Terbacil Solution	5902-51-2	5ML	S-13800B4-1ML	Simazine Solution	122-34-9	1ML
S-13510M1-1ML	Terbufos Solution	13071-79-9	1ML	S-13800B4-5ML	Simazine Solution	122-34-9	5ML
S-13511M1-1ML	Terbumeton Solution	33693-04-8	1ML	S-13800M1-1ML	Simazine Solution	122-34-9	1ML
S-13512A1-1ML	Terbutylazine Solution	5915-41-3	1ML	S-13800M1-5ML	Simazine Solution	122-34-9	5ML
S-13512T4-1ML	Terbutylazine Solution	5915-41-3	1ML	S-13802A1-1ML	Simetryn Solution	1014-70-6	1ML
S-13514A1-1ML	Terbutryne Solution	886-50-0	1ML	S-13802T1-1ML	Simetryn Solution	1014-70-6	1ML
S-13514T1-1ML	Terbutryne Solution	886-50-0	1ML	S-13802T1-5ML	Simetryn Solution	1014-70-6	5ML
S-13514T1-5ML	Terbutryne Solution	886-50-0	5ML	S-13818A1-1ML	Imazosulfuron	122548-33-8	1ML
S-13517A1-1ML	Etridiazole Solution	2593-15-9	1ML	S-13920B1-1ML	Fipronil desulfinyl Solution	205650-63-3	1ML
S-13517T1-1ML	Etridiazole Solution	2593-15-9	1ML	S-13921M1-1ML	Butoxycarboxim Solution	34681-23-7	1ML
S-13537A1-1ML	Tetrachlorvinphos Solution	22248-79-9	1ML	S-13926A1-1ML	Dimoxystrobin Solution	149961-52-4	1ML
S-13537X1-1ML	Tetrachlorvinphos Solution	22248-79-9	1ML	S-13927A1-1ML	Picoxystrobin Solution	117428-22-5	1ML
S-13537X1-5ML	Tetrachlorvinphos Solution	22248-79-9	5ML	S-13928A1-1ML	Acrinathrin Solution	101007-06-1	1ML
S-13538A1-1ML	Tetraconazole Solution	112281-77-3	1ML	S-13929K4-1ML	o,p'-Dicofol Solution		1ML
S-13538U1-1ML	Tetraconazole Solution	112281-77-3	1ML	S-13931A1-1ML	Simeconazole Solution	149508-90-7	1ML
S-13539A1-1ML	Tetradifon Solution	116-29-0	1ML	S-13933A1-1ML	Isocarboxophos Solution	24353-61-5	1ML
S-13539B3-1ML	Tetradifon Solution	116-29-0	1ML	S-13939M1-1ML	Flamprop-M-isopropyl Solution	57973-67-8	1ML
S-13539B3-5ML	Tetradifon Solution	116-29-0	5ML	S-13940M1-1ML	Fluxofenim Solution	88485-37-4	1ML
S-13540A1-1ML	Tetraethylthiopyrophosphate Solution	3689-24-5	1ML	S-13941A1-1ML	Flupyradifurone Solution	951659-40-8	1ML
S-13540J1-1ML	Tetraethylthiopyrophosphate Solution	3689-24-5	1ML	S-13942M1-1ML	Nuarimol Solution	63284-71-9	1ML
S-13552A4-1ML	Tetramethrin Solution	7696-12-0	1ML	S-13944B1-1ML	Orthosulfamuron Solution	213464-77-8	1ML
S-13552A4-5ML	Tetramethrin Solution	7696-12-0	5ML	S-13946M0-1ML	sec-Acetochlor OA Solution	152019-74-4	1ML
S-13552U1-1ML	Tetramethrin Solution	7696-12-0	1ML	S-13955A1-1ML	Daimuron Solution	42609-52-9	1ML
S-13560M4-1ML	Thiabendazole Solution	148-79-8	1ML	S-13956K1-1ML	Triazamate Solution		1ML
S-13561A1-1ML	Thiacloprid Solution	111988-49-9	1ML	S-13960M1-1ML	Amisulbrom Solution	348635-87-0	1ML
S-13561U1-1ML	Thiacloprid Solution	111988-49-9	1ML	S-13962A1-1ML	Fenothicarb Solution	62850-32-2	1ML
S-13562M1-1ML	Thiamethoxam Solution	153719-23-4	1ML	S-13963A1-1ML	Pyrimidifen Solution	105779-78-0	1ML
S-13563A1-1ML	Thiazopyr Solution	117718-60-2	1ML	S-13964A1-1ML	Dioxabenzofos Solution	3811-49-2	1ML
S-13563U1-1ML	Thiazopyr Solution	117718-60-2	1ML	S-13965A1-1ML	Tridiphane Solution	58138-08-2	1ML
S-13564A1-1ML	Thidiazuron Solution	51707-55-2	1ML	S-13966A1-1ML	Thenylchlor Solution	96491-05-3	1ML
S-13564U1-1ML	Thidiazuron Solution	51707-55-2	1ML	S-13967A1-1ML	Propaphos Solution	7292-16-2	1ML
S-13568A1-1ML	Thiodicarb Solution	59689-26-0	1ML	S-13968A1-1ML	Imibenconazole Solution	86598-92-7	1ML
S-13570U1-1ML	Thiometon Solution	640-15-3	1ML	S-13970A1-1ML	Tepraloxidim Solution	149979-41-9	1ML
S-13571U1-1ML	Thiophanate Solution	23564-06-9	1ML	S-13971A1-1ML	Isouron Solution	55861-78-4	1ML
S-13572A1-1ML	Thiophanate-methyl Solution	23564-05-8	1ML	S-13973A1-1ML	Chromafenozide Solution	143807-66-3	1ML
S-13576A1-1ML	Propiconazole Solution	60207-90-1	1ML	S-13974A1-1ML	Carproamid Solution	104030-54-8	1ML
S-13576E1-1ML	Propiconazole Solution	60207-90-1	1ML	S-13975A1-1ML	Azimsulfuron Solution	102162-55-2	1ML
S-13577A1-1ML	Prothiophos Solution	34643-46-4	1ML	S-13976A1-1ML	Azaconazole Solution	60207-31-0	1ML
S-13577U1-1ML	Prothiophos Solution	34643-46-4	1ML	S-13977M1-1ML	Bromobutide Solution	74712-19-9	1ML
S-13578A1-1ML	Tolclofos-Methyl Solution	57018-04-9	1ML	S-13978A1-1ML	Dimepiperate Solution	61432-55-1	1ML
S-13578U1-1ML	Tolclofos-Methyl Solution	57018-04-9	1ML	S-13979A1-1ML	Diofenolan Solution	63837-33-2	1ML
S-13586B6-1ML	Camphechlor Solution	8001-35-2	1ML	S-13980M1-1ML	Fluacrypyrim Solution	229977-93-9	1ML
S-13586J4-1ML	Camphechlor Solution	8001-35-2	1ML	S-13981A1-1ML	Flurtamone Solution	96525-23-4	1ML
S-13586J4-5ML	Camphechlor Solution	8001-35-2	5ML	S-13982M1-1ML	Picolinafen Solution	137641-05-5	1ML
S-13586M1-1ML	Camphechlor Solution	8001-35-2	1ML	S-13983A1-1ML	Proquinazid Solution		1ML
S-13586M1-5ML	Camphechlor Solution	8001-35-2	5ML	S-13984J1-1ML	2,3,5,6-Tetrachloroanisole Solution	6936-40-9	1ML
S-13587A1-1ML	Tralkoxydim Solution	87820-88-0	1ML	S-13985F1-1ML	2,4-D Choline salt Solution	1048373-72-3	1ML
S-13587U1-1ML	Tralkoxydim Solution	87820-88-0	1ML	S-13988A1-1ML	Fluopyram Solution	658066-35-4	1ML
S-13588A1-1ML	Talomechlin Solution	66841-25-6	1ML	S-13990A1-1ML	2,3,4,5-Tetrachloroaniline Solution	634-83-3	1ML
S-13588T1-1ML	Talomechlin Solution	66841-25-6	1ML	S-13991A1-1ML	Acionifen Solution	74070-46-5	1ML
S-13615M1-1ML	trans-Chlordane Solution	5103-74-2	1ML	S-13994A1-1ML	(1R)-trans-Tetramethrin Solution	1166-46-7	1ML
S-13615M1-5ML	trans-Chlordane Solution	5103-74-2	5ML	S-13995X1-1ML	Tetramine Solution	80-12-6	1ML
S-13618M1-1ML	trans-Heptachlor epoxide Solution	28044-83-9	1ML	S-13996A1-1ML	Isoprazam Solution	881685-58-1	1ML
S-13620A1-1ML	trans-Permethrin Solution	61949-77-7	1ML	S-14002M1-1ML	Perfluorooctanesulfonic acid Solution	11763-23-1	1ML
S-13620T1-1ML	trans-Permethrin Solution	61949-77-7	1ML	S-14006A1-1ML	Imiprothrin Solution	72963-72-5	1ML
S-13628A1-1ML	Tri-allate Solution	2303-17-5	1ML	S-14007A1-1ML	Propaquizafop Solution	111479-05-1	1ML
S-13628K1-1ML	Tri-allate Solution	2303-17-5	1ML	S-14010X1-1ML	Benzovindiflupyr Solution	1072957-71-1	1ML
S-13636A1-1ML	Triadimefon Solution	43121-43-3	1ML	S-14011A1-1ML	Oxadiazyl Solution	39807-15-3	1ML
S-13636T1-1ML	Triadimefon Solution	43121-43-3	1ML	S-14012X1-1ML	Octhilonone Solution	26530-20-1	1ML
S-13636T1-5ML	Triadimefon Solution	43121-43-3	5ML	S-14013M1-1ML	Fenfuram Solution	24691-80-3	1ML
S-13638A1-1ML	Triasulfuron Solution	82097-50-5	1ML	S-14029A1-1ML	Flusulfamide Solution	106917-52-6	1ML
S-13639A1-1ML	Triazophos Solution	24017-47-8	1ML	S-14030A1-1ML	Ferimzone Solution	89269-64-7	1ML
S-13639U1-1ML	Triazophos Solution	24017-47-8	1ML	S-14031A1-1ML	(E)-Ferimzone Solution	77359-18-3	1ML
S-13639U1-5ML	Triazophos Solution	24017-47-8	5ML	S-14032A1-1ML	Uniconazole Solution	83657-22-1	1ML
S-13640A1-1ML	Tribenuron methyl Solution	101200-48-0	1ML	S-14035A1-1ML	Penflufen Solution	494793-67-8	1ML
S-13652T1-1ML	Trichloroacetic acid Solution	76-03-9	1ML	S-14038A1-1ML	Cyhalofop Solution	122008-78-0	1ML
S-13652T1-5ML	Trichloroacetic acid Solution	76-03-9	5ML	S-14042A1-1ML	Fluazaindolizine Solution	1254304-22-7	1ML
S-13657A1-1ML	Trichloronate Solution	327-98-0	1ML	S-14050M1-1ML	Cinmethylin Solution	87818-31-3	1ML
S-13657L1-1ML	Trichloronate Solution	327-98-0	1ML	S-14051A1-1ML	Cymiazole Solution	61676-87-7	1ML
S-13657L1-5ML	Trichloronate Solution	327-98-0	5ML	S-14052U1-1ML	Dimethylvinphos Solution	2274-67-1	1ML
S-13659A1-1ML	Triclopyr Solution	55335-06-3	1ML	S-14054J1-1ML	Oxy-chlordane Solution	33386-84-4	1ML
S-13659T1-1ML	Triclopyr Solution	55335-06-3	1ML	S-14058A1-1ML	Thiofanox Solution	39196-18-4	1ML
S-13661A1-1ML	Triclopyr methyl ester Solution	60825-26-5	1ML	S-14059A1-1ML	Fentrazamide Solution	158237-07-1	1ML
S-13661T1-1ML	Triclopyr methyl ester Solution	60825-26-5	1ML	S-FC100S-1.2ML	Heptachlor (13C10) Solution		1.2ML
S-13662M1-1ML	Tricyclohexylin hydroxide Solution	13121-70-5	1ML	S-FC101S-1.2ML	Heptachlor epoxide (Isomer B) (13C1) Solution		1.2ML
S-13665A1-1ML	Tridemorph (Mix of Isomers) Solution	81412-43-3	1ML	S-FC102S-1.2ML	BHC (alpha isomer) (13C6) Solution		1.2ML
S-13665U1-1ML	Tridemorph (Mix of Isomers) Solution	81412-43-3	1ML	S-FC103S-2X1.2ML	BHC (beta isomer) (13C6) Solution		2X 1.2ML
S-13666A1-1ML	Trietazine Solution	1912-26-1	1ML	S-FC104S-1.2ML	Lindane (BHC gamma isomer) (13C6) Solution		1.2ML
S-13666X1-1ML	Trietazine Solution	1912-26-1	1ML	S-FC105S-1.2ML	BHC (delta isomer) (13C6) Solution		1.2ML
S-13684A1-1ML	Trifloxystrobin Solution	141517-21-7	1ML	S-FC2005S-1.2ML	Carbaryl (ring-13C6) Solution		1.2ML
S-13687A1-1ML	Triflurumuron Solution	64628-44-0	1ML	S-FC2006S-1.2ML	Carbofuran (ring 13C6) Solution		1.2ML
S-13687U1-1ML	Triflurumuron Solution	64628-44-0	1ML	S-FC202S-1.2ML	a-Endosulfan (13C9) Solution		1.2ML
S-13689A1-1ML	Trifluralin Solution	1582-09-8	1ML	S-FC203S-1.2ML	b-Endosulfan (13C9) Solution		1.2ML

Part #	Description	CAS #	Size	Part #	Description	CAS #	Size
S-FC2072S-1.2ML	Alachlor (ring-13C6) Solution		1.2ML				
S-FC2073S-1.2ML	Simazine (ring-13C3) Solution		1.2ML				
S-FC2093S-1.2ML	Bromoxynil (ring-13C6) Solution		1.2ML				
S-FC2121S-1ML	Mirex (13C8) Solution		1ML				
S-FC2208S-1.2ML	Atrazine (13C3) Solution		1.2ML				
S-FC2216AS-1.2ML	cis-Permethrin (phenoxy-13C6) Solution		1.2ML				
S-FC2216BS-1.2ML	trans-Permethrin (phenoxy-13C6) Solution		1.2ML				
S-FC2236S-1.2ML	Propazine (ring-13C3) Solution		1.2ML				
S-FC2241S-1.2ML	Metolachlor (ring-13C6) Solution		1.2ML				
S-FC2306S-1.2ML	o,p'-DDE (13C12) Solution		1.2ML				
S-FC2307S-1.2ML	o,p'-DDT (13C12) Solution		1.2ML				
S-FC2441S-1.2ML	Kepone (TM) (13C10) Solution		1.2ML				
S-FC89S-1.2ML	Aldrin (TM) (13C12, 99%) Solution		1.2ML				
S-FC908S-1.2ML	trans-Chlordane (13C10) Solution		1.2ML				
S-FC90AS-1.2ML	Dieldrin (13C12, 99%) Solution		1.2ML				
S-FC91S-1.2ML	Chlordane (random-13C1) Solution		1.2ML				
S-FC92S-1.2ML	4,4'-DDT (13C12) Solution		1.2ML				
S-FC93S-1.2ML	4,4'-DDE (13C12) Solution		1.2ML				
S-FC962S-1.2ML	Dichlorprop (ring-13C6) Solution		1.2ML				
S-FC972S-1.2ML	2,4,5-T (ring-13C6) Solution		1.2ML				
S-FCD104S-1.2ML	Lindane (BHC gamma isomer) (13C6, D6) Solution		1.2ML				
S-FCN1104S-1.2ML	Glyphosate (2-13C, 15N) Solution		1.2ML				
S-FD202S-1.2ML	a-Endosulfan (d4) Solution		1.2ML				
S-FD203S-1.2ML	b-Endosulfan (d4) Solution		1.2ML				
S-FD2057S-1.2ML	Chlorpyrifos (diethyl-d10) Solution	285138-81-0	1.2ML				
S-FD2060S-1.2ML	Diazinon (diethyl-d10) Solution		1.2ML				
S-FD2118S-1.2ML	Malathion (d10) Solution		1.2ML				
S-FD2129S-1.2ML	Phosmet (dimethyl-d6) (TM) Solution		1.2ML				
S-FD2140S-1.2ML	Trifluralin (di-n-propyl-d14) Solution		1.2ML				
S-FD2208S-1.2ML	Atrazine (ethylamine-d5) Solution		1.2ML				
S-FD94S-1.2ML	4,4'-DDD (ring-d8) Solution		1.2ML				
S-FD998S-1.2ML	Parathion (TM) (diethyl-d10) Solution	350820-04-1	1.2ML				



Pesticide & Metabolite Standards Listing by CAS Number

a



50-00-0

b



Formaldehyde (in water)

c



N-12012-1G

a Chemical Abstracts Registry Number

c Part Number

b Description

CAS#	Description	Part #	CAS#	Description	Part #
50-00-0	Formaldehyde (in water)	N-12012-1G	78-48-8	S,S,S-Tributyl phosphorothriothioate	N-13194-250MG
50-29-3	4,4'-DDT	N-10876-100MG	78-48-8	S,S,S-Tributyl phosphorothriothioate Solution	S-13194A1-1ML
50-29-3	4,4'-DDT Solution	S-10876J1-1ML	78-48-8	S,S,S-Tributyl phosphorothriothioate Solution	S-13194T1-1ML
50-30-6	2,6-Dichlorobenzoic acid	MET-10688A-250MG	78-87-5	1,2-Dichloropropane	N-10156-1G
50-31-7	2,3,6-Trichlorobenzoic acid	N-10605-100MG	78-87-5	1,2-Dichloropropane Solution	N-10156M1-1ML
50-65-7	Niclosamide	N-12653-250MG	78-87-5	1,2-Dichloropropane Solution	N-10156M1-5ML
51-03-6	Piperonyl butoxide	N-13061-100MG	79-11-8	Chloroacetic acid	N-11438-250MG
51-03-6	Piperonyl butoxide Solution	S-13061A1-1ML	79-11-8	Chloroacetic acid Solution	S-11438T1-1ML
51-03-6	Piperonyl butoxide Solution	S-13061K1-1ML	80-12-6	Tetramine Solution	S-13995X1-1ML
51-36-5	3,5-Dichlorobenzoic acid	N-10770-250MG	80-32-0	Sulfachloropyridazine	N-13241-250MG
51-36-5	3,5-Dichlorobenzoic acid Solution	S-10770A1-1ML	80-33-1	Ovex	N-12730-250MG
51-36-5	3,5-Dichlorobenzoic acid Solution	S-10770T5-1ML	80-33-1	Ovex Solution	S-12730A1-1ML
51-36-5	3,5-Dichlorobenzoic acid Solution	S-10770T5-5ML	80-33-1	Ovex Solution	S-12730J1-1ML
52-51-7	Bronopol	N-11348-250MG	80-38-6	Fenson	N-13812-250MG
52-51-7	Bronopol Solution	S-11348M1-1ML	81-81-2	Warfarin	N-13750-250MG
52-68-6	Trichlorophen	N-11843-250MG	81-81-2	Warfarin Solution	S-13750A1-1ML
52-85-7	Famphur	N-11944-250MG	81-81-2	Warfarin Solution	S-13750U1-1ML
52-85-7	Famphur Solution	S-11944K1-1ML	81-84-5	1-8-Naphthalic anhydride	N-10012-500MG
53-19-0	o,p'-DDD	N-12706-250MG	81-84-5	1-8-Naphthalic anhydride Solution	S-10012A1-1ML
53-19-0	o,p'-DDD Solution	S-12706M1-1ML	81-84-5	1-8-Naphthalic anhydride Solution	S-10012U1-1ML
54-11-5	(S)-Nicotine	N-12655-1G	82-66-6	Diphacinone	N-11793-100MG
54-11-5	Nicotine Solution	S-12655M1-1ML	82-66-6	Diphacinone Solution	S-11793A1-1ML
55-38-9	Fenthion	N-11964-250MG	82-68-8	Pentachloronitrobenzene	N-12830-250MG
55-38-9	Fenthion Solution	S-11964S1-1ML	82-68-8	Pentachloronitrobenzene Solution	S-12830H5-1ML
55-38-9	Fenthion Solution	S-11964S1-5ML	82-68-8	Pentachloronitrobenzene Solution	S-12830H5-5ML
55-68-5	Phenylmercury nitrate(Technical)	N-13032-100MG	82-68-8	Pentachloronitrobenzene Solution	S-12830K4-1ML
56-35-9	Bis[tri-n-butyltin]oxide(Technical)	N-11322-250MG	82-68-8	Pentachloronitrobenzene Solution	S-12830K4-5ML
56-38-2	Parathion	N-12819-100MG	82-68-8	Pentachloronitrobenzene Solution	S-12830K7-1ML
56-38-2	Parathion Solution	S-12819A1-1ML	82-68-8	Pentachloronitrobenzene Solution	S-12830K7-5ML
56-38-2	Parathion Solution	S-12819U1-1ML	82-68-8	Pentachloronitrobenzene Solution	S-12830M1-1ML
56-72-4	Coumaphos	N-11507-100MG	82-68-8	Pentachloronitrobenzene Solution	S-12830M1-5ML
56-72-4	Coumaphos Solution	S-11507A1-1ML	83-05-6	p,p'-DDA	MET-12708A-100MG
56-72-4	Coumaphos Solution	S-11507U1-1ML	83-26-1	2-Pivaloyl-1,3-indandione	N-10512-250MG
56-72-4	Coumaphos Solution	S-11507U1-5ML	83-26-1	2-Pivaloyl-1,3-indandione Solution	S-10512A1-1ML
57-00-1	Creatine	MET-11840A-250MG	83-26-1	2-Pivaloyl-1,3-indandione Solution	S-10512U1-1ML
57-68-1	Sulfamethazine	N-13245-250MG	83-28-3	2-Isovaleryl-1,3-indanedione	N-10384-250MG
57-74-9	Chlordane	N-11425-250MG	83-28-3	2-Isovaleryl-1,3-indanedione Solution	S-10384A1-1ML
57-74-9	Chlordane Solution	S-11425J4-1ML	83-28-3	2-Isovaleryl-1,3-indanedione Solution	S-10384U1-1ML
57-74-9	Chlordane Solution	S-11425J4-5ML	83-79-4	Rotenone	N-13184-250MG
57-74-9	Chlordane Solution	S-11425M1-1ML	83-79-4	Rotenone Solution	S-13184A1-1ML
58-89-9	Lindane (BHC gamma isomer)	N-12319-250MG	84-65-1	9,10-Anthraquinone	N-10970-250MG
58-89-9	Lindane (BHC gamma isomer) Solution	S-12319A1-1ML	84-65-1	9,10-Anthraquinone Solution	S-10970M1-1ML
58-90-2	2,3,4,6-Tetrachlorophenol	N-10599-50MG	84-70-8	Phenyl mercuric phthalate	N-13016-100MG
58-90-2	2,3,4,6-Tetrachlorophenol Solution	S-10599M1-1ML	84-70-8	Phenyl mercuric phthalate Solution	S-13016A1-1ML
59-40-5	Sulfaquinoxaline	N-13251-10MG	84-70-8	Phenyl mercuric phthalate Solution	S-13016U1-1ML
59-40-5	Sulfaquinoxaline Solution	S-13251M1-1ML	84-74-2	Di-n-butyl phthalate Solution	S-11589A1-1ML
60-51-5	Dimethoate	N-11758-100MG	84-74-2	Di-n-butyl phthalate Solution	S-11589J1-1ML
60-51-5	Dimethoate Solution	S-11758J1-1ML	84-74-2	Di-n-butyl phthalate Solution	S-11589J1-5ML
60-51-5	Dimethoate Solution	S-11758J1-5ML	85-29-0	2,4'-Dichlorobenzophenone	MET-11948-100MG
60-57-1	Dieldrin	N-11688-100MG	85-34-7	Fenac	N-11945-250MG
60-57-1	Dieldrin Solution	S-11688M1-1ML	85-34-7	Fenac Solution	S-11945A1-1ML
60-57-1	Dieldrin Solution	S-11688M1-5ML	85-34-7	Fenac Solution	S-11945T1-1ML
61-82-5	Aminotriazole	N-11067-250MG	85-72-3	N-m-Tolylphthalamic acid	N-12557-100MG
61-82-5	Aminotriazole Solution	S-11067A1-1ML	85-72-3	N-m-Tolylphthalamic acid Solution	S-12557A1-1ML
61-82-5	Aminotriazole Solution	S-11067B1-1ML	85-72-3	N-m-Tolylphthalamic acid Solution	S-12557T1-1ML
62-38-4	Phenyl mercuric acetate(Technical)	N-13011-250MG	86-29-3	Diphenylacetone	N-11800-100MG
62-73-7	Dichlorvos	N-11675-250MG	86-29-3	Diphenylacetone Solution	S-11800M1-1ML
62-73-7	Dichlorvos Solution	S-11675J4-1ML	86-50-0	Azinphos-methyl	N-12139-250MG
62-73-7	Dichlorvos Solution	S-11675U1-1ML	86-50-0	Azinphos-methyl Solution	S-12139A1-1ML
62-74-8	Sodium fluoroacetate(Technical)	N-13216-250MG	86-50-0	Azinphos-methyl Solution	S-12139U1-1ML
63-25-2	Carbaryl	N-11402-250MG	86-50-0	Azinphos-methyl Solution	S-12139U1-5ML
63-25-2	Carbaryl Solution	S-11402A4-1ML	86-86-2	1-Naphthaleneacetamide	N-10074-1G
63-25-2	Carbaryl Solution	S-11402A4-5ML	86-86-2	1-Naphthaleneacetamide Solution	S-10074M1-1ML
66-81-9	Cycloheximide	N-11534-100MG	86-86-2	1-Naphthaleneacetamide Solution	S-10074U1-1ML
66-81-9	Cycloheximide Solution	S-11534A1-1ML	86-87-3	1-Naphthaleneacetic acid	N-10075-1G
66-81-9	Cycloheximide Solution	S-11534U1-1ML	86-87-3	1-Naphthaleneacetic acid Solution	S-10075M1-1ML
67-97-0	Cholecalciferol	N-11463-100MG	86-87-3	1-Naphthaleneacetic acid Solution	S-10075T1-1ML
67-97-0	Cholecalciferol Solution	S-11463A1-1ML	86-88-4	a-Naphthylthiourea	N-10983-250MG
67-97-0	Cholecalciferol Solution	S-11463U1-1ML	86-88-4	a-Naphthylthiourea Solution	S-10983A1-1ML
70-30-4	Hexachlorophene	N-12163-250MG	86-88-4	a-Naphthylthiourea Solution	S-10983U1-1ML
70-30-4	Hexachlorophene Solution	S-12163M1-1ML	87-51-4	3-Indoleacetic acid	N-10727-250MG
70-30-4	Hexachlorophene Solution	S-12163M1-5ML	87-51-4	3-Indoleacetic acid Solution	S-10727A1-1ML
72-14-0	Sulfathiazole	N-13252-250MG	87-51-4	3-Indoleacetic acid Solution	S-10727T1-1ML
72-20-8	Endrin	N-11854-100MG	87-86-5	Pentachlorophenol	N-12831-500MG
72-20-8	Endrin Solution	S-11854M1-1ML	87-86-5	Pentachlorophenol Solution	S-12831M1-1ML
72-20-8	Endrin Solution	S-11854M1-5ML	87-86-5	Pentachlorophenol Solution	S-12831M1-5ML
72-43-5	o,p'-Methoxychlor Solution	MET-12403AM1-1ML	87-87-6	Tetrachlorohydroquinone	MET-12159A-250MG
72-43-5	Methoxychlor	N-12403-100MG	88-06-2	2,4,6-Trichlorophenol	N-10657-250MG
72-43-5	Methoxychlor Solution	S-12403M1-1ML	88-06-2	2,4,6-Trichlorophenol Solution	S-10657M1-1ML
72-43-5	Methoxychlor Solution	S-12403M1-5ML	88-06-2	2,4,6-Trichlorophenol Solution	S-10657M1-5ML
72-54-8	4,4'-DDD	N-10874-250MG	88-82-4	2,3,5-Triiodobenzoic acid	N-10601-250MG
72-54-8	4,4'-DDD Solution	S-10874A1-1ML	88-82-4	2,3,5-Triiodobenzoic acid Solution	S-10601A1-1ML
72-54-8	4,4'-DDD Solution	S-10874K1-1ML	88-82-4	2,3,5-Triiodobenzoic acid Solution	S-10601T1-1ML
72-54-8	4,4'-DDD Solution	S-10874K1-5ML	88-85-7	Dinoseb	N-11786-100MG
72-55-9	4,4'-DDE	N-10875-100MG	88-85-7	Dinoseb Solution	S-11786A1-1ML
72-55-9	4,4'-DDE Solution	S-10875M1-1ML	88-85-7	Dinoseb Solution	S-11786B1-1ML
72-55-9	4,4'-DDE Solution	S-10875M1-5ML	88-97-1	Phthalamic acid	MET-12007A-100MG
72-56-0	Ethylan	N-12849-250MG	89-61-2	2,5-Dichloronitrobenzene	N-10670-250MG
72-56-0	Ethylan Solution	S-12849A1-1ML	89-61-2	2,5-Dichloronitrobenzene Solution	S-10670A1-1ML
72-56-0	Ethylan Solution	S-12849K4-1ML	89-61-2	2,5-Dichloronitrobenzene Solution	S-10670U1-1ML
72-56-0	Ethylan Solution	S-12849K4-5ML	90-43-7	2-Biphenylol	N-10280-250MG
75-60-5	Dimethylarsinic acid	N-11779-500MG	90-43-7	2-Biphenylol Solution	S-10280A1-1ML
75-60-5	Dimethylarsinic acid Solution	S-11779F4-1ML	90-43-7	2-Biphenylol Solution	S-10280M5-1ML
75-99-0	Dalapon	N-11562-250MG	90-97-1	4,4'-Dichlorobenzhydrol	MET-11442A-250MG
75-99-0	Dalapon Solution	S-11562B1-1ML	90-98-2	4,4'-Dichlorobenzophenone	MET-11442B-100MG
75-99-0	Dalapon Solution	S-11562M1-1ML	91-53-2	Ethoxyquin	N-11877-100MG
76-03-9	Trichloroacetic acid	N-13652-250MG	92-52-4	Biphenyl	N-11207-500MG
76-03-9	Trichloroacetic acid Solution	S-13652T1-1ML	92-52-4	Biphenyl Solution	S-11207J1-1ML
76-03-9	Trichloroacetic acid Solution	S-13652T1-5ML	92-52-4	Biphenyl Solution	S-11207J1-5ML
76-06-2	Trichloronitromethane	N-11452-100MG	92-52-4	Biphenyl Solution	S-11207M4-1ML
76-06-2	Trichloronitromethane Solution	S-11452B1-1ML	92-52-4	Biphenyl Solution	S-11207M4-5ML
76-06-2	Trichloronitromethane Solution	S-11452B1-5ML	93-35-6	7-Hydroxycoumarin	MET-13750A-250MG
76-44-8	Heptachlor	N-12147-100MG	93-71-0	N,N-Diallyl-2-chloroacetamide	N-12617-250MG
76-44-8	Heptachlor Solution	S-12147M1-1ML	93-71-0	N,N-Diallyl-2-chloroacetamide Solution	S-12617M1-1ML
76-44-8	Heptachlor Solution	S-12147M1-5ML	93-72-1	Silvex	N-13212-250MG
76-87-9	Fentin hydroxide	N-11966-250MG	93-72-1	Silvex Solution	S-13212A1-1ML
76-87-9	Fentin hydroxide Solution	S-11966A1-1ML	93-72-1	Silvex Solution	S-13212B1-1ML
77-06-5	Gibberellic acid	N-12108-100MG	93-72-1	Silvex Solution	S-13212B1-5ML
77-06-5	Gibberellic acid Solution	S-12108A1-1ML	93-76-5	2,4,5-T	N-10646-250MG
77-06-5	Gibberellic acid Solution	S-12108T1-1ML	93-76-5	2,4,5-T Solution	S-10646A1-1ML
77-73-6	Dicyclopentadiene	N-11686-500MG	93-76-5	2,4,5-T Solution	S-10646B1-1ML
77-73-6	Dicyclopentadiene Solution	S-11686K1-1ML	93-76-5	2,4,5-T Solution	S-10646B1-5ML
78-34-2	Dioxathion	N-11790-100MG	93-78-7	2,4,5-T isopropyl ester	N-10650-100MG
78-34-2	Dioxathion Solution	S-11790B1-1ML	93-78-7	2,4,5-T isopropyl ester Solution	S-10650A1-1ML

CAS#	Description	Part #
93-78-7	2,4,5-T isopropyl ester Solution	S-10650T1-1ML
93-79-8	2,4,5-T n-butyl ester	N-10651-250MG
93-79-8	2,4,5-T n-butyl ester Solution	S-10651A1-1ML
93-79-8	2,4,5-T n-butyl ester Solution	S-10651J1-1ML
94-11-1	2,4-D isopropyl ester	N-10617-250MG
94-11-1	2,4-D isopropyl ester Solution	S-10617A1-1ML
94-11-1	2,4-D isopropyl ester Solution	S-10617U1-1ML
94-43-9	Phenyl mercuric benzoate	N-13012-100MG
94-74-6	4-Chloro-o-tolylxyacetic acid	N-10818-250MG
94-74-6	4-Chloro-o-tolylxyacetic acid Solution	S-10818B1-1ML
94-74-6	4-Chloro-o-tolylxyacetic acid Solution	S-10818M1-1ML
94-74-6	4-Chloro-o-tolylxyacetic acid Solution	S-10818B1-1ML
94-75-7	2,4-D	N-10609-250MG
94-75-7	2,4-D Solution	S-10609B1-1ML
94-75-7	2,4-D Solution	S-10609M1-1ML
94-75-7	2,4-D sec-butyl ester Solution	S-10621A4-1ML
94-75-7	2,4-D sec-butyl ester Solution	S-10621A4-1ML
94-75-7	2,4-D sec-butyl ester Solution	S-10621U1-1ML
94-79-1	2,4-D sec-butyl ester	N-10621-250MG
94-80-4	2,4-D butyl ester	N-10611-250MG
94-80-4	2,4-D butyl ester Solution	S-10611A1-1ML
94-80-4	2,4-D butyl ester Solution	S-10611U1-1ML
94-81-5	MCPB	N-12357-100MG
94-81-5	MCPB Solution	S-12357A1-1ML
94-81-5	MCPB Solution	S-12357T1-1ML
94-82-6	2,4-DB	N-10622-250MG
94-82-6	2,4-DB Solution	S-10622A1-1ML
94-82-6	2,4-DB Solution	S-10622B1-1ML
94-82-6	2,4-DB Solution	S-10622B1-5ML
94-96-2	2-Ethyl-1,3-hexanediol	N-10346-500MG
94-96-2	2-Ethyl-1,3-hexanediol Solution	S-10346B6-1ML
94-96-2	2-Ethyl-1,3-hexanediol Solution	S-10346B6-5ML
95-06-7	Sulfallate	N-13243-250MG
95-06-7	Sulfallate Solution	S-13243M1-1ML
95-06-7	Sulfallate Solution	S-13243M1-5ML
95-50-1	1,2-Dichlorobenzene	N-10152-1G
95-50-1	1,2-Dichlorobenzene Solution	S-10152M1-1ML
95-50-1	1,2-Dichlorobenzene Solution	S-10152M1-5ML
95-69-2	4-Chloro-2-methylaniline	MET-11426A-250MG
95-69-2	4-Chloro-2-methylaniline Solution	MET-11426AM1-1ML
96-12-8	1,2-Dibromo-3-chloropropane	N-10149-250MG
96-12-8	1,2-Dibromo-3-chloropropane Solution	S-10149M1-1ML
96-12-8	1,2-Dibromo-3-chloropropane Solution	S-10149M1-5ML
96-45-7	2-Imidazolidinethione	N-10381-250MG
96-45-7	2-Imidazolidinethione Solution	S-10381A7-1ML
96-45-7	2-Imidazolidinethione Solution	S-10381A7-5ML
96-45-7	2-Imidazolidinethione Solution	S-10381H1-1ML
96-45-7	2-Imidazolidinethione Solution	S-10381H1-5ML
97-16-5	2,4-Dichlorophenyl benzenesulfonate	N-10628-250MG
97-16-5	2,4-Dichlorophenyl benzenesulfonate Solution	S-10628A1-1ML
97-16-5	2,4-Dichlorophenyl benzenesulfonate Solution	S-10628U1-1ML
97-17-6	Dichlofenthion	N-11659-250MG
97-17-6	Dichlofenthion Solution	S-11659A1-1ML
97-17-6	Dichlofenthion Solution	S-11659J4-1ML
97-17-6	Dichlofenthion Solution	S-11659J4-5ML
97-23-4	Dichlorophen	N-11667-250MG
97-23-4	Dichlorophen Solution	S-11667M4-1ML
97-23-4	Dichlorophen Solution	S-11667M4-5ML
97-53-0	4-Allyl-2-methoxyphenol	N-10803-250MG
97-53-0	4-Allyl-2-methoxyphenol Solution	S-10803M1-1ML
98-28-2	4-tert-Butyl-2-chlorophenol	MET-13186A-10MG
99-30-9	Dicloran	N-11678-250MG
99-30-9	Dicloran Solution	S-11678A1-1ML
99-30-9	Dicloran Solution	S-11678U4-1ML
99-54-7	3,4-Dichloronitrobenzene	N-10785-250MG
99-54-7	3,4-Dichloronitrobenzene Solution	S-10785A1-1ML
99-54-7	3,4-Dichloronitrobenzene Solution	S-10785U1-1ML
100-56-1	Phenyl mercuric chloride(Technical)	N-13013-250MG
100-57-2	Phenyl mercuric hydroxide(Technical)	N-13014-250MG
101-05-3	Anilazine	N-11075-250MG
101-05-3	Anilazine Solution	S-11075A1-1ML
101-05-3	Anilazine Solution	S-11075U1-1ML
101-05-3	Anilazine Solution	S-11075U1-5ML
101-10-0	2-(3-Chlorophenoxy)propionic acid	N-10254-100MG
101-10-0	2-(3-Chlorophenoxy)propionic acid Solution	S-10254A1-1ML
101-10-0	2-(3-Chlorophenoxy)propionic acid Solution	S-10254T1-1ML
101-21-3	Chlorpropham	N-11458-250MG
101-21-3	Chlorpropham Solution	S-11458A1-1ML
101-21-3	Chlorpropham Solution	S-11458T1-1ML
101-21-3	Chlorpropham Solution	S-11458T1-5ML
101-27-9	Barban	N-11123-100MG
101-27-9	Barban Solution	S-11123M1-1ML
101-27-9	Barban Solution	S-11123M1-5ML
101-42-8	Fenuron	N-11967-250MG
101-42-8	Fenuron Solution	S-11967A4-1ML
101-42-8	Fenuron Solution	S-11967A4-5ML
101-42-8	Fenuron Solution	S-11967U1-1ML
102-98-7	Phenylmercury borate(Technical)	N-13031-100MG
103-17-3	Chlorbendside	N-11423-100MG
103-17-3	Chlorbendside Solution	S-11423A1-1ML
103-17-3	Chlorbendside Solution	S-11423U1-1ML
103-26-4	2,3-D methyl ester	N-12861-100MG
103-33-3	Azobenzene	N-11110-500MG
103-33-3	Azobenzene Solution	S-11110M1-1ML
103-33-3	Azobenzene Solution	S-11110M1-5ML
103-85-5	Phenyl-2-thiourea	N-13020-250MG
103-85-5	Phenyl-2-thiourea	S-13020A1-1ML
103-85-5	Phenyl-2-thiourea	S-13020U1-1ML
104-55-2	Cinnamaldehyde	N-11468-500MG
104-55-2	Cinnamaldehyde Solution	S-11468A1-1ML
104-55-2	Cinnamaldehyde Solution	S-11468U1-1ML
104-60-9	Phenylmercuric oleate(Technical)	N-13030-100MG
104-60-9	Phenylmercuric oleate Solution	S-13030A1-1ML
104-60-9	Phenylmercuric oleate Solution	S-13030U1-1ML
106-46-7	1,4-Dichlorobenzene	N-10216-1G
106-46-7	1,4-Dichlorobenzene Solution	S-10216M1-1ML
106-46-7	1,4-Dichlorobenzene Solution	S-10216M1-5ML
106-93-4	1,2-Dibromoethane	N-10150-1G
106-93-4	1,2-Dibromoethane Solution	S-10150M1-1ML
107-27-7	Ethylmercuric chloride(Technical)	N-11940-100MG
107-49-3	Tetraethylpyrophosphate	N-13543-500MG
112-12-9	2-Undecanone	N-10520-1G
112-12-9	2-Undecanone Solution	S-10520N1-1ML
112-56-1	2-[2-Butoxyethoxy]ethyl thiocyanate	N-10260-250MG

CAS#	Description	Part #
112-56-1	2-[2-Butoxyethoxy]ethyl thiocyanate Solution	S-10260A1-1ML
112-56-1	2-[2-Butoxyethoxy]ethyl thiocyanate Solution	S-10260J4-1ML
112-56-1	2-[2-Butoxyethoxy]ethyl thiocyanate Solution	S-10260J4-5ML
113-48-4	MGK 264 (TM)	N-12483-250MG
113-48-4	MGK 264 (TM) Solution	S-12483A1-1ML
113-48-4	MGK 264 (TM) Solution	S-12483T1-1ML
113-48-4	MGK 264 (TM) Solution	S-12483T1-5ML
114-26-1	Propoxur	N-11128-250MG
114-26-1	Propoxur Solution	S-11128A1-1ML
115-29-7	a,b-Endosulfan	N-10992-500MG
115-29-7	a,b-Endosulfan Solution	S-10992A1-1ML
115-29-7	a,b-Endosulfan Solution	S-10992J1-1ML
115-29-7	a,b-Endosulfan Solution	S-10992J1-5ML
115-31-1	Isobornyl thiocanoacetate	N-12232-250MG
115-31-1	Isobornyl thiocanoacetate Solution	S-12232A1-1ML
115-31-1	Isobornyl thiocanoacetate Solution	S-12232U1-1ML
115-32-2	Dicofol	N-12290-100MG
115-78-6	Phosphan	N-13040-50MG
115-78-6	Phosphan Solution	S-13040M1-1ML
115-84-4	2-Butyl-2-ethyl-1,3-propanediol	N-10306-1G
115-84-4	2-Butyl-2-ethyl-1,3-propanediol Solution	S-10306M1-1ML
115-90-2	Fensulfthion	N-11963-100MG
115-90-2	Fensulfthion Solution	S-11963M1-1ML
115-90-2	Fensulfthion Solution	S-11963M1-5ML
115-93-5	Cythioate	N-11551-50MG
115-93-5	Cythioate Solution	S-11551A1-1ML
115-93-5	Cythioate Solution	S-11551U1-1ML
116-06-3	Aldicarb	N-11044-100MG
116-06-3	Aldicarb Solution	S-11044M1-1ML
116-16-5	Hexachloroacetone	N-12158-500MG
116-16-5	Hexachloroacetone Solution	S-12158A1-1ML
116-16-5	Hexachloroacetone Solution	S-12158J1-1ML
116-29-0	Tetradifon	N-13539-250MG
116-29-0	Tetradifon Solution	S-13539A1-1ML
116-29-0	Tetradifon Solution	S-13539B3-1ML
116-29-0	Tetradifon Solution	S-13539B3-5ML
117-18-0	1,2,4,5-Tetrachloro-3-nitrobenzene	N-10111-250MG
117-18-0	1,2,4,5-Tetrachloro-3-nitrobenzene Solution	S-10111M1-1ML
117-80-6	Dichlone	N-11661-250MG
117-80-6	Dichlone Solution	S-11661M1-1ML
117-80-6	Dichlone Solution	S-11661U1-1ML
117-80-6	Dichlone Solution	S-11661U1-5ML
118-74-1	Hexachlorobenzene	N-12159-250MG
118-74-1	Hexachlorobenzene Solution	S-12159K4-1ML
118-74-1	Hexachlorobenzene Solution	S-12159M1-1ML
118-74-1	Hexachlorobenzene Solution	S-12159M1-5ML
118-75-2	Chloranil	N-11420-250MG
118-75-2	Chloranil Solution	S-11420M1-1ML
118-75-2	Chloranil Solution	S-11420U1-1ML
119-12-0	Pyridaphenthion	N-13154-250MG
119-12-0	Pyridaphenthion Solution	S-13154A1-1ML
119-12-0	Pyridaphenthion Solution	S-13154U1-1ML
120-23-0	2-Naphthoxyacetic acid	N-10437-250MG
120-23-0	2-Naphthoxyacetic acid Solution	S-10437M1-1ML
120-23-0	2-Naphthoxyacetic acid Solution	S-10437T1-1ML
120-32-1	o-Benzyl-p-chlorophenol	N-12676-100MG
120-32-1	o-Benzyl-p-chlorophenol Solution	S-12676M1-1ML
120-36-5	Dichlorprop	N-11671-250MG
120-36-5	Dichlorprop Solution	S-11671B1-1ML
120-36-5	Dichlorprop Solution	S-11671B1-5ML
120-36-5	Dichlorprop Solution	S-11671M1-1ML
120-51-4	Benzyl benzoate	N-11182-100MG
120-51-4	Benzyl benzoate Solution	S-11182A1-1ML
120-51-4	Benzyl benzoate Solution	S-11182J5-1ML
120-51-4	Benzyl benzoate Solution	S-11182J7-1ML
120-51-4	Benzyl benzoate Solution	S-11182J7-5ML
120-62-7	Sulfoxide	N-13257-100MG
120-62-7	Sulfoxide Solution	S-13257M1-1ML
120-62-7	Sulfoxide Solution	S-13257M1-5ML
121-57-3	Sulfanilic acid	N-13249-250MG
121-75-5	Malathion	N-12346-100MG
121-75-5	Malathion Solution	S-12346A1-1ML
121-75-5	Malathion Solution	S-12346U1-1ML
121-75-5	Malathion Solution	S-12346U1-5ML
122-11-2	Sulfadimethoxine	N-13242-250MG
122-11-2	Sulfadimethoxine Solution	S-13242A1-1ML
122-14-5	Fenitrothion	N-11955-250MG
122-14-5	Fenitrothion Solution	S-11955A1-1ML
122-14-5	Fenitrothion Solution	S-11955J4-1ML
122-34-9	Simazine	N-13800-100MG
122-34-9	Simazine Solution	S-13800B4-1ML
122-34-9	Simazine Solution	S-13800B4-5ML
122-34-9	Simazine Solution	S-13800M1-1ML
122-34-9	Simazine Solution	S-13800M1-5ML
122-42-9	Propham	N-13109-250MG
122-42-9	Propham Solution	S-13109A4-1ML
122-59-8	Phenoxyacetic acid	N-13004-250MG
122-59-8	Phenoxyacetic acid Solution	S-13004M1-1ML
122-59-8	Phenoxyacetic acid Solution	S-13004T1-1ML
122-64-5	Phenyl mercuric lactate	N-13015-250MG
122-88-3	p-Chlorophenoxyacetic acid	N-12763-250MG
122-88-3	p-Chlorophenoxyacetic acid Solution	S-12763M1-1ML
122-88-3	p-Chlorophenoxyacetic acid Solution	S-12763T1-1ML
123-33-1	Maleic hydrazide	N-12350-250MG
123-33-1	Maleic hydrazide Solution	S-12350M1-1ML
124-58-3	Methylarsonic acid	N-12922-100MG
126-75-0	Demeton S	N-11582-100MG
126-75-0	Demeton S Solution	S-11582A1-1ML
126-75-0	Demeton S Solution	S-11582U1-1ML
127-20-8	Dalapon Sodium	N-11564-250MG
127-71-9	Sulfabenzamide	N-13240-250MG
127-79-7	Sulfamerazine	N-13244-250MG
127-90-2	S 421	N-13919-10MG
128-03-0	Busan 85	N-11352-100MG
128-03-0	Busan 85 Solution	S-11352A1-1ML
128-03-0	Busan 85 Solution	S-11352F1-1ML
128-03-0	Busan 85 Solution	S-11352M1-1ML
128-04-1	Dibam	N-12955-10MG
131-11-3	Dimethyl phthalate	N-11770-500MG
131-11-3	Dimethyl phthalate Solution	S-11770H1-1ML
131-11-3	Dimethyl phthalate Solution	S-11770H1-5ML
131-11-3	Dimethyl phthalate Solution	S-11770M1-1ML
131-11-3	Dimethyl phthalate Solution	S-11770M1-5ML
131-72-6	Meptyldinocap	N-13288-25MG

CAS#	Description	Part #	CAS#	Description	Part #
131-72-6	Meptyldinocap Solution	S-13288M1-1ML	327-98-0	Trichloronate Solution	S-13657L1-1ML
132-66-1	N-1-Naphthylphthalamic acid	N-12507-250MG	327-98-0	Trichloronate Solution	S-13657L1-5ML
132-66-1	N-1-Naphthylphthalamic acid Solution	S-12507A1-1ML	330-54-1	Diuron	N-11827-250MG
132-66-1	N-1-Naphthylphthalamic acid Solution	S-12507T1-1ML	330-54-1	Diuron Solution	S-11827A4-1ML
133-06-2	Captan	N-11400-250MG	330-54-1	Diuron Solution	S-11827U1-1ML
133-06-2	Captan Solution	S-11400U1-1ML	330-55-2	Linuron	N-12322-250MG
133-06-2	Captan Solution	S-11400U1-5ML	330-55-2	Linuron Solution	S-12322A6-1ML
133-07-3	Folpet	N-12007-250MG	330-55-2	Linuron Solution	S-12322A6-5ML
133-07-3	Folpet Solution	S-12007M1-1ML	330-95-0	Nicarbazin	N-12651-100MG
133-07-3	Folpet Solution	S-12007U1-1ML	333-41-5	Diazinon	N-11621-250MG
133-32-4	3-Indolebutyric acid	N-10728-100MG	333-41-5	Diazinon Solution	S-11621U1-1ML
133-32-4	3-Indolebutyric acid Solution	S-10728A1-1ML	338-45-4	Z-Mevinphos	N-12906-50MG
133-32-4	3-Indolebutyric acid Solution	S-10728T1-1ML	344-96-7	Methyl o-trifluoromethylbenzoate	MET-12004A-100MG
133-49-3	Pentachlorothiophenol	MET-12159B-100MG	404-86-4	Capsaicin, Natural	N-11398-100MG
133-90-4	Chloramben	N-11418-100MG	404-86-4	Capsaicin, Natural Solution	S-11398A1-1ML
133-90-4	Chloramben Solution	S-11418A1-1ML	404-86-4	Capsaicin, Natural Solution	S-11398U1-1ML
133-90-4	Chloramben Solution	S-11418B1-1ML	465-73-6	Isodrin	N-12248-100MG
133-90-4	Chloramben Solution	S-11418B1-5ML	465-73-6	Isodrin Solution	S-12248M1-1ML
134-20-3	Methyl anthranilate	N-12414-100MG	465-73-6	Isodrin Solution	S-12248M1-5ML
134-20-3	Methyl anthranilate Solution	S-12414M1-1ML	467-69-6	Flurenol	N-13815-250MG
134-62-3	N,N-Diethyl-m-toluamide	N-12618-250MG	470-90-6	Chlorfenvinphos	N-11429-250MG
134-62-3	N,N-Diethyl-m-toluamide Solution	S-12618A1-1ML	470-90-6	Chlorfenvinphos Solution	S-11429K1-1ML
134-62-3	N,N-Diethyl-m-toluamide Solution	S-12618B5-1ML	485-31-4	Binapacryl	N-11204-100MG
134-62-3	N,N-Diethyl-m-toluamide Solution	S-12618B5-5ML	485-31-4	Binapacryl Solution	S-11204A1-1ML
136-45-8	Di-n-propyl isocinchomerate	N-12484-250MG	485-31-4	Binapacryl Solution	S-11204T1-1ML
137-30-4	Ziram	N-13761-250MG	492-86-4	p-Chloromandelic acid	N-12761-250MG
139-40-2	Propazine	N-13108-100MG	492-86-4	p-Chloromandelic acid Solution	S-12761A1-1ML
139-40-2	Propazine Solution	S-13108A1-1ML	492-86-4	p-Chloromandelic acid Solution	S-12761T1-1ML
139-40-2	Propazine Solution	S-13108T1-1ML	502-55-6	Diethyl dithiobis(thionoformate)	N-11696-250MG
139-40-2	Propazine Solution	S-13108T1-5ML	502-55-6	Diethyl dithiobis(thionoformate) Solution	S-11696A1-1ML
140-38-5	p-Chlorophenylurea	MET-11722A-100MG	502-55-6	Diethyl dithiobis(thionoformate) Solution	S-11696K1-1ML
140-41-0	Monuron TCA	N-12498-100MG	504-24-5	4-Aminopyridine	N-10806-1G
140-56-7	Fenaminosulf	N-11586-250MG	504-24-5	4-Aminopyridine Solution	S-10806M1-1ML
140-56-7	Fenaminosulf Solution	S-11586M1-1ML	510-15-6	Chlorobenzilate	N-11442-100MG
140-57-8	Aramite	N-11085-10MG	510-15-6	Chlorobenzilate Solution	S-11442A1-1ML
140-57-8	Aramite Solution	S-11085J1-1ML	510-15-6	Chlorobenzilate Solution	S-11442J1-1ML
141-00-4	Cadmium succinate(Technical)	N-11391-250MG	510-15-6	Chlorobenzilate Solution	S-11442J1-5ML
141-03-7	Dibutyl succinate	N-11644-500MG	525-79-0	Kinetin	N-12294-250MG
141-03-7	Dibutyl succinate Solution	S-11644A1-1ML	527-20-8	Pentachloroaniline	MET-12383A-100MG
141-03-7	Dibutyl succinate Solution	S-11644U1-1ML	532-34-3	Butylmesityl oxide oxalate	N-11385-250MG
141-66-2	Dicrotophos	N-11680-100MG	533-23-3	2,4-D ethyl ester	N-10614-250MG
141-66-2	Dicrotophos Solution	S-11680U1-1ML	533-23-3	2,4-D ethyl ester Solution	S-10614A1-1ML
141-66-2	Dicrotophos Solution	S-11680U1-5ML	533-23-3	2,4-D ethyl ester Solution	S-10614J1-1ML
142-59-6	Nabam	N-12639-250MG	533-74-4	Dazomet	N-11566-250MG
142-59-6	Nabam Solution	S-12639F1-1ML	533-74-4	Dazomet Solution	S-11566M1-1ML
142-59-6	Nabam Solution	S-12639U1-1ML	534-52-1	4,6-Dinitro-o-cresol (contains ~10% water)	N-10889-250MG
143-50-0	Chlordecone	N-12291-100MG	534-52-1	4,6-Dinitro-o-cresol Solution	S-10889M4-1ML
143-50-0	Chlordecone Solution	S-12291A1-1ML	539-47-9	2-Furanacrylic acid	N-10363-100MG
143-50-0	Chlordecone Solution	S-12291J1-1ML	539-47-9	2-Furanacrylic acid Solution	S-10363A1-1ML
143-50-0	Chlordecone Solution	S-12291J1-5ML	539-47-9	2-Furanacrylic acid Solution	S-10363T1-1ML
143-50-0	Chlordecone Solution	S-12291T4-1ML	542-75-6	1,3-Dichloropropene	N-10193-500MG
143-50-0	Chlordecone Solution	S-12291T4-5ML	542-75-6	1,3-Dichloropropene Solution	S-10193M1-1ML
144-21-8	Disodium methyl arsonate hexahydrate	N-11817-500MG	554-00-7	2,4-Dichloroaniline	MET-11521A-1G
144-21-8	Disodium methyl arsonate hexahydrate Solution	S-11817F1-1ML	555-37-3	Neburon	N-12647-250MG
144-82-1	Sulfamethizole	N-13246-250MG	555-37-3	Neburon Solution	S-12647A4-1ML
144-83-2	Sulfapyridine	N-13250-250MG	555-37-3	Neburon Solution	S-12647A4-5ML
148-24-3	8-Hydroxyquinoline	N-10965-250MG	556-22-9	Glyodin(Technical)	N-12131-250MG
148-24-3	8-Hydroxyquinoline Solution	S-10965M1-1ML	556-61-6	Methyl isothiocyanate	MET-12392A-250MG
148-79-8	Thiabendazole	N-13560-250MG	557-30-2	Ethanedial dioxime	N-11867-100MG
148-79-8	Thiabendazole Solution	S-13560M4-1ML	557-30-2	Ethanedial dioxime	N-11867-250MG
150-50-5	Tributylphosphoro-trithioite	N-13648-250MG	557-30-2	Ethanedial dioxime Solution	S-11867M1-1ML
150-68-5	Monuron	N-12497-250MG	563-12-2	Ethion	N-11873-250MG
150-68-5	Monuron Solution	S-12497A4-1ML	563-12-2	Ethion Solution	S-11873M1-1ML
150-68-5	Monuron Solution	S-12497M4-1ML	563-12-2	Ethion Solution	S-11873M1-5ML
152-16-9	Octamethylpyrophosphoramide	N-12715-50MG	575-36-0	N-1-Naphthylacetamide	N-12506-500MG
288-88-0	1,2,4-Triazole	MET-13576A-250MG	575-36-0	N-1-Naphthylacetamide Solution	S-12506M1-1ML
297-78-9	Isobenzan	N-12231-10MG	575-90-6	2,6-D Acid	N-12996-10MG
297-78-9	Isobenzan Solution	S-12231M1-1ML	579-66-8	2,6-Diethylaniline	MET-11043C-500MG
297-97-2	Zinophos (TM)	N-13760-50MG	584-79-2	Allethrin	N-11051-100MG
297-97-2	Zinophos (TM) Solution	S-13760A1-1ML	584-79-2	Allethrin Solution	S-11051X1-1ML
297-97-2	Zinophos (TM) Solution	S-13760J1-1ML	584-79-2	Allethrin Solution	S-11051X1-5ML
297-97-2	Zinophos (TM) Solution	S-13760J1-5ML	590-28-3	Potassium cyanate(Technical)	N-13087-250MG
298-00-0	Methyl parathion	N-12452-100MG	594-42-3	Perchloromethyl mercaptan	N-12845-250MG
298-00-0	Methyl parathion Solution	S-12452A1-1ML	594-42-3	Perchloromethyl mercaptan Solution	S-12845A1-1ML
298-00-0	Methyl parathion Solution	S-12452U1-1ML	594-42-3	Perchloromethyl mercaptan Solution	S-12845U1-1ML
298-02-2	Phorate	N-13035-100MG	594-72-9	1,1-Dichloro-1-nitroethane	N-10122-1G
298-02-2	Phorate Solution	S-13035A1-1ML	594-72-9	1,1-Dichloro-1-nitroethane Solution	S-10122M1-1ML
298-02-2	Phorate Solution	S-13035J1-1ML	594-72-9	1,1-Dichloro-1-nitroethane Solution	S-10122X1-1ML
298-03-3	Demeton O	N-11580-100MG	598-02-7	Diethyl phosphate	MET-11621C-100MG
298-03-3	Demeton O Solution	S-11580U1-1ML	598-21-0	Bromoacetyl bromide	N-11333-500MG
298-04-4	Disulfoton	N-11819-250MG	598-21-0	Bromoacetyl bromide Solution	S-11333A1-1ML
298-04-4	Disulfoton Solution	S-11819M1-1ML	598-21-0	Bromoacetyl bromide Solution	S-11333U1-1ML
299-84-3	Fenchlorphos	N-11952-100MG	608-73-1	BHC (mixed isomers)	N-11198-250MG
299-84-3	Fenchlorphos Solution	S-11952U1-1ML	608-73-1	BHC (mixed isomers) Solution	S-11198A1-1ML
299-84-3	Fenchlorphos Solution	S-11952U1-5ML	608-73-1	BHC (mixed isomers) Solution	S-11198U1-1ML
299-86-5	Ruelene (TM)	N-13186-50MG	611-06-3	2,4-Dichloronitrobenzene	N-10625-250MG
299-86-5	Ruelene (TM) Solution	S-13186A1-1ML	611-06-3	2,4-Dichloronitrobenzene Solution	S-10625A1-1ML
299-86-5	Ruelene (TM) Solution	S-13186U1-1ML	611-06-3	2,4-Dichloronitrobenzene Solution	S-10625U1-1ML
300-76-5	Naled	N-12640-250MG	614-61-9	o-Chlorophenoxy acetic acid	N-12679-250MG
300-76-5	Naled Solution	S-12640U1-1ML	614-61-9	o-Chlorophenoxy acetic acid Solution	S-12679M1-1ML
301-12-2	Oxydemeton-methyl	N-12741-50MG	614-61-9	o-Chlorophenoxy acetic acid Solution	S-12679T1-1ML
309-00-2	Aldrin	N-11049-250MG	618-62-2	3,5-Dichloronitrobenzene	N-10792-250MG
309-00-2	Aldrin Solution	S-11049M1-1ML	618-62-2	3,5-Dichloronitrobenzene Solution	S-10792A1-1ML
309-00-2	Aldrin Solution	S-11049M1-5ML	618-62-2	3,5-Dichloronitrobenzene Solution	S-10792U1-1ML
311-45-5	Paraoxon	N-12816-100MG	619-08-9	2-Chloro-4-nitrophenol	MET-11658A-250MG
311-45-5	Paraoxon Solution	S-12816A1-1ML	621-00-1	1,3-di-p-tolylurea	N-14078-100MG
314-40-9	Bromacil	S-12816U1-1ML	626-43-7	3,5-Dichloroaniline	MET-12220B-250MG
314-40-9	Bromacil Solution	N-11330-250MG	626-67-5	N-Methylpiperidine	MET-12371A-250MG
314-40-9	Bromacil Solution	S-11330A1-1ML	634-83-3	2,3,4,5-Tetrachloroaniline Solution	S-13990A1-1ML
315-18-4	Zectran	S-11330T1-1ML	640-15-3	Thiometon	N-13570-50MG
315-18-4	Zectran Solution	N-13753-500MG	640-15-3	Thiometon Solution	S-13570U1-1ML
315-18-4	Zectran Solution	S-13753A4-1ML	640-19-7	2-Fluoroacetamide	N-10358-100MG
319-84-6	BHC (alpha isomer)	S-11194-100MG	645-78-3	2-Fluoroacetamide Solution	S-10358M1-1ML
319-84-6	BHC (alpha isomer) Solution	S-11194U1-1ML	645-78-3	Monothiono TEPP	N-12496-50MG
319-85-7	BHC (beta isomer)	N-11195-100MG	645-92-1	Monothiono TEPP Solution	S-12496U1-1ML
319-85-7	BHC (beta isomer) Solution	S-11195U1-1ML	671-04-5	Carbanolate	MET-11106E-50MG
319-85-7	BHC (beta isomer) Solution	S-11195U1-5ML	671-04-5	Carbanolate Solution	N-11122-100MG
319-86-8	BHC (delta isomer)	N-11196-50MG	671-04-5	Carbanolate Solution	S-11122A1-1ML
319-86-8	BHC (delta isomer) Solution	S-11196U1-1ML	672-99-1	BDMC	S-11122U1-1ML
319-86-8	BHC (delta isomer) Solution	S-11196U1-5ML	672-99-1	BDMC Solution	N-11131-100MG
321-54-0	Coumaphos-O-analog	MET-11507D-50MG	672-99-1	BDMC Solution	S-11131M1-1ML
327-98-0	Trichloronate	N-13657-100MG	673-04-1	Simeton	S-11131M1-5ML
327-98-0	Trichloronate Solution	S-13657A1-1ML	709-98-8	Propanil	N-12960-50MG
					N-10800-500MG

CAS#	Description	Part #	CAS#	Description	Part #
709-98-8	Propanil Solution	S-10800A1-1ML	1194-65-6	2,6-Dichlorobenzonitrile Solution	S-10688M1-1ML
709-98-8	Propanil Solution	S-10800Z4-1ML	1214-39-7	6-Benzylaminopurine	N-10956-250MG
709-98-8	Propanil Solution	S-10800Z4-5ML	1214-39-7	6-Benzylaminopurine Solution	N-10956A1-1ML
723-46-6	Sulfamethoxazole	N-13247-250MG	1216-44-0	Flurenol-methyl	N-13813-250MG
731-27-1	Tolyfluamid	N-13584-250MG	1219-99-4	N,N'-Bis(p-chlorophenyl)urea	MET-12497A-500MG
732-11-6	Phosmet	N-12207-100MG	1332-65-6	Copper oxychloride(Technical)	N-11502-250MG
732-11-6	Phosmet Solution	S-12207A1-1ML	1338-02-9	Copper Naphthenates(Technical)	N-11501-250MG
732-11-6	Phosmet Solution	S-12207J1-1ML	1405-10-3	Neomycin trisulfate hydrate(Technical)	N-12649-250MG
732-11-6	Phosmet Solution	S-12207J1-5ML	1441-02-7	Pentachlorophenyl acetate	N-13824-100MG
741-58-2	Bensulide	N-11140-250MG	1469-48-3	cis-1,2,3,6-Tetrahydrophthalimide	MET-11399A-250MG
741-58-2	Bensulide Solution	S-11140M4-1ML	1563-38-8	2,3-Dihydro-2,2-dimethylbenzofuran-7-ol	MET-11405B-250MG
741-58-2	Bensulide Solution	S-11140M4-5ML	1563-66-2	Carbofuran	N-11405-250MG
759-94-4	S-Ethyl dipropylthiocarbamate	N-13188-500MG	1563-66-2	Carbofuran Solution	S-11405M1-1ML
759-94-4	S-Ethyl dipropylthiocarbamate Solution	S-13188A1-1ML	1570-64-5	4-Chloro-2-methylphenol	MET-10818A-1G
759-94-4	S-Ethyl dipropylthiocarbamate Solution	S-13188T1-1ML	1582-09-8	Trifluralin	N-13689-250MG
759-94-4	S-Ethyl dipropylthiocarbamate Solution	S-13188T1-5ML	1582-09-8	Trifluralin Solution	S-13689A1-1ML
786-19-6	Carbophenothion	N-11408-250MG	1582-09-8	Trifluralin Solution	S-13689U1-1ML
786-19-6	Carbophenothion Solution	S-11408M1-1ML	1582-09-8	Trifluralin Solution	S-13689U1-5ML
789-02-6	o,p'-DDT	N-12708-50MG	1593-77-7	Dodermorph	N-11839-250MG
789-02-6	o,p'-DDT Solution	S-12708M1-1ML	1593-77-7	Dodermorph Solution	S-11839M1-1ML
789-02-6	o,p'-DDT Solution	S-12708M1-5ML	1596-84-5	Daminozide	N-11565-250MG
830-96-6	3-Indolepropionic acid	N-10729-250MG	1596-84-5	Daminozide Solution	S-11565A1-1ML
830-96-6	3-Indolepropionic acid Solution	S-10729A1-1ML	1596-84-5	Daminozide Solution	S-11565T1-1ML
830-96-6	3-Indolepropionic acid Solution	S-10729T1-1ML	1610-17-9	Atraton	N-11105-250MG
841-06-5	Methoprotyrine	N-12401-100MG	1610-17-9	Atraton Solution	S-11105A1-1ML
841-06-5	Methoprotyrine Solution	S-12401A1-1ML	1610-17-9	Atraton Solution	S-11105T1-1ML
841-06-5	Methoprotyrine Solution	S-12401T1-1ML	1610-17-9	Atraton Solution	S-11105T1-5ML
867-27-6	Demeton-O-methyl	N-11583-50MG	1610-18-0	Prometon	N-13102-250MG
873-32-5	2-Chlorobenzonitrile	MET-11490A-250MG	1610-18-0	Prometon Solution	S-13102A1-1ML
886-50-0	Terbutyrine	N-13514-250MG	1610-18-0	Prometon Solution	S-13102T1-1ML
886-50-0	Terbutyrine Solution	S-13514A1-1ML	1610-18-0	Prometon Solution	S-13102T1-5ML
886-50-0	Terbutyrine Solution	S-13514T1-1ML	1634-78-2	Malathion-O-analog	MET-12346C-100MG
886-50-0	Terbutyrine Solution	S-13514T1-5ML	1646-75-9	Aldicarb oxime	MET-11044A-10MG
887-54-7	Monomethyl tetrachloroterephthalate	MET-11462A-100MG	1646-87-3	Aldicarb sulfoxide	N-11046-10MG
900-95-8	Phentoinacetate	N-13007-250MG	1646-87-3	Aldicarb sulfoxide Solution	S-11046A1-1ML
900-95-8	Phentoinacetate Solution	S-13007A1-1ML	1646-88-4	Aldicarb sulfone	N-11045-10MG
900-95-8	Phentoinacetate Solution	S-13007U1-1ML	1689-83-4	Aldicarb sulfone Solution	S-11045M1-1ML
919-86-8	Metasystox I (TM)	N-12384-50MG	1689-83-4	3,5-Diiodo-4-hydroxybenzonitrile	N-10794-100MG
919-86-8	Metasystox I (TM) Solution	S-12384A1-1ML	1689-83-4	3,5-Diiodo-4-hydroxybenzonitrile Solution	S-10794M1-1ML
919-86-8	Metasystox I (TM) Solution	S-12384U1-1ML	1689-84-5	Bromoxynil	N-11345-250MG
931-17-9	1,2-Cyclohexanediol (cis & trans mixture)	MET-12727B-250MG	1689-84-5	Bromoxynil Solution	S-11345M1-1ML
934-32-7	2-Aminobenzimidazole	MET-11138A-250MG	1689-99-2	Bromoxynil octanoic acid ester	N-11346-250MG
939-05-9	3-Hydroxy-5-phenylisoxazole	MET-13932A-10MG	1689-99-2	Bromoxynil octanoic acid ester Solution	S-11346A1-1ML
940-31-8	2-Phenoxypropionic acid	N-10507-500MG	1689-99-2	Bromoxynil octanoic acid ester Solution	S-11346K4-1ML
940-31-8	2-Phenoxypropionic acid Solution	S-10507M1-1ML	1689-99-2	Bromoxynil octanoic acid ester Solution	S-11346K4-5ML
940-31-8	2-Phenoxypropionic acid Solution	S-10507T1-1ML	1698-60-8	Pyrazon	N-13145-250MG
944-22-9	Fonofos	N-11842-100MG	1698-60-8	Pyrazon Solution	S-13145M1-1ML
944-22-9	Fonofos Solution	S-11842A1-1ML	1702-17-6	Clopyralid	N-12324-250MG
944-22-9	Fonofos Solution	S-11842J4-1ML	1702-17-6	Clopyralid Solution	S-12324A1-1ML
944-22-9	Fonofos Solution	S-11842J4-5ML	1702-17-6	Clopyralid Solution	S-12324T1-1ML
947-02-4	Phosfolan	N-13038-100MG	1713-15-1	2,4-D isobutyl ester	N-10532-250MG
947-02-4	Phosfolan Solution	S-13038A1-1ML	1713-15-1	2,4-D isobutyl ester Solution	S-10532A1-1ML
947-02-4	Phosfolan Solution	S-13038U1-1ML	1713-15-1	2,4-D isobutyl ester Solution	S-10532U1-1ML
950-10-7	Mephosfolan	N-12370-100MG	1716-09-2	Fenthion-ethyl	N-11965-500MG
950-10-7	Mephosfolan Solution	S-12370A1-1ML	1716-09-2	Fenthion-ethyl Solution	S-11965A1-1ML
950-10-7	Mephosfolan Solution	S-12370U1-1ML	1716-09-2	Fenthion-ethyl Solution	S-11965U1-1ML
950-35-6	Dimethyl-p-nitrophenylphosphate	N-11775-100MG	1746-81-2	Monolinuron	N-12494-250MG
950-35-6	Dimethyl-p-nitrophenylphosphate Solution	S-11775A1-1ML	1746-81-2	Monolinuron Solution	S-12494A1-1ML
950-35-6	Dimethyl-p-nitrophenylphosphate Solution	S-11775T1-1ML	1746-81-2	Monolinuron Solution	S-12494U1-1ML
950-35-6	Dimethyl-p-nitrophenylphosphate Solution	S-11775T1-5ML	1757-18-2	Akton (TM)	N-11042-100MG
950-37-8	Methidathion	N-12397-250MG	1757-18-2	Akton (TM) Solution	S-11042A1-1ML
950-37-8	Methidathion Solution	S-12397A1-1ML	1757-18-2	Akton (TM) Solution	S-11042U1-1ML
950-37-8	Methidathion Solution	S-12397U1-1ML	1763-23-1	Perfluorooctanesulfonic acid Solution	S-14002M1-1ML
953-17-3	Carbophenothion methyl Solution	MET-11408A1-1ML	1824-09-5	Mesoprazine	N-13296-50MG
957-51-7	Diphenamid	N-11794-250MG	1824-09-5	Mesoprazine Solution	S-13296A1-1ML
957-51-7	Diphenamid Solution	S-11794A1-1ML	1825-21-4	Pentachloroanisole	N-12826-100MG
957-51-7	Diphenamid Solution	S-11794T1-1ML	1825-21-4	Pentachloroanisole Solution	S-12826M1-1ML
957-51-7	Diphenamid Solution	S-11794T1-5ML	1836-75-5	Nitrofen	N-12663-100MG
959-98-8	a-Endosulfan	N-10979-100MG	1836-75-5	Nitrofen Solution	S-12663M1-1ML
959-98-8	a-Endosulfan Solution	S-10979A1-1ML	1861-32-1	Chlorthal-dimethyl	N-11462-250MG
959-98-8	a-Endosulfan Solution	S-10979J1-1ML	1861-32-1	Chlorthal-dimethyl Solution	S-11462B1-1ML
959-98-8	a-Endosulfan Solution	S-10979J1-5ML	1861-32-1	Chlorthal-dimethyl Solution	S-11462B1-5ML
959-98-8	a-Endosulfan Solution	S-10979T1-1ML	1861-32-1	Chlorthal-dimethyl Solution	S-11462M1-1ML
959-98-8	a-Endosulfan Solution	S-10979T1-5ML	1861-40-1	Benfluralin	N-11136-250MG
960-25-8	Famphur-O-analog	MET-11944A-50MG	1861-40-1	Benfluralin Solution	S-11136A1-1ML
961-22-8	Azinphos-methyl oxon Solution	MET-12139AU1-1ML	1861-40-1	Benfluralin Solution	S-11136J4-1ML
962-58-3	Diazinon-O-analog	MET-11621A-100MG	1861-40-1	Benfluralin Solution	S-11136J4-5ML
992-20-1	Salannin	N-12911-5MG	1878-91-7	p-Bromophenoxy acetic acid	N-12757-250MG
999-81-5	Chlormequat chloride	N-11435-250MG	1878-91-7	p-Bromophenoxy acetic acid Solution	S-12757A1-1ML
999-81-5	Chlormequat chloride Solution	S-11435F1-1ML	1878-91-7	p-Bromophenoxy acetic acid Solution	S-12757T1-1ML
1007-28-9	Atrazine desisopropyl	MET-13213A-50MG	1891-95-8	Chloroxynil	N-11454-500MG
1014-69-3	Desmetryn	N-12973-10MG	1891-95-8	Chloroxynil Solution	S-11457M1-1ML
1014-70-6	Simetryn	N-13802-50MG	1897-45-6	Chlorothalonil	N-11454-250MG
1014-70-6	Simetryn Solution	S-13802A1-1ML	1897-45-6	Chlorothalonil Solution	S-11454A1-1ML
1014-70-6	Simetryn Solution	S-13802T1-1ML	1897-45-6	Chlorothalonil Solution	S-11454T1-1ML
1014-70-6	Simetryn Solution	S-13802T1-5ML	1910-42-5	Paraquat dichloride hydrate	N-12818-500MG
1022-22-6	p,p'-DDMU	MET-12810B-10MG	1910-42-5	Paraquat dichloride Solution	S-12818F1-1ML
1024-57-3	Heptachlor epoxide (Isomer B)	N-12148-50MG	1910-42-5	Paraquat dichloride hydrate Solution	S-12818F17-1ML
1024-57-3	Heptachlor epoxide (Isomer B) Solution	S-12148M1-1ML	1910-42-5	Paraquat dichloride hydrate Solution	S-12818F17-5ML
1024-57-3	Heptachlor epoxide (Isomer B) Solution	S-12148M1-5ML	1912-24-9	Atrazine	N-11106-250MG
1031-07-8	Endosulfan sulfate	N-11851-50MG	1912-26-1	Trietazine	N-13666-100MG
1031-07-8	Endosulfan sulfate Solution	S-11851A1-1ML	1912-26-1	Trietazine Solution	S-13666A1-1ML
1031-07-8	Endosulfan sulfate Solution	S-11851J1-1ML	1912-26-1	Trietazine Solution	S-13666X1-1ML
1031-07-8	Endosulfan sulfate Solution	S-11851J1-5ML	1918-00-9	Dicamba	N-11656-250MG
1066-51-9	Aminomethyl phosphonic acid	MET-12133A-100MG	1918-00-9	Dicamba Solution	S-11656A1-1ML
1071-83-6	Glyphosate	N-12133-250MG	1918-00-9	Dicamba Solution	S-11656B1-1ML
1076-38-6	4-Hydroxycoumarin	MET-13750B-250MG	1918-00-9	Dicamba Solution	S-11656B1-5ML
1085-98-9	Dichlofluanid	N-11660-250MG	1918-02-1	Picloram	N-13050-250MG
1085-98-9	Dichlofluanid Solution	S-11660U1-1ML	1918-02-1	Picloram Solution	S-13050A1-1ML
1113-01-5	(Dimethoxythiophosphorylthio)acetic acid	MET-12017A-100MG	1918-02-1	Picloram Solution	S-13050A4-1ML
1113-02-6	Omethoate	N-12726-100MG	1918-02-1	Picloram Solution	S-13050A4-5ML
1113-02-6	Omethoate Solution	S-12726U1-1ML	1918-02-1	Picloram Solution	S-13050B1-1ML
1114-71-2	S-Propyl butylethylthiocarbamate Solution	S-13192A1-1ML	1918-02-1	Picloram Solution	S-13050B1-5ML
1114-71-2	S-Propyl butylethylthiocarbamate Solution	S-13192T1-1ML	1918-11-2	Terbutol	N-13513-100MG
1129-41-5	Metolcarb	N-12479-100MG	1918-13-4	Chlorthiamid	N-12992-10MG
1129-41-5	Metolcarb Solution	S-12479A1-1ML	1918-16-7	Propachlor	N-13104-250MG
1129-41-5	Metolcarb Solution	S-12479U1-1ML	1918-16-7	Propachlor Solution	S-13104A1-1ML
1134-23-2	Cycloate	N-11522-250MG	1918-16-7	Propachlor Solution	S-13104T1-1ML
1134-23-2	Cycloate Solution	S-11522A1-1ML	1918-16-7	Propachlor Solution	S-13104T1-5ML
1134-23-2	Cycloate Solution	S-11522T1-1ML	1918-18-9	SWEP	N-13260-250MG
1134-23-2	Cycloate Solution	S-11522T1-5ML	1918-18-9	SWEP Solution	S-13260A4-1ML
1166-46-7	(1R)-trans-Tetramethrin	N-13994-100MG	1918-18-9	SWEP Solution	S-13260A4-5ML
1166-46-7	(1R)-trans-Tetramethrin Solution	S-13994A1-1ML	1918-18-9	SWEP Solution	S-13260U1-1ML
1194-65-6	2,6-Dichlorobenzonitrile	N-10688-250MG	1928-37-6	(2,4,5-Trichlorophenoxy)acetic acid methyl ester	N-10005-250MG

CAS#	Description	Part #	CAS#	Description	Part #
1928-37-6	(2,4,5-Trichlorophenoxy)acetic acid methyl ester	S-10005M1-1ML	2524-04-1	O,O-Diethyl phosphorochloridothioate	N-12703-250MG
1928-37-6	(2,4,5-Trichlorophenoxy)acetic acid methyl ester	S-10005M1-5ML	2524-04-1	O,O-Diethyl phosphorochloridothioate Solution	S-12703A1-1ML
1928-38-7	2,4-D methyl ester	N-10619-250MG	2524-04-1	O,O-Diethyl phosphorochloridothioate Solution	S-12703U1-1ML
1928-38-7	2,4-D methyl ester Solution	S-10619M1-1ML	2536-31-4	Chlorfluorecol-methyl ester	N-11431-100MG
1928-38-7	2,4-D methyl ester Solution	S-10619M1-5ML	2536-31-4	Chlorfluorecol-methyl ester Solution	S-11431M1-1ML
1928-43-4	2-Ethylhexyl-2,4-dichlorophenoxy acetate	N-10356-1G	2540-82-1	Formothion	N-12017-100MG
1928-43-4	2-Ethylhexyl-2,4-dichlorophenoxy acetate Solution	S-10356A1-1ML	2540-82-1	Formothion Solution	S-12017A1-1ML
1928-43-4	2-Ethylhexyl-2,4-dichlorophenoxy acetate Solution	S-10356J1-1ML	2540-82-1	Formothion Solution	S-12017J1-1ML
1928-43-4	2,4-D 2-ethylhexyl ester Solution	S-10531A4-1ML	2545-59-7	2,4,5-T butoxyethyl ester	N-10648-100MG
1928-43-4	2,4-D 2-ethylhexyl ester Solution	S-10531A4-5ML	2545-59-7	2,4,5-T butoxyethyl ester Solution	S-10648A4-1ML
1928-47-8	2,4,5-T 2-ethylhexyl ester	N-10647-250MG	2545-59-7	2,4,5-T butoxyethyl ester Solution	S-10648A4-5ML
1928-47-8	2,4,5-T 2-ethylhexyl ester Solution	S-10647A1-1ML	2545-59-7	2,4,5-T butoxyethyl ester Solution	S-10648U1-1ML
1928-47-8	2,4,5-T 2-ethylhexyl ester Solution	S-10647T1-1ML	2564-95-6	N,N-Dimethylsuccinamic acid	N-12630-250MG
1929-73-3	2,4-D butoxyethyl ester	N-10610-250MG	2564-95-6	N,N-Dimethylsuccinamic acid Solution	S-12630A1-1ML
1929-73-3	2,4-D butoxyethyl ester Solution	S-10610A1-1ML	2564-95-6	N,N-Dimethylsuccinamic acid Solution	S-12630T1-1ML
1929-73-3	2,4-D butoxyethyl ester Solution	S-10610A4-1ML	2581-34-2	3-Methyl-4-nitrophenol	MET-11955B-100MG
1929-73-3	2,4-D butoxyethyl ester Solution	S-10610A4-5ML	2588-03-6	Phorate sulfoxide	MET-13035C-100MG
1929-77-7	Vernolate	N-13744-250MG	2588-04-7	Phorate sulfone	MET-13035B-100MG
1929-77-7	Vernolate Solution	S-13744A1-1ML	2588-04-7	Phorate sulfone Solution	MET-13035BJ1-1ML
1929-77-7	Vernolate Solution	S-13744T1-1ML	2593-15-9	Etridiazole	N-13517-250MG
1929-77-7	Vernolate Solution	S-13744T1-5ML	2593-15-9	Etridiazole Solution	S-13517A1-1ML
1929-82-4	Nitrapyrin	N-12659-100MG	2593-15-9	Etridiazole Solution	S-13517T1-1ML
1929-82-4	Nitrapyrin Solution	S-12659A1-1ML	2595-54-2	Mecarbam	N-12359-50MG
1942-71-8	2-(4-tert-Butylphenoxy)cyclohexanol	MET-12727A-100MG	2595-54-2	Mecarbam Solution	S-12359A1-1ML
1967-16-4	Chlorbufam	N-12970-10MG	2595-54-2	Mecarbam Solution	S-12359U1-1ML
1982-47-4	Chloroxuron	N-11456-250MG	2597-03-7	Phenthoate	N-13006-100MG
1982-47-4	Chloroxuron Solution	S-11456A1-1ML	2597-03-7	Phenthoate Solution	S-13006A1-1ML
1982-49-6	Siduron	N-13211-100MG	2597-03-7	Phenthoate Solution	S-13006U1-1ML
1982-49-6	Siduron Solution	S-13211A4-1ML	2597-93-5	Emmi (Technical)	N-10224-100MG
1982-49-6	Siduron Solution	S-13211A4-5ML	2599-11-3	Methiocarb sulfoxide Solution	MET-12398BU1-1ML
2008-39-1	2,4-D dimethylamine salt	N-10612-250MG	2599-11-3	Simazine hydroxy	MET-13800B-25MG
2008-39-1	2,4-D dimethylamine salt Solution	S-10612M1-1ML	2600-69-3	Phorate-oxon Solution	MET-13035AU1-1ML
2008-41-5	Butylate	N-11383-250MG	2631-37-0	Promecarb	N-13101-100MG
2008-41-5	Butylate Solution	S-11383A1-1ML	2631-37-0	Promecarb Solution	S-13101M4-1ML
2008-41-5	Butylate Solution	S-11383T1-1ML	2631-40-5	Isoprocab	N-12264-250MG
2008-41-5	Butylate Solution	S-11383T1-5ML	2631-40-5	Isoprocab Solution	S-12264A1-1ML
2008-58-4	2,6-Dichlorobenzamide	MET-10688B-1G	2631-40-5	Isoprocab Solution	S-12264U1-1ML
2012-74-0	2-(p-Chlorophenyl)-3-methyl butyric acid	MET-13201A-100MG	2636-26-2	Cyanophos	N-11517-100MG
2032-59-9	Aminocarb	N-11065-250MG	2636-26-2	Cyanophos Solution	S-11517A1-1ML
2032-59-9	Aminocarb Solution	S-11065A4-1ML	2636-26-2	Cyanophos Solution	S-11517U1-1ML
2032-59-9	Aminocarb Solution	S-11065A4-5ML	2642-71-9	Azinphos-ethyl	N-12140-100MG
2032-59-9	Aminocarb Solution	S-11065U1-1ML	2642-71-9	Guthion Ethyl Solution	S-12140A1-1ML
2032-65-7	Methiocarb	N-12398-100MG	2655-14-3	XMC	S-13275-100MG
2032-65-7	Methiocarb Solution	S-12398M1-1ML	2655-14-3	XMC Solution	S-13275A1-1ML
2058-46-0	Oxytetracycline hydrochloride	N-12743-250MG	2655-15-4	2,3,5-Trimethylphenyl methyl carbamate	N-10603-50MG
2104-64-5	EPN	N-11858-250MG	2655-15-4	2,3,5-Trimethylphenyl methyl carbamate Solution	S-10603A1-1ML
2104-64-5	EPN Solution	S-11858A1-1ML	2675-77-6	Chloroneb	N-11449-250MG
2104-64-5	EPN Solution	S-11858J1-1ML	2675-77-6	Chloroneb Solution	S-11449A1-1ML
2104-64-5	EPN Solution	S-11858J1-5ML	2675-77-6	Chloroneb Solution	S-11449T1-1ML
2104-64-5	EPN Solution	S-11858J4-1ML	2675-77-6	Chloroneb Solution	S-11449T1-5ML
2104-64-5	EPN Solution	S-11858J4-5ML	2686-99-9	3,4,5-Trimethylphenyl methyl carbamate	N-10769-10MG
2104-96-3	Bromophos methyl	N-11343-50MG	2758-06-7	(2-Bromoethyl)trimethyl ammonium bromide	N-10001-1G
2104-96-3	Bromophos methyl Solution	S-11343A1-1ML	2758-42-1	2,4-DB dimethylamine salt	N-10623-250MG
2104-96-3	Bromophos methyl Solution	S-11343J1-1ML	2758-42-1	2,4-DB dimethylamine salt Solution	S-10623M1-1ML
2122-70-5	Ethyl-1-naphthalene acetate	N-11914-250MG	2797-51-5	Quinoclamine	N-13825-100MG
2122-70-5	Ethyl-1-naphthalene acetate Solution	S-11914A1-1ML	2814-20-2	2-Isopropyl-6-methyl-4-pyrimidinol	MET-11621B-250MG
2122-70-5	Ethyl-1-naphthalene acetate Solution	S-11914T1-1ML	2876-78-0	Methyl-1-naphthalene acetate	N-12464-500MG
2132-70-9	p,p'-Methoxychlor-olefin	MET-12403B-50MG	2876-78-0	Methyl-1-naphthalene acetate Solution	S-12464A1-1ML
2132-70-9	p,p'-Methoxychlor-olefin Solution	MET-12403B1-1ML	2876-78-0	Methyl-1-naphthalene acetate Solution	S-12464U1-1ML
2157-19-9	Endosulfan alcohol	MET-11848A-50MG	2905-67-1	3,5-Dichlorobenzoic acid methyl ester	N-10791-100MG
2163-68-0	Atrazine-2-hydroxy	MET-11106F-100MG	2905-67-1	3,5-Dichlorobenzoic acid methyl ester Solution	S-10791A1-1ML
2163-80-6	Monosodium acid methane arsonate sesquihydrate	N-12495-100MG	2905-67-1	3,5-Dichlorobenzoic acid methyl ester Solution	S-10791T5-1ML
2163-80-6	Monosodium acid methane arsonate sesquihydrate	S-12495F1-1ML	2905-67-1	3,5-Dichlorobenzoic acid methyl ester Solution	S-10791T5-5ML
2164-08-1	Lenacil	N-12314-100MG	2921-88-2	Chlorpyrifos	N-11459-250MG
2164-08-1	Lenacil Solution	S-12314M1-1ML	2921-88-2	Chlorpyrifos Solution	S-11459A1-1ML
2164-17-2	Fluometuron	N-11992-250MG	2921-88-2	Chlorpyrifos Solution	S-11459K1-1ML
2164-17-2	Fluometuron Solution	S-11992A4-1ML	2921-88-2	Chlorpyrifos Solution	S-11459K1-5ML
2179-25-1	Methiocarb sulfone	MET-12398A-50MG	2976-74-1	2,3-Dichlorophenoxyacetic acid	N-12953-10MG
2212-67-1	Molinate	N-12487-250MG	3004-51-1	Allyl trimethylammonium bromide	N-11059-250MG
2212-67-1	Molinate Solution	S-12487A1-1ML	3004-51-1	Allyl trimethylammonium bromide Solution	S-11059F1-1ML
2212-67-1	Molinate Solution	S-12487T1-1ML	3004-51-1	Allyl trimethylammonium bromide Solution	S-11059U1-1ML
2212-67-1	Molinate Solution	S-12487T1-5ML	3060-89-7	Metobromuron	N-12477-250MG
2227-17-0	Dienochlor	N-12824-250MG	3060-89-7	Metobromuron Solution	S-12477A1-1ML
2227-17-0	Dienochlor Solution	S-12824M1-1ML	3060-89-7	Metobromuron Solution	S-12477U1-1ML
2255-17-6	Fenitrothion-O-analog Solution	MET-11955AJ1-1ML	3064-70-8	Hexachlorodimethyl sulfone	N-12161-250MG
2274-67-1	Dimethylvinphos Solution	S-14052U1-1ML	3064-70-8	Hexachlorodimethyl sulfone Solution	S-12161A1-1ML
2275-23-2	Vamidothion	N-13741-100MG	3064-70-8	Hexachlorodimethyl sulfone Solution	S-12161U1-1ML
2275-23-2	Vamidothion Solution	S-13741A1-1ML	3070-15-3	Fensulfthion-sulfide	MET-11963B-100MG
2303-16-4	Di-allate	N-11587-100MG	3209-22-1	2,3-Dichloronitrobenzene	N-10584-250MG
2303-16-4	Di-allate Solution	S-11587M1-1ML	3209-22-1	2,3-Dichloronitrobenzene Solution	S-10584A1-1ML
2303-16-4	Di-allate Solution	S-11587M1-5ML	3209-22-1	2,3-Dichloronitrobenzene Solution	S-10584U1-1ML
2303-17-5	Tri-allate	N-13628-250MG	3244-90-4	Aspon (TM)	N-11103-10MG
2303-17-5	Tri-allate Solution	S-13628A1-1ML	3244-90-4	Aspon (TM) Solution	S-11103A1-1ML
2303-17-5	Tri-allate Solution	S-13628K1-1ML	3244-90-4	Aspon (TM) Solution	S-11103J4-1ML
2307-68-8	Pentanochlor	N-12844-100MG	3244-90-4	Aspon (TM) Solution	S-11103J4-5ML
2307-68-8	Pentanochlor Solution	S-12844A1-1ML	3336-39-8	Bromoxynil-methyl ether	MET-11345B-100MG
2307-68-8	Pentanochlor Solution	S-12844U1-1ML	3336-40-1	loxynil-methyl	MET-10794A-100MG
2310-17-0	Phosalone	N-13036-100MG	3337-71-1	Asulam	N-11104-250MG
2310-17-0	Phosalone Solution	S-13036M1-1ML	3337-71-1	Asulam Solution	S-11104M1-1ML
2312-35-8	Propargite	N-12727-100MG	3337-71-1	Asulam Solution	S-11104M1-5ML
2312-35-8	Propargite Solution	S-12727A1-1ML	3344-12-5	Isomalathion Solution	MET-12346AJ1-1ML
2312-35-8	Propargite Solution	S-12727U1-1ML	3347-22-6	Dithianon	N-11821-250MG
2385-85-5	Mirex	N-12486-100MG	3369-52-6	Endosulfan ether	MET-11848B-50MG
2385-85-5	Mirex Solution	S-12486M1-1ML	3380-34-5	Triclosan	N-13076-100MG
2416-30-0	1-(4-tert-Butylphenoxy)propan-2-ol	MET-11085A-100MG	3380-34-5	Triclosan Solution	S-13076M1-1ML
2425-06-1	Captafol	N-11399-250MG	3383-96-8	Temephos	N-10996-100MG
2425-10-7	MPMC	N-13274-10MG	3383-96-8	Temephos Solution	S-10996A1-1ML
2436-73-9	4-Chloro-o-tolxyloxyacetic acid methyl ester	N-10819-100MG	3383-96-8	Temephos Solution	S-10996U1-1ML
2436-73-9	4-Chloro-o-tolxyloxyacetic acid methyl ester Solution	S-10819M1-1ML	3397-62-4	2-Chloro-4,6-diamino-1,3,5-triazine	MET-13213C-100MG
2436-73-9	4-Chloro-o-tolxyloxyacetic acid methyl ester Solution	S-10819M1-5ML	3401-80-7	3,6-Dichloro-2-hydroxy benzoic acid	MET-11656B-50MG
2439-01-2	Quinomethionate	N-13172-100MG	3424-82-6	o,p'-DDD olefin	MET-12706A-10MG
2439-10-3	Dodine (Technical)	N-11840-250MG	3424-82-6	o,p'-DDE	N-12707-50MG
2439-99-8	Glyphosine	N-12135-100MG	3424-82-6	o,p'-DDE Solution	S-12707M1-1ML
2439-99-8	Glyphosine Solution	S-12135F1-1ML	3424-82-6	o,p'-DDE Solution	S-12707M1-5ML
2439-99-8	Glyphosine Solution	S-12135U1-1ML	3478-94-2	Piperalin	S-13058-100MG
2463-84-5	Dicaphon	N-11658-100MG	3478-94-2	Piperalin Solution	S-13058A1-1ML
2463-84-5	Dicaphon Solution	S-11658A1-1ML	3478-94-2	Piperalin Solution	S-13058U1-1ML
2463-84-5	Dicaphon Solution	S-11658U1-1ML	3481-20-7	2,3,5,6-Tetrachloroaniline	N-10530-250MG
2497-06-5	Disulfoton sulfone	MET-11819A-100MG	3481-20-7	2,3,5,6-Tetrachloroaniline Solution	S-10530M1-1ML
2497-07-6	Disulfoton sulfoxide	MET-11819B-100MG	3547-33-9	2-Hydroxyethyl-n-octyl sulphide	N-10378-100MG
2514-53-6	EGT	N-11845-100MG	3547-33-9	2-Hydroxyethyl-n-octyl sulphide Solution	S-10378M1-1ML
2524-03-0	O,O-Dimethyl phosphorochloridothioate	N-12704-250MG	3566-10-7	Amobam	N-11071-100MG
2524-03-0	O,O-Dimethyl phosphorochloridothioate Solution	S-12704A1-1ML	3566-10-7	Amobam Solution	S-11071M1-1ML
2524-03-0	O,O-Dimethyl phosphorochloridothioate Solution	S-12704U1-1ML	3567-62-2	1-(3,4-Dichlorophenyl)-3-methylurea	MET-11827C-250MG

CAS#	Description	Part #	CAS#	Description	Part #
3586-60-5	Asomate	N-13904-50MG	6420-47-9	2-sec-Butyl-4,6-dinitrophenol triethanolamine salt	S-10517F1-1ML
3599-58-4	2,4-D ethanolamine salt	N-10613-250MG	6420-47-9	2-sec-Butyl-4,6-dinitrophenol triethanolamine salt	S-10517U1-1ML
3599-58-4	2,4-D ethanolamine salt Solution	S-10613M1-1ML	6515-38-4	3,5,6-Trichloro-2-pyridinol	MET-11459A-250MG
3653-48-3	MCPA sodium salt	MET-12348-25050MG	6552-12-1	Fenthion-O-analog Solution	MET-11964CK1-1ML
3689-24-5	Tetraethyldithiopyrophosphate	N-13540-50MG	6552-13-2	Fenthion oxon sulfoxide Solution	MET-11964FA1-1ML
3689-24-5	Tetraethyldithiopyrophosphate Solution	S-13540A1-1ML	6552-21-2	Fensulfthion Oxon Solution	MET-11963DA1-1ML
3689-24-5	Tetraethyldithiopyrophosphate Solution	S-13540J1-1ML	6597-78-0	Dicamba methyl ester	N-11657-100MG
3691-35-8	Chlorophacinone	N-11451-250MG	6597-78-0	Dicamba methyl ester Solution	S-11657M1-1ML
3691-35-8	Chlorophacinone Solution	S-11451A1-1ML	6597-78-0	Dicamba methyl ester Solution	S-11657M1-5ML
3734-33-6	Bitrex (TM)	N-11324-100MG	6636-55-1	6-Chloro-2-picolinic acid methyl ester Solution	MET-12659B1-1ML
3734-33-6	Bitrex (TM) Solution	S-11324F1-1ML	6711-48-4	Iminobis (methylphosphonic acid)	MET-12135A-500MG
3734-33-6	Bitrex (TM) Solution	S-11324U1-1ML	6734-80-1	Metam sodium	N-12382-250MG
3734-48-3	Chlordane	MET-11425D-250MG	6923-22-4	Monocrotophos	N-12493-100MG
3735-33-9	Phosmet-oxon Solution	MET-12207AU1-1ML	6923-22-4	Monocrotophos Solution	S-12493A1-1ML
3739-38-6	3-Phenoxybenzoic acid	MET-12848C-100MG	6923-22-4	Monocrotophos Solution	S-12493U1-1ML
3740-92-9	Fenchlorim	N-13265-250MG	6923-22-4	Monocrotophos Solution	S-12493U1-5ML
3761-41-9	Fenthion sulfoxide Solution	MET-11964BX1-1ML	6936-40-9	2,3,5,6-Tetrachloroanisole Solution	S-13984J1-1ML
3761-41-9	Fenthion sulfoxide	MET-11964D-100MG	6947-94-0	2-Methyl-3-furic acid	MET-14013A-100MG
3761-42-0	Fenthion-sulfone	MET-11964A-100MG	6967-29-9	2-Chloro-2',6'-diethylacetanilide	MET-11043A-100MG
3766-60-7	Buturon	N-12944-10MG	6981-18-6	Ormetoprim	N-12728-250MG
3766-81-2	Baycarb (TM)	N-11126-250MG	6981-18-6	Ormetoprim Solution	S-12728A1-1ML
3766-81-2	Baycarb (TM) Solution	S-11126M1-1ML	6981-18-6	Ormetoprim Solution	S-12728U1-1ML
3810-74-0	Streptomycin sulfate	N-13228-250MG	6988-21-2	Dioxacarb	N-11789-250MG
3811-49-2	Dioxabenzofos Solution	S-13964A1-1ML	6988-21-2	Dioxacarb Solution	S-11789MA-1ML
3813-05-6	Benazolin	N-12901-100MG	7047-28-1	Chlorbenside sulfoxide	MET-11423B-100MG
3861-47-0	Ioxynil octanoate	N-12218-100MG	7047-28-1	Chlorbenside sulfoxide Solution	MET-11423BX1-1ML
3861-47-0	Ioxynil octanoate Solution	S-12218A1-1ML	7082-99-7	Chlorbenside sulfone	MET-11423A-25MG
3861-47-0	Ioxynil octanoate Solution	S-12218U1-1ML	7082-99-7	Chlorbenside sulfone Solution	MET-11423AJ1-1ML
3868-61-9	Endosulfan lactone	MET-11848C-50MG	7085-19-0	Mecoprop	N-12360-100MG
3878-19-1	Fuberidazole	N-12021-100MG	7085-19-0	Mecoprop Solution	S-12360A1-1ML
3878-19-1	Fuberidazole Solution	S-12021M1-1ML	7085-19-0	Mecoprop Solution	S-12360B1-1ML
3964-56-5	4-Bromo-2-chlorophenol	MET-13097A-500MG	7173-84-4	Carbophenothion methyl-O-analog Solution	MET-11408BH10-1ML
3983-45-7	Fenchlorphos-oxon	MET-11952A-100MG	7286-69-3	Sebuthylazin	N-12950-250MG
4147-51-7	Dipropetryn	N-11808-100MG	7286-84-2	Chloramben methyl ester	N-11419-100MG
4147-51-7	Dipropetryn Solution	S-11808A1-1ML	7286-84-2	Chloramben methyl ester Solution	S-11419Y1-1ML
4147-51-7	Dipropetryn Solution	S-11808T1-1ML	7286-84-2	Chloramben methyl ester Solution	S-11419Y1-5ML
4300-43-0	N,N'-Bis(3,4-dichlorophenyl)urea	MET-11827A-100MG	7287-19-6	Prometryne	N-13103-250MG
4482-55-7	Fenuron TCA	N-11968-50MG	7287-19-6	Prometryne Solution	S-13103A1-1ML
4726-14-1	Nitralin	N-12658-250MG	7287-19-6	Prometryne Solution	S-13103T1-1ML
4726-14-1	Nitralin Solution	S-12658A1-1ML	7287-19-6	Prometryne Solution	S-13103T1-5ML
4726-14-1	Nitralin Solution	S-12658J1-1ML	7287-36-7	Monalide	N-13917-10MG
4727-29-1	Phthalanilic acid	N-13043-250MG	7292-16-2	Propaphos Solution	S-13967A1-1ML
4727-29-1	Phthalanilic acid Solution	S-13043A1-1ML	7332-32-3	Carbophenothion methyl-O-analog	MET-11408B-50MG
4727-29-1	Phthalanilic acid Solution	S-13043T1-1ML	7359-55-9	Thionazon-O-analog	MET-13760A-50MG
4812-20-8	2-Isopropoxyphenol	MET-11128A-250MG	7374-53-0	Propazine-2-hydroxy Solution	MET-13108AAA1-1ML
4824-78-6	Bromophos ethyl	N-11342-100MG	7421-93-4	Endrin aldehyde	N-11855-10MG
4824-78-6	Bromophos ethyl Solution	S-11342J1-1ML	7421-93-4	Endrin aldehyde Solution	S-11855M1-1ML
4841-20-7	Silvex methyl ester	N-13214-250MG	7421-93-4	Endrin aldehyde Solution	S-11855M1-5ML
4841-20-7	Silvex methyl ester Solution	S-13214M1-1ML	7446-18-6	Thallium sulfate(Technical)	N-13559-250MG
4841-20-7	Silvex methyl ester Solution	S-13214M1-5ML	7696-12-0	Tetramethrin	N-13552-100MG
4849-32-5	Karbutilate	N-12945-250MG	7696-12-0	Tetramethrin Solution	S-13552A4-1ML
5103-71-9	cis-Chlordane	N-11480-10MG	7696-12-0	Tetramethrin Solution	S-13552A4-5ML
5103-71-9	cis-Chlordane Solution	S-11480M1-1ML	7696-12-0	Tetramethrin Solution	S-13552U1-1ML
5103-71-9	cis-Chlordane Solution	S-11480M1-5ML	7700-17-6	Crotoxypfos	N-11513-50MG
5103-74-2	trans-Chlordane	N-13615-10MG	7700-17-6	Crotoxypfos Solution	S-11513K1-1ML
5103-74-2	trans-Chlordane Solution	S-13615M1-1ML	7700-17-6	Crotoxypfos Solution	S-11513K1-5ML
5103-74-2	trans-Chlordane Solution	S-13615M1-5ML	7773-06-0	Ammonium sulfamate(Technical)	N-11070-500MG
5131-24-8	Ditalimfos	N-11820-100MG	7778-43-0	Sodium arsenate	N-13215-250MG
5131-24-8	Ditalimfos Solution	S-11820A1-1ML	7778-43-0	Sodium arsenate Solution	S-13215F1-1ML
5131-24-8	Ditalimfos Solution	S-11820J1-1ML	7778-43-0	Sodium arsenate Solution	S-13215U1-1ML
5131-24-8	Ditalimfos Solution	S-11820J1-5ML	7779-90-0	Zinc phosphate(Technical)	N-13758-250MG
5131-24-8	Ditalimfos Solution	S-11820U1-1ML	7784-40-9	Lead arsenate(Technical)	N-12310-100MG
5234-68-4	Carboxin	N-11410-250MG	7786-34-7	Mevinphos	N-13037-250MG
5234-68-4	Carboxin Solution	S-11410A1-1ML	7786-34-7	Mevinphos Solution	S-13037A1-1ML
5259-88-1	Oxycarboxin	N-12740-100MG	7786-34-7	Mevinphos Solution	S-13037U1-1ML
5259-88-1	Oxycarboxin Solution	S-12740A1-1ML	7786-34-7	Mevinphos Solution	S-13037U1-5ML
5259-88-1	Oxycarboxin Solution	S-12740U1-1ML	8001-35-2	Camphechlor	N-13586-250MG
5259-88-1	Oxycarboxin Solution	S-12740U1-1ML	8001-35-2	Camphechlor Solution	S-13586B6-1ML
5326-23-8	6-Chloronicotinic acid	MET-12066B-100MG	8001-35-2	Camphechlor Solution	S-13586J4-1ML
5326-23-8	6-Chloronicotinic acid Solution	MET-12066BM1-1ML	8001-35-2	Camphechlor Solution	S-13586J4-5ML
5333-99-3	N-(1-Naphthyl)phthalimide	MET-12507A-250MG	8001-35-2	Camphechlor Solution	S-13586M1-1ML
5352-88-5	3-(4-Chlorophenyl)methyl urea	MET-11827B-100MG	8001-35-2	Camphechlor Solution	S-13586M1-5ML
5459-93-8	N-Ethylcyclohexylamine	MET-11522A-250MG	8001-50-1	Strobane (TM)	N-13229-100MG
5598-13-0	Chlorpyrifos Methyl	N-11460-250MG	8001-50-1	Strobane (TM) Solution	S-13229A1-1ML
5598-13-0	Chlorpyrifos Methyl Solution	S-11460A1-1ML	8001-50-1	Strobane (TM) Solution	S-13229K4-1ML
5598-13-0	Chlorpyrifos Methyl Solution	S-11460J4-1ML	8003-34-7	Pyrethrum	N-13151-100MG
5598-15-2	Chlorpyrifos Oxon	MET-11459B-100MG	8003-34-7	Pyrethrum Solution	S-13151A1-1ML
5598-52-7	Chlorpyrifos methyl-O-analog	MET-11460B-50MG	8003-34-7	Pyrethrum Solution	S-13151T1-1ML
5734-69-0	2-Ethylamino-6-methyl-4(3H)-Pyrimidinone	MET-13064C-100MG	8017-34-3	DDT (mixture p.p' & o,p')	N-11567-250MG
5742-17-6	2,4-D isopropylamine salt	N-10618-250MG	8017-34-3	DDT (mixture p.p' & o,p') Solution	S-11567M1-1ML
5742-17-6	2,4-D isopropylamine salt Solution	S-10618M1-1ML	8018-01-7	Mancozeb (Technical)	N-12353-250MG
5825-87-6	2-[3-Chlorophenoxy]propionamide	N-10261-1G	8065-48-3	Demeton O&S	N-11581-100MG
5825-87-6	2-[3-Chlorophenoxy]propionamide Solution	S-10261M1-1ML	8065-48-3	Demeton O&S Solution	S-11581U1-1ML
5836-10-2	Chloropropylate	N-11453-50MG	9002-91-9	Metaldehyde(Technical)	N-12381-1G
5836-29-3	Coumatetralyl	N-11509-250MG	9003-13-8	Butoxypolypropylene glycol(Technical)	N-11354-250MG
5836-29-3	Coumatetralyl Solution	S-11509M1-1ML	9006-42-2	Metiram(Technical)	N-13086-250MG
5902-51-2	Terbacil	N-13509-250MG	9016-72-2	Propineb(Technical)	N-13111-250MG
5902-51-2	Terbacil Solution	S-13509A1-1ML	10004-44-1	Hymexazole	N-12196-25MG
5902-51-2	Terbacil Solution	S-13509T1-1ML	10004-44-1	Hymexazole Solution	S-12196M1-1ML
5902-51-2	Terbacil Solution	S-13509T1-5ML	10265-92-6	Methamidophos	N-12393-100MG
5902-85-2	Vanicide-20S (TM)	N-13742-250MG	10265-92-6	Methamidophos Solution	S-12393A1-1ML
5902-85-2	Vanicide-20S (TM) Solution	S-13742A1-1ML	10265-92-6	Methamidophos Solution	S-12393U1-1ML
5902-85-2	Vanicide-20S (TM) Solution	S-13742U1-1ML	10380-28-6	Copper-8-Quinolinate	N-11503-500MG
5903-13-9	N-Methyl-N-1-naphthyl acetamide	N-12559-100MG	10453-86-8	Resmethrin	N-13176-250MG
5903-13-9	N-Methyl-N-1-naphthyl acetamide Solution	S-12559M1-1ML	10453-86-8	Resmethrin Solution	S-13176A1-1ML
5915-41-3	Terbutylazine	N-13512-50MG	10453-86-8	Resmethrin Solution	S-13176T4-1ML
5915-41-3	Terbutylazine Solution	S-13512A1-1ML	10453-86-8	Resmethrin Solution	S-13176T4-5ML
5915-41-3	Terbutylazine Solution	S-13512T4-1ML	10548-10-4	Terbufos sulfoxide	MET-13510B-100MG
5945-86-8	Nimbin	N-12912-5MG	10552-74-6	Nitrothal isopropyl	N-12665-100MG
6099-79-2	Dinoseb methyl ether	N-11787-100MG	10605-21-7	Carbendazim	N-11404-100MG
6099-79-2	Dinoseb methyl ether Solution	S-11787M1-1ML	10605-21-7	Carbendazim Solution	S-11404M1-1ML
6099-79-2	Dinoseb methyl ether Solution	S-11787M1-5ML	10605-21-7	Carbendazim Solution	S-11404M1-5ML
6132-17-8	Fensulfthion oxon sulfone Solution	MET-11963CA1-1ML	11141-17-6	Azadirachtin(Technical)	N-11107-10MG
6164-98-3	Chlordimeform	N-11426-250MG	12122-67-7	Zineb(Technical)	N-13759-250MG
6164-98-3	Chlordimeform Solution	S-11426A1-1ML	12771-68-5	Ancymidol	N-11073-100MG
6164-98-3	Chlordimeform Solution	S-11426U1-1ML	12771-68-5	Ancymidol Solution	S-11073M1-1ML
6190-65-4	Atrazine desethyl	MET-11106B-250MG	13067-93-1	Cyanofenphos	N-12881-50MG
6270-63-9	2-Hydroxypyrazine Solution	MET-13760BM1-1ML	13067-93-1	Cyanofenphos Solution	S-12881M1-1ML
6365-72-6	2,4-D isopropanol amine salt	N-10616-250MG	13071-79-9	Terbufos	N-13510-250MG
6365-72-6	2,4-D isopropanol amine salt Solution	S-10616M1-1ML	13071-79-9	Terbufos Solution	S-13510M1-1ML
6385-62-2	Diquat Dibromide monohydrate	N-11816-500MG	13121-70-5	Tricyclohexyltin hydroxide	N-13662-250MG
6385-62-2	Diquat Dibromide monohydrate Solution	S-11816F1-1ML	13121-70-5	Tricyclohexyltin hydroxide Solution	S-13662M1-1ML
6385-62-2	Diquat Dibromide monohydrate Solution	S-11816F13-5ML	13171-21-6	Phosphamidon	N-13039-250MG
6420-47-9	2-sec-Butyl-4,6-dinitrophenol triethanolamine salt	N-10517-500MG	13171-21-6	Phosphamidon Solution	S-13039A1-1ML

CAS#	Description	Part #	CAS#	Description	Part #
13171-21-6	Phosphamidon Solution	S-13039U1-1ML	21087-64-9	Metribuzin Solution	S-12481A1-1ML
13171-21-6	Phosphamidon Solution	S-13039U1-5ML	21087-64-9	Metribuzin Solution	S-12481T1-1ML
13181-17-4	Bromofenoxim	N-12961-10MG	21087-64-9	Metribuzin Solution	S-12481T1-5ML
13194-48-4	Prophos	N-13110-250MG	21564-17-0	TCMTB	N-13500-10MG
13194-48-4	Prophos Solution	S-13110A1-1ML	21564-17-0	TCMTB Solution	S-13500A1-1ML
13194-48-4	Prophos Solution	S-13110U1-1ML	21564-17-0	TCMTB Solution	S-13500X4-1ML
13194-48-4	Prophos Solution	S-13110U1-5ML	21564-17-0	TCMTB Solution	S-13500X4-5ML
13356-08-6	Fenbutatin-oxide(Technical)	N-11951-250MG	21609-90-5	Leptophos	N-12315-100MG
13360-45-7	Chlorbromuron	N-11424-250MG	21609-90-5	Leptophos Solution	S-12315A1-1ML
13360-45-7	Chlorbromuron Solution	S-11424A1-1ML	21609-90-5	Leptophos Solution	S-12315U1-1ML
13360-45-7	Chlorbromuron Solution	S-11424U1-1ML	21725-46-2	Bladex	N-11327-250MG
13457-18-6	Pyrazophos	N-13146-100MG	21725-46-2	Bladex Solution	S-11327A1-1ML
13457-18-6	Pyrazophos Solution	S-13146A1-1ML	21725-46-2	Bladex Solution	S-11327T4-1ML
13457-18-6	Pyrazophos Solution	S-13146U1-1ML	21725-46-2	Bladex Solution	S-11327T4-5ML
13532-25-7	Trietazine-2-hydroxy	MET-13666B-10MG	22175-22-0	1-(3-Chloro-4-methylphenyl)-3-methylurea	MET-11455A-100MG
13532-25-7	Trietazine-2-hydroxy Solution	MET-13666BA1-1ML	22212-55-1	Benzoylprop ethyl	N-11177-100MG
13593-03-8	Quinalphos	N-13168-250MG	22212-55-1	Benzoylprop ethyl Solution	S-11177A4-1ML
13593-03-8	Quinalphos Solution	S-13168A1-1ML	22212-55-1	Benzoylprop ethyl Solution	S-11177A4-5ML
13593-03-8	Quinalphos Solution	S-13168U1-1ML	22212-55-1	Benzoylprop ethyl Solution	S-11177U1-1ML
13684-56-5	Desmedipham	N-11584-250MG	22224-92-6	Fenamiphos	N-12854-250MG
13684-56-5	Desmedipham Solution	S-11584M1-1ML	22224-92-6	Fenamiphos Solution	S-12854A1-1ML
13684-63-4	Phenmedipham	N-12859-250MG	22224-92-6	Fenamiphos Solution	S-12854T1-1ML
13684-63-4	Phenmedipham Solution	S-12859A1-1ML	22224-92-6	Fenamiphos Solution	S-12854T1-5ML
13684-63-4	Phenmedipham Solution	S-12859U1-1ML	22248-79-9	Tetrachlorvinphos	N-13537-250MG
13826-35-2	m-Phenoxybenzyl alcohol	MET-12848B-100MG	22248-79-9	Tetrachlorvinphos Solution	S-13537A1-1ML
14086-35-2	Fenthion-oxon-sulfone	MET-11964E-100MG	22248-79-9	Tetrachlorvinphos Solution	S-13537X1-1ML
14143-25-0	Aminocarb-phenol	MET-11065-100MG	22248-79-9	Tetrachlorvinphos Solution	S-13537X1-5ML
14143-55-6	Picloram methyl ester	N-13051-100MG	22781-23-3	Bendiocarb	N-11135-250MG
14143-55-6	Picloram methyl ester Solution	S-13051M1-1ML	22781-23-3	Bendiocarb Solution	N-11135U1-1ML
14214-32-5	Difenoxuron	N-12977-10MG	22936-75-0	Dimepax	N-12984-10MG
14255-72-2	Fensulfthion-sulfone	MET-11963A-50MG	22936-86-3	Cyprazine	N-11546-50MG
14437-17-3	Chlorfenprop-methyl	N-12898-50MG	22936-86-3	Cyprazine Solution	S-11546A1-1ML
14484-64-1	Ferbam	N-11970-250MG	22936-86-3	Cyprazine Solution	S-11546B1-1ML
14484-64-1	Ferbam Solution	S-11970M1-1ML	23031-36-9	Prallethrin	N-13088-10MG
14484-64-1	Ferbam Solution	S-11970U1-1ML	23031-36-9	Prallethrin Solution	S-13088A1-1ML
14698-29-4	Oxolinic acid	N-12738-250MG	23031-36-9	Prallethrin Solution	S-13088K1-1ML
14698-29-4	Oxolinic acid Solution	S-12738A1-1ML	23103-98-2	Pirimicarb	N-13062-250MG
14763-60-1	4-Methylsulfonylphenol	MET-11964G-50MG	23103-98-2	Pirimicarb Solution	S-13062A1-1ML
14816-18-3	Phoxim	N-13041-100MG	23103-98-2	Pirimicarb Solution	S-13062U1-1ML
14816-18-3	Phoxim Solution	S-13041A1-1ML	23135-22-0	Oxamyl	N-12735-100MG
14816-18-3	Phoxim Solution	S-13041U1-1ML	23135-22-0	Oxamyl Solution	S-12735A1-1ML
14920-87-7	2,6-Dichlorobenzoic acid methyl ester	MET-10688C-250MG	23184-66-9	Butachlor	N-11353-100MG
14958-90-7	Simeconazole	N-13931-25MG	23184-66-9	Butachlor Solution	S-11353A1-1ML
15020-57-2	p-Hydroxy-N,N-dimethylbenzene sulfonamide	MET-11944B-10MG	23184-66-9	Butachlor Solution	S-11353T1-1ML
15165-67-0	Dichlorprop-P	N-11674-250MG	23184-66-9	Butachlor Solution	S-11353T1-5ML
15165-67-0	Dichlorprop-P Solution	S-11674A1-1ML	23422-53-9	Formetanate hydrochloride	N-12015-250MG
15165-67-0	Dichlorprop-P Solution	S-11674U1-1ML	23422-53-9	Formetanate hydrochloride Solution	S-12015A1-1ML
15263-52-2	Cartap hydrochloride	N-11412-250MG	23505-41-1	Pirimiphos-ethyl	N-13063-250MG
15299-99-7	Napropamide	N-11585-250MG	23505-41-1	Pirimiphos-ethyl Solution	S-13063U1-1ML
15299-99-7	Napropamide Solution	S-11585A1-1ML	23560-59-0	Heptenophos	N-12153-100MG
15299-99-7	Napropamide Solution	S-11585T1-1ML	23560-59-0	Heptenophos Solution	S-12153A1-1ML
15299-99-7	Napropamide Solution	S-11585U1-5ML	23560-59-0	Heptenophos Solution	S-12153U1-1ML
15310-01-7	Benodanil	N-12991-10MG	23564-05-8	Thiophanate-methyl	N-13572-250MG
15457-05-3	Fluorodifen	N-11998-250MG	23564-05-8	Thiophanate-methyl Solution	S-13572A1-1ML
15457-05-3	Fluorodifen Solution	S-11998M1-1ML	23564-06-9	Thiophanate	N-13571-100MG
15454-48-9	Chlorotoluron	N-11455-250MG	23564-06-9	Thiophanate Solution	S-13571U1-1ML
15454-48-9	Chlorotoluron Solution	S-11455A1-1ML	23844-56-6	Mecoprop methyl ester	N-12361-100MG
15454-48-9	Chlorotoluron Solution	S-11455U1-1ML	23844-56-6	Mecoprop methyl ester Solution	S-12361M1-1ML
15879-93-3	a-Chloralose(Technical)	N-10978-250MG	23844-56-6	Mecoprop methyl ester Solution	S-12361M1-5ML
15972-60-8	Alachlor	N-11043-250MG	23844-56-6	Mecoprop-2,4,4-trimethylpentyl ester Solution	S-12362B1-1ML
15972-60-8	Alachlor Solution	S-11043M1-1ML	23844-56-6	Mecoprop-2,4,4-trimethylpentyl ester Solution	S-12362B1-5ML
15972-60-8	Alachlor Solution	S-11043M1-5ML	23906-13-0	2-Hydrazino-4,6-dimethylpyrimidine	MET-14030A-100MG
16079-88-2	Agribrom	N-11041-500MG	23906-13-0	2-Hydrazino-4,6-dimethylpyrimidine Solution	MET-14030AM1-1ML
16118-49-3	Carbetamide	N-12971-10MG	23947-60-6	Ethirimol	N-12986-10MG
16376-36-6	b-Chloralose(Technical)	N-11117-100MG	23950-58-5	Propyzamide	N-13136-250MG
16484-77-8	Mecoprop-P	N-12363-250MG	23950-58-5	Propyzamide Solution	S-13136M1-1ML
16484-77-8	Mecoprop-P Solution	S-12363A1-1ML	23950-58-5	Propyzamide Solution	S-13136M1-5ML
16484-77-8	Mecoprop-P Solution	S-12363T1-1ML	24017-47-8	Triazophos	N-13639-100MG
16662-85-4	Carbophenothion sulfone	MET-11408D-25MG	24017-47-8	Triazophos Solution	S-13639A1-1ML
16672-87-0	(2-Chloroethyl)phosphonic acid	N-10002-250MG	24017-47-8	Triazophos Solution	S-13639U1-1ML
16709-30-1	Carbofuran-3-keto	MET-11405E-100MG	24017-47-8	Triazophos Solution	S-13639U1-5ML
16752-77-5	Methomyl	N-12399-100MG	24096-53-5	Dimethachlon	N-11753-500MG
16752-77-5	Methomyl Solution	S-12399A1-1ML	24096-53-5	Dimethachlon Solution	S-11753M1-1ML
16919-19-0	Ammonium fluosilicate(Technical)	N-11069-500MG	24151-93-7	Piperophos	N-12989-250MG
17040-19-6	Dioxydemeton-S-methyl	N-11791-100MG	24157-81-1	2,6-Diisopropyl-naphthalene	N-10550-250MG
17040-19-6	Dioxydemeton-S-methyl Solution	S-11791A1-1ML	24157-81-1	2,6-Diisopropyl-naphthalene Solution	S-10550A1-1ML
17109-49-8	Edifenphos	N-11844-100MG	24157-81-1	2,6-Diisopropyl-naphthalene Solution	S-10550U1-1ML
17109-49-8	Edifenphos Solution	S-11844A1-1ML	24307-26-4	Mepiquat chloride	N-12371-100MG
17109-49-8	Edifenphos Solution	S-11844U1-1ML	24307-26-4	Mepiquat chloride Solution	S-12371F1-1ML
17297-40-4	Carbophenothion sulfoxide	MET-11408C-100MG	24353-61-5	Isocarboxophos Solution	S-13933A1-1ML
17606-31-4	Bensultap	N-11141-250MG	24579-73-5	Propamocarb	N-13105-250MG
17606-31-4	Bensultap Solution	S-11141M1-1ML	24691-80-3	Fenfuram	N-14013-250MG
17640-02-7	Dalapon methyl ester	N-11563-100MG	24691-80-3	Fenfuram Solution	S-14013M1-1ML
17640-02-7	Dalapon methyl ester Solution	S-11563M1-1ML	24934-91-6	Chlormephos	N-11434-250MG
17781-15-6	2,3-Dihydro-2,2-dimethyl benzofuran-3,7-diol Solution	MET-11405DA1-1ML	24934-91-6	Chlormephos Solution	S-11434A1-1ML
17781-16-7	Carbofuranphenol-3-keto	MET-11405F-50MG	24934-91-6	Chlormephos Solution	S-11434U1-1ML
17804-35-2	Benomyl	N-11138-100MG	25006-32-0	Leptophos oxon	MET-12315A-50MG
18181-80-1	Bromopropylate	N-10997-250MG	25057-89-0	Bentazon	N-11142-250MG
18181-80-1	Bromopropylate Solution	S-10997A1-1ML	25057-89-0	Bentazon Solution	S-11142A4-1ML
18467-77-1	Dikegulac acid	N-11750-250MG	25057-89-0	Bentazon Solution	S-11142A4-5ML
18467-77-1	Dikegulac acid Solution	S-11750T1-1ML	25057-89-0	Bentazon Solution	S-11142B1-1ML
18530-56-8	Noruron	N-12669-250MG	25059-80-7	Benazolin-ethyl	N-11134-250MG
18530-56-8	Noruron Solution	S-12669A1-1ML	25059-80-7	Benazolin-ethyl Solution	S-11134A1-1ML
18530-56-8	Noruron Solution	S-12669U1-1ML	25059-80-7	Benazolin-ethyl Solution	S-11134T1-1ML
18625-12-2	2,4-DB methyl ester	N-10533-100MG	25168-15-4	2,4,5-T isooctyl ester	N-10649-100MG
18625-12-2	2,4-DB methyl ester Solution	S-10533M1-1ML	25168-15-4	2,4,5-T isooctyl ester Solution	S-10649A1-1ML
18625-12-2	2,4-DB methyl ester Solution	S-10533M1-5ML	25168-15-4	2,4,5-T isooctyl ester Solution	S-10649T1-1ML
18691-97-9	Methabenzthiazuron	N-12387-100MG	25168-26-7	2,4-D isooctyl ester	N-10615-250MG
18691-97-9	Methabenzthiazuron Solution	S-12387A1-1ML	25168-26-7	2,4-D isooctyl ester Solution	S-10615A1-1ML
18691-97-9	Methabenzthiazuron Solution	S-12387U1-1ML	25205-08-7	Isafenphos-des-N-isopropyl Solution	S-10615T1-1ML
18854-01-8	Isoxathion	N-13932-25MG	25205-08-7	Isafenphos-des-N-isopropyl Solution	MET-12249AJ1-1ML
19044-88-3	Oryzalin	N-12729-100MG	25311-71-1	Isafenphos	MET-12249C-100MG
19044-88-3	Oryzalin Solution	S-12729A4-1ML	25311-71-1	Isafenphos Solution	N-12249-100MG
19666-30-9	Oxadiazon	N-12731-250MG	25311-71-1	Isafenphos Solution	S-12249A1-1ML
19666-30-9	Oxadiazon Solution	S-12731A1-1ML	25311-71-1	Isafenphos Solution	S-12249U1-1ML
19666-30-9	Oxadiazon Solution	S-12731U1-1ML	25606-41-1	Propamocarb HCl	N-13106-250MG
19932-84-4	6-Chloro-2,3-dihydrobenzoxazol-2-one	MET-11957A-100MG	25606-41-1	Propamocarb HCl Solution	S-13106F1-1ML
19937-59-8	Metoxuron	N-12480-250MG	25713-60-4	2,4,6-Tris-(2,4,6-tribromophenoxy)-1,3,5-triazine	N-12888-50MG
19937-59-8	Metoxuron Solution	S-12480A1-1ML	25713-60-4	2,4,6-Tris-(2,4,6-tribromophenoxy)-1,3,5-triazine Solution	S-12888U1-1ML
19937-59-8	Metoxuron Solution	S-12480U1-1ML	25954-13-6	Fosamine Ammonium Solution	S-12018F1-1ML
19988-24-0	Atrazine desethyl-2-hydroxy	MET-11106D-50MG	26002-80-2	d-(cis-trans)-Phenothrin	N-11554-50MG
19988-24-0	Atrazine-desethyl-2-hydroxy Solution	MET-11106DAM1-1ML	26002-80-2	d-(cis-trans)-Phenothrin Solution	S-11554A4-1ML
20416-09-5	5-(Phenylmethyl)-3-furanmethanol	MET-13176A-100MG	26002-80-2	d-(cis-trans)-Phenothrin Solution	S-11554U1-1ML
21087-64-9	Metribuzin	N-12481-250MG	26087-47-8	Iprobenfos	N-12219-250MG

CAS#	Description
26087-47-8	Iprobenfos Solution
26087-47-8	Iprobenfos Solution
26225-79-6	Ethofumesate
26225-79-6	Ethofumesate Solution
26225-79-6	Ethofumesate Solution
26399-36-0	Profuralin
26399-36-0	Profuralin Solution
26399-36-0	Profuralin Solution
26399-36-0	Profuralin Solution
26530-20-1	Octhilinone
26530-20-1	Octhilinone Solution
26644-46-2	Triforine
26644-46-2	Triforine Solution
26718-65-0	E-Mevinphos
26718-65-0	E-Mevinphos Solution
26718-65-0	E-Mevinphos Solution
27304-13-8	Oxychlorane Solution
27314-13-2	Norflurazon
27314-13-2	Norflurazon Solution
27314-13-2	Norflurazon Solution
27314-13-2	Norflurazon Solution
27355-22-2	Fthalide
27605-76-1	Probenazole
27605-76-1	Probenazole Solution
28044-83-9	trans-Heptachlor epoxide
28044-83-9	trans-Heptachlor epoxide Solution
28057-48-9	Bioallethrin
28057-48-9	Bioallethrin Solution
28057-48-9	Bioallethrin Solution
28086-13-7	Phenylmercury salicylate(Technical)
28086-13-7	Phenylmercury salicylate Solution
28086-13-7	Phenylmercury salicylate Solution
28159-98-0	Cybutryne
28159-98-0	Cybutryne Solution
28159-98-0	Cybutryne Solution
28167-49-9	N-Methyl-methamidophos
28167-49-9	N-Methyl-methamidophos Solution
28167-49-9	N-Methyl-methamidophos Solution
28249-77-6	Benthiocarb
28249-77-6	Benthiocarb Solution
28434-00-6	s-Bioallethrin
28434-00-6	s-Bioallethrin Solution
28434-00-6	s-Bioallethrin Solution
28434-01-7	Bioresmethrin
28434-01-7	Bioresmethrin Solution
28434-01-7	Bioresmethrin Solution
28772-56-7	Super-Rozol
28772-56-7	Super-Rozol Solution
28772-56-7	Super-Rozol Solution
29091-05-2	Dinitramine
29091-21-2	Prodiamine
29091-21-2	Prodiamine Solution
29091-21-2	Prodiamine Solution
29232-93-7	Pirimiphos-methyl
29232-93-7	Pirimiphos-methyl Solution
29232-93-7	Pirimiphos-methyl Solution
29973-13-5	Ethiofencarb
29973-13-5	Ethiofencarb Solution
29973-13-5	Ethiofencarb Solution
30043-49-3	Ethidimuron
30043-49-3	Ethidimuron Solution
30043-49-3	Ethidimuron Solution
30125-63-4	Terbutylazine-desethyl
30125-64-5	Terbutmeton-desethyl
30391-89-0	2-Amino-N-isopropyl benzamide
30558-43-1	Oxamyl oxime
30560-19-1	Acephate
30560-19-1	Acephate Solution
30560-19-1	Acephate Solution
30560-19-1	Acephate Solution
30614-22-3	Pirimicarb-desmethyl
30667-99-3	o,p'-Methoxychlor
30864-28-9	Methacrifos
30864-28-9	Methacrifos Solution
30864-28-9	Methacrifos Solution
31120-83-9	Isofenphos-des-N-isopropyl-O-analogue Solution
31218-83-4	Propetamphos
31218-83-4	Propetamphos Solution
31218-83-4	Propetamphos Solution
31717-87-0	Dodemorph acetate
31972-43-7	Fenamiphos sulfoxide
31972-44-8	Fenamiphos sulfone
32428-71-0	2-Chloro-N-(2-ethyl-6-methylphenyl) acetamide
32809-16-8	Procymidone
32809-16-8	Procymidone Solution
32809-16-8	Procymidone Solution
32889-48-8	Procyazine
32889-48-8	Procyazine Solution
32889-48-8	Procyazine Solution
33089-61-1	Amirtraz
33089-61-1	Amirtraz Solution
33089-74-6	Amirtraz Metabolite BTS 27271
33213-65-9	b-Endosulfan
33213-65-9	b-Endosulfan Solution
33213-65-9	b-Endosulfan Solution
33213-65-9	b-Endosulfan Solution
33245-39-5	Fluchloralin
33245-39-5	Fluchloralin Solution
33245-39-5	Fluchloralin Solution
33245-39-5	Fluchloralin Solution
33386-84-4	Oxy-chlordene Solution
33629-47-9	Butralin
33629-47-9	Butralin Solution
33629-47-9	Butralin Solution
33693-04-8	Terbutmeton
33693-04-8	Terbutmeton Solution
33820-53-0	Isopropalin
33820-53-0	Isopropalin Solution
33820-53-0	Isopropalin Solution
33820-53-0	Isopropalin Solution
34014-18-1	Tebuthiuron
34014-18-1	Tebuthiuron Solution
34014-18-1	Tebuthiuron Solution
34014-18-1	Tebuthiuron Solution

Part #
S-12219A1-1ML
S-12219U1-1ML
N-11874-250MG
S-11874A1-1ML
S-11874U1-1ML
N-13098-250MG
S-13098A1-1ML
S-13098J4-1ML
S-13098J4-5ML
N-14012-250MG
S-14012X1-1ML
N-13691-250MG
S-13691U1-1ML
N-11481-100MG
S-11481A1-1ML
S-11481T1-1ML
MET-11425AM1-1ML
N-12668-100MG
S-12668A1-1ML
S-12668T1-1ML
S-12668T1-5ML
N-13085-10MG
N-13092-500MG
S-13092M1-1ML
N-13618-50MG
S-13618M1-1ML
N-11205-250MG
S-11205A1-1ML
S-11205T1-1ML
N-13033-250MG
S-13033A1-1ML
S-13033U1-1ML
N-11520-100MG
S-11520A1-1ML
S-11520U1-1ML
N-12558-50MG
S-12558A1-1ML
S-12558U1-1ML
N-11144-250MG
S-11144M1-1ML
N-13187-250MG
S-13187A1-1ML
S-13187T1-1ML
N-11206-250MG
S-11206A1-1ML
S-11206T1-1ML
N-13259-100MG
S-13259A1-1ML
S-13259U1-1ML
N-12976-10MG
N-13096-100MG
S-13096A1-1ML
S-13096K1-1ML
N-13064-250MG
S-13064A1-1ML
S-13064U1-1ML
N-11872-100MG
S-11872A1-1ML
S-11872U1-1ML
N-11871-100MG
S-11871A1-1ML
S-11871U1-1ML
MET-13512A-50MG
MET-13511-50MG
MET-11142C-50MG
MET-12735-10MG
N-11002-250MG
S-11002A1-1ML
S-11002B4-1ML
S-11002B4-5ML
MET-13801A-50MG
MET-12403A-100MG
N-12388-100MG
S-12388A1-1ML
S-12388U1-1ML
MET-12249BJ1-1ML
N-13197-100MG
S-13197A1-1ML
S-13197U1-1ML
N-12948-100MG
MET-11947A-50MG
MET-11947B-100MG
MET-11013-50MG
N-13095-250MG
S-13095A1-1ML
S-13095U1-1ML
N-13094-100MG
S-13094A1-1ML
S-13094T1-1ML
N-11068-250MG
S-11068A1-1ML
MET-11068B-10MG
N-11118-100MG
S-11118A1-1ML
S-11118J1-1ML
S-11118J1-5ML
N-11983-50MG
S-11983A1-1ML
S-11983J4-1ML
S-11983J4-5ML
S-14054J1-1ML
N-11355-500MG
S-11355A1-1ML
S-11355K1-1ML
N-13511-250MG
S-13511M1-1ML
N-12265-250MG
S-12265A1-1ML
S-12265J4-1ML
S-12265J4-5ML
N-13505-250MG
S-13505A1-1ML
S-13505T1-1ML
S-13505T1-5ML

CAS#	Description
34123-57-4	Monodesmethyl Isoproturon
34123-59-6	Isoproturon
34123-59-6	Isoproturon Solution
34205-21-5	Dimefuron
34256-82-1	Acetochlor
34256-82-1	Acetochlor Solution
34388-29-9	Methothrin
34622-58-7	Orbencarb
34643-46-4	Prothiophos
34643-46-4	Prothiophos Solution
34643-46-4	Prothiophos Solution
34643-47-5	Sulprofos-sulfoxide
34681-10-2	Butocarboxim
34681-10-2	Butocarboxim Solution
34681-23-7	Butoxycarboxim
34681-23-7	Butoxycarboxim Solution
35256-85-0	Tebutam
35256-85-0	Tebutam Solution
35256-85-0	Tebutam Solution
35367-38-5	Diflubenzuron
35367-38-5	Diflubenzuron Solution
35400-43-2	Sulprofos
35400-43-2	Sulprofos Solution
35400-43-2	Sulprofos Solution
35554-44-0	Imazali
35554-44-0	Imazali Solution
35554-44-0	Imazali Solution
35884-76-5	Malathion monocarboxylic acid Solution
35884-76-5	Malathion monocarboxylic acid
36378-61-7	Pirimiphos-ethyl-oxon
36734-19-7	Iprodione
36734-19-7	Iprodione Solution
37019-18-4	Sebutylazine-desethyl
37764-25-3	Dichlorimid
37764-25-3	Dichlorimid Solution
38260-54-7	Etrifos
38641-94-0	Glyphosate-isopropylammonium
38641-94-0	Glyphosate-isopropylammonium Solution
38641-94-0	Glyphosate-isopropylammonium Solution
38727-55-8	Diethyl ethyl
38727-55-8	Diethyl ethyl Solution
38902-68-0	Trietazine-desethyl
38902-68-0	Trietazine-desethyl Solution
39148-24-8	Fosetyl-Aluminum(Technical)
39196-18-4	Thiofanox
39196-18-4	Thiofanox Solution
39222-73-6	Tebuthiuron TP 108
39300-45-3	Dinocap
39300-45-3	Dinocap Solution
39300-45-3	Dinocap Solution
39515-40-7	d-trans-Cyphenothrin
39515-40-7	d-trans-Cyphenothrin Solution
39515-40-7	d-trans-Cyphenothrin Solution
39515-41-8	Fenprothrin
39515-41-8	Fenprothrin Solution
39515-41-8	Fenprothrin Solution
39807-15-3	Oxadiazon-hydroxy
39807-15-3	Oxadiazon-hydroxy Solution
39807-19-7	Oxadiazon-hydroxy
39807-19-7	Oxadiazon-hydroxy Solution
39807-20-0	Pivaloyl 2-(2,4-dichloro-5-hydroxyphenyl)hydrazide Solution
40487-42-1	Pendimethalin
40487-42-1	Pendimethalin Solution
40487-42-1	Pendimethalin Solution
40596-69-8	Methoprene
40596-69-8	Methoprene Solution
40596-69-8	Methoprene Solution
40843-25-2	Diclofop acid
41083-11-8	Azocyclotin (Technical)
41083-11-8	Azocyclotin (Technical) Solution
41198-08-7	Profenofos
41198-08-7	Profenofos Solution
41198-08-7	Profenofos Solution
41394-05-2	Metamitron
41394-05-2	Metamitron Solution
41483-43-6	Bupirimate
41814-78-2	Tricyclazole
41814-78-2	Tricyclazole Solution
41814-78-2	Tricyclazole Solution
41814-78-2	Tricyclazole Solution
42487-72-9	2-Diethylamino-6-methyl pyrimidin-4-ol
42509-80-8	Isazophos
42509-80-8	Isazophos Solution
42576-02-3	Bifenox
42576-02-3	Bifenox Solution
42588-37-4	Kinoprene
42588-37-4	Kinoprene Solution
42588-37-4	Kinoprene Solution
42588-37-4	Kinoprene Solution
42609-52-9	Daimuron
42609-52-9	Daimuron Solution
42874-03-3	Oxyfluorfen
42874-03-3	Oxyfluorfen Solution
43121-43-3	Triadimefon
43121-43-3	Triadimefon Solution
43121-43-3	Triadimefon Solution
43121-43-3	Triadimefon Solution
43171-49-9	2,4-Dichlorobenzhydrol
43222-48-6	Difenoquat methylsulfate
43222-48-6	Difenoquat methylsulfate Solution
43222-48-6	Difenoquat methylsulfate Solution
50471-44-8	Vinclozolin
50471-44-8	Vinclozolin Solution
50471-44-8	Vinclozolin Solution
50512-35-1	Isoprothiolane
50512-35-1	Isoprothiolane Solution
50512-35-1	Isoprothiolane Solution
50563-36-5	Dimethachlor
50594-66-6	Acifluorfen
50594-66-6	Acifluorfen Solution
50594-66-6	Acifluorfen Solution
50594-66-6	Acifluorfen Solution

Part #
MET-12279A-250MG
N-12279-100MG
S-12279A1-1ML
N-12972-10MG
N-11013-100MG
S-11013M1-1ML
N-14037-100MG
N-12951-10MG
N-13577-100MG
S-13577A1-1ML
S-13577U1-1ML
MET-13258C-100MG
N-12963-100MG
S-12963M1-1ML
N-13921-100MG
S-13921M1-1ML
N-13504-250MG
S-13504A1-1ML
S-13504U1-1ML
N-11722-250MG
S-11722A1-1ML
N-13258-100MG
S-13258A1-1ML
S-13258J1-1ML
N-12198-100MG
S-12198A1-1ML
S-12198J1-1ML
MET-12346BT1-1ML
MET-12346D-100MG
MET-13063A-250MG
N-12220-100MG
S-12220A1-1ML
MET-12950A-100MG
N-11662-100MG
S-11662A1-1ML
S-11662U1-1ML
N-12954-10MG
N-12134-250MG
S-12134F1-1ML
S-12134U1-1ML
N-11690-100MG
S-11690M1-1ML
MET-13666A-100MG
MET-13666AA1-1ML
N-12019-100MG
N-14058-25MG
S-14058A1-1ML
N-13948-100MG
N-11783-100MG
S-11783M1-1ML
S-11783M1-5ML
N-11561-100MG
S-11561A1-1ML
S-11561T1-1ML
N-11960-250MG
S-11960A1-1ML
S-11960T1-1ML
N-14011-100MG
S-14011A1-1ML
MET-12731A-50MG
MET-12731AU10-1ML
MET-12731BM1-1ML
N-13140-100MG
S-13140A1-1ML
S-13140K4-1ML
N-12400-100MG
S-12400T4-1ML
S-12400T4-5ML
MET-11677B-100MG
N-11111-250MG
S-11111A1-1ML
N-13097-250MG
S-13097A1-1ML
S-13097J1-1ML
N-12383-250MG
S-12383M1-1ML
N-11349-100MG
S-11349A1-1ML
S-11349T1-1ML
S-11349T1-5ML
MET-13064A-100MG
N-12222-100MG
S-12222U1-1ML
N-11202-250MG
S-11202M1-1ML
N-12295-100MG
S-12295M1-1ML
S-12295T4-1ML
S-12295T4-5ML
N-13955-100MG
S-13955A1-1ML
N-12742-250MG
S-12742M1-1ML
N-13636-500MG
S-13636A1-1ML
S-13636T1-1ML
S-13636T1-5ML
MET-11948A-100MG
N-11721-250MG
S-11721F1-1ML
S-11721U1-1ML
N-13745-250MG
S-13745A1-1ML
S-13745K1-1ML
N-12278-250MG
S-12278A1-1ML
S-12278T1-1ML
N-11754-250MG
N-11027-250MG
S-11027A1-1ML
S-11027B1-1ML
S-11027B1-5ML

CAS#	Description	Part #	CAS#	Description	Part #
50594-67-7	Acifluorfen methyl ester	N-11028-100MG	57153-17-0	Dichlorprop methyl ester Solution	S-11672M1-1ML
50594-67-7	Acifluorfen methyl ester Solution	S-11028Y1-1ML	57153-17-0	Dichlorprop methyl ester Solution	S-11672M1-5ML
50594-67-7	Acifluorfen methyl ester Solution	S-11028Y1-5ML	57153-18-1	MCPB methyl ester	N-12358-250MG
51026-28-9	Busan 40	N-11351-100MG	57153-18-1	MCPB methyl ester Solution	S-12358A1-1ML
51026-28-9	Busan 40 Solution	S-11351A1-1ML	57153-18-1	MCPB methyl ester Solution	S-12358T1-1ML
51026-28-9	Busan 40 Solution	S-11351F1-1ML	57213-69-1	2,4,5-T-triethylamine salt	N-10652-250MG
51218-45-2	Metolachlor	N-12478-100MG	57213-69-1	2,4,5-T-triethylamine salt Solution	S-10652M1-1ML
51218-45-2	Metolachlor Solution	S-12478A1-1ML	57369-32-1	Pyroquilon	N-12998-10MG
51218-45-2	Metolachlor Solution	S-12478T1-1ML	57646-30-7	Furalaxyl	N-12978-10MG
51218-45-2	Metolachlor Solution	S-12478T1-5ML	57837-19-1	Metalaxyl	N-12380-100MG
51218-49-6	Pretlilachlor	N-13089-250MG	57837-19-1	Metalaxyl Solution	S-12380A1-1ML
51218-49-6	Pretlilachlor Solution	S-13089A1-1ML	57837-19-1	Metalaxyl Solution	S-12380U1-1ML
51218-49-6	Pretlilachlor Solution	S-13089U1-1ML	57960-19-7	Acequinocyl	N-11003-250MG
51229-78-8	1-cis-3-Chloroallyl-3,5,7-triaza-1-azonia-adamantane chlorid	N-10045-100MG	57960-19-7	Acequinocyl Solution	S-11003A1-1ML
51229-78-8	1-cis-3-Chloroallyl-357-triaza-1-azonia-adamantane chloride	S-10045M1-1ML	57960-19-7	Acequinocyl Solution	S-11003U1-1ML
51235-04-2	Velpar	N-13743-250MG	57966-95-7	Cymoxanil	N-11544-100MG
51235-04-2	Velpar Solution	S-13743A1-1ML	57966-95-7	Cymoxanil Solution	S-11544A1-1ML
51235-04-2	Velpar Solution	S-13743T1-1ML	57966-95-7	Cymoxanil Solution	S-11544U1-1ML
51235-04-2	Velpar Solution	S-13743T1-5ML	57973-67-8	Flamprop-M-isopropyl	N-13939-250MG
51338-27-3	Diclofop methyl	N-11677-250MG	57973-67-8	Flamprop-M-isopropyl Solution	S-13939M1-1ML
51338-27-3	Diclofop methyl Solution	S-11677A1-1ML	58011-68-0	Pyrazolynate	N-13083-50MG
51338-27-3	Diclofop methyl Solution	S-11677T1-1ML	58011-68-0	Pyrazolynate Solution	S-13083A1-1ML
51630-58-1	Fenvalerate	N-13201-250MG	58138-08-2	Triphane Solution	S-13965A1-1ML
51630-58-1	Fenvalerate Solution	S-13201A4-1ML	58594-72-2	Imazalil sulfate	N-12199-100MG
51630-58-1	Fenvalerate Solution	S-13201A4-5ML	58594-72-2	Imazalil sulfate Solution	S-12199F1-1ML
51630-58-1	Fenvalerate Solution	S-13201U1-1ML	58769-20-3	Kadethrine	N-12947-10MG
51707-55-2	Thidiazuron	N-13564-250MG	58810-48-3	Ofurace Solution	S-12988K10-1ML
51707-55-2	Thidiazuron Solution	S-13564A1-1ML	58905-16-1	1-(2,4-Dichlorophenyl)-2-(1H-1,2,4-triazole-1-yl)-ethanone	MET-13576B-100MG
51707-55-2	Thidiazuron Solution	S-13564U1-1ML	59669-26-0	Thiodicarb	N-13568-250MG
52207-48-4	Dimehypo	N-11752-250MG	59669-26-0	Thiodicarb Solution	S-13568A1-1ML
52207-48-4	Dimehypo Solution	S-11752F1-1ML	59756-60-4	Fluridone	N-13217-100MG
52315-07-8	Cypermethrin	N-11545-100MG	59756-60-4	Fluridone Solution	S-13217T1-1ML
52315-07-8	zeta-Cypermethrin	N-13754-250MG	59756-60-4	Fluridone Solution	S-13217T1-5ML
52315-07-8	Cypermethrin Solution	S-11545A1-1ML	60168-88-9	Fenarimol	N-11948-100MG
52315-07-8	Cypermethrin Solution	S-11545U1-1ML	60168-88-9	Fenarimol Solution	S-11948A1-1ML
52315-07-8	zeta-Cypermethrin Solution	S-11737A1-1ML	60168-88-9	Fenarimol Solution	S-11948T1-1ML
52315-07-8	zeta-Cypermethrin Solution	S-13754U1-1ML	60168-88-9	Fenarimol Solution	S-11948T1-5ML
52645-53-1	Permethrin	N-12848-250MG	60207-31-0	Azaconazole Solution	S-13976A1-1ML
52645-53-1	Permethrin Solution	S-12848A1-1ML	60207-90-1	Propiconazole	N-13576-250MG
52645-53-1	Permethrin Solution	S-12848T1-1ML	60207-90-1	Propiconazole Solution	S-13576A1-1ML
52756-22-6	Flamprop isopropyl	N-11972-100MG	60207-90-1	Propiconazole Solution	S-13576E1-1ML
52756-22-6	Flamprop isopropyl Solution	S-11972A1-1ML	60207-93-4	Etaconazole	N-12985-100MG
52756-22-6	Flamprop isopropyl Solution	S-11972T1-1ML	60238-56-4	Chlorthiophos	N-12969-100MG
52756-25-9	Flamprop-methyl	N-11973-100MG	60397-77-5	N-(2,4-Dimethylphenyl)formamide	MET-11068A-50MG
52756-25-9	Flamprop-methyl Solution	S-11973A1-1ML	60541-86-8	Ethyl-3-amino-2,5-dichlorobenzoate	N-11916-100MG
52756-25-9	Flamprop-methyl Solution	S-11973T1-1ML	60541-86-8	Ethyl-3-amino-2,5-dichlorobenzoate Solution	S-11916A1-1ML
52756-25-9	Flamprop-methyl Solution	S-11973A1-1ML	60541-86-8	Ethyl-3-amino-2,5-dichlorobenzoate Solution	S-11916T1-1ML
52888-80-9	Prosulfocarb	N-13918-10MG	60568-05-0	Furmecyclox	N-12979-100MG
52918-63-5	Deltamethrin	N-11579-250MG	60568-05-0	Furmecyclox Solution	S-12979M1-1ML
52918-63-5	Deltamethrin Solution	S-11579A1-1ML	60825-26-5	Triclopyr methyl ester	N-13661-100MG
52918-63-5	Deltamethrin Solution	S-11579T1-1ML	60825-26-5	Triclopyr methyl ester Solution	S-13661A1-1ML
53112-28-0	Pyrimethanil	N-13159-100MG	60825-26-5	Triclopyr methyl ester Solution	S-13661T1-1ML
53112-28-0	Pyrimethanil Solution	S-13159A1-1ML	61213-25-0	Flurochloridone	N-12987-100MG
53112-28-0	Pyrimethanil Solution	S-13159U1-1ML	61213-25-0	Flurochloridone Solution	S-12987M1-1ML
53179-78-5	3-(2,2-Dibromovinyl)-2,2-dimethyl-(1-cycloprop.)carb.acid(cis)	MET-11579A-10MG	61432-55-1	Dimepiperate Solution	S-13978A1-1ML
53179-78-5	3-(2,2-Dibromovinyl)-2,2-dimethyl-(1-cycloprop.)carb.acid(cis)	MET-11579B1-1ML	61520-53-4	2-[(2-Ethyl-6-methylphenyl)amino]-1-propanol	MET-12478A-50MG
53250-83-2	2-Chloro-4-(methylsulfonyl)benzoic acid	MET-12890A-100MG	61592-45-8	Bentazon methyl derivative	N-11143-50MG
53380-22-6	Ethiofencarb sulfonide	MET-11872A-10MG	61592-45-8	Bentazon methyl derivative Solution	S-11143Y1-1ML
53380-22-6	Ethiofencarb sulfone	MET-11872B1-1ML	61592-45-8	Bentazon methyl derivative Solution	S-11143Y1-5ML
53380-23-7	Ethiofencarb sulfone Solution	MET-11872BA1-1ML	61676-87-7	Cymiazole Solution	S-14051A1-1ML
53404-31-2	Dichlorprop-butoxyethyl ester	N-11673-50MG	61949-76-6	cis-Permethrin	N-11483-50MG
53494-70-5	Endrin ketone	N-11856-10MG	61949-76-6	cis-Permethrin Solution	S-11483A1-1ML
53494-70-5	Endrin ketone Solution	S-11856A1-1ML	61949-76-6	cis-Permethrin Solution	S-11483T1-1ML
53494-70-5	Endrin ketone Solution	S-11856U1-1ML	61949-77-7	trans-Permethrin	N-13620-50MG
53558-25-1	Pyrinuron	N-13738-100MG	61949-77-7	trans-Permethrin Solution	S-13620A1-1ML
53558-25-1	Pyrinuron Solution	S-13738A1-1ML	61949-77-7	trans-Permethrin Solution	S-13620T1-1ML
53558-25-1	Pyrinuron Solution	S-13738U1-1ML	62059-43-2	Endothal monohydrate	N-11852-100MG
53780-34-0	Mefluidide	N-12367-100MG	62059-43-2	Endothal monohydrate Solution	S-11852A1-1ML
53780-34-0	Mefluidide Solution	S-12367M1-1ML	62059-43-2	Endothal monohydrate Solution	S-11852A1-5ML
54406-48-3	Empenthrin	N-13913-10MG	62059-43-2	Endothal monohydrate Solution	S-11852F0-1ML
54460-46-7	Cycloprate	N-11541-500MG	62059-43-2	Endothal monohydrate Solution	S-11852F10-1ML
54460-46-7	Cycloprate Solution	S-11541A1-1ML	62059-43-2	Endothal monohydrate Solution	S-11852F10-5ML
54460-46-7	Cycloprate Solution	S-11541T4-1ML	62059-43-2	Endothal monohydrate Solution	S-11852U1-1ML
54460-46-7	Cycloprate Solution	S-11541T4-5ML	62625-13-2	N-[(2,4-Dichlorophenoxy)acetyl]-DL-methionine	N-12505-10MG
54593-83-8	Chlorethoxyfos	N-11427-250MG	62625-13-2	N-[(2,4-Dichlorophenoxy)acetyl]-DL-methionine Solution	S-12505A1-1ML
54593-83-8	Chlorethoxyfos Solution	S-11427A1-1ML	62625-13-2	N-[(2,4-Dichlorophenoxy)acetyl]-DL-methionine Solution	S-12505T1-1ML
54593-83-8	Chlorethoxyfos Solution	S-11427U1-1ML	62850-32-2	Fenothiocarb Solution	S-13962A1-1ML
54965-21-8	Albendazole	N-12959-10MG	62924-70-3	Flumetralin	N-11988-250MG
55179-31-2	Bitertanol	N-11127-100MG	62924-70-3	Flumetralin Solution	S-11988A1-1ML
55179-31-2	Bitertanol Solution	S-11127A1-1ML	62924-70-3	Flumetralin Solution	S-11988B1-1ML
55179-31-2	Bitertanol Solution	S-11127U1-1ML	62924-70-3	Flumetralin Solution	S-11988B1-5ML
55219-65-3	Triadimenol	N-11129-250MG	62924-70-3	Flumetralin Solution	S-11988K1-1ML
55219-65-3	Triadimenol Solution	S-11129A1-1ML	63284-71-9	Nuarimol	N-13942-100MG
55219-65-3	Triadimenol Solution	S-11129U1-1ML	63284-71-9	Nuarimol Solution	S-13942M1-1ML
55283-68-6	Ethalfuralin	N-11865-250MG	63837-33-2	Diofenolan Solution	S-13979A1-1ML
55283-68-6	Ethalfuralin Solution	S-11865A1-1ML	64249-01-0	Anilofos	N-13131-100MG
55283-68-6	Ethalfuralin Solution	S-11865J4-1ML	64249-01-0	Anilofos Solution	S-13131K1-1ML
55283-68-6	Ethalfuralin Solution	S-11865J4-5ML	64529-56-2	Ethiozin	N-13731-100MG
55285-14-8	Carbosulfan	N-11409-250MG	64529-56-2	Ethiozin Solution	S-13731M1-1ML
55285-14-8	Carbosulfan Solution	S-11409A1-1ML	64628-44-0	Triflumuron	N-13687-100MG
55290-64-7	Dimethipin	N-11757-100MG	64628-44-0	Triflumuron Solution	S-13687A1-1ML
55290-64-7	Dimethipin Solution	S-11757M1-1ML	64628-44-0	Triflumuron Solution	S-13687U1-1ML
55335-06-3	Triclopyr	N-13659-500MG	64700-56-7	Triclopyr 2-butoxyethyl ester	S-13660-250MG
55335-06-3	Triclopyr Solution	S-13659A1-1ML	64709-45-1	Pyrimiphos-methyl-oxon	MET-13064B-250MG
55335-06-3	Triclopyr Solution	S-13659T1-1ML	64902-72-3	Chlorsulfuron	N-11461-100MG
55512-33-9	Pyridate	N-13155-500MG	65731-84-2	beta-Cypermethrin	N-11192-250MG
55512-33-9	Pyridate Solution	S-13155A1-1ML	65731-84-2	beta-Cypermethrin Solution	S-11192A1-1ML
55512-33-9	Pyridate Solution	S-13155T1-1ML	65731-84-2	beta-Cypermethrin Solution	S-11192T1-1ML
55701-05-8	3-(2,2-Dichlorovinyl)-2,2-dimethyl-(1-cycloprop.)carb. acid	MET-12848A1-10MG	65733-16-6	(S)-Methoprene	N-13832-100MG
55701-05-8	3-(2,2-Dichlorovinyl)-2,2-dimethyl-(1-cycloprop.)carb. acid	MET-12848AM1-1ML	65733-18-8	S-Hydroprene	N-13189-50MG
55814-41-0	Mepronil	N-12372-100MG	65733-18-8	S-Hydroprene Solution	S-13189A1-1ML
55814-41-0	Mepronil Solution	S-12372M1-1ML	65733-18-8	S-Hydroprene Solution	S-13189J1-1ML
55861-78-4	Isouron Solution	S-13971A1-1ML	66003-55-2	Alloxydim-Sodium	N-11052-50MG
56070-16-7	Terbufos sulfone	MET-13510A-50MG	66063-05-6	Pencycuron	N-12488-100MG
56073-07-5	Difenacoum	N-11719-100MG	66063-05-6	Pencycuron Solution	S-12488A1-1ML
56073-07-5	Difenacoum Solution	S-11719M1-1ML	66063-05-6	Pencycuron Solution	S-12488U1-1ML
56073-10-0	Brodifacoum	N-11329-50MG	66215-27-8	Cyromazine	N-11550-250MG
56073-10-0	Brodifacoum Solution	S-11329A1-1ML	66215-27-8	Cyromazine Solution	S-11550A1-1ML
56425-91-3	Flurprimidol	N-12002-100MG	66215-27-8	Cyromazine Solution	S-11550B1-1ML
56425-91-3	Flurprimidol Solution	S-12002M1-1ML	66230-04-4	Esfenvalerate	N-11102-100MG
56634-95-8	Bromoxynil-heptanoate	MET-11345A-50MG	66230-04-4	Esfenvalerate Solution	S-11102T1-1ML
57018-04-9	Tolclofos-Methyl	N-13578-500MG	66230-04-4	Esfenvalerate Solution	S-11102T1-1ML
57018-04-9	Tolclofos-Methyl Solution	S-13578A1-1ML	66246-88-6	Penconazole	N-12822-100MG
57018-04-9	Tolclofos-Methyl Solution	S-13578U1-1ML	66246-88-6	Penconazole Solution	S-12822A1-1ML
57153-17-0	Dichlorprop methyl ester	N-11672-100MG			

CAS#	Description	Part #
66246-88-6	Penconazole Solution	S-12822U1-1ML
66332-96-5	Flutolanil	N-12004-250MG
66332-96-5	Flutolanil Solution	S-12004M1-1ML
66441-23-4	Fenoxaprop ethyl	N-11957-100MG
66441-23-4	Fenoxaprop ethyl Solution	S-11957M1-1ML
66840-71-9	DMST Solution	MET-13584AM1-1ML
66841-25-6	Tralomethrin	N-13588-100MG
66841-25-6	Tralomethrin Solution	S-13588A1-1ML
66841-25-6	Tralomethrin Solution	S-13588T1-1ML
66952-49-6	Methasulfocarb	N-13283-25MG
66952-49-6	Methasulfocarb Solution	S-13283A1-1ML
67018-59-1	Pirimiphos-methyl-N-desethyl	MET-13064D-100MG
67129-08-2	Metazachlor	N-12385-100MG
67375-30-8	alpha-Cypermethrin	N-11061-100MG
67375-30-8	alpha-Cypermethrin Solution	S-11061A1-1ML
67375-30-8	alpha-Cypermethrin Solution	S-11061J1-1ML
67485-29-4	Amdro	N-11063-100MG
67485-29-4	Amdro Solution	S-11063A1-1ML
67485-29-4	Amdro Solution	S-11063J1-1ML
67747-09-5	Prochloraz	N-13093-250MG
67747-09-5	Prochloraz Solution	S-13093A1-1ML
67747-09-5	Prochloraz Solution	S-13093U1-1ML
68157-60-8	Forchlorfenuron	N-12010-100MG
68157-60-8	Forchlorfenuron Solution	S-12010A1-1ML
68359-37-5	Cyfluthrin	N-11130-250MG
68359-37-5	beta-Cyfluthrin	N-11191-250MG
68359-37-5	Cyfluthrin Solution	S-11130A4-1ML
68359-37-5	Cyfluthrin Solution	S-11130A4-5ML
68359-37-5	Cyfluthrin Solution	S-11130U1-1ML
68359-37-5	beta-Cyfluthrin Solution	S-11191A1-1ML
68359-37-5	beta-Cyfluthrin Solution	S-11191T1-1ML
68694-11-1	Triflumizole	N-13686-250MG
69327-76-0	Buprofezin	N-11350-100MG
69327-76-0	Buprofezin Solution	S-11350A1-1ML
69327-76-0	Buprofezin Solution	S-11350U1-1ML
69335-90-6	Fluazifop-methyl	N-13906-50MG
69335-90-6	Fluazifop-P-methyl	N-13908-50MG
69335-91-7	Fluazifop	N-11977-10MG
69335-91-7	Fluazifop Solution	S-11977A1-1ML
69335-91-7	Fluazifop Solution	S-11977A1-5ML
69377-81-7	Fluroxypyr	N-12000-100MG
69377-81-7	Fluroxypyr Solution	S-12000M1-1ML
69770-45-2	Flumethrin	N-13139-100MG
69770-45-2	Flumethrin Solution	S-13139A1-1ML
69806-34-4	Haloxifyop (free acid)	N-12143-10MG
69806-50-4	Fluazifop-butyl	N-11978-250MG
69806-50-4	Fluazifop-butyl Solution	S-11978A1-1ML
69806-50-4	Fluazifop-butyl Solution	S-11978U1-1ML
70124-77-5	Flucythrinate	N-11984-100MG
70124-77-5	Flucythrinate Solution	S-11984A1-1ML
70124-77-5	Flucythrinate Solution	S-11984U1-1ML
70288-86-7	Ivermectin	N-12289-100MG
70288-86-7	Ivermectin Solution	S-12289A1-1ML
70288-86-7	Ivermectin Solution	S-12289U1-1ML
70628-36-3	Propachlor OA	MET-13104B-10MG
70628-36-3	Propachlor OA Solution	MET-13104BM1-1ML
70630-17-0	Mefenoxam	N-12365-100MG
70630-17-0	Mefenoxam Solution	S-12365A1-1ML
70630-17-0	Mefenoxam Solution	S-12365T1-1ML
71283-80-2	Fenoxaprop-P-ethyl	N-11958-100MG
71283-80-2	Fenoxaprop-P-ethyl Solution	S-11958A1-1ML
71283-80-2	Fenoxaprop-P-ethyl Solution	S-11958U1-1ML
71422-67-8	Chlorfluazuron	N-11430-50MG
71422-67-8	Chlorfluazuron Solution	S-11430A1-1ML
71561-11-0	Pyrazoxyfen	N-13148-100MG
71561-11-0	Pyrazoxyfen Solution	S-13148M1-1ML
71626-11-4	Benalaxyl	N-11133-250MG
71626-11-4	Benalaxyl Solution	S-11133M1-1ML
71751-41-2	Abamectin	N-10995-100MG
71751-41-2	Abamectin Solution	S-10995A1-1ML
72178-02-0	Fomesafen	N-12008-500MG
72178-02-0	Fomesafen Solution	S-12008M1-1ML
72490-01-8	Fenoxycarb	N-11959-100MG
72490-01-8	Fenoxycarb Solution	S-11959A1-1ML
72619-32-0	Haloxifyop-P methyl	N-12146-50MG
72963-72-5	Imiprothrin	N-14006-50MG
72963-72-5	Imiprothrin Solution	S-14006A1-1ML
73250-68-7	Mefenacet	N-12364-250MG
73250-68-7	Mefenacet Solution	S-12364M1-1ML
74051-80-2	Sethoxydim	N-13210-100MG
74051-80-2	Sethoxydim Solution	S-13210U1-1ML
74070-46-5	Aclofenif	N-13991-250MG
74070-46-5	Aclofenif Solution	S-13991A1-1ML
74115-24-5	Clofentezine	N-11490-100MG
74115-24-5	Clofentezine Solution	S-11490M1-1ML
74222-97-2	Sulfometuron methyl	N-13254-100MG
74222-97-2	Sulfometuron methyl Solution	S-13254U1-1ML
74223-64-6	Metsulfuron methyl	N-12482-100MG
74223-64-6	Metsulfuron methyl Solution	S-12482A1-1ML
74712-19-9	Bromobutide Solution	S-13977M1-1ML
76280-91-6	Tecloftalam	N-14068-10MG
76578-12-6	Quizalofop	MET-13174A-50MG
76578-13-7	Quizalofop methyl	MET-13174B-50MG
76578-14-8	Quizalofop ethyl	N-13174-100MG
76578-14-8	Quizalofop ethyl Solution	S-13174A1-1ML
76578-14-8	Quizalofop ethyl Solution	S-13174U1-1ML
76674-21-0	Flutriaol	N-12005-100MG
76674-21-0	Flutriaol Solution	S-12005A1-1ML
76674-21-0	Flutriaol Solution	S-12005U1-1ML
76703-62-3	gamma-Cyhalothrin	N-12106-100MG
76738-62-0	Paclobutrazol	N-12811-500MG
76738-62-0	Paclobutrazol Solution	S-12811A1-1ML
76738-62-0	Paclobutrazol Solution	S-12811U1-1ML
77182-82-2	Glufosinate-ammonium	N-12111-100MG
77359-18-3	(E)-Fenimzone Solution	S-14031A1-1ML
77458-30-1	1-(4-Chlorophenyl)-1H-pyrazol-4-ol	MET-13142A-100MG
77501-63-4	Lactofen	N-12304-100MG
77501-63-4	Lactofen Solution	S-12304A1-1ML
77501-63-4	Lactofen Solution	S-12304U1-1ML
77501-90-7	Fluoroglycofen-ethyl	N-11999-100MG
77732-09-3	Oxadixyl	N-12732-100MG
77732-09-3	Oxadixyl Solution	S-12732A1-1ML
77732-09-3	Oxadixyl Solution	S-12732U1-1ML
77738-92-2	4-(Ethylthio)-6-methyl-2-(1-methylethyl)-pyrimidine	MET-11621D-10MG

CAS#	Description	Part #
78587-05-0	Hexythiazox	N-12184-100MG
78587-05-0	Hexythiazox Solution	S-12184A1-1ML
79241-46-6	Fluazifop-P-butyl	N-11979-100MG
79241-46-6	Fluazifop-P-butyl Solution	S-11979A1-1ML
79241-46-6	Fluazifop-P-butyl Solution	S-11979U1-1ML
79270-78-3	Dichlorprop 2-ethylhexyl ester	N-13821-250MG
79277-27-3	Thifensulfuron-methyl	N-13056-100MG
79277-27-3	Thifensulfuron-methyl Solution	S-13056A1-1ML
79622-59-6	Fluazinam	N-11980-100MG
79622-59-6	Fluazinam Solution	S-11980A1-1ML
79622-59-6	Fluazinam Solution	S-11980U1-1ML
79983-71-4	Hexaconazole	N-12165-100MG
79983-71-4	Hexaconazole Solution	S-12165A1-1ML
79983-71-4	Hexaconazole Solution	S-12165U1-1ML
80060-09-9	Diafenthion	N-11612-250MG
80060-09-9	Diafenthion Solution	S-11612A1-1ML
80060-09-9	Diafenthion Solution	S-11612U1-1ML
80844-07-1	Etofenprox	N-11941-100MG
80844-07-1	Etofenprox Solution	S-11941A1-1ML
80844-07-1	Etofenprox Solution	S-11941U1-1ML
80854-21-3	Etofenprox-desethyl Solution	MET-11941AM1-1ML
81334-34-1	Imazapyr	N-12203-100MG
81334-34-1	Imazapyr Solution	S-12203A1-1ML
81335-37-7	Imazaquin	N-12204-100MG
81335-37-7	Imazaquin Solution	S-12204A1-1ML
81335-37-7	Imazaquin Solution	S-12204T1-1ML
81335-77-5	Imazethapyr	N-12205-100MG
81335-77-5	Imazethapyr Solution	S-12205M1-1ML
81405-85-8	Imazamethabenz-methyl	N-12200-100MG
81405-85-8	Imazamethabenz-methyl Solution	S-12200A1-1ML
81405-85-8	Imazamethabenz-methyl Solution	S-12200U1-1ML
81406-37-3	Fluroxypyr-meptyl	N-12001-100MG
81406-37-3	Fluroxypyr-meptyl Solution	S-12001M1-1ML
81412-43-3	Tridemorph (Mix of Isomers)	N-13665-100MG
81412-43-3	Tridemorph (Mix of Isomers) Solution	S-13665A1-1ML
81412-43-3	Tridemorph (Mix of Isomers) Solution	S-13665U1-1ML
81777-89-1	Command	N-11495-100MG
81777-89-1	Command Solution	S-11495M1-1ML
82097-50-5	Triasulfuron	N-13638-100MG
82097-50-5	Triasulfuron Solution	S-13638A1-1ML
82558-50-7	Isoxaben	N-12285-250MG
82558-50-7	Isoxaben Solution	S-12285M1-1ML
82560-54-1	Benfuracarb	N-11137-100MG
82657-04-3	Bifenthrin	N-11203-100MG
82657-04-3	Bifenthrin Solution	S-11203A1-1ML
82657-04-3	Bifenthrin Solution	S-11203T1-1ML
83055-99-6	Bensulfuron-methyl	N-11139-100MG
83055-99-6	Bensulfuron-methyl Solution	S-11139A1-1ML
83066-88-0	Fluazifop-P	N-13907-50MG
83121-18-0	Teflubenzuron	N-13506-250MG
83121-18-0	Teflubenzuron Solution	S-13506A1-1ML
83164-33-4	Diffenican	N-11723-100MG
83657-17-4	Uniconazole-P	N-13735-10MG
83657-17-4	Uniconazole-P Solution	S-13735M1-1ML
83657-22-1	Uniconazole	N-14032-50MG
83657-22-1	Uniconazole Solution	S-14032A1-1ML
83657-24-3	Diniconazole	N-11782-100MG
83657-24-3	Diniconazole Solution	S-11782M1-1ML
84030-86-4	Esbiothrin	N-11864-100MG
84087-01-4	Quinclorac	N-13169-500MG
84087-01-4	Quinclorac Solution	S-13169A1-1ML
84087-01-4	Quinclorac Solution	S-13169T1-1ML
84625-61-6	Itraconazole	N-13958-100MG
85509-19-9	Flusilazole	N-12003-100MG
85509-19-9	Flusilazole Solution	S-12003A1-1ML
85509-19-9	Flusilazole Solution	S-12003U1-1ML
86209-51-0	Primisulfuron-methyl	N-13090-100MG
86209-51-0	Primisulfuron-methyl Solution	S-13090A1-1ML
86479-06-3	Hexaflumuron	N-12167-250MG
86479-06-3	Hexaflumuron Solution	S-12167A1-1ML
86479-06-3	Hexaflumuron Solution	S-12167U1-1ML
86598-92-7	Imibenconazole Solution	S-13968A1-1ML
86763-47-5	Propisochlor	N-13823-50MG
87130-20-9	Diethofencarb	N-11691-250MG
87130-20-9	Diethofencarb Solution	S-11691M1-1ML
87392-12-9	S-Metolachlor	N-13191-100MG
87392-12-9	S-Metolachlor Solution	S-13191A1-1ML
87392-12-9	S-Metolachlor Solution	S-13191T1-1ML
87546-18-7	Flumiclorac-pentyl	N-11990-500MG
87546-18-7	Flumiclorac-pentyl Solution	S-11990A1-1ML
87546-18-7	Flumiclorac-pentyl Solution	S-11990T1-1ML
87674-68-8	Dimethenamid	N-11755-100MG
87674-68-8	Dimethenamid Solution	S-11755A1-1ML
87674-68-8	Dimethenamid Solution	S-11755U1-1ML
87818-31-3	Cinmethylin Solution	S-14050M1-1ML
87820-88-0	Tralkoxydim	N-13587-500MG
87820-88-0	Tralkoxydim Solution	S-13587A1-1ML
87820-88-0	Tralkoxydim Solution	S-13587U1-1ML
88283-41-4	Pyriphenox	N-12990-10MG
88283-41-4	Pyriphenox Solution	S-12990M1-1ML
88349-88-6	Quinclotocet	N-13820-50MG
88485-37-4	Fluxofenim	N-13940-250MG
88485-37-4	Fluxofenim Solution	S-13940M1-1ML
88671-89-0	Myclobutanil	N-13261-100MG
88671-89-0	Myclobutanil Solution	S-13261A1-1ML
88671-89-0	Myclobutanil Solution	S-13261U1-1ML
88805-35-0	Prohexadione	N-13099-100MG
89269-64-7	Ferimzone	N-14030-100MG
89269-64-7	Ferimzone Solution	S-14030A1-1ML
89784-60-1	Pyraclofos	N-13142-10MG
89784-60-1	Pyraclofos Solution	S-13142A1-1ML
90035-08-8	Flocoumafen	N-11975-10MG
90035-08-8	Flocoumafen Solution	S-11975A1-1ML
90717-03-6	Quinmerac	N-13925-10MG
90982-32-4	Chlorimuron ethyl	N-11432-100MG
90982-32-4	Chlorimuron ethyl Solution	S-11432A1-1ML
91465-08-6	lambda-Cyhalothrin	N-12307-100MG
91465-08-6	lambda-Cyhalothrin Solution	S-12307A1-1ML
91465-08-6	lambda-Cyhalothrin Solution	S-12307J1-1ML
92642-35-8	Sulprofos sulfone Solution	MET-13258AU1-1ML
92642-35-8	Sulprofos sulfone Solution	MET-13258AU1-5ML
93567-85-2	Trietazine-ethylamine	N-13989-100MG
93697-74-6	Pyrazosulfuron-ethyl	N-13147-100MG
93697-74-6	Pyrazosulfuron-ethyl Solution	S-13147A1-1ML

CAS#	Description	Part #	CAS#	Description	Part #
94125-34-5	Prosulfuron	N-13137-100MG	115086-54-9	Imidacloprid olefin	MET-12206-10MG
94125-34-5	Prosulfuron Solution	S-13137A1-1ML	115852-48-7	Fenoxanil	N-11956-100MG
94361-06-5	Cyproconazole	N-11547-50MG	116255-48-2	Bromuconazole	N-11347-100MG
94361-06-5	Cyproconazole Solution	S-11547A1-1ML	116714-46-6	Novaluron	N-12670-100MG
94361-06-5	Cyproconazole Solution	S-11547E1-1ML	116714-46-6	Novaluron Solution	S-12670A1-1ML
94593-91-6	Cinosulfuron	N-11470-100MG	116714-46-6	Novaluron Solution	S-12670U1-1ML
95266-40-3	Trinexapac-ethyl	N-13706-250MG	117337-19-6	Fluthiacet-methyl	N-13129-25MG
95266-40-3	Trinexapac-ethyl Solution	S-13706A1-1ML	117337-19-6	Fluthiacet-methyl Solution	S-13129A1-1ML
95266-40-3	Trinexapac-ethyl Solution	S-13706T1-1ML	117428-22-5	Picoxystrobin	N-13927-100MG
95465-99-9	Cadusafos	N-11392-100MG	117428-22-5	Picoxystrobin Solution	S-13927A1-1ML
95465-99-9	Cadusafos Solution	S-11392A1-1ML	117718-60-2	Thiazopyr	N-13563-500MG
95465-99-9	Cadusafos Solution	S-11392U1-1ML	117718-60-2	Thiazopyr Solution	S-13563A1-1ML
95617-09-7	Fenoxaprop	N-13830-50MG	117718-60-2	Thiazopyr Solution	S-13563U1-1ML
95737-68-1	Pyriproxyfen	N-13160-100MG	118134-30-8	Spiroxamine	N-13226-100MG
95737-68-1	Pyriproxyfen Solution	S-13160M1-1ML	118134-30-8	Spiroxamine Solution	S-13226U1-1ML
95977-29-0	Haloxypyr-P (free acid)	N-12145-50MG	118712-89-3	Transfluthrin	N-13626-100MG
96182-53-5	Tebupirimfos	N-13503-100MG	119168-77-3	Tebufenpyrad	N-13502-100MG
96182-53-5	Tebupirimfos Solution	S-13503A1-1ML	119446-68-3	Difenoconazole	N-11720-250MG
96182-53-5	Tebupirimfos Solution	S-13503U1-1ML	119446-68-3	Difenoconazole Solution	S-11720M1-1ML
96489-71-3	Pyridaben	N-13152-25MG	119515-38-7	Icaridin	N-12197-100MG
96489-71-3	Pyridaben Solution	S-13152M1-1ML	119515-38-7	Icaridin Solution	S-12197A1-1ML
96491-05-3	Thenylchlor Solution	S-13966A1-1ML	119791-41-2	Emamectin	N-11847-100MG
96525-23-4	Flurtamone Solution	S-13981A1-1ML	119791-41-2	Emamectin Solution	S-11847A1-1ML
97780-06-8	Ethametsulfuron-methyl	N-11866-100MG	120067-83-6	Fipronil sulfide	MET-11971B-25MG
97780-06-8	Ethametsulfuron-methyl Solution	S-11866A1-1ML	120067-83-6	Fipronil sulfide Solution	MET-11971BA1-1ML
97886-45-8	Dithiopyr	N-11822-100MG	120068-36-2	Fipronil sulfone	MET-11971A-50MG
97886-45-8	Dithiopyr Solution	S-11822A1-1ML	120068-36-2	Fipronil sulfone Solution	MET-11971AA1-1ML
97886-45-8	Dithiopyr Solution	S-11822U1-1ML	120068-37-3	Fipronil	N-11971-100MG
98730-04-2	Benoxacor	N-12957-10MG	120068-37-3	Fipronil Solution	S-11971A1-1ML
98886-44-3	Fosthiazate	N-12020-30MG	120068-37-3	Fipronil Solution	S-11971U1-1ML
98886-44-3	Fosthiazate Solution	S-12020A1-1ML	120116-88-3	Cyazofamid	N-11519-100MG
98886-44-3	Fosthiazate Solution	S-12020U1-1ML	120162-55-2	Azimsulfuron Solution	S-13975A1-1ML
98967-40-9	Flumetsulam	N-11989-100MG	120868-66-8	Imidacloprid urea	MET-12206A-50MG
98967-40-9	Flumetsulam Solution	S-11989A1-1ML	120923-37-7	Amidosulfuron	N-11949-100MG
98967-40-9	Flumetsulam Solution	S-11989U1-1ML	120928-09-8	Fenazaquin	S-11949A1-1ML
99105-77-8	Sulcotriane	N-12890-100MG	120928-09-8	Fenazaquin Solution	S-11949A1-1ML
99129-21-2	Clethodim	N-11488-100MG	121552-61-2	Cyprodinil	N-11548-250MG
99607-70-2	Clointocet-mexyl	N-11491-100MG	121552-61-2	Cyprodinil Solution	S-11548A1-1ML
99675-03-3	Isofenphos-Methyl	N-12250-50MG	121552-61-2	Cyprodinil Solution	S-11548U1-1ML
100646-51-3	Quizalofop-P-ethyl	N-13175-100MG	121776-33-8	Furilazole	N-13817-10MG
100646-51-3	Quizalofop-P-ethyl Solution	S-13175A1-1ML	122008-78-0	Cyhalofop	N-14038-25MG
100646-51-3	Quizalofop-P-ethyl Solution	S-13175U1-1ML	122008-78-0	Cyhalofop Solution	S-14038A1-1ML
100728-84-5	Imazamethabenz (free acid) Solution	MET-12200AA0-1ML	122008-85-9	Cyhalofop-butyl	N-11543-50MG
100784-20-1	Halosulfuron-methyl	N-12142-100MG	122008-85-9	Cyhalofop-butyl Solution	S-11543A1-1ML
100784-20-1	Halosulfuron-methyl Solution	S-12142A1-1ML	122008-85-9	Cyhalofop-butyl Solution	S-11543U1-1ML
101007-06-1	Acrinathrin	N-13928-100MG	122453-73-0	Chlorfenapyr	N-11428-100MG
101007-06-1	Acrinathrin Solution	S-13928A1-1ML	122453-73-0	Chlorfenapyr Solution	S-11428A1-1ML
101200-48-0	Tribenuron methyl	N-13640-250MG	122453-73-0	Chlorfenapyr Solution	S-11428U1-1ML
101200-48-0	Tribenuron methyl Solution	S-13640A1-1ML	122548-33-8	Imazosulfuron	N-13818-100MG
101205-02-1	Cycloxydime	N-13916-10MG	122548-33-8	Imazosulfuron	S-13818A1-1ML
101463-69-8	Flufenoxuron	N-11987-50MG	122836-35-5	Sulfentrazone	N-13253-250MG
101463-69-8	Flufenoxuron Solution	S-11987A1-1ML	122836-35-5	Sulfentrazone Solution	S-13253M1-1ML
102851-06-9	tau-Fluvalinate	N-13263-100MG	122931-48-0	Rimsulfuron	N-13180-100MG
102851-06-9	tau-Fluvalinate Solution	S-13263A1-1ML	123312-89-0	Pymetrozine	N-13141-250MG
102851-06-9	tau-Fluvalinate Solution	S-13263K1-1ML	123312-89-0	Pymetrozine Solution	S-13141A1-1ML
103055-07-8	Lufenuron	N-12325-100MG	123312-89-0	Pymetrozine Solution	S-13141U1-1ML
103055-07-8	Lufenuron Solution	S-12325A1-1ML	123343-16-8	Pyrithiobac-Sodium	N-13161-500MG
103361-09-7	Flumioxazin	N-11991-100MG	123343-16-8	Pyrithiobac-Sodium Solution	S-13161F1-1ML
103361-09-7	Flumioxazin Solution	S-11991A1-1ML	124427-60-7	4-(2-phenylethoxy)-quinazoline	N-12908-10MG
103361-09-7	Flumioxazin Solution	S-11991U1-1ML	124495-18-7	Quinoxifen	N-13173-100MG
103564-78-9	2,4-D propylene glycol ester	N-10620-250MG	124495-18-7	Quinoxifen Solution	S-13173A1-1ML
103564-78-9	2,4-D propylene glycol ester Solution	S-10620A1-1ML	124495-18-7	Quinoxifen Solution	S-13173U1-1ML
103564-78-9	2,4-D propylene glycol ester Solution	S-10620J1-1ML	125116-23-6	Metconazole	N-12386-100MG
104030-54-8	Carproamid Solution	S-13974A1-1ML	125401-92-5	Bispyribac Sodium	N-11323-100MG
104040-78-0	Flazasulfuron	N-11974-50MG	126535-15-7	Triflusalufuron-methyl	N-13690-50MG
104098-48-8	Imazapic	N-12202-250MG	126535-15-7	Triflusalufuron-methyl Solution	S-13690A1-1ML
104098-48-8	Imazapic Solution	S-12202A1-1ML	126801-58-9	Ethoxysulfuron	N-11878-100MG
104206-82-8	Mesotrione	N-12377-100MG	126833-17-8	Fenhexamid	N-11954-250MG
104273-73-6	Trinexapac	MET-13706A-100MG	126833-17-8	Fenhexamid Solution	S-11954A1-1ML
105024-66-6	Silaflofen	N-13914-10MG	127277-53-6	Prohexadione-Calcium	N-13100-100MG
105512-06-9	Clodinafop-propargyl	N-11489-100MG	128639-02-1	Carfentrazone-ethyl	N-11411-10MG
105512-06-9	Clodinafop-propargyl Solution	S-11489A1-1ML	128639-02-1	Carfentrazone-ethyl Solution	S-11411A1-1ML
105512-06-9	Clodinafop-propargyl Solution	S-11489U1-1ML	128639-02-1	Carfentrazone-ethyl Solution	S-11411U1-1ML
105779-78-0	Pyrimidifen Solution	S-13963A1-1ML	129558-76-5	Tolpyrpyrad	N-13579-100MG
106917-52-6	Flusulfamide	N-14029-100MG	129630-17-7	Pyraflufen	MET-13144A-10MG
106917-52-6	Flusulfamide Solution	S-14029A1-1ML	129630-19-9	Pyraflufen-ethyl	N-13144-100MG
107534-96-3	Tebuconazole	N-12006-250MG	129909-90-6	Amicarbazone	N-13075-100MG
107534-96-3	Tebuconazole Solution	S-12006A1-1ML	129909-90-6	Amicarbazone Solution	S-13075A1-1ML
107534-96-3	Tebuconazole Solution	S-12006K1-1ML	130000-40-7	Thifluzamide	N-12884-100MG
109293-97-2	Diflufenzopyr	N-11724-100MG	130000-40-7	Thifluzamide Solution	S-12884M1-1ML
109293-98-3	Diflufenzopyr sodium salt	N-13822-100MG	131341-86-1	Fludioxonil	N-11985-100MG
110235-47-7	Mepanipyrim	N-12369-100MG	131341-86-1	Fludioxonil Solution	S-11985A1-1ML
110488-70-5	Dimethomorph	N-11759-100MG	131341-86-1	Fludioxonil Solution	S-11985U1-1ML
110488-70-5	Dimethomorph Solution	S-11759A1-1ML	131807-57-3	Famoxadone	N-11943-100MG
111479-05-1	Propaquizafop	N-14007-250MG	131807-57-3	Famoxadone Solution	S-11943A1-1ML
111479-05-1	Propaquizafop Solution	S-14007A1-1ML	131807-57-3	Famoxadone Solution	S-11943U1-1ML
111697-06-4	N,N-diethyl-4-methoxy-6-methyl-2-Pyrimidinamine	MET-13064E-100MG	131860-33-8	Azoxystrobin	N-11113-100MG
111872-58-3	Halfenprox	N-13079-10MG	131860-33-8	Azoxystrobin Solution	S-11113A1-1ML
111872-58-3	Halfenprox Solution	S-13079M1-1ML	131860-33-8	Azoxystrobin Solution	S-11113U1-1ML
111988-49-9	Thiacloprid	N-13561-100MG	131983-72-7	Triticonazole	N-13729-100MG
111988-49-9	Thiacloprid Solution	S-13561A1-1ML	132827-25-9	STB	MET-11138B-500MG
111988-49-9	Thiacloprid Solution	S-13561U1-1ML	133220-30-1	Indanofan	N-12209-100MG
111991-09-4	Nicosulfuron	N-10998-100MG	133855-98-8	Epoxiconazole	N-11859-100MG
112143-82-5	Triazamate	N-13956-100MG	134098-61-6	Fenpyroximate	N-11962-100MG
112226-61-6	Halofenozide	N-12141-500MG	134098-61-6	Fenpyroximate Solution	S-11962A1-1ML
112226-61-6	Halofenozide Solution	S-12141M1-1ML	134098-61-6	Fenpyroximate Solution	S-11962U1-1ML
112281-77-3	Tetraconazole	N-13538-100MG	135158-54-2	Acibenzolar-S-methyl	N-11026-250MG
112281-77-3	Tetraconazole Solution	S-13538A1-1ML	135158-54-2	Acibenzolar-S-methyl Solution	S-11026A1-1ML
112281-77-3	Tetraconazole Solution	S-13538U1-1ML	135158-54-2	Acibenzolar-S-methyl Solution	S-11026U1-1ML
112410-23-8	Tebufenozide	N-13501-100MG	135410-20-7	Acetamidiprid	N-11009-100MG
112410-23-8	Tebufenozide Solution	S-13501A1-1ML	135410-20-7	Acetamidiprid Solution	S-11009A1-1ML
113136-77-9	Cyclanilide	N-11521-100MG	135410-20-7	Acetamidiprid Solution	S-11009U1-1ML
113136-77-9	Cyclanilide Solution	S-11521A1-1ML	135590-91-9	Mefenpyr-diethyl	N-12366-100MG
113136-77-9	Cyclanilide Solution	S-11521U1-1ML	135590-91-9	Mefenpyr-diethyl Solution	S-12366A1-1ML
113158-40-3	Fenoxaprop-P	MET-13802-100MG	135590-91-9	Mefenpyr-diethyl Solution	S-12366U1-1ML
113210-98-3	Z-Dimethomorph	N-13752-100MG	136849-15-5	Cyclosulfamuron	N-11542-250MG
113210-98-3	Z-Dimethomorph Solution	S-13752M1-1ML	136849-15-5	Cyclosulfamuron Solution	S-11542A1-1ML
113561-03-8	Sethoxydim-desethoxy	MET-13210A-50MG	136849-15-5	Cyclosulfamuron Solution	S-11542U1-1ML
113830-50-5	cis-Decamethrinic acid methyl ester	MET-11579C-10MG	137641-05-5	Picolinafen Solution	S-13982M1-1ML
114311-32-9	Imazamox	N-12201-250MG	138261-41-3	Imidacloprid	N-12206-100MG
114311-32-9	Imazamox Solution	S-12201A1-1ML	138261-41-3	Imidacloprid Solution	S-12206M1-1ML
114311-32-9	Imazamox Solution	S-12201T1-1ML	139968-49-3	Metaflumizone	N-12379-100MG
114369-43-6	Fenbuconazole	N-11950-100MG	140923-17-7	Iprovalicarb	N-13924-10MG
114369-43-6	Fenbuconazole Solution	S-11950M1-1ML	140939-15-7	Alachlor ESA sodium salt	MET-11043D-25MG

CAS#	Description	Part #	CAS#	Description	Part #
140939-15-7	Alachlor ESA sodium salt Solution	MET-11043DM1-1ML	211867-47-9	Flumorph Solution	S-13280A1-1ML
141112-29-0	Isoxaflutole	N-12287-100MG	213464-77-8	Orthosulfamuron Solution	S-13944B1-1ML
141112-29-0	Isoxaflutole Solution	S-12287U1-1ML	219714-96-2	Penoxsulam	N-12823-50MG
141517-21-7	Trifloxystrobin	N-13684-100MG	220899-03-6	Metrafenone	N-12897-100MG
141517-21-7	Trifloxystrobin Solution	S-13684A1-1ML	221667-31-8	Cyprosulfamide	N-11549-100MG
141776-32-1	Sulfosulfuron	N-13255-100MG	222966-66-7	BHC (alpha isomer) (13C6) Solution	S-FC102S-1.2ML
141776-32-1	Sulfosulfuron Solution	S-13255A1-1ML	223580-51-6	Tiadinil	N-13282-25MG
141776-32-1	Sulfosulfuron Solution	S-13255U1-1ML	223580-51-6	Tiadinil Solution	S-13282A1-1ML
142459-58-3	Flufenacet	N-11986-100MG	229977-93-9	Fluacrypyrim Solution	S-13980M1-1ML
142459-58-3	Flufenacet Solution	S-11986M1-1ML	239110-15-7	Fluopicolide	N-11993-100MG
142994-06-7	2-Methylsulfonyl-4-trifluoromethyl benzoic acid	MET-12287B-100MG	243973-20-8	Pinoxaden	N-13057-100MG
143390-89-0	Kresoxim-methyl	N-12296-250MG	252564-94-6	Cyahalofop-diacid	MET-14038A-100MG
143390-89-0	Kresoxim-methyl Solution	S-12296A1-1ML	256412-89-2	Metamifop	N-13084-25MG
143390-89-0	Kresoxim-methyl Solution	S-12296U1-1ML	272451-65-7	Flubendiamide	N-11981-100MG
143701-75-1	1-(Methylsulfonyl-4-trifluoromethyl phenyl)-2-cyano-3-cyclop	MET-12287A-100MG	283159-90-0	Valifenalate	N-13281-25MG
143807-66-3	Chromafenozide Solution	S-13973A1-1ML	283159-90-0	Valifenalate Solution	S-13281A1-1ML
144550-36-7	Iodosulfuron-methyl-sodium	N-13923-100MG	283594-90-1	Spiromesifen	N-13224-100MG
145701-21-9	Diclosulam	N-11679-250MG	285138-81-0	Chlorpyrifos (diethyl-d10) Solution	S-FD2057S-1.2ML
145701-21-9	Diclosulam Solution	S-11679A1-1ML	317815-83-1	Thiencarbazone-methyl	N-13565-50MG
145701-23-1	Florasulam	N-13922-100MG	318290-98-1	Fluensulfone	N-14084-100MG
147150-35-4	Cloransulam-methyl	N-11492-100MG	332855-88-6	Fenamidone Metabolite	MET-11946-100MG
147150-35-4	Cloransulam-methyl Solution	S-11492A1-1ML	335104-84-2	Tembotrione	N-13508-100MG
147150-35-4	Cloransulam-methyl Solution	S-11492U1-1ML	337458-27-2	Pyrifluquinazon	N-13158-100MG
148476-30-6	Spiromesifen alcohol	MET-13224-100MG	348635-87-0	Amisulbrom Solution	S-13960M1-1ML
148477-71-8	Spirodiclofen	N-13223-100MG	350820-04-1	Parathion (TM) (diethyl-d10) Solution	S-FD998S-1.2ML
149508-90-7	Simeconazole Solution	S-13931A1-1ML	361377-29-9	Fluoxastrobin	N-13816-100MG
149877-41-8	Bifenazate	N-11201-50MG	365400-11-9	Pyrasulfotole	N-14001-100MG
149877-41-8	Bifenazate Solution	S-11201A1-1ML	372137-35-4	Saflufenacil	N-13195-100MG
149877-41-8	Bifenazate Solution	S-11201U1-1ML	374726-62-2	Mandipropamid	N-12354-100MG
149961-52-4	Dimoxystrobin Solution	S-13926A1-1ML	400882-07-7	Cyflumetofen	N-13972-50MG
149979-41-9	Tepraloxymid Solution	S-13970A1-1ML	400882-57-7	2-Methoxyethyl 2-[4-(tert-butyl)phenyl]cyanoacetate	MET-13987A-50MG
150114-71-9	Aminopyralid	N-11066-100MG	413594-06-6	Sethoxydim-oxazole	MET-13210B-50MG
150824-47-8	Nitenpyram	N-12657-100MG	422556-08-9	Pyroxulam	N-13164-100MG
150824-47-8	Nitenpyram Solution	S-12657M1-1ML	447399-55-5	Pyroxasulfone	N-13070-100MG
152019-73-3	Metolachlor OA	MET-12478D-10MG	447399-55-5	Pyroxasulfone Solution	S-13070M1-1ML
152019-73-3	Metolachlor OA Solution	MET-12478DM1-1ML	494793-67-8	Penflufen	N-14035-100MG
152019-74-4	sec-Acetochlor OA Solution	S-13946M0-1ML	494793-67-8	Penflufen Solution	S-14035A1-1ML
153233-91-1	Etoxazole	N-11942-50MG	500008-45-7	Chlorantraniliprole	N-11422-100MG
153233-91-1	Etoxazole Solution	S-11942A1-1ML	658066-35-4	Fuopyram	N-13988-50MG
153233-91-1	Etoxazole Solution	S-11942U1-1ML	658066-35-4	Fuopyram Solution	S-13988A1-1ML
153719-23-4	Thiamethoxam	N-13562-100MG	676228-91-4	Thiacloprid-amide	MET-13561-100MG
153719-23-4	Thiamethoxam Solution	S-13562M1-1ML	736994-63-1	Cyantraniliprole	N-12886-10MG
155569-91-8	Emamectin benzoate	N-12920-100MG	858956-08-8	Aminocyclopyrachlor	N-13826-250MG
156052-68-5	Zoxamide	N-13762-100MG	874967-67-6	Sedaxane	N-13068-100MG
158062-67-0	Flonicamid	N-11976-50MG	874967-67-6	Sedaxane Solution	S-13068M1-1ML
158062-67-0	Flonicamid Solution	S-11976A1-1ML	875915-78-9	Isofetamid	N-14079-100MG
158062-67-0	Flonicamid Solution	S-11976U1-1ML	881685-58-1	Isoprazam Solution	S-13996A1-1ML
158063-66-2	4-(Trifluoromethyl)nicotinic acid	MET-11976A-100MG	907204-31-3	Fluxapyroxad	N-12921-100MG
158237-07-1	Fentrazamide	N-14059-100MG	925689-10-7	3-(difluoromethyl)-1-methyl-1H-Pyrazole-4-carboxamide Solution	MET-13068AA1-1ML
158237-07-1	Fentrazamide Solution	S-14059A1-1ML	946578-00-3	Sulfoxaflor	N-12883-10MG
161050-58-4	Methoxyfenozide	N-12406-100MG	946578-00-3	Sulfoxaflor Solution	S-12883M1-1ML
161050-58-4	Methoxyfenozide Solution	S-12406A1-1ML	947601-85-6	Metolachlor ESA sodium salt	MET-12478C-25MG
161326-34-7	Fenamidone	N-11946-100MG	947601-85-6	Metolachlor ESA sodium salt Solution	MET-12478CM1-1ML
161326-34-7	Fenamidone Solution	S-11946A1-1ML	947601-87-8	Flufenacet ESA sodium salt	MET-11986A-25MG
161326-34-7	Fenamidone Solution	S-11946U1-1ML	947601-87-8	Flufenacet ESA sodium salt Solution	MET-11986AM1-1ML
163165-75-1	Atrazine (ethylamine-d5)	N-FD2208-5-5MG	947601-88-9	Propachlor ESA sodium salt	MET-13104A-25MG
163515-14-8	Dimethenamid-P	N-11756-100MG	947601-88-9	Propachlor ESA sodium salt Solution	MET-13104AM1-1ML
163515-14-8	Dimethenamid-P Solution	S-11756A1-1ML	950782-86-2	Indaziflam	N-13074-100MG
163515-14-8	Dimethenamid-P Solution	S-11756U1-1ML	950782-86-2	Indaziflam Solution	S-13074A1-1ML
163520-33-0	Isoxadifen-ethyl	N-12286-100MG	951659-40-8	Flupyradifurone	N-13941-100MG
165252-70-0	Dinotefuran	N-11788-50MG	951659-40-8	Flupyradifurone Solution	S-13941A1-1ML
165252-70-0	Dinotefuran Solution	S-11788M1-1ML	1048373-72-3	2,4-D Choline salt	N-13985-100MG
165252-70-0	Dinotefuran Solution	S-11788U1-1ML	1048373-72-3	2,4-D Choline salt Solution	S-13985F1-1ML
166407-15-4	Alachlor acetate	MET-11043B-100MG	1072957-71-1	Benzovindiflupyr Solution	S-14010X1-1ML
168088-61-7	Pyribenzoxim	N-14070-100MG	1173022-75-7	Butachlor ESA sodium salt	MET-11353A-25MG
168316-95-8	Spinosad	N-13222-100MG	1254304-22-7	Fluazaindolizine Solution	S-14042A1-1ML
168316-95-8	Spinosad Solution	S-13222A1-1ML	205939-58-8	Dimethenamid ESA sodium salt	MET-11755A-25MG
171103-04-1	N-Desmethyldimethoxam	MET-13562A-50MG	205939-58-8	Dimethenamid ESA sodium salt Solution	MET-11755AM1-1ML
171103-04-1	N-Desmethyldimethoxam Solution	MET-13562AA1-1ML			
171262-17-2	Alachlor OA	MET-11043E-10MG			
171262-17-2	Alachlor OA Solution	MET-11043EM1-1ML			
173159-57-4	Foramsulfuron	N-12009-100MG			
173159-57-4	Foramsulfuron Solution	S-12009A1-1ML			
173584-44-6	S-Indoxacarb	N-12213-100MG			
173584-44-6	S-Indoxacarb Solution	S-12213A1-1ML			
173584-44-6	S-Indoxacarb Solution	S-12213U1-1ML			
175013-18-0	Pyraclostrobin	N-13143-100MG			
175013-18-0	Pyraclostrobin Solution	S-13143A1-1ML			
176969-34-9	3-(Difluoromethyl)-1-methyl-1H-pyrazole-4-carboxylic acid	MET-13068BA1-1ML			
177406-68-7	Benthiavalicarb-isopropyl	N-14076-50MG			
178928-70-6	Prothioconazole	N-13138-100MG			
178928-70-6	Prothioconazole Solution	S-13138B1-1ML			
179101-81-6	Pyridalyl	N-13153-50MG			
180409-60-3	Cyflufenamid	N-13132-25MG			
181274-15-7	Propoxycarbazone-sodium	N-13119-100MG			
181274-15-7	Propoxycarbazone-sodium Solution	S-13119A1-1ML			
181274-17-9	Flucarbazone-sodium	N-11982-250MG			
181274-17-9	Flucarbazone-sodium Solution	S-11982F1-1ML			
181587-01-9	Ethiprole	N-12870-100MG			
183675-82-3	Penthiopyrad	N-17917-50MG			
187022-11-3	Acetochlor ESA sodium salt	MET-11013A-25MG			
187022-11-3	Acetochlor ESA sodium salt Solution	MET-11013AM1-1ML			
187166-40-1	Spinetoram	N-13221-100MG			
188425-85-6	Boscalid	N-11328-100MG			
188425-85-6	Boscalid Solution	S-11328A1-1ML			
189278-12-4	Proquinazid	N-14080-50MG			
190604-92-3	Acetamidrid-N-desmethyl	MET-11009A-10MG			
194992-44-4	Acetochlor OA	MET-11013B-10MG			
194992-44-4	Acetochlor OA Solution	MET-11013BM1-1ML			
199119-58-9	Trifloxysulfuron-Sodium	N-13685-100MG			
201668-31-7	Flufenacet OA	MET-11986B-10MG			
201668-31-7	Flufenacet OA Solution	MET-11986BM1-1ML			
203313-25-1	Spirotetramat	N-13225-100MG			
203313-25-1	Spirotetramat Solution	S-13225B1-1ML			
203448-76-4	2-(1,3-dimethylbutyl)-Benzenamine	MET-14035A-10MG			
203645-53-8	Dichlorvos (dimethyl-d6)	N-FD2061-E-0.01G			
205650-65-3	Fipronil desulfanyl Solution	S-13920B1-1ML			
207233-95-2	Sodium dimethyldithiocarbamate hydrate	N-12863-100MG			
208465-21-8	Mesosulfuron-methyl	N-12376-100MG			
208465-21-8	Mesosulfuron-methyl Solution	S-12376A1-1ML			
210631-68-8	Topramezone	N-13585-100MG			
210880-92-5	Clothianidin	N-11493-100MG			
210880-92-5	Clothianidin Solution	S-11493A1-1ML			
211867-47-9	Flumorph	N-13280-25MG			

Product Documentation

SAFETY DATA SHEETS (SDS) - Federal Regulation (OSHA) 29 CFR 1910.1200 and other applicable state laws require us to furnish a Safety Data Sheet with the initial delivery of each hazardous chemical. However, as a service to you, we supply a SDS with every order. Where hazards are known, appropriate warnings are noted on the label. Always READ your Safety Data Sheet before handling chemicals.

CERTIFICATE OF YOUR ANALYSIS (COA) - With the purchase of any Chem Service "Standards Grade" environmental, pesticide, organic chemicals, solutions or mixtures, a COA will be provided free of charge, where available. Please note that "Standards Grade" products assigned a purity of tech, will not receive a COA. COAs are not available for "Tech Grade" and "No Grade" products, but, where available, these products will come with an IR.

SAFETY DATA SHEETS (SDS) - English language SDS's are shipped in the box with the chemicals, exceeding OSHA's requirements. Additional copies, replacement copies, updated versions, SDS in advance of an order, or additional languages, can be downloaded from our website www.chemservice.com.

Replacement Policy

"If it's on the label, we guarantee it."

Prior to the expiration date shown on the label and exclusive of any customer contamination, we will replace any products showing purity degradation. You must advise our Customer Service Department immediately to receive approval and instruction on either disposal or return of the product.

Shelf Life Information

Shelf life at Chem Service is part of our warranty. Our normally assigned shelf life is eighteen (18) months or longer. However, a standard may have a **short shelf life** (less than 12 months) assigned due to a purity or a stability problem. If Chem Service offers a standard with a short shelf life, we notify the buyer of the material's short shelf life at the time of purchase.

Contact us at:



660 Tower Lane • PO Box 599
West Chester, PA 19381-0599

INTERNET and E-MAIL

Web Site

www.chemservice.com

General

info@chemservice.com

Marketing

marketing@chemservice.com

Orders, Pricing and Availability

quote@chemservice.com

Technical Inquiries & Customer Service:

tech@chemservice.com

ORDERS

Phone:

1-610-692-3026 or 1-800-452-9994

Fax: 1-610-692-8729

TECHNICAL INQUIRIES & CUSTOMER SERVICE

Phone: 1-610-692-3026 or 1-800-452-9994

Fax: 1-610-692-8729



Return Policy

Goods may not be returned for credit except with Seller's written permission and then only in strict compliance with Seller's return shipment instructions, since most of Seller's products are classified as Hazardous Materials (HazMats) by the U.S. D.O.T. or Dangerous Goods under the applicable international conventions (ICAO/IATA or IMDG). It is imperative that Buyer not return anything without specific instructions from Seller to prevent Buyer's violation of applicable laws and regulations. Buyer will at times handle and dispose of goods in accordance with Seller's instructions, the current Safety Data Sheet, and all applicable law and government regulations.

Buyer has 15 days to notify Seller of any damage or shortage, and 30 days to notify Seller of item ordered incorrectly. If a return is necessary, a return authorization number (RO#) must be given by the Customer Service Manager, and the item(s) must be returned within 10 days after issuance of the RO#.

Unauthorized returns may result in an assessment of a charge to the Buyer, equal to Chem Service's actual expense in handling the unauthorized return. In addition, Buyer may incur regulatory or legal sanctions for illegal shipments and will be held liable by Chem Service for any expenses, including legal fees, fines, or other penalties incurred by Chem Service as a result of Buyer's unauthorized shipment.

The restocking fee is 30% or a \$10.00 minimum and is decided upon by Seller.

Grade Information

"STANDARDS GRADE" (S.G.) - Material has been analyzed and assigned a certified purity and valid expiration date. The analysis is traceable by lot number. Minimum purity requirement is 98% for neat materials, except in a few cases where it is not physically or economically possible (a few of these items may be assigned a purity of "Tech" due to the nature of the compound - not to be confused with "Technical Grade"). Acceptable as a standard by various government agencies.

"NO Grade" (S.G.) - Material is 95% or less. Neither an analysis or certified purity is available. However, many "No Grade" products will be provided with an Infrared Spectrum (IR). When an adequate purity can be assigned to "No Grade" material, they will be changed to "Standard Grade" and will come with the appropriate documentation.

Products furnished by Chem Service, Incorporated are **FOR LABORATORY USE ONLY!** Our products may **NOT BE USED** as drugs, cosmetics, agricultural or pesticidal products, food additives or as household chemicals

WARNING - Always read your Safety Data Sheet before handling or using any chemicals.

General Information

Extended Business Hours - Chem Service's business hours are Monday - Thursday 8:00 AM to 6:00 PM ET and Friday 8:00am - 5pm ET for ordering and support of all Chem Service products.

Rush Hours - Place any rush order by 4:00 PM. ET and have guaranteed same day shipment of any in-stock Chem Service product. With our extensive 16,000 plus inventory of reference materials, Chem Service has the products you need when you need them.

Purity Minimums - Chem Service's policy is to provide all neat Standards Grade materials that meet the 98% purity minimum. There will be a few exceptions for a small number of compounds where the chemical stability or the nature of the compound makes this purity unobtainable or impractical, due to chemical stability, isomer separation, purification, synthesis, cost or other factors, but nearly all Standards Grade materials will ship as 98% or higher purity.

Certified Neat Products and Certified Custom Solutions - Our certified in-house neat chemicals are the basis of our quality custom solutions. We can tailor a product to meet your specification, certification, budget and delivery requirements. Use our Custom Solution & Mixture Request form on page 8.

Prices - Are Subject to Change Without Notice.



NOTICE: THIS LANGUAGE FORMS YOUR CONTRACT WITH CHEM SERVICE, INC.

1. **ACCEPTANCE GOVERNING PROVISIONS** - Buyer places any order subject to the terms and conditions herein contained, whether additional to or different from those contained in Buyer's purchase order or any other form or document heretofore or hereafter supplied by Buyer to Seller. Seller's failure to object to provisions contained in any purchase order or any other form or document from Buyer shall not be construed as a waiver of these terms and conditions and not an acceptance of any such provision. If Buyer places an order from Seller's catalogue or catalogue supplement, Buyer is deemed to have assented to these terms. If Buyer places an order without reference to the catalogue, Buyer is deemed to have assented to these terms within ten (10) days of receipt of Seller's invoice or delivery of all or any part of the order, whichever is first, unless Seller receives written notice of any objection within 10 days after Buyer's receipt of this form. This contract and the terms and conditions shall be governed by and construed in accordance with the laws of the Commonwealth of Pennsylvania.
2. **CHANGES, CANCELLATION** - No change by Buyer of any term or condition of this contract or any of Seller's rights or remedies hereunder shall be binding on Seller, nor shall the order thereby acknowledged be cancelled or changed by Buyer, unless Seller shall expressly consent thereto in writing by Seller's authorized officer. There are no representations, agreements, promises or understandings between Buyer and Seller that are not expressed herein.
3. **DELIVERY, CLAIMS, DELAYS, RETURNS** - Unless specified differently in writing, all sales are F.O.B. Seller's shipping point. Delivery of goods to the carrier at Seller's plant or other loading point shall constitute delivery to Buyer and, regardless of shipping terms, all risk of loss or damage in transit shall be borne by Buyer. Seller reserves the right to make delivery in installments, all such installments to be separately invoiced and paid for when due per invoice, without regard to subsequent deliveries. Delay in delivery of any installment shall not relieve Buyer of its obligation to accept remaining deliveries.

Immediately upon Buyer's receipt of any goods shipped hereunder, Buyer shall inspect the same and shall notify Seller, in writing, of any claims for shortages, defects or damages and shall hold the goods for Seller's written instructions concerning disposition. If Buyer shall fail to so notify Seller, within five days after the goods have been received by Buyer, such goods shall conclusively be deemed to conform to the terms and conditions thereof.

Seller shall not be liable for any loss, damage or penalty as a result of any delay in or failure to manufacture, deliver or otherwise perform hereunder due to any cause beyond Seller's reasonable control, including, without limitation, unsuccessful reactions, act of Buyer, embargo or other governmental act, regulation or request, including but not limited to reinterpretation of any regulation, or change in enforcement or policy, affecting the conduct of Seller's business, fire, explosion, accident, theft, vandalism, riot, acts of war, strikes or other labor difficulties, lightning, flood, windstorm or other act of God, delay in transportation or inability to obtain necessary labor, fuel, materials, supplies or power at current prices.

4. **ALLOCATION OF GOODS** - If Seller is unable for any reason to supply the total demands for goods specified in Buyer's order, Seller may allocate its available supply among any or all buyers on such basis as Seller may deem fair and practical, without liability for any failure of performance which may result therefrom.
5. **PAYMENT** - Terms of sale are net 30 days of date of invoice, unless otherwise stated. Materials will be billed at the price in effect at the time shipment is made. If the financial condition of Buyer does not justify the terms of payment specified, Seller may demand full or partial payment in advance before proceeding with the contract. If Buyer defaults in any payment when due, under this or any other order, Seller at its option without prejudice to its other lawful remedies, may defer delivery or cancel this contract.
6. **TAXES AND OTHER CHARGES** - Any use tax, sales tax, excise tax, duty, custom, inspection or testing fee, or any other tax, fee or charge of any nature whatsoever imposed by any governmental authority, on or measured by the transaction between Seller and Buyer shall be paid by Buyer in addition to price quoted or invoiced. In the event Seller is required to pay any such tax, fee or charge, Buyer shall reimburse Seller therefore or, in lieu of such payment, Buyer shall provide Seller at the time the order is submitted an exemption certificate or other document acceptable to the authority imposing the tax, fee or charge.
7. **WARRANTIES** - Unless specified differently in writing, Seller warrants that its products shall conform to the description of such products as provided in the Seller's catalog. THIS WARRANTY IS EXCLUSIVE, AND SELLER MAKES NO OTHER WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURCHASE. Seller's warranties made in connection with this sale shall not be effective if Seller has determined in its sole discretion, that Buyer has misused the products in any manner or has failed to use the products in accordance with instructions, if any, furnished by the Seller and/or in accordance with any Material Safety Data Sheet.

Seller's sole and exclusive liability and Buyer's exclusive remedy with respect to products proved to Seller's satisfaction to be defective or nonconforming shall be the replacement of such product without charge or refund of the purchase price, at Seller's sole discretion, upon the disposition of such products in accordance with Seller's instructions. SELLER SHALL NOT BE LIABLE FOR ANY INCIDENTAL, CONSEQUENTIAL OR CONTINGENT DAMAGES.

8. **COMPLIANCE WITH LAWS, REGULATIONS** - Seller certifies that to the best of its knowledge its products are produced in compliance with applicable requirements of the Fair Labor Standards Act, as amended, and the Occupational Safety and Health Standards Act of 1970 and regulations, rules and orders issued pursuant thereto. Seller also certifies to the best of its knowledge, it is in compliance with Executive Order 11246 and its regulations, rules and orders issued pursuant thereto, that it does not discriminate against any employee or prospective employee because of race, creed, color, national origin, sex, age or handicap, nor permit discrimination in any form nor maintain segregated facilities for its employees; that it actively pursues employment of minorities, females, handicapped, disabled veterans and veterans of the Vietnam era; and that it uses its best efforts to award contracts to and place purchase orders with minority business enterprises and with labor surplus area concerns and small business concerns.



NOTICE: THIS LANGUAGE FORMS YOUR CONTRACT WITH CHEM SERVICE, INC. (cont'd)

9. **BUYER'S USE OF PRODUCTS** - Seller's products are intended only for laboratory use and, unless otherwise indicated on Seller's invoice, other writings, or on product labels, are not to be used for other purposes, including but not limited to as, or a component in a food, drug, medical device (including in vitro diagnostic reagents) or cosmetics as defined in the Federal Food, Drug and Cosmetic Act, as amended, nor as a pesticide unless otherwise stated by Seller in writing furnished to Buyer. Buyer expressly represents and warrants to Seller that Buyer will properly test, use, manufacture and market any products purchased from Seller and any final articles made from them in accordance with the practices of a reasonable person who is an expert in the field and in strict compliance with all applicable federal and state food, drug, device and cosmetic and other relevant laws and regulations, now and hereinafter enacted. Buyer further warrants to Seller that any final articles manufactured from Seller's products shall not be adulterated or misbranded within the meaning of the Federal Food, Drug, and Cosmetic Act and shall not be articles which may not be introduced into interstate commerce. Buyer agrees to comply with instructions, if any, furnished by Seller relating to the use of the products and not misuse the products in any manner.

Buyer acknowledges (1) that Seller's products are small quantities solely for purposes of scientific experimentation, or analysis, or chemical research and intended only for laboratory use by technically qualified individuals, (2) that the products will be used only for such purposes, and by such qualified individuals and, (3) that some of Seller's products may not be on the Inventory of Existing Chemical Substances defined pursuant to the Toxic Substances Control Act (TSCA). Seller assumes no responsibility to assure that products purchased hereunder may be used in a commercial application under the provisions of TSCA and regulations promulgated under TSCA. Seller will provide a Material Safety Data Sheet (MSDS) to Buyer which will inform persons using the product of the health risks which Seller has reason to believe may be associated with the product. Buyer agrees that it will inform its employees (1) that some of Seller's products may not be on the Inventory of Existing Chemical Substances defined pursuant to the Toxic Substances Control Act (TSCA) and that (2) a current MSDS is available to any and all employees or other personnel working, handling and/or in contact with or proximity to Seller's products for which an MSDS is required. Buyer agrees that it will comply with all provisions of 29 C.F.R. Part 1910. If the products purchased from Seller are to be repackaged, relabeled, or used as starting material or components of other products, Buyer will verify Seller's assay of the product and accept its obligations to conform to the requirements of TSCA.

Buyer agrees that it has the responsibility to verify the hazards and to conduct any further research necessary to learn the hazards involved in using products purchased from Seller. Buyer also has the duty to warn Buyer's employees, customers and any auxiliary personnel (such as freight handlers, etc.) of any risks involved in using or handling the products.

10. **INDEMNIFICATION** - Buyer agrees to indemnify and hold Seller harmless from and against any and all losses, damages and expenses, (including attorney's fees and other costs of defending any action) that Seller may sustain or incur as a result of any claim for damages for personal injury, death, damage to property or other injury by reason of any claim of negligence, breach of express or implied warranty, strict liability in tort or other theory of law made by Buyer, its officers, agents or employees, its successors and assigns, and its customers, whether direct or indirect, in connection with the use of Seller's products, or by reason of Buyer's failure to perform the obligations herein contained. This obligation of indemnification shall exist notwithstanding any ordinary negligence or fault by the Seller.

In the event of any claim against Seller by an employee of Buyer or an employee of a subcontractor of Buyer or other person who might or would benefit from an immunity from making indemnity, contained in any state or federal workers' compensation law or similar law, the Buyer specifically waives the immunity and agrees to indemnify Seller and hold Seller harmless from any claim for any type of injury or damages.

Buyer shall notify Seller within 15 days of Buyer's receipt of knowledge of any claim involving Seller's products resulting in personal injury, death, or damage to property and Buyer shall fully cooperate with Seller in the investigation and determination of the cause of such accident and shall make available to Seller all statements, reports and tests made by Buyer or made available to Buyer by others. The furnishing of such information to Seller and any investigation by Seller shall not constitute an assumption of any liability by Seller.

11. **PATENT DISCLAIMER** - Seller does not warrant that the use or sale of the products delivered hereunder will not infringe the claims of any United States or other patents covering the product itself or the use thereof in combination with other products or in the operation of any process.
12. **RETURNS AND DISPOSAL** - Goods may not be returned for credit except with Seller's written permission and then only in strict compliance with Seller's return shipment instructions. Buyer agrees that it will, at all times, handle and dispose of goods in accordance with Seller's instructions, the current Material Safety Data Sheet, and all applicable law and government regulations.
13. **TECHNICAL ASSISTANCE** - At Buyer's request, Seller may furnish technical assistance and information with respect to Seller's products. Unless otherwise agreed, all such technical assistance and information will be provided gratis, and Buyer assumes sole responsibility for results obtained in reliance thereon. SELLER MAKES NO WARRANTIES OF ANY KIND OR NATURE WITH RESPECT TO TECHNICAL ASSISTANCE OR INFORMATION PROVIDED BY IT. ANY SUGGESTIONS BY SELLER REGARDING USE, APPLICATION OR SUITABILITY OF THE PRODUCTS SHALL NOT BE CONSTRUED A WARRANTY UNLESS EXPRESSLY DESIGNATED AS SUCH, IN WRITING, SIGNED BY SELLER.
14. **RUSH ORDERS (SAME DAY SHIPMENT)** - Seller will accept rush orders for any in-stock items until 4pm EST, for next day delivery by Federal Express or UPS. This service is provided without additional charge to all customers, but Buyer will be responsible for any additional charge for shipping made by the carrier. Seller will use reasonable efforts to ship by next day service but Seller assumes no liability for failure to deliver a rush order the next day.



CHEM-SERVICE^{INC}

Pesticide and Metabolite Standards Catalog

Orders, Technical Inquiries, and Customer Service
www.chemservice.com
1-610-692-3026 or 1-800-452-9994

General: info@chemservice.com
Marketing: marketing@chemservice.com
Orders, Pricing, Availability: quote@chemservice.com
Technical Inquiries and Customer Service: tech@chemservice.com

Hours
Monday - Thursday: 8am - 6pm EST
Friday: 8am - 5pm EST

660 Tower Lane P.O. Box 599
West Chester, PA 19381-0599



Certificate Number:31610



Accreditation Number:63520



Reference Material Producer
ISO GUIDE 34

Accreditation Number:63520